

Cell Wash Solution II/Acid Wash

cobas®

CellIC2/Detergent 2

Order information

REF	CONTENT	Analyzer(s) on which kit(s) can be used
04880307 190	Cell Wash Solution II/Acid Wash (2 x 1.8 L)	Roche/Hitachi MODULAR Roche/Hitachi cobas c

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche Diagnostics representative.

English

Intended use

Cell Wash Solution II and Acid Wash are mainly used as acid wash solutions for sample and reagent probes.

Summary

Reagent and sample probe washes may be required due to potential interference from other reagents or samples. These special washes maintain reagent and sample integrity.

Extra wash cycles are required with specific reagent type (sample type) combinations, e.g. if a preceding test interferes with an assay through carry over by reagent (sample) probe.

Reagents - working solutions

Buffer; detergent: 310 mmol/L citric acid monohydrate

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Hazardous component: Citric acid monohydrate



Danger

H314 Causes severe skin burns and eye damage.

Prevention:

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

Response:

P301 + P330 IF SWALLOWED: rinse mouth. Do NOT induce vomiting. + P331

P303 + P361 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. + P353

P305 + P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. + P338

P310 Immediately call a POISON CENTER or doctor/physician

Disposal:

P501 Dispose of contents/container to an approved waste disposal plant.

Contact phone: all countries: +49-621-7590, USA: +1-800-428-2336

Reagent handling

Ready for use

Storage and stability

Shelf life at 15-25 °C: See expiration date on reagent label.

On-board stability: 12 weeks

Materials provided

See "Reagents – working solutions" section for reagents.

Materials required (but not provided)

- See "Order information" section

Roche/Hitachi analyzers and reagents.

See appropriate package inserts and operator's manuals for additional required materials.

Assay

Use Cell Wash Solution II/Acid Wash as specified in the respective instructions for use of the system reagents.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.

CONTENT	Contents of kit
REAGENT	Reagent
CALIBRATOR	Calibrator
→	Volume after reconstitution or mixing

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS C and MODULAR are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

Distribution in USA by:
Roche Diagnostics, Indianapolis, IN
US Customer Technical Support 1-800-428-2336

