REF	CONTENT	Analizatorius (-iai), su kuriuo (-iais) galima naudoti cobas c kasetę (-es)
04525299 190	Creatine Kinase-MB 100 tyrimų	Sistemos-ID 07 5924 4
11447394 216	Calibrator f.a.s. CK-MB (3 x 1 mL)	Kodas 402
11447378 122	Precinorm CK-MB (4 x 3 mL)	Kodas 320
04358210 190	Precipath CK-MB (4 x 3 mL, neprieinama JAV)	Kodas 356
05117003 190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Kodas 391
05947626 190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Kodas 391
05947626 160	PreciControl ClinChem Multi 1 (4 x 5 mL, skirta JAV)	Kodas 391
05117216 190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Kodas 392
05947774 190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Kodas 392
05947774 160	PreciControl ClinChem Multi 2 (4 x 5 mL, skirta JAV)	Kodas 392
04489357 190	Diluent NaCl 9 % (50 mL)	Sistemos-ID 07 6869 3

Lietuvių

Sistemos informacija

Skirta **cobas c** 311/501 analizatoriams:**CKMBL**: ACN 060Skirta **cobas c** 502 analizatoriui:**CKMB**: ACN 8060

Paskirtis

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- Lott JA, Stang JM. Serum enzymes and isoenzymes in the diagnosis and differential diagnosis of myocardial ischemia and necrosis. Clin Chem 1980;26:1241-1250.
- Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders 1987;346-421.
- Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. Clin Chem Lab Med 2002;40(6):635-642.
- Mathieu M, Bretaudiere JP, Galteau MM, et al. Recommendations for measuring the catalytic concentration of creatine kinase in human serum at 30°C. Ann Biol Clin 1982;40:138-149.
- Hørder M, Magid E, Pitkänen E, et al. Recommended method for the determination of creatine kinase in blood modified by the inclusion of EDTA. Scand J Clin Lab Invest 1979;39:1-5.
- Bergmeyer HU, Breuer H, Büttner H, et al. Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Standard-Methode zur Bestimmung der Aktivität der Creatin-Kinase. J Clin Chem Clin Biochem 1977;15:249-254.
- Würzburg U, Hennrich N, Lang H, et al. Determination of creatine kinase-MB in serum using inhibiting antibodies. Klin Wschr 1976;54(8):357-360.
- Hørder M, Elser RC, Gerhardt W, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Provisional recommendation IFCC method for creatine kinase Appendix A. J Int Fed Clin Chem 1990;2:26-35.
- Braun S, Röschenhaler F, Jarausch J, et al. Analyte Stability of CK-MB Activity and cTnT in ICU Patient Serum and Heparin Plasma. Poster presented at Medica 2004, Düsseldorf. (Roche Diagnostics GmbH No. 04587979990).
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- Wu AHB, Bowers GN. Evaluation and comparison of immunoinhibition and immunoprecipitation methods for differentiating MB from BB and macro forms of creatine kinase isoenzymes in patients and healthy individuals. Clin Chem 1982;28:2017-2021.
- Remaley AT, Wilding P. Macroenzymes: Biochemical Characterization, Clinical Significance, and Laboratory Detection. Clin Chem 1989;35:2261-2270.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- Klein G, Berger A, Bertholf R, et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. Clin Chem 2001;47:Suppl. A30.
- Thomas L, Müller M, Schumann G, et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29(5):301-308.
- Stein W. Strategie der klinischen-chemischen Diagnostik des frischen Myokardinfarktes. Med Welt 1985;36:572-577.
- Myocardial Infarction Redefined - A Consensus Document of the Joint European Society of Cardiology/ American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Eur Heart J 2000;21:1502-1513.

- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

JAV platina:

Roche Diagnostics, Indianapolis, IN

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© 2014, Roche Diagnostics