REF	CONTENT		Analizatorius (-iai), su kuriuo (-iais) galima naudoti <b>cobas c</b> kasetę (-es)
06600239 190	Tina-quant Cystatin C Gen.2 (225 tyrimai)	Sistemos-ID 07 7550 9	Roche/Hitachi <b>cobas c</b> 311 Roche/Hitachi <b>cobas c</b> 501/502
04975901 191	C.f.a.s. Cystatin C (4 x 1 mL)	Kodas 407	
04975936 190*	Cystatin C Control Set Control I (low) (4 x 1 mL) Control II (high) (4 x 1 mL)	Kodas 121 Kodas122	
06729371 190	Cystatin C Control Set Gen.2 Control 1 (3 x 1 mL) Control 2 (3 x 1 mL) Control 3 (3 x 1 mL)	Kodas 139 Kodas 140 Kodas 141	

\*Neskiarta JAV

## Lietuvių

## Sistemos informacija

Skirta **cobas c** 311/501 analizatoriams:**CYSC2:** ACN 109Skirta **cobas c** 502 analizatoriui:**CYSC2:** ACN 8109

## Paskirtis

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**Nuorodos**

- 1 Levey AS, Coresh J, Balk E, et al. National Kidney Foundation practice guidelines for chronic kidney disease: evaluation, classification, and stratification. Ann Intern Med 2003;139:137-147.
- 2 Rule AD, Larson TS, Bergstralh EJ, et al. Using serum creatinine to estimate glomerular filtration rate: accuracy in good health and in chronic kidney disease. Ann Intern Med 2004;141:929-937.
- 3 Wasen E, Isoaho R, Mattila K, et al. Estimation of glomerular filtration rate in the elderly: a comparison of creatinine based formulae with serum Cystatin C. J Intern Med 2004;256:70-78.

- 4 Levey AS, Coresh JK. DOQI clinical practice guidelines on chronic kidney disease. Guideline 4. Estimation of GFR. Am J Kidney Dis 2002;39 (suppl1):76-92.
  - 5 Kyhse-Andersen J, Schmidt C, Nordin G, et al. Serum Cystatin C, determined by a rapid, automated particle-enhanced turbidimetric method, is a better marker than serum creatinine for glomerular filtration rate. Clin Chem 1994;40(10):1921-1926.
  - 6 Grubb A, Simonsen O, Sturfelt G, et al. Serum concentration of Cystatin C, factor D and  $\beta$ 2-microglobulin as a measure of glomerular filtration rate. Acta Med Scand 1985;218(5):499-503.
  - 7 Jung K, Jung M. Cystatin C: a promising marker of glomerular filtration rate to replace creatinine. Nephron 1995;70(3):370-371.
  - 8 Newman DJ, Thakkar H, Edwards RG, et al. Serum Cystatin C measured by automated immunoassay: a more sensitive marker of changes in GFR than serum creatinine. Kidney Int 1995;47(1):312-318.
  - 9 Mussap M, Dalla Vestra M, Fioretto P, et al. Cystatin C is a more sensitive marker than creatinine for the estimation of GFR in type 2 diabetic patients. Kidney Int. 2002;61:1453-1461.
  - 10 Dhamidharka VR, Kwon C, Stevens G. Serum cystatin C is superior to serum creatinine as a marker of kidney function: a meta analysis. Am J Kidney Dis 2002 Aug;40(2):221-226.
  - 11 Risch L, Drexel H, Huber AR. Differences in glomerular filtration rate estimates by two cystatin C-based equations. Clinical Chemistry 2005;51:2211-2212.
  - 12 Grubb A, Nyman U, Bjork J, et al. Simple Cystatin C-Based prediction equations for glomerular filtration rate compared with the modification of diet in renal disease prediction equation for adults and the Schwartz and the Counahan-Barratt prediction equations for children. Clinical Chemistry 2005;51:1420-1431.
  - 13 Hoek FJ, Kemperman FA, Krediet RT. A comparison between Cystatin C, plasma creatinine and the Cockcroft and Gault formula for the estimation of glomerular filtration rate. Nephrol Dial Transplant 2003;18:2034-2031.
  - 14 Lassus J, Harjola VP, Sund R, et al. Prognostic value of Cystatin C in acute heart failure in relation to other markers of renal function and NT-proBNP. European Heart Journal (February) doi:10.1093/eurheart/ehl507.
  - 15 Strevens H, Wide-Svensson D, Torffvit O, et al. Serum Cystatin C for assessment of glomerular filtration rate in pregnant and non-pregnant women. Indications of altered filtration process in pregnancy. Scand J Clin Lab Invest 2002;62:141-148.
  - 16 Andersen TB, Eskild-Jensen A, Frokiaer J, et al. Measuring glomerular filtration rate in children; can Cystatin C replace established methods? Pediatr Nephrol 2009;24:929-41.
  - 17 Grubb A, Nyman U, Björk J. Improved estimation of glomerular filtration rate (GFR) by comparison of eGFR (cystatin C) and eGFR (creatinine). Scand J Clin Lab Invest 2012;72:73-77.
  - 18 Filler G, Huang S-HS, Yasin A.. The usefulness of Cystatin C and related formulae in pediatrics. Clin Chem Lab Med 2012;50(12):2081-2091.
  - 19 Grubb A. Cystatin C- and creatinine-based GFR-prediction equations for children and adults. Clinical Biochemistry 2011;44:501-502.
  - 20 Rollins G. Predicting risk of ESRD in Diabetics, Clinical Laboratory Strategies, Jan.2013,1-4.
  - 21 Leseley A, Inker M.D., Christopher H, et al, Estimating Glomerular Filtration Rate from Serum Creatinine and Cystatin C. N Engl J Med 2012;367:20-29.
  - 22 Horio M, Imai E, Yasuda Y, et al. GFR Estimation Using Standardized Serum Cystatin C in Japan. Am J Kidney Dis 2013;61(2):197-203
  - 23 Grubb A, Horio M, Hansson LO, et al. Generation of a New Cystatin C-Based Estimating Equation for Glomerular Filtration Rate Using Seven Assays Standardized to the International Calibrator.Clin Chem 2014, in press.
  - 24 van Deutekom AW, Zur B, van Wijk JAE, et al. Measurement of cystatin C in capillary blood samples in pediatric patients. Clin Biochem 2010;43:335-337.
  - 25 Erlandsen EJ, Randers E, Kristensen JH. Evaluation of the Dade Behring N Latex Cystatin C assay on the Dade Behring Nephelometer II System. Scand J Clin Lab Invest 1999;59:1-8.
  - 26 Mussap M, Ruzzante N, Varagnolo M, et al. Quantitative Automated Particle-Enhanced Immunonephelometric Assay for the Routinary Measurement of Human Cystatin C. Clin Chem Lab Med 1998;36(11):859-865.
  - 27 Gislefoss RE, Grimsrud TK, Morkrid L. Stability of selected serum proteins after long-term storage in the Janus Serum Bank. Clin Chem Lab Med 2009;47(5):596-603.
  - 28 Data on file at Roche Diagnostics.
  - 29 Bökenkamp A, van Wijk JAE, Lentze ML, et.al. Effect of Corticosteroid Therapy on Serum Cystatin C and  $\beta$ 2-Microglobulin Concentrations. Clin Chem 2002;48:1123-1126.
  - 30 Wiesli P, Schwegler B, Spinaz AS, et al. Serum Cystatin C is sensitive to small changes in thyroid function Clin Chim Acta 2003;338(1):87-90.
  - 31 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
  - 32 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
  - 33 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
  - 34 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
  - 35 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Šiuose informaciniuose lapeliuose kaip dešimtainės trupmenos skyrklis visada naudojamas taškas, skiriantis sveikąjį skaičių nuo dešimtainės trupmenos skaitmenų. Tūkstančių skyrkliai nenaudojami.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
www.roche.com

JAV platina:

Roche Diagnostics, Indianapolis, IN

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© 2014, Roche Diagnostics