
GE Medical Systems
 Information Technologies, Inc.
 8200 West Tower Avenue
 Milwaukee, WI 53223 USA
 Tel: + 1 414 355 5000
 1 800 558 7044 (US Only)
 1 800 668 0732 (Canada Only)
 Fax: + 1 414 355 3790

GE Medical Systems
 Information Technologies GmbH
 Munzinger Straße 5
 79111 Freiburg
 GERMANY
 Tel: +49 761 4543-0
 Fax: +49 761 4543-233

Asia Headquarters
 GE (China) Co., Ltd.
 No.1 Huatuo Road,
 Zhangjiang Hi-Tech Park, Pudong,
 Shanghai, P.R.China 201203
 Tel: +86 21 38777888
 Fax: +86 21 38777402

通用电气医疗系统信息技术公司
 8200 West Tower Avenue
 Milwaukee, WI 53223 USA 美国
 Tel: + 1 414 355 5000
 1 800 558 7044 (US Only)
 1 800 668 0732 (Canada Only)
 Fax: + 1 414 355 3790

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Electrode Contact Spray

Description

Electrode Contact Spray for ECG application systems employing suction electrodes.

200 ml 21730701
 10 x 200 ml 21730702
 2 L 21730705

Intended Use

The Electrode Contact Spray is intended to support the quality of electrical connection between devices and the human body for ECG acquisition with reusable suction electrodes.
 The Electrode Contact Spray is intended for use on adult and pediatric patients.

The Electrode Contact Spray is intended for use by qualified medical personnel under the direct supervision of a licensed health-care practitioner.

Purpose of Document

This document contains the instructions how to use the Electrode Contact Spray. For further instructions on skin preparation, see the device operator's manual.

For external use only! Not suitable for long-term use!

Instructions for Use

- Apply an adequate amount of spray to the **skin**, not onto electrodes.
- Avoid direct inhalation of fumes.
- After use, remove remaining spray with water.
- Avoid contact with open wounds.
- Avoid eye contact.
- If spray has entered the eyes, rinse thoroughly with water.

/ See bottom of bottle.

Responsibility of the Manufacturer

GE Healthcare is responsible for the effects of safety, reliability, and performance only if the Electrode Contact Spray is used in accordance with the instruction for use.

Device Classification


This product is assigned to a Class I medical device according to 93/42/EEC for Europe.

CE Marking Information

The product bears the CE marking, indicating it is in compliance with the regulatory requirements of the European Directive 93/42/EEC concerning medical devices.

Reordering Information

To reorder the Electrode Contact Spray, contact your local GE Healthcare representative.

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Electrode Contact Spray

Safety Data Sheet according to 91/155/EEC and 2001/58/EC

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1	Product Identification
1.1	Commercial Name: Electrode Contact Spray
1.2	Commercial Part No: 217 307 01 217 307 02 217 307 05
1.3	Manufacturer: GE Medical Systems Information Technologies GmbH Munzinger Straße 5 79111 Freiburg GERMANY Tel.: +49 761 4543 - 0
2	Composition (according to INCI)
2.1	Aqua, Isopropyl Alcohol, Sodium Chloride, Cetylpyridinium Chloride, Sodium Borate.
3	Health and environmental Hazards
3.1	Health: This product does not originate any special health hazard.
3.2	Environment: This product does not originate any special environmental hazard.
4	First Aid and Emergency Actions
4.1	Inhalation: Avoid direct inhalation of the fumes. Airing of closed rooms.
4.2	Skin contact: No precautions required.
4.3	Eye contact: Rinse the eyes with water.
4.4	Ingestion: Drink plenty of water. Consult physician.
5	Measures in Case of Fire
5.1	No particular measures required. Put out in case of fire with water or CO ₂ .
5.2	Fumes heavier than air. Under unfavourable conditions formation of explosive mixture possible.

6	Measures in case of accidental Leaks
6.1	No inhalation of the fumes. Airing of closed rooms. Wash with clear water.
6.2	Residues do not dispose in sewage system.
7	Handling and Storage
7.1	Handling: No precautions required
7.2	Storage: Keep in a cool and dry and well aired place far from source of heat and fire. Do not expose to ultraviolet rays for a long time.
8	Particular Protection Measures
8.1	Exposition limit: Not applicable
8.2	Personal protective measures: Inhalation protection: No protective measures required. Skin protection: No protective measures required. Eye protection: No protective measures required. Other protections: Clean the hands after finished work
9	Physical and chemical Data
9.1	Physical aspects: water based propanol solution
9.2	Color: colorless, clear
9.3	Odor: alcohol
9.4	Boiling point: 82-100 °C
9.5	Inflammation point: > 450 °C
9.6	Flash point: 32 °C
9.7	pH: 6.5 – 7.5
9.8	Weight density: 0.94 – 0.99 g/cm ³
9.9	Viscosity: Not applicable
9.10	Consistency: Not applicable
9.11	Electrical Conductivity: >5.0 mS/cm
10	Stability and Reactivity Data
10.1	Conditions to be avoided: enforced heating
10.2	Materials to be avoided: organic nitro compounds.
10.3	Dangerous products in case of decomposition: Peroxides.
11	Toxic Data
11.1	No toxic potential. No skin irritation.
12	Ecological Data
12.1	In case of proper use and handling no adverse environmental influence to expect.
12.2	Biodegradable, COD = 200 mgO ₂ /g

13 Disposal Data

13.1 According statutory and regulatory requirements of the country where used.
Please contact the appropriate agency of your country.

14 Transport

14.1 No transport requirements.

15 Regulatory and normative Data

15.1 Product is in compliance with the Medical Device Directive 93/42/EEC.
15.2 Product meets requirements for biocompatibility according ISO 10993.

16 Other Data

None

Information contained in this document are based on the momentary state of knowledge of the product and are provided to describe appropriate means of measures for safety. It is referred to the product indicated and does not give any particular quality guarantees.

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© 2011 General Electric Company - All rights reserved.