
Safe, effective, simple

Regional citrate anticoagulation during CRRT

Prismocitrate® 18/0 anticoagulant solution
Prismocal® B22 dialysis solution

 **GAMBRO®**

Leveraging on experience to improve therapy

Clinical use of dilute citrate solutions makes treatment safe, effective and simple.¹⁻⁶

Over the past decade, Gambro has partnered with the medical community to develop safe solutions for regional citrate anticoagulation (RCA) in continuous renal replacement therapy (CRRT).

One of the key findings is the advantages of dilute citrate solutions over concentrated solutions. While the original descriptions of citrate anticoagulation in CRRT used concentrated solutions (ACD-A, 4% TSC), these protocols required pharmacy-made, hypotonic solutions and complex monitoring, and have been associated with serious adverse events.⁷⁻⁹

In contrast, the clinical use of dilute citrate solutions provides a simpler and safer alternative for RCA during CRRT.¹⁻⁶

Introducing new solutions for regional citrate anticoagulation

Prismocitrate 18/0 anticoagulant solution

| | |
|----------|------------|
| Citrate | 18 mmol/l |
| Sodium | 140 mmol/l |
| Chloride | 86 mmol/l |

Prism0cal B22 dialysis solution

| | |
|-------------|--------------|
| Bicarbonate | 22 mmol/l |
| Lactate | 3 mmol/l |
| Sodium | 140 mmol/l |
| Potassium | 4 mmol/l |
| Calcium | 0 mmol/l |
| Magnesium | 0.75 mmol/l |
| Chloride | 120.5 mmol/l |
| Glucose | 6.1 mmol/l |

Safe

- **Low risk of adverse events**, such as hypernatremia or alkalosis, which are often reported when using high concentrated citrate solutions.¹⁻⁶
- **Safe management of acid-base balance** through bicarbonate-buffered, calcium-free dialysate with key electrolytes at physiological concentrations.^{2,4}
- **Low risk of metabolic complications** in case solutions are inadvertently interchanged.

Effective


- **Optimal anticoagulation.** Prismocitrate 18/0 solution easily provides the blood citrate concentration needed to effectively lower ionized calcium level in the circuit and inhibit coagulation.
- **Improved treatment quality.** Longer circuit patency with citrate yields reduced treatment downtime, higher cumulative dose and thus improved treatment quality. This property ensures that each patient receives the prescription of fluid exchange prescribed.
- **Easy achievement of targets.** The fluid and electrolyte balance is achieved predominantly by the addition of concomitant dialysate and/or post replacement fluid at a rate to achieve either a weight-based dose or a specific end-point such as the normalization of acid-base status.

Simple

- **Pre-formulated.** Literature reports multiple citrate protocols, some quite complex requiring pharmacy-made fluids. With pre-formulated solutions from Gambro, citrate protocols become much simpler and could be standardized.
- **Easier citrate implementation.** Prismocitrate 18/0 solution with the Prismaflex eXeed™ II system significantly simplifies the implementation of citrate protocols. The system allows convective and diffusive therapies with fully integrated citrate management and calcium substitution to the patient.

Simple implementation with the Prismaflex eXeed™ II system

Typical set-up for regional citrate anticoagulation

Citrate dose: 3 mmol/l

Suggested CRRT dose: 35 ml/kg/h

References

1. Palsson, R., Niles, J.L. Regional citrate anticoagulation in continuous venovenous hemofiltration in critically ill patients with a high risk of bleeding. *Kidney International* (1999), 55:1991-1997.
2. Tolwani, J.A. et al. A practical citrate anticoagulation continuous venovenous hemodiafiltration protocol for metabolic control and high solute clearance. *Clin J Am Soc Nephrol* (2006) 1:79-87.
3. Wybo, R. et al. Comparison of two citrate solutions for regional anticoagulation during hemofiltration. *Intensive Care Med* (2006), 32 (Suppl. 1):S160.
4. Davies, H. et al. A regional citrate anticoagulation protocol for pre-dilutional CVVHDF: the Modified Alabama Protocol. *Aust Crit Care* (2008) Aug;21(3):154-65.
5. Schmitz, M. et al. Continuous venovenous haemofiltration using a citrate buffered substitution fluid. *Anaesth Intensive Care* (2007); 35:529-535.
6. Egi, M. et al. A comparison of two citrate anticoagulation regimens for continuous veno-venous hemofiltration. *Int J Artif Organs* (2005); 28(12):1211-1218.
7. Mehta, R.L. et al. Regional citrate anticoagulation for continuous arteriovenous hemodialysis in critically ill patients. *Kidney International* (1990); 38:976-981.
8. Monchi, M. et al. Citrate vs. Heparin for anticoagulation in continuous venovenous hemofiltration: a prospective randomized study. *Intensive Care Med* (2004); 30:260-265.
9. Cointault, O. et al. Regional citrate anticoagulation in continuous venovenous haemodiafiltration using commercial solutions. *Nephrol Dial Transplant*. (2004); 19(1):171-8.

Gambro, Prismocitrate, Prismocal and Prismaflex eXeed are trademarks of the Gambro Group.

These products are CE-marked in accordance with the requirements in EC Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Gambro Lundia AB
PO Box 10101, SE-22010 Lund
Sweden
Phone + 46 46 16 90-00
partner@gambro.com
www.gambro.com

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