INSTRUCTIONS FOR USE

LYFO DISK®

- **LYFO DISK® Microorganisms**
- **KWIK-STIK™ Microorganisms**
- **KWIK-STIK™ Plus Microorganisms**

INTENDED USE

LYFO DISK®, KWIK-STIK™ and KWIK-STIK™ Plus Microorganisms are lyophilized, reference stock culture preparations containing a single strain of a microorganism. These microorganism preparations are intended to be used for quality control of culture media, educational/instructional programs and industrial applications. The microorganism preparations are traceable to the American Type Culture Collection (ATCC®) or other authentic reference culture collections.

SUMMARY AND HISTORY

A reliable source of reference stock cultures for use in microbiology quality assurance programs is essential. Microorganisms with known and predictable characteristics are used in quality control, education and proficiency programs. Lyophilization is a well-documented and recommended method for long-term preservation of microorganisms. The use of this lyophilized material provides equivalent results to traditional methods used in preparing, storing and maintaining reference stock culture collections.

PRINCIPLE

LYFO DISK®, KWIK-STIK™ and KWIK-STIK™ Plus Microorganisms incorporate a lyophilization method reported by Obara et.al. which uses a suspending medium consisting of gelatin, skim milk, ascorbic acid, dextrose, and charcoal. The gelatin serves as a carrier for the microorganism. Skim milk, ascorbic acid, and dextrose protect the microorganism by preserving the integrity of the cell wall during freeze-drying and storage. The charcoal is included to neutralize any toxic substances formed during the lyophilization process.

A safer, healthier world.

KWIK-STIK™

KWIK-STIK™ PLUS

PRODUCT DESCRIPTION

- A. KWIK-STIK™ Microorganisms** Each KWIK-STIK™ unit contains a lyophilized pellet of a single microorganism strain, a reservoir of hydrating fluid and an inoculating swab. Each device is sealed within a laminated pouch that contains a desiccant to prevent adverse moisture accumulation. The KWIK-STIK™ Microorganisms offer an additional feature.
- Each lyophilized microorganism preparation is less than or equal to four (4) passages from a reference culture.
- B. KWIK-STIK™ Plus Microorganisms** The packaging of the KWIK-STIK™ Plus Microorganism is identical to the KWIK-STIK™ Microorganisms. The KWIK-STIK™ Plus Microorganisms offer two additional features.
- Each lyophilized microorganism preparation is two (2) passages from a reference culture.
 - A Certificate of Assay is provided. It offers documentation regarding the identity and traceability of the microorganism preparation to a reference culture and the number of passages the microorganism preparation has been removed from the reference culture.
- C. LYFO DISK® Microorganisms** LYFO DISK® Microorganisms are packaged in a resealable vial that contains ten (10) lyophilized pellets of a single microorganism strain and a desiccant to prevent adverse accumulations of moisture. The LYFO DISK® Microorganisms offer an additional feature.
- Each lyophilized microorganism preparation is less than or equal to four (4) passages from a reference culture.

PRECAUTIONS AND LIMITATIONS

- These products are for in-vitro use only.
- Refer to the MSDS for more detailed information. The MSDS can be found on our website at www.microbiologics.com in the support center document library.
- These devices, and growth of these microorganisms, are considered biohazard material.
- These devices contain viable microorganisms that may produce disease. Proper techniques must be employed to avoid exposure and contact with any microorganism growth.
- The microbiology laboratory must be equipped, and have the facilities to receive, process, maintain, store and dispose of biohazard material.
- Only trained laboratory personnel should use these devices.
- Agencies and statutes regulate the disposal of all biohazard materials. Each laboratory must be aware of, and comply with, the proper disposal of biohazard materials.
- All of Microbiologics products and packaging are latex free.

STORAGE AND EXPIRATION

Store the LYFO DISK®, KWIK-STIK™ and KWIK-STIK™ Plus Microorganisms at 2°C to 8°C in the original, sealed vial or pouch containing the desiccant. Stored as directed, the lyophilized microorganism preparation will retain, until the expiration date stated on the device label, its specifications and performance within the stated limits.

The LYFO DISK®, KWIK-STIK™ and KWIK-STIK™ Plus Microorganisms should not be used if:

- Stored improperly;
- There is evidence of excessive exposure to heat or moisture; or,
- The expiration date has passed.

INSTRUCTIONS FOR USE ---

A. KWIK-STIK™ and KWIK-STIK Plus™ Microorganism Procedure

1. Allow the unopened KWIK-STIK™ pouch to equilibrate to room temperature. Tear open pouch at notch and remove the KWIK-STIK™ unit.
2. Tear off Pull-Tab portion on the label and attach it to the primary culture plate or QC record. Do not disassemble the device during hydration.
3. Pinch (once only) the ampoule at the top of the KWIK-STIK™ (just below the fluid meniscus of the ampoule) found in the cap to release the hydrating fluid.
4. Hold vertically and tap on a hard surface to facilitate flow of fluid through shaft into bottom of unit containing pellet. Allow the hydrating fluid to flow through the swab shaft and into the bottom portion of the unit containing the pellet.
5. Using a pinching action on the bottom portion of the unit, crush the pellet in the fluid until the pellet suspension is homogenous.
6. **IMMEDIATELY** heavily saturate the swab with the hydrated material and transfer to agar medium.
7. Inoculate the primary culture plate(s) by gently rolling the swab over one-third of the plate.
8. Using a sterile loop, streak to facilitate colony isolation.
9. Using proper biohazard disposal, discard the KWIK-STIK™.
10. **IMMEDIATELY** incubate the inoculated primary culture plate(s) at temperature and conditions appropriate to the microorganism.

B. LYFO DISK® Microorganism Procedure

1. Remove the unopened LYFO DISK® vial from 2°C to 8°C storage and allow the unopened vial to equilibrate to room temperature.
2. Aseptically remove one (1) pellet with sterile forceps from the vial. Do not remove desiccant.
3. Place the pellet in 0.5 mL of sterile fluid (water, saline, TSB, or BHIB). **IMMEDIATELY** stopper and recap vial and return the resealed vial to 2°C to 8°C storage.
4. Crush the pellet with a sterile swab until the suspension is homogenous. **IMMEDIATELY** heavily saturate the same swab with the hydrated material and transfer to agar medium.
5. Inoculate the primary culture plate(s) by gently rolling the swab over one-third of the plate.
6. Using a sterile loop, streak to facilitate colony isolation.
7. Using proper biohazard disposal, discard the remaining hydrated material.
8. **IMMEDIATELY** incubate the inoculated media at temperature and conditions appropriate to the microorganism.

MATERIALS REQUIRED BUT NOT PROVIDED ---

- LYFO DISK® Microorganisms require sterile tubes and 0.5 ml of sterile liquid such as, Tryptic Soy Broth, Brain Heart Infusion Broth, saline, or deionized water to hydrate the lyophilized preparation. Sterile swabs or inoculating loops are needed to transfer the hydrated preparation to an agar plate.
- LYFO DISK®, KWIK-STIK™ and KWIK-STIK Plus Microorganisms require non-selective, nutrient or enriched agar media and specific incubation times and conditions to optimize growth and recovery.

The Technical Information Bulletin (TIB.081) *“Recommended Growth Requirements”* lists the recommended media and incubation requirements. This bulletin is available from our website support center document library at www.microbiologics.com.

KEY OF SYMBOLS

Authorized Representative in the European Community

Batch Code (Lot)

Biological Hazards Biological Risk

CE Mark

Catalog Number

Caution consult accompanying documents Attention, see instructions for use

In Vitro Diagnostic Medical Device

Manufacturer

Temperature Limitation

Use By

QUALITY CONTROL

This product is developed, manufactured, and distributed:

- in compliance with the mandates of FDA: Quality System Regulation (QSR), 21CFR Part 820;
- in conformance with the elements of ISO 9001; and,
- in conformance with CE Mark requirements.

Quality control functions may include, but are not limited to:

- purity and growth characteristics;
- morphological features;
- biochemical activity;
- the identity and traceability of the microorganism preparation to a reference culture; and,
- the number of passages the microorganism preparation has been removed from the reference culture.

The decision to perform additional quality control is the responsibility of each individual laboratory.

PRODUCT WARRANTY

These products are covered under warranty to meet the specifications and performance printed and illustrated in product inserts, instructions, and supportive literature. The warranty, expressed or implied, is limited when:

- the procedures employed in the laboratory are contrary to printed and illustrated directions and instructions or
- the products are employed for applications other than the intended use cited in product inserts, instructions, and supportive literature.

REFERENCES

The following reference cites the basis for the lyophilization method employed on these microorganism preparations.

1. Y. Obara, S. Yamai, T. Nikkawa, Y. Shimoda, and Y. Miyamoto. 1981. J. Clin. Microbiol. 14:61-66.

The selection of reference stock cultures is only one integral part of the overall scheme for QC challenge procedures and techniques. Reference to guidelines for each laboratory's applications is essential. Examples might include:

1. AOAC Compendium of Microbiological Methods.
2. Clinical Microbiology Procedures Handbook. ASM. Washington, D.C.
3. FDA Bacteriological Analytical Manual.
4. Manual of Clinical Microbiology, ASM, Washington, D.C.
5. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically. CLSI.
6. Official Methods of Analysis of the Association of Official Analytical Chemists.
7. Performance Standards for Antimicrobial Disk Susceptibility Tests. CLSI.
8. Quality Assurance for Commercially Prepared Microbiological Culture Media. CLSI.
9. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. CLSI.
10. Standard Methods for the Examination of Dairy Products.
11. Standard Methods for the Examination of Water and Wastewater.
12. US Pharmacopoeia and National Formulary.

WEBSITE

Visit our website for current technical information and product availability.

www.microbiologics.com

ACKNOWLEDGEMENTS

Microbiologics, Inc.
217 Osseo Avenue North
St. Cloud, MN 56303 USA
Tel. 320-253-1640
Fax. 320-253-6250
Email. info@microbiologics.com

MediMark® Europe
11, rue Emile Zola B.P. 2332
38033 Grenoble Cedex 2, France
Tel. 33 (0)4 76 86 43 22
Fax. 33 (0)4 76 17 19 82
Email. info@medimark-europe.com

* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

ILLUSTRATED INSTRUCTIONS

The LYFO DISK® vial contains 10 lyophilized pellets of an individual microorganism strain.

1

Remove the unopened LYFO DISK® vial from 2°C to 8°C storage and allow the unopened vial to equilibrate to room temperature.

2


Aseptically remove one (1) pellet with sterile forceps from the vial. **Do not remove desiccant.**

3

Place the pellet in 0.5 mL of sterile fluid (water, saline, TSB, or BHI(B)). **IMMEDIATELY** stopper and recap vial and return the resealed vial to 2° to 8° storage.

4

Crush the pellet with a sterile swab until the suspension is homogenous. **IMMEDIATELY** heavily saturate the same swab with the hydrated material and transfer to agar medium.

5

Inoculate the primary culture plate(s) by gently rolling the swab over one-third of the plate.

6


Using a sterile loop, streak to facilitate colony isolation.

7

Using proper biohazard disposal, discard the remaining hydrated material.

8

IMMEDIATELY incubate the inoculated media at temperature and conditions appropriate to the microorganism.

A safer, healthier world.

KWIK-STIK™ & KWIK-STIK™ PLUS

ILLUSTRATED INSTRUCTIONS


KWIK-STIK™ and KWIK-STIK™ Plus microorganism preparations contain a lyophilized pellet of a single strain of microorganism.

1

Allow the unopened KWIK-STIK™ pouch to equilibrate to room temperature. Tear open pouch at notch and remove the KWIK-STIK™ unit.

2


Tear off Pull-Tab portion on the label and attach it to the primary culture plate or QC record. Do not disassemble the device during hydration.

3


Pinch (once only) the ampoule at the top of the KWIK-STIK™ (just below the fluid meniscus of the ampoule) found in the cap to release the hydrating fluid.

4

Hold vertically and tap on a hard surface to facilitate flow of fluid through shaft into bottom of unit containing pellet. Allow the hydrating fluid to flow through the swab shaft and into the bottom portion of the unit containing the pellet.

5

Using a pinching action on the bottom portion of the unit, crush the pellet in the fluid until the pellet suspension is homogenous.

6

IMMEDIATELY heavily saturate the swab with the hydrated material and transfer to agar medium.

7

Inoculate the primary culture plate(s) by gently rolling the swab over one-third of the plate.

8

Using a sterile loop, streak to facilitate colony isolation.

9

Using proper biohazard disposal, discard the KWIK-STIK™.

10

IMMEDIATELY incubate the inoculated primary culture plate(s) at temperature and conditions appropriate to the microorganism.

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