
CHROMOGENIC URINE AGAR II

Powdered medium and ready to use plates for the isolation, counting and direct identification of urinary tract pathogens.

Chromogenic Urine Agar II: *E.coli* (pink colonies), *K.pneumoniae* blue colonies, *Enterococcus* (small green colonies)

Typical formula (g/l)

Peptones	18.7
Tryptophan	2.0
Chromogenic mix	13.8
MUG	0.1
Agar	12.0

Directions for powdered medium

Suspend 46,6 in 1000 ml of cold distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C and distribute into sterile Petri dishes. Use within 48 hours from the preparation.

pH 6.8 ± 0.2.

Description

Chromogenic Urine Agar II is a diagnostic medium useful for the isolation, counting and direct presumptive rapid identification of urinary tract pathogens: *E.coli*, *Enterobacter-Klebsiella-Serratia* (KES), *Proteus-Morganella-Providencia*, Enterococci, Staphylococci.

The differentiation between the different bacterial species or genus is achieved by:

- * a chromogenic substrate for β -galactosidase (GAL) which is split with the liberation of an insoluble pink dye.

- * a chromogenic glucopyranoside derivative which is split by β -glucosidase (GLU) with the formation of an insoluble blue-green dye.

- * MUG for the fluorogenic detection of β -glucuronidase enzyme (GUR).

- * Tryptophan for the detection of tryptophan deaminase (TDA) of *Proteus* spp. *Morganella* spp., *Providencia* spp. and for indole test of *E.coli* colonies.

The opaque white background of the medium enhances the colour of the colonies.

Technique

Chromogenic Urine Agar II (CUA II) can be used according to the usual laboratory practices for urine bacterial count, by spreading the specimen on the agar surface and incubating at 37°C for 18-24 hours. The cultivated colonies can be identified with the following scheme:

Colonies on CUA II	Enzyme activity	Additional characterisation on CUA II	Enzyme activity	Results
Large pink colonies	GAL +	fluorescent under Wood lamp indole positive	GUR + IND +	<i>E.coli</i>
Pink colonies	GAL +	not fluorescent under Wood lamp		GUR - supplementary identification needed
Small blue-green colonies	GLU + GAL -			Enterococci <i>St.agalactiae</i>
Blue to purple colonies	GLU + GAL +			<i>Klebsiella, Enterobacter, Serratia</i>
Brown colonies with brown halo	TDA + GLU -			<i>Proteus-Morganella-Providencia</i> glucosidase negative strains
Blue colonies with brown halo	TDA + GLU +			<i>Proteus-Providencia</i> glucosidase positive strains
Colourless colonies				supplementary identification tests needed

Notes

- * It is reported that some strains of the bacterial genus reported above have abnormal biochemical patterns .
- * Between the *Proteus-Morganella-Providencia* group, *P.mirabilis* is indole negative and can be easily recognised
- * Biochemical identification is needed for species identification within *Klebsiella Enterobacter Serratia* (KES) group
- * A Pyroglutamate (PYR) test might be necessary to differentiate enterococci from *S.agalactiae*
- * *S.saprophyticus* and *S.xylosus* produces small pink colonies and are MUG negative.
- * Biochemical tests are necessary for species identification of microorganisms producing colourless colonies.
- * Gram staining is recommended to confirm any doubtful colour reactions.

User quality assurance (24 h-37°C)

Productivity control

E.coli ATCC 25922: growth, large pink colonies fluorescent under Wood's lamp, indole positive

E.aerogenes ATCC 13048: growth, blue colonies not fluorescent under Wood's lamp, indole negative

S.aureus ATCC 25923: growth, colourless colonies

Storage

Dehydrated medium: 2-8°C

Ready to use plates: 2-8°C

User prepared plates: up to 48 hours at 2-8°C

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

E.coli : blue colonies; *Citrobacter freundii*: red colonies

Typical formulas

Dehydrated and ready to use media (g/L)

Enzymatic digest of casein	1,00
Yeast extract	2,00
Sodium chloride	5,00
Sodium dihydrogen phosphate x 2H ₂ O	2,20
Di-sodium hydrogen phosphate	2,70
Sodium pyruvate	1,00
Sorbitol	1,00
Tryptophane	1,00
Tergitol® 7	0,15
Salmon-β-D-galactoside	0,20
X-β-glucuronide CHX salt	0,10
Isopropyl-β-D-thiogalactopyranoside (IPTG)	0,10
Bacteriological agar	10,60

Intended use

Powdered medium and ready to use plates for the enumeration of *Escherichia coli* and coliform bacteria with membrane filtration method for waters with low bacterial background, according to ISO 9308-1:2014

Principle of the method and explanation

Chromogenic Coliform Agar ISO Formulation is a selective and differential medium for the simultaneous enumeration of *E. coli* and coliform bacteria in waters with low bacterial background with the membrane filtration technique, according to the ISO Standard 9308-1:2014.

The medium is made selective by the presence of Tergitol 7 an inhibitor of Gram positive bacteria. The differentiation between coliform bacteria, different from *E.coli* and *E. coli* is given by the presence of Salmon-β-D-galactoside, a chromogenic substrate for the detection of β-galactosidase and X-β-Glucuronide, a chromogenic substrate for the detection of β-glucuronidase.

Salmon-β-D-galactoside is hydrolysed by coliforms releasing a pink-red colour; this reaction is strengthened in the medium by the presence of IPTG (isopropyl-β-D-thiogalactopyranoside).

X-β-glucuronide is hydrolysed, among *Enterobacteriaceae*, by *E. coli*, and by a few strains of *Salmonella* and *Shigella* releasing a blue pigment.

Directions for preparation from dehydrated medium

Suspend 27,1 g in 1000 ml of cold purified water. Heat to boiling with frequent agitation until completely dissolved. Do not autoclave, do not overheat. Some turbidity may occur after boiling; it normally disappears when the temperature decreases to 45-50 °C and doesn't affect the medium performances


Dispense in petri dishes to a depth of at least 4 mm. If not for immediate use, the plates can be stored at 5 ± 3 °C in the dark and protected against evaporation for not more than one month. There should be no visible moisture on the plates before use. When moisture is present, the plates should be dried for the minimum time required to remove visible moisture.

Physical characteristics

Dehydrated medium appearance: yellowish, fine, homogeneous, free-flowing powder.

Prepared medium appearance: light yellow, limpid or slightly opalescent.

Final pH: 6,8 ± 0,2

Biolife

Technical Sheet

N° 401297-CCA ISO EN-0 Nov-14 page 2/2

Technique

Filter 100 ml (or other volumes, e.g. 250 ml for bottled water) of the sample using a membrane filter usually about 47 mm or 50 mm in diameter, with filtration characteristics equivalent to a rated nominal pore diameter of 0,45 µm and, preferentially, with grid lines. The minimum volume for filtration is 10 ml of sample or dilutions thereof to ensure even distribution of the bacteria on the membrane filter.

After filtration place the membrane filter on the Chromogenic Coliform Agar ISO Form., ensuring that no air is trapped underneath, invert petri dish, and incubate at 36 ± 2 °C for 18-24 h.

Examine the membrane filters and count:

- All pink to red colonies (positive for β -D-galactosidase reaction) as presumptive coliform bacteria that are not *E. coli*.
- All dark-blue to violet colonies (positive for β -D-galactosidase and β -D-glucuronidase reactions) as *E. coli*.

To avoid false-positive results, caused by oxidase positive bacteria, for example, *Aeromonas* spp, the presumptive colonies shall be confirmed by a negative oxidase reaction (Oxidase test Strips, cat. N° 191040ST)

Perform the oxidase test preferentially on all, or at least 10 selected pink to red colonies (coliform bacteria, different from *E.coli*), the test should be negative.

The count of total coliform bacteria is the sum of all oxidase negative pink to red colonies plus all dark-blue to violet colonies.

Limitations

Some strains of *E. coli* which are β -D-glucuronidase negative, such as *Escherichia coli* O157, will not be detected as *E. coli*. As they are β -D-galactosidase positive, they will appear as coliform bacteria on this chromogenic agar.

Quality control: microbiological characteristics

CONTROL STRAINS			INOCULATION INOCUBATION	SPECIFICATIONS	
<i>E. coli</i>	ATCC	25922	Q 37°C /24H / A	BLUE COLONIES	A/C $\geq 0,7$
<i>E. coli</i>	ATCC	8739	Q 37°C /24H / A	BLUE COLONIES	A/C $\geq 0,7$
<i>E.aerogenes</i>	ATCC	13048	Q 37°C /24H / A	PINK COLONIES	A/C $\geq 0,7$
<i>C.freundii</i>	ATCC	43864	Q 37°C /24H / A	RED COLONIES	A/C $\geq 0,7$
<i>P.aeruginosa</i>	ATCC	10145	E 37°C /24H / A	GOOD GROWTH COLOURLESS COLONIES	
<i>E. faecalis</i>	ATCC	19433	MM 37°C /24H / A	GROWTH INHIBITED	

A: AEROBIC INCUBATION

A/C (Productivity Ratio): CFU OBTAINED ON TEST BATCH / CFU OBTAINED ON TRYPTIC SOY AGAR

Q: QUANTITATIVE MF (MEMBRANE FILTRATION) METHOD / E: ECOMETRY / MM: MODIFIED MILES MISRA

Storage

Dehydrated medium REF 4012972: keep tightly closed, away from bright light, at +2°C to +8 °C

User prepared plates : up to 1 month at +2°C to +8 °C, in the dark, protected against evaporation.

Ready to use plates (REF 491297): keep at +2°C to +8 °C

Precautions

The products described are for Laboratory use only.

References

ISO 9308-1:2014 Water quality - Enumeration of *Escherichia coli* and coliform bacteria - Part 1: Membrane filtration method for waters with low bacterial background flora.

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PEPTONE BACTERIOLOGICAL

Peptone obtained by enzymatic hydrolysis of animal tissue.

USE

PEPTONE BACTERIOLOGICAL is an enzymatic hydrolysate of meat that supplies a limpid, colorless and very stable watery solution. It is used in the preparation of culture media as a nitrogen source readily available for bacterial growth. It is a general use very nutritive peptone, with neutral pH.

PHYSICO-CHEMICAL CHARACTERISTICS

	Standard
Solubility in water at 2%	Complete
pH of 2% solution	7.0+/-0.5
Loss on drying	≤ 6.0%
Total nitrogen	>12.5%
α-amino nitrogen AN	3-4.5%
Ash	5.0%

TECHNIQUE

Peptone Bacteriological can be used as an ingredient of dehydrated culture media and need dissolution in distilled or deionized water and sterilization by autoclaving.

QUALITY CONTROL

Dehydrated powder

Appearance: free-flowing, homogeneous.

Color: cream.

STORAGE

The powder is very hygroscopic; store the powder at 10-30 °C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident.

REFERENCES

1. Standard Methods for Examination of Water and Sewage, 15th ed., (1980).
2. J. Dairy Science, 16: 277: (1933).

PRESENTATION


Product	REF	
PEPTONE BACTERIOLOGICAL	611701	500 g
PEPTONE BACTERIOLOGICAL	621701	100 g
PEPTONE BACTERIOLOGICAL	6117015	5 Kg

TABLE OF SYMBOLS

LOT Batch code	 Caution, consult accompanying documents	 Manufacturer	 Contains sufficient for <n> tests	IVD In Vitro Diagnostic Medical Device
REF Catalogue number	 Fragile, handle with care	 Use by	 Temperature limitation	 Keep away from heat source

 **LIOFILCHEM s.r.l.**
Via Scozia, Zona Industriale - 64026, Roseto D.A. (TE) - ITALY
Phone +390858930745 Fax +390858930330
Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

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TRYPTONE

Enzymatic digest of casein for preparing microbiological culture media.

PHYSIC-CHEMICAL CHARACTERISTICS

	Standard
Solubility in water at 5%	Complete
pH (5% solution)	7.2 ± 0.5
Loss on drying	≤ 6.0%
Total nitrogen	12.1-13.1%
α-amino nitrogen AN	2-3%
AN/TN x 100	26-37
Chloride (as NaCl)	≤ 1.0%
Residue on ignition	≤ 17.0%

DESCRIPTION

TRYPTONE is an enzymatic digest of casein for use in preparing microbiological culture media.

PRINCIPLE

TRYPTONE is a nitrogen source rich in tryptophan and carbohydrate free. Hence, it is suitable for use in detecting indole production and in differentiating bacteria on the basis of their ability to ferment various carbohydrates.

PREPARATION

Refer to appropriate references for specific preparations.

TECHNIQUE

TRYPTONE can be used as an ingredient of dehydrated culture media and need dissolution in distilled or deionized water and sterilization by autoclaving.

INTERPRETATION OF RESULTS

Refer to appropriate references for test results.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous, it is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used by properly trained operators only.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

1. U.S. Department of Agriculture (1998) Microbiology laboratory guidebook, 3rd ed. Food safety and Inspection service, USDA, Washington, DC.
2. Horowitz (ed) (2000) Official method of analysis of AOAC International, Gaithersburg, Md.

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

PRODUCT SPECIFICATIONS

NAME
TRYPTONE

PRESENTATION
Dehydrated powdered

STORAGE
10-30°C

PACKAGE

Ref.	Content	Packaging
611004	500 g	500 g of powder in plastic bottle
621004	100 g	100 g of powder in plastic bottle
6110045	5000 g	5 kg of powder in plastic container

pH OF THE MEDIUM
7.2 ± 0.5

USE
TRYPTONE is an enzymatic digest of casein for use in preparing microbiological culture media

TECHNIQUE
Refer to technical sheet of the product

SHELF LIFE
4 years

QUALITY CONTROL
Dehydrated powder
Appearance: free-flowing, homogeneous
Colour: cream

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instructions for use	 Keep away from heat sources

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NUORODOS

1. U.S. Departamento f Agriculture (1998) Microbiology laboratory guidebook, 3rd ed. Food safety and Inspection service, USDA, Washington, DC.
2. Horowitz (ed) (2000) Official methods of analysis of AOAC International, Gqathersburg, Md.

PAKUOTĖS

Code 611004TRYPTONE.....500 g

Code 621004TRYPTONE.....500 g

Code 6110045TRYPTONE.....5Kg