INDEX, ÍNDICE, ΕΥΡΕΤΗΡΙΟ, INDEKSAS	
1- 9	ENGLISH
10-19	FRANÇAIS
20- 21	SYMBOLS SIMBOLOS SIMBOLI SYMBOL SYMBOLES ΣΥΜΒΟΛΟ SIMBOLIS SYMBOLY SYMBOLE

Instructions for Use

Thermo Scientific™ Sensititre™ Susceptibility Plates for
Anaerobes (*B. fragilis* group)

007- Anaerobe ROW Only –CID9248
Revision Date: 4 th December 2018

Sensititre Susceptibility Plates for Anaerobes (B.fragilis group)

Sensititre 46-48 hours

For *in vitro* diagnostic use

For full plate information, including plate layout, QC information, interpretative criteria, performance data and references, please refer to www.trekds.com/techinfo. The plate code and batch number will be required.

INTENDED USE

The Sensititre Anaerobe Susceptibility Plate is designed for *in vitro* diagnostic use in determining antimicrobial susceptibility of *B. fragilis* isolates.

SUMMARY AND PRINCIPLES OF USE

The Sensititre susceptibility plate is a micro broth dilution method in a dried format and provides qualitative (susceptible, intermediate, resistant) and semi- quantitative (minimum inhibitory concentration, MIC) results. Each susceptibility plate is dosed with antimicrobics at appropriate dilutions and dried to extend plate storage periods. Plates are inoculated with organisms suspended in Thermo Scientific Sensititre broth and incubated at 34-36°C in an anaerobic atmosphere for 46 – 48 hours. Results are determined by visually observing the lowest antimicrobial concentration showing inhibition of growth.

PRECAUTIONS

Results should be used as an aid in selecting the drug of choice for treatment. This product is for *in vitro* diagnostic use and should be used by properly trained personnel. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test plates after use. Directions should be read and followed carefully.

STORAGE AND SHELF LIFE

Store Sensititre susceptibility plates in original containers at 15 – 25°C and away from direct sunlight and heat until used. Each plate is packaged in foil with a desiccant. Plates should not be used beyond their expiration date, if the foil pouch is damaged, or the desiccant color is not orange. Plates should be inoculated within 5 hours of removing from pouch.

SPECIMEN COLLECTION AND PREPARATION

Specimens should be collected and handled following recommended guidelines.

PROCEDURE

Materials supplied

- 10 Sensititre™ anaerobe susceptibility plates
- 10 Perforated adhesive seals

Materials required but not supplied [Product code]:

Broth:

- Sensititre™ Cation Adjusted Mueller-Hinton Broth with TES (CAMHBT) [T3462-05]
- Sensititre Supplemented Brucella Broth for Anaerobes [YT3450]

Turbidity:

- Thermo Scientific 0.5 McFarland Standard [E1041]
- Sensititre Nephelometer™ [V3011]

Plate inoculation:

For manual plate inoculation

- 100µL pipettor and disposable tips

For automated plate inoculation

- Sensititre AIM™ Automated Inoculation Delivery System [V3020]
- Doseheads for Sensititre AIM for plate inoculation [E3010]

Plate reading:

For manual plate reading

- Sensititre Manual Viewbox [V4007]

For optimized manual plate reading

- Sensititre Vizion Digital MIC Viewing System™ [V2021]

Other supplies:

- Sterile inoculum reservoir
- Quality control strains (see Table 1)
- Incubator 34-36°C, anaerobic jar or chamber
- Vortex mixer
- Current CLSI, EUCAST, or local guideline documents

SELECTION OF SUSCEPTIBILITY TEST BROTH

Thermo Scientific Sensititre broth is validated for use with Sensititre susceptibility plates. Sensititre susceptibility plates for testing anaerobic organisms have been validated for use with Sensititre Supplemented Brucella broth.

INOCULATION PROCEDURE

Allow broth to equilibrate to room temperature before use.

1. Test organisms must be grown anaerobically in pure culture and examined by Gram stain prior to use in the system.
2. Test organisms may be removed from Brucella agar with sheep's blood. Cultures used for inoculum preparation should be less than 72 hours old. The agar plate should not remain in an aerobic atmosphere for more than 30 minutes prior to making the suspension.
3. Using a cotton swab or inoculating loop, suspend 3 – 5 colonies from the agar plate in 5mL Sensititre Mueller-Hinton Broth and adjust to achieve a turbidity equivalent to 0.5 McFarland turbidity standard. Mix well.
4. Invert the tube 8-10 times
5. Using a calibrated pipette transfer 100µL of the suspension into a tube of Sensititre Supplemented Brucella broth to achieve an inoculum of 1×10^6 cfu/mL. The final inoculum should be inoculated into the plate within 15 minutes of preparation.
6. Transfer 100µl to each well of the plate by either:
 - a. **Sensititre AIM.** Replace the tube cap with a Sensititre single-use dosehead and inoculate the plate according to the AutoInoculator/AIM instructions.

Remove the test tube/dosehead combination from the AutoInoculator /AIM and store inverted in a rack, or discard.

- b. **Manual pipette.** Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette (e.g. 8-channel pipette). Carefully inoculate 100µL into

each well of the Sensititre® plate. Avoid drawing up any air bubbles as they will interfere with the accurate delivery of the inoculum. Since the tips are used to inoculate a series of wells, avoid touching the tip to the bottom of any well.

7. A periodic check of the colony count of the positive control well is recommended. (See Appendix 1). Isolates should have an inoculum of approximately 1×10^6 cfu/mL.
8. Cover all wells with the perforated adhesive seal. Avoid creases as these can lead to skips.
9. Plates should be inoculated within 5 hours of removal from their pouch

INCUBATION

Incubate inoculated plates at 34-36°C in an anaerobic atmosphere for 46-48 hours. No more than three plates should be stacked while incubating.

READING TEST RESULTS

After incubation, results can be manually read using the Sensititre Manual Viewbox or the Vizion for computer assisted reading. Refer to the Sensititre Vizion User Manual for instructions on operating the Vizion. Plate adhesive seals do not need to be removed prior to reading results. Growth appears as turbidity or as a deposit of cells at the bottom of a well. The MIC is usually recorded as the lowest concentration of the antimicrobial that inhibits visible growth. The positive growth control wells should be read first. If there is no growth in any of the positive growth control wells, the results are invalid.

The following points should be noted:

a. **Trailing End Points**

Most organism / antimicrobial combinations give distinct end points. With some combinations there may be a gradual fading of growth over 2 to 3 wells. End points should be determined as the first well that inhibits visible growth with the exception of Sulphonamides where the end point is determined as an 80% decrease in growth compared to the positive control well.

Examples for reading and determining endpoints are shown in the CLSI, M11-A8 document (2).

b. **Contamination**

Single well contamination appears as growth in one well that is bordered by wells showing no growth. Single well contamination can be ignored, however, if multiple well contaminants are suspected, the test should be repeated.

c. **Skips**

Occasionally a “skip” may be seen - a well showing no growth bordered by wells showing growth. There are a variety of explanations for skips including contamination, mutation, creased seal and misaligned dosing. A single skip can be ignored. However, in order to ensure effective antimicrobial therapy, **never** read the skip well as the MIC; always read the lowest well concentration above which there is consistently no growth.

d. **Mixed Cultures**

With the exception of trailing end points as described in section (a) above, if two end points are seen as a distinct “button” of cells followed by several wells of diffuse growth with the “button” no longer visible (or seen as smaller buttons), there may be a mixed

bacterial population. Purity can be checked by sub-culturing the growth onto a suitable agar plate. Test results are invalid if a mixed culture is detected.

INTERPRETATION OF RESULTS

For interpretation of results, refer to the MIC interpretive guidelines provided by CLSI, EUCAST, or local standard guidelines.

MIC Results

The MIC is recorded as the lowest concentration of antimicrobial that inhibits visible growth or shows a significant reduction in growth. If no growth occurs in any well containing the dilutions of antimicrobial, other than the positive control well, the MIC should be recorded as less than or equal to the lowest concentration of antimicrobial on the plate.

QUALITY CONTROL – Refer to www.trekds.com/techinfo

Frequency of quality control testing should be conducted in accordance to local guidelines. Inoculum should be cultured onto a suitable medium to check for purity. Test results are invalid if a mixed culture is detected.

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The inoculum suspension should be cultured onto a suitable medium to check for purity, and incubated both aerobically and anaerobically. Test results are invalid if a mixed culture is detected.

All Sensititre® plates include positive control wells. Tests are invalid unless there is distinct growth in all positive control wells.

A number of factors influence MICs including organism state, inoculum density, temperature, and broth. In practice replicate MICs form a normal distribution with the majority of results lying within one dilution of the modal value. The test procedure can be considered satisfactory if control organism MICs are within range. Results should **not** be reported if QC results are not in range.

b. Selection of Broth

Each lot of Sensititre® supplemented Brucella broth is performance tested for use in Sensititre anaerobe susceptibility products.

Contact Thermo Fisher Scientific Microbiology Technical Support for assistance in the event that quality control discrepancies cannot be resolved.

EXPECTED VALUES

For expected values of quality control organisms, refer to Table 1. For expected values of clinically occurring organisms, it is recommended that each testing site generate a comprehensive antibiogram to determine the percentage of susceptibility to each antimicrobial tested. National based antibiograms are available to provide general guidelines of expected values (2). Antibiotic resistance may vary due to nosocomial infections and geographic location. It has been reported for certain antimicrobics that increased resistance may occur over the life span of the antimicrobial

Table 1: Acceptable Quality Control ranges* for QC organisms

Organism ATCC Number	<i>B. fragilis</i> ATCC® 25285™	<i>B. thetaiotaomicron</i> ATCC® 29741™
ANTIMICROBIC		
AMPICILLIN	<u>8->16</u>	<u>8-32</u>
AMOXICILLIN / CLAVULANIC ACID	0.25/0.12-1/0.5	0.25/0.12-1/0.5
AMPICILLIN/SULBACTAM	0.5/0.25-2/1	0.5/0.25-2/1
CEFOTETAN	1-8	16-128
CEFOXITIN	2-8	8-64
CEFTIZOXIME	-	-
CHLORAMPHENICOL	4-16	<u>4-16</u>
CLINDAMYCIN	0.5-2	2-8
DOXYCYCLINE	-	2-8
ERTAPENEM	0.06-0.5	0.5-2
GARENOXACIN	0.06-0.25	0.25-2
IMIPENEM	0.03-0.25	<u>0.12-0.5</u>
LINEZOLID	2-8	2-8
MEROPENEM	0.03-0.25	0.06-0.5
METRONIDAZOLE	0.25-2	0.5-4
MEZLOCILLIN	<u>8-32</u>	<u>8-32</u>
MOXIFLOXACIN	0.12-.5	1-8
PENICILLIN G	8-32	8-32
PIPERACILLIN	4-16	8-64
PIPERACILLIN/TAZOBACTAM	0.03/4-0.25/4	2/4-16/4
TETRACYCLINE	<u><0.25-1</u>	<u>8-32</u>
TICARCILLIAN/CLAVULANIC ACID	0.06/2-0.5/2	0.5/2-2/2

Note: *Ranges different or additional to Table 2J-1 of CLSI M100 S26E are underlined.

Contact Thermo Scientific Sensititre Technical Support in the event that quality control discrepancies cannot be resolved .

LIMITATIONS

1. This product is designed to conform to recommendations of the CLSI anaerobe reference procedure (2).
2. In common with all other methods of antimicrobial susceptibility testing, the results generated by Sensititre susceptibility plates are *in vitro* results.
3. Trained personnel are necessary to set up and interpret results.

PERFORMANCE

Plates are designed to give comparable performance to the CLSI and EUCAST reference micro broth procedure. Comparable performance is defined as >90% essential agreement with the reference MIC on a typical mix of organisms (see References).

For further information, please contact Thermo Scientific Sensititre Technical Support.

APPENDIX 1: Colony Count Procedure.

Perform immediately after inoculating a plate.

1. Transfer 10µL of the inoculated broth into 10 ml of reduced 0.85% saline or reduced broth.
2. Mix and transfer 100µL of this dilution onto an appropriate agar plate.
3. Spread the inoculum over the entire surface of the plate with the aid of a sterile loop or bent glass rod.
4. Incubate at 34-36° C for 48 hours in an anaerobic atmosphere.
5. Count the colonies, 100 colonies indicates an inoculum of 1×10^6 cfu/mL.

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DISCLAIMER

The information provided in this technical insert is current at the time of printing and may change without notice.

The latest information can be downloaded from the www.TREKDS.com/techinfo or by contacting Thermo Fisher Scientific Microbiology Technical Support

Legal Manufacturer:
TREK Diagnostic Systems Ltd
Units 17-19 Willard Way, Birches Industrial Estate,
East Grinstead,
West Sussex.
RH19 1XZ,
United Kingdom
Tel: +44-1342-318777

The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC cultures.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Technical Support email addresses

Austria	microbiology.techsupport.de@thermofisher.com
Belgium	microbiology.techsupport.benelux@thermofisher.com
Czech Republic	mikrobiologie.tech.podpora.cz@thermofisher.com
Denmark	microbiology.techsupport.nordic@thermofisher.com
France	microbiologie.techsupport.fr@thermofisher.com
Finland	mikrobiologia.tekninentuki.fi@thermofisher.com
Germany	microbiology.techsupport.de@thermofisher.com
Italy	microbiologia.supportotecnico.it@thermofisher.com
Netherlands	microbiology.techsupport.benelux@thermofisher.com
Norway	microbiology.techsupport.nordic@thermofisher.com
Spain	microbiologia.soporte.es@thermofisher.com
Sweden	microbiology.techsupport.nordic@thermofisher.com
UK	microbiology.techsupport.uk@thermofisher.com

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