INDEX, ÍNDICE, ΕΥΡΕΤΗΡΙΟ, INDEKSAS	
1-15	ENGLISH
16-31	ESPAÑOL
32-46	ITALIANO
47-62	DEUTSCHE
63-79	ΕΛΛΗΝΙΚΟ
80-95	FRANÇAIS
96-109	CZECH
110-124	<i>LIETUVIŲ KALBA</i>
125-140	POLSKI
141-156	ROMÂNĂ
157-158	SYMBOLS SIMBOLOS SIMBOLI SYMBOL SYMBOLES ΣΥΜΒΟΛΟ SIMBOLIS SYMBOLY SYMBOLE

Instructions for Use

Thermo Scientific™ Sensititre™ Susceptibility Plates
for testing non-fastidious Gram negative and Gram
positive isolates

039-NFAST ROW Only -CID9494
Revision Date:13th December 2018

Sensititre Susceptibility Plates for testing non-fastidious Gram negative and Gram positive isolates

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For *in vitro* Diagnostic use

For full plate information, including plate layout, QC information, Interpretative criteria, performance data and references please refer to www.trekds.com/techinfo. The plate type and plate lot number will be required.

INTENDED USE

The Sensititre MIC and Breakpoint Susceptibility system is an *in vitro* diagnostic product for clinical susceptibility testing of non-fastidious Gram negative isolates, comprising of *Enterobacteriaceae*, *Pseudomonas aeruginosa*, and other non-*Enterobacteriaceae* and of non-fastidious Gram positive isolates, comprising of *Staphylococcus* sp., *Enterococcus* sp., and Beta haemolytic *Streptococci* other than *S. pneumoniae*. The Sensititre ESBL confirmatory test plate is an *in vitro* diagnostic product for detection of ESBLs in clinical isolates of *Klebsiella pneumoniae*, *Klebsiella oxytoca* and *Escherichia coli*. MIC and ESBL plates can either be read manually or automatically on the Sensititre Autoreader / OptiRead and/or ARIS. Thermo Scientific manufactured broths have only been validated with Sensititre products .

SUMMARY AND PRINCIPLES OF USE

Each plate is dosed with antimicrobial agents at appropriate dilutions. Results can be read manually by visual reading of growth or automatically on an ARIS / OptiRead /Autoreader using fluorescence.

The Sensititre Autoreader /OptiRead systems utilize fluorescence technology. The technology involves the detection of bacterial growth by monitoring the activity of specific surface enzymes produced by the test organism. Growth is determined by generating a fluorescent product from a non-fluorescent (fluorogenic) substrate. The non-fluorescent substrate is prepared by conjugating a fluorescent compound to the specific enzyme substrates with a bond, which prevents fluorescence. The fluorophore is then said to be quenched. Enzymatic action of the bacterial surface enzymes on the specific substrates cleaves this bond releasing the fluorophore which is now capable of fluorescence. The amount of fluorescence detected is directly related to the activity of the bacterial surface enzymes and, therefore, to the bacterial growth.

Only plates which have the format name suffixed with F can be read automatically.

PRECAUTIONS

Only instruments supported by Sensititre (i.e. a simple manual viewer, Sensitouch, Vizion, Sensititre Autoreader / OptiRead and ARIS) must be used to report results with CE IVD and FDA cleared Sensititre products. Any other system used will not be supported.

Results should be used as an aid in selecting the drug of choice for treatment. This product is for *in vitro* diagnostic use and should be used by properly trained personnel. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test plates after use. Directions should be read and followed carefully.

STORAGE AND SHELF LIFE

Store Sensititre susceptibility plates in original containers at 15 – 25°C and away from direct

sunlight and heat until used. Each plate is packaged in foil with a desiccant. Plates should not be used beyond their expiration date, if the foil pouch is damaged, or the desiccant color is not orange. Plates should be inoculated within 5 hours of removing from pouch.

SPECIMEN COLLECTION AND PREPARATION

Specimens should be collected and handled following recommended guidelines.

Materials included

- 10 Sensititre™ plate with substrate in wells
- 10 adhesive seals

Materials not included [Product Code]:

Broth:

- Sensititre™ demineralised water [T3339]
- Sensititre™ Cation adjusted Mueller-Hinton broth with TES (CAMHBT) [T3462]

Turbidity:

- Sensititre Nephelometer™ [V3011]
- 0.5 McFarland polymer turbidity standard [E1041]

Plate inoculation:

- Doseheads for Sensititre AIM for plate inoculation [E3010]
- Sensititre Automated Inoculation Delivery System™ [V3020] / Sensititre Autoinoculator*

Plate reading:

- Sensititre AutoReader™* OptiRead™ [V3030] or ARIS™ [V3090]
- Sensititre Vizion Digital MIC Viewing System™ [V2021]
- Sensititre Manual Viewbox [V4007]

Other supplies:

- Sterile inoculum reservoir
- Quality control strains (see Table 1)
- Incubator 34-36°C
- Vortex mixer
- Current CLSI, EUCAST, or local guideline documents

**No longer sold or supported*

SELECTION OF SUSCEPTIBILITY TEST BROTH

Thermo Scientific Sensititre broth is validated for use with Sensititre susceptibility plates.

INOCULATION PROCEDURES FOR SENSITITRE PLATES

Allow all broths equilibrate to room temperature before use.

Standard Inoculum

For all plates pick 3-5 colonies of non-fastidious isolates from fresh primary agar plate, emulsify in sterile water, and adjust to a 0.5 McFarland Standard. Mix well.

Transfer 10µl of the McFarland suspension into 11 ml tube of cation adjusted Muller Hinton broth with TES buffer. Proceed to step 2.

Increased inoculum method

Depending on the strains tested, an increase to 30 µl may aid in detecting resistance mechanisms.

For both Gram positive and negative isolates, the transfer of 30µl of the suspension into an 11 ml Sensititre Mueller-Hinton broth tube results in colony counts which fall within the cleared Sensititre range of 5.0×10^4 and 5.0×10^5 cfu/ml. Proceed to step 2.

Proteus, Providencia and Morganella

For tribe *Proteae* (including *Proteus* species, *Providencia* species and *Morganella morganii*), transfer 1µl into the broth.

2. Transfer 50µl of the broth suspension into the plate within 30 min of preparation by either:
 1. **Sensititre AutoInoculator / AIM.** Replace the tube cap with a Sensititre single-use dosehead and inoculate the plate according to the AutoInoculator / AIM instructions. Remove the test tube/dosehead combination from the AutoInoculator / AIM within 30 seconds of dosing a plate.
 2. **Manual pipette.** Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette.
3. Perform a purity check on all final inocula after plate inoculation.
4. A periodic check of the colony count should be performed (see Appendix 1). Isolates should have an inoculum of 1×10^5 cfu/ml, (range 5.0×10^4 – 5.0×10^5) with the exception of tribe *Proteae* isolates which should be 1×10^4 cfu/ml (range 5.0×10^3 – 5.0×10^4).
5. Cover all wells of each plate with the adhesive seal provided. Avoid creases as these can lead to skips.
6. Plates should be inoculated within 5 hours of removal from their pouch.

INCUBATION

All non-fastidious aerobic organisms should be incubated at 34-36°C either in the ARIS or in a non-CO₂ incubator for 18 –24 hours.

- In order to ensure detection of vancomycin-resistant *Enterococci* and oxacillin- resistant *Staphylococci*, incubate for 24 hours.
- *Acinetobacter* spp., *Burkholderia cepacia*, and *Stenotrophomonas maltophilia* should be incubated for 20 -24 hours.

Sealed plates can be stacked up to 3 high in an incubator if not incubated in the ARIS.

Summary of set up and incubation

Organism group	McFarland Suspension medium	Intended Final Inoculum* (cfu/ml)	Inoculum transfer	Broth	Plate Reconstitution	Incubation conditions	Incubation Hours
Non-fastidious Gram negative	Water	1x10 ⁵	10µL	MHB	50µL	35°C O ₂	18-24
Non-fastidious Gram positive	Water	1x10 ⁵	10µL	MHB	50µL	35°C O ₂	18-24
Increased inoculum for both Gram negative and Gram positive organisms**	Water	3x10 ⁵	30µL	MHB	50µL	35°C O ₂	18-24
Tribe <i>Proteae</i>	Water	1x10 ⁴	1µL	MHB	50µL	35°C O ₂	18-24

* See appendix #1 for inocula ranges

** Depending on strains being tested, an increase to 30ul may aid in the detecting resistant mechanisms

READING TEST RESULTS

Check purity test plates. Results are invalid if a mixed culture is present

1. Automatically

Plates may be read on the AutoReader / OptiRead or the ARIS according to the instructions in the appropriate User Manual or by using the SWIN HELP files.

Sensititre plates with intended use of Manual, Sensitouch or Vizion only cannot be read on the AutoReader / OptiRead or ARIS.

2. Manually

After incubation, results can be read using the Sensititre manual viewer or Vizion. See User Manuals. It is not necessary to remove the adhesive seal. Growth appears as turbidity or as a deposit of cells at the bottom of a well. The MIC is recorded as the lowest concentration of antimicrobial that inhibits visible growth. Reading faint growth on the Vizion can be improved by adjusting the lighting. The positive growth control wells should be read first. If any show no growth, results are invalid.

The following points should be noted:

a. Fading End Points

Most organism / antimicrobial combinations give distinct end points. With some combinations there may be a gradual fading of growth over 2 to 3 wells. The end points should be taken as the first well that inhibits visible growth, except with sulphonamides and linezolid when the MIC must be read as an 80-90% decrease in growth compared to

the Positive Control well.

b. Contamination

Contamination may result in growth in a well bordered by wells showing no growth. Such a single well contamination can be ignored, but if multiple well contaminants are suspected, the test should be repeated.

c. Skips

Occasionally a “skip” may be seen - a well showing no growth bordered by wells showing growth. There are a variety of explanations for this including contamination, mutation, creased seal and misaligned dosing. A single skip can be ignored. However, in order to ensure effective antimicrobial therapy, NEVER read the skip well as the MIC; always read the lowest well concentration above which there is consistently no growth.

d. Mixed Cultures

Except as referred to in (a) above, if two end points are seen as a distinct “button” of cells followed by several wells of diffuse growth with the “button” no longer visible (or seen as smaller buttons), there may be a mixed bacterial population. Purity should be checked by sub-culturing growth onto suitable agar. Test results are invalid if a mixed culture is detected.

3. Additional Information

MIC Results for ESBL Confirmation

To confirm ESBLs a ≥ 3 twofold concentration decrease in an MIC for the antimicrobial agent tested in combination with clavulanic acid versus its MIC when tested alone = ESBL. Example: ceftazidime MIC = 8 µg/ml, ceftazidime/clavulanic acid MIC = 1 µg/ml. For more details refer to CLSI M7.

Cefoxitin Screen

The “Cefoxitin Screen (6 µg/ml)” can be used to predict the presence of *mecA*-mediated resistance in *Staphylococcus aureus*. Isolates for which the cefoxitin MIC is > 6 (positive growth) should be reported as oxacillin resistant. Those which have cefoxitin MICs of ≤ 6 (no growth) should be reported as oxacillin susceptible.

D Test

The D test for Broth Microdilution test is for *Staphylococcus* spp. resistant to (MIC's ≥ 8 µg/ml) erythromycin and susceptible or intermediate (MIC's ≤ 2 µg/ml) to clindamycin. With these isolates, a D test 1 (4/0.5 µg/ml) or D Test 2 (8/1.5 µg/ml) or both 1 or 2 are positive (growth), the test should be reported as inducible clindamycin resistance. If with these isolates there is no growth in both D test 1 and D test 2, the isolates should be reported as negative for inducible resistance. If erythromycin is resistant, clindamycin sensitive and the D test is negative after 18 hours incubation, reread after 24 hours incubation.

INTERPRETATION OF RESULTS

For interpretation of results, refer to the MIC interpretive guidelines provided by CLSI, EUCAST, or local standard guidelines.

There may be antimicrobial agents included in the plate that have not been proven to be

effective for treating infections for all organisms tested. Refer to the individual FDA approved pharmaceutical antimicrobial agent package insert for interpreting and reporting results of antimicrobial agents that have shown to be active against organism groups both *in vitro* and in clinical infections.

For automatically read results, the software applies interpretive criteria based on CLSI, EUCAST or FDA guidelines, depending on software configuration

When reporting these ESBL Confirmatory Test antimicrobial results, please refer to the CLSI M100: Screening and Confirmatory Tests for ESBL's in *Klebsiella pneumoniae*, *Klebsiella oxytoca* and *Escherichia coli*, for the appropriate interpretation of results.

Breakpoint Results

Breakpoint testing is a broth dilution method for qualitative susceptibility testing. The Sensititre breakpoint system has been developed to provide a simple standardized method for sensitivity testing based on the concept of breakpoint concentrations. A breakpoint is defined as the concentration of an antibiotic that inhibits the growth of sensitive, but not resistant, organisms. For most antimicrobics, two concentrations are used - a lower concentration, which represents the upper limit of the susceptible category; and a higher concentration, which represents the upper limit of the intermediate category.

Breakpoint interpretations for Antimicrobics present in 1 or 2 well concentrations.

Result	Interpretation One well	Interpretation Two well
No growth	Susceptible	-
Growth	Resistant	-
No growth in both wells	-	Susceptible
Growth in lower concentration only.	-	Intermediate
Growth in both wells	-	Resistant
Growth in higher concentration only.	-	Invalid - Repeat

QUALITY CONTROL

1. The frequency of quality control testing should be set according to laboratory standard guidelines.
2. The inoculum should be cultured onto a suitable medium to check for purity. Test results are invalid if a mixed culture is detected.
3. All Sensititre plates include positive control wells. Tests are invalid unless there is distinct growth in all positive control wells. Some plate formats also include a "negative growth" well. This well may, in some cases, be used for calibration of the AutoReader /Optiread and is not required for manual reading.
4. Several factors contribute to achieving the correct MICs including using 3-5 colonies from a pure 16-24 hour agar plate as well as a standardized inoculum density made using our Mcfarland standard and nephelometer, temperature and broth. In practice, replicate MICs will form a normal distribution with the majority of results lying within one dilution of the modal value.

5. The control organisms listed in the quality control tables are recommended for monitoring the susceptibility test procedure and checking the potency of antimicrobics. Refer to the quality control tables for expected quality control results. The test procedure can be considered satisfactory if the susceptibility results obtained with the control organisms are within the stated ranges. Results should **not** be reported if QC results are outside the stated ranges for a drug/bug combination.
6. The following organism/antimicrobial quality control ranges are different from the CLSI recommended ranges: *E. coli* (ATCC® 25922™) versus cefotetan , *Ps. aeruginosa* (ATCC® 27853™) versus mezlocillin, *S. aureus* (ATCC® 29213™) versus ampicillin, and penicillin.

Contact Thermo Fisher Microbiology Technical Support for assistance in the event that quality control discrepancies cannot be resolved. See page below for contact information.

LIMITATIONS

- Trained clinical personnel are necessary to make proper interpretations of test results.
- In common with all other methods of antimicrobial susceptibility testing, the results generated by Sensititre susceptibility plates are *in vitro* results.
- The testing of fastidious organisms requires specialized procedures and they should be tested by other methods. The Sensititre *Haemophilus influenzae* / *Streptococcus pneumoniae* MIC plate can be used for the testing of *Haemophilus influenzae* and *Streptococcus pneumoniae*. Refer to the Sensititre fastidious technical insert or CLSI M7 for appropriate modifications to susceptibility procedures.
- Studies were performed using the Autoinoculator/ AIM. When testing Minocycline, Doxycycline, Telavancin, Dalbavancin, Oritavancin, or Tedizolid if an alternative inoculation system is used, validate the performance of the alternative inoculation system with Quality Control testing and perform colony counts.

Gram Negative Limitations

1. Results obtained with *C. freundii* and cefotetan have shown discrepant MICs when compared with an overnight reference method, so the antimicrobial agent should not be reported. If the antimicrobial agent is critical to patient care, an alternative procedure should be used
2. Intra-laboratory reproducibility of amoxicillin/clavulanic acid results on Sensititre and a reference procedure was less than 95%. As the MIC is dependent on the level of β -lactamase produced, variability in MIC results may be attributed to variation in enzyme production. Intra-laboratory reproducibility of cefotetan and Gram negative organisms using Sensititre and the CLSI reference procedures was also less than 95 %.
3. Performance has not been established for *Proteus* spp., tested with piperacillin/tazobactam and cefepime. Results should not be reported for patient isolates.
4. *Klebsiella* spp., and non-fermenting organisms tested with piperacillin / tazobactam should be read manually due to lack of performance data on the Sensititre Autoreader / OptiRead system.
5. The Autoreader / OptiRead system should not be used to predict resistance to cefpodoxime with *Serratia* spp., since clinical trials produced an unacceptable error rate when reading this organism /antimicrobial combination.
6. There are no approved primary “Indications for Use” and the clinical significance of cefdinir is unknown. For reporting results which have been shown to be against groups of organisms both in vitro and in clinical infections, refer to CLSI M100 and FDA approved drug labelling.
7. Do not report resistance for the following antimicrobial/organism combinations due to lack

of resistant isolates at the time of testing

Organism	Antimicrobial Agent
<i>Citrobacter</i> species	Piperacillin/tazobactam
<i>Serratia</i> species	Piperacillin/tazobactam
<i>Proteus</i> species	Lomefloxacin and cefpodoxime
Non- <i>Enterobacteriaceae</i>	Levofloxacin

8. The ability of the Sensititre system to detect resistance to moxifloxacin and gatifloxacin is unknown because resistant strains were not available at the time of comparative testing.
9. Tigecycline has decreased *in vitro* activity against *Morganella* spp., *Proteus* spp and *Providencia* spp.
10. The current absence of resistant isolates to tigecycline and gemifloxacin precludes defining any results other than “susceptible”. Isolates yielding MIC results suggestive of “non-susceptible” category should be submitted to a reference laboratory for further testing.
11. Gemifloxacin exhibits *in vitro* MICs of 0.25µg/ml or less against most (≥ 90%) strains of the following micro-organisms: *Acinetobacter lwoffii*, *Klebsiella oxytoca* and *Proteus vulgaris*. However, the safety and effectiveness of gemifloxacin in treating clinical infections due to these micro-organisms has not been established in adequate and well-controlled clinical trials
12. The current absence of resistant isolates to Doripenem precludes defining any results other than “susceptible”. Isolates yielding MIC results suggestive of “Nonsusceptible” category, should be submitted to a reference laboratory for further testing.
13. The ability of the Sensititre system to detect resistance of *Acinetobacter* spp. to minocycline is unknown because resistant strains were not available at the time of comparative testing.
14. Studies were performed using the AutoInoculator / AIM. Use of an alternative inoculation system when testing Ceftaroline has not been evaluated.
15. The ability of the Sensititre system to detect resistance of *Citrobacter koseri* to Ceftaroline is unknown because resistance strains were not available at the time of comparative testing.
16. Ceftaroline is not active against Gram negative bacteria producing extended spectrum beta-lactamases (ESBLs) from the TEM, SHV or CTX-M families, serine carbapenemases (such as KPC), class B metallo-beta-lactamases or class C (AmpC cephalosporins).
17. Performance has not been established for *Moraxella catarrhalis* with Ceftaroline. Results should not be reported for patient isolates.

Gram Positive Limitations

1. The Autoreader / OptiRead system can be used to read coagulase-negative *Staphylococci* and *Enterococci* providing the panel label states ‘For reading gram positives isolates’ and the wells contain fluorogenic substrates
2. The Autoreader / OptiRead should not be used to read nitrofurantoin with *Enterococcus*

- spp. Nitrofurantoin should be read manually. A nitrocefin β -lactamase test should be performed to detect β -lactamase producing strains of *Enterococci*.
3. Sensititre susceptibility plates are configured to meet the CLSI recommendations to detect methicillin-resistant *Staphylococci*. A two per cent salt supplement is therefore included with oxacillin. Inoculum must be prepared directly from an overnight agar plate and not from a fresh broth culture. Be aware that most hetero-resistant *Staphylococci* are usually resistant to multiple antimicrobics including β -lactams, aminoglycosides, macrolides, clindamycin, chloramphenicol and tetracycline and this should be used as a clue to detecting cross-resistance among the penicillinase-resistant penicillins.
 4. *Staphylococci* tested against penicillin G should also be tested for β -lactamase production, especially in strains with borderline MICs (0.06 to 0.25 $\mu\text{g/mL}$). If these strains test β -lactamase positive, a repeat test should be done to ensure resistance is reported.
 5. Poor growth of non-enterococcal strains of *Streptococci* in Mueller-Hinton broth may give unreliable results with aminoglycosides.
 6. Results obtained for coagulase negative oxacillin susceptible *Staphylococci* with clarithromycin have shown discrepant MICs when compared with an overnight reference method. If the antimicrobial agent is critical to patient care, an alternative procedure should be used or the antimicrobial agent should not be reported.
 7. Cefdinir has been shown to be active against most strains of methicillin sensitive *S. aureus* (MSSA) and *Streptococcus pyogenes* in clinical infections. *In vitro* correlation was established for *S. epidermidis* (MSSE) and *Streptococcus agalactiae*. Such activity has not been established for any other organisms appropriate for the Gram positive plate.
 8. The Autoreader / OptiRead system should not be used to read Trimethoprim/sulphamethoxazole (SXT) with CNS. SXT should be read manually
 9. The clinical significance of cefdinir is unknown and there are no approved primary "Indications for Use". For reporting results which have been shown to be against groups of organisms both *in vitro* and in clinical infections, refer to CLSI M100.
 10. The current absence of resistant isolates to Tigecycline and Gemifloxacin precludes defining any results other than "susceptible". Isolates yielding MIC results suggestive of "Non-susceptible" category, should be submitted to a reference laboratory for further testing.
 11. Do not report resistance for the following organism/antimicrobial combinations due to lack of resistant isolates at the time of testing

Organism	Antimicrobial Agent
Gram positive	Linezolid, Moxifloxacin, Gatifloxacin, Tigecycline, Gemifloxacin
<i>Staphylococcus</i> species	Ampicillin/sulbactam
<i>Enterococcus</i> species	Amoxicillin/clavulanic acid
<i>Enterococcus</i> species <i>Staphylococcus</i> species <i>Streptococcus</i> species	Daptomycin

The ability of Sensititre system to detect resistance or non-susceptibility to the following antimicrobial is unknown because an insufficient number of resistant or non-susceptible strains were available at the time of comparative testing. If such a strain is observed, the strain should be sent to a reference laboratory for further testing.

Organism	Antimicrobial Agent
Gram positive	Moxifloxacin, Gatifloxacin,
<i>Staphylococcus aureus</i> (methicillin and erythromycin susceptible only isolates)	Telithromycin
<i>Staphylococcus aureus</i>	Doxycycline
<i>Streptococcus</i> species	Ertapenem
<i>S. pyogenes</i>	Doxycycline
<i>Staphylococcus aureus</i> (both methicillin- resistant and methicillin-susceptible isolates)	Dalbavancin, Oritavancin, Telavancin
<i>Enterococcus faecalis</i> (VSE only)	
<i>Staphylococcus aureus</i> (methicillin-resistant and methicillin-susceptible isolates) <i>Enterococcus faecalis</i>	Tedizolid

12. Sensititre 18-24 hour MIC susceptibility plates are capable of detecting vancomycin resistance in the VRSA *S. aureus* strains available at the time of comparative testing. The ability of the Sensititre 18-24 hour MIC susceptibility plates to detect vancomycin resistance in other *S. aureus* strains is unknown due to the limited number of resistant strains available for comparative testing.
13. Gemifloxacin exhibits *in vitro* MICs of 0.25µg/mL or less against most (≥ 90%) strains of the following micro-organisms: *Staphylococcus aureus* (Methicillin-susceptible strains only) and *Streptococcus pyogenes*. However, the safety and effectiveness of gemifloxacin in treating clinical infections due to these micro-organisms has not been established in adequate and well-controlled clinical trials.
14. Use an alternative testing method with lysed horse blood to test *S. pyogenes* when no growth or poor growth is observed in the positive growth control well.
15. The ability of the Sensititre system to detect resistance of *S. aureus* to minocycline is unknown because resistant strains were not available at the time of comparative testing.
16. Use the Sensititre fastidious plate as a testing method with lysed horse blood to test *Streptococcus spp.* when no growth or poor growth is observed in the positive growth control well on the non-fastidious Sensititre plate.
17. Studies were performed using the Autolnoculator / AIM. When testing telavancin or

minocycline, if an alternative inoculation system is used, validate the performance of the alternative inoculation system with Quality Control testing and perform viable counts.

18. Studies were performed using the AutoInoculator / AIM. When testing Ceftaroline, if an alternative inoculation system is used, validate the performance of the alternate inoculation system with Quality Control testing and perform colony counts.
19. The performance of oritavancin with *S.aureus* and *E. faecalis* was performed using the AIM autoinocultuor. The use of an alternative inoculation system when testing oritavancin has not been evaluated.
20. Due to the lack of intermediate and resistant interpretations for Oritavancin, there is a potential very major error rate. There were 3 isolates out of 9 non susceptible isolates that reported one doubling dilution lower than the reference. Use an alternative testing method prior to reporting results for *S.aureus* and *E. faecalis* with Oritavancin when the Sensititre MIC is 0.12µg/mL (breakpoint) if critical to patient care.

PERFORMANCE

APPENDIX 1: Colony Count Procedure

The 2 methods outlined below can both be used to determine viable counts.

Direct Method

1. Immediately following inoculation of a plate, using a 1µl loop, sample from the positive growth control well and inoculate onto a blood agar.
2. Take another loop (1µl), sample from the same growth well and mix with 50µl sterile deionised water. Inoculate 1µl of this dilution onto a blood agar plate to obtain countable colonies.
3. Incubate both plates at 34–36°C overnight under appropriate conditions.
4. Read as follows:

Intended Inoculum dilution	Colony Count	Number of colonies	
		0.001 plate	0.001 of 1/50
1×10^4	$5 \times 10^3 - 5 \times 10^4$	5-50	0
1×10^5	$5 \times 10^4 - 5 \times 10^5$	50 - 500	1 - 10
5×10^5	$2.5 \times 10^5 - 2.5 \times 10^6$	250 - 2500	5 - 50

In-Direct Method

1. A. For an Inoculum of 1×10^5 cfu/ml

Using a sterile pipette tip, transfer 100µl of broth culture into a fresh 10ml tube of broth and label as tube A.

B. For an Inoculum of 5×10^5 cfu/ml

Using a sterile pipette tip, transfer 10µl of broth culture into a fresh 10ml tube of broth and label as tube A.

2. Using a sterile pipette tip, transfer 100µl from tube A to the surface of a suitable labelled

and dated blood plate.

3. Using a sterile inoculating loop, spread this 100µl of culture across the entire plate surface.
4. Incubate the plate at 34-36°C for 18-24 hours.
5. Count the colonies on the plate and multiply by 10^3 (1×10^5 cfu/ml) or 10^4 (5×10^5 cfu/ml) to obtain the viable count of the original suspension.

For an Inoculum of 1×10^5 cfu/ml	For an Inoculum of 5×10^5 cfu/ml
Viable count range = 5×10^4 - 5×10^5 cfu/ml	Viable count range = 2.5×10^5 - 2.5×10^6 cfu/ml

REFERENCES

1. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7 Approved Standard. The Clinical and Laboratory Standards Institute.
2. Performance Standards for Antimicrobial Susceptibility Testing: Informational Supplement M100. The Clinical and Laboratory Standards Institute.

DISCLAIMER

The information provided in this technical insert is current at the time of printing and may change without notice.

The latest information can be downloaded from the www.TREKDS.com/techinfo or by contacting Thermo Fisher Scientific Microbiology Technical Support.

Legal Manufacturer:
 TREK Diagnostic Systems Ltd
 Units 17-19 Willard Way, Birches Industrial Estate,
 East Grinstead,
 West Sussex.
 RH19 1XZ,
 United Kingdom
 Tel: +44-1342-318777

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. ATCC is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Technical Support email addresses

Austria	microbiology.techsupport.de@thermofisher.com
Belgium	microbiology.techsupport.benelux@thermofisher.com
Czech Republic	mikrobiologie.tech.podpora.cz@thermofisher.com
Denmark	microbiology.techsupport.nordic@thermofisher.com
France	microbiologie.techsupport.fr@thermofisher.com
Finland	mikrobiologia.tekninentuki.fi@thermofisher.com
Germany	microbiology.techsupport.de@thermofisher.com
Italy	microbiologia.supportotecnico.it@thermofisher.com
Netherlands	microbiology.techsupport.benelux@thermofisher.com
Norway	microbiology.techsupport.nordic@thermofisher.com
Spain	microbiologia.soporte.es@thermofisher.com
Sweden	microbiology.techsupport.nordic@thermofisher.com
UK Technical Support	microbiology.techsupport.uk@thermofisher.com

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