# STERILE AIRWAY NEEDLE/SUBCULTURE UNITS

For use with blood culture bottles.

100 Units/20 Gauge

For Laboratory Use.

Disposable: Use once and destroy.

Sterility assured unless individual seal is broken.

## Preliminary comments and precautions

1. Remove culture bottles from instrument as soon as possible after being flagged positive.
2. Place positive culture bottles under a safety hood and wear protective equipment.
3. When venting, general caution should be taken as positive culture bottles could have been overfilled or contain high gas-producing organisms. **Positive culture bottle contents may be under increased internal pressure.**
4. If overfilling or increased gas pressure is suspected (e.g. by bulging bottle septum) place the bottle at room temperature for at least 10 minutes before continuing.

## Directions for Venting Procedure

1. Disinfect stopper of positive culture bottle.
2. Break seal of sterile airway/subculture unit. Remove white cap from unit, exposing short beveled needle. Do not disturb the clear plastic sheath.
3. Always direct the positive culture bottle to be vented away from the user to minimize the risk of exposure should liquid under pressure escape from the bottle.
4. Insert beveled needle through center of disinfected stopper. To ensure adequate venting, partially lift the plastic sheath off the venting needle base (approximately 0.5 cm). Hold for two or three seconds before reseating the sheath.
5. Excess gas will escape from the positive culture bottle.

## Directions for Subculture Procedure

1. If dispersion of the bottle contents is required, keep the previously vented bottle in an upright position and gently swirl the bottle in a circular motion.
2. Aseptically remove the clear plastic sheath, exposing the longer, blunt-ended subculture needle.
3. Tip bottle and place 1-2 drops of broth onto surface of plated medium and 1 drop onto a glass slide for staining.
4. After subculture, carefully remove airway/subculture unit and dispose into sharps container. Please follow your institution's guidelines for proper disposal of sharps. **DO NOT REUSE AIRWAY/SUBCULTURE UNIT**

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To open Airway/Subculture Unit, break seal with a twisting motion.

Dist. by  
bioMérieux, Inc.  
Box 15969  
Durham, North Carolina 27704-0969  
© BIOMÉRIEUX 2004, 2006

January 2006

## INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Caution, consult accompanying documents
	Fragile, handle with care
	Authorized Representative in the European Community
	Do not reuse
	Sterile

**WARRANTY**

*bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.*

bioMérieux® sa  
au capital de 11 879 045 €  
673 620 399 RCS LYON  
69280 Marcy-l'Etoile / France  
Tel. 33 (0)4 78 87 20 00  
Fax 33 (0)4 78 87 20 90  
<http://www.biomerieux.com>

The logo is a registered and protected trademark of bioMérieux sa or one of its subsidiaries.

bioMérieux, Inc  
Box 15969,  
Durham, NC 27704-0969 / USA  
Tel. (1) 919 620 20 00  
Fax (1) 919 620 22 11

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Dist. by  
bioMérieux, Inc.  
Box 15969  
Durham, North Carolina 27704-0969

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## DENOTES REVISED SECTION

### INTENDED USE

BacT/ALERT® FA Plus Culture Bottles are used with the BacT/ALERT® Microbial Detection System in qualitative procedures for recovery and detection of aerobic and facultative anaerobic microorganisms (bacteria and yeast) from blood and other normally sterile body fluids.

### SUMMARY AND EXPLANATION

The BacT/ALERT Microbial Detection System is used to determine if microorganisms are present in blood or other normally sterile body fluid samples taken from a patient suspected of having bacteremia/fungemia. The BacT/ALERT System and culture bottles provide both a microbial detection system and a culture medium with suitable nutritional and environmental conditions for organisms commonly encountered in blood infections and other normally sterile body fluid infections. An inoculated bottle is placed into the instrument where it is incubated and continuously monitored for the presence of microorganisms that will grow in the BacT/ALERT FA Plus Culture Bottle.

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**BacT/ALERT® FA Plus** (color-coded pale green) – BacT/ALERT FA Plus disposable culture bottles contain 30 ml of complex medium and 1.6 g adsorbent polymeric beads. The medium consists of casein peptone (1.0% w/v), yeast extract (0.45% w/v), soybean peptone (0.3% w/v), meat peptone (0.1% w/v), sodium polyanethol sulfonate (SPS) (0.083% w/v), menadione (0.00005% w/v), hemin (0.0005% w/v), L-cysteine (0.03% w/v), pyruvic acid (0.1% w/v), pyridoxine HCl (0.001% w/v), nicotinic acid (0.0002% w/v), pantothenic acid (0.0002% w/v), thiamine HCl (0.0001% w/v), and other complex amino acid and carbohydrate substrates in purified water. Bottles contain an atmosphere of N<sub>2</sub>, O<sub>2</sub>, and CO<sub>2</sub>, under vacuum. The composition of the medium may be adjusted to meet specific performance requirements.

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**CAUTION: BacT/ALERT FA Plus Culture Bottles used to culture non-blood specimens (normally sterile body fluids) or very small blood specimen volumes (0.5 ml or less) will require added blood such as sterile defibrinated horse blood (10.0% v/v) to support growth, particularly for the recovery of fastidious organisms such as *Haemophilus influenzae*, *Streptococcus pneumoniae*, and *Neisseria gonorrhoeae*.<sup>3</sup>**

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### Storage instructions

BacT/ALERT FA Plus Culture Bottles are ready for use. Store in an upright position protected from direct light at room temperature (15-30°C). An expiration date is printed on each bottle label. Do not use the culture bottles beyond the expiration date indicated. If the bottles are exposed to temperatures less than 15°C, precipitates may form that will disappear when the bottles are warmed to room temperature. Bottles must be at room temperature before use.

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Prior to use, the BacT/ALERT FA Plus Culture Bottles should be examined for evidence of damage or deterioration (discoloration). Bottles exhibiting evidence of damage, leakage, or deterioration should be discarded. The medium in undisturbed bottles should be clear, but there may be a slight opalescence or a trace of precipitate due to the anticoagulant SPS or the presence of adsorbent polymeric beads; do not confuse this with turbidity indicative of microbial growth. Do not use a bottle which contains medium exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination.

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### General Considerations

1. BacT/ALERT FA Plus Culture Bottles should be utilized by trained healthcare personnel. Correct specimen collection is extremely important when obtaining blood culture specimens. Refer to Cumitech 1C<sup>4</sup> for the proper specimen collection procedure.
2. Take care to prevent contamination during both bottle preparation and inoculation of the patient sample. Proper skin disinfection is an essential requirement to reduce the incidence of contamination.
3. Although not recommended by bioMérieux, blood may be drawn directly into collection tubes containing SPS. Tubes containing other anticoagulants should never be used for blood culture.<sup>5</sup>

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**Bottle preparation**

1. Label the culture bottle with patient information. The icons on the bottle label (☹, #, ☺) can be defined by the user.
2. Remove plastic flip-top from the culture bottle. Prior to inoculation, disinfect the culture bottle top with an alcohol swab or equivalent. Allow to air dry.
3. Clean the selected venipuncture site as recommended by your institution's approved procedure.

**Direct draw inoculation procedure**

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**CAUTION:** A contaminated culture bottle could contain positive pressure, and if used for direct draw, may cause reflux into the patient's vein. Culture bottle contamination may not be readily apparent. Monitor the direct draw process closely to avoid reflux. Do not use a bottle that contains medium exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination.

**Syringe draw inoculation procedure**

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- a. Perform venipuncture and blood transfer to the BacT/ALERT culture bottle according to your institution's established procedures.
4. Transfer the inoculated culture bottle promptly to the testing laboratory.

**BacT/ALERT FA PLUS CULTURE BOTTLE TEST PROCEDURE****Preliminary comments and precautions**

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**Laboratory procedure**

**CAUTION:** General caution should be taken when subculturing positive culture bottles as they could have been overfilled or contain high gas-producing organisms. Positive culture bottle contents may be under increased internal pressure. Positive culture bottles should be transiently vented before staining or disposal to release any gas produced during microbial metabolism.

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### QUALITY CONTROL

A Certificate of Conformance is available for each lot of culture bottles. If desired, individual laboratories can perform quality control testing of BacT/ALERT FA Plus Culture Bottles. Refer to the BacT/ALERT User Manual and CLSI® document M22-A3.<sup>6</sup>

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1. Patient specimens determined positive by BacT/ALERT may contain organisms that are positive by smear that will not grow on routine subculturing media. When this is suspected, specimens should be subcultured on special media. Also, BacT/ALERT positive specimens may contain organisms that are not seen with routine smear staining methods and may require both specialized staining and subculturing media for detection and recovery.
2. It is possible that certain rare, fastidious microorganisms will not grow or may grow slowly in the BacT/ALERT FA Plus Culture Bottle growth medium. In addition, on rare occasions, organisms may be encountered that grow in the BacT/ALERT FA Plus Culture Bottle growth medium but do not produce sufficient carbon dioxide to be determined positive. If rare, fastidious organisms requiring specialized media and culture conditions are suspected, alternative methods or extended incubation time should be considered for recovery.
3. Certain strains of *Haemophilus influenzae*, *Neisseria meningitidis*, and *Neisseria gonorrhoeae* may be sensitive to the anticoagulant SPS, which may result in lack of growth or low production of CO<sub>2</sub> by these strains if an insufficient amount of sample is inoculated into the culture bottles.
4. Infrequently, if there is a very high number of white blood cells present in the sample, the BacT/ALERT may indicate a culture bottle positive. In this case, the smear and subculture results may be negative.
5. Organisms are often few in numbers and may appear intermittently in the blood stream; therefore, consecutive blood samples should be collected from each patient.

6. Promptly remove positive culture bottles when they are signaled by BacT/ALERT to avoid possible non-viable cultures due to autolysis or other reasons. Certain strains of *Streptococcus pneumoniae* may be particularly prone to autolysis if they are not removed promptly after being signaled positive.
7. A Gram-stained smear from a negative bottle may sometimes contain a small number of non-viable organisms that were derived from culture medium components, staining reagents, immersion oil, or glass slides, resulting in a false positive smear.
8. bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. But, in the unavoidable cases when there is a delay in bottle receipt by the laboratory, delayed entry information is provided from seeded studies in the "Performance Characteristics of the Test" section.
9. Antimicrobial neutralization was not achieved for ceftazidime or cefepime.
10. Extended times to detection or negative results may be observed for *Stenotrophomonas maltophilia* in cerebrospinal fluid for 12 month old bottles.

### EXPECTED VALUES

1. Percent positive cultures were observed to be 12.3% (range: 7.9%-14.3%) overall and 9.4% (range: 5.0%-11.2%) for significant isolates from three clinical trial sites in BacT/ALERT FA Plus Culture Bottles that received 6-10 ml of blood.
2. Percent positive cultures were observed to be 20.2% (range: 16.4%-24.3%) overall and 15.4% (range: 8.2%-21.6%) for significant isolates from three clinical trial sites in BacT/ALERT FA Plus Culture Bottles that received sterile body fluids.
3. Expected percent positives will vary based on factors such as patient population, prevalence of significant organisms, site location, and contamination rates. The expected values provided are based on clinical study data.

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Neutralization of antimicrobials by adsorbent polymeric beads varies depending upon dosage level and timing of specimen collection. Internal studies have demonstrated that antimicrobials are effectively neutralized by the BacT/ALERT FA Plus medium based on 100% recovery of the organisms tested. In these tests, antimicrobials were added in clinically relevant concentrations directly to culture bottles during inoculation with susceptible strains. The effectiveness of the antimicrobials was confirmed by parallel testing using a non-neutralizing medium as a control. Antimicrobials from the following categories were neutralized by the medium: penicillins, glycolcyclines, polyenes, macrolides, triazoles, echinocandins, cefazolin, cefoxitin, ceftaroline, aminoglycosides, fluoroquinolones, lincosamides, glycopeptides, and oxazolidinones.

Antimicrobial neutralization was not achieved for ceftazidime or cefepime. Less than complete neutralization was observed for cefotaxime and ceftriaxone. Cefotaxime was neutralized at ranges of 50% peak serum level (PSL) to 2% PSL depending on the microorganism. Ceftriaxone was neutralized at ranges of 50% PSL to 1% PSL depending on the organism.

For additional information on antimicrobial agents neutralized by BacT/ALERT FA Plus Culture Bottles, contact your local bioMérieux representative.

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Data in Table 1 represent results from in-house seeded studies. A minimum of 30 replicates were tested per species. Data in Table 1 were generated using bottles at the end of shelf life. Bottles inoculated with *H. influenzae* received 4 ml pooled human blood supplementation. At least 95% detection was achieved at LoD.

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**NOTE:** 96.7% of the bottles were subcultured within 30 minutes of being declared positive. STL 104016 was sourced from bioMérieux's internal culture collection.

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Data in Table 2 represent results from in-house seeded studies with and without blood obtained from healthy human volunteers. Multiple strains were tested for each species at target inoculum levels of 125 CFU/bottle. The actual inoculum levels ranged from 3-298 CFU/bottle. In this seeded study, the BacT/ALERT FA Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. The species listed are representatives of clinically prevalent organisms in blood cultures and sterile body fluids.

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Less than 100% detection was observed for some species, including *Capnocytophaga ochracea*, *Cardiobacterium hominis*, *Eikenella corrodens*, *Haemophilus parainfluenzae*, *Granulicatella adiacens*, and *Helicobacter cinaedi*.

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Results compare BacT/ALERT FA Plus to BacT/ALERT FA blood cultures (for all compliant pairs).

A multi-center clinical study was conducted at three different geographic sites in the U.S. comparing the performance of the BacT/ALERT FA Plus and BacT/ALERT FA blood culture bottles for aerobic culture pairs that received blood volumes between 6 ml and 10 ml (compliant pairs). A total of 1656 bottle pairs were obtained from 728 adult patients suspected of blood stream bacterial/yeast infections. Subcultures of both bottles were performed when either bottle in the set was determined to be positive by the BacT/ALERT System. A pair of bottles was determined to have a positive status if the subculture of either the BacT/ALERT FA Plus or BacT/ALERT FA Culture Bottle was positive. A culture bottle was determined to be a "True Positive" if the culture was flagged positive by the BacT/ALERT System and resulted in growth of the isolate upon subculture of this bottle. True positive rates were calculated for the BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles, and the ratio of BacT/ALERT FA Plus true positives to BacT/ALERT FA true positives was calculated to compare performance. Clinical isolates recovered were classified as significant, contaminant, or unknown based on determination by the clinical trial sites.

A total of 267 isolates were recovered from all compliant aerobic blood culture pairs with a positive status. There were a total of 238 bottle pairs that recovered at least one isolate by subculture of BacT/ALERT FA Plus or BacT/ALERT FA Culture Bottles. A total of 214 bottle pairs recovered a single isolate, 19 bottle pairs recovered two isolates, and five bottle pairs recovered three isolates. The total population reported in Table 3 comprises the 267 isolates recovered from positive bottle pairs and 1418 negative bottle pairs for a total of 1685 results. The BacT/ALERT FA Plus Culture Bottle detected a total of 208 isolates compared to the BacT/ALERT FA Culture Bottle that detected 194 isolates. Of the significant isolates, the BacT/ALERT FA Plus Culture Bottle detected a total of 159 isolates compared to the BacT/ALERT FA Culture Bottle that detected 135 isolates. Five false positives were identified by subculture of positive BacT/ALERT FA Plus Culture Bottles and comprised 0.30% (5/1685) of the study population.

The following tables compare results of the BacT/ALERT FA Plus to BacT/ALERT FA blood cultures for all compliant blood culture bottles that yielded any number of isolates on subculture (Table 3), a single isolate alone on subculture (Table 4), and multiple isolates on subculture (Table 5).

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\*One hundred thirty five (135) isolates were detected by both BacT/ALERT FA Plus and BacT/ALERT FA, 73 isolates were detected only by BacT/ALERT FA Plus, and 59 isolates were detected only by BacT/ALERT FA. The ratio of true positive rates for overall isolates was 1.072 (208/194) with a 95% CI (0.952, 1.193).<sup>7</sup>

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\*One hundred eleven (111) isolates were detected by both BacT/ALERT FA Plus and BacT/ALERT FA, 61 isolates were detected only by BacT/ALERT FA Plus, and 42 isolates were detected only by BacT/ALERT FA. The ratio of true positive rates for overall single isolates was 1.124 (172/153) with a 95% CI (0.986, 1.262).<sup>7</sup>

**Table 5. All Compliant Pairs With Multiple Isolates (Blood Cultures)**

Clinical Determination	BacT/ALERT FA Plus True Positives	BacT/ALERT FA True Positives	Ratio of True Positives*
Significant	21	24	0.875
Contaminant	10	14	0.714
Unknown	5	3	1.667
Total	36	41	0.878

\*Twenty four (24) isolates were detected by both BacT/ALERT FA Plus and BacT/ALERT FA, 12 isolates were detected only by BacT/ALERT FA Plus, and 17 isolates were detected only by BacT/ALERT FA. The ratio of true positive rates for overall multiple isolates was 0.878 (36/41) with a 95% CI (0.637, 1.119).<sup>7</sup>

In this clinical study, there were 1413 pairs of BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles with negative instrument results for both bottles after 5 days of incubation. Among these pairs, terminal subcultures on both bottles were performed for 95 pairs, and two false negative results by both BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles were observed; subculture on BacT/ALERT FA Plus Culture Bottles alone was performed for 1312 pairs, and one false negative result was observed; both subcultures were not performed for six pairs of bottles. Results are summarized in Table 6.

**Table 6. Summary Of Percent False Negatives From Compliant Aerobic Blood Culture Pairs That Were Flagged Negative By The Instrument For Both Bottles**

Subculture Performed BacT/ALERT FA Plus	Subculture Performed BacT/ALERT FA	% False Negative BacT/ALERT FA Plus	% False Negative BacT/ALERT FA
Yes	Yes	2.11 (2/95)	2.11 (2/95)
Yes	No	0.08 (1/1312)	-

Overall false negative rate for BacT/ALERT FA Plus based on a subset of terminal subcultures was 0.2% (3/1407).

A comparative yield of microorganisms (number of isolates) recovered on subculture of BacT/ALERT FA Plus and BacT/ALERT FA cultures is presented in Table 7.

**Table 7. Comparative Yield Of Microorganisms From BacT/ALERT FA Plus And BacT/ALERT FA (Number Of Isolates) – Blood Cultures**

Group	BacT/ALERT FA Plus	BacT/ALERT FA
Enterobacteriaceae	34	32
<i>Enterococcus</i> spp.	33	28
Yeasts	17	18
Non-fermentative Gram-Negative Bacilli	8	11
Other Gram-Negative	4	3
Other Gram-Positive	6	8
Coagulase-Negative <i>Staphylococcus</i>	37	53
<i>Staphylococcus aureus</i>	52	30
<i>Streptococcus</i> spp.	17	11

NOTE: Isolate table includes polymicrobial cultures.

#### Clinical Study Results (Sterile Body Fluid Cultures)

A multi-center clinical study was conducted at four different geographic sites in the U.S. and Canada comparing the performance of the BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles with sterile body fluid specimens. A total of 404 bottle pairs were obtained from 369 adult patients suspected of sterile body fluid bacterial/yeast infections. Sterile body fluid types evaluated were amniotic fluid, continuous ambulatory peritoneal dialysis (CAPD) fluid, cerebrospinal fluid (CSF), peritoneal fluid, pleural fluid, and synovial fluid. Clinical isolates recovered were classified as significant, contaminant, or unknown based on determination by the clinical trial sites.

A total of 92 isolates were recovered from all aerobic sterile body fluid culture pairs with a positive status. There were a total of 75 bottle pairs that recovered at least one isolate by subculture of BacT/ALERT FA Plus or BacT/ALERT FA Culture Bottles. A total of 62 bottle pairs recovered a single isolate, nine bottle pairs recovered two isolates, and four bottle pairs recovered three isolates. The total population reported in Table 8 comprises the 92 isolates recovered from positive bottle pairs and 329 negative bottle pairs for a total of 421 results. The BacT/ALERT FA Plus Culture Bottle detected a total of 85 isolates compared to the BacT/ALERT FA Culture Bottle that detected 67 isolates. Of the significant isolates, the BacT/ALERT FA Plus Culture Bottle detected a total of 65 isolates compared to the BacT/ALERT FA Culture Bottle that detected 59 isolates. No false positives were observed for the BacT/ALERT FA Plus Culture Bottle from the study population (0/421).

Table 8 compares results of the BacT/ALERT FA Plus to BacT/ALERT FA sterile body fluid cultures that yielded single or multiple isolates on subculture.

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\*Sixty (60) isolates were detected by both BacT/ALERT FA Plus and BacT/ALERT FA, 25 isolates were detected only by BacT/ALERT FA Plus, and seven isolates were detected by only BacT/ALERT FA. The ratio of true positive rates for overall isolates was 1.269 (85/67) with a 95% CI (1.083, 1.455).<sup>7</sup>

NOTE: A limited number of amniotic fluid (n=2) and cerebrospinal fluid specimens (n=38) were obtained during the clinical study.

Table 9 summarizes the minimum specimen volume achieved in aerobic sterile body fluid clinical trials.

**Table 9. BacT/ALERT FA Plus Sterile Body Fluids Fill Volume (ml) – Positive Status**

Specimen Type	Total # of Specimen	# of Positives	Minimum Specimen Volume (ml)
Amniotic Fluid	2	1	1.0
CAPD Fluid	94	26	1.0
CSF	38	4	0.1
Peritoneal Fluid	116	19	1.0
Pleural Fluid	106	17	0.5
Synovial Fluid	48	8	0.5
Total	404	75	-

In this clinical study, there were 329 pairs of BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles with negative instrument results for both bottles after 5 days of incubation. Among these pairs, terminal subcultures on both bottles were performed for 297 pairs and no false negative results by both BacT/ALERT FA Plus and BacT/ALERT FA Culture Bottles were observed. Subculture on BacT/ALERT FA Plus Culture Bottles alone was performed for 32 pairs and no false negative results were observed. Results are summarized in Table 10.

**Table 10. Summary Of Percent False Negatives From Aerobic Sterile Body Fluid Culture Pairs That Were Flagged Negative By The Instrument For Both Bottles**

Subculture Performed BacT/ALERT FA Plus	Subculture Performed BacT/ALERT FA	% False Negative BacT/ALERT FA Plus	% False Negative BacT/ALERT FA
Yes	Yes	0.0 (0/297)	0.0 (0/297)
Yes	No	0.0 (0/32)	-

A comparative yield of microorganisms (number of isolates) recovered on subculture of BacT/ALERT FA Plus and BacT/ALERT FA cultures is presented in Table 11.

**Table 11. Comparative Yield Of Microorganisms From BacT/ALERT FA Plus And BacT/ALERT FA (Number Of Isolates) – Sterile Body Fluid Cultures**

Group	BacT/ALERT FA Plus	BacT/ALERT FA
Enterobacteriaceae	8	7
<i>Enterococcus</i> spp.	12	10
Yeasts	11	11
Non-fermentative Gram-Negative Bacilli	8	3
Other Gram-Positive	5	3
Other Gram-Negative	-	-
Coagulase-Negative <i>Staphylococcus</i>	21	19
<i>Staphylococcus aureus</i>	12	9
<i>Streptococcus</i> spp.	8	5

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Quality control was performed during the clinical study on each of 13 organisms (*Candida albicans*, *Candida krusei*, *Enterococcus faecalis*, *Escherichia coli*, *Haemophilus influenzae*, *Neisseria meningitidis*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus agalactiae*, *Streptococcus pneumoniae*, and *Streptococcus pyogenes*), which were prepared using serial dilution and seeded into the BacT/ALERT FA Plus Culture Bottle at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle. Overall quality control results were found to be acceptable. Instances where unacceptable quality control results were observed were found to be due to technical errors (i.e., colony counts out of range, contaminated and mislabeled bottles). Repeat testing resulted in acceptable results.

#### Delayed Entry

Table 12 includes results from seeded studies using 11 species (*Staphylococcus aureus*, *Candida albicans*, *Candida krusei*, *Escherichia coli*, *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Streptococcus pneumoniae*, *Enterococcus faecium*, *Haemophilus influenzae*, and *Neisseria meningitidis*), at target concentrations 100 CFU/bottle (acceptable range of 30-300 CFU/bottle) were generated at three sites. Actual inoculum levels ranged from 35-290 CFU/bottle. All bottles contained human blood from healthy volunteers and were held at specified temperatures and times prior to loading into the BacT/ALERT instrument. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 12. Delayed Entry**

Sample Input	Incubation Temperature (°C)	Hold Time (hours)	% Recovery	Time to Detection from Sample Inoculation (Hold Time + Instrument TTD in hours)	
				Mean	Range
Inoculated Test Bottles	Control	No delay	100.0 (459/459)	14.3	8.5-84.0
	2-8	48	98.6 (292/296)	63.7	57.5-103.2
	20-25	24	98.0 (291/297)	31.8	26.2-74.4
	20-25	36	91.9 (272/296)	41.8	38.0-70.5
	35-37	8	98.9 (454/459)	16.1	10.2-53.8
	35-37	24	56.6 (259/458)*	28.3	26.0-74.4
Negative Controls	All conditions		0.5 (1/221)**	-	-

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\*\*False positive observed during seeded study (1/221).

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Data in Table 13 represent results from in-house seeded studies conducted on 12 days on multiple instruments by multiple operators. Organisms were grown in the presence of clinically relevant concentrations of antimicrobials to which they are susceptible. In this seeded study BacT/ALERT FA Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. A minimum of 108 replicates were tested for each organism/antimicrobial combination.

**Table 13. Within-Laboratory Precision (Repeatability)**

Sample Input		Range (CFU/bottle)	% Recovery				Time to Detection (hours)	
Organism	Antimicrobial		Lot 1	Lot 2	Lot 3	Overall	Mean	Range
<i>C. albicans</i>	Fluconazole	140-364	100.0	100.0	100.0	100.0	26.0	22.8-31.3
<i>E. coli</i>	Amikacin	26-156	100.0	100.0	100.0	100.0	12.0	11.2-13.0
<i>K. pneumoniae</i>	Levofloxacin	108-170	100.0	100.0	100.0	100.0	13.4	11.7-15.2
<i>P. aeruginosa</i>	Piperacillin	80-148	100.0	97.2	100.0	99.1	19.2	17.4-24.1
<i>S. pneumoniae</i>	Penicillin G	9-505	100.0	100.0	100.0	100.0	13.2	11.6-15.5
<i>S. aureus</i>	Vancomycin	94-158	100.0	100.0	100.0	100.0	16.9	14.6-20.3

**Reproducibility**

Data in Table 14 represent results from seeded studies conducted at three sites using a target of 162 replicates per site on 3 days with a minimum of two operators per site. Reproducibility was evaluated on each of nine organisms. Two organisms (*C. albicans* and *S. pneumoniae*) were prepared using serial dilution and the other seven organisms were prepared using BioBall® products. *C. albicans* and *S. pneumoniae* were seeded into the BacT/ALERT FA Plus Culture Bottle, at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle and the other seven organisms at a target range of 1-17 CFU/bottle. The actual inoculum ranged from 6-700 CFU/bottle for the 30-300 CFU/bottle range, and from 1-270 CFU/bottle for the 1-17 CFU/bottle range. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 14. Reproducibility**

Sample Input	% Recovery				Time to Detection (hours)		Inoculum Ranges (CFU/Bottle)
	Site 1	Site 2	Site 3	Overall	Mean	Range	
<i>S. aureus</i>	100.0% (18/18)	87.5% (21/24)	100.0% (30/30)	95.8% (69/72)	15.6	14.6-16.7	2-11
<i>C. albicans</i>	100.0% (18/18)	83.3% (30/36)	100.0% (33/33)	93.1% (81/87)	36.6	24.6-76.8	14-700
<i>E. coli</i>	100.0% (27/27)	77.8% (21/27)	100.0% (30/30)	92.9% (78/84)	12.8	11.8-14.1	1-38
<i>P. aeruginosa</i>	100.0% (24/24)	75.0% (18/24)	97.0% (32/33)	91.4% (74/81)	18.4	17.1-21.1	1-11
<i>E. faecalis</i>	100.0% (18/18)	79.2% (19/24)	96.7% (29/30)	91.7% (66/72)	13.9	12.6-15.3	1-15
<i>E. aerogenes</i>	74.4% (29/39)	72.2% (26/36)	85.4% (41/48)	78.1% (96/123)	14.9	11.7-20.8	1-270*
<i>L. monocytogenes</i>	100.0% (18/18)	100.0% (24/24)	100.0% (30/30)	100.0% (72/72)	24.1	20.4-36.4	1-14
<i>S. enterica</i>	100.0% (24/24)	75.0% (18/24)	100.0% (33/33)	92.6% (75/81)	13.5	2.3-14.8	1-13
<i>S. pneumoniae</i>	100.0% (30/30)	100.0% (36/36)	100.0% (21/21)	100.0% (87/87)	14.2	11.6-18.9	6-500
Overall	95.4% (206/216) 95% CI: 91.7%, 97.8%	83.5% (213/255) 95% CI: 78.4%, 87.9%	96.9% (279/288) 95% CI: 94.2%, 98.6%	92.0% (698/759) 95% CI: 89.8%, 93.8%	-		

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These data include repeat testing performed as a result of laboratory errors at a single site (i.e., contaminated bottles/reagents, colony counts out of range, and site failure to change bottle status after positive instrument signal and positive subculture). Data excluding the laboratory errors demonstrated 100% recovery with the exception of *E. aerogenes*, which exhibited 85% recovery for all sites combined.

## REFERENCES

1. Thorpe TC, Wilson ML, Turner JE, DiGuseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. Principles and Procedures for Blood Cultures; Approved Guideline. CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI. Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.
7. Kondratovich M. Comparing Two Medical Tests When Results of Reference Standard Are Unavailable for Those Negative via Both Tests. *J. Biopharm Stat* 2008;18:1,145-166.

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## INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Authorized representative in the European Community
	This way up
	In Vitro Diagnostic Medical Device
	Do not reuse
	Latex-free

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## REVISION TABLE

This section contains a summary of changes made to the BacT/ALERT FA Plus Instructions for Use (9305048 C).

Revision Date	Revision Number	Change Type	Change Summary
2013-04	9305048 C	Requirement:	<b>Intended Use, Specimen Collection and Preparation, Expected Values, Neutralization of Antimicrobials, Performance Characteristics of the Test</b> – Revised text to include additional information on product performance based on clinical studies
		Administrative:	<b>Limitations of the Test</b> - Added Limitations 2, 8, 9, and 10 following FDA review <b>Reagents</b> - Moved Cautions to Limitations of the Test section

**NOTE:** Minor typographical, grammar, and formatting changes are not included in the revision history.

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**LITERATŪRA**

1. Thorpe TC, Wilson ML, Turner JE, DiGuiseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. Principles and Procedures for Blood Cultures; Approved Guideline. CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI. Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.
7. Kondratovich M. Comparing Two Medical Tests When Results of Reference Standard Are Unavailable for Those Negative via Both Tests. *J. Biopharm Stat* 2008;18:1,145-166.

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## DENOTES REVISED SECTION

### INTENDED USE

BacT/ALERT® FN Plus Culture Bottles are used with the BacT/ALERT® Microbial Detection System in qualitative procedures for recovery and detection of anaerobic and facultative anaerobic microorganisms from blood and other normally sterile body fluids.

### SUMMARY AND EXPLANATION

The BacT/ALERT Microbial Detection System is used to determine if microorganisms are present in blood or other normally sterile body fluid samples taken from a patient suspected of having septicemia. The BacT/ALERT System and culture bottles provide both a microbial detection system and a culture medium with suitable nutritional and environmental conditions for organisms commonly encountered in blood infections and other normally sterile body fluid infections. An inoculated bottle is placed into the instrument where it is incubated and continuously monitored for the presence of microorganisms that will grow in the BacT/ALERT FN Plus Culture Bottle.

### PRINCIPLE OF THE TEST

The BacT/ALERT Microbial Detection System utilizes a colorimetric sensor and reflected light to monitor the presence and production of carbon dioxide (CO<sub>2</sub>) dissolved in the culture medium. If microorganisms are present in the test sample, carbon dioxide is produced as the organisms metabolize the substrates in the culture medium. When growth of the microorganisms produces CO<sub>2</sub>, the color of the gas-permeable sensor installed in the bottom of each culture bottle changes from blue-green to yellow.<sup>1</sup> The lighter color results in an increase of reflectance units monitored by the system. Bottle reflectance is monitored and recorded by the instrument every 10 minutes.

### REAGENTS

For *in vitro* diagnostic use.

**CAUTION: Handle specimens and inoculated culture bottles as though capable of transmitting infectious agents. All inoculated culture bottles, specimen collection needles, and blood-drawing devices should be decontaminated according to your institution's procedures.<sup>2</sup>**

**BacT/ALERT® FN Plus** (color-coded orange) – BacT/ALERT FN Plus disposable culture bottles contain 40 ml of complex medium and 1.6 g adsorbent polymeric beads. The medium consists of a combination of peptones (1.48% w/v), yeast extract (0.5% w/v), sodium polyanethol sulfonate (SPS) (0.083% w/v), menadione (0.00005% w/v), hemin (0.001% w/v), pyridoxine HCl (0.0008% w/v), pyruvic acid (0.1% w/v), reducing agents (0.38% w/v), and other complex amino acid and carbohydrate substrates in purified water. Bottles contain an atmosphere of N<sub>2</sub> and CO<sub>2</sub>, under vacuum. The composition of the medium may be adjusted to meet specific performance requirements.

**CAUTION: BacT/ALERT culture bottles contain polycarbonate. Not all disinfectants are intended for use with polycarbonate surfaces, and may cause bottle deterioration. Verify disinfectant compatibility with polycarbonate before use on BacT/ALERT culture bottle surfaces.**

**CAUTION: BacT/ALERT FN Plus Culture Bottles used to culture non-blood specimens (normally sterile body fluids) or very small blood specimen volumes (0.5 ml or less) will require added blood such as sterile defibrinated horse blood (10.0% v/v) to support growth, particularly for the recovery of fastidious organisms such as *Haemophilus influenzae* and *Streptococcus pneumoniae*.<sup>3</sup>**

### Additional materials required

BacT/ALERT® Microbial Detection Systems

Blood-drawing device

Sterile Airway Needles/Subculture Units

Disposable gloves

Appropriate biohazard waste containers for materials potentially contaminated with infectious agents

### Materials available from bioMérieux

Blood Collection Adapter Cap

BacT/ALERT® Microbial Detection Systems

Sterile Airway Needles/Subculture Units

### Storage instructions

BacT/ALERT FN Plus Culture Bottles are ready for use. Store in an upright position protected from direct light at room temperature (15-30°C). An expiration date is printed on each bottle label. Do not use the culture bottles beyond the expiration date indicated. If the bottles are exposed to temperatures less than 15°C, precipitates may form that will disappear when the bottles are warmed to room temperature. Bottles must be at room temperature before use.

### Chemical or physical indications of instability

Prior to use, the BacT/ALERT FN Plus Culture Bottles should be examined for evidence of damage or deterioration (discoloration). Bottles exhibiting evidence of damage, leakage, or deterioration should be discarded. The medium in undisturbed bottles should be clear, but there may be a slight opalescence or a trace of precipitate due to the anticoagulant SPS or the presence of adsorbent polymeric beads; do not confuse this with turbidity indicative of microbial growth. Do not use a bottle which contains medium exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination.

### INSTRUMENTS

Review the appropriate BacT/ALERT Microbial Detection System User Manual before use.

## SPECIMEN COLLECTION AND PREPARATION

### General Considerations

1. BacT/ALERT FN Plus Culture Bottles should be utilized by trained healthcare personnel. Correct specimen collection is extremely important when obtaining blood culture specimens. Refer to Cumitech 1C<sup>4</sup> for the proper specimen collection procedure.
2. Take care to prevent contamination during both bottle preparation and inoculation of the patient sample. Proper skin disinfection is an essential requirement to reduce the incidence of contamination.
3. Although not recommended by bioMérieux, blood may be drawn directly into collection tubes containing SPS. Tubes containing other anticoagulants should never be used for blood culture.<sup>5</sup>

4. bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. If there is an unavoidable delay, inoculated bottles may be maintained at room temperature up to 24 hours before loading into the instrument.
5. Optimal recovery of isolates will be achieved by adding maximum amounts of specimen. Use of lower volumes may adversely affect recovery and/or detection times of some organisms. Do not fill above the bottle's maximum specimen volume of 10 ml. The vacuum in the bottle will usually exceed 10 ml; monitor the volume collected by means of the 5 ml incremental markings on the bottle label.

#### Bottle preparation

1. Label the culture bottle with patient information. The icons on the bottle label (☹, #, ☺) can be defined by the user.
2. Remove plastic flip-top from the culture bottle. Prior to inoculation, disinfect the culture bottle top with an alcohol swab or equivalent. Allow to air dry.
3. Clean the selected venipuncture site as recommended by your institution's approved procedure.
4. For inoculation, use needles with a gauge size of 21 or higher [for example 22, 23, or 25]. Use of needles with less than a 21 gauge size [for example 16, 18, 19, or 20] may vent the bottle's atmosphere due to larger bore size.

#### Direct draw inoculation procedure

**NOTE:** If inoculating more than one type of BacT/ALERT blood culture bottle using a butterfly blood collection set and direct draw adapter cap, inoculate first the aerobic culture bottle and then the anaerobic culture bottle so that any oxygen trapped in the tubing will not be transferred to the anaerobic bottle.

**NOTE:** Monitor the direct draw process closely at all times during collection to assure proper flow is obtained and to avoid flow of the bottle contents into the adapter tubing. Due to the presence of chemical additives in the culture bottle, it is important to prevent possible backflow and subsequent adverse reactions by following all steps below.

- a. Hold the culture bottle at a position below the patient's arm with the bottle in an upright position (stopper uppermost).
- b. Collect the blood using a butterfly blood collection set and a blood collection adapter cap as recommended by your institution's approved procedure and inoculate directly into the culture bottle at the patient's bedside. Although lower sample volumes can be used, recovery may be improved using a sample volume closer to the recommended 10 ml. To prevent over inoculation, monitor the blood volume intake into the culture bottle, using the 5 ml incremental markings on the bottle label.
- c. Release the tourniquet as soon as the blood starts to flow into the culture bottle, or within 2 minutes of application.
- d. Do not allow the culture bottle contents to touch the stopper or the end of the needle during the collection procedure.

**CAUTION: A contaminated culture bottle could contain positive pressure, and if used for direct draw, may cause reflux into the patient's vein. Culture bottle contamination may not be readily apparent. Monitor the direct draw process closely to avoid reflux. Do not use a bottle that contains medium exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination.**

#### Syringe draw inoculation procedure

**NOTE:** If inoculating more than one type of BacT/ALERT blood culture bottle using syringe draw, inoculate first the anaerobic culture bottle and then the aerobic culture bottle so that any oxygen trapped in the syringe will not be transferred to the anaerobic bottle. Line demarcations on the bottle label should be used to assist in estimating the sample volume.

- a. Perform venipuncture and blood transfer to the BacT/ALERT culture bottle according to your institution's established procedures.
5. Transfer the inoculated culture bottle promptly to the testing laboratory.

### BacT/ALERT FN PLUS CULTURE BOTTLE TEST PROCEDURE

#### Preliminary comments and precautions

1. Use disposable gloves and handle inoculated bottles cautiously as though capable of transmitting infectious agents. Consult a physician immediately if contaminated materials are ingested or come in contact with open lacerations, lesions, or other breaks in skin.
2. Immediately clean up any spillage of contaminated material using a 1:10 dilution of 5% sodium hypochlorite. Dispose of the cleaning material by an acceptable method.
3. All inoculated culture bottles, specimen collection needles, and blood-drawing devices should be decontaminated according to your institution's procedures.<sup>2</sup>
4. These bottles should be utilized by trained healthcare personnel.

#### Procedural notes and precautions

1. Great care must be taken to prevent contamination of the patient sample during venipuncture and during inoculation into the culture bottle since contamination could lead to a specimen being determined positive when a clinically relevant isolate is not actually present.
2. Obtain blood samples prior to initiating antibiotic therapy. If this is not possible, draw blood immediately before administering the next antibiotic dose.
3. If inoculated culture bottles have been delayed in their receipt into the laboratory or have been incubated prior to entry into the BacT/ALERT instrument, visually inspect for indications of microbial growth. If microbial growth is evident, treat the bottles as positive and do not place in the BacT/ALERT Microbial Detection System for monitoring.

#### Laboratory procedure

**CAUTION: General caution should be taken when subculturing positive culture bottles as they could have been overfilled or contain high gas-producing organisms. Positive culture bottle contents may be under increased internal pressure. Positive culture bottles should be transiently vented before staining or disposal to release any gas produced during microbial metabolism.**

1. Visually inspect bottles before testing. Do not use bottles with evidence of damage, leakage, or deterioration. Consider bottles with hemolysis, turbidity, excess gas pressure, yellow sensors, and/or evidence of growth as positive. Smear and subculture. Do not incubate unless smear is negative.
2. After culture bottles have been loaded into the instrument, incubate 5 days or until designated positive.
3. Smear and subculture all positive bottles. If the smear is negative, indicating a possible false positive, the bottle should be reloaded into the instrument until growth of the subculture or redesignation as positive. Bottles that were initially determined false positive and were redesignated positive should be smeared and subcultured.

4. Negative bottles may be checked by smear and/or subculture at some point prior to discarding as negative.
5. Procedures for loading and unloading culture bottles into the BacT/ALERT instrument are given in the User Manual.
6. **Do not reuse BacT/ALERT culture bottles.** Dispose of inoculated BacT/ALERT culture bottles according to your laboratory protocol. Autoclaving and/or incinerating inoculated BacT/ALERT culture bottles is appropriate.<sup>2</sup>
7. Utilization of coring devices (i.e., blunt needle) to puncture the septum may result in bottle leakage.

### QUALITY CONTROL

A Certificate of Conformance is available for each lot of culture bottles. If desired, individual laboratories can perform quality control testing of BacT/ALERT FN Plus Culture Bottles. Refer to the BacT/ALERT User Manual and CLSI® document M22-A3.<sup>6</sup>

### Instrument

A BacT/ALERT Reflectance Standards kit is provided with each instrument for the QC and Calibration procedures. All quality control should be part of normal system maintenance. Refer to the BacT/ALERT User Manual for more information.

### RESULTS

Positive or negative culture bottles are determined by decision-making software contained in the BacT/ALERT Microbial Detection System. No action is required until the BacT/ALERT instrument signals a culture bottle either positive or negative.

### LIMITATIONS OF THE TEST

Many variables involved in blood culture testing cannot be practically controlled to provide total confidence that results obtained are due solely to proper or improper performance of any culture medium or detection system.

1. Patient specimens determined positive by BacT/ALERT may contain organisms that are positive by smear that will not grow on routine subculturing media. When this is suspected, specimens should be subcultured on special media. Also, BacT/ALERT positive specimens may contain organisms that are not seen with routine smear staining methods and may require both specialized staining and subculturing media for detection and recovery.
2. It is possible that certain rare, fastidious microorganisms will not grow or may grow slowly in the BacT/ALERT FN Plus Culture Bottle growth medium. In addition, on rare occasions, organisms may be encountered that grow in the BacT/ALERT FN Plus Culture Bottle growth medium but do not produce sufficient carbon dioxide to be determined positive. If rare, fastidious organisms requiring specialized media and culture conditions are suspected, alternative methods or extended incubation time should be considered for recovery.
3. Certain strains of *Peptostreptococcus anaerobius* may be sensitive to the anticoagulant SPS, which may result in a lack of growth or low production of CO<sub>2</sub> by these strains if an insufficient amount of sample is inoculated into the culture bottles.
4. Infrequently, if there is a very high number of white blood cells present in the sample, the BacT/ALERT may indicate a culture bottle positive. In this case, the smear and subculture results may be negative.
5. Organisms are often few in numbers and may appear intermittently in the blood stream; therefore, consecutive blood samples should be collected from each patient.

6. Promptly remove positive culture bottles when they are signaled by BacT/ALERT to avoid possible non-viable cultures due to autolysis or other reasons. Certain strains of *Streptococcus pneumoniae* may be particularly prone to autolysis if they are not removed promptly after being signaled positive.
7. A Gram-stained smear from a negative bottle may sometimes contain a small number of non-viable organisms that were derived from culture medium components, staining reagents, immersion oil, or glass slides, resulting in a false positive smear.
8. bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. But, in the unavoidable cases when there is a delay in bottle receipt by the laboratory, delayed entry information is provided from seeded studies in the "Performance Characteristics of the Test" section.
9. Antimicrobial neutralization was not achieved for ceftazidime, ceftriaxone, or cefepime.

### EXPECTED VALUES

1. Percent positive cultures were observed to be 11.1% (range: 8.7%-14.3%) overall and 8.0% (range: 6.4%-9.3%) for significant isolates from three clinical trial sites in BacT/ALERT FN Plus Culture Bottles that received 6-10 ml of blood.
2. Percent positive cultures were observed to be 20.3% (range: 17.0%-25.3%) overall and 14.6% (range: 6.1%-24.1%) for significant isolates from three clinical trial sites in BacT/ALERT FN Plus Culture Bottles that received sterile body fluids.
3. Expected percent positives will vary based on factors such as patient population, prevalence of significant organisms, site location, and contamination rates. The expected values provided are based on clinical study data.

### NEUTRALIZATION OF ANTIMICROBIALS

Neutralization of antimicrobials by adsorbent polymeric beads varies depending upon dosage level and timing of specimen collection. Internal studies have demonstrated that antimicrobials are effectively neutralized by the BacT/ALERT FN Plus medium based on 100% recovery of the organisms tested. In these tests, antimicrobials were added in clinically relevant concentrations directly to culture bottles during inoculation with susceptible strains of obligate and facultative anaerobic microorganisms. The effectiveness of the antimicrobials was confirmed by parallel testing using a non-neutralizing medium as a control. Studies demonstrated that the following antimicrobials/antimicrobial categories were neutralized by the medium: imipenem, meropenem, oxacillin, glycolylcyclines, macrolides, cefoxitin, ceftaroline, aminoglycosides, fluoroquinolones, lincosamides, ketolides, and glycopeptides.

Antimicrobial neutralization was not achieved for ceftazidime, ceftriaxone, or cefepime. Less than complete neutralization was observed for cefotaxime, cefazolin, ampicillin, penicillin, and ertapenem. Cefotaxime was neutralized at ranges of 40% peak serum level (PSL) to 3% PSL depending on the microorganism. Cefazolin was neutralized at ranges of 25% PSL to 5% PSL depending on the microorganism. Ertapenem was neutralized at 5% PSL. Ampicillin was neutralized at 75% PSL for *E. faecalis*. Penicillin was neutralized at 120% PSL for *S. pneumoniae*. No neutralization was observed for *C. perfringens* at 100% PSL of either ampicillin or penicillin.

For additional information on antimicrobial agents neutralized by BacT/ALERT FN Plus Culture Bottles, contact your local bioMérieux representative.

## PERFORMANCE CHARACTERISTICS OF THE TEST

### Potentially Interfering Substances

In-house seeded studies were conducted with cerebrospinal fluid, pleural fluid, synovial fluid, plasma, blood, and blood clots. Aliquots of each of these fluids also received white blood cells at concentrations relevant to bacteremia in each given body fluid. Testing was conducted with and without microorganisms. These substances neither interfered with recovery and detection of organisms, nor did they generate false positive results in the absence of organisms.

### Analytical Sensitivity: Limit of Detection (LoD)

Data in Table 1 represent results from in-house seeded studies. A minimum of 30 replicates were tested per species. Data in Table 1 were generated using bottles at the end of shelf life. Bottles inoculated with *B. fragilis* and *S. pneumoniae* received 1 ml pooled human blood supplementation. At least 95% detection was achieved at LoD.

**Table 1. Analytical Sensitivity: Limit of Detection (LoD)**

Microorganism	Strain ID	LoD (CFU/bottle)
<i>Bacteroides fragilis</i>	ATCC® 25285™	5
<i>Clostridium perfringens</i>	NCTC 8798	4
<i>Enterobacter aerogenes</i>	ATCC® 13048™	8
<i>Enterococcus faecalis</i>	NCTC 12697	4
<i>Escherichia coli</i>	NCTC 12923	4
<i>Listeria monocytogenes</i>	ATCC® 15313™	6
<i>Salmonella enterica</i>	ATCC® 14028™	5
<i>Staphylococcus aureus</i>	NCTC 10788	4
<i>Streptococcus pneumoniae</i>	ATCC® 6305™	6

**NOTE:** 86.6% of the bottles were subcultured within 30 minutes of being declared positive.

### Analytical Sensitivity: Growth Performance

Data in Table 2 represent results from in-house seeded studies with and without blood obtained from healthy human volunteers. Multiple strains were tested for each species at target inoculum levels of 125 CFU/bottle. The actual inoculum levels ranged from 3-299 CFU/bottle. In this seeded study, the BacT/ALERT FN Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. The species listed are representatives of clinically prevalent organisms in blood cultures and sterile body fluids.

**Table 2. Analytical Sensitivity: Growth Performance**

Microorganism	Blood				No Blood			
	% Recovery (n)	Range (CFU/bottle)	Time to Detection (hours)		% Recovery* (n=3)	Range (CFU/bottle)	Time to Detection (hours)	
			Mean	Range			Mean	Range
<i>Staphylococcus aureus</i>	100.0 (15/15)	54-150	14.6	12.9-16.7	100.0	116-150	21.8	21.3-22.0
<i>Escherichia coli</i>	100.0 (15/15)	73-254	10.9	10.4-12.4	100.0	73-176	11.6	10.4-12.9
<i>Bacteroides fragilis</i>	100.0 (18/18)	9-154	29.7	24.3-43.6	66.7 (2/3)	19-154	97.2	79.2-115.2
<i>Streptococcus pneumoniae</i>	100.0 (15/15)	4-260	16.5	11.5-43.8	100.0	4-25	17.5	16.0-19.3
<i>Clostridium perfringens</i>	100.0 (18/18)	58-210	17.3	12.0-40.1	100.0 (8/8)	76-210	27.1	14.1-35.7
<i>Klebsiella pneumoniae</i>	100.0 (15/15)	89-123	11.2	10.4-13.1	100.0	95-123	12.6	12.1-13.4
<i>Fusobacterium nucleatum</i>	70.6 (12/17)	19-204	74.4	36.0-108.0	75.0 (3/4)	19-116	56.1	41.2-64.8
<i>Streptococcus agalactiae</i>	100.0 (15/15)	14-194	16.9	12.7-28.9	100.0	21-34	25.9	20.5-33.4
<i>Enterococcus faecalis</i>	100.0 (15/15)	63-259	13.7	11.9-19.4	100.0	71-169	22.4	17.8-24.8
<i>Parvimonas micra</i>	80.0 (16/20)	46-154	51.4	37.3-69.6	0.0 (0/4)	46-154	-	-
<i>Enterobacter cloacae</i>	100.0 (15/15)	111-200	11.9	11.1-12.5	100.0	111-185	13.0	11.9-14.7
<i>Proteus mirabilis</i>	100.0 (15/15)	36-213	11.4	10.9-12.5	100.0	36-213	11.9	11.5-12.7
<i>Eggerthella lenta</i>	86.7 (13/15)	83-175	41.0	34.8-60.0	66.7 (2/3)	83-151	46.0	44.0-48.0
<i>Staphylococcus epidermidis</i>	100.0 (15/15)	44-135	21.0	17.4-25.3	100.0	44-105	29.3	24.5-36.8
<i>Listeria monocytogenes</i>	100.0 (15/15)	121-251	17.1	15.5-19.3	100.0	121-251	19.2	17.7-20.3
<i>Clostridium tertium</i>	100.0 (15/15)	24	12.5	11.4-13.5	100.0	24	14.8	14.1-16.0
<i>Clostridium septicum</i>	50.0 (20/40)	25-146	31.7	13.9-62.4	40.0 (2/5)	90-146	43.4	17.1-69.6
<i>Peptoniphilus asaccharolyticus</i>	60.0 (12/20)	49-296	50.9	34.4-79.2	25.0 (1/4)	81-296	44.7	-

\*In case of less than 100.0% recovery, it is recommended to add blood such as sterile defibrinated horse blood (10.0% v/v).<sup>3</sup>

Less than 100% detection was observed for some species, including *Capnocytophaga ochracea*, *Cardiobacterium hominis*, *Haemophilus parainfluenzae*, and *Granulicatella adiacens*.

### Clinical Study Results (Blood Cultures)

Results compare BacT/ALERT FN Plus to BacT/ALERT FN blood cultures (for all compliant pairs).

A multi-center clinical study was conducted at three different geographic sites in the U.S. comparing the performance of the BacT/ALERT FN Plus and BacT/ALERT FN blood culture bottles for anaerobic culture pairs that received blood volumes between 6 ml and 10 ml (compliant pairs). A total of 2514 anaerobic bottle pairs were obtained from 1080 adult patients suspected of blood stream bacterial/yeast infections. Subcultures of both bottles were performed when either bottle in the set was determined to be positive by the BacT/ALERT System. A pair of bottles was determined to have a positive status if subculture of either the BacT/ALERT FN Plus or BacT/ALERT FN Culture Bottle was positive. A culture bottle was determined to be a "True Positive" if the culture was flagged positive by the BacT/ALERT System and resulted in growth of the isolate upon subculture of this bottle. True positive rates were calculated for the BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles, and the ratio of BacT/ALERT FN Plus true positives to BacT/ALERT FN true positives was calculated to compare performance. Clinical isolates recovered were classified as significant, contaminant, or unknown based on determination by the clinical trial sites.

A total of 312 isolates were recovered from all compliant anaerobic blood culture pairs with a positive status. There were a total of 289 bottle pairs that recovered at least one isolate by subculture of BacT/ALERT FN Plus or BacT/ALERT FN Culture Bottles. A total of 266 bottle pairs recovered a single isolate and 23 bottle pairs recovered two isolates. The total population reported in Table 3 comprises the 312 isolates recovered from positive bottle pairs and 2225 negative bottle pairs for a total of 2537 results. The BacT/ALERT FN Plus Culture Bottle detected a total of 282 isolates compared to the BacT/ALERT FN Culture Bottle that detected 192 isolates. Of the significant isolates, the BacT/ALERT FN Plus Culture Bottle detected a total of 202 isolates compared to the BacT/ALERT FN Culture Bottle that detected 150 isolates. Three false positives were identified by subculture of positive BacT/ALERT FN Plus Culture Bottles and comprised 0.12% (3/2537) of the study population.

The following tables compare results of the BacT/ALERT FN Plus to BacT/ALERT FN blood cultures for all compliant blood culture bottles that yielded any number of isolates on subculture (Table 3), a single isolate alone on subculture (Table 4), and multiple isolates on subculture (Table 5).

**Table 3. All Compliant Pairs With Single And Multiple Isolates Combined (Blood Cultures)**

Clinical Isolate Determination	BacT/ALERT FN Plus True Positives	% of BacT/ALERT FN Plus True Positives in Population	BacT/ALERT FN True Positives	% of BacT/ALERT FN True Positives in Population	Ratio of True Positives*
Significant	202	8.0 (202/2537)	150	5.9 (150/2537)	1.347
Contaminant	58	2.3 (58/2537)	30	1.2 (30/2537)	1.933
Unknown	22	0.9 (22/2537)	12	0.5 (12/2537)	1.833
Total	282	11.1 (282/2537)	192	7.6 (192/2537)	1.469

\*One hundred sixty two (162) isolates were detected by both BacT/ALERT FN Plus and BacT/ALERT FN, 120 isolates were detected only by BacT/ALERT FN Plus, and 30 isolates were detected only by BacT/ALERT FN. The ratio of true positive rates for overall isolates was 1.469 (282/192) with a 95% CI (1.317, 1.621).<sup>7</sup>

**Table 4. All Compliant Pairs With Single Isolates (Blood Cultures)**

Clinical Determination	BacT/ALERT FN Plus True Positives	BacT/ALERT FN True Positives	Ratio of True Positives*
Significant	177	133	1.331
Contaminant	51	27	1.889
Unknown	15	9	1.667
Total	243	169	1.438

\*One hundred forty six (146) isolates were detected by both BacT/ALERT FN Plus and BacT/ALERT FN, 97 isolates were detected only by BacT/ALERT FN Plus and 23 isolates were detected only by BacT/ALERT FN. The ratio of true positive rates for overall single isolates was 1.438 (243/169) with a 95% CI (1.286, 1.590).<sup>7</sup>

**Table 5. All Compliant Pairs With Multiple Isolates (Blood Cultures)**

Clinical Determination	BacT/ALERT FN Plus True Positives	BacT/ALERT FN True Positives	Ratio of True Positives*
Significant	25	17	1.471
Contaminant	7	3	2.333
Unknown	7	3	2.333
Total	39	23	1.696

\*Sixteen (16) isolates were detected by both BacT/ALERT FN Plus and BacT/ALERT FN, 23 isolates were detected only by BacT/ALERT FN Plus, and seven isolates were detected only by BacT/ALERT FN. The ratio of true positive rates for overall multiple isolates was 1.696 (39/23) with a 95% CI (1.088, 2.304).<sup>7</sup>

In this clinical study, there were a total of 2218 compliant anaerobic blood culture pairs of BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles with negative instrument results for both bottles after 5 days of incubation. Among these pairs, terminal subcultures on both bottles were performed for 63 pairs, and seven false negative results by both BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles were observed; subculture on BacT/ALERT FN Plus Culture Bottles alone was performed for 844 pairs, and four false negative results were observed; both subcultures were not performed for 1307 pairs of bottles. Results are summarized in Table 6.

**Table 6. Summary Of Percent False Negatives From Compliant Anaerobic Blood Culture Pairs That Were Flagged Negative By The Instrument For Both Bottles**

Subculture Performed BacT/ALERT FN Plus	Subculture Performed BacT/ALERT FN	% False Negative BacT/ALERT FN Plus	% False Negative BacT/ALERT FN
Yes	Yes	11.1 (7/63)	11.1 (7/63)
Yes	No	0.5 (4/844)	-

Of these seven positive subcultures, two subcultures yielded isolates that are strict aerobes (*Acinetobacter baumannii* and *Pseudomonas aeruginosa*). The BacT/ALERT FN Plus Culture Bottle is not intended to detect strict aerobes from blood or other normally sterile body fluids. Overall false negative rate for BacT/ALERT FN Plus based on a subset of terminal subcultures was 1.2% (11/907) and excluding strict aerobes was 1.0% (9/907).

A comparative yield of microorganisms (number of isolates) recovered on subculture of BacT/ALERT FN Plus and BacT/ALERT FN cultures is presented in Table 7.

**Table 7. Comparative Yield Of Microorganisms From BacT/ALERT FN Plus And BacT/ALERT FN (Number Of Isolates) – Blood Cultures**

Group	BacT/ALERT FN Plus	BacT/ALERT FN
Anaerobes*	10	4
Enterobacteriaceae	50	43
<i>Enterococcus</i> spp.	28	19
Yeasts	0	0
Other Gram-Negative	6	3
Other Gram-Positive	9	6
Coagulase-Negative <i>Staphylococcus</i>	93	53
<i>Staphylococcus aureus</i>	65	42
<i>Streptococcus</i> spp.	21	22

**NOTE:** Isolate table includes polymicrobial cultures.

\*Anaerobic isolates recovered in clinical trial: *Bacteroides caccae* (4), *Bacteroides fragilis* (2), *Lactobacillus* spp. (1), *Peptostreptococcus micros* (1), *Peptostreptococcus* spp. (1), *Prevotella denticola* (1), *Propionibacterium acnes* (1), *Propionibacterium* spp. (3)

### Clinical Study Results (Sterile Body Fluid Cultures)

A multi-center clinical study was conducted at four different geographic sites in the U.S. and Canada comparing the performance of the BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles with sterile body fluid specimens. A total of 339 bottle pairs were obtained from 310 adult patients suspected of sterile body fluid bacterial/yeast infections. Sterile body fluid types evaluated were continuous ambulatory peritoneal dialysis (CAPD) fluid, cerebrospinal fluid (CSF), peritoneal fluid, pleural fluid, and synovial fluid. Clinical isolates recovered were classified as significant, contaminant, or unknown based on determination by the clinical trial sites.

A total of 77 isolates were recovered from all anaerobic sterile body fluid culture pairs with a positive status. There were a total of 61 bottle pairs that recovered at least one isolate by subculture of BacT/ALERT FN Plus or BacT/ALERT FN Culture Bottles. A total of 50 bottle pairs recovered a single isolate, seven bottle pairs recovered two isolates, three bottle pairs recovered three isolates, and one bottle pair recovered four isolates. The total population reported in Table 8 comprises the 77 isolates recovered from positive bottle pairs and 278 negative bottle pairs for a total of 355 results. The BacT/ALERT FN Plus Culture Bottle detected a total of 72 isolates compared to the BacT/ALERT FN Culture Bottle that detected 59 isolates. Of the significant isolates, the BacT/ALERT FN Plus Culture Bottle detected a total of 52 isolates compared to the BacT/ALERT FN Culture Bottle that detected 50 isolates. No false positives were observed for the BacT/ALERT FN Plus Culture Bottle from the study population (0/355).

Table 8 compares results of the BacT/ALERT FN Plus to BacT/ALERT FN sterile body fluid cultures that yielded single or multiple isolates on subculture.

**Table 8. All Pairs With Single And Multiple Isolates Combined (Sterile Body Fluid Cultures)**

Clinical Determination	BacT/ALERT FN Plus True Positives	% of BacT/ALERT FN Plus True Positives in Population	BacT/ALERT FN True Positives	% of BacT/ALERT FN True Positives in Population	Ratio of True Positives*
Significant	52	14.6 (52/355)	50	14.1 (50/355)	1.040
Contaminant	12	3.4 (12/355)	2	0.6 (2/355)	6.000
Unknown	8	2.2 (8/355)	7	2.0 (7/355)	1.143
Total	72	20.3 (72/355)	59	16.6 (59/355)	1.220

\*Fifty four (54) isolates were detected by both BacT/ALERT FN Plus and BacT/ALERT FN, 18 isolates were detected only by BacT/ALERT FN Plus, and five isolates were detected only by BacT/ALERT FN. The ratio of true positive rates for overall isolates was 1.220 (72/59) with a 95% CI (1.044, 1.396).<sup>7</sup>

Table 9 summarizes the minimum specimen volume achieved in anaerobic sterile body fluid clinical trials.

**Table 9. BacT/ALERT FN Plus Sterile Body Fluids Fill Volume (ml) – Positive Status**

Specimen Type	Total # of Specimen	# of Positives	Minimum Specimen Volume (ml)
CAPD Fluid	75	22	1.0
CSF	24	2	0.1
Peritoneal Fluid	106	16	0.3
Pleural Fluid	93	22	0.2
Synovial Fluid	41	10	0.2
Total	339	72	-

In this clinical study, there were 278 pairs of BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles with negative instrument results for both bottles after 5 days of incubation. Among these pairs, terminal subcultures on both bottles were performed for 166 pairs and 11 false negative results by both BacT/ALERT FN Plus and BacT/ALERT FN Culture Bottles were observed. Subculture on BacT/ALERT FN Plus Culture Bottles alone was performed for 112 pairs and no false negative results were observed. Results are summarized in Table 10.

**Table 10. Summary Of Percent False Negatives From Anaerobic Sterile Body Fluid Culture Pairs That Were Flagged Negative By The Instrument For Both Bottles**

Subculture Performed BacT/ALERT FN Plus	Subculture Performed BacT/ALERT FN	% False Negative BacT/ALERT FN Plus	% False Negative BacT/ALERT FN
Yes	Yes	6.6 (11/166)*	6.6 (11/166)*
Yes	No	0.0 (0/112)	-

\*Of these 11 positive subcultures, 8 subcultures yielded isolates that were strict aerobes: *Acinetobacter baumannii* (2), *Candida albicans* (1), *Candida parapsilosis* (2), *Cryptococcus neoformans* (1), *Micrococcus* spp. (1), *Pseudomonas aeruginosa* (1).

**NOTE:** The BacT/ALERT FN Plus Culture Bottle is not intended to detect strict aerobes from blood or other normally sterile body fluids.

Overall false negative rate for BacT/ALERT FN Plus based on a subset of terminal subcultures is 4.0% (11/278) and excluding strict aerobes is 1.1% (3/278).

A comparative yield of microorganisms (number of isolates) recovered on subculture of BacT/ALERT FN Plus and BacT/ALERT FN cultures is presented in Table 11.

**Table 11. Comparative Yield Of Microorganisms From BacT/ALERT FN Plus And BacT/ALERT FN (Number Of Isolates) – Sterile Body Fluid Cultures**

Group	BacT/ALERT FN Plus	BacT/ALERT FN
Anaerobes*	6	6
Enterobacteriaceae	9	8
<i>Enterococcus</i> spp.	14	11
Yeasts	4	2
Other Gram-Negative	2	1
Other Gram-Positive	1	0
Coagulase-Negative <i>Staphylococcus</i>	19	16
<i>Staphylococcus aureus</i>	9	9
<i>Streptococcus</i> spp.	8	6

**NOTE:** Isolate table includes polymicrobial cultures.

\*Anaerobes recovered in clinical trial: *Bacteroides fragilis* (2: 1 BacT/ALERT FN only; 1 BacT/ALERT FN Plus and BacT/ALERT FN), *Clostridium ramosum* (1 BacT/ALERT FN Plus and BacT/ALERT FN), *Clostridium* spp. (1 BacT/ALERT FN Plus and BacT/ALERT FN), *Prevotella loeschii* (1 BacT/ALERT FN Plus and BacT/ALERT FN), *Prevotella melaninogenica* (1 BacT/ALERT FN Plus only), *Veillonella* spp. (1 BacT/ALERT FN Plus and BacT/ALERT FN)

Quality control was performed during the clinical study on each of three organisms (*Bacteroides fragilis*, *Bacteroides vulgatus*, and *Clostridium perfringens*), which were prepared using serial dilution and seeded into the BacT/ALERT FN Plus Culture Bottle at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle. Overall quality control results were found to be acceptable. Instances where unacceptable quality control results were observed were found to be due to technical errors (i.e., supplement not added, colony counts out of range, incorrect diluent, site failure to change bottle status after positive instrument signal and positive subculture, failure to use pre-reduced media/reagents). Repeat testing resulted in acceptable results.

#### Delayed Entry

Table 12 includes results from seeded studies using six species (*Staphylococcus aureus*, *Escherichia coli*, *Streptococcus pneumoniae*, *Enterococcus faecium*, *Bacteroides fragilis*, and *Clostridium perfringens*) at target concentrations of 100 CFU/bottle (acceptable range of 30-300 CFU/bottle) were generated at three sites. Actual inoculum levels ranged from 41-253 CFU/bottle. All bottles contained human blood from healthy volunteers and were held at specified temperatures and times prior to loading into the BacT/ALERT instrument. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 12. Delayed Entry**

Sample Input	Incubation Temperature (°C)	Hold Time (hours)	% Recovery	Time to Detection from Sample Inoculation (Hold Time + Instrument TTD in hours)	
				Mean	Range
Inoculated Test Bottles	Control	No delay	100.0 (89/89)	15.9	9.5-52.8
	2-8	48	100.0 (65/65)	63.3	50.1-90.4
	20-25	24	100.0 (62/62)	34.7	26.0-79.2
	20-25	36	100.0 (62/62)	43.6	38.0-78.6
	35-37	8	100.0 (72/72)	17.7	10.0-53.4
	35-37	24	80.0 (64/80)*	28.6	26.0-52.3
Negative Controls	All conditions		0.0 (0/51)	-	-

\*CAUTION: Culture bottles held at 35-37°C for 24 hours or longer before loading may not detect microorganisms and should be subcultured.

**Within-Laboratory Precision (Repeatability)**

Data in Table 13 represent results from in-house seeded studies conducted on 12 days on multiple instruments by multiple operators. Organisms were grown in the presence of clinically relevant concentrations of antimicrobials to which they are susceptible. In this seeded study BacT/ALERT FN Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. A minimum of 108 replicates were tested for each organism/antimicrobial combination.

**Table 13. Within-Laboratory Precision (Repeatability)**

Sample Input		Range (CFU/bottle)	% Recovery				Time to Detection (hours)	
Organism	Antimicrobial		Lot 1	Lot 2	Lot 3	Overall	Mean	Range
<i>B. fragilis</i>	Imipenem	136-406	100.0	100.0	100.0	100.0	36.9	30.2-55.2
<i>C. perfringens</i>	Vancomycin	75-204	100.0	94.4	100.0	98.2	14.5	11.1-22.0
<i>S. aureus</i>	Oxacillin	94-158	100.0	100.0	100.0	100.0	17.7	15.1-24.3

**Reproducibility**

Data in Table 14 represent results from seeded studies conducted at three sites using a target of 144 replicates per site on 3 days with a minimum of two operators per site. Reproducibility was evaluated on each of eight organisms. One organism (*S. pneumoniae*) was prepared using serial dilution and the other seven organisms were prepared using BioBall® products. *S. pneumoniae* was seeded into the BacT/ALERT FN Plus Culture Bottle at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle and the other seven organisms at a target range of 1-17 CFU/bottle. The actual inoculum ranged from 5-500 CFU/bottle for the 30-300 CFU/bottle range, and from 1-270 CFU/bottle for the 1-17 CFU/bottle range. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 14. Reproducibility**

Sample Input	% Recovery				Time to Detection (hours)		Inoculum Ranges (CFU/Bottle)
	Site 1	Site 2	Site 3	Overall	Mean	Range	
<i>S. aureus</i>	96.3% (26/27)	79.2% (19/24)	100.0% (33/33)	92.9% (78/84)	20.2	18.5-35.7	2-12
<i>E. coli</i>	100.0% (18/18)	79.2% (19/24)	100.0% (33/33)	93.3% (70/75)	12.8	11.4-20.8	2-11
<i>E. faecalis</i>	100.0% (30/30)	83.3% (20/24)	97.0% (32/33)	94.3% (82/87)	24.6	17.9-30.4	2-15
<i>C. perfringens</i>	100.0% (18/18)	96.8% (61/63)	100.0% (30/30)	98.2% (109/111)	12.2	10.2-17.3	3-122
<i>E. aerogenes</i>	90.0% (27/30)	75.0% (18/24)	90.5% (38/42)	86.5% (83/96)	14.6	11.9-16.7	1-270*
<i>L. monocytogenes</i>	100.0% (21/21)	100.0% (24/24)	100.0% (33/33)	100.0% (78/78)	22.8	20.6-37.0	1-13
<i>S. enterica</i>	100.0% (24/24)	79.2% (19/24)	100.0% (30/30)	93.6% (73/78)	13.3	12.4-14.4	1-16
<i>S. pneumoniae</i>	100.0% (30/30)	100.0% (36/36)	100.0% (18/18)	100.0% (84/84)	17.5	13.3-23.1	5-500
Overall	98.0% (194/198) 95% CI: 94.9%, 99.5%	88.9% (216/243) 95% CI: 84.3%, 92.6%	98.0% (247/252) 95% CI: 95.4%, 99.4%	94.8% (657/693) 95% CI: 92.9%, 96.3%	-		

\*Plate count of 270 CFU/bottle was arrived at by serial dilution.

These data include repeat testing performed as a result of laboratory errors at a single site (i.e., contaminated bottles/reagents, colony counts out of range and site failure to change bottle status after positive instrument signal and positive subculture). Data excluding the laboratory errors demonstrated 100% recovery with the exception of *E. aerogenes*, which exhibited 96.3% recovery for all sites combined.

## REFERENCES

1. Thorpe TC, Wilson ML, Turner JE, DiGiuseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. Principles and Procedures for Blood Cultures; Approved Guideline. CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI. Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.
7. Kondratovich M. Comparing Two Medical Tests When Results of Reference Standard Are Unavailable for Those Negative via Both Tests. *J. Biopharm Stat* 2008;18:1,145-166.

## AVAILABILITY

**bioMérieux**  
**BacT/ALERT® FN Plus**

100/case

**REF** 410852

For technical assistance in the USA, contact bioMérieux Customer Service at 1-800-682-2666. Outside the USA, contact your local bioMérieux Representative.

## INDEX OF SYMBOLS

Symbol	Meaning
<b>REF</b>	Catalogue number
	Manufacturer
	Temperature limitation
	Use by
<b>LOT</b>	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
<b>ECREP</b>	Authorized representative in the European Community
	This way up
<b>IVD</b>	In Vitro Diagnostic Medical Device
	Do not reuse
	Latex-free

Instructions for use provided in the kit or downloadable from [www.biomerieux.com/techlib](http://www.biomerieux.com/techlib)

## REVISION TABLE

This section contains a summary of changes made to the BacT/ALERT FN Plus Instructions for Use (9305049 C).

Revision Date	Revision Number	Change Type	Change Summary
2013-04	9305049 C	Requirement:	<b>Specimen Collection and Preparation, Expected Values, Neutralization of Antimicrobials, Performance Characteristics of the Test</b> – Revised text to include additional information on product performance based on clinical studies
			<b>Limitations of the Test</b> - Added Limitations 2, 8, and 9 following FDA review
		Administrative:	<b>Reagents</b> - Moved Cautions to Limitations of the Test section

**NOTE: Minor typographical, grammar, and formatting changes are not included in the revision history.**

### Change Type categories:

- **Correction** = Correction of documentation anomalies.
- **Requirement** = Implementation of new and modified (updated) intended use and performance characteristics.
- **Administrative** = Implementation of non-technical changes noticeable to the user.

bioMérieux, the blue logo, BacT/ALERT, BioBall, and FAN are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a registered trademark belonging to Clinical and Laboratory Standards Institutes, Inc.

Any other name or trademark is the property of its respective owner.

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## DENOTES REVISED SECTION

### INTENDED USE

BacT/ALERT® PF Plus Culture Bottles are used with the BacT/ALERT® Microbial Detection System in qualitative procedures for recovery and detection of aerobic and facultative anaerobic microorganisms (bacteria and yeast) from blood.

### SUMMARY AND EXPLANATION

The BacT/ALERT Microbial Detection System is used to determine if microorganisms are present in blood taken from a patient suspected of having bacteremia/fungemia. The BacT/ALERT System and culture bottles provide both a microbial detection system and a culture medium with suitable nutritional and environmental conditions for organisms commonly encountered in blood infections. The BacT/ALERT PF Plus Culture Bottle provides for detection of microorganisms when a small volume of blood is available. An inoculated bottle is placed into the instrument where it is incubated and continuously monitored for the presence of microorganisms that will grow in the BacT/ALERT PF Plus Culture Bottle.

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**BacT/ALERT® PF Plus** (color-coded yellow) – BacT/ALERT PF Plus disposable culture bottles contain 30 ml of complex medium and 1.6 g adsorbent polymeric beads. The medium consists of casein peptone (1.0% w/v), yeast extract (0.45% w/v), soybean peptone (0.3% w/v), meat peptone (0.1% w/v), sodium polyanethol sulfonate (SPS) (0.083% w/v), menadione (0.00005% w/v), hemin (0.0005% w/v), L-cysteine (0.03% w/v), pyruvic acid (0.1% w/v), pyridoxine HCl (0.001% w/v), nicotinic acid (0.0002% w/v), pantothenic acid (0.0002% w/v), thiamine HCl (0.0001% w/v), and other complex amino acid and carbohydrate substrates in purified water. Bottles contain an atmosphere of N<sub>2</sub>, O<sub>2</sub>, and CO<sub>2</sub>, under vacuum. The composition of the medium may be adjusted to meet specific performance requirements.

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**CAUTION: BacT/ALERT PF Plus Culture Bottles used to culture very small blood specimen volumes (0.5 ml or less) will require added blood such as sterile defibrinated horse blood (10.0% v/v) to support growth, particularly for the recovery of fastidious organisms such as *Haemophilus influenzae*, *Streptococcus pneumoniae*, and *Neisseria gonorrhoeae*.<sup>3</sup>**

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### Storage instructions

BacT/ALERT PF Plus Culture Bottles are ready for use. Store in an upright position protected from direct light at room temperature (15-30°C). An expiration date is printed on each bottle label. Do not use the culture bottles beyond the expiration date indicated. If the bottles are exposed to temperatures less than 15°C, precipitates may form that will disappear when the bottles are warmed to room temperature. Bottles must be at room temperature before use.

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Prior to use, the BacT/ALERT PF Plus Culture Bottles should be examined for evidence of damage or deterioration (discoloration). Bottles exhibiting evidence of damage, leakage, or deterioration should be discarded. The medium in undisturbed bottles should be clear, but there may be a slight opalescence or a trace of precipitate due to the anticoagulant SPS or the presence of adsorbent polymeric beads; do not confuse this with turbidity indicative of microbial growth. Do not use a bottle which contains medium exhibiting turbidity, a yellow sensor, or excess gas pressure; these are signs of possible contamination.

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### General Considerations

1. BacT/ALERT PF Plus Culture Bottles should be utilized by trained healthcare personnel. Correct specimen collection is extremely important when obtaining blood culture specimens. Refer to Cumitech 1C<sup>4</sup> for the proper specimen collection procedure.
2. Take care to prevent contamination during both bottle preparation and inoculation of the patient sample. Proper skin disinfection is an essential requirement to reduce the incidence of contamination.
3. Although not recommended by bioMérieux, blood may be drawn directly into collection tubes containing SPS. Tubes containing other anticoagulants should never be used for blood culture.<sup>5</sup>

4. bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. If there is an unavoidable delay, inoculated bottles may be maintained at room temperature up to 24 hours before loading into the instrument.
5. Optimal recovery of isolates will be achieved by adding maximum amounts of specimen. Use of lower volumes may adversely affect recovery and/or detection times of some organisms. The bottle's recommended specimen volume is up to 4 ml and the volume collected should be monitored by means of the 4 ml incremental markings on the bottle label.
6. In general, culturing higher volumes of blood containing small numbers of bacteria will improve recovery of bacteria by culture.<sup>5</sup>
7. The clinical studies were conducted with blood volumes as low as 0.1 ml. It is recommended, however, that up to 4 ml of blood be inoculated in the BacT/ALERT PF Plus Culture Bottle. To prevent over inoculation, monitor the blood volume intake into the culture bottle, using the 4 ml incremental markings on the bottle label.

**CAUTION: Direct vacuum draw using the BacT/ALERT PF Plus Culture Bottle can result in the collection of greater than 4 ml blood sample volume. As a result, caution is advised when performing direct draw collection of blood from young children, infants, and neonates where total blood volume is a concern.**

#### Bottle preparation

1. Label the culture bottle with patient information. The icons on the bottle label (☺, #, ☹) can be defined by the user.
2. Remove plastic flip-top from the culture bottle. Prior to inoculation, disinfect the culture bottle top with an alcohol swab or equivalent. Allow to air dry.
3. Clean the selected venipuncture site as recommended by your institution's approved procedure.

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**NOTE:** If inoculating more than one type of BacT/ALERT blood culture bottle using a butterfly blood collection set and direct draw adapter cap, inoculate first the aerobic culture bottle and then the adult anaerobic culture bottle so that any oxygen trapped in the tubing will not be transferred to the anaerobic bottle.

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- a. Hold the culture bottle at a position below the patient's arm with the bottle in an upright position (stopper uppermost).
- b. Collect the blood using a butterfly blood collection set and a blood collection adapter cap as recommended by your institution's approved procedure and inoculate directly into the culture bottle at the patient's bedside. Although lower sample volumes can be used, recovery may be improved using a sample volume closer to the recommended 4 ml.
- c. Release the tourniquet as soon as the blood starts to flow into the culture bottle, or within 2 minutes of application.
- d. Do not allow the culture bottle contents to touch the stopper or the end of the needle during the collection procedure.

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**NOTE:** If inoculating more than one type of BacT/ALERT blood culture bottle using syringe draw, inoculate first the adult anaerobic culture bottle and then the aerobic culture bottle so that any oxygen trapped in the syringe will not be transferred to the adult anaerobic bottle. Line demarcations on the bottle label should be used to assist in estimating the sample volume.

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### BacT/ALERT PF PLUS CULTURE BOTTLE TEST PROCEDURE

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1. Visually inspect bottles before testing. Do not use bottles with evidence of damage, leakage, or deterioration. Consider bottles with hemolysis, turbidity, excess gas pressure, yellow sensors, and/or evidence of growth as positive. Smear and subculture. Do not incubate unless smear is negative.

- After culture bottles have been loaded into the instrument, incubate 5 days or until designated positive.
- Smear and subculture all positive bottles. If the smear is negative, indicating a possible false positive, the bottle should be reloaded into the instrument until growth of the subculture or redesignation as positive. Bottles that were initially determined false positive and were redesignated positive should be smeared and subcultured.
- Negative bottles may be checked by smear and/or subculture at some point prior to discarding as negative.
- Procedures for loading and unloading culture bottles into the BacT/ALERT instrument are given in the User Manual.
- Do not reuse BacT/ALERT culture bottles.** Dispose of inoculated BacT/ALERT culture bottles according to your laboratory protocol. Autoclaving and/or incinerating inoculated BacT/ALERT culture bottles is appropriate.<sup>2</sup>
- Utilization of coring devices (i.e., blunt needle) to puncture the septum may result in bottle leakage.

### QUALITY CONTROL

A Certificate of Conformance is available for each lot of culture bottles. If desired, individual laboratories can perform quality control testing of BacT/ALERT PF Plus Culture Bottles. Refer to the BacT/ALERT User Manual and CLSI® document M22-A3.<sup>6</sup>

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- Patient specimens determined positive by BacT/ALERT may contain organisms that are positive by smear that will not grow on routine subculturing media. When this is suspected, specimens should be subcultured on special media. Also, BacT/ALERT positive specimens may contain organisms that are not seen with routine smear staining methods and may require both specialized staining and subculturing media for detection and recovery.
- It is possible that certain rare, fastidious microorganisms will not grow or may grow slowly in the BacT/ALERT PF Plus Culture Bottle growth medium. In addition, on rare occasions, organisms may be encountered that grow in the BacT/ALERT PF Plus Culture Bottle growth medium, but do not produce sufficient carbon dioxide to be determined positive. If rare, fastidious organisms requiring specialized media and culture conditions are suspected, alternative methods or extended incubation time should be considered for recovery.
- Certain strains of *Haemophilus influenzae*, *Neisseria meningitidis*, and *Neisseria gonorrhoeae* may be sensitive to the anticoagulant SPS, which may result in a lack of growth or low production of CO<sub>2</sub> by these strains if an insufficient amount of sample is inoculated into the culture bottles.
- Infrequently, if there is a very high number of white blood cells present in the sample, the BacT/ALERT may indicate a culture bottle positive. In this case, the smear and subculture results may be negative.

- Organisms are often few in numbers and may appear intermittently in the blood stream; therefore, consecutive blood samples should be collected from each patient.
- Promptly remove positive culture bottles when they are signaled by BacT/ALERT to avoid possible non-viable cultures due to autolysis or other reasons. Certain strains of *Streptococcus pneumoniae* may be particularly prone to autolysis if they are not removed promptly after being signaled positive.
- A Gram-stained smear from a negative bottle may sometimes contain a small number of non-viable organisms that were derived from culture medium components, staining reagents, immersion oil, or glass slides, resulting in a false positive smear.
- bioMérieux recommends that inoculated culture bottles be placed into the BacT/ALERT Microbial Detection System as soon as possible after collection. But, in the unavoidable cases when there is a delay in bottle receipt by the laboratory, delayed entry information is provided from seeded studies in the "Performance Characteristics of the Test" section.
- Antimicrobial neutralization was not achieved for ceftazidime or cefepime.

### EXPECTED VALUES

- Percent positive cultures were observed to be 6.3% (range: 4.9%-8.1%) overall and 4.1% (range: 2.5%-6.4%) for significant isolates from three clinical trial sites in BacT/ALERT PF Plus Culture Bottles that received 0.1-4 ml of blood.
- Expected percent positives will vary based on factors such as patient population, prevalence of significant organisms, site location, and contamination rates. The expected values provided are based on clinical study data.

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Neutralization of antimicrobials by adsorbent polymeric beads varies depending upon dosage level and timing of specimen collection. Internal studies have demonstrated that antimicrobials are effectively neutralized by the BacT/ALERT PF Plus medium based on 100% recovery of the organisms tested. In these tests, antimicrobials were added in clinically relevant concentrations directly to culture bottles during inoculation with susceptible strains. The effectiveness of the antimicrobials was confirmed by parallel testing using a non-neutralizing medium as a control. Antimicrobials from the following categories were neutralized by the medium: penicillins, glycolcyclines, polyenes, macrolides, triazoles, echinocandins, cefazolin, cefoxitin, ceftaroline, aminoglycosides, fluoroquinolones, lincosamides, glycopeptides, and oxazolidinones.

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For additional information on antimicrobial agents neutralized by BacT/ALERT PF Plus Culture Bottles, contact your local bioMérieux representative.

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In-house seeded studies were conducted with plasma, blood, and blood clots. Aliquots of each of these fluids also received white blood cells at concentrations relevant to bacteremia in blood. Testing was conducted with and without microorganisms. These substances neither interfered with recovery and detection of organisms, nor did they generate false positive results in the absence of organisms.

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Data in Table 1 represent results from in-house seeded studies. A minimum of 30 replicates were tested per species. At least 95% detection was achieved at LoD. Data in Table 1 were generated using bottles at the end of shelf life. Bottles inoculated with *H. influenzae* received 4 ml pooled human blood supplementation.

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Data in Table 2 represent results from in-house seeded studies with blood obtained from healthy human volunteers. Multiple strains were tested for each species at target inoculum levels of 125 CFU/bottle. The actual inoculum levels ranged from 3-298 CFU/bottle. In this seeded study, the BacT/ALERT PF Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. The species listed are representatives of clinically prevalent organisms in blood cultures.

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Results compare BacT/ALERT PF Plus to BacT/ALERT PF blood cultures (for all compliant pairs).

A multi-center clinical study was conducted at three different geographic sites in the U.S. comparing the performance of the BacT/ALERT PF Plus and BacT/ALERT PF blood culture bottles for pediatric culture pairs that received blood volumes between 0.1 ml and 4 ml (compliant pairs). A total of 2188 bottle pairs were obtained from 1086 pediatric patients suspected of blood stream bacterial/yeast infections. Subcultures of both bottles were performed when either bottle in the set was determined to be positive by the BacT/ALERT System. A pair of bottles was determined to have a positive status if the subculture of either the BacT/ALERT PF Plus or BacT/ALERT PF Culture Bottle was positive. A culture bottle was determined to be a "True Positive" if the culture was flagged positive by the BacT/ALERT System and resulted in growth of the isolate upon subculture of this bottle. True positive rates were calculated for the BacT/ALERT PF Plus and BacT/ALERT PF Culture Bottles, and the ratio of BacT/ALERT PF Plus true positives to BacT/ALERT PF true positives was calculated to compare performance. Clinical isolates recovered were classified as significant, contaminant, or unknown based on determination by the clinical trial sites.

A total of 172 isolates were recovered from all compliant pediatric blood culture pairs with a positive status. There were a total of 145 bottle pairs that recovered at least one isolate by subculture of BacT/ALERT PF Plus or BacT/ALERT PF Culture Bottles. A total of 126 bottle pairs recovered a single isolate, 12 bottle pairs recovered two isolates, six bottle pairs recovered three isolates, and one bottle pair recovered four isolates. The total population reported in Tables 3, 4, and 5 comprises the 172 isolates recovered from positive bottle pairs and 2043 negative bottle pairs for a total of 2215 results. The BacT/ALERT PF Plus Culture Bottle detected a total of 140 isolates compared to the BacT/ALERT PF Culture Bottle that detected 128 isolates. Of the significant isolates, the BacT/ALERT PF Plus Culture Bottle detected a total of 91 isolates compared to the BacT/ALERT PF Culture Bottle that detected 77 isolates. One false positive was identified by subculture of a positive BacT/ALERT PF Plus Culture Bottle and comprised 0.05% (1/2215) of the study population.

The following tables compare results of the BacT/ALERT PF Plus to BacT/ALERT PF blood cultures for all compliant blood culture bottles that yielded any number of isolates on subculture (Table 3), a single isolate alone on subculture (Table 4), and multiple isolates on subculture (Table 5).

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Clinical Isolate Determination	BacT/ALERT PF Plus True Positives	% of BacT/ALERT PF Plus True Positives in Population	BacT/ALERT PF True Positives	% of BacT/ALERT PF True Positives in Population	Ratio of True Positives*
Significant	91	4.1 (91/2215)	77	3.5 (77/2215)	1.182
Contaminant	24	1.1 (24/2215)	29	1.3 (29/2215)	0.828
Unknown	25	1.1 (25/2215)	22	1.0 (22/2215)	1.136
Total	140	6.3 (140/2215)	128	5.8 (128/2215)	1.094

\*Ninety six (96) isolates were detected by both BacT/ALERT PF Plus and BacT/ALERT PF, 44 isolates were detected only by BacT/ALERT PF Plus, and 32 isolates were detected only by BacT/ALERT PF. The ratio of true positive rates for overall isolates was 1.094 (140/128) with a 95% CI (0.954, 1.234).<sup>7</sup>

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\*Seventy three (73) isolates were detected by both BacT/ALERT PF Plus and BacT/ALERT PF, 32 isolates were detected only by BacT/ALERT PF Plus, and 21 isolates were detected only by BacT/ALERT PF. The ratio of true positive rates for overall single isolates was 1.117 (105/94) with a 95% CI (0.957, 1.277).<sup>7</sup>

**Table 5. All Compliant Pairs With Multiple Isolates (Blood Cultures)**

Clinical Determination	BacT/ALERT PF Plus True Positives	BacT/ALERT PF True Positives	Ratio of True Positives*
Significant	22	16	1.375
Contaminant	7	12	0.583
Unknown	6	6	1.000
Total	35	34	1.029

\*Twenty three (23) isolates were detected by both BacT/ALERT PF Plus and BacT/ALERT PF, 12 isolates were detected only by BacT/ALERT PF Plus, and 11 isolates were detected only by BacT/ALERT PF. The ratio of true positive rates for overall multiple isolates was 1.029 (35/34) with a 95% CI (0.748, 1.310).<sup>7</sup>

In this clinical study, there were 2041 pairs of BacT/ALERT PF Plus and BacT/ALERT PF Culture Bottles with negative instrument results for both bottles after 5 days of incubation. Among these pairs, terminal subcultures on both bottles were performed for three pairs, and no false negative results by both BacT/ALERT PF Plus and BacT/ALERT PF Culture Bottles were observed; subculture on BacT/ALERT PF Plus Culture Bottles alone was performed for 2034 pairs, and one false negative result was observed; both subcultures were not performed for four pairs of bottles. Results are summarized in Table 6.

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Subculture Performed BacT/ALERT PF Plus	Subculture Performed BacT/ALERT PF	% False Negative BacT/ALERT PF Plus	% False Negative BacT/ALERT PF
Yes	Yes	0.00 (0/3)	0.00 (0/3)
Yes	No	0.05 (1/2034)	-

Overall false negative rate for BacT/ALERT PF Plus based on a subset of terminal subcultures was 0.05% (1/2037).

A comparative yield of microorganisms (number of isolates) recovered on subculture of BacT/ALERT PF Plus and BacT/ALERT PF cultures is presented in Table 7.

**Table 7. Comparative Yield Of Microorganisms From BacT/ALERT PF Plus And BacT/ALERT PF**

Organism Group	Pediatric Subgroup	BacT/ALERT PF Plus	BacT/ALERT PF Plus Fill Range (ml)	BacT/ALERT PF	BacT/ALERT PF Fill Range (ml)
Enterobacteriaceae	Newborn (< 1 mo)	6	0.1-1.4	7	0.1-1.6
	Infant (> 1 mo–2 yrs)	19	0.1-3.7	13	0.1-3.2
	Child (> 2 yrs–12 yrs)	9	0.6-3.9	7	0.3-3.2
	Adolescent (>12–21 yrs)	0	1.4	1	1.6
Fastidious ( <i>Neisseria meningitidis</i> , <i>Neisseria sicca</i> )	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	1	0.2	1	0.5
	Child (> 2 yrs–12 yrs)	0	-	0	-
	Adolescent (>12–21 yrs)	1	1.1	1	0.6
Yeast ( <i>Candida albicans</i> , <i>C. guilliermondii</i> , <i>C. krusei</i> , <i>C. lusitanae</i> )	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	0	0.5	1	0.9
	Child (> 2 yrs–12 yrs)	5	0.9-3.7	6	1.0-3.4
	Adolescent (>12–21 yrs)	1	0.2-3.0	2	2.1-2.5
Non-fermentative Gram-Negative Bacilli	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	5	1.7-3.5	6	1.5-2.5
	Child (> 2 yrs–12 yrs)	2	0.9-2.2	3	1.0-2.8
	Adolescent (>12–21 yrs)	0	-	0	-
Coagulase-Negative <i>Staphylococcus</i>	Newborn (< 1 mo)	5	0.1-0.5	5	0.1-0.9
	Infant (> 1 mo–2 yrs)	12	0.1-3.0	10	0.1-3.4
	Child (> 2 yrs–12 yrs)	15	0.1-3.8	12	0.5-3.6
	Adolescent (>12–21 yrs)	6	0.5-3.5	7	0.5-3.2
<i>Staphylococcus aureus</i>	Newborn (< 1 mo)	0	0.3	1	0.1
	Infant (> 1 mo–2 yrs)	5	0.5-1.5	5	0.6-1.6
	Child (> 2 yrs–12 yrs)	7	0.8-4.0	3	0.1-3.6
	Adolescent (>12–21 yrs)	2	1.5-1.7	2	1.3-1.4
<i>Enterococcus</i> spp.	Newborn (< 1 mo)	1	0.1	1	0.1
	Infant (> 1 mo–2 yrs)	9	0.2-2.9	10	0.1-3.2
	Child (> 2 yrs–12 yrs)	2	0.2-1.0	1	0.8-1.8
	Adolescent (>12–21 yrs)	8	1.5-3.1	7	1.9-2.9
<i>Streptococcus pneumoniae</i>	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	2	1.0-2.7	1	1.1-1.6
	Child (> 2 yrs–12 yrs)	0	-	0	-
	Adolescent (>12–21 yrs)	0	-	0	-
<i>Streptococcus</i> spp., Group A, B	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	1	0.5	0	1.0
	Child (> 2 yrs–12 yrs)	0	-	0	-
	Adolescent (>12–21 yrs)	0	-	0	-
Other <i>Streptococcus</i> spp.	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	4	0.1-1.6	5	0.1-1.5
	Child (> 2 yrs–12 yrs)	3	0.7-1.1	2	0.6-2.6
	Adolescent (>12–21 yrs)	1	0.7-2.4	2	1.1-2.4
Other Gram-negative*	Newborn (< 1 mo)	0	-	0	-
	Infant (> 1 mo–2 yrs)	1	0.6	0	0.4
	Child (> 2 yrs–12 yrs)	0	-	0	-
	Adolescent (>12–21 yrs)	0	-	0	-
Other Gram-positive**	Newborn (< 1 mo)	1	0.1-0.7	2	0.1-0.4
	Infant (> 1 mo–2 yrs)	3	0.5-3.0	2	0.9-2.3
	Child (> 2 yrs–12 yrs)	3	0.1-3.4	1	0.1-3.2
	Adolescent (>12–21 yrs)	0	3.3	1	3.2

\*Other Gram-Negative Organisms recovered in clinical study: Unidentified Gram-Negative Rods (1)

\*\*Other Gram-Positive Organisms recovered in clinical study: *Bacillus* spp. (7), *Corynebacterium* spp. (2), Diphtheroids (1), *Micrococcus* spp. (2), *Stomatococcus* spp. (1).

Organisms were recovered in the clinical studies at blood volumes greater than 0.1 ml. They are:

- ≥1.0 ml for *Streptococcus pneumoniae*
- ≥0.9 ml for non-fermenting Gram-negative rods; ≥ 0.6 ml for unidentified Gram-negative rods
- ≥0.5 ml for Group A, B *Streptococcus* spp.
- ≥0.3 ml for *S. aureus*
- ≥0.2 ml for fastidious organism (*N. meningitidis* and *N. sicca*)
- ≥0.2 ml for Yeast (*Candida albicans*, *C. guilliermondii*, *C. krusei*, and *C. lusitanae*)

Quality control was performed during the clinical study on each of 13 organisms (*Candida albicans*, *Candida krusei*, *Enterococcus faecalis*, *Escherichia coli*, *Haemophilus influenzae*, *Neisseria meningitidis*, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus agalactiae*, *Streptococcus pneumoniae*, and *Streptococcus pyogenes*), which were prepared using serial dilution and seeded into the BacT/ALERT PF Plus Culture Bottle at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle. Overall quality control results were found to be acceptable. Instances where unacceptable quality control results were observed were found to be due to technical errors (i.e., colony counts out of range, site failure to change bottle status after positive instrument signal and positive subculture, and no supplement added). Repeat testing resulted in acceptable results.

#### Delayed Entry

Table 8 includes results from seeded studies using 11 species (*Staphylococcus aureus*, *Candida albicans*, *Candida krusei*, *Escherichia coli*, *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Pseudomonas aeruginosa*, *Streptococcus pneumoniae*, *Enterococcus faecium*, *Haemophilus influenzae*, and *Neisseria meningitidis*), at target concentrations of 100 CFU/bottle (acceptable range of 30-300 CFU/bottle) were generated at three sites. Actual inoculum levels ranged from 35-290 CFU/bottle, although there was an instance of a colony count that was >300 CFU/bottle which was not included in the final data analysis. All bottles contained human blood from healthy volunteers and were held at specified temperatures and times prior to loading into the BacT/ALERT instrument. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 8. Delayed Entry**

Sample Input	Incubation Temperature (°C)	Hold Time (hours)	% Recovery	Time to Detection from Sample Inoculation (Hold Time + Instrument TTD in hours)		Inoculum Range (CFU/bottle)
				Mean	Range	
Inoculated Test Bottles	Control	No delay	100.0 (459/459)	14.3	8.5-84.0	35-288
	2-8	48	98.6 (292/296)	63.7	57.5-103.2	48-288
	20-25	24	98.0 (291/297)	31.8	26.2-74.4	50-288
	20-25	36	91.9 (272/296)	41.8	38.0-70.5	50-290
	35-37	8	98.9 (454/459)	16.1	10.2-53.8	35-288
	35-37	24	56.6 (259/458)*	28.3	26.0-74.4	35-288
Negative Controls	All conditions		0.5 (1/221)**	-	-	-

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\*\*One false positive result observed during seeded study. Negative confirmation by Gram stain/subculture.

#### Within-Laboratory Precision (Repeatability)

Data in Table 9 represent results from in-house seeded studies conducted on 12 days on multiple instruments by multiple operators. Organisms were grown in the presence of clinically relevant concentrations of antimicrobials to which they are susceptible. In this seeded study BacT/ALERT PF Plus Culture Bottles were subcultured at least 24 hours after being flagged positive by the instrument. A minimum of 108 replicates were tested for each organism/antimicrobial combination.

**Table 9. Within-Laboratory Precision (Repeatability)**

Sample Input		Range (CFU/bottle)	% Recovery				Time to Detection (hours)	
Organism	Antimicrobial		Lot 1	Lot 2	Lot 3	Overall	Mean	Range
<i>C. albicans</i>	Fluconazole	140-364	100.0	100.0	100.0	100.0	26.0	22.8-31.3
<i>E. coli</i>	Amikacin	26-156	100.0	100.0	100.0	100.0	12.0	11.2-13.0
<i>K. pneumoniae</i>	Levofloxacin	108-170	100.0	100.0	100.0	100.0	13.4	11.7-15.2
<i>P. aeruginosa</i>	Piperacillin	80-148	100.0	97.2	100.0	99.1	19.2	17.4-24.1
<i>S. pneumoniae</i>	Penicillin G	9-505	100.0	100.0	100.0	100.0	13.2	11.6-15.5
<i>S. aureus</i>	Vancomycin	94-158	100.0	100.0	100.0	100.0	16.9	14.6-20.3

**Reproducibility**

Data in Table 10 represent results from seeded studies conducted at three sites using a target of 162 replicates per site on 3 days with a minimum of two operators per site. Reproducibility was evaluated on each of nine organisms. Two organisms (*C. albicans* and *S. pneumoniae*) were prepared using serial dilution and the other seven organisms were prepared using BioBall® products. *C. albicans* and *S. pneumoniae* were seeded into the BacT/ALERT PF Plus Culture Bottle at a target inoculum of 100 CFU/bottle, with an acceptable range of 30-300 CFU/bottle and the other seven organisms at a target range of 1-17 CFU/bottle. The actual inoculum ranged from 6-700 CFU/bottle for the 30-300 CFU/bottle range, and from 1-270 CFU/bottle for the 1-17 CFU/bottle range. Percent recovery reflects positive flag by the instrument and Gram-stain/subculture consistent with the seeded organism.

**Table 10. Reproducibility**

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These data include repeat testing performed as a result of laboratory errors at a single site (i.e., contaminated subculture). Data excluding the laboratory errors demonstrated 100% recovery with the exception of *E. aerogenes*, which exhibited 85% recovery for all sites combined.

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## REVISION TABLE

This section contains a summary of changes made to the BacT/ALERT PF Plus Instructions for Use (9305130 C).

Revision Date	Revision Number	Change Type	Change Summary
2013-04	9305130 C	Requirement:	<b>Intended Use, Specimen Collection and Preparation, Expected Values, Neutralization of Antimicrobials, Performance Characteristics of the Test</b> – Revised text to include additional information on product performance based on clinical studies
			<b>Limitations of the Test</b> - Added Limitations 2, 8, and 9 following FDA review
		Administrative:	<b>Reagents</b> - Moved Cautions to Limitations of the Test section

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CLSI yra registruotasis „Clinical and Laboratory Standards Institutes, Inc.“ prekės ženklas.