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▪ Introduction

▪ Main features

- Automated process from blood tube
- 10 different tests in parallel
- **STAT* testing for emergency**
- **Bidirectional LIS/HIS**

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▪ Advanced features

- Pressure sensor for the accurate sampling
- Air blower
- Result analysis with AI
- **Extra Tipstation**
- Online update
- 1:1 matching with cartridge and test parameter
- Auto door close
- Exact reading of barcode

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
▪ Unique selling proposition

* **STAT** laboratory tests and services are those that are needed immediately in order to manage medical emergencies.

Introduction

Applicable field of AFIAS system

COVID-19 testing

Neonatal Intensive Care Unit (NICU)

Emergency Room


Clinical laboratories

Introduction


Improvement for customer needs

User convenience

Accurate diagnosis

Automatic process

High throughput

Main features

- **Automated process from blood tube**
- **10 different tests in parallel**
- **STAT* testing for emergency**
- **Bidirectional LIS/HIS**


* **STAT** laboratory tests and services are those that are needed immediately in order to manage medical emergencies.

Automated process from blood tube

Operation procedure –blood tube to result


10 blood tubes in a tube rack

Insert a tube rack into AFIAS-10

Auto sampling process


Test result

Automated process from blood tube

Auto sampler


- **Various types of collection tube available**
 - Whole blood tube (7 cm & 10 cm, 1.5 mL micro blood tube)
 - Microcentrifuge tube (1.5 mL & 0.5 mL)
- **Dedicated rack provided**

10 different tests in parallel

- AFIAS-10 has 10 cartridge slots that enable 10 different parameters testing in parallel
- Each cartridge can be identified with its own color
- Multiple testing available from 1 sample


Automated process from blood tube

Fully random access with 10 slots

- Maximum usage of cartridge slots
 - AFIAS-6 : 2 bay random access
 - An empty slot in a working bay cannot be accessed
 - AFIAS-10 : 10 slots random access
 - You can insert an unplanned cartridge to be tested next



STAT testing for emergency

- **AFIAS-10** can implement **STAT testing** for medical emergency priority while regular routine processing
- **STAT** laboratory tests and services are those that are needed immediately in order to manage medical emergencies. STAT test requests are given the highest priority by the Clinical Laboratory for processing, analysis and reporting

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Emergency

Insert cartridge and select 'Emergency test'

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■ Bidirectional LIS function using the PID (barcode) of the sample tube

- Fully automated sample and test cartridge assignment
- Necessary to apply the dedicated protocol to the hospital network

■ Usage Scenarios

- After inserting the cartridge, insert the sample rack (after the sample rack)
- Insert the sample rack first, then insert the cartridge (the sample rack first)

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Loading sample rack first (1)


1. Insert the tube rack

2. The types of cartridge and the amount required

3. Insert cartridges


4. Identify the cartridges inserted

Bidirectional LIS/HIS

Loading sample rack first (2)

5. The cartridges and samples are automatically matched based on the test order through the bidirectional LIS/HIS

Bidirectional LIS/HIS

Loading cartridges first (1)

1. Insert and identify the cartridges

2. Insert the tube rack

Bidirectional LIS/HIS

Loading cartridges first (2)

3. The cartridges and samples are automatically matched


Advanced features

- Pressure sensor for the accurate sampling
- Air blower
- Result analysis with AI
- Extra tipstation
- Online update
- 1:1 matching with cartridge and test parameter
- Auto door close
- Exact reading of barcode

Pressure sensor for the accurate sampling


- Blood tube – capable of detecting the level/bubble of a whole blood sample with intrinsic foam
 - Differentiates water level/bubble by tracking the pressure change of the material sucked into the tip by a pressure sensor
- Detect bubble-caused errors to avoid misleading results

Pressure sensor for the accurate sampling


- Sensing and removal of bubble in the sample
 - Repeat 4 times (max.) by checking pressure between tip and sample surface
- When the foam is in the sample;
 - Marked as (B)
 - Wrong results or invalid error possible
 - Retest suggested

No	Patient ID	Item	Sample Type	Result	Date	Lot
1385	210728 008	Tn-I Plus BloodTube	S/P (B)	<0.01 ng/mL	Jul.28.2021 14:19:55	TNRCB04G
1384	210728 007	CK-MB BloodTube	S/P (B)	Invalid	Jul.28.2021 14:09:52	CKQGA27F

Air Blower

- Pressure Sensor, Blower


- When the tip is ejected from the adapter, the air blower enables to remove the residual solution
- Washing-free tip adapter


Result analysis with AI

- Wrong result error
 - Error patterns that RawCalculation logic (previous analysis) cannot detect

Result analysis with AI



- Deep learning process of accumulated error signal patterns

Extra TipStation

- AFIAS-1 / 6 : Putting the tip into the tip hole in the cartridge
- AFIAS-10 :
 - Putting the tip into the tip hole in the cartridge
 - Putting 10 tips in the Extra TipStation at a time.
 - Enhanced user convenience
- TipStation UI
 - Real-time display of the number of tips mounted on the instrument

Online update – S/W (TBD)

- Check the version of existed S/W
- Request for update online with the latest version

Online update – IDChip information (TBD)


- Check the version of IDChip for selected items
- Request for update online with the latest version

1:1 matching with cartridge and test parameter

- Automatic assign of cartridges requested
 - Applicable only for single test parameter
 - User convenience
 - Suitable for the mass testing with single parameter



Auto Door Close

- A function to automatically close the door by moving the bay into the device when the door is left open for a certain period
 - Avoiding the contamination from dust
 - Constant temperature

A countdown timer is displayed at the top right of the bay drive button


Emphasize to the user if the remaining time is less than 10 seconds

- Cancel the countdown when inserting a new cartridge
- When the countdown is completed, Bay Insertion & Door Close automatically
- Waiting time can be specified

Exact reading of barcode

- Using two or more sample barcodes in one tube
- Improved function so that auto sampler of AFIAS-10 can recognize even in multi-barcode situation (Korean hospital requested)

(X)

(O)

AFIAS-10

Specification

▪ Interface

- 10.1" TFT LCD touch screen
- LIS/HIS compatible (Bidirectional)
- Built-in thermal printer
- RS 232 / USB port / Ethernet / SD Card slot

▪ Memory

- 10,000 patient results / 5,000 QC results
- 100 ID chips /100 User IDs

▪ Throughput

- Up to 50 tests/hour

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• Analyzer

- 426 x 443 x 395 mm / 20 kg
- 100-240 V AC, 50/60 Hz, 2.6-1.0 A

AFIAS-10

*Available parameters**

Category	Item	Category	Item	Category	Item	Category	Item
Cardiac	Tn-I Plus	Diabetes	HbA1c	Infection	COVID-19 Ag	Autoimmune	Anti-CCP Plus
	CK-MB		Microalbumin		COVID-19 Ab		Total IgE
	D-Dimer	Hormone	TSH		COVID-19/Flu Ag Combo	Others	Ferritin
	NT-proBNP		TSH Plus		COVID-19 nAb		Vitamin D
	Myoglobin		T3		CRP		Infliximab
	Cardiac Triple		T4		PCT		Total Anti-Infliximab
	ST2		FSH		PCT Plus		
	BNP		Total β hCG		IL-6		
	Tn-T		β -HCG Plus		HBsAg		
Cancer	PSA		LH		Anti-HBs		
	PSA Plus		PRL		Anti-HCV		
	AFP		Progesterone		ROTA		
	AFP Plus		Testosterone		NORO		
	CEA		Cortisol		Rota/Adeno		
	CEA Plus		AMH				

*In-house evaluation has been completed and documentation has been processing

AFIAS-10

Target field

- **For medium-large hospitals & laboratories**
 - **More than 500** tests in a month
 - Fully automated diagnostic system
 - 10 different tests in parallel
 - Bidirectional LIS/HIS

- **For emergency testing**
 - **Fast result** in a single test
 - Cardiac / Inflammatory marker (PCT / IL-6 / CRP)
 - **STAT** testing
 - Quick test order with bidirectional LIS/HIS

AFIAS-10

Comparison with Competitors

	AFIAS-10	MINI VIDAS®	VIDAS® 3	AQT90 Flex
Manufacturer	Boditech Med Inc.	bioMerieux SA	bioMerieux SA	Radiometer Medical ApS
Dimensions (mm)	426 x 443 x 395	450 x 550 x 575	450 x 550 x 575	460 x 480 x 450
Weight (kg)	20	40	65	35
Principle	FIA Eu-TRF (Time-Resolved Fluorescence)	ELFA (Enzyme Linked Fluorescent Assay)	ELFA (Enzyme Linked Fluorescent Assay)	Eu-TRF (Time-Resolved Fluorescence)
Display	10.1" TFT LCD	7" LCD	10.1" LCD	8.4" color TFT-LCD
Input	Touch screen	Manual keypad	Touch screen	Touch screen
Cartridge slot	10 (5 x 2 bays)	12 (6 x 2 bays)	12 (3 x 4 bays)	2
Testing time (min)	12	20	20	18
Throughput (test/hr)	Up to 50	36	36	30
Shelf life (month)	20	12	12	12
Storage temperature	2-30°C	2-8°C	2-8°C	2-8°C
Blood tube	O	X	O	O
Whole blood test	O	X	X	O
Temperature control	O	X	O	O
LIS/HIS	bidirectional	Unidirectional	bidirectional	bidirectional
Calibration	Unnecessary	14 days	14 days	30 days

Advantages for users

- **Precise test result**
 - Accurate sampling
 - Air blower
 - Result analysis with AI
- **Enhanced user convenience**
 - Random Access
 - Auto sampler
 - Extra TipStation
 - STAT testing
 - Bidirectional LIS/HIS
 - Online update
- **Advanced functions that reflect user needs**
 - 1:1 matching with cartridge and test parameter
 - Auto door close
 - Exact reading of barcode

Unique selling proposition

Automated process from blood tube

- Providing convenience to users
- Higher performance

10 different tests in parallel

- High throughput result with various parameters
- Multi testing performance

STAT testing for emergency

- Fast result derived from medical decision
- Useful in ER/ICU

Bidirectional LIS/HIS

- Fast and quick order from clinicians
- Suitable for medical center information system

Expansion in market field!

Thank you for your attention!


Automatic Fluorescence Immuno Assay System

AFIAS-10

1 p.

- Automated process from blood collection tube
- 10 different tests in parallel
- STAT testing for emergency
- Bidirectional LIS/HIS

AFIAS-10 is a unique integrated Sample-to-Answer system from blood collection tube featuring 10 parallel tests of differing parameters

► **Fully automated diagnostic system**

- Sample loading to result
- PID barcode reading
- Various types of sample tube available
 - Blood collection tube (0.5 & 3 & 5 mL)
 - Microcentrifuge tube (0.5 & 1.5 mL)

► **Emergency test function**

- Emergency test first during regular test routine
- Processing the rest of the test routine after finishing emergency test

► **10 channels with 2 bays**


- Each bay can be operated separately
- 10 different tests in parallel
- Multi tests with one sample

► **Convenient data connectivity**

- Bidirectional LIS/HIS function
- Getting the medical order from clinicians directly

Simple procedure

1 Insert cartridge

2 Insert sample tube into tube rack

3 Insert tube rack into AFIAS-10

Automatic Fluorescence Immuno Assay System

AFIAS-10

All-in-One cartridge for AFIAS platform

- All the reagents for assay contained
- Over 50 test parameters
- Ready-to-use cartridge

AFIAS-10 can provide

- Convenient testing procedure to users
- High throughput results with various parameters
- Multi testing performance
- Prompt result derived from medical decision
- Quick transmission of medical information

Specifications

Model	AFIAS-10	Data Output	Internal Printer/RS232 port
Ordering Information	FPRR038	Operating temperature	15 – 35 °C
Dimensions	426 x 443 x 395 mm	Connectivity	• RS232C • 3 USB port • 1 mini USB port • LAN port • SD card slot
Weight	20 kg		
Analyzer input	100~240V AC, 50/60Hz, 2.6-1.0 A		
User Interface	10.1" TFT LCD touch screen		

AFIAS Parameters

Category	Item	Sample type	Working range	Test time (min)
Cardiac	Tn-I Plus	WB/S/P	0.01-15.00 ng/mL	12
	CK-MB	WB/S/P	3-100 ng/mL	12
	D-Dimer	WB/P	50-10,000 ng/mL	12
	NT-proBNP	WB/S/P	10-30,000 pg/mL	12
	Myoglobin	WB/S/P	2-500 ng/mL	12
	Cardiac Triple	WB/S/P	Tn-I: 0.01-15 ng/mL CK-MB: 3-100 ng/mL Myoglobin: 5-500 ng/mL	12
	ST2	WB/S/P	3.1-200 ng/mL	12
	BNP	WB/S/P	10-20,000 pg/mL	12
Cancer	Troponin T	WB/S/P	5-5,000 pg/mL	12
	PSA	WB/S/P	0.5-100.0 ng/mL	15
	PSA Plus	WB/S/P	0.07-50 ng/mL	12
	AFP	WB/S/P	5-350 ng/mL	15
	AFP Plus	WB/S/P	0.5-350 ng/mL	12
	CEA	S/P	1-500 ng/mL	12
Diabetes	CEA Plus	WB/S/P	0.5-200 ng/mL	12
	HbA1c	WB	(NGSP) 4-15%	10
	Microalbumin	Urine	2-300 mg/L	12
Hormone	TSH	WB/S/P	0.09-80.0 µIU/mL (WB: 0.4-80.0 µIU/mL)	15
	TSH Plus	WB/S/P	0.08-50 µIU/mL	12
	T3	S/P	0.77-7.70 nmol/L	10
	T4	S/P	10.23-300.00 nmol/L	10
	FSH	WB/S/P	1-100 mIU/mL	12
	Total βhCG	WB/S/P	5-50,000 mIU/mL	15
	β-HCG Plus	WB/S/P	2-1,500 mIU/mL	12
	LH	WB/S/P	1-100 mIU/mL	12
	Progesterone	S/P	4.45-127.2 nmol/L (1.4-40 ng/mL)	15
	PRL	WB/S/P	1-100 ng/mL	10
	Testosterone	S/P	1-10 ng/mL	15
	Cortisol	WB/S/P	80-800 nmol/L	10
Infection	AMH	S/P	0.02-10 ng/mL	12
	COVID-19 Ab	WB/S/P	Qualitative	10
	COVID-19 Ag	Nasopharyngeal swab	Qualitative	12
	COVID-19/Flu Ag Combo	Nasopharyngeal swab	Qualitative	20
	COVID-19 nAb	WB/S/P	Qualitative	20
	CRP	WB/S/P	0.5-200 mg/L	3
	PCT	WB/S/P	0.1-100.0 ng/mL	12
	PCT Plus	WB/S/P	0.02-50 ng/mL	12
	IL-6	WB/S/P	2-2,500 pg/mL	12
	HBsAg	WB/S/P	Qualitative	12
	Anti-HBs	WB/S/P	Qualitative	15
	Anti-HCV	WB/S/P	Qualitative	12
	NORO	Feces	Qualitative	12
	ROTA	Feces	Qualitative	12
Autoimmune	Rota/Adeno	Feces	Qualitative	12
	Anti-CCP Plus	WB/S/P	3.2-300 U/mL (Semi-quantitative)	12
Other test items	Total IgE	WB/S/P	1-1,000 IU/mL	10
	Ferritin	S/P	10-1,000 ng/mL	10
	Vitamin D	S/P	8-70 ng/mL	28
	Infliximab	WB/S/P	0.2-50 µg/mL	10
	Total Anti-Infliximab	S/P	8-250 AU/mL	22

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8 p.


AFIAS parametrai

Kategorija	Elementas	2 p.	Mėginio tipas	3 p.	Darbinis intervalas	4 p.	7 p.	Testo laikas (min.)
Širdis	Tn-I Plus		WB/S/P		0.01-15.00 ng/mL			12
	CK-MB		WB/S/P		3-100 ng/mL			12
	D-Dimer		WB/P		50-10,000 ng/mL			12
	NT-proBNP		WB/S/P		10-30,000 pg/mL			12
	Mioglobinas		WB/S/P		2-500 ng/mL			12
	Širdies trigubas		WB/S/P		Tn-I: 0.01-15 ng/mL CK-MB: 3-100 ng/mL Mioglobinas: 5-500 ng/mL			12
	ST2		WB/S/P		3.1-200 ng/mL			12
	BNP		WB/S/P		10-20,000 pg/mL			12
	Troponinas T		WB/S/P		5-5,000 pg/mL			12
Vėžys	PSA		WB/S/P		0.5-100.0 ng/mL			15
	PSA Plus		WB/S/P		0.07-50 ng/mL			12
	AFP		WB/S/P		5-350 ng/mL			15
	AFP Plus		WB/S/P		0.5-350 ng/mL			12
	CEA		S/P		1-500 ng/mL			12
	CEA Plus		WB/S/P		0.5-200 ng/mL			12
Diabetas	HbA1c		WB		(NGSP) 4-15%			10
	Mikroalbuminas		Šlapimas		2-300 mg/L			12
Hormonas	TSH		WB/S/P		0.09-80.0 μIU/mL (WB: 0.4-80.0 μIU/mL)			15
	TSH Plus		WB/S/P		0.08-50 μIU/mL			12
	T3		S/P		0.77-7.70 nmol/L			10
	T4		S/P		10.23-300.00 nmol/L			10
	FSH		WB/S/P		1-100 mIU/mL			12
	Total βhCG		WB/S/P		5-50,000 mIU/mL			15
	β-HCG Plus		WB/S/P		2-1,500 mIU/mL			12
	LH		WB/S/P		1-100 mIU/mL			12
	Progesteronas		S/P		4.45-127.2 nmol/L (1.4-40 ng/mL)			15
	PRL		WB/S/P		1-100 ng/mL			10
	Testosteronas		S/P		1-10 ng/mL			15
	Kortizolis		WB/S/P		80-800 nmol/L			10
	AMH		S/P		0.02-10 ng/mL			12
Infekcija	COVID-19 Ab		WB/S/P		Kokybinis			10
	COVID-19 Ag		Niosiaryklės tepinėlis		Kokybinis			12
	COVID-19/Flu Ag Combo		Niosiaryklės tepinėlis		Kokybinis			20
	COVID-19 nAb		WB/S/P		Kokybinis			20
	CRP		WB/S/P		0.5-200 mg/L			3
	PCT		WB/S/P		0.1-100.0 ng/mL			12
	PCT Plus		WB/S/P		0.02-50 ng/mL			12
	IL-6		WB/S/P		2-2,500 pg/mL			12
	HBsAg		WB/S/P		Kokybinis			12
	Anti-HBs		WB/S/P		Kokybinis			15
	Anti-HCV		WB/S/P		Kokybinis			12
	NORO		Išmatos		Kokybinis			12
	ROTA		Išmatos		Kokybinis			12
Autoimuninės ligos	Rota/Adeno		Išmatos		Kokybinis			12
	Anti-CCP Plus		WB/S/P		3.2-300 U/mL (pusiau kiekybinis)			12
	Bendras IgE		WB/S/P		1-1,000 IU/mL			10
Kiti testo elementai	Feritinas		S/P		10-1,000 ng/mL			10
	Vitaminas D		S/P		8-70 ng/mL			28
	Infliksimumas		WB/S/P		0.2-50 μg/mL			10
	Bendras anti-infliksimumas		S/P		8-250 AU/mL			22

1 p.

AFIAS-10

AFIAS-10

The warranty for Boditech Med Inc. applies only if the product is used in accordance with the method provided by Boditech Med Inc.

In no event shall Boditech Med Inc. be liable for any direct or indirect damages incurred.

For technical assistance, contact:

Boditech Med Inc. Technical and Service Support Team

Tel: +82 33 243 1400, E-mail: TS@boditech.co.kr

Boditech Med Inc.

43, Geodudanji 1gil, Dongnae-myeon,

Chuncheon-si, Gang-won-do 24398, Republic of Korea

Phone: +82 33 243 1400 / Fax: +82 33 243 9373

www.boditech.co.kr

Obelis s.a

Bd. Général Wahis 53, 1030 Brussels, BELGIUM

Phone: +(32) 2.732.59.54

Fax: +(32) 2.732.60.03

E-Mail : mail@obelis.net

The actual content on the screen may vary depending on the specifications and customer requirements.

This product is an in vitro medical device.

contents

1. General warnings and important guidelines	4
2. Works	7
3. Components	8
4. Technical specifications	10
5. Functional and behavior elements	11
6. installation	14
7. Perform AFIAS-10testing	16
8. Other support features	52
9. AFIAS-10 Preferences	63
10. Error Message	67
11. Maintenance, repair and disposal	69
12. Warranty	72

1. General warnings and important guidelines

Graphic symbols

The following list describes the AFIAS-10 operating manual, AFIAS-10 product labels and packaging boxes, and the various markings available in the packaging design.

representation	explanation
	Manufacturer
	Please refer to the user's service order.
	In vitro eddyesting device
	Remove and dispose of electrical and electronic components.
	Biological risks
	USB devices
	Temperature limit
	Serial Number
	Laser radiation Don't look at the light.
	Potential biological risks
	This mark is used when special warnings or precautions are required related to the use of devices that are not marked on the label. This indicator is also used to mean "Please, see the user's note"
	This indicator is used to indicate situations where you are not careful or do not follow instructions, or if you are displaying to users of the product or people around you that you may be injured by electrical shock or UV exposure.

Electrical hazards!

If you use AFIAS-10 in an inappropriate way or do not pay attention to it, you may experience electric shock, burns, fire or other dangers. To reduce these risks, please follow these steps:

- Use only the supplied powercord, such as AFIAS-10.
- The AFIAS-10 does not have any parts that can be repaired by the user. Never disassemble the AFIAS-10.
- In any case, do not let water or liquids touch the AFIAS-10.
- If cleaning is required, wipe the cartridge holder and exterior with a soft cloth. Remove the power cord before cleaning.
- If the AFIAS-10 is stored unused for a long time, remove and store the power cord.

Potential biological risks

Patient samples and used cartridges should be treated as potentially infectious substances. To reduce biological risk, follow these steps:

- Dispose of the cartridges used in accordance with biological waste disposal regulations.

If your AFIAS-10 has a full life, please send it to Boditech Med Inc. for proper disposal or recycling. If you are disposing of the AFIAS-10 on its own, please dispose of it in accordance with the electrical/electronics and component disposal regulations.

Always dispose of biological waste in accordance with local ordinances

Laser radiation

The AFIAS-10 is a Class 1 laser device or an internal laser module is Class 3R. Since the AFIAS-10 has a built-in laser module, laser radiation, emission can pose a risk; to reduce this risk, please follow the following:

- Never retrofit or disassemble the AFIAS-10. This instrument may be exposed to laser radiation due to the built-in laser module and laser window. Use the procedures in the user's manual when using and operating.

- The AFIAS-10 cannot see the racer outside the product during normal use, so no eye protection is required.
- Do not look inside the machine during testing. Looking at the laser beams inside through optical equipment (e.g. microscopes, telescopes) can be dangerous.
- All exterior cases of the product protect the user from laser rays, so never open the exterior in case of the product.
- Laser radiation intensity cannot be adjusted other than the manufacturer.

Attention

The AFIAS-10 is used in laboratory environments.

To reduce AFIAS-10 damage, follow these steps:

If not follow the below guideline, it will not be valid the warranty.

- To prevent damage to internal components, please take care when using the AFIAS-10.
- Use in a clean, dust-free place.
- Do not expose AFIAS-10 to excessive humidity and high temperatures.
- Do not place any objects on top of the AFIAS-10.
- Do not use anything other than the products and accessories we offer.
- If you do not have any special instructions, do not insert anything into the sockets, drawers, etc. openings of the AFIAS-10.
- When transporting or storing, store AFIAS-10 in a dry place between -20 and 60°C to prevent mechanical defects or damage.

To avoid errors in test results, please follow these instructions:

Failure to follow these precautions will invalidate the warranty.

- Before testing, make sure your AFIAS-10 is in the proper environment.
- Tests should be prepared to be conducted at 15-35°C.
- Cartridges for testing must be stored under the storage conditions recommended by us and must not be used past their expiration date.
- Before testing, make sure that the test name to be measured, the lot number of the cartridge, and the lot number of the ID Chip match exactly.
- Do not remove the test cartridge from the aluminum pouch until the test is carried out, leave at room temperature for at least 30 minutes before use.
- Before testing, be sure to read the enclosed instructions carefully and test according to the instructions.

2. Works

Purpose of use:

The AFIAS-10 is a medical immunofluorescent photometer that can be used with a variety of fluorescent immuno-immunospatio test cartridges provided by Boditech Med Inc. The sample used can be blood, urinary and other samples, and the type of sample depends on the type of test.

Principle

AFIAS-10 can be used with a variety of fluorescence immune test cartridges based on antigen-antibody response and fluorescence immunoassay.

AFIAS-10 uses a semiconductor diode laser as a light source to investigate lasers in cartridges where samples are placed; lasers are investigated on fluorescent labeled analytical materials gathered in the cartridge's test line, and fluorescence light is emitted from these materials.

The fluorescence light is collected with the laser and only the fluorescent light is filtered using a filter. The intensity of the collected fluorescent light is converted into electrical signals.

The microprocessor mounted on the circuit board uses the formula information stored in the ID-Chip to calculate the concentration of the analysis.

The test results, converted and calculated, are quantitatively displayed on the screen of the AFIAS-10.

AFIAS-10 testing is used for screening, monitoring, and routine testing in laboratories such as hospitals and medical institutions.

3. Components

The AFIAS-10 and accessories are provided in one box. For proper use of the product, the user must familiarly know the technical terminology and how to use each component.

The user must verify that the following components are included in the box: If any components are missing, please call our customer support representative at 033-243-1400.

	AFIAS-10	1EA
	Power code	1EA
	Thermal Printer paper	1EA
	BLOOD TUBE Rack	1EA/RACK
	Waste bin	1EA
	User's Manual	1EA

The following components are available as an option (supplied on demand): For the TUBE Rack, one basic BLOOD TUBE Rack is provided and can be exchanged for another TUBE Rack when purchasing the product. If additional TUBE Rack required, it can be for purchase.

	Eppendorf (1.5ml) tube mount RACK	1EA/RACK
	Eppendorf (0.5ml) tube mount RACK	1EA/RACK
	mini blood tube mount RACK	1EA/RACK
	Barcode reader	1EA
	USB KEY BOARD	1EA
	SD CARD	1EA

4. Technical specifications

Exterior specifications

- ☐ **size** 426 mm (W) x 443 mm (D) x 395 mm (H)
- ☐ **weight** 20 Kg
- ☐ **Power** 100-240 V AC, 50/60 Hz, 2.6-1.0 A
- ☐ **Data output** LCD / Built-in Printer / RS-232 Port

Experience

- ☐ **Use temperature** 15 ~ 35 °C
- ☐ **Storage temperature** -20 ~ 50 °C
- ☐ **humidity** 10 ~ 70 %
- ☐ **environment** Avoid direct sunlight, mechanical vibration and magnetic fields, and install on dry and clean, flat horizontal surfaces

Optical system

- ☐ **illuminant** Laser Diode
- ☐ **receiver** Silicone photo diodes

Other specifications


- ☐ **interface** USB 3 Ports
LAN 1 Port
LCD / Built-in Printer / RS-232 Port
SD Card Slot 1 Port
Mini USB 1 Port
- ☐ **display** 10.1" TFT LCD(touchscreen)
- ☐ **fuse** T4L-250V
- ☐ **laser** Laser rating: Class 1
Radiation wave:635 nm
Radiate Angle: Equilibrium Light
Impulse cycle: -
Maximum radial power: < 1 mW

*****AFIAS-10 meets the standards EN 61326-2-6, CHINA: GB4793. 1-2007 & YY 0648-2008 & GB7247. 1-2012 & GB/T18268.1-2010 & GB/T18268.26-2010 & GB/T14710-2009.**

***** AFIAS-10 is a raser analyzer following IEC 60825-1 : 2007**


5. Functional and behavior elements

Front


①	LCD	Display equipment operation and measurement results, error messages, etc.
②	Front Power Button	Power button located on the front, rear power switch on when equipment On/Off function
③	ID Chip Socket	ID Chip insert port with built-in information from item that you want to check
④,⑤	Front cover	Cover for Bay protection when equipment is not used
⑥,⑩	Sample Rack Bay	Bay with Rack with Sample Tube to inspect
⑦	Waste bins	Cartridges and tip outlets discarded after test
⑧	Left Cartridge Bay	Left Bay with cartridges and consumables for up to 5 tests
⑨	Right Cartridge Bay	Cartridges and consumables for up to 5 tests can be fitted right Bay

Left side


Rear


①	Power switch	Equipment Mains Power Switch
②	Power code	Power Cable Connections
③	LAN port	AFIAS-10 and LAN connection
④	SD card	Storage of Data and equipment upgrades (for manufacturers, up to 32 G)
⑤	USB port	Connect barcode readers, etc.
⑥	Mini USB	Developer port (for manufacturers)
⑦	RS232	AFIAS-10 and computer integration

Top

Printer	Output test results of equipment
---------	----------------------------------

Ba y internal structure

Cartridge Holder	Where to mount cartridges for test (up to 5 in onebay)
Tip Holder	Where to mount the Tip for test

6. installation

Before using AFIAS-10, please read and use the warnings and precautions in this manual chapter 1.

Powerrequirements:

The AFIAS-10 works on commercial 100-240 V AC, 50to 60 Hz, 2.2-1.1 A. If power information is uncertain, contact your administrator.

The useof the affected:

The AFIAS-10 is designed for indoor use and must be used in a dry, clean, flat, horizontal place; in a place unaffected by direct sunlight, vibration, and electronic waves.

Adjust the room temperature to a range of 15 to 35 °C.

Installationprocess:

Remove the component from the packaging and check for physical damage or missing components.

1. Connect the power cord to the power connector on the back of the appliance.
2. Connect the power cord to the power point and press and turn on the power switch at the bottom of the rear of the machine.
3. After turning on the rear switch, check the front switch for red light.
4. Pressing the front switch will turn on the equipment with a blue light.


(Becareful! If the front switch is switched on, the screen does not fail for 9-10 seconds.)


5. When powered on, the AFIAS-10 briefly displays the Boditech Med logo and begins system startup and self-check.


6. After completing the self-check, the login window is displayed.

7. Once you're signed in, go through the device reset process.

8. The main menu will appear.

7. Perform AFIAS-10 testing

	Before testing, read the cartridge's package inserts in fine print.
	Before testing, make sure that the lot number printed on the test cartridge and ID Chip matches exactly.
	Cartridges must be stored under the storage conditions required by us. If you are using a product that is out of date, the test results may be error.
	If the test ends abruptly or fails while scanning the cartridge, an associated error code is displayed. For more information on the causes and solutions of error codes that occur, please refer to chapter 9 troubleshooting.

1) Preparing for testing

- To conduct the test, please prepare the following products:
- XXX ID Chip
- XXX Cartridge
- Patient sample or Boditech XXX control

(Note: XXX refers to the name of the Test to use.)

THE AFIAS-10 CAN PERFORM UP TO 10 TESTS OF THE SAME OR DIFFERENT (UP TO 10 TYPES) AT A TIME.

- Before pulling the cartridge out of the cartridge pouch, make sure that the lot number printed on the cartridge pouch matches the test name printed on the ID Chip.
- Leave patient samples, controls, and cartridges at room temperature for at least 30 minutes to ensure the appropriate temperature for testing.


2) Login /enter required information:

- As mentioned in Chapter 6, AFIAS-10 displays the login window.

(1) Select a username, and enter a password. The initial username and password are "Admin" and "0000". Users can change the password for "Administrator" and add other users. For more information, see section 8.2, User Management.

If the login is successful, the main menu is displayed.


(2) The description of each feature shown in the main menu is as follows:


LCD display control window: guidance, operation, status display and measurement result display to set the desired measurement environment (input via touch screen)

- Menu bar: button area to go Analysis/History/Sting
- Command Button: Button area where you can start checking and view stored ID Chip information
- Bay Info Area: Area showing settings and test progress information for the test consumables Cartridge and Tip
- Sample Rack Area : An area that displays the settings and information of The Sample to be inspected
- System Info Area: Information from the operator. Displays date time information, currently running S/W information of AFIAS-10

(3) Bay Info Area & Sample Rack Area

Each corresponding menu is configured to match the consumables on the machine.

- The bay can be inserted and released via a button below each bay screen, and the status of the current bay is displayed.

	Unloaded	The bay is ejected, Click to insert
	Loading	Bay being inserted
	Loaded	The bay is inserted, can be released by Long Click
	Unloading	Bay being released
	ReactionWaiting	Cartridge reaction waiting, Long Click to cancel
	Running	Test in progress

- Cartridge Bay: Cartridge information and status are displayed.
- Sinclude Rack Bay: Sample information and status are displayed.

(4) System Infor Area


3) Enter test information

(1) ID-Chip identification


Insert the ID Chip of the items you want to test into the ID Chip port.

Insert the ID chip into the ID Chip port as pictured above and touch the ID-Chip Button on screen.

When touched, the currently stored ID-Chip List is displayed and the Load button is pressed to load and save the mounted ID-Chip information.

For the same Item, up to 10 are stored per lot, and if the same lot is in the same Item, it is overwritten with the most recently inserted ID Chip information. Up to 100 ID Chip information is stored, and after 100, it is erased in order of oldest history. The ID chip list is displayed in the latest order and does not display expired ID-Chip information.

If you are checking multiple items at the same time, make a case of whether the ID Chip information exists in the ID Chip list. If it does not exist, the cartridge will not be identified, so please insert the ID-Chip to store the ID Chip information.

(2) Insert the cartridges

① Please insert the cartridges

As pictured above, align the cartridge in the direction and slide it all the way in.
If you are inspecting multiple items at the same time, it is efficient to equip items with similar reaction times in the same bay.

When inserting the cartridge, it must be inserted in the direction of the arrow.
Otherwise, the cartridge barcode error will occur.


When inserting the cartridge, push it all the way in.
If you do not push it all the way, you will experience operational problems with the equipment.


(O) (O) (O) (X) (X)

When the cartridge is fitted, the cartridge mounted on the screen changes to activated.

- ② Put the tip to the cartridge you have inserted.

	When using tip, you should always use the tip we provided; otherwise, it may happen the problem on the operation. Be sure to put the tip to the cartridge. If the tip is not fitted, a 'tip empty error' will occur during the test and the next test will take place.
---	--


※ Additional Tip


Generally, you will need 1 tip per item. However, some items (ex: Vitamin D, T3, T4, etc.) may require additional tip. You can equip up to 10 additional tips per Bay. After additional Tip mounting, match the number of tips fitted on the screen.

(For additional tips, check the Insert of Item.)


Insert additional tips(max 10 per bay)



③ Bay Load

Press the Load button on the Bay with cartridge and tip to read the Barcode on the cartridge and display the Item name of the cartridge you are wearing.

Identify the Barcode

Show the identified item

Complete Loading

	If the cartridge is not identified, please check the item name and Lot number in the ID-Chip information and please re-verify it. 
---	--


(3) Preparing samples

Four types of RACK(samples pickup type) are available for the AFIAS-10.




	Blood tube rack  3ml (75mm) 5ml (100mm)
	 1.5ml 1.5ml ependorf tube rack
	 mini blood tube rack
	 0.5ml 0.5ml ependorf tube rack


① Blood collection tube

Open the lid of the blood collection tube with the blood blood below.

	Use the blood collection tube rack we provided. If the popped side is positioned towards the user's body, it is number 1 from the right; it must be filled in order. A total of 10 blood collection tube can be used for one rack.
---	---

	
	When inserting into the machine, insert it with the side with the barcode positioned to the right of the machine. If the insertion direction is not correct, it may be a problem to drive the machine. 
	The size of the blood collection tube you use is only available for commercial 3ml and 5ml. • 3ml  • 5ml  Larger or smaller than the baseline can cause problems with machine behavior.
	When using the B blood collection tube, samples with a sample volume of 3ml tube can only be used for samples 1 ml or 3ml or less, and for 5 ml tube, samples of 1 ml or less can be used only; samples of 1 ml or less may not result in quantitative sample collection, which may affect the test results.
	If the internal sample has bubbles when using the Blood collection tube, the sample collection of the quantification may not be able to take effect, which may affect the test results.


	 (X) (O)
	If you use your own barcode on the Blood collection tube, Barcodes must be properly present on the blood collection tube for recognition. 

② Eppendorf tube (1.5ml)

Place a 1.5 ml Eppendorf tube with sample on a 1.5ml Eppendorf tube rack after the lid is fully lidded.

	If you are using an Eppendorf tube with a lid (1.5 ml/0.5 ml), place the lid in the groove as follows:
---	---

	   (O) (X) (X)
	Please open the lid of eppendorf tube completely.    (O) (X) (X)
	If you do not mount rack exactly when using Eppendorf tube, you may not be able to obtain samples of quantitation, which may affect the results of the test.
	Within a 1.5 ml Eppendorf tube, the sample titration is 40 ul + sample quantity. If it is more or less than that amount, the resulting value may be error.

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Remove the lidon the Mini blood tube containing the sample and mount it on the Mini blood tube rack.

	Mini blood tube rack and 1.5ml Eppendorf tube rack have the same appearance. Please be careful not to mix (you can use a 1.5ml Eppendorf tube on a mini blood tube rack, or a mix of mini blood tube and a 1.5 mlEppendorf tube on the same tub rack).
---	--

	 (X)	If you place mini blood tube on rack, be sure to remove the lid.
		The amount of sample in mini blood tube is 40 ul +sample quantity. If it is more or less than that amount, the resulting value may be error.

④ Eppendorf tube (0.5ml)

Place a 0.5 ml Eppendorf tube with sample on a 0.5ml Eppendorf tube rack after the lid is fully lidded.


	How to place a container within a 0.5 ml Eppendorf tube rack is the same as how to place a 1.5 ml Eppendorf tube. Please read the cautions in section 7.3.3.2(② Eppendorf tube(1.5ml)) once again.
	Within 0.5 ml Eppendorf tube, the sample titration is 40 ul +sample quantity. If it is more or less than that amount, the resulting value may be error.
	When using Eppendorf tube and mini blood tube, pay attention to the bubbles in the sample. Bubbles inside the sample may cause errors in sampling and errors in the resulting values.
Caution	The first tube rack is one Blood tube rack. For the 5ml Eppendorf tube rack, 0.5ml Eppendorf tube rack, mini blood tube rack, it is offered instead of the Blood tube rack or requires additional purchase.

- (4) Insert the prepared Tube rack into the Loadingpart of the Symbling station as shown.



It automatically loads when mounted, reads the mounted Rack Type, detects and displays the mounted Tube.

(5) Specify a test

① Automatic unspecified

If the Item type of cartridge is one type, it is automatically assigned to the cartridge in the tube order mounted on the Tube Rack.


② Manually selection


If you have more than one item type in Cartridge, or if you change it, touch Sample Rack Area on the screen to activate the sample assignments window.


- The Item List shows the item list of the inserted cartridges and the quantity that can be set for test.


- Select sample type, and then touch the item you want to specify in the item list.
 - The process can be set automatically when bidirectional (Bi-LIS) is used.
 - You can assign multiple items to a sample for test. (Up to 10 items can be specified)
 - When you assign an item to a sample, it is displayed in cartridge color.
- Example of assigning an item to a sample


e. x. 1) Match the blood and 'hsCRP' to the sample at no.1


e. x. 2) Full blood andhsCRP in the sample 1-3times, full blood and Tn-I Plus matching in the sample4-6 times


e. x. 3) HS CRP and Tn-I Plus matching to sample 7

- Delete the designated item
 - Select the item in the item list to change it to the 'remove' state, and then touch the sample to removeit. (If you add it again later, touch the item again to change it to the 'Add' state.)

- Insert Patient ID(Sample ID)
 - Product ID (Barcode): Reads and automatically displays the barcode attached to the Blood tube.
 - Patient ID (User): Enter any ID by touching the text box as needed.


4) Start testing

(1) Press the 'Start' button.

(2) In the message confirming the start of the scan, press the 'OK' button to start the scan.


(3) Test progress


- ① Once the test begins, the machine will automatically take samples and collect reaction cartridges. Tests are conducted by bay, and if both are set, the test is conducted from the left bay.


The Test starts from the left bay

- ② As the scan progresses, progress is displayed on the screen.
 - Cartridge color background: Waiting for test and waiting for reaction
 - Blue border: Sample acquisition and dosing
 - Red border: Optical system scanning and counting
 - White background: Test done


③ The information configuration displayed in the cartridge is as follows:

5) Restart of test (1)

The test of the AFIAS-10 is conducted by bay, and the sample rack is released at the end of all dosing for the test set in the sample rack.


If you have a bay that can start a scan, you can set up and start a new check in that bay.


6) Interruption of test


- If you want to stop the scan in progress, touch the Stop icon.


- If you see the following pop-up, press the 'OK' button to stop it:


- When stopping, stop all checks except those that are in response.

- To stop the scan being responsive, touch and hold the "Reaction Waiting" button.

- In the pop-up, touch the 'OK' button to stop the scan being responded to.

7) Urgent test

(1) Insert the cartridge for the required test as an emergency check during drive.


(2) The sample amount required for testing is put into the 'Sample Well' of the cartridge.


Please refer to the AFIAS insert for detailed testing methods, such as sample preparation and sample volume.


(3) After inserting the tip on the cartridge, load the cartridge into the screen.

(4) Clip the keyboard button to the identified cartridge.


(5) Enter test information and click the 'set' button.

	When taking samples with C-tip, check the correct sample volume; low or large volumes can lead to errors in tests.  (○) (X) (X)
	Do not shock the C-tip when the sample is taken and placed in the tip holder. Samples within the C-tip may be ejected.
	C-When using tip, it should be used within 2 minutes of taking the C-tip sample. If the sample in C-tip hardens and fails to pump after 2 minutes, it affects the result.

(6) After the current sample acquisition and dosing test, an urgent test follows.


Never put your hands inside the machine during operation. There is a risk of injury.


Do not place debris inside the machine. It may cause a failure.

8) AFIAS-10 Shutdown:

- (1) Before exiting the AFIAS-10, make sure that the IDchip, cartridge, tip, and tube rack are all removed.
- (2) Touch the lower left corner of the Main screen.


- (3) Touch the power off button on the pop up window that appears, and then touch the OK button.

- log out: Use if you want to change the operators of the equipment.
- Reboot: Use to turn the device off and on for a while.

- (4) Click 'OK' will automatically shut down the machine after the shutdown operation.

- (5) When exiting the machine, please cover the front area and store it.
- (6) If the equipment is not used for more than 2days, turn off the power switch on the rear and remove the power cord from the power point and store it.
- (7) Disconnect the power cord from the AFIAS-10.

Please exit in the order suggested. In addition to the suggested shutdown method, internal data memory may be erased if you exit with the power button located on the back.
If the internal data memory is cleared, please contact your nearest dealer or our technician for action.

Before testing, make sure that the name of the item you want to test matches the lot number of the cartridge with the lot number of the ID Chip.

The AFIAS-10 differs from the rack type and seal type available depending on the sample carrier being operated. Below shows the rack type and sample type according to the operating instrument.

Operating vessels	Operating rack type	Sample type	Remarks
Blead collection tube	Blood collection tube rack	Whole blood (Venus blood) Serum,Control	3 ml(or 5ml) general blood collection tube rack Minimum volume 1 ml Maximum volume 3ml(or 5ml)
1.5 ml Eppendorf tube	1.5 ml Eppendorf tube rack	Whole blood (Venus blood) Serum,Control	Maximum volume 1.5ml Minimum volume 40 μl
0.5 ml Eppendorf tube	0.5 ml Eppendorf tube rack	Whole blood (Venus blood) Serum,Control	Maximum volume 500 μl Minimum volume 40 μl
mini bleed tube	Mini BLead Tube Rack	Whole blood (Venus blood) Serum,Control	Maximum volume 500 μl Minimum volume 40 μl


8. Other support features

In addition to conducting tests to ensure efficient and reliable performance, as well as proper reporting and recording of results, The AFIAS-10 includes features such as QC, system check, and view/management of stored test results.


1) QC

- QC(Degree Management)
 - To proceed with QC, a certain degree of management material (control) is required. Degree management material (control) is offered according to the user's needs.
- QC conducts the following cases:
 - Cartridge lot changed.
 - The test results are thought to be abnormal or measured differently from the reference value.
 - The test results do not match the symptoms.
 - System training purposes


(1) Press the Analysis button to switch to QC mode.

(2) Press the Control Material button to open the Control Material window.

(3) If the control material is not registered, press add to add the control (2 level each).

LotAACQK16_LV1Exp.2021-11

HbA1cHbA1cMean5.30Low4.51High6.10

RegisterClose

LotAACQK16_LV2Exp.2021-11

HbA1cHbA1cMean10.30Low8.76High11.85

RegisterClose

QCControl MaterialHistorySetting

Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACQK16_LV2	2021.11.30		9.9	8.42	11.39
HbA1c	AACQK16_LV1	2021.11.30		5.3	4.51	6.1


AddClose

Operator: admin
L : 35-13°C, 0.009
R : 35-00°C, 0.009
2021.03.12 11:31:22
A10: Boditech
AA10: 30, 30, 30, 30
PA10: 30, 30, 30, 30


(4) Insert the control material and tip to be inspected in the cartridge to be used.


(5) Press the load button to identify the cartridge.


(6) If there are more than one cartridge, Lv. 1 and Lv. 2 are automatically granted.


(7) Press the keyboard button on the cartridge to set the control material.


(8) Press the Start button to proceed to QC.


(9) QC results are available in History and ID-Chip information.

- History

- ID-Chip

2) View/print/manage saved test results:


- You can store up to 10,000 AFIAS-10 test results in internal memory.
- When you touch the History tab in the main menu, you will see the results of the currently saved tests. Each test result is displayed, including patient ID, user, test type, results, test date, and LOT information.
- If you want to print any of the test results displayed, select the result and touch the Print button.
- If you want specific test results, you can tap the Search button to enter the patient ID, user, test type, test date, and so on to view the results you want.

(1) History screen

Analysis							
No	Patient ID	Item	Sample Type	Result	Date	Lot	Print
24	210311 024 U 000 B 00010	hsCRP BloodTube	WB	5.65 mg/L	2021.03.11 12:29:12	CAPLR01A	Search UIS
23	210311 023 U 000 B 00009	hsCRP BloodTube	WB	5.68 mg/L	2021.03.11 12:28:22	CAPLR01A	
22	210311 022 U 000 B 00008	hsCRP BloodTube	WB	5.58 mg/L	2021.03.11 12:27:25	CAPLR01A	
21	210311 021 U 000 B 00007	hsCRP BloodTube	WB	5.58 mg/L	2021.03.11 12:26:29	CAPLR01A	
20	210311 020 U 000 B 00006	hsCRP BloodTube	WB	5.77 mg/L	2021.03.11 12:25:41	CAPLR01A	
19	210311 019 U 000 B 00005	hsCRP BloodTube	WB	5.65 mg/L	2021.03.11 12:24:42	CAPLR01A	
18	210311 018 U 000 B 00004	hsCRP BloodTube	WB	5.44 mg/L	2021.03.11 12:23:49	CAPLR01A	
17	210311 017 U 000 B 00003	hsCRP BloodTube	WB	5.62 mg/L	2021.03.11 12:22:47	CAPLR01A	
16	210311 016 U 000 B 00002	hsCRP BloodTube	WB	5.78 mg/L	2021.03.11 12:21:51	CAPLR01A	
15	210311 015 U 000 B 00001	hsCRP BloodTube	WB	5.30 mg/L	2021.03.11 12:21:02	CAPLR01A	
Operator: admin				L : 35.00°C, Good R : 35.00°C, Good		2021.03.11 13:41:36	A10: Bedi tech AA10: XX.XX.XX.XX FA10: XX.XX.XX.XX

- Check details: You can check the details by 1-2 seconds long (1-2seconds) to see the details.

Analysis							
No	Patient ID	Item	Sample Type	Result	Date	Lot	Print
24	210311 024 U 010 B 00010	hsCRP BloodTube	WB	5.65 mg/L	2021.03.11 12:29:12	CAPLR01A	Search UIS
23	210311 023 U 000 B 00009	hsCRP BloodTube	WB	5.68 mg/L	2021.03.11 12:28:22	CAPLR01A	
22	210311 022 U 000 B 00008	hsCRP BloodTube	WB	5.58 mg/L	2021.03.11 12:27:25	CAPLR01A	
21	210311 021 U 000 B 00007	hsCRP BloodTube	WB	5.58 mg/L	2021.03.11 12:26:29	CAPLR01A	
20	210311 020 U 000 B 00006	hsCRP BloodTube	WB	5.77 mg/L	2021.03.11 12:25:41	CAPLR01A	
19	210311 019 U 000 B 00005	hsCRP BloodTube	WB	5.65 mg/L	2021.03.11 12:24:42	CAPLR01A	
18	210311 018 U 000 B 00004	hsCRP BloodTube	WB	5.44 mg/L	2021.03.11 12:23:49	CAPLR01A	
17	210311 017 U 000 B 00003	hsCRP BloodTube	WB	5.62 mg/L	2021.03.11 12:22:47	CAPLR01A	
16	210311 016 U 000 B 00002	hsCRP BloodTube	WB	5.78 mg/L	2021.03.11 12:21:51	CAPLR01A	
15	210311 015 U 000 B 00001	hsCRP BloodTube	WB	5.30 mg/L	2021.03.11 12:21:02	CAPLR01A	
Operator: admin				L : 35.00°C, Good R : 35.00°C, Good		2021.03.11 13:34:02	A10: Bedi tech AA10: XX.XX.XX.XX FA10: XX.XX.XX.XX


- Print of save results

Touch the result you want to output (v: activate the print mark red) and click 'Print'.


Click 'OK' in the pop-up that appears to print the selected results.


- Search for saved results

If you have saved results that you want to search for,click 'Search'on the right.

When you type a search terms in the pop-up window that appears, they display the results, including the corresponding strings.


- Transfer of save results

Touch to select the result of the LIS (v: activate the print mark red), and then click 'LIS'.

Click 'OK' in the pop-up that appears to send LIS the selected results.

9. AFIAS-10 Preferences

The AFIAS-10 can set the environment using the following settings menu.

1) System Settings


The screenshot shows the 'System' settings page of the AFIAS-10. The left sidebar has tabs for 'System', 'User management', 'Version', and 'LIS'. The 'System' tab is active. The main content area contains the following settings:

- Date:** 2021.03.11, with a 'Set' button.
- Format:** MM.DD.YYYY, DD.MM.YYYY, YYYY.MM.DD (radio buttons).
- Time:** 14:09:24, with a 'Set' button.
- Time Zone:** Asia/Seoul, with a 'Set' button.
- Language:** Korean, English, Chinese (radio buttons).
- Print:** Auto, Manual only (radio buttons).
- Log-out Time:** 00:00, 30:00, 60:00, 90:00, 120:00, 180:00 (radio buttons).

The bottom status bar displays: Operator admin, L: 35.00°C, 60% H: 34.80°C, 60%, 2021.03.11 14:09:24, and version information: A10: Boditech, A10.XX.XX.XX.XX, FA10.XX.XX.XX.XX.


- From the main menu, touch the Setting tab, and then touch the System tab to go to the System Settings menu.
- To set the date, tap the 'Set' button on date. You will see a screen for setting the date, and select the date you want to set.
- To change the date display format, select the format you want in Date Format.
- To set the time, touch the Set button on time. The screen for setting the time is displayed, and set the time you want to set.
- To change the language, select the language type that you want to change to Language.
- Touch the Printer button to set up the printer. You can select Auto or Manual only to set whether to print automatically during test.
- You can set the automatic log-out time for long periods of unused non-entry. Selecting '00:00' does not use automatic logout.

2) Manage users


- From the main menu, touch the Settings tab, and then touch the User manage tab to go to the User Management menu.
- Users can only be added and deleted if they are logged in with Administrator ("Admin"), and if you are logged in as a regular user, the Add, Delete button will not appear.
- To add a user, touch the Add button and enter the user ID and password to register a new user.
- To delete a registered user, select the user you want to delete and tap Delete to delete it. If you want to delete the entire user, select the full selection box at the top and touch Delete to proceed with the entire deletion. (In this case, administrators are not deleted.)
- If you are logged in with a regular user ID, you can change your own password.

3) Versioning

- From the main menu, touch the Settings tab, and then touch the User manage tab to go to the versioning menu.
- You can check the operating environment and version of AFIAS-10, such as Firmware and Software.
- In the event of a future update, you can proceed with the update via our SD-card or USB.


4) LIS Management

- From the main menu, touch the Settings tab, and then touch the LIS tab to go to the versioning menu.
- You can select the LIS type, Protocol, and interface that are sent.

- If you are using Serial communication, you can set the Baud rate. (Base 115,200)
- You can set the address of the LIS server when using TCP/IP communication, and disable DHCP if you are using static IP for your device.
- When checking, you can select the LIS automatic transport function active in Export.

5) Setting menu tree

10. Error Message

1) Dialog Error Message

Error Message	Cause	Solution
Printer Error No Paper	Print in no printer paper or printer cover not to be closed well	Check printer paper or printer cover
Printer Error No Reponse	Problem with printer communication	Contact on the service center if there is continuous error.
Abnormal Temperature, Test banned	Start test not to reach the appropriate temperature in the analyzer	Restart the test when the analyzer reach in appropriate temperature
Can't start	Limit the test in case there is a severe error in Self Test	Contact on the service center
Unsupported Sample Type	Use a sample type that is not supported by the Rack you have put it(ex) Use Swab sample in Blood Tube Rack)	Use the sample type that is supported by Rack
Expired ID chip	Insert the expired ID Chip	Do not use the expired ID Chip
Invalid ID chip	Insert the ID Chip that is not matched	Do not use the ID Chip
All cartridges are incubating. Please wait	Press the Button(Bay Load/Unload) when there is no empty slot and all of the cartridges are tested	Do not Load/Unload the bay when all of the bays are tested
Extra Tip 1~X	Insufficient Extra Tip when assigning test items	Replenish the Extra tips and try again
The test is still in progress.	Attempt a firmware update while test is progressed	Update a firmware after all of tests are completed
No rack found	load the rack without inserting it	Try again after insert the rack
Duplicated User Name	Create the same user ID	Do not create the same user ID
Failed to remove the tip.	Do not remove used tips	If the waste bin is full, empty it and try again. If there is a continuous error, please contact the service center

2) Cartridge Error Message

Error Message	Cause	Solution
Invalid 25~32	Cartridge Error	Please check with the new cartridge. If you have continuous errors, please contact our service center.
Invalid	Cartridge error or poor dosing position error	
Value Invalid 1	Cartridge Error	
Flow Invalid	Poor dosing position or pump failure	
Strip Invalid	Cartridge Loading Line Error	
Cancel	Show at user shutdown Under test	The cartridge that the user forced to cancel the scan.
Tip Required	No Tip	Replenish the Extra tips and try again
Sample Invalid	Sample Detect Error	Check the sample and Try again
Scan Invalid	Optical system drive error	If you have continuous errors, please contact our service center.
Unknown	Do not match with Cartridge Barcode and ID Chip Lot	Register the ID chip first


11. Maintenance, repair and disposal

1) Maintenance:

- Mounting print media on the built-in printer
- Emptying waste bins

To install and replace the print media on the built-inprinter:

- Open the printer cover as shown below, insert the printed paper properly, and close the printer cover.

- Place the appropriate print paper.

Recommended paper

Paper Type : Infection Roll

Paper width: 57 mm

Paper thickness: 60±5µm

	When the printer cover is closed, the printer reset is automatically run when it is functioning normally, at which time the printable paper is discharged.
	When opening and closing the printer cover, do not use any force when replacing the print media. This may cause failure.
	Make sure the print paper is of the appropriate specifications. Do not use paper with the central part of the printing paper glued together. The built-in printer does not have the ability to check if the paper is used up.
	While mounting the print paper, make sure that the printer is not receiving any data. Otherwise, the data will be lost.

	If the previous test results are printed and the print paper is not printed due to lack of print paper, the old results will be printed automatically when a new print paper is fitted. However, if you turn off power during printing, unprinted data is lost.
	If the printer cover is not closed normally, no data will be printed.

Emptying wastebins:

- As the test progresses and completes, the used tip is discharged into the waste bin.
- If the tip used for the waste bin is full, remove the waste bin and discard the used tip.

	When handling biological waste, you must follow the relevant regulations.
---	--

2) Repairs:

- The AFIAS-10 does not has parts that can be repaired by the user.
- If you need technical repair or maintenance, please send it to us or contact your local distributor.

3) Disposal:

- If the AFIAS-10 isrequired to be disposed of, it must be disposed of inaccordance with the regulations on the disposal of electrical and electronic products in each region.

 	Important! For proper disposal and recycling of the life-long AFIAS-10, please send us an AFIAS-10 to dispose of. If you dispose of it on its own, it must be disposed of in accordance with the relevant regulations.
---	--

4) Transportation and storage:

- Use the original packaging box when shipping the product. Also, please use the original packaging box even if it is not used for a long period of time.

When transporting or storing AFIAS-10, store it in a dry place and avoid physical impact.

12. Warranty

AFIAS-10 has passed strict quality assurance and testing procedures.

Boditech Med Inc.'s expressed and implied warranties are conditional upon full observance of Boditech Med Inc.'s published directions with respect to the use of Boditech Med Inc.'s products.

To avail the warranty service, the consumer must return the defective AFIAS-10 unit or any component thereof; along with proof of purchase to Boditech Med Inc. or the nearest Authorized Service Agent.

Warranty Information

In the event of a manufacturer's defect or failure within the warranty period (24 months from the date of purchase), we will compensate you in accordance with the free repair or consumer damagecompensation regulations for such products or parts.

In the following cases, the warranty period is also covered by the warranty.

- Failure due to user's mishandling and mishandling
- Failure due to consumer will or negligence
- Failure and damage caused by repairs and modifications other than our services
- Loss of product warranty and arbitrary change of mention
- Failure due to natural problems

For after-sales inquiry, please contact:

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, REPUBLIC OF KOREA

Phone: +82 243-1400 / Fax: +82 243-9373

Email: TS@boditech.co.kr / Website: www.boditech.co.kr

Product Certificate

Thank you for purchasing AFIAS-10.

Please fill the required information and send a copy of this warranty card within 30 days of date of purchase to the above address.

Product Name		AFIAS-10	
Serial No.			
Date of Manufacture			
Warranty Period		24 months from the date of purchase	
Date of Purchase	Year	Month	Day
Manufacturer		Boditech Med Inc.	
Purchase Location			

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

The top two screenshots show individual test result screens. The first screen is for Lot AACQK16_LV.1 with a Mean of 5.30, Low of 4.51, and High of 6.10. The second screen is for Lot AACQK16_LV.2 with a Mean of 10.30, Low of 8.75, and High of 11.85. Both screens have 'Register' and 'Close' buttons.

The bottom screenshot shows the 'Control Material' history table with the following data:

Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACQK16_LV.2	2021.11.30		9.9	8.42	11.39
HbA1c	AACQK16_LV.1	2021.11.30		5.3	4.51	6.1

The bottom screenshot also shows an 'Add' button and a 'Close' button. At the bottom of the screen, there is a status bar with the text 'Operating status: OK', '25.12°C', '10.00', '2021.03.12 11:31:22', and 'ATQ: 3000 TACTY', 'AA10: XX, XX, XX, XX', 'FA10: XX, XX, XX, XX'.

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Quick Operation Guide


AFIAS-10 Start-up


- 1)  Turn on the switch 'Power On/OFF'.
(The power switch is located at the bottom left of AFIAS-10 back)
- 2)  Press the power button on the front of the equipment to turn it on.
- 3) When turned on, AFIAS-10 briefly displays the Boditech Med logo and enters the login window.

① For the first log in, please insert the password **0000**.

The instrument will request to change the password after the first login.*

- 4) When the login is completed, after going through the device initialisation process the main screen is displayed.

- 5) On the main screen, check the temperature of the equipment at the bottom.


If it is confirmed that the equipment temperature is green (**35 ± 2 °C**), the equipment is ready for testing.

- ▶ To conduct the AFIAS test, prepare the following items:
- ▶ AFIAS ID Chip.
: If the LOT to be tested has not been registered before, proceed with registering the ID Chip.
Registration is needed only once per LOT.
- ▶ AFIAS cartridge.
- ▶ Patient sample or Boditech control.

* refer to separate manual for more detailed instructions.


Main Screen Features



ID Chip Registration

- ① Insert the ID Chip you want to register into the ID Chip port and touch the 'ID Chip' button on the main screen. When the ID Chip List window appears, the ID Chip information is automatically saved.

* Insert the ID Chip into the ID Chip port (side with the item name visible).


- ② When registering more than one item, insert the second ID chip while the ID chip List window is open and press the Load button to register any additional item.

* AFIAS-10 can register 100 ID chip information.

* When changing Lot, the Lot specific ID Chip needs to be registered.

Cartridge Structure


Pre-inspection Precautions

- ▶ Verify that the lot numbers printed on the test cartridge and ID Chip match.
- ▶ Verify the expiry date.
- ▶ Use cartridges immediately after opening packaging.
- ▶ If the cartridge was refrigerated, leave it at room temperature for at least 15 minutes before proceeding with the test.
- ▶ Use a separate tip for dispensing samples onto the cartridge well (inside the instrument only insert tips from the box).
- ▶ Do not use tips if they are curved or damaged.



Pre-inspection Preparations

[Tip Fill]


1) Tip Fill & Cartridge placing on Tip Station

- ① Insert the tip in the tip stations on the side of the holder.

- ② Touch the tip station part on the main screen and correctly insert the number of tips installed.

[Cartridge Inserting]

* When a cartridge is installed, the slot changes to activated.

① Insert cartridges.

② Insert the tips into the tip holder of the inserted cartridge.


③ Click the 'Load' button on the holder where the cartridges are inserted.
(AFIAS-10 will read the cartridge barcode and display the cartridge item name).

After checking the barcode, the corresponding item will be shown.

2) Specimen and Rack Preparation

[Sample preparation]

* Start from the right, where the rack barcode is located.

① Blood tube rack.


(A) Prepare the blood tube rack and specimen.

* Blood tube available size: 13mm x 75~100 mm

* Minimum sample volume in the tube: 1mL (Dead volume) + a

(B) Open the lid of the blood tube and place the tubes on the rack.

* In the case of a whole blood sample, before taking off the lid slowly invert twice up and down then proceed. (Be careful not to create bubbles in the tube when inverting)

* When attaching a patient barcode, position it vertically in a similar way as the container's own label.

* Sample barcode recognition area: 1cm lower than the smaller tube edge to 1cm higher than the smaller tube bottom.
* smaller tube size: 75mm.

② Eppendorf tube rack

(A) Prepare the Eppendorf tube rack and the Eppendorf tube containing the specimen.

Minimum sample volume in the tube: 40 uL (Dead volume) + α (refer to the product insert for the volume)


(B) After fully opening the lid of the Eppendorf tube, place it in the rack following the instruction below.

* Start from the right, where the rack barcode is located.


[Installing the Specimen Tube Rack]

① Press the rack 'unload' button on the bottom right corner, in order for the rack holder to be ejected.

② Insert the tube rack.

*When loading the rack, it automatically detects the presence of the tubes and displays them.

Test Method 1) Automatic order designation (bidirectional LIS)

[Place the cartridges first]

[place the sample tube first]

① After placing the cartridge in the holder, select "unloaded". **Unloaded**

② Place the tip into the tip position.

③ Open the lid of the sample tube and place in the sample tube rack.

→ The sample rack enters the equipment, automatically recognises the tube barcodes, select the order, and matches them directly to the cartridges.

① Open the cap of the sample tube and install the sample tube rack.


→ The sample rack enters the equipment, automatically recognises the tube barcodes, select the order, and the required quantity of cartridge per item appears in the pop-up window.

② After placing the requested cartridges, press "unloaded". **Unloaded**

→ The cartridges and the sample tubes will match automatically.

③ Insert the tip into the tip position.


[In case of changing the sample type is needed]



(A) Touch the sample rack area to open the sample window.
(B) Select sample type.

Blood tube **WB** **S/P**



(C) Select the sample tube you want to change the sample type.


D) Touch the 'x' button to save the entered information.


④ Touch the 'Start' button.

⑤ Select 'OK'.


Test Method 2) Manual Ordering


If not all cartridges are of the same item, it is required to manually match sample tubes and cartridges.

- ① When the cartridges have been loaded, select. Unloaded
- ② Insert tips in the tip positions.
- ③ Load the sample tube rack.
- ④ Touch the sample rack area on the screen to designate sample tubes to the cartridges.

* The item list shows the item name of the inserted cartridge (with available quantity).

⑤ Select a sample type. WB S/P


⑥ Touch the item from the list (select add).

⑦ Select a sample (tube) which corresponds to the selected parameter that needs to be tested.

* It is possible to match up to 10 items with one sample.


e.g.1) Match S/P (serum/plasma) and Tn-I to samples 1 and 2.


e.g.2) Match whole blood and CRP for samples 1 to 3 and whole blood for samples 4 to 6 and Tn-I Plus.



e.g.3) Matching S/P and NT-proBNP to CK-MB in samples 3 and 4.

* Remove designated items.


- ① Select the corresponding item from the item list part.
- ② Change to 'Remove'.
- ③ Touch and remove the sample (if you add an item again later, change to 'Add' and proceed).

- ⑧ If you wish to add a patient ID, it is possible to type it using the keyboard on the screen. When finished, touch the X button to save the information.

- ⑨ Select the button 'Start'.

- ⑩ Press 'OK'.



Test Method 3) Emergency Test

Emergency test allows to manually insert the sample, this can be done in case of small sample volume or emergency examination.


- ① Insert the cartridge of the test needed for an emergency check.
- ② Put the required sample amount into the cartridge sample well.



- ③ Insert the tip in the tip position.

- ④ Click the 'Load' button on the corresponding bay.

- ⑤ Click the keyboard button on the bottom of the cartridge.

- ⑥ Enter Patient ID (optional).

- ⑦ Select Emergency Test.


- ⑧ Select General Tip.

- ⑨ Select the correct sample type.


- ⑩ Press Set.


- ⑪ The emergency test will then acquire the priority and will be processed immediately.

Delete the Test

1) How to cancel while a test is processing


- ① Touch the icon 'Stop' → Select OK.


- ② AFIAS-10 will cancel all tests, with the exception of those that are already at the reaction time step.

2) Canceling an individual test

- ① Touch the center of the cartridge slot for at least 3 seconds → Select OK.


- ② Only the selected cartridge will be cancelled.

3) Canceling all bay tests

- ① Touch for at least 3 seconds the button under the Bay that you want to cancel → press OK.

- ② All tests will be cancelled.

Maintenance

1) Clean Tip removal area

* Slide on the left side of the instrument.

2) Cleaning and emptying the tip waste box

3) System check

- ① Remove any remains from the tip removal slide on the left side of the equipment using an alcohol solution.

- ① Remove the waste box attached to the side of the equipment, discard the tip, and clean the inside of the tip box.
- ② Attach the waste box in place using the magnet.

- ① Prepare the System check cartridge, System check ID chip
When to use it: When AFIAS-10 is first installed, when AFIAS-10 is moved to another location, when an error occurs, when the test result is considered abnormal or different from the reference value.


- ② How to use.

(A) Register the System Check ID chip.

(B) Insert the System check cartridge into the cartridge slot and press the 'Unloaded' button.


C) Press the 'Start' button to proceed with the inspection and check the result.

Register and Test Quality Control (QC)

- ① Press the **Analyzer** button at the top of the main screen to switch to **QC** mode.
- ② Press the button **Control Material** to open the control material window.
- ③ If the control substance is not registered or the control Lot has changed, press **Add** to add controls (repeat for each level) and press the 'Register' button.

↓


Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACQK16_LV2	2021.11.30		9.9	8.42	11.39
HbA1c	AACQK16_LV1	2021.11.30		5.3	4.51	6.1

- ④ Insert the cartridge to be used and press the **Unloaded** button to recognise the cartridge.
- ⑤ If 2 or more cartridges have been inserted, the values Lv1 and Lv2 will be assigned automatically.
If manual application of the control material is required, press,  to set the control material.
- ⑤ **(a)** Control Lot (Lv1 / Lv2) **(b)** Tip type **(c)** Sample type.

Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACQK16_LV2	2021.11.30		9.9	8.42	11.39
HbA1c	AACQK16_LV1	2021.11.30		5.3	4.51	6.1


- ⑥ Add the control material to be inspected into the sample well.
- ⑦ Add the tip into the tip hole.
- ⑧ Press the 'Start' button to proceed with the QC.
- ⑨ QC results can be found both in History and in ID-Chip information.

- History

No	Patient ID	Item	Sample Type	Result	Date	Lot	QC
✓ 6	210225 805 ZZZ	HbA1c	WB	4.6 %	2021.02.25 11:25:24	AAQGB18F	Pass
✓ 5	210225 805 111	HbA1c	WB	4.2 %	2021.02.25 11:12:42	AAQGB18F	Pass
✓ 4	210225 804 ZZZ	HbA1c	WB	Invalid	2021.02.25 10:29:52	AAQGB18F	Fail
✓ 3	210225 803 111	HbA1c	WB	Invalid	2021.02.25 10:28:16	AAQGB18F	Fail
✓ 2	210225 802 ZZZ	HbA1c	WB	5.2 %	2021.02.25 09:50:53	AAQGB18F	Pass
✓ 1	210225 801 111	HbA1c	WB	6.3 %	2021.02.25 09:48:18	AAQGB18F	Pass

- ID Chip

Item	Lot	Exp	Rev	QC
HbA1c	AAQGB18F	12.28.2022	0	Pass
Tn-I Plus	TNRD810G	12.28.2022	0	Required
hsCRP	CAPLR01A	12.28.2022	0	Required
PSA	PSRDA39F	12.02.2022	0	Required
Tn-I Plus	TNRD805G	11.29.2022	0	Required

Item	Lot	L1	L2	Status	Date
HbA1c	AAQGB18F	AAQGB18F.LV1 (4.31 ~ 6.1) 5.0	AAQGB18F.LV2 (8.76 ~ 11.80) 10.0	Pass	05.28.2020
HbA1c	AAQGB18F	AAQGB18F.LV1 (4.31 ~ 6.1) 10.0	AAQGB18F.LV2 (8.76 ~ 11.80) 5.0	Fail	05.28.2020

AFIAS-10 Shut Down

1) Sleep mode

- ① Press briefly the power button  to enter sleep mode.

2) Reboot

- ① Press and hold the button  Operator: admin for at least 3 seconds and touch 'Reboot' in the active window.
- ② The equipment reboots again.

3) Shut down

- ① -1  Press and hold the power button for at least 4 seconds to enter the 'Power Off' window.
- ① -2 Press and hold the button  Operator: admin button or at least 3 seconds and touch 'Power off' in the active window.
- ② When you touch the 'OK' button, the instrument will proceed with the shutting down procedure.
- ③ Turn off the main power on/off switch on the back of the instrument.


Frequently Asked Questions

- The below issues can appear in the pop-up window or on the screen.

Issue	Cause	Solution
No ID Chip 	If the inserted cartridge and the lot on the ID chip do not match.	Please register the ID Chip and proceed with the test again.
	If the inserted cartridge and the lot on the ID chip do not match.	Insert the cartridges of Lot that match the ID Chip and proceed with the test.
Unknown 	Barcode printed on the cartridge: If the barcode is obscured by other material or the printing is damaged and the equipment does not recognise the cartridge barcode.	1. Remove any other material that obscures the barcode and reinsert the cartridge. 2. Recognition failure: Press the center of the cartridge slot picture on the screen to select manually the correct ID chip information (ID chip and cartridge LOT must match) 3. If the same error code continues to occur, test using a new cartridge.
	Barcode module failure.	Contact Boditech Med or the local representative for technical support.
Tip required 	If the test was carried out without a tip.	Check whether the tip has been inserted, add tip(s) to the cartridge or Extra Tip station, and touch the center of the slot in question to re-examine.
Sample Empty 	Sample detection error.	Check the sample and proceed with the examination again.
	If sample amount is small.	
Sample Error 	If the sample tube lid is not removed when inserted in the sample rack.	Check the tube lid removal status and proceed with the inspection again.
	Sample Error. (Bubble Detection or Blood Clot)	If there is too much foam in the sample tube, remove the bubbles in the tube and re-examine the sample. If there are clots, take the sample again and re-examine.
WB (B) S/P (B)	Sampling attempted 2-3 times before processing (possible problem in the sample, such as bubble).	Please consider examining again and report the final result.
Invalid TF / Invalid 25~ 34	In case of abnormal graph, the instrument shows one of these errors.	Please test with the new cartridge. If the same error occurs again, call Boditech Med or the local representative for technical support.
Extra Tip 1~X	If the Extra Tip quantity is insufficient to process specific items. (T3, Testosterone, Total Anti-infliximab ...)	After replenishing the Extra Tip station, proceed with the inspection.
Abnormal Temperature, Test banned	If the test is started with the temperature inside the device not reaching the appropriate temperature.	When the temperature inside the instrument reaches the appropriate temperature, the test will resume.

Trouble shooting

Issue	Cause	Solution
Tip is not removed	After tip removing step, the device still recognises that there is a tip.	After the technician cleans the tip sensor contaminated area, proceed with the test again.
	The instrument attempted to remove the Tip, but it could not actually remove it.	Clean the contaminated area of the tip removal slide.
If the sample rack goes in but comes out immediately	Sample rack type not recognised.	Attach a new sample rack barcode and proceed with the inspection again.
If the sample rack does not enter automatically	If the sample rack itself is not recognised.	Contact Boditech Med or the local representative for technical support.
If the equipment does not recognize the sample tube barcode	The sample barcode location is outside the barcode recognition range.	If the barcode is outside the recognition area, reattach it to the correct position and proceed with the inspection again. -13 mm x 75 mm tube : attach it at least 1 cm above the bottom of the tube -13 mm x 100 mm Tube : 1cm above the bottom of the tube
	If the sample barcode is damaged and cannot be read from the equipment.	Re-print the barcode, attach it, and proceed with the testing it again.
	If the equipment fails to rotate the sample barcode.	1. If the sample barcodes are pasted in layers, remove the existing barcodes, attach a single barcodes, and proceed with the re-examination. 2. If the sample barcode gets detached, secure it and conduct again the inspection. If the problem rises when the sample barcode is not loosen or detached, contact Boditech Med or the local representative for technical support.
If the cartridge is not present but the instrument recognizes it as present	Possible cartridge sensor malfunctions due to a light source entering the equipment.	Please change the position of the equipment in a direction where no light source enters directly.
If the cartridge itself is not recognized (inactivated)	Cartridge sensor failure.	Contact Boditech Med or the local representative for technical support.
Initialization fail	If the equipment moving parts are not in the correct position.	Please proceed with the equipment reboot.

Your life, Our inspiration
www.boditech.co.kr

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Lot	AACOQK16_LV1	Exp.	2021-11
HbA1c	HbA1c	Mean	5.30
		Low	4.51
		High	6.10
Register		Close	

Lot	AACOQK16_LV2	Exp.	2021-11
HbA1c	HbA1c	Mean	10.30
		Low	8.76
		High	11.85
Register		Close	

Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACOQK16_LV2	2021.11.30		9.9	8.42	11.39
HbA1c	AACOQK16_LV1	2021.11.30		5.3	4.51	6.1

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


European Representative: **Obelis s.a.**
Bd. Général Wahis, 53, 1030 Brussels, Belgium
(SRN: BE-AR-000000106)

Product: **AFIAS-10**

Basic UDI-DI	880613301031BM
Reference number (Catalog number)	FPRR038

Intended purpose: AFIAS-10 is an analyzer intended for use in conjunction with fluorescence immunoassay (FIA) kits for quantitative, semi-quantitative and qualitative measurements of various analytes.
For *in vitro* diagnostic use only.

Classification: Class A (Rule 5)
- According to the Annex VIII of the REGULATION (EU) 2017/746

Conformity Assessment Route: According to the Annex II and Annex III of the REGULATION (EU) 2017/746

We, Boditech Med Inc., herewith declare under our sole responsibility that the above-mentioned product is in conformity with the following European Union harmonisation legislation.

- REGULATION (EU) 2017/746 OF EUROPEAN PARLIAMENT AND OF THE COUNCIL OF 5 April 2017 on *in vitro* diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU
- DIRECTIVE 2011/65/EU OF EUROPEAN PARLIAMENT AND OF THE COUNCIL OF 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- COMMISSION DELEGATED DIRECTIVE (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances

All supporting documentation is retained under the premises of the manufacturer.

Standards applied: EN ISO 13485:2016, EN ISO 14971:2019, EN 13612:2002,
EN ISO 20417:2021, EN ISO 15223-1:2021, ISO 18113-1:2022,
ISO 18113-3:2022, ISO/TR 20416:2020, EN IEC 61326-1:2021,
EN IEC 61326-2-6:2021, EN 61010-1:2010/A1:2019,
EN IEC 61010-2-101:2022, EN 62366-1:2015, EN 62304:2006,
EN 62321-4:2014/A1:2017, EN 62321-5:2014, EN 62321-6:2015,
EN 62321-7-1:2015, EN 62321-7-2:2017, EN 62321-8:2017

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November 26, 2024

To whom it may concern,

We hereby authorize INTERLUX Ltd. With its principal place of business at **Aviečių str. 16, LT-08418, Vilnius, Lithuania** to act as the exclusive distributor for the below-stated products manufactured by Boditech Med Inc. with its principal place of business at 43 Geodudanji 1-gil, Dongnae-myeon, Chuncheon, Gangwon-do, Republic of Korea, Postal Code 24398 (the "Company").

INTERLUX Ltd. Is authorized to register, sell, market, provide full guarantees and technical service in Lithuania, Latvia, and Estonia market.

- **ichroma™ analyzers and reagents**
- **AFIAS analyzers and reagents**

This Letter of Authorization shall become effective as of January 1, 2025 and expires on December 31, 2025. The terms of this Letter of Authorization may be renewed for an additional year period, subject to the terms of any definitive agreement which shall be entered into by and between the Distributor and the Company. No change or amendment to this Letter of Authorization shall be valid unless agreed in writing by the Company.

Yours faithfully,

Representative Director
Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gangwon-do, South Korea
Boditech Med Inc.

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(MB „Beikeris“, įm .k. 304539005)

Certificate

This is to certify that _____ at Interlux has successfully completed the basic training requirement in accordance with the provisions of Boditech Med, Inc. to use and service AFIAS series.

This entitles the recipient to perform limited services for the above products.

Certificate No. BT240328-3

Date of Issuance : March 28th, 2024

Authorized by _____

Boditech Med, INC.

Boditech Med Inc.
43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gangwon-do, South Korea 200-883
T +82-33-243-1400 F +82-33-243-9373

Certificate

This is to certify that _____ at Interlux has successfully completed the basic training requirement in accordance with the provisions of Boditech Med, Inc. to use and service AFIAS series.

This entitles the recipient to perform limited services for the above products.

Certificate No. BT240328-2

Date of Issuance : March 28th, 2024

Authorized by _____

Boditech Med, INC.

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, South Korea 200.883
T: +82-33-243-1400 F: +82-33-243-9373

Certificate

This is to certify that at Interlux has successfully completed
the basic training requirement in accordance with the provisions
of Boditech Med, Inc. to use and service
AFIAS series.

This entitles the recipient to perform limited services
for the above products.

Certificate No. BT240328-1

Date of Issuance : March 28th, 2024

Authorized by

Boditech Med, INC.

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gangwon-do, South Korea 200-883
T: +82-33-243-1400 F: +82-33-243-9373

Boditech AFIAS Quality Control Guidelines

9 p.

5 p.

Instrument Quality Control

1. To perform system check, you need to use a system check cartridge and system check ID chip.
2. AFIAS Analyzers measure the intensity of fluorescent signal from the system check cartridge and compare the measured intensity with preset value to confirm the system's reliability.
3. System checks are needed in the following cases:
 - When an analyzer is first installed or moved to another location.
 - If there are frequent errors during tests.
 - When there are some doubts about the test results.
4. Boditech Med Inc. recommends performing the system check test 1 time/month to secure analyzer's functional integrity.

Cartridge Calibration

1. A master calibration curve for lot specific cartridges is saved in an ID chip included in the cartridge box. Therefore, there is no need to perform a separate calibration process.
2. Regular calibration is not required/recommended for AFIAS as the system is designed for on-site (point of care), emergency and small-volume testing.

Cartridge Quality Control

1. Each cartridge contains a separate internal control band to prevent and correct technical errors.
2. Boditech Med Inc. recommends performing quality control to ensure the validity and reliability of test results.
3. To perform quality control testing, you need to use control materials provided by Boditech Med Inc. on demand or 3rd party control materials (not applicable to every test).
4. The quality control testing is needed as below.
 - When a new cartridge lot is used.
 - When a test result is considered abnormal or differs considerably from the reference values.
 - When a test result is not consistent with symptoms.
5. Boditech Med Inc. suggests recommends quality control testing 1 time/month to secure analyzer's functional integrity.
6. Quality Control testing should be in accordance with any relevant local regulations.

Name:

Signed:

Date: October 11, 2023

Position: Quality Management Representative

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STATEMENT

To whom it may concern,

We **Boditech Med Inc.**, a corporation duly organized and existing under the laws of the Republic of Korea, having its principal office at 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si Gangwon-do, 24398 Korea, hereby certify that:

The **AFIAS-10** analyzer incorporates a feature allowing users to fill in the User ID field during login through either manual input or barcode scanning. Patient ID field can also be populated by either manual input or barcode scanning.

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1. MSH segment structure

9 p.

- This segment is the HL7 MSH message segment, indicating the Message Header.
- Segment Field

- MSH-1: Field Separator (ST)
- MSH-2: Encoding Characters (ST)
- MSH-3: Sending Application (HD) optional
- MSH-4: Sending Facility (HD) optional
- MSH-5: Receiving Application (HD) optional
- MSH-6: Receiving Facility (HD) optional
- MSH-7: Date/Time of Message (DTM)
- MSH-8: Security (ST) optional
- MSH-9: Message Type (MSG)
- MSH-10: Message Control ID (ST)
- MSH-11: Processing ID (PT)
- MSH-12: Version ID (VID)
- MSH-13: Sequence Number (NM) optional
- MSH-14: Continuation Pointer (ST) optional
- MSH-15: Accept Acknowledgment Type (ID) optional
- MSH-16: Application Acknowledgment Type (ID) optional
- MSH-17: Country Code (ID) optional
- MSH-18: Character Set (ID) optional repeating
- MSH-19: Principal Language Of Message (CWE) optional
- MSH-20: Alternate Character Set Handling Scheme (ID) optional
- MSH-21: Message Profile Identifier (EI) optional repeating
- MSH-22: Sending Responsible Organization (XON) optional
- MSH-23: Receiving Responsible Organization (XON) optional
- MSH-24: Sending Network Address (HD) optional
- MSH-25: Receiving Network Address (HD) optional

2. PID segment structure

- This segment is the HL7 PID segment, indicating the Patient Identification
- Segment Field

- PID-1: Set ID - PID (SI) optional
- PID-2: Patient ID (CX) optional
- PID-3: Patient Identifier List (CX) repeating
- PID-4: Alternate Patient ID - PID (CX) optional repeating
- PID-5: Patient Name (XPN) repeating
- PID-6: Mother's Maiden Name (XPN) optional repeating
- PID-7: Date/Time of Birth (DTM) optional
- PID-8: Administrative Sex (IS) optional
- PID-9: Patient Alias (XPN) optional repeating
- PID-10: Race (CWE) optional repeating

3. OBR segment structure

- This segment is the HL7 OBR message segment, indicating the Observation Request
- Segment Field

- OBR-1: Set ID - OBR (SI) optional
- OBR-2: Placer Order Number (EI) optional
- OBR-3: Filler Order Number (EI) optional
- OBR-4: Universal Service Identifier (CWE)
- OBR-5: Priority (ID) optional
- OBR-6: Requested Date/Time (DTM) optional
- OBR-7: Observation Date/Time # (DTM) optional
- OBR-8: Observation End Date/Time # (DTM) optional
- OBR-9: Collection Volume * (CQ) optional
- OBR-10: Collector Identifier * (XCN) optional repeating
- OBR-11: Specimen Action Code * (ID) optional

Document ID	Name	LIS Communication Protocol Specifications	Rev.	0.2
R-A1R-010	Product	For POCcelerator	Page	2 / 6

- OBR-12: Danger Code (CWE) optional

4. ORC segment structure

-. This segment is the HL7 ORC message segment, indicating the Common Order

-. Segment Field

- ORC-1: Order Control (ID)
- ORC-2: Placer Order Number (EI) optional
- ORC-3: Filler Order Number (EI) optional
- ORC-4: Placer Group Number (EI) optional
- ORC-5: Order Status (ID) optional
- ORC-6: Response Flag (ID) optional
- ORC-7: Quantity/Timing (TQ) optional repeating
- ORC-8: Parent (EIP) optional
- ORC-9: Date/Time of Transaction (DTM) optional
- ORC-10: Entered By (XCN) optional repeating
- ORC-11: Verified By (XCN) optional repeating
- ORC-12: Ordering Provider (XCN) optional repeating
- ORC-13: Enterer's Location (PL) optional
- ORC-14: Call Back Phone Number (XTN) optional repeating
- ORC-15: Order Effective Date/Time (DTM) optional
- ORC-16: Order Control Code Reason (CWE) optional
- ORC-17: Entering Organization (CWE) optional
- ORC-18: Entering Device (CWE) optional
- ORC-19: Action By (XCN) optional repeating

5. SPM segment structure

-. This segment is the HL7 SPM message segment, indicating the specimen

-. Segment Field

- SPM-1: Set ID - SPM (SI) optional
- SPM-2: Specimen ID (EIP) optional
- SPM-3: Specimen Parent IDs (EIP) optional repeating
- SPM-4: Specimen Type (CWE)
- SPM-5: Specimen Type Modifier (CWE) optional repeating
- SPM-6: Specimen Additives (CWE) optional repeating
- SPM-7: Specimen Collection Method (CWE) optional
- SPM-8: Specimen Source Site (CWE) optional
- SPM-9: Specimen Source Site Modifier (CWE) optional repeating
- SPM-10: Specimen Collection Site (CWE) optional
- SPM-11: Specimen Role (CWE) optional repeating
- SPM-12: Specimen Collection Amount (CQ) optional
- SPM-13: Grouped Specimen Count (NM) optional
- SPM-14: Specimen Description (ST) optional repeating
- SPM-15: Specimen Handling Code (CWE) optional repeating
- SPM-16: Specimen Risk Code (CWE) optional repeating
- SPM-17: Specimen Collection Date/Time (DR) optional
- SPM-18: Specimen Received Date/Time (DTM) optional
- SPM-19: Specimen Expiration Date/Time (DTM) optional

6. OBX segment structure

-. This segment is the HL7 OBX message segment, indicating the Observation/Result

-. Segment Field.

- OBX-1: Set ID - OBX (SI) optional
- OBX-2: Value Type (ID) optional
- OBX-3: Observation Identifier (CWE)
- OBX-4: Observation Sub-ID (ST) optional
- OBX-5: Observation Value (Varies) optional repeating

Document ID	Name	LIS Communication Protocol Specifications	Rev.	0.2
R-A1R-010	Product	For POCcelerator	Page	3 / 6

- OBX-6: Units (CWE) optional
- OBX-7: References Range (ST) optional
- OBX-8: Abnormal Flags (IS) optional repeating
- OBX-9: Probability (NM) optional
- OBX-10: Nature of Abnormal Test (ID) optional repeating
- OBX-11: Observation Result Status (ID)

7. Message example

-. Message formulation from Analyzer to server

-. <SB>: 0x0B

-. <EB>: 0x1C

-. <CR>: 0x0D

```

<SB>
MSH|^~\&|Protocol|ModelName|DeviceNum||Date||OUL^R24^OUL_R24|1|T|2.6<CR>
PID||PatientID1|PatientID2|PatientID3|||Age|Gender<CR>
OBR||TestName|TestType|SampleType|||Date|||UserID<CR>
ORC|OK|||||||DeviceNum|SlotNum<CR>
SPM|1|CartridgeLotName|||||||CartridgeExpiryDate<CR>
SPM|2|ControlLotName|ControlName|||||||ControlExpiryDate <CR>
OBX|1|TX|Result1 TestName||Result1 Value|Unit1|Result1 Flag|Result1 Type||R<CR>
OBX|2|TX|Result2 TestName||Result2 Value|Unit2|Result2 Flag|Result2 Type||R<CR>
OBX|3|TX|Result3 TestName||Result3 Value|Unit3|Result3 Flag|Result3 Type||R<CR>
OBX|4|TX|Result4 TestName||Result4 Value|Unit4|Result4 Flag|Result4 Type||R<CR>
OBX|5|TX|Result5 TestName||Result5 Value|Unit5|Result5 Flag|Result5 Type||R<CR>
<EB><CR>

```

-. The meaning of segment

Segment example	Each parameter segment meaning
<SB>	Start Block character
MSH ^~\& Protocol ModelName DeviceNum Date OUL^R24^OUL_R24 1 T 2.6<CR>	Message header
PID PatientID1 PatientID2 PatientID3 Age Gender<CR>	Patient information
OBR TestName TestType SampleType Date UserID<CR>	Test information
ORC OK DeviceNum SlotNum<CR>	Device information
SPM 1 CartridgeLotName CartridgeExpiryDate<CR>	Cartridge information
SPM 2 ControlLotName ControlName ControlExpiryDate <CR>	Control information
OBX 1 TX Result1 TestName Result1 Value Unit1 Result1 Flag Result1 Type R<CR>	Test result 1
OBX 2 TX Result2 TestName Result2 Value Unit2 Result2 Flag Result2 Type R<CR>	Test result 2
OBX 3 TX Result3 TestName Result3 Value Unit3 Result3 Flag Result3 Type R<CR>	Test result 3
OBX 4 TX Result4 TestName Result4 Value Unit4 Result4 Flag Result4 Type R<CR>	Test result 4
OBX 5 TX Result5 TestName Result5 Value Unit5 Result5 Flag Result5 Type R<CR>	Test result 5
<EB><CR>	End Block character, Carriage return

Document ID	Name	LIS Communication Protocol Specifications	Rev.	0.2
R-A1R-010	Product	For POCcelerator	Page	4 / 6

-. each segment location of the results in HL7 field

Number	Parameter	Description	Valid value	Segment location
1	Protocol	Protocol version	Digit	MSH-3
2	ModelName	Device model name	alphanumeric	MSH-4
3	DeviceNum	Device serial number	alphanumeric	MSH-5, ORC-18
4	Date	Test date (YYYYMMDDHHMMSS)	alphanumeric	MSH-7, OBR-7
5	UserID	Operator ID	alphanumeric	OBR-10
6	Patient ID1	User input Patient ID	alphanumeric	PID-2
7	Patient ID2	Specimen barcode ID	alphanumeric	PID-3
8	Patient ID3	Default Patient ID	alphanumeric	PID-4
9	Age	Patient age	Digit	PID-7
10	Gender	Patient gender	Male / Female	PID-8
11	TestName	Test name Ex) CRP, PCT	alphanumeric	OBR-2
12	TestType	0: General mode 1: QC mode 2: QC calibration mode	Digit	OBR-3
13	SampleType	0: Whole blood 1: Serum/Plasma 2: Urine 3: Feces 4: Swab 5: Etc	Digit	OBR-4
14	CartridgeLotName	Test cartridge lot name	alphanumeric	1 st SPM-2
15	CartridgeExpiry Date	Test cartridge expiry date (YYYYMMDD)	alphanumeric	1 st SPM-19
16	ControlLotName	Control lot name	alphanumeric	2 nd SPM-2
17	ControlName	Control name	alphanumeric	2 nd SPM-4
18	ControlExpiryDate	Control expiry date (YYYYMM)	alphanumeric	2 nd SPM-19
19	Result1TestName	Test name of result 1 Ex) CRP, PCT	alphanumeric	1 st OBX-3
20	Result2TestName	Test name of result 2 Ex) CRP, PCT	alphanumeric	2 nd OBX-3
21	Result3TestName	Test name of result 3 Ex) CRP, PCT	alphanumeric	3 rd OBX-3

Document ID	Name	LIS Communication Protocol Specifications	Rev.	0.2
R-A1R-010	Product	For POCcelerator	Page	5 / 6

22	Result4TestName	Test name of result 4 Ex) CRP, PCT	alphanumeric	4 th OBX-3
23	Result5TestName	Test name of result 5 Ex) CRP, PCT	alphanumeric	5 th OBX-3
24	Result1Value	Quantitative result of Result 1	alphanumeric	1 st OBX-5
25	Result1Flag	Qualitative result of Result 1 0: Positive 1: Negative 2: Indeterminate	Digit	1 st OBX-7
26	Result1Type	Result type of Result 1 0: Quantitative 1: Qualitative 2: Semi quantitative	Digit	1 st OBX-8
27	Result2Value	Quantitative result of Result 2	alphanumeric	2 nd OBX-5
28	Result2Flag	Qualitative result of Result 2 0: Positive 1: Negative 2: Indeterminate	Digit	2 nd OBX-7
29	Result2Type	Result type of Result 2 0: Quantitative 1: Qualitative 2: Semi quantitative	Digit	2 nd OBX-8
30	Result3Value	Quantitative result of Result 3	alphanumeric	3 rd OBX-5
31	Result3Flag	Qualitative result of Result 3 0: Positive 1: Negative 2: Indeterminate	Digit	3 rd OBX-7
32	Result3Type	Result type of Result 3 0: Quantitative 1: Qualitative 2: Semi quantitative	Digit	3 rd OBX-8
33	Result4Value	Quantitative result of Result 4	alphanumeric	4 th OBX-5
34	Result4Flag	Qualitative result of Result 4 0: Positive 1: Negative 2: Indeterminate	Digit	4 th OBX-7
35	Result4Type	Result type of Result 4 0: Quantitative 1: Qualitative 2: Semi quantitative	Digit	4 th OBX-8
36	Result5Value	Quantitative result of Result 5	alphanumeric	5 th OBX-5
37	Result5Flag	Qualitative result of Result 5 0: Positive 1: Negative 2: Indeterminate	Digit	5 th OBX-7
38	Result5Type	Result type of Result 5 0: Quantitative 1: Qualitative 2: Semi quantitative	Digit	5 th OBX-8

Document ID	Name	LIS Communication Protocol Specifications	Rev.	0.2
R-A1R-010	Product	For POCcelerator	Page	6 / 6

39	Unit1	Unit of result 1	alphanumeric	1 st OBX-6
40	Unit2	Unit of result 2	alphanumeric	2 nd OBX-6
41	Unit3	Unit of result 3	alphanumeric	3 rd OBX-6
42	Unit4	Unit of result 4	alphanumeric	4 th OBX-6
43	Unit5	Unit of result 5	alphanumeric	5 th OBX-6
44	SlotNum	Device slot number	alphanumeric	ORC-19

- . Message formulation from Server to analyzer (Only 3 byte)
- . <SB>: 0x0B
- . <EB>: 0x1C
- . <CR>: 0x0D

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Product: Boditech CRP Control
Cat. No. : CFPO-100

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European Representative: OBELIS S.A.
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech Tn-I Plus Control
Cat. No. : CFPO-212

Classification: Others
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002, EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019, E

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European Representative: OBELIS S.A
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: AFIAS MxA/CRP
Cat. No. : SMFP-102

Classification: Others
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, March 20, 2022

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European Representative: OBELIS S.A.
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech MxA/CRP Control
Cat. No. : CFPO-382

Classification: Others
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

Boditech Med Inc. www.boditech.co.kr

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

RA-DOC-III-199 (Rev. 00)

AFIAS CRP

INTENDED USE

AFIAS CRP is a fluorescence immunoassay (FIA) for the quantitative determination of C-Reactive Protein (CRP) in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of autoimmune diseases and infectious processes, such as rheumatoid arthritis.

For *in vitro* diagnostic use only.

INTRODUCTION

C-reactive protein (CRP) is a protein found in blood; the level of which rises in response to inflammation. CRP is first acute-phase protein to be described and is an exquisitely sensitive systemic marker of inflammation and tissue damage. The serum CRP level may rise from a normal level of < 5 mg/L to 500 mg/L during the body's general, nonspecific response to infections and other acute inflammatory conditions. Measurement of CRP concentration has been widely used as a clinical tool for monitoring the status of inflammation, effectiveness of treatment of various infections and autoimmune diseases such as rheumatoid arthritis. And high sensitive CRP (hsCRP) has been suggested that along with serum cholesterol, may be marker of diagnosis in cardiovascular diseases (CVD). And hsCRP is emerging as the strongest and most independent predictive risk factor for CVD.

PRINCIPLE

The test uses a sandwich immunodetection method.

The detector antibodies in detector bind to antigens in the sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized-antibodies on a test strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by the instrument for AFIAS tests to show CRP concentration in the sample.

COMPONENTS

AFIAS CRP consists of 'cartridges'.

- Each sealed aluminum pouch contains two cartridges.
- Each sealed cartridge packaged in an aluminum pouch has three components including a cartridge part, a detector part and detector diluent part.
- The cartridge part contains the membrane called a test strip which has mouse monoclonal anti-human CRP at the test line, human C-reactive protein at the antigen line, and streptavidin at the control line.
- The detector part has a granule containing mouse monoclonal anti-human CRP-fluorescence conjugator, mouse monoclonal anti-human CRP, Biotin-BSA fluorescence conjugator, bovine serum albumin (BSA) and sucrose as a stabilizer and sodium azide as a preservative in phosphate buffered saline (PBS).
- The detector diluent part contains antifoam and tween 20 as a detergent, sodium azide as a preservative in phosphate buffered saline (PBS).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instructions for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges, pipette tips, and C-tips. It should be used for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use the cartridge, if pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- If test components and/or sample are stored in refrigerator, then allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges and pipette tips, and C-tip should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- The cartridge contains sodium azide (NaN₃), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.
- Biotin interference was not observed in **AFIAS CRP** when biotin concentration in the sample was up to 3,500 ng/mL. In reference, if biotin is taken at a level of 300 mg/day, the biotin concentration in the body is about 1,200 ng/mL. However, if a patient has been taking biotin at dosage of more than 300mg a day, it is recommended to collect blood again 24 hours after discontinuation of biotin intake.
- AFIAS CRP** will provide accurate and reliable results subject to the below conditions.
 - AFIAS CRP** should be used only in conjunction with the instrument for AFIAS tests.
 - Have to use recommended anticoagulant.

Recommended anticoagulant
K ₂ EDTA, K ₃ EDTA, Sodium heparin, Lithium heparin, Sodium citrate
- C-tip should be used when the following conditions are met.**
 - C-tip provided with the kit is recommended to obtain correct test result.
 - Whole blood should be immediately tested after collection.
 - Do not perform a test with C-tip on General Mode. It might cause an erroneous result.
 - Excess whole blood around the C-tip should be wiped off.
 - In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.
 - AFIAS cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.
 - While collecting blood, be careful not to create air bubbles in the C-tip.

LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 - 30 °C	20 months	Unopened
		1 month	Resealed

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Following items can be purchased separately from **AFIAS CRP**. Please contact our sales division for more information.

Instrument for AFIAS tests	REF FPRR019
- AFIAS-1	REF FPRR040
- AFIAS-3	REF FPRR020
- AFIAS-6	REF FPRR038
- AFIAS-10	REF CFPO-100
Boditech CRP Control	REF CFPO-112
Boditech CRP Calibrator	

SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS CRP** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The samples (serum, plasma) should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.

- The samples (whole blood, serum, plasma) may be stored for a week at 2-8 °C prior to being tested. If testing will be delayed more than a week, samples should be frozen at below -20 °C.
- The samples stored frozen at below -20 °C for 3 months showed no performance difference.
- However, the whole blood sample should not be kept in a freezer in any case.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.
- Collection of whole blood sample using C-tip
 - Hold the C-tip horizontally and touch the surface of the blood with the tip of the C-tip.
 - Capillary action will automatically draw the blood sample to C-tip and stop.
 - Wipe off any excess blood around the tip.
 - Double-check if whole blood is filled accurately in the C-tip and AFIAS reader is ready for a test on the 'C-tip mode'.

TEST SETUP

- Check the components of the **AFIAS CRP** as described below.: Cartridges, pipette tips, C-tips, an ID chip, a spare cartridge zipper bag and an instructions for use.
 - Ensure that the lot number of the cartridge matches that of the ID chip.
 - If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
 - Turn on the instrument for AFIAS tests.
 - Empty the tip box.
 - Insert the ID chip into the 'ID chip port'.
- ※ **Please refer to the instrument for AFIAS tests operation manual for complete information and operating instructions.**

TEST PROCEDURE

► AFIAS-1, AFIAS-3, AFIAS-6,

General mode

- Insert the cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Select the 'General mode' in the instrument for AFIAS tests.
- Take 100 µL of the sample (whole blood/serum/plasma/control) using a pipette and dispense it into the sample well of the cartridge.
- Tap the 'Start' button on the screen.
- The test result will be displayed on the screen after 3 minutes.

C-tip Mode

- Insert a cartridge into the cartridge holder.
- Take 10 µL of whole blood using a C-tip.
- Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- Select the 'C-tip mode' in the instrument for AFIAS tests.
- Tap the 'Start' button on the screen.
- The test result will be displayed on the screen after 3 minutes.

► AFIAS-10

Normal mode

- Insert a cartridge into the cartridge holder.

- 2) Insert a tip into the tip hole of the cartridge.
- 3) Tap the 'Load' button of the bay that holds the cartridge with the tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 4) Insert the sample tube into the tube rack.
- 5) Insert the tube rack into the loading part of the sampling station.
- 6) Tap the 'Start' button on the screen.
- 7) The test result will be displayed on the screen after 3 minutes.

Emergency mode (General tip)

- 1) The test procedure is same with the 'Normal mode 1) – 3)'.
- 2) Convert the 'Emergency mode' in AFIAS-10.
- 3) Select the tip type (General tip) on the screen.
- 4) Select the sample type (whole blood/serum/plasma) on the screen.
- 5) Take 100 µL of the sample using a pipette and dispense it into the sample well of the cartridge.
- 6) Tap the 'Start' button on the screen.
- 7) The test result will be displayed on the screen after 3 minutes.

Emergency mode (C-tip)

- 1) Insert a cartridge into the cartridge holder.
- 2) Take 10 µL of whole blood using a C-tip.
- 3) Insert the C-tip with sample into the tip hole of the cartridge.
- 4) Tap the 'Load' button of the bay that holds the cartridge with a tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 5) Convert the 'Emergency mode' in AFIAS-10.
- 6) Select the tip type (C-tip) on the screen.
- 7) Tap the 'Start' button on the screen.
- 8) The test result will be displayed on the screen after 3 minutes.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays CRP concentration of the test sample in terms of mg/L.
- Cut-off : 10 mg/L
- Working range : 0.5-200 mg/L
- Effect of Hematocrit

The CRP Whole Blood of the instrument for AFIAS tests is calibrated to read the CRP serum concentration of a blood sample with a hematocrit of 40 %. If the actual hematocrit value deviates from 40 %, the result should be corrected by multiplying with the respective factor in the table: deviates from 40 %, the result should be corrected by multiplying with the respective factor in the table:

Hct (%)	Factor	Hct (%)	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6
43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

Reference range, HCT:

- Woman: 35-44 %
- Men: 39-48 %

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **AFIAS CRP**. For more information regarding obtaining the control materials, contact **Boditech Med Inc.'s Sales Division for assistance**.
(Please refer to the instructions for use of control material.)

PERFORMANCE CHARACTERISTICS

- **Analytical sensitivity**
 - Limit of Blank (LoB) 0.19 mg/L
 - Limit of Detection (LoD) 0.32 mg/L
 - Limit of Quantitation (LoQ) 0.50 mg/L

- **Analytical Specificity**
 - **Cross-reactivity**
Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS CRP** test results did not show any significant cross-reactivity with these biomolecules.

Cross reactants	Concentration
IL-6	50 µg/mL
PCT	50 µg/mL
Serum amyloid P	500 µg/mL
Serum amyloid A	500 µg/mL
Ferritin	500 µg/mL

- **Interference**
Interferents listed in the following table were added to the test sample at the concentration mentioned below. **AFIAS CRP** test results did not show any significant interference with these materials.

Interferents	Concentration
Ascorbic acid	350 µmol/L
Bilirubin (conjugated)	475 µmol/L
Albumin	60 g/L
Glucose	1,000 mg/dL
Triglyceride mixture (Triglyceride, total)	1,500 mg/dL
Hemoglobin	10 g/L
Biotin	3,500 ng/ml

- **Precision**
 - **Single site study**
3 Lots of **AFIAS CRP** were tested for 20 days. Each standard material was tested 2 times per day. For each test, each material was duplicated.
Repeatability (within-run precision)
Within-laboratory precision (Total precision)
Lot to lot precision
 - **Between persons**
Three different persons tested one lot of **AFIAS CRP**, ten times at each concentration of the control standard.
 - **Between sites**
One person tested **AFIAS CRP** at three different sites, ten times at each concentration of the control standard.
 - **Between readers**
One person tested **AFIAS CRP** with three different readers, ten times at each concentration of the control standard.

Expected Value [mg/L]	Repeatability			Within-laboratory precision		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	1.97	0.16	8.0	2.02	0.17	8.5
10.0	9.96	0.63	6.3	10.01	0.61	6.1
150.0	148.76	13.95	9.4	149.33	13.79	9.2
Expected Value [mg/L]	Lot to lot precision			Between site		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	2.00	0.18	8.9	2.04	0.16	7.6
10.0	9.99	0.6	6.0	10.02	0.71	7.1
150.0	149.04	14.39	9.65	152.32	9.09	6.0
Expected Value [mg/L]	Between person			Between reader		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	2.02	0.13	6.5	2.00	0.18	9.0
10.0	10.03	0.81	8.1	9.93	0.70	7.0
150.0	155.47	11.92	7.7	152.49	11.79	7.7


■ **Accuracy**

The accuracy was confirmed by testing 3 different lots, three times with international standard material (ERM-DA474/IFCC, human serum).

	Lo1		Lot 2		Lot 3	
	Measured Value [mg/L]	Bias (%)	Measured Value [mg/L]	Bias (%)	Measured Value [mg/L]	Bias (%)
Run 1	41.82	1.5	39.88	-3.2	42.68	3.6
Run 2	40.09	-2.7	41.42	0.5	40.52	-1.7
Run 3	42.06	2.1	41.96	1.8	40.95	-0.6

■ **Comparability**

CRP concentration of 100 clinical samples were quantified independently with **AFIAS CRP (AFIAS-6)** and **comparator A** as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and correlation coefficient (R). The regression equation and correlation coefficient are as follows.

REFERENCES

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. Mol Immunol 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. Circulation 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C-reactive protein and Lipid Screening. Clin. Chem. 2001; 47:28-30
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. Clin. Chem.

2001; 47(3): 403-411.

5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. Circulation 1999; 99:855-860
6. Taubes G. Does inflammation cut to the heart of the matter Science 2002; 296:242-245
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. N Engl J Med 2000;342(12): 836-843
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. Clin Chem 1999; 45:1676-1678
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of -care testing. Clin Chim Acta 2005; 356:172-177
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J. Lab. Clin Med 1976;87:120-128
11. Kindmark CO. The concentration of C-reactive protein in sera from healthy individuals. Scand J Clin Lab Invest 1972;29:407-411

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
LITERATŪROS ŠARŠAS

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. Mol Immunol 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. Circulation 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C reactive protein and Lipid Screening. Clin. Chem. 2001; 47:28-30
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. Clin. Chem. 2001; 47(3): 403-411.
5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. Circulation 1999; 99:855-860
6. Taubes G. Does inflammation cut to the heart of the matter Science 2002; 296:242-245
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. N Engl J Med 2000;342(12): 836-843
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. Clin Chem 1999; 45:1676-1678
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of-care testing. Clin Chim Acta 2005; 356:172-177
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J. Lab. Clin Med 1976;87:120-128
11. Kindmark CO. The concentration of C-reactive protein in sera from healthy individuals. Scand J Clin Lab Invest 1972;29:407-411

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
Boditech CRP Control

INTENDED USE

Boditech CRP Control is intended for the quality control of CRP assay kits manufactured/provided by Boditech Med Inc.
For *in vitro* diagnostic use only.

COMPONENTS

Boditech CRP Control consists of 'Boditech CRP Control Level 1', 'Boditech CRP Control Level 2', 'Instructions for use' and 'Control value & Barcode sheet'.

- Boditech CRP Control is provided in liquid form.
- The control contains CRP antigen and horse serum.
- The control materials are contained in vials, and the vials are further packaged in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech CRP Control should not be used past the expiration date.
- Boditech CRP Control is solely designed for the quality control of CRP cartridges manufactured/provided by Boditech Med Inc.
- Human-derived materials in Boditech CRP Control were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be non-reactive. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human-derived materials and clinical samples should be handled as though they are capable of transmitting infectious diseases and should be disposed of as hazardous wastes.
- All waste materials should be disposed of in accordance with the requirements of your local waste management authorities.
- Boditech CRP Control contains sodium azide (NaN₃), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.

STORAGE AND STABILITY

- Storage and stability condition of Boditech CRP Control.

	Unopened	Opened (After reconstitution)
Temperature	2 ~ 8 °C	2 ~ 8 °C
Expiration date	Until expiration date on the label	28 days

- Close the opened vial tightly after use.
- After use, any residual substance should not be returned to the original vial.
- Bacterial contamination of opened Boditech CRP Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be used.

TEST PROCEDURE

Boditech CRP Control is supplied in liquid form.

1. Thoroughly mix the contents of the vial before each use by gently inverting for several times.
2. When removed from the refrigerator, allow the vial to stand for about 30 minutes before using the control.

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In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Sales.**

MATERIALS SUPPLIED

REF CFPO-100

Boditech CRP Control Box (2 vials)

- | | |
|---------------------------------|---|
| - Boditech CRP Control Level 1 | 1 |
| - Boditech CRP Control Level 2 | 1 |
| - Instructions for use | 1 |
| - Control value & Barcode sheet | 1 |

INTERPRETATION OF THE RESULT

The test result of the 'Boditech CRP Control' should be consistent with the expected result of control value sheet.

If the test results fall outside the expected result, check the following for potential sources of error, and retest after resolving these matters. If the error still persists, contact the **Boditech Med Inc.'s Technical Sales.**

- ※ Potential sources of error
 - Errors in a testing process
 - Incorrect storage condition of Boditech CRP Control
 - Use of expired or contaminated Boditech CRP Control.
 - Faulty Boditech's CRP assay kits.
 - Faulty Boditech's instruments.

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	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

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
5. Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. Clin Chem 1986; 32:B64-70.

6. Goldstein DE, Little RR, England JD, Wiedmeyer H-M, McKenzie E. Methods of glycosylated hemoglobins: high performance liquid chromatography and thiobarbituric acid colorimetric methods. In: Clarke WL, Larner J, Pohl SL, eds. Methods in diabetes research, Vol. 2. New York: John Wiley, 1986:475-504.

7. Tahara Y, Shima K. The response of GHb to stepwise plasma glucose change over time in diabetic patients. Diabetes Care 1993; 16:1313-4.

8. Svendsen PA, Lauritzen T, Soegaard U, Nerup J. Glycosylated haemoglobin and steady-state mean blood glucose concentration in type 1 (insulin-dependent) diabetes. Diabetologia 1982; 23:403-5.

9. Cefalu WT, Wang ZQ, Bell-Farrow A, Kiger FD, Izlar C. Glycohemoglobin measured by automated affinity HPLC correlates with both short-term and long-term antecedent glycemia. Clin Chem 1994; 40:1317-21.

10. Singer DE, Coley CM, Samet JH, Nathan DM. Tests of glycemia in diabetes mellitus. Their use in establishing a diagnosis and in treatment. Ann Intern Med 1989; 110:125-37.

11. Molnar GD. Clinical evaluation of metabolic control in diabetes. Diabetes 1978; 27:216-25.

12. UK Prospective Diabetes Study. Reduction in HbA1c with basal insulin supplement, sulfonylurea or biguanide therapy in maturity-onset diabetes. Diabetes 1985; 34:793-8.

13. Baker JR, Johnson RN, Scott DJ. Serum fructosamine concentrations in patients with type II (non-insulin-dependent) diabetes mellitus during changes in management. BMJ (Clin Res Ed) 1984; 288:1484-6.

14. Tahara Y, Shima K. Kinetics of HbA1c, glycated albumin, and fructosamine and analysis of their weight functions against preceding plasma glucose level. Diabetes Care 1995; 18:440-7.

15. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point of care testing. Clin Chem 1999; 45:1676-1678.

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INTENDED USE

AFIAS HbA1c is a fluorescence Immunoassay (FIA) for the quantitative determination of HbA1c (Hemoglobin A1c) in human whole blood. It is useful as an aid in management and monitoring of the long-term glycemic status in patients with diabetes mellitus.

For *in vitro* diagnostic use only.

INTRODUCTION

Glycated protein is formed post-translationally through the slow, nonenzymatic reaction between glucose and amino groups on proteins. HbA1c is a clinically useful index of mean glycemia during the preceding 120 days, the average life span of erythrocytes. Carefully controlled studies have documented a close relationship between the concentrations of HbA1c and mean glycemia. HbA1c is considered as a more reliable parameter in monitoring glycemia over the glycemic reading with the conventional glucometer.

PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibody in buffer binds to antigen in sample, forming antigen-antibody complexes, and migrates onto nitrocellulose matrix to be captured by the other immobilized-antibody on test strip.

The more antigen in sample forms the more antigen-antibody complex and leads to stronger intensity of fluorescence signal on detector antibody, which is Instrument for the instrument for AFIAS tests displays the content of glycated hemoglobin in terms of percent of the total hemoglobin in blood.

COMPONENTS

AFIAS HbA1c consists of 'Cartridge', 'C- tip', 'ID chip' and 'Instruction for use'.

- Each sealed aluminum pouch contains two cartridges.
- Each cartridge packaged in an aluminum pouch has three components, detection buffer part, hemolysis buffer part and cartridge part.
- Cartridge part contains a test strip, the membrane which has anti human HbA1c at the test line, while rabbit IgG at the control line.
- Detection buffer part contains anti human HbA1c-fluorescence conjugate, anti rabbit IgG-fluorescence conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide in phosphate buffered saline (PBS) as a preservative.
- Hemolysis buffer part contains non-ionic detergent and sodium azide in phosphate buffered saline (PBS) as a preservative.
- C-tip (Capillary tip) is a useful tool for point of care testing which requires a small volume of capillary blood from

fingertip, heel site (in infants) or ear-lobe.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instruction for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect of test result(s).
- Do not reuse cartridges. A cartridge should be used for testing one sample only.
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge, if pouch is damaged.
- For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used and should be recollected.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used C-tips, and cartridges should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- AFIAS HbA1c** will provide accurate and reliable results subject to the below conditions.
 - AFIAS HbA1c** should be used only in conjunction with the instrument for AFIAS tests.
 - Have to use recommended anticoagulant sample.

Sample type	Recommended anticoagulant
Whole blood	Na ₂ EDTA, K ₂ EDTA, K ₃ EDTA, Sodium citrate, Lithium heparin

- C-tip should be used when the following conditions are met.**
 - C-tip provided with the kit is recommended to obtain correct test result.
 - Capillary blood should be immediately tested after collection.
 - Do not leave C-tip after collection of blood, test immediately.
 - Excess capillary blood around the C-tip should be wiped off.
 - In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.
 - AFIAS cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.

LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample

components to the capture/detector antibodies.

- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

STORAGE AND STABILITY

- The cartridge is stable for 20 months (while sealed in an aluminum pouch) if stored at 2-8 °C.
- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.

MATERIALS SUPPLIED

REF SMFP-28

Components of **AFIAS HbA1c**

- Cartridge Box Contains
 - Cartridge 24
 - C- Tip (Zipper bag, 10 μL) 24
 - ID Chip 1
 - Instruction for Use 1
 - Spare Cartridge Zipper bag 1

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS HbA1c**

Please contact our sales division for more information.

- AFIAS-1** **REF** FPRR019
- AFIAS-6** **REF** FPRR020
- Boditech HbA1c Control** **REF** CFPO-96
- Boditech HbA1c Calibrator** **REF** CFPO-108

SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS HbA1c** is human whole blood.

- It is recommended to test the sample within 12 hours after collection.
- Samples may be stored for up to a week at 2-8 °C prior to being tested.
- If testing will be delayed more than a week, samples should be frozen at -70 °C or below. Samples stored frozen at -70 °C or below for 3 months showed no performance difference.
- Once the sample was frozen, it should be used one time only for test, because repeated freezing and thawing can

result in the change of test values.

- Collection of capillary blood sample using C-tip (Fingertip blood)
 - Clean the area with a pre-injection swab.
 - Pierce with a sterile lancet.
 - Wipe away first drop of blood.
 - Gently massage the surrounding area towards a C-tip for a second drop.
 - Hold a C-tip horizontally and touch the tip of C-tip to the blood drop.
 - Capillary action will automatically draw the blood sample to C-tip and stop
 - Wipe off any excess blood around the tip.
 - Double-check if capillary blood is fully filled in the C-tip and AFIAS reader is ready for a test on the 'C-tip mode'.

TEST SETUP

- Check the components of the **AFIAS HbA1c** as described below.: Cartridge, C-tip, ID chip and instruction for use.
- Keep the sealed cartridge (if stored in refrigerator) at room temperature for at least 30 minutes just prior to the test. Place the cartridge on a clean, dust-free and flat surface.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the "ID chip port".
- Please refer to the instrument for AFIAS tests 'Operation Manual' for complete information and operating instructions.

TEST PROCEDURE

- Select "Test Mode" in the instrument for AFIAS tests.
 - General Mode : Tube blood (with anticoagulant)
 - C-tip Mode : Fingertips blood (non-anticoagulant blood)
 - You must use C-tip only to collect test samples in both test mode (General / C-tip mode).**
 - Do not use a general pipette tip to collect test samples.**
- Insert the cartridge into the cartridge holder.
- Take 10 μL of whole blood (tube blood or fingertip blood) with a C-tip. (Please refer to "SAMPLE COLLECTION AND PROCESSING")
- Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test result will be displayed on the screen after 10 minutes.

※ Note: Refer to the instrument for AFIAS tests Operation Manual to select sample type.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculate the test result automatically and displays HbA1c concentration of the test sample in terms of percentage (% , NGSP) and mmol/mol (IFCC).
- The reference value
 - NGSP (%): 4.5-6.5 %
 - IFCC (mmol/mol): 26-48 mmol/mol

- The working range
 - NGSP (%): 4-15 %
 - IFCC (mmol/mol): 20.2-140.4 mmol/mol
 - eAG (mg/dL): 68.1-383.8 mg/dL

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with **AFIAS HbA1c**. For more information regarding obtaining the control materials, contact **Boditech Med Inc.'s Sales Division** for assistance. (Please refer to the instruction for use of control material.)

PERFORMANCE CHARACTERISTICS

- **Analytical sensitivity**
 - Limit of Blank (LoB) 2.26%
 - Limit of Detection (LoD) 2.39%
 - Limit of Quantification (LoQ) 3.0 %
- **Analytical Specificity**
 - Cross-reactivity
There was no significant cross-reactivity from these materials with the **AFIAS HbA1c** test measurements.

cross-reactivity material	Standard material conc. (%)		
	5.3	8.1	13.5
Recovery (%)			
HbA0 (2 mg/mL)	98.9	102.0	101.0
HbA1a&A1b mixture (5 mg/mL)	100.9	101.9	100.1
Acetylated hemoglobin (100 mg/mL)	100.7	100.2	100.5
Carbamylated hemoglobin (100 mg/mL)	100.2	101.2	100.6
Glycated h-Albumin (2.2 mg/mL)	100.1	100.6	102.0
HbA1d (100 mg/mL)	101.6	101.3	100.5
Acetylaldehyde hemoglobin (100 mg/mL)	99.8	101.4	100.7


- Cross-reactivity
There was no significant interference from these materials with the **AFIAS HbA1c** test measurements.

Interference material	Standard material conc. (%)		
	5.3	8.1	13.5
Recovery (%)			
Non-interference	100.0	99.8	99.5
Acetaminophen (1.5 mmol/L)	100.3	100.6	99.6
L-Ascorbic acid (5 mg/dL)	100.1	100.6	100.8
Bilirubin (20 mg/mL)	100.3	100.3	100.0
D-Glucose (500 mg/dL)	99.8	100.6	102.0
Intralipid (1 mg/mL)	101.0	99.7	100.5
Triglyceride (3000 mg/mL)	100.4	99.6	99.9

Urea (100 mmol/L)		100.3	100.0	100.8				
■ Precision								
- Between lot One person tested three different lots of AFIAS HbA1c , five times at each concentration of the control standard.								
- Between reader Three different persons tested AFIAS HbA1c at three different reader; five times at each concentration of the control standard.								
- Between day One person tested AFIAS HbA1c for five days; five times at each concentration of the control standard.								
- Between site One person tested AFIAS HbA1c at three different sites; five times at each concentration of the control standard.								
HbA1c (%)	Between lot		Between reader		Between day		Between site	
	Avg.	CV(%)	Avg.	CV(%)	Avg.	CV(%)	Avg.	CV(%)
5.3	5.21	2.7	5.25	3.6	5.31	2.5	5.27	2.3
8.1	8.15	3.1	8.12	2.7	8.20	2.9	8.23	3.0
13.5	13.48	2.8	13.38	3.3	13.60	2.8	13.57	2.8

■ Accuracy The accuracy was confirmed each different concentration.				
HbA1c (%)	Between lot		Between reader	
	Avg.	Bias(%)	Avg.	Bias(%)
5.3	5.21	-0.02	5.25	-0.01
8.1	8.15	0.01	8.12	0
13.5	13.48	0	13.38	-0.01
HbA1c (%)	Between day		Between site	
	Avg.	Bias(%)	Avg.	Bias(%)
5.3	5.31	0	5.27	-0.01
8.1	8.2	0.01	8.23	0.02
13.5	13.6	0.01	13.57	0.01

- **Comparability**
HbA1c concentrations of 165 clinical samples were quantified independently with **AFIAS HbA1c** and Tosoh G7 HPLC (Tosoh Bioscience Inc. Japan) as per prescribed test procedures. Test results were compared and their comparability was investigated with linear regression and coefficient of correlation (R). Linear regression and coefficient of correlation between the two tests were $Y = 0.997X + 0.0309$ and $R = 0.9982$ respectively.

REFERENCES

1. Goldstein DE, Little RR, Lorenz RA, Malone JI, Nathan D, Peterson CM. Tests of glycemia in diabetes. Diabetes Care 1995; 18:896-909.
2. Bunn HF. Nonenzymatic glycosylation of protein: relevance to diabetes. Am J Med 1981; 70:325-30.
3. Jovanovic L, Peterson CM. The clinical utility of glycosylated hemoglobin. Am J Med 1981; 70:331-8.
4. Nathan DM, Singer DE, Hurxthal K, Goodson JD. The clinical information value of the glycosylated hemoglobin assay. N Engl J Med 1984; 310:341-6.
5. Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. Clin Chem 1986; 32:B64-70.
6. Goldstein DE, Little RR, England JD, Wiedmeyer H-M, McKenzie E. Methods of glycosylated hemoglobins: high performance liquid chromatography and thiobarbituric acid colorimetric methods. In: Clarke WL, Larner J, Pohl SL, eds. Methods in diabetes research, Vol. 2. New York: John Wiley, 1986:475-504.
7. Tahara Y, Shima K. The response of GHb to stepwise plasma glucose change over time in diabetic patients. Diabetes Care 1993; 16:1313-4.
8. Svendsen PA, Lauritzen T, Soegaard U, Nerup J. Glycosylated haemoglobin and steady-state mean blood glucose concentration in type 1 (insulin-dependent) diabetes. Diabetologia 1982; 23:403-5.
9. Cefalu WT, Wang ZQ, Bell-Farrow A, Kiger FD, Izlar C. Glycohemoglobin measured by automated affinity HPLC correlates with both short-term and long-term antecedent glycemia. Clin Chem 1994; 40:1317-21.
10. Singer DE, Coley CM, Samet JH, Nathan DM. Tests of glycemia in diabetes mellitus. Their use in establishing a diagnosis and in treatment. Ann Intern Med 1989; 110:125-37.
11. Molnar GD. Clinical evaluation of metabolic control in diabetes. Diabetes 1978; 27:216-25.
12. UK Prospective Diabetes Study. Reduction in HbA1c with basal insulin supplement, sulfonylurea or biguanide therapy in maturity-onset diabetes. Diabetes 1985; 34:793-8.
13. Baker JR, Johnson RN, Scott DJ. Serum fructosamine concentrations in patients with type II (non-insulin-dependent) diabetes mellitus during changes in management. BMJ (Clin Res Ed) 1984; 288:1484-6.
14. Tahara Y, Shima K. Kinetics of HbA1c, glycated albumin, and fructosamine and analysis of their weight functions against preceding plasma glucose level. Diabetes Care 1995; 18:440-7.
15. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point of care testing. Clin Chem 1999; 45:1676-1678.

Note: Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance; please contact:
Boditech Med Inc.'s Technical Services
Tel: +82 33 243-1400
E-mail: sales@boditech.co.kr

Boditech Med Incorporated
43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gang-won-do, 24398
Republic of Korea
Tel: +(82) -33-243-1400
Fax: +(82) -33-243-9373
www.boditech.co.kr

Obelis s.a
Bd. Général Wahis 53,
1030 Brussels, BELGIUM
Tel: +(32) -2-732-59-54
Fax: +(32) -2-732-60-03
E-Mail: mail@obelis.net

Boditech HbA1c Control

INTENDED USE

Boditech HbA1c Control is intended for *in vitro* diagnostic use in the quality control of HbA1c Assay Kit. **For *in vitro* diagnostic use only.**

INTRODUCTION

The use of Boditech HbA1c Control may be considered as an objective assessment of the precision of HbA1c Assay Kits and is an integral part of Good Laboratory Practices. Boditech HbA1c Control is provided in lyophilized form.

COMPONENTS

Boditech HbA1c Control consists of 'Boditech HbA1c Control level 1', 'Boditech HbA1c Control level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The control contains HbA1c stock and Horse serum
- Each control vial packed in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech HbA1c Control should not be used past the expiration date.
- Boditech HbA1c Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and HbA1c Assay Kits.
- Human source materials from which Boditech HbA1c Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

STORAGE AND STABILITY

- Storage and stability condition of Boditech HbA1c Control.

	Unopened	Opened (after reconstitution)	
Temperature	+2 to +8 °C	+2 to +8 °C	-20 to -80 °C
Expiration date	Until expiration date on the label.	7 days	20 days

- Close the opened Boditech HbA1c Control bottle tightly after use.
- Once the Boditech HbA1c Control was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech HbA1c Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

INSTRUCTIONS FOR USE

Boditech HbA1c Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 0.5 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently.
Avoid formation of foam. Do not shake.

Please refer to package inserts of the test cartridges for detailed test procedure.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

In the event of damage to the package, contact the
Boditech Med Inc.'s Technical Services.

MATERIALS SUPPLIED

REF CFPO-96

Boditech HbA1c Control Box (2 vials)	
Boditech HbA1c Control level 1 (0.5 mL)	1
Boditech HbA1c Control level 2 (0.5 mL)	1
Instruction For Use	1
Control value & Barcode Sheet	1

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

For Technical Assistance

Boditech Med Inc.'s Technical Services at

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AFIAS
Tn-I Plus

INTENDED USE

AFIAS Tn-I Plus is a fluorescence Immunoassay (FIA) for the quantitative determination of cardiac troponin-I (Tn-I) in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of acute myocardial infarction (AMI).

For *in vitro* diagnostic use only.

INTRODUCTION

Cardiac troponins are currently the most sensitive and specific biochemical markers of myocardial necrosis. There are three types of troponin in heart muscle fibers. Those are troponin-C, -I, and -T. Together they contribute to make cardiac muscle fibers contract. The clinical measurement of serum Tn-I has become an important tool in the diagnosis of acute myocardial infarction. Serum Tn-I is a more reliable than creatine kinase as a prognostic marker in people with ischemic chest pain. National and international scientific organizations have suggested the use of troponins, Tn-I and Tn-T, when implementing new diagnostic strategies in patients with acute coronary syndrome.

PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized antibodies on the strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS tests to show Tn-I concentration in the sample.

COMPONENTS

AFIAS Tn-I Plus consists of a ‘cartridges’.

- Each cartridge packaged in an aluminum pouch has three components including cartridge part, detector part and diluent part.
- Cartridge part contains the membrane called a test strip, which has streptavidin at the test line, and chicken IgY at the control line.
- Detector part contains anti Tn-I-fluorescence conjugate, anti Tn-I-biotin conjugate , anti chicken IgY-fluorescence conjugate, bovine serum albumin (BSA), sucrose, bromophenol blue, and sodium azide as a preservative.
- Diluent part contains tween 20, CA-630,and sodium azide as a preservative.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow the instructions and procedures described in this ‘Instruction for use’.
- Use fresh samples only and avoid direct sunlight.
- It is possible to use frozen samples. Please refer to “SAMPLE COLLECTION AND PROCESSING.”

- Lot numbers of all the test components (cartridge and ID chip) must match each other.
 - Do not interchange test components between different lots or use test components after the expiration date, either of which might yield incorrect test result(s).
 - The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge if the pouch found damaged or has already been opened.
 - Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local requirement. Sample with severe hemolysis and/or hyperlipidemia must not be used.
 - Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
 - The instrument for AFIAS tests may generate slight vibration during use.
 - Used cartridges and pipette tips should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
 - An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
 - No Biotin interference was observed in **AFIAS Tn-I Plus** when biotin concentration in the sample was below 5 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
 - AFIAS Tn-I Plus** will provide accurate and reliable results subject to the below conditions.
 - AFIAS Tn-I Plus** should be used only in conjunction with the instrument for AFIAS tests.
 - Have to use recommended anticoagulant sample.
- | Recommended anticoagulant | | | |
|---|--|--|--|
| Sodium heparin, Lithium heparin, Sodium citrate | | | |

STORAGE AND STABILITY			
Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 - 8 °C	20 months	Unopened
		1 month	Resealed
■ Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.			

LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause the false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors,

degradation of the test components/reagents or presence of interfering substances in the test samples.

- Any clinical diagnosis based on the test result should be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

MATERIALS SUPPLIED	
REF SMFP-35	
Components of AFIAS Tn-I Plus	
Cartridge Box Contains	
- Cartridge	24
- Pipette tip (Zipper bag)	24
- ID chip	1
- Instruction for use	1
- Spare cartridge zipperbag	1

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Following items can be purchased separately from **AFIAS Tn-I Plus**.

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SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS Tn-I Plus** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood. If longer storage is required, e.g. if the test could not be performed within 24 hours, serum or plasma should be immediately frozen below -20 °C. The freezing storage of sample up to 3 months does not affect the quality of results.
- However, the whole blood sample should not be kept in a freezer in any case.
- Once the sample was frozen, it should be thawed one time and only for test, because repeated freezing and thawing can result in the changed test values.

TEST SETUP

- Check the components of the **AFIAS Tn-I Plus** as described below: Cartridges, Pipette tips, ID chip, Cartridge spare zipper bag and Instruction for use
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests
- Empty the tip box.
- Insert the ID chip into the “ID chip port”.

(Please refer to the instrument for AFIAS tests ‘Operation Manual’ for complete information and operating instructions.

TEST PROCEDURE

- Take 100 µL of sample with a pipette and dispense it into the sample well in the cartridge.
 - Insert the cartridge into the cartridge holder
 - Insert a tip into the tip hole of the cartridge.
 - Tap the ‘START’ icon on the screen.
 - The test result will be displayed on the screen after 12 minutes.

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INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays Tn-I concentration of the test sample in terms of ng/mL.
- The working range of the AFIAS Tn-I Plus is 0.01-15.00 ng/mL.
- Expected Values**
 - In studies performed with the **AFIAS Tn-I Plus** assay involving 70 healthy volunteers in Korea, the upper reference limit (99th percentile) for Tn-I was 0.04 ng/mL. The lowest concentration with a CV less than or equal to 10 % with the **AFIAS Tn-I Plus** assay was 0.04 ng/mL.
 - Due to the release kinetics of Tn-I, a result below the decision limit within the first hours of the onset of symptoms does not rule out myocardial infarction with certainty. If myocardial infarction is still suspected, repeat the test at appropriate intervals.
 - A cut-off of 0.3 ng/mL Tn-I is recommended for diagnosis of AMI, as this yields optimal performance of 91 % of sensitivity and 92.1 % of specificity. However, laboratories should establish their own diagnostic cut-off concentration based on the clinical practice at their perspective institutions.

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with AFIAS Tn-I Plus. For more information regarding obtaining the control materials, contact Boditech Med Inc.’s Sales Division for assistance. (Please refer to the instruction for use of control material.)

PERFORMANCE CHARACTERISTICS	
Analytical sensitivity	
Limit of Blank (LoB):	0.004 ng/mL
Limit of Detection (LoD):	0.01 ng/mL
Limit of Quantitation (LoQ):	0.03 ng/mL

■ Analytical specificity

- Cross-reactivity

Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS Tn-I Plus** test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactivity material	Conc. (ng/mL)
CK-MB	60
NT-proBNP	1,000
Myoglobin	1,000
D-Dimer	1,000

- Interference

Interference materials such as below the ones in the table were added to the test sample(s) the same as the below concentrations. **AFIAS Tn-I Plus** test results did not show any significant interference with these materials except for EDTA.

Interference material	Conc.
Bilirubin	350 µmol/L
Cholesterol	13 mmol/L
D-glucose	1,000 mg/dL
Hemoglobin	2 g/L
L-Ascorbic acid	350 µmol/L
Triglyceride mixture	500 mg/dL
Heparin	3000 U/L
Sodium citrate	2 mg/mL
EDTA	3.4 µmol/L
Biotin	5 ng/mL

■ Precision

- Intra-assay (Between lot)

One person tested three different lots of **AFIAS Tn-I Plus**, ten times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.94	0.95	0.92
7.50	7.65	7.32	7.51
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.5
0.94	0.93	0.057	6.1
7.50	7.50	0.444	5.9

- Inter-assay (Between person)

Three different persons tested **AFIAS Tn-I Plus**; five times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.96	0.94
7.50	7.75	7.50	7.48
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.4
0.94	0.95	0.056	5.9
7.50	7.58	0.400	5.3



■ Accuracy

The accuracy was confirmed by 3 different lots testing ten times each different concentration.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.91	0.93
7.50	7.56	7.35	7.37
Conc. (ng/mL)	Mean	CV (%)	Recovery(%)
0.23	0.23	5.8	101.0
0.94	0.93	5.3	99.0
7.50	7.43	5.7	99.0

■ Comparability

Using Access2 (Beckman Coulter) as a comparison reference device for **AFIAS Tn-I Plus**, 100 clinical samples were independently tested for its Tn-I concentration following each instrument's procedure. Results of both the test methods were analyzed and their comparability was investigated with linear regression and coefficient of correlation (R).

REFERENCES

1. Mauro Panteghini, Franca Pacani, Kiang-Teck J.Yeo, Fred S. Apple, Robert H. Christenson. Francesco Dati, Johannes Mair, Jan Ravkilde, and Alan H.B. We. Evaluation of Imprecision for Cardiac Troponin Assays at Low-Range Concentrations. 2004;50:2:327-332.
2. Alan McNeil, PhD, FRACP, FRCPA. The Trouble with Troponin. Heart, Lung and Circulation 2007;16:S13-S16.
3. David M. Bunk and Micahel J. Welch. Characterization of a New Certified Reference Material for Human Cardiac Troponin I. Clinical Chemistry 2002;52:2:212-219
4. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, Katus H. It's time for a change to a troponin standard. Circulation 2000;102:1216 –1220.
5. Jillan R. Tate, David Heathcote, Gus Koerbin, Gary Thean,

David Andriske, Jone Bonar, Janice Gill. The harmonization of cardiac troponin I measurement is independent of sample time collection but is dependent on the source of calibrator. Clinica Chimica Acta 324:2002:13-23.

6. Ohman EM, Armstrong PW, Christenson RH, et al. Cardiac troponin T levels for risk stratification in acute myocardial ischemia. N Engl J Med 1996;335:1333– 41.
7. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996;335:1342 –9.

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For technical assistance; please contact:

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Boditech Tn-I Plus Control

INTENDED USE

Boditech Tn-I Plus Control is intended for *in vitro* diagnostic use in the quality control of Tn-I Assay Kit. **For *in vitro* diagnostic use only.**

SUMMARY AND PRINCIPLE

The use of Boditech Tn-I Plus Control may be considered as an objective assessment of the precision of Tn-I Assay Kits and is an integral part of Good Laboratory Practices. Boditech Tn-I Plus Control is provided in lyophilized form.

COMPONENTS

Boditech Tn-I Plus Control consists of 'Boditech Tn-I Plus Control level 1', 'Boditech Tn-I Plus Control level 2', 'Instruction for Use' and 'Control value & Barcode Sheet'.

- The control contains Human Cardiac Troponin I/C Complex Antigen and horse serum
- Each control vial packed in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech Tn-I Plus Control should not be used past the expiration date.
- Boditech Tn-I Plus Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and Tn-I Assay Kits.
- Human source materials from which Boditech Tn-I Plus Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

STORAGE AND STABILITY

- Storage and stability condition of Boditech Tn-I Plus Control.

Unopened	
Temperature	Expiration date
+2 to +8 °C	Until expiration date on the label.

Opened (After reconstitution)	
Temperature	Expiration date
+2 to +8 °C	1 day
-20 to -80 °C	7 days

- Close the opened Boditech Tn-I Plus Control bottle tightly after use.
- Once the Boditech Tn-I Plus Control was frozen, it should be used **ONE TIME ONLY** for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should **NOT BE RETURNED** to the original vial.

- Bacterial contamination of reconstituted Boditech Tn-I Plus Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

INSTRUCTIONS FOR USE

Boditech Tn-I Plus Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 1 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently. Avoid formation of foam. Do not shake.

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In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services.**

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INTENDED USE

AFIAS MxA/CRP is a fluorescence immunoassay (FIA) for the quantitative determination of Myxovirus resistance protein A (MxA) and C-Reactive Protein (CRP) in human whole blood. This test aids in identification of viral and/or bacterial infection in patients who have symptoms of acute respiratory infection within 7 days.
For *in vitro* diagnostic use only.

INTRODUCTION

Viral and bacterial respiratory infections represent a major source of morbidity, mortality, and healthcare costs. Approximately 80% of all antimicrobials are prescribed in primary care, and up to 80% of these are for respiratory tract indications⁸. But sensitive and specific diagnostic tools to aid in the diagnosis of ARIs (particularly in differentiating bacterial and viral infections) in primary and urgent healthcare settings have been lacking⁹.

Mx proteins are large GTPases and belong to a group of IFN-induced GTPases involved in the control of intracellular pathogens⁵. In humans, two Mx homologs (MxA and MxB, also called MX1 and MX2, respectively) mediate antiviral activity against a broad range of viruses include Covid-196. Elevated levels of MxA protein could be an indicator of endogenous interferon production mediated by an unknown viral activation and so, the MxA protein levels could be used as a general marker of viral infection⁷.

CRP is one of the cytokine-induced acute-phase proteins¹ whose blood levels rise during a general, unspecific response to infections and non-infectious inflammatory processes². CRP tests provide information for the diagnosis, therapy, and monitoring of inflammatory diseases.¹⁰ During infectious or inflammatory disease states, CRP levels rise rapidly within the first 6 to 8 hours and peak at levels of up to 350–400 mg/L after 48 hours. Measurement of CRP concentration has been widely used as a clinical tool for monitoring the status of inflammation, effectiveness of treatment of various infections and autoimmune diseases such as rheumatoid arthritis.

MxA/CRP test will aid in the differentiation of viral and bacterial acute febrile respiratory infections.

PRINCIPLE

The test uses a sandwich immunodetection method.
The detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes, and migrate onto nitrocellulose matrix to be captured by the other immobilized-antibodies on a test strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by the instrument for AFIAS tests to show MxA and CRP concentration in the sample.

COMPONENTS

- AFIAS MxA/CRP** consists of a ‘cartridges’.
- Each sealed aluminum pouch contains two cartridges.
 - Each cartridge packaged in an aluminum pouch has three components including a cartridge part, a detector part and a diluent part.

- The cartridge part contains the membrane called a test strip which has anti-MxA and anti-CRP at the test line, CRP Ag at the antigen line, and chicken IgY at the control line.
- The detector part has a granule containing the anti-MxA-fluorescence complex, anti-CRP-fluorescence conjugate, anti-chicken IgY-fluorescence conjugate, and sodium azide as a preservative in phosphate buffered saline (PBS).
- The diluent part contains sodium azide as a preservative in phosphate buffered saline (PBS).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this ‘Instructions for use’.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges. A cartridge should be used for processing of one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use a cartridge, if the pouch is damaged or has already been opened.
- If test components and/or sample are stored in refrigerator, then allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges, C-tips and pipette tips should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- The cartridge contains sodium azide (NaN₃), and it may cause certain health issues like convulsions, low blood pressure, low heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.
- No biotin interference was observed in **AFIAS MxA/CRP** when biotin concentration in the sample was below 3,500 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
- **AFIAS MxA/CRP** will provide accurate and reliable results subject to the below conditions.

- AFIAS MxA/CRP should be used only in conjunction with instrument for AFIAS tests.
- Have to use recommended anticoagulant.
Recommended anticoagulant
K ₂ EDTA, K ₃ EDTA, Sodium heparin, Lithium heparin, Sodium citrate

- **C-tip should be used when the following conditions are met.**
 - C-tip provided with the kit is recommended to obtain correct test result.
 - Whole blood should be immediately tested after collection.
 - Do not perform a test with C-tip on General Mode. It might cause an erroneous result.
 - Excess whole blood around the C-tip should be wiped off.
 - In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.
 - AFIAS cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.
 - While collecting blood, be careful not to create air bubbles in the C-tip.

LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigens to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigens with time and/or temperature may also cause false negative result as it makes antigens unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

✖ The following cases may affect the measurement results.

- Receiving interferon therapy (e.g. MS, HIV, HBV, HCV) in the last 30 days.
- Immunocompromised state (e.g. HIV) or taking immunosuppressive or chemotherapeutic medications in the last 30 days (e.g. oral steroids, Methotrexate, Cyclosporine, Antimetabolite chemotherapy, interferon therapy).
- Taking antibiotics or antiviral therapy in the last 14 days.
- Received a live viral immunization in the last 14 days.
- Significant trauma or burns (> 5% total body surface area or full thickness (3rd°)) in the last 30 days.
- Major surgery (requiring intravenous anesthesia and/or respiratory assistance) in the last 30 days- - History of a myocardial infarction or stroke in the last 30 days.
- Taking high biotin-containing drugs within 24 hours.

STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 ~ 30°C	20 months	Unopened
		1 month	Re-sealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

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Following items can be purchased separately with **AFIAS MxA/CRP**. Please contact our sales division for more information.

▪ Instrument for AFIAS tests	
- AFIAS-1	REF FPRR019
- AFIAS-3	REF FPRR040
- AFIAS-6	REF FPRR020
- AFIAS-10	REF FPRR038
▪ Boditech MxA/CRP Control	REF CFPO-382
▪ Boditech MxA/CRP Calibrator	REF CFPO-383

SAMPLE COLLECTION AND PROCESSING

- The sample type for **AFIAS MxA/CRP** is human whole blood.
- It is recommended to test the sample within 1 hour after collection.
 - The samples (whole blood) may be stored for 12 hours at 2-8°C prior to being tested.
 - However, the whole blood sample should not be kept in a freezer in any case.
 - Collection of whole blood sample using C-tip

- ① Hold the C-tip horizontally and touch the surface of the blood with the tip of the C-tip.
- ② Capillary action will automatically draw the blood sample to C-tip and stop.
- ③ Wipe off any excess blood around the tip.
- ④ Double-check if whole blood is filled accurately in the C-tip and the instrument for AFIAS tests is ready for a test on the ‘C-tip mode’.

TEST SETUP

- Check the components of the **AFIAS MxA/CRP** as described below: Cartridges, pipette tips, C-tips, an ID chip, a spare cartridge zipper bag and an instructions for use.
- Ensure that the lot number of the cartridge matches that of the ID chip.
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the “ID chip port”.

✖ Please refer to the instrument for AFIAS tests operation manual for complete information and operating instructions.

TEST PROCEDURE

► AFIAS-1, AFIAS-3, AFIAS-6
General mode
1) Insert a cartridge into the cartridge holder.
2) Insert a tip into the tip hole of the cartridge.
3) Select the ‘General mode’ in the instrument for AFIAS tests.
4) Take 150 µL of the sample (whole blood/control) using pipette and dispense it into the sample well of the cartridge.
5) Tap the ‘Start’ button on the screen.
6) The test results will be displayed on the screen after 12 minutes.
C-tip mode
1) Insert a cartridge into the cartridge holder.

- 2) Take 10 µL of whole blood using a C-tip.
- 3) Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- 4) Select the ‘C-tip mode’ in the instrument for AFIAS tests.
- 5) Tap the ‘Start’ button on the screen.
- 6) The test result will be displayed on the screen after 12 minutes.

► AFIAS-10

Normal mode

- 1) Insert a cartridge into the cartridge holder.
- 2) Insert a tip into the tip hole of the cartridge.
- 3) Tap the “load” button of the bay that holds the cartridge with the tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 4) Insert the sample tube into the tube rack.
- 5) Insert the tube rack into the loading part of the sampling station.
- 6) Tap the “Start” button on the screen.
- 7) The test result will be displayed on the screen after 12 minutes.

Emergency mode – General tip

- 1) The test procedure is same with the ‘Normal mode 1) – 3)’.
- 2) Convert the ‘Emergency mode’ in AFIAS-10.
- 3) Select the tip type (general tip) on the screen.
- 4) Select the sample type (whole blood) on the screen.
- 5) Take 150 µL of the sample using a pipette and dispense it into the sample well of the cartridge.
- 6) Tap the ‘Start’ button on the screen.
- 7) The test result will be displayed on the screen after 12 minutes.

Emergency mode – C-tip


- 1) Insert a cartridge into the cartridge holder.
- 2) Take 10 µL of whole blood using a C-tip.
- 3) Insert the C-tip with sample into the tip hole of the cartridge.
- 4) Tap the ‘Load’ button of the bay that holds the cartridge with a tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 5) Convert the ‘Emergency mode’ in AFIAS-10.
- 6) Select the tip type (C-tip) on the screen.
- 7) Tap the ‘Start’ button on the screen.
- 8) The test result will be displayed on the screen after 12 minutes.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays MxA, CRP concentration of the test sample in terms of ng/mL. and mg/L, respectively.

Unit	MxA [ng/mL]	CRP [mg/L]
Clinical Cut-off	15.0	10.0
measuring range	10.0 – 300.0	1.0-200.0

MxA [ng/mL]	CRP [mg/L]	Interpretation
MxA ≥ 15.00	CRP < 10.00	viral infection
	10.00 ≤ CRP < 20.00	viral infection and low likelihood of bacterial infection
	20.00 ≤ CRP < 100.00	viral infection and bacterial infection likely
	CRP ≥ 100.00	bacterial infection very likely
MxA < 15.00	10.00 ≤ CRP < 20.00	likelihood of bacterial infection
	20.00 ≤ CRP < 100.00	bacterial infection likely
	CRP ≥ 100.00	bacterial infection very likely
	CRP < 10.00	non-infection

- This product is intended for screening tool only. A negative result does not preclude a viral or bacterial infection. Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.
- NICE(National Institute for Health and Care Excellence) recommend CRP test for antibiotics prescription as follows.
Pneumonia in adults: diagnosis and management (CG191)
[Presentation with lower respiratory tract infection¹⁰]
For people presenting with symptoms of lower respiratory tract infection in primary care, consider a point care C-reactive protein test if after clinical assessment a diagnosis of pneumonia has not been made and it is not clear whether antibiotics should be prescribed. Use the results of the C-reactive protein test to guide antibiotic prescribing in people without a clinical diagnosis of pneumonia as follows:
- Do not routinely offer antibiotic therapy if the C-reactive protein concentration is less than 20 mg/litre.
- Consider a delayed antibiotic prescription (a prescription for use at a later date if symptoms worsen) if the C-reactive protein concentration is between 20 mg/litre and 100 mg/litre.
- Offer antibiotic therapy if the C-reactive protein concentration is greater than 100 mg/litre.

Effect of Hematocrit

The CRP Whole Blood of the instrument for AFIAS tests is calibrated to read the CRP serum concentration of a blood sample with a hematocrit of 40%. If the actual hematocrit value deviates from 40%, the result should be corrected by multiplying with the respective factor in the table: deviates from 40%, the result should be corrected by multiplying with the respective factor in the table:

Hct (%)	Factor	Hct (%)	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6

43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

Reference range, HCT:

- Woman: 35-44%
- Men: 39-48%

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **AFIAS MxA/CRP**. For more information regarding obtaining the control materials, contact **Boditech Med Inc.’s Sales Division for assistance**. (Please refer to the instructions for use of control material.)

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

	MxA [ng/mL]	CRP [mg/L]
LOB	4.50	0.30
LOD	5.70	0.40
LOQ	10.00	1.00

Analytical specificity

Cross-reactivity

Biomolecules listed in the following table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS MxA/CRP** test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactants	Concentration
MxB	100 ng/mL
Human IRGM (Immunity-related GTPase family M protein)	100 ng/mL
Human GBP1 (guanylate-binding protein 1)	100 ng/mL
Human GBP2 (guanylate-binding protein 2)	100 ng/mL
Human Gvin1 (Interferon-induced very large GTPase 1)	100 ng/mL
Human Interferon-γ-Inducible Protein-10	100 ng/mL
Interferon α	100 pg/mL
Interferon β 1a	100 pg/mL
Interferon λ	100 pg/mL
TNF-α	100 pg/mL
IL-10	100 pg/mL
IL-6	100 pg/mL
IL-2	100 pg/mL
PCT	50 ng/mL

Interference

Interferents listed in the following table were added to the test sample at the concentration mentioned below. **AFIAS MxA/CRP** test results did not show any significant interference with these materials.

Interferents	Concentration
Ascorbic acid	350 µmol/L
Bilirubin (conjugated)	475 µmol/L
Albumin	60 g/L

Glucose	1,000 mg/dL
Triglycerides mixture	1,500 mg/dL
Hemoglobin	10 g/L
K ₂ EDTA	5.45mg/mL
K ₃ EDTA	5.45mg/mL
Na-Heparin	10.41mg/mL
Li-Heparin	10.44mg/mL
Sodium Citrate	32 mg/mL

Precision

Single-site study

Repeatability (within-run precision)

within-laboratory precision (Total precision)

Lot to lot precision

3 Lots of **AFIAS MxA/CRP** were tested for 20 days. Each standard material was tested 2 times per day. For each test, each material was duplicated.

CRP Conc. [mg/L]	Repeatability (within-run)		within-laboratory		Lot to lot precision	
	Mean [mg/L]	CV (%)	Mean [mg/L]	CV (%)	Mean [mg/L]	CV (%)
10	10.15	5.7	10.05	6.4	10.00	6.4
40	39.90	6.5	40.13	6.4	39.75	6.3
150	151.24	6.3	150.63	6.2	149.78	6.3

MxA Conc. [ng/mL]	Repeatability (within-run)		within-laboratory		Lot to lot precision	
	Mean [ng/mL]	CV (%)	Mean [ng/mL]	CV (%)	Mean [ng/mL]	CV (%)
20	20.24	5.7	20.22	5.8	19.99	6.2
50	48.56	6.1	49.68	6.6	50.39	6.6
100	99.26	6.5	99.99	6.2	99.87	6.3

Multi-site study

Reproducibility

1 Lot of **AFIAS MxA/CRP** was tested for 5 days in 3 different sites (1 person per 1 site, 1 instrument per 1 site). Each standard material was tested 1 time per and 5 replicates per day.

Multi-site study (Reproducibility)			
CRP Conc. [mg/L]	Mean [mg/L]	SD	CV (%)
10	9.95	0.57	5.7
40	39.54	2.56	6.5
150	151.54	10.41	6.9
MxA Conc. [ng/mL]	Mean [ng/mL]	SD	CV (%)
20	20.42	1.26	6.2
50	49.53	3.46	7.0
100	100.88	6.70	6.6

Accuracy

The accuracy was confirmed by testing with 3 different lots of **AFIAS MxA/CRP**. Each MxA sample was tested in 10 replicates, and recovery (%) calculated. Each CRP sample was tested in 10 replicates, and bias (%) was calculated.

Expected value CRP [mg/L]	Mean [mg/L]			3 LOTs of AVG [mg/L]	Bias (%)
	Lot 1	Lot 2	Lot 3		
41.20	38.84	40.47	41.60	40.30	-2.2%

Expected value MxA	Mean [ng/mL]			3 LOTs of AVG [ng/mL]	Recovery (%)
	Lot 1	Lot 2	Lot 3		

[ng/mL]					
20.00	20.10	19.39	19.05	19.51	97.6
100.00	97.09	99.92	95.24	97.42	97.4

▪ **Comparability**

AFIAS MxA/CRP has demonstrated the following clinical performance results.
The results were evaluated as 39 samples infected with bacteria, 25 samples infected with the virus.

Clinical Diagnosis	Sensitivity (95% CI)	Specificity (95% CI)
Viral infection	88.0 (75.3 - 100.0)	94.9 (87.9 – 100.0)
Bacterial infection	87.2 (76.7 - 97.7)	88.0 (75.3 – 100.0)

REFERENCES

1. Kushner I. The Phenomenon of the Acute Phase Response. Ann NY Acad Sci 1982, 39:389.
2. James V. Callaghan. Review Criteria for Assessment of C-Reactive Protein (CRP), High Sensitivity C-Reactive Protein (hsCRP) and Cardiac C-Reactive Protein (cCRP) Assays. 2005, CDRH.
3. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003, 111:1805-1812.
4. Nicola R. Sproston. Role of C-Reactive Protein at Sites of inflammation and infection. Immunol. 2018, 9:754.
5. Song Gao. Structure of Myxovirus Resistance Protein A Reveals Intra- and Intermolecular Domain Interactions Required for the Antiviral Function. Immunity 35, 2011, 514–525.
6. Alexej Dick. Role of Nucleotide Binding and GTPase Domain Dimerization in Dynamin-like Myxovirus Resistance Protein A for GTPase Activation and Antiviral Activity. THE JOURNAL OF BIOLOGICAL CHEMISTRY VOL. 290, NO. 20, 2015, 12779–12792.
7. Mizuho Kawamura. New Sandwich-Type Enzyme-Linked Immunosorbent Assay for Human MxA Protein in a Whole Blood Using Monoclonal Antibodies Against GTP-Binding Domain for Recognition of Viral Infection. Journal of Clinical Laboratory Analysis 26, 2012, 174–183
8. Robert Sambursky. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. European Clinical Respiratory Journal. 2015.
9. Matt Shirley. A Rapid Diagnostic Test for Differentiating Bacterial and Viral Aetiologies in Acute Respiratory Infections. Molecular Diagnosis & Therapy 23, 2019, 803–809.
10. Woodhead, M. et al. Pneumonia in adults: diagnosis and management (CG191), National Institute for Health and Care Excellence_Clinical guideline. Published: 3 December 2014, Last updated 7 July 2022. Available at <https://www.nice.org.uk/guidance/cg191>.

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
For technical assistance, please contact:

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
NUORODOS

1. Kushner I. The Phenomenon of the Acute Phase Response. Ann NY Acad Sci 1982, 39:389.
2. James V. Callaghan. Review Criteria for Assessment of C-Reactive Protein (CRP), High Sensitivity C-Reactive Protein (hsCRP) and Cardiac C-Reactive Protein (cCRP) Assays. 2005, CDRH.
3. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003, 111:1805-1812.
4. Nicola R. Sproston. Role of C-Reactive Protein at Sites of inflammation and infection. Immunol. 2018, 9:754.
5. Song Gao. Structure of Myxovirus Resistance Protein A Reveals Intra- and Intermolecular Domain Interactions Required for the Antiviral Function. Immunity 35, 2011, 514–525.
6. Alexej Dick. Role of Nucleotide Binding and GTPase Domain Dimerization in Dynamin-like Myxovirus Resistance Protein A for GTPase Activation and Antiviral Activity. THE JOURNAL OF BIOLOGICAL CHEMISTRY VOL. 290, NO. 20, 2015, 12779–12792.
7. Mizuho Kawamura. New Sandwich-Type Enzyme-Linked Immunosorbent Assay for Human MxA Protein in a Whole Blood Using Monoclonal Antibodies Against GTP-Binding Domain for Recognition of Viral Infection. Journal of Clinical Laboratory Analysis 26, 2012, 174–183.
8. Robert Sambursky. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. European Clinical Respiratory Journal. 2015.
9. Matt Shirley. A Rapid Diagnostic Test for Differentiating Bacterial and Viral Aetiologies in Acute Respiratory Infections. Molecular Diagnosis & Therapy 23, 2019, 803–809.
10. Woodhead, M. et al. Pneumonia in adults: diagnosis and management (CG191), National Institute for Health and Care Excellence Clinical guideline. Published: 3 December 2014, Last updated 7 July 2022. Available at <https://www.nice.org.uk/guidance/cg191>.

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Boditech MxA/CRP Control

INTENDED USE

Boditech MxA/CRP Control is intended for the quality control of MxA/CRP assay kit manufactured/provided by Boditech Med Inc. For *in vitro* diagnostic use only.

COMPONENTS

Boditech MxA/CRP Control consists of 'Boditech MxA/CRP Control Level 1', 'Boditech MxA/CRP Control Level 2', 'Instructions for use' and 'Control value & Barcode sheet'.

- **Boditech MxA/CRP Control** is provided in lyophilized form.
- The control contains MxA recombinant protein, Human C-reactive protein and sodium azide (NaN₃) as a preservative in horse serum.
- The control materials are contained in vials, and the vials are further packaged in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- **Boditech MxA/CRP Control** should not be used past the expiration date.
- **Boditech MxA/CRP Control** is solely designed for the quality control of MxA/CRP cartridges manufactured/provided by Boditech Med Inc.
- Human-derived materials in **Boditech MxA/CRP Control** were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be non-reactive. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human-derived materials and clinical samples should be handled as though they are capable of transmitting infectious diseases and should be disposed of as hazardous wastes.
- All waste materials should be disposed of in accordance with the requirements of your local waste management authorities.
- **Boditech MxA/CRP Control** contains sodium azide (NaN₃), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.

STORAGE AND STABILITY

- Storage and stability condition of **Boditech MxA/CRP Control**.

	Unopened	Opened (After reconstitution)	
Temperature	2 ~ 8 °C	2 ~ 8 °C	-80 ~ -20 °C
Expiration date	Until expiration date on the label.	7 days	30 days

- Close the opened vial tightly after use.
- Once the **Boditech MxA/CRP Control** was frozen, it should be used one time only for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual substance should not be returned to the original vial.
- Bacterial contamination of reconstituted **Boditech MxA/CRP Control** will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

TEST PROCEDURE

Boditech MxA/CRP Control is supplied in lyophilized form.

- 1) Carefully reconstitute each vial of lyophilized control material with exactly 1.0mL of sterilized distilled water.
- 2) Close the vial and allow it to stand for 30 minutes before use. Ensure the contents are completely dissolved by swirling the vial gently.
(To avoid the formation of foam, do not shake the vial.)

Please refer to the instructions for use of the test cartridges for detailed test procedure.

In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services**.

MATERIALS SUPPLIED

REF CFPO-382

Boditech MxA/CRP Control box (2 vials):

- | | |
|------------------------------------|---|
| - Boditech MxA/CRP Control Level 1 | 1 |
| - Boditech MxA/CRP Control Level 2 | 1 |
| - Instructions for use | 1 |
| - Control value & Barcode sheet | 1 |

INTERPRETATION OF THE RESULT

The test result of the '**Boditech MxA/CRP Control**' should be consistent with the expected result of control value sheet. If the test results fall outside the expected result, repeat the test.

- ※ Possible causes of test errors
 - Errors in a testing process
 - Incorrect storage condition of **Boditech MxA/CRP Control**
 - Use of expired or contaminated **Boditech MxA/CRP Control**
 - Faulty Boditech's MxA/CRP assay kits
 - Faulty Boditech's instruments

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

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For technical assistance, please contact:
Boditech Med Inc.'s Technical Services at
Tel: +(82) -33-243-1400
E-mail: TS@boditech.co.kr

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I. General Information

- A. Product Name/Catalogue Number
: AFIAS CRP / SMFP-2
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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[Cartridge]

	Component	Synonyms	Classification No.	Amount (%)
Housing	2-Methyl-2-propenoic acid methyl ester homopolymer	PMMA	9011-14-7	57 ~ 67
	2-Propenenitrile polymer with 1,3-butadiene and ethenylbenzene	ABS	9003-56-9	33 ~ 43

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I. General Information

- A. Product Name/Catalogue Number
: AFIAS HbA1c / SMFP-28
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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Other hazards

		Cartridge
	Test strip	Housing
		Health: 1
NEPA	No data available	Flammability: 1
		Physical hazard: 0

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[Cartridge]

	Component	Synonyms	Classification No.	Amount
Housing	2-Methyl-2-propenoic acid methyl ester homopolymer	PMMA	9011-14-7	57~67 (%)
	2-Propenenitrile polymer with 1,3-butadiene and ethenylbenzene	ABS	9003-56-9	33~43 (%)

IV. First Aid Measures

- In case of eye contact
: Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
: Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
: Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed
: If irritation or signs of toxicity occur, seek medical attention.
- * The following symptoms may occur: Irreversible eye damage.
- * Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- * This material should be considered as being potentially infectious.

V. Fire Fighting Measures

- Recommended extinguishing media
: Suitable extinguishing agents including CO₂, WATER SPRAY or regular form.
- Specific hazard from the chemical
: Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
: Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
: Chemical resisted protective equipment for fire-fighter.

MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS415
Rev. No.	07
Rev. Date	2022. 05. 09

VI. Accidental Release Measures

- Personal precautions
 - * If not danger, stop to leak
 - * Do not contact with exposed material without protective equipment
 - * Note the avoid conditions and materials
 - * Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions
 - : No data available
- Methods and materials for purification and cleaning up
 - : After absorbed exposed materials, clean up with cleaner and water

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- Engineering management
 - : Good insolation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes
- Standard of exposure	
Chemical	No data available
Biological	No data available

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X. Stability and Reactivity

	Cartridge		Diluent Buffer
	Test strip	Housing	
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C	Heat, high temperature
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk	No data available
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency conditions		

XI. Toxicological Information

	Cartridge	Diluent buffer
Information on the likely routes of exposure	No data available	No data available
Acute toxicity	No data available	No data available
Skin corrosion/Irritation	No data available	No data available
Serious eye damage/Eye irritation	No data available	No data available
Respiratory/Skin sensitization	No data available	No data available
Germ cell mutagenicity	No data available	No data available
Carcinogenicity	No data available	No data available
Reproductive toxicity	No data available	No data available
Specific target organ toxicity(single exposure)	No data available	No data available
Specific target organ toxicity(repeated exposure)	No data available	No data available
Aspiration hazard	No data available	No data available

XII. Ecological Information

	Cartridge	Diluent buffer
Toxicity	No data available	No data available
Persistence and degradability	No data available	No data available

MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS415
Rev. No.	07
Rev. Date	2022. 05. 09

Bio-accumulative potential	No data available	No data available
Mobility in soil	No data available	No data available
Other adverse effects	No data available	No data available

XIII. Disposal Considerations

	Test strip	Not available
Cartridge	Granule	Not available
	Housing	Consider the notice on Waste management act
Diluent buffer		Not available

XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous good
- Not a dangerous good under DOT, IATA, ADR, IMDG, or RID	

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I. General Information

- A. Product Name/Catalogue Number
: Boditech HbA1c Control / CFPO-96
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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- * Avoid skin friction.
- * Note the avoid conditions and materials from fire or flame.
- * Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
- * Note the avoid conditions and materials.
- * Store container tightly closed in a well-ventilated area.
- * Storage temperature: low temperature

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I. General Information

- A. Product Name/Catalogue Number
: AFIAS Tn-I Plus / SMFP-35
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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Other hazards

	Cartridge	
	Test strip	Housing
		Health: 1
NEPA	No data available	Flammability: 1
		Physical hazard: 0

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X. Stability and Reactivity

	Cartridge		Diluent Buffer
	Test strip	Housing	
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C	Heat, high temperature
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk	No data available
Hazardous decomposition products formed under fire conditions Irritation or toxic gases and fume may be occurred in fire-emergency conditions			

XI. Toxicological Information

	Cartridge	Diluent buffer
Information on the likely routes of exposure	No data available	No data available
Acute toxicity	No data available	No data available
Skin corrosion/Irritation	No data available	No data available
Serious eye damage/Eye irritation	No data available	No data available
Respiratory/Skin sensitization	No data available	No data available
Germ cell mutagenicity	No data available	No data available
Carcinogenicity	No data available	No data available
Reproductive toxicity	No data available	No data available
Specific target organ toxicity(single exposure)	No data available	No data available
Specific target organ toxicity(repeated exposure)	No data available	No data available
Aspiration hazard	No data available	No data available

XII. Ecological Information

	Cartridge	Diluent buffer
--	-----------	----------------

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	05
	Rev. Date	2022. 05. 09

Toxicity	No data available	No data available
Persistence and degradability	No data available	No data available
Bio-accumulative potential	No data available	No data available
Mobility in soil	No data available	No data available
Other adverse effects	No data available	No data available

XIII. Disposal Considerations

Cartridge	Test strip	Not available
	Granule	Not available
	Housing	Consider the notice on Waste management act
Diluent buffer		Not available

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I. General Information

- A. Product Name/Catalogue Number
: Boditech Tn-I Plus Control / CFPO-212
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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IV. First Aid Measures

- In case of eye contact
 - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
 - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
 - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed
 - : If irritation or signs of toxicity occur, seek medical attention.
- * The following symptoms may occur: Irreversible eye damage.
- * Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- * This material should be considered as being potentially infectious.

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- Engineering management
 - : Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

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XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

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I. General Information

- A. Product Name/Catalogue Number
: AFIAS MxA/CRP SMFP-102
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

II. Hazard Identification

GHS classification of the product and national or regional

	Cartridge
GHS classification	Not applicable
IMDG Code	Not applicable

GHS label elements including precautionary

	Cartridge
Symbols	No data available
Signal word	No data available
Hazard statements	No data available

Precautionary statements

	Cartridge
Prevention	No data available
Reaction	No data available
Storage	No data available
Disposal	No data available

Other hazards

	Cartridge
	Housing
	Health: 1
NEPA	Flammability: 1
	Physical hazard: 0

MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS469
Rev. No.	00
Rev. Date	2021. 11. 19

III. Composition/Information on Ingredients

[Cartridge]

	Component	Synonyms	Classification No.	Amount (%)
Housing	Poly(methyl methacrylate)	PMMA	9011-14-7	57 ~ 67
	Acrylonitrile butadiene styrene	ABS	9003-56-9	33 ~ 43

IV. First Aid Measures

- In case of eye contact
 - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
 - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
 - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring.
- If swallowed
 - : If irritation or signs of toxicity occur, seek medical attention.
- * The following symptoms may occur: Irreversible eye damage.
- * Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- * This material should be considered as being potentially infectious.

V. Fire Fighting Measures

- Recommended extinguishing media
 - : Suitable extinguishing agents including CO₂, WATER SPRAY or regular form.
- Specific hazard from the chemical
 - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
 - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
 - : Chemical resisted protective equipment for fire-fighter.

VI. Accidental Release Measures

- Personal precautions
 - * If not danger, stop to leak
 - * Do not contact with exposed material without protective equipment
 - * Note the avoid conditions and materials
 - * Do not breathe dust / fume / gas / mist / vapours / spray

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS469
	Rev. No.	00
	Rev. Date	2021. 11. 19

- Environmental precautions
- : No data available
- Methods and materials for purification and cleaning up
- : After absorbed exposed materials, clean up with cleaner and water

VII. Handling and Storage

- Precaution for safe handling
 - * Avoid skin friction.
 - * Note the avoid conditions and materials from fire or flame.
- * Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
 - * Note the avoid conditions and materials.
 - * Store container tightly closed in a well-ventilated area.
- * Storage temperature: low temperature

VIII. Exposure controls and Personal protection

- Engineering management
- : Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes

- Standard of exposure

Chemical	No data available
Biological	No data available

IX. Physical and Chemical properties

	Cartridge
	Housing
State	Natural color or pigment processed color
Odor	Little odorless
Odor threshold value	No data available
pH	Natural
Melting/	180 ~ 200°C/

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS469
	Rev. No.	00
	Rev. Date	2021. 11. 19

Freezing point	No data available
Initial boiling point or boiling range	No data available
Flash point	No data available
Evaporating rate	No data available
Lower explosion limit	No data available
Upper explosion limit	No data available
Vapor pressure	No data available
Water solubility	Not soluble in water, soluble in acetone
Density	1.11
Vapor density	No data available
Specific gravity	No data available (water=1)
N-octan/water partition coefficient	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Molecular weight	50,000 ~ 250,000

X. Stability and Reactivity

	Cartridge
	Housing
Chemical stability	Stable at high temperature and pressure
Conditions to avoid	Heat, high temperature above 200°C
Materials to avoid	Oxidizing agent: fire or explosive risk
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency

XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS469
	Rev. No.	00
	Rev. Date	2021. 11. 19

Specific target organ toxicity (single exposure)	No data available
Specific target organ toxicity (repeated exposure)	No data available
Aspiration hazard	No data available

XII. Ecological Information

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

XIII. Disposal Considerations

Cartridge	Housing	Consider the notice on Waste management act
-----------	---------	---

XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods

Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

XV. Regulatory Information

Korea Industrial Safety and Health Act	No data available
Korea Hazardous Materials Safety Control Act	No data available
Korea Toxic Chemical Control Act	No data available
Korea Wastes Control Act	No data available
Other internal and foreign acts	No data available

XVI. Other Notes

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.

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Dokumento Nr.	BT-MSDS469
Rev. Nr.	00
Rev. data	2021 m. lapkričio

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I. General Information

- A. Product Name/Catalogue Number
: Boditech MxA control / CFPO-352
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

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GHS classification	Not applicable
IMDG Code	Not applicable

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We recommend handling all chemicals with caution.

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X. Stability and Reactivity

Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature
Materials to avoid	No data available
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDSC104
	Rev. No.	00
	Rev. Date	2022. 01. 20

XI. Toxicological Information

	Detection buffer
Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

XII. Ecological Information

	Detection buffer
Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

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