VIDAS[®] TOXO IgM (TXM)

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VIDAS[®] CMV IgG (CMVG)

IVD

VIDAS CMV IgG is an automated quantitative enzyme immunoassay for use on the VIDAS family instruments for the quantitative measurement, of anti-cytomegalovirus IgG (CMVG) in human serum, using the technique ELFA (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Cytomegalovirus (CMV) is a herpes virus which can cause serious disease in infants and adults. CMV can persist in the human body for several years and can cause recurrent infections or be transmitted to other individuals.

CMV infections are very common; 60 to 85% of the population have been infected but most cases are asymptomatic.

1 to 3% of women are infected during pregnancy and in 1 out of every 2 cases the infection is passed on to the fetus. The infection is most often asymptomatic but in about 5% of cases, the consequences can be very severe: hepatosplenomegaly, hydrocephalus, microcephalism, prematurity. Fetal death is frequent (1,2). Even in cases of asymptomatic infection, approximately 10% of infants show sensorineural after-effects, such as deafness or partial or total blindness.

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The detection of IgG antibodies to CMV can be useful in the diagnosis of current primary infection, particularly in pregnant women (4, 5, 6).

The detection of anti-cytomegalovirus IgG is useful in assessing the immune status of patients.

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After dilution, the sample is incubated with the SPR. Anti-CMV IgG antibodies present in the specimen will bind to the CMV antigen coating the interior of the SPR. Unbound components are eliminated during the washing steps.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibodies present in the sample.

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CONTENT OF THE KIT (60 TESTS):

60 CMVG strips	STR	Ready-to-use.
60 CMVG SPRs 2 x 30	SPR	Ready-to-use. SPRs coated with CMV antigen (strain AD169).
CMVG Positive control 1 x 1.5 mL (liquid)	C1	Ready-to-use. Human serum* containing anti-CMV IgG + 1 g/L sodium azide. MLE data indicate the confidence interval in aU/mL (arbitrary unit/mL) ("Control C1 (+) Dose Value Range").
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin + preservatives.
CMVG Calibrator 1 x 2 mL (liquid)	S1	Ready-to-use. Human serum* containing anti-CMV IgG + 1 g/L sodium azide. MLE data indicate the concentration in aU/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value ("Calibrator (S1) RFV Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR[®]

The interior of the SPR is coated during production with purified CMV antigen. Each SPR is identified by the CMVG code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the CMVG Strip

Wells	Reagents
1	Sample well.
2	Sample diluent: Phosphate buffer (10 mmol/L) pH 7.2 - Tween + protein and chemical stabilizers + 1 g/L sodium azide (300 µL).
3	Pre-wash solution: Phosphate (10 mmol/L) pH 7.2 - Tween + protein and chemical stabilizers + 1 g/L sodium azide (600 µL).
4 - 5 - 7 - 8	Wash solution: TRIS (50 mmol/L) pH 7.4 + 0.9 g/L sodium azide (600 µL).
6	Conjugate: Alkaline phosphatase-labeled monoclonal anti-human IgG antibodies (mouse) + 1 g/L sodium azide (400 µL).
9	Empty well.
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine DEA* (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

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- Pipette with disposable tip to dispense 100 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instrument of the VIDAS[®] family.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The optical cuvette with substrate (well 10) contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

STORAGE CONDITIONS

- Store the VIDAS CMV IgG kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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The use of clearly hemolyzed, lipemic, icteric or inactivated samples has not been validated.

Specimen Stability

Samples can be stored at 2-8°C in stoppered tubes for up to 5 days; if longer storage is required, freeze the sera at -25 ± 6°C.

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The calibrator, identified by S1, must be tested **in duplicate** (see User Manual). The calibrator value must be within the set RFV (Relative Fluorescence Value) range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "CMVG" strip and one "CMVG" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**

3. The test is identified by the code "CMVG" on the instrument (to do so, refer to the Instrument User Manual). The calibrator must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by C1. If the negative control is to be tested, it should be identified by C2.

4. Mix the calibrator, controls and samples using a Vortex-type mixer (for serum separated from the pellet).

5. **For this test, the calibrator, control, and sample test portion is 100 µL.**

6. Insert the "CMVG" SPRs and "CMVG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.

7. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.

8. Restopper the vials and return them to 2-8°C after pipetting.

9. The assay will be completed within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.

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The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistic model). The concentrations are expressed in aU/mL.

Samples with anti-CMV IgG concentrations > 400 aU/mL should be reassayed after dilution by 1/4 in saline solution. If the dilution factor has not been entered when the Work List was created (see User Manual), multiply the result by the dilution factor to obtain the sample concentration. For serological monitoring, to obtain an accurate serum titer, dilute the samples with concentrations > 150 aU/mL in order to obtain a result ranging between 50-150 aU/mL.

Since no international standard is available for the determination of anti-CMV IgG, the VIDAS CMV IgG reagent is calibrated against collection sera.

Thresholds and Interpretation of Results

Value (aU/mL)	Interpretation
< 4	Negative
from ≥ 4 to < 6	Equivocal
≥ 6	Positive

It is recommended to confirm equivocal results using a second sample.

Critical Ratio

The critical ratio determines a significant increase in antibody level. This increase is indicative of an active infection and is calculated as follows:

$$\text{Critical ratio} = \frac{2^{\text{nd}} \text{ collection of sample aU/mL}}{1^{\text{st}} \text{ collection of sample aU/mL}}$$

Critical ratio	Interpretation
< 2	No significant change in antibody level in 25% of cases
2 - 4	Suggestive of an increase in antibody level in 15% of cases.
> 4	Highly significant increase in antibody level in 60% of cases.

If the critical ratio increases by a factor of less than 4, it is recommended to collect a new specimen within 2 to 3 weeks.

This critical ratio must in no case be used to exclude an active infection.

For serological monitoring, it is recommended to test the samples from the different collections in the same run using the same lot.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS CMV IgG kit.

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LIMITATIONS OF THE METHOD

1. The VIDAS CMVG assay results should be interpreted taking into account the clinical symptoms and the medical analysis results from other tests (e.g. detection of anti-CMV IgM in case of primary infection).
2. Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient history, and the results of any other tests performed
3. VIDAS CMV IgG has not been validated for use with specimens collected post mortem.

PERFORMANCE

Studies performed using VIDAS CMVG gave the following results:

Study 1: During evaluation at an external site, 284 samples were tested in comparison with another commercially available EIA.

Discrepant samples were checked using another method. Equivocal samples were excluded from the performance calculation.

N = 284		VIDAS CMV IgG		
		+	Equivocal	-
EIA results	+	169	10	3
	Equivocal	0	1	0
	-	0	1	100

Sensitivity: 98.26%. Range at 95%: 94.90 - 99.42%

Specificity: 100%. Range at 95%: 96.15 - 100%

Agreement: 98.9%. Range at 95%: 96.75 - 99.63%

Study 2: During evaluation at an external site, 396 samples were tested in comparison with another commercially available EIA.

Discrepant samples were checked using another method. Equivocal samples were excluded from the performance calculation.

N = 396		VIDAS CMV IgG		
		+	Equivocal	-
EIA results	+	241	3	2
	Equivocal	0	0	2
	-	0	10	138

Sensitivity: 99.18%. Range at 95%: 96.98 - 99.78%

Specificity: 100%. Range at 95%: 97.18 - 100%

Agreement: 99.5%. Range at 95%: 98.06 - 99.86%

Precision

Within-Run Reproducibility:

3 samples were tested 30 times in a same run.

Sera	Mean RFV	Intra-Run CV %
Negative	60	4
Weakly positive	1151	2.5
Strongly positive	2784	2.1

Between-Run Reproducibility:

2 samples were tested singly in 8 different runs on the same VIDAS instrument.

Sera	Average Titer aU/mL	Between-Run CV %
Weakly positive	34.5	4.1
Strongly	296	15

Sera	Mean RFV	Between-Run CV %
Weakly positive	2926	2.8
Strongly positive	5434	1.6

RANGE OF EXPECTED VALUES

The incidence of serum anti-CMV IgG differs according to the populations studied. In developed countries, the percentage of CMV positive women between 15-30 years of age does not exceed 40-50%.

The curve showing IgG acquisition as a function of age reveals 2 peaks: one in the first year of life, the other after puberty, proving that the infection is sexually transmissible.

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LITERATURE REFERENCES

1. STARR S.E., and FRIEDMAN H.M. "Human Cytomegalovirus", Chapter 65 in Manual of Clinical Microbiology, fourth edition, edited by E.H. Lennette, A. Balows, W.J. Hausler, Jr., H.J. Shadomy, American Society for Microbiology, 1985, 711-719.
2. PECKHAM C.S. "Cytomegalovirus in the neonate". Journal of Antimicrobial Chemotherapy. 1989, 23, 17-21.
3. CHOU S., et al. "Immunoglobulin M to cytomegalovirus in primary and reactivation infections in renal transplant recipients". Journal of Clinical Microbiology. 1987, 25, 52-55.
4. ANDERSSON J. Cytomegalovirus Infection in Pregnancy, Scand. J. Infect. Dis., 1990, Suppl. 71, 67-70.
5. HORODNICEANU F., MICHELSON S. "Assessment of human cytomegalovirus antibody detection techniques". Archives of Virology. 1980, 64, 287-301.
6. JOASSIN L., REGINSTER M. "Elimination of Nonspecific Cytomegalovirus Immunoglobulin M, Activities in the Enzyme-Linked Immunosorbent Assay by Using Anti-Human Immunoglobulin G", Journal of Clinical Microbiology, 1986, 23, 576-581.

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Release date	Part Number	Change Type	Change Summary
2015/01	05936Q	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 tests) WARNINGS AND PRECAUTIONS
2015/06	05936R	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE
2016/10	05936S	Technical	CONTENT OF THE KIT (60 TESTS) PERFORMANCE

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1. STARR S.E., and FRIEDMAN H.M. "Human Cytomegalovirus", Chapter 65 in Manual of Clinical Microbiology, fourth edition, edited by E.H. Lennette, A. Balows, W.J. Hausler, Jr., H.J. Shadomy, American Society for Microbiology, 1985, 711-719.
2. PECKHAM C.S. "Cytomegalovirus in the neonate". Journal of Antimicrobial Chemotherapy.1989, 23, 17-21.
3. CHOU S., et al. "Immunoglobulin M to cytomegalovirus in primary and reactivation infections in renal transplant recipients". Journal of Clinical Microbiology.1987, 25, 52-55.
4. ANDERSSON J. Cytomegalovirus Infection in Pregnancy, Scand. J. Infect. Dis., 1990, Suppl.71, 67-70.
5. HORODNICEANU F., MICHELSON S. "Assessment of human cytomegalovirus antibody detection techniques". Archives of Virology.1980, 64, 287-301.
6. JOASSIN L., REGINSTER M. "Elimination of Nonspecific Cytomegalovirus Immunoglobulin M, Activities in the Enzyme-Linked Immunosorbent Assay by Using Anti-Human Immunoglobulin G", Journal of Clinical Microbiology, 1986, 23, 576-581.

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VIDAS[®] CMV IgM (CMVM)

IVD

VIDAS CMV IgM is an automated qualitative enzyme immunoassay for use on the VIDAS family instruments, for the detection of anti-cytomegalovirus IgM (CMVM) in human serum, using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Cytomegalovirus (CMV) is a herpesvirus which can cause serious diseases in infants and adults. CMV can persist in the human body for several years and can cause recurrent infections or be transmitted to other individuals. CMV infections are very common; 60 to 85% of the population have been infected but most cases are asymptomatic.

1 to 3% of women are infected during pregnancy and in 1 out of every 2 cases the infection is passed on to the fetus. The infection is most often asymptomatic but in about 5% of cases, the consequences can be very severe: hepatosplenomegaly, hydrocephalus, microcephalism, prematurity. Fetal death is frequent (1, 2).

Even in cases of asymptomatic infection, approximately 10% of infants show sensorineural after-effects, such as deafness or partial or total blindness.

In immunocompromised patients (HIV positive or organ transplanted), CMV infections can be severe (3).

The detection of anti-CMV IgM can be useful in the diagnosis of current primary infection, particularly in pregnant women (4, 5, 6).

The anti-CMV IgM present in approximately 70% of primary infections, generally persist 16 to 20 weeks after infection, and may reappear irregularly during reactivation (7).

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The Solid Phase Receptacle (SPR[®]) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

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After IgG and rheumatoid factor adsorption, the sample is cycled in and out of the SPR for a specified length of time. Anti-CMV IgM antibodies present in the specimen will bind to the CMV antigen coating the interior of the SPR. Unbound components are eliminated during the washing steps.

An Alkaline phosphatase-labeled mouse monoclonal anti-human IgM antibody is cycled in and out of the SPR and binds to the human anti-CMV IgM coated on the SPR. A final wash step removes unbound components.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (30 TESTS):

30 CMVM Strips	STR	Ready-to-use.
30 CMVM SPRs	SPR [®]	Ready-to-use. Interior of SPRs coated with CMV antigen (cell culture of the virus of the AD169 strain).
CMVM Positive Control 1 x 0.5 mL (liquid)	C1	Ready-to-use. Human serum* containing anti-CMV IgM + 1 g/L sodium azide. MLE data indicate the index: confidence interval ("Control C1 (+)Test Value Range).
Negative Control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin + preservatives.
CMVM Standard 1 x 1 mL (liquid)	S1	Ready-to-use. Human serum* containing anti-CMV IgM + 1 g/L sodium azide.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

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The SPR

The interior of the SPR is coated during production with purified CMV antigen. Each SPR is identified by the CMVM code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

The Strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay

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- Pipette with disposable tip to dispense 100 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory Biosafety Manual - WHO - Geneva - latest Edition).

- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.

- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

STORAGE CONDITIONS

- Store the VIDAS CMV IgM kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen Type and Collection

Human serum.

It is recommended that each laboratory checks the compatibility of used collection tubes.

Since the use of hemolyzed, lipemic, icteric or inactivated samples has not been validated, it is recommended to collect a new specimen.

Specimen Stability

Samples can be stored at 2-8°C in stoppered tubes for up to 5 days; if longer storage is required, freeze the sera at $-25 \pm 6^\circ\text{C}$.

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INSTRUCTIONS FOR USE

For complete instructions, see the User Manual.

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Note: the master lot data need to be entered only once for each lot.

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The standard, identified by S1, must be tested **in duplicate** (see User Manual). The standard value must be within the set RFV (Relative Fluorescence Value) range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "CMVM" strip and one "CMVM" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "CMVM" code on the instrument. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by C1. If the negative control is to be tested, it should be identified by C2.
4. Mix the standard, controls and samples using a vortex-type mixer (for serum separated from the pellet).
5. **For this test, the standard, control, and sample test portion is 100 µL.**
6. Insert the "CMVM" SPRs and "CMVM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
8. Restopper the vials and return them to 2-8°C after pipetting.
9. The assay will be completed within approximately 60 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate receptacle.

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The index (ratio of the fluorescent signal found for the serum to be tested, over the standard signal stored in the memory) is calculated by the instrument.

This index and the interpretation are indicated on the result sheet.

Threshold and Interpretation of Results

Index	Interpretation
$i < 0.70$	Negative
$0.70 \leq i < 0.90$	Equivocal
$i \geq 0.90$	Positive

Repeat the assay for indexes between 0.70 and 0.90. If after verification interpretation is still equivocal, repeat the assay on a new sample collected 10 to 15 days later.

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As no international standard is available for the determination of anti-CMV IgM, the VIDAS CMV IgM reagent is calibrated against collection sera.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS CMVM kit.

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LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted as part of a complete clinical profile.

The detection of anti-CMV IgM may be difficult in immunocompromised patients. The diagnosis of cytomegalovirus infection or non-infection in immunocompromised patients should not be based solely on the determination of anti-CMV IgM. Interpretation of VIDAS CMV IgM test results should be made taking into consideration the patient's history and the results of any other tests performed (viral culture, viral DNA, anti-CMV IgG, etc.).

VIDAS CMV IgM has not been validated for use with specimens collected post mortem or neonatal specimens (cord blood, etc.).

PERFORMANCE

The VIDAS CMV IgM reagent was evaluated at an external site in comparison with 3 other commercially available EIA methods.

The presence or absence of anti-CMV IgM in each sample was determined by taking into account the concordant results obtained with at least 2 or 3 of the commercially available tests.

Equivocal samples were excluded from the performance calculation.

Sensitivity Study of VIDAS CMV IgM

During the evaluation, 43 typically anti-CMV IgM positive samples were tested. The following results were obtained with VIDAS CMV IgM:

N = 43		VIDAS CMV IgM		
		Positive	Equivocal	Negative
Status	Positive	37	2	4*

* The 4 specimens that were found to be negative were all collected from transplant patients.

Sensitivity: 90.24%

95% Confidence Interval: 76.87 – 97.28%.

Specificity Study of VIDAS CMV IgM

During the evaluation, 171 typically anti-CMV IgM-negative samples were tested. The following results were obtained with VIDAS CMV IgM:

N = 171		VIDAS CMV IgM		
		Positive	Equivocal	Negative
Status	Negative	1	1	169

Specificity: 99.41%

95% Confidence Interval: 96.65 – 99.90%.

Specificity Study of VIDAS CMV IgM on potentially interfering samples

This study was carried out using 171 serum samples with the following characteristics:

- 40 serum samples from patients with HIV1 virus,
- 11 serum samples from infants less than 3 months old,
- 57 serum samples from pregnant women,
- 17 serum samples from carriers of Epstein-Barr virus,
- 5 serum samples from carriers of varicella virus,
- 5 serum samples from carriers of herpes virus,
- 16 serum samples from subjects possessing antinuclear antibodies,
- 20 serum samples from subjects possessing rheumatoid factor.

158 out of the 171 serum samples were diagnosed as anti-CMV IgM-negative. The following results were obtained with VIDAS CMV IgM:

N = 158		VIDAS CMV IgM		
		Positive	Equivocal	Negative
Status	Negative	2*	2	154

*One sample belongs to the group of subjects with rheumatoid factor, the second to the group of patients with antinuclear antibodies.

Specificity: 98.72% 95%

Confidence Interval: 95.34 – 99.66%.

Within-Run Reproducibility

3 samples were tested 30 times in the same run.

Sera	Mean Index	Standard Deviation	Intra-Assay CV %
High positive	1.63	0.09	5.52
Low positive	1.00	0.06	6.00
Negative	0.08	0.01	12.50

Between-Run Reproducibility

3 samples were tested singly in 10 different runs on the same VIDAS instrument.

Sera	Mean Index	Standard Deviation	Inter-Assay CV %
High positive	1.51	0.06	3.97
Low positive	0.96	0.05	5.21
Negative	0.14	0.06	42.86

The figures in the above tables are given for information only.

RANGE OF EXPECTED VALUES

The incidence of cytomegalovirus infection differs according to the populations studied. In developed countries, the percentage of CMV-positive women between 15-30 years of age does not exceed 40-50%.

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LITERATURE REFERENCES

- STARR S.E., and FRIEDMAN H.M. "Human Cytomegalovirus", Chapter 65 in Manual of Clinical Microbiology, fourth edition, edited by E.H. Lennette, A. Balows, W.J. Hausler, Jr., H.J. Shadomy, American Society for Microbiology, 1985, 711-719.
- PECKHAM C.S. "Cytomegalovirus in the neonate". Journal of Antimicrobial Chemotherapy. 1989, 23, 17-21.
- CHOU S., et al. "Immunoglobulin M to cytomegalovirus in primary and reactivation infections in renal transplant recipients". Journal of Clinical Microbiology. 1987, 25, 52-55.
- ANDERSSON J. Cytomegalovirus Infection in Pregnancy, Scand. J. Infect. Dis., 1990, Suppl. 71, 67-70.
- HORODNICEANU F., MICHELSON S. "Assessment of human cytomegalovirus antibody detection techniques". Archives of Virology. 1980, 64, 287-301.
- JOASSIN L., REGINSTER M. "Elimination of Nonspecific Cytomegalovirus Immunoglobulin M, Activities in the Enzyme-Linked Immunosorbent Assay by Using Anti-Human Immunoglobulin G", Journal of Clinical Microbiology, 1986, 23, 576-581.
- RANGER-ROGEZ S., VENOT C., AUBARD Y. et al. CytomégaloVirus in: Les virus transmissibles de la mère à l'enfant DENIS F. Ed, 1999, 214-239.

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Release date	Part Number	Change Type	Change Summary
2015/01	06637K	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (30 TESTS) WARNINGS AND PRECAUTIONS
2015/06	06637L	Technical	CONTENT OF THE KIT (30 TESTS) INSTRUCTIONS FOR USE
2016/10	06637M	Technical	CONTENT OF THE KIT (30 TESTS)

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1. STARR S.E., and FRIEDMAN H.M. "Human Cytomegalovirus", Chapter 65 in Manual of Clinical Microbiology, fourth edition, edited by E.H. Lennette, A. Balows, W.J. Hausler, Jr., H.J. Shadomy, American Society for Microbiology, 1985, 711-719.
2. PECKHAM C.S. "Cytomegalovirus in the neonate". Journal of Antimicrobial Chemotherapy. 1989, 23, 17-21.
3. CHOU S., et al. "Immunoglobulin M to cytomegalovirus in primary and reactivation infections in renal transplant recipients". Journal of Clinical Microbiology. 1987, 25, 52-55.
4. ANDERSSON J. Cytomegalovirus Infection in Pregnancy, Scand. J. Infect. Dis., 1990, Suppl. 71, 67-70.
5. HORODNICEANU F., MICHELSON S. "Assessment of human cytomegalovirus antibody detection techniques". Archives of Virology. 1980, 64, 287-301.
6. JOASSIN L., REGINSTER M. "Elimination of Nonspecific Cytomegalovirus Immunoglobulin M, Activities in the Enzyme-Linked Immunosorbent Assay by Using Anti-Human Immunoglobulin G", Journal of Clinical Microbiology, 1986, 23, 576-581.
7. RANGER-ROGEZ S., VENOT C., AUBARD Y. et al. Cytomégalo virus in: Les virus transmissibles de la mère à l'enfant DENIS F. Ed, 1999, 214-239.

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VIDAS[®] TOXO IgG II (TXG)

VIDAS TOXO IgG II is an automated quantitative test for use on the VIDAS family instruments for the quantitative measurement of anti-toxoplasma IgG in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Toxoplasma gondii, (1) an obligate intracellular protozoan parasite, is a significant pathogen among humans. The parasite, whose primary host is the Felidae, is scattered in nature and invades all orders of mammals.

Toxoplasmosis is usually benign or asymptomatic, but can have severe consequences if it occurs in immunodeficient subjects or fetuses (2). Women who are seropositive before they become pregnant are essentially protected from transmitting the infection to their unborn child (3, 4, 5, 6, 7). Women who are seronegative are at risk of becoming infected during gestation. Transmission to the fetus occurs during the acute stage of the infection in the mother. The frequency with which *Toxoplasma* crosses the placental barrier to cause congenital infection will depend on a number of factors including the virulence of the infecting strain of the parasite, the size of the original inoculum, the immune response of the patient and when, during gestation, the mother became infected (8).

The prevalence of toxoplasmosis varies depending upon geographical location, age and gender of the population studied and other factors. In Europe, the prevalence rate ranges from 20% to 85%. In the United States, the prevalence is lower: 12% to 41%. Prevalence in other countries can vary from 18% to 65%.

The diagnosis of *T.gondii* infection is most commonly made by biological examination: specific immunoglobulin detection (IgM and IgG)(9, 10).

The diagnosis of the acute acquired infection during pregnancy is established by demonstration of a seroconversion, or by a significant rise in antibody titer in sequential sera assayed concomitantly (6, 11).

VIDAS TOXO IgG II is a screening test intended for use as an aid in the determination of the immunity status of patients.

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After a sample dilution step, the sample is cycled in and out of the SPR. Anti-*T.gondii* IgG antibodies present in the specimen will bind to the *T.gondii* antigen coating the interior of the SPR. Unbound components are eliminated during the washing steps. Mouse monoclonal anti-human IgG conjugated with alkaline phosphatase is cycled through the SPR and will attach to any human IgG bound to the SPR wall. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibodies present in the sample. At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (60 TESTS):

60 TXG strips	STR	Ready-to-use.
60 TXG SPRs 2 x 30	SPR [®]	Ready-to-use. SPRs coated with membrane and cytoplasmic Toxoplasma antigen, RH Sabin strain grown in mice (12).
TXG positive control 1 x 2 mL (liquid)	C1	Ready-to-use. Human serum* containing anti-Toxoplasma IgG + protein stabilizer + 1 g/L sodium azide. MLE data indicate the confidence interval in IU/mL (international unit per milliliter) ("Control C1 (+) Dose Value Range").
TXG negative control 1 x 3 mL (liquid)	C2	Ready-to-use. Human serum* negative for anti-Toxoplasma IgG + protein stabilizer + 1 g/L sodium azide.
TXG calibrator 1 x 1 mL (liquid)	S1	Ready-to-use. Human serum* containing anti-Toxoplasma IgG and calibrated against the 2 nd WHO International Standard + protein stabilizer + 1 g/L sodium azide. MLE data indicate the concentration in IU/mL ("Calibrator (S1) Dose Value") and the confidence interval in "Relative Fluorescence Value ("Calibrator (S1) RFV Range").
Specifications for the factory master data required to calibrate the test:		
• MLE data (Master Lot Entry) provided in the kit, - or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR®

The interior of the SPR is coated during production with membrane and cytoplasmic Toxoplasma antigen (RH Sabin strain). Each SPR is identified by the TXG code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the TXG strip

Wells	Reagents
1	Sample well.
2	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
3	Pre-wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
4 - 5 - 7 - 8	Wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
6	Conjugate: Alkaline phosphatase-labeled monoclonal anti-human IgG antibodies (mouse) + 0.9 g/l sodium azide (400 µl).
9	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (400 µl).
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

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- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS® family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of

transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The optical cuvette with substrate (well 10) contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS TOXO IgG II kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection:

Serum or plasma (lithium heparin or EDTA). It is recommended that each laboratory checks the compatibility of collection tubes used.

Since the use of hemolyzed, lipemic, or icteric samples has not been validated, it is recommended to collect a new specimen.

Sera inactivated at 56°C for 30 minutes can be tested with VIDAS TOXO IgG II.

Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 5 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C.

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A study performed on frozen samples over a period of one year, showed that the quality of results is not affected.

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Calibration

Calibration, using the calibrator provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested in **duplicate** (see User's Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "TXG" strip and one "TXG" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "TXG" code on the instrument. The calibrator must be identified by "S1", and tested in duplicate. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.
4. Mix the calibrator, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

5. **For this test, the calibrator, control, and sample test portion is 100 µl.**

6. Insert the "TXG" SPRs and "TXG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Restopper the vials and return them to 2-8°C after pipetting.
9. The assay will be completed within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

NOTE: If a comparison needs to be made with a previously collected specimen, test both sera in the same run.

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The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistic model) and are expressed in « IU/ml » (2nd WHO International Standard). Equivocal samples should be retested. If the interpretation remains equivocal, a new sample must be collected.

Samples with an IgG concentration greater than 300 IU/ml should be reassayed after diluting 1/4 with the negative control.

If the dilution factor was not entered when the Work List was created (see User's Manual), multiply the result by the dilution factor to obtain the sample concentration.

Samples with concentrations greater than 300 IU/ml will be reported as "greater than 300 IU/ml". Samples with high anti-toxoplasma IgM titers may affect dilution test results.

For serological monitoring, it is recommended to test the samples from the different collections, in the same run, using the same lot.

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QUALITY CONTROL

A positive and a negative control are included in each VIDAS TOXO IgG II kit.

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LIMITATIONS OF THE METHOD

1. Sera collected very early in the acute stage of disease may have IgG levels < 4 IU/ml.
2. Positive test results from cord blood should be interpreted with caution. The presence of total or IgG anti-*T. gondii* antibodies in cord blood is usually the result of passive transfer from the mother to the fetus.
3. Positive test results may not be valid in persons who have received blood transfusions or other blood products within the past several months.

4. VIDAS TOXO IgG II results should be used in conjunction with clinical symptoms and laboratory test results, such as anti-*T. gondii* IgM results.

5. Test results from immunosuppressed patients may be difficult to interpret due to diminished immune response.

6. The VIDAS TOXO IgG II assay may only be used with serum or plasma (EDTA, heparin) samples. The use of amniotic fluid or other body fluids has not been established.

7. Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

PERFORMANCE

Studies performed using VIDAS TOXO IgG II gave the following results:

1. Sensitivity - Specificity

1940 samples were evaluated on 3 sites in comparison with another immunoenzymatic reference technique.

49 samples were found to be equivocal with one and/or two other techniques, and were not included in the analysis.

Consolidated results for the 3 sites are as follows:

		EIA Technique		
		Positive	Negative	Total
VIDAS	Positive	572	3	575
	Negative	9	1307	1316
Total		581	1310	1891

The 12 discrepant samples were analyzed using a third reference technique: the Sabin-Feldmann test (Dye Test). The definitive interpretation for each of these samples is that obtained with at least two out of the three techniques used (2/3 rule).

Of the 3 VIDAS TOXO IgG II-positive/EIA-negative samples, 2 were positive according to the 2/3 rule.

Of the 9 VIDAS TOXO IgG II-negative/EIA-positive samples, 6 were negative according to the 2/3 rule and 1 could not be tested using the reference technique due to insufficient quantity.

Consolidated results after resolving discrepant samples are as follows:

		2/3 Rule		
		Positive	Negative	Total
VIDAS	Positive	574	1	575
	Negative	2	1313	1315
Total		576	1314	1890

Sensitivity: 99.65 %

(Range at 95 %: 98.75 – 99.96 %)

Specificity: 99.92 %

(Range at 95 %: 99.58 % - 100 %)

Precision

Precision was tested with a negative sample, a low positive sample and a high positive sample.

Each sample was tested in duplicate two times per day for twenty days.

Within-run reproducibility (intra-assay precision) and between-run reproducibility (total precision) were calculated according to the recommendations of NCCLS Document EP5-T2, volume 12 number 4.

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No cross-reactivity or interference in the VIDAS TOXO IgG II assay was seen with any of the samples tested positive for rheumatoid factor, antinuclear antibodies, or Epstein Barr virus.

RANGE OF EXPECTED VALUES

Toxoplasma gondii is a strict pathogen whose prevalence differs from one country to another or even one region to another. Contamination by *T. gondii* can vary according to cultural customs and eating habits, resulting in a prevalence ranging from less than 10% in certain regions of Northern Europe to more than 90% in Africa.

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Release date	Part Number	Change Type	Change Summary
2015/01	09065I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS
2015/06	09065J	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE

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VIDAS[®] Varicella-Zoster IgG (VZG)

VIDAS Varicella-Zoster IgG (VZG) is an automated qualitative test for use on the instruments of the VIDAS family, for the detection of IgG antibodies to varicella zoster virus in human serum using the ELFA technique (Enzyme Linked Fluorescent Assay). VIDAS VZG is intended as an aid in the determination of immunological experience with varicella-zoster virus.

SUMMARY AND EXPLANATION

Varicella (chickenpox) and zoster (shingles) are two different clinical manifestations of infection with the varicella-zoster virus (VZV). Chickenpox, the primary infection, is a highly contagious disease which spreads rapidly in susceptible populations (e.g. children or immunocompromised adults) (1, 2, 3). It is spread by direct contact or by airborne droplets. The contagious period is from 24 hours before until six or seven days after the appearance of vesicles (4). In most cases, chickenpox is a mild illness occurring in childhood, characterized by a mild fever and malaise accompanied by a vesicular rash (4). However, adults with chickenpox usually have a high fever and severe constitutional symptoms. Approximately 20% develop pneumonia (2).

Shingles generally occurs in adults, resulting from a reactivation of latent VZV. Although there is often no apparent reason for zoster to occur, it often is linked to stress, trauma, or certain types of drug therapy.

VZV infection in immunosuppressed patients can cause serious or fatal disease. Progressive, generalized varicella occurs in up to 30% of children who acquire chickenpox while on cancer therapy. Children with deficient cell-mediated immunity also are at risk for progressive varicella. In older immunodeficient patients there is an increased risk of disseminated zoster (1,3).

Laboratory diagnosis of VZV infection can be accomplished via viral isolation by culture, viral antigen detection, or by serological methods. Serological methods are preferred over viral isolation techniques because serologies give results rapidly and because the virus is extremely labile. Serological methods include neutralizing antibody, complement fixation (CF), fluorescent antibody to membrane antigen (FAMA), anticomplement immunofluorescence (ACIF) (5) and enzyme-linked immunosorbent assay (ELISA).

Measurement of FAMA is rapid and more sensitive than CF (6). ELISA methods are easier to perform than FAMA, are extremely sensitive, and do not require propagation of the virus (7).

In most cases, VZV infection can be diagnosed by observing clinical symptoms. However ELISA is valuable for determining immunological experience particularly in adults whose memory of previous varicella may be inaccurate (8).

The VIDAS Varicella-Zoster IgG (VZG) assay is intended as an aid in the determination of immunological experience with varicella-zoster virus.

PRINCIPLE

The assay principle combines a 2-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR[®]), serves as the solid phase as well as the pipetting device for the assay.

Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

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After a preliminary wash step and a sample dilution step, anti-varicella IgG antibodies present in the specimen will bind to the varicella antigen coating the interior of the SPR. Unbound sample components are washed away.

Anti-human IgG conjugated with alkaline phosphatase bind to any human IgG bound to the SPR wall. A final wash step removes unbound conjugate.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of anti-varicella IgG antibodies present in the sample.

At the end of the assay, results are automatically analyzed by the instrument, a test value is generated and a report is printed for each sample.

CONTENT OF THE KIT (60 TESTS):

60 VZG strips	STR	Ready-to-use.
60 VZG SPRs (2 x 30)	SPR	Ready-to-use. Interior of SPRs coated with inactivated varicella-zoster viral antigen (Sissy strain).
VZG standard 1 x 2 ml (liquid)	S1	Ready-to-use. Human* defibrinated plasma containing antibodies to varicella + 1 g/l sodium azide. MLE data indicate the confidence interval in "Relative Fluorescence Value (RFV)" ("Standard (S1) RFV Range").
VZG positive control 1 x 2 ml (liquid)	C1	Ready-to-use. Human* defibrinated plasma containing antibodies to varicella + 1 g/l sodium azide. MLE data indicate the Test Value (TV) range: Control C1 (+) Test Value Range.
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives. MLE data indicate the Test Value (TV) range: Control C2 (-) Test Value Range.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar codes printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR

The interior of the SPR is coated during production with inactivated varicella-zoster viral antigen (Sissy strain). Each SPR is identified by the VZG code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the VZG strip:

Wells	Reagents
1	Sample well: a minimum of 100 µl of serum is placed into the well.
2	Sample diluent: 400 µl of phosphate buffered saline (0.01 mol/l, pH 7.2) + bovine albumin + 1 g/l sodium azide.
3	Prewash buffer: Phosphate buffered saline (0.01 mol/l ; pH = 7.2) + bovine albumin + 1 g/l sodium azide. (400 µl).
4 - 5	Wash buffer: casein (bovine) in TRIS buffered saline (0.05 mol/l, pH 7.4) + 1 g/l sodium azide (600 µl).
6	Conjugate: prediluted mixture of anti-human IgG (mouse) conjugated to alkaline phosphatase + bovine albumin + 1 g/l sodium azide (400 µl).
7 - 8	Wash buffer: casein (bovine) in TRIS buffered saline (0.05 mol/l; pH = 7.4) + 1 g/l sodium azide (600 µl).
9	Empty well.
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6 %) pH 9.2 + 1 g/l sodium azide (300 µl).

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MATERIALS REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials, please refer to the Instrument Operator's Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Consider all patient specimens potentially infectious and observe routine biosafety precautions. Though the varicella-zoster virus coating the SPRs has been inactivated, handle the SPRs as if they were infectious. Dispose of all used components and other contaminated materials by acceptable procedures for potentially biohazardous human blood products.
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests (9).
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the Operator's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the Operator's Manual).

STORAGE CONDITIONS

- Store the VIDAS VZG kit at 2-8° C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection:

Whole blood should be collected and the serum separated by standard procedures. Do not heat the serum. Clarify samples containing particulate matter by centrifugation or filtration prior to testing. Lipemic serum may interfere with the assay. Do not use samples that are clearly hemolyzed, lipemic or icteric. Collect a new sample.

The use of plasma has not been established for this test. Do not test specimens with obvious microbial contamination.

Specimen stability:

If the sera cannot be tested on the day of collection, store them at 2 - 8°C in stoppered tubes for up to five days. If longer storage is required, freeze the sera at -25 ± 6°C. Avoid successive freezing and thawing.

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Calibration

Calibration, using the standard provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration data and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by "S1", must be tested in **duplicate** (see Operator's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.
2. Use one "VZG" strip and one "VZG" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "VZG" code on the instrument. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
4. **If needed, label the Reagents Strips with the appropriate sample identification in the space provided.**
5. Mix the standard, controls and sample using a vortex-type mixer (for serum separated from the pellet).
6. **For this test, the standard, control, and sample test portion is 100 µl.**
NOTE: When manual pipetting is performed, check the sample wells for bubbles. If necessary, gently tap the strips to remove any bubbles present.
7. Insert the "VZG" SPRs and "VZG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the Operator's Manual. All the assay steps are performed automatically by the instrument.
9. Reclose the vials and return them to 2–8°C after pipetting.
10. The assay will be completed within approximately 35–40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result.

A test value is generated for each sample by forming a ratio from the RFV of the sample to that of a standard. Test values from patient and control samples are compared to a set of thresholds stored in the computer representing a low varicella antibody titer.

The following table shows the thresholds and the interpreted results.

Test value	Interpretation
$TV < 0.60$	Negative
$0.60 \leq TV < 0.90$	Equivocal
$TV \geq 0.90$	Positive

Results with test values less than the negative threshold (0.60) indicate patients with undetectable concentrations of anti-VZV IgG antibodies. Patients with a current primary infection of varicella may not begin producing measurable concentrations of anti-VZV IgG until several days after infection.

Samples with test values greater than or equal to the threshold (0.90) are reported as positive.

A report is printed which records

- the type of test performed,
- the sample identification,
- the date and time,
- the lot number and the expiration date of the reagent being used,
- each sample's RFV, test value and interpreted result.

The imprecision inherent in any method implies a lack of confidence in samples with test values very close to the thresholds. Consequently, an equivocal zone is established between the thresholds based on a statistical understanding of this imprecision.

Equivocal samples should be repeated with a fresh specimen. If it is impractical to obtain a fresh specimen, repeat the assay with the original specimen (provided it is not microbially contaminated or lipemic). If the sample repeats as an equivocal, test it using another method.

Invalid results are reported when the background reading is above a predetermined cutoff (indicating low-level substrate contamination). In this case, repeat the assay with the original specimen.

An invalid result is also seen if there is no standard available for the lot number of the patient test strip. In this case, run a standard in duplicate in VZG strips with the same lot number as the invalid patient test. The patient test result can then be recalculated using the new stored standard, (See the Operator's Manual for complete information).

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS VZG kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2. The expected values for the VIDAS VZG assay are entered into the instrument using the master lot entry (MLE) data. If the results from testing the controls do not fall within these ranges, the result will be flagged – do not report patient results.

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LIMITATIONS OF THE METHOD

1. The VIDAS Varicella-Zoster IgG assay should be used as a means of determining the immunological experience of an individual through the detection of anti-varicella IgG antibodies.
2. The results of this test are qualitative and should be considered as either positive or negative for the presence of antibodies to VZV.
3. The presence of anti-varicella IgG antibodies indicates prior exposure to the virus but does not exclude the possibility of an ongoing infection. Quantitative testing on acute and convalescent phase paired sera or an assay for the presence of IgM varicella antibody should be used to diagnose recent infection.
4. Test results for IgG/IgM antibodies to VZV should be evaluated in relation to patient symptoms and clinical history to establish a diagnosis.
5. A positive test result in immunocompromised patients receiving prophylactic treatment with zoster immune serum globulin or zoster immune plasma may not be indicative of prior infection with VZV (10).
6. The correlation of negative test results determined with this test as they relate to protection from varicella infection has not yet been definitely established. If a definitive determination of immunity is required, all persons with negative results should be retested by a more sensitive assay, e.g. FAMA (6).
7. Although negative results by FAMA, ACIF, and EIA are reliable indicators of susceptibility to VZV infection, the occurrence of clinical VZV infections in a few individuals with low titers of antibody demonstrable by FAMA and by ACIF indicates that low levels of antibody are not always a reliable indication of protection against clinical illness. However, the detection of serum antibody to VZV has correlated with protection in up to 96% of individuals in one study (1).
8. Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

RANGE OF EXPECTED VALUES

In urban populations of the United States, an estimated 95% of the adult population have antibody indicating prior exposure and immunity to varicella-zoster virus (2). The incidence of reactivation in older adults of a latent varicella infection is 0-2% - 2% (2).

The expected values will vary according to the age of the population, geographical location, and seasonal trends.

The VIDAS Varicella-Zoster IgG assay was evaluated at a clinical reference laboratory using samples from a normal patient population. The positivity rate was 88%.

PERFORMANCE

Sensitivity / specificity:

364 serum samples were tested at two sites - a clinical reference laboratory and bioMérieux Vitek. Each sample was tested using the VIDAS VZG assay, a commercially available IFA, and a commercially available latex agglutination test. In the final interpretation, a sample was considered to be positive if at least two of the three tests were positive. A sample was considered to be negative if at least two of the three tests were negative. A sample was considered unresolved if one of the tests was equivocal and the other two were discrepant. There were three unresolved samples. There were originally 8 (2.2%) equivocal results. Upon retesting, 2 were true positives, 4 were true negatives, and 2 (0.5%) repeated as equivocal. The following table shows the VIDAS Varicella-Zoster IgG results compared to the final interpretation as described above.

Combined sites - VIDAS VZG assay compared to Latex agglutination / IFA / Final interpretation

364 samples tested.

2 equivocal results (excluded from calculations).

3 unresolved samples (excluded from calculations).

		Latex agglutination	
		+	-
VIDAS	+	315	4
VZG	-	6	37

Sensitivity = 98.1%

(95% confidence interval: 95.9% – 99.2%)

Specificity = 90.2%

(95% confidence interval: 76.9% – 97.3%)

Overall agreement: 97.2%

(95% confidence interval: 94.9% – 98.5%)

		IFA *	
		+	-
VIDAS	+	302	2
VZG	-	5	31

Sensitivity = 98.4%

(95% confidence interval: 96.2% – 99.3%)

Specificity = 93.9%

(95% confidence interval: 95%: 79.8% – 99.3%)

Overall agreement: 97.9%

(95% confidence interval: 95.8% – 99.0%)

* 22 samples were IFA "equivocal" IFA – indeterminate due to nonspecific fluorescence.

		VIDAS / LA / IFA	
		+	-
VIDAS	+	317	1
VZG	-	1	40

Sensitivity = 99.7%
(95% confidence interval: 98.2% – 99.9%)

Specificity = 97.6%
(95% confidence interval: 87.1% – 99.9%)

Relative positive predictive value 99.7%
(95% confidence interval: 98.2% – 99.9%)

Relative negative predictive value 97.6%
(95% confidence interval: 87.1% – 99.9%)

Overall agreement: 99.4%
(95% confidence interval: 98.0% – 99.9%)

Additional testing – negative samples

Forty-nine serum samples were obtained from several clinical laboratories. These samples had been tested with immunoenzyme (VZV-IgG) and indirect immunofluorescence (IFA) tests and had been determined to be negative. At bioMérieux Vitek, the samples were tested using the VIDAS VZG assay and two other commercially available tests: one latex agglutination test and an IFA test. Of the 49 samples, 43 were found to be negative. Of these 43 samples, there were 40 VIDAS VZG true negatives, 2 false positives and 1 equivocal (which repeated as equivocal). Three of the remaining 6 samples were positive by both of the commercially available tests. One sample was positive by both the IFA and VIDAS VZG, the 2 others were unresolved (discrepancies between the 3 techniques).

Intra-assay reproducibility:

Intra-assay reproducibility was evaluated at two sites. High positive, low positive, and negative controls were tested «n» times in a single run.

SITE A	n =	Mean RFV	% CV
- High positive	22	1816	7.3 %
- Low positive	22	1073	6.7 %
- Negative	21	107	7.9 %

SITE B	n =	Mean RFV	% CV
- High positive	28	1961	7.0 %
- Low positive	27	926	11.2 %
- Negative	28	195	12.2 %

Inter-assay reproducibility:

Inter-assay reproducibility was evaluated at two sites. At each site, high positive, low positive, and negative controls were each tested in triplicate in a single work list for 9 consecutive days. The mean RFV for each day was used in the calculations.

SITE A	n =	Mean RFV	% CV
- High positive	9	1722	6.4%
- Low positive	9	1052	3.6%
- Negative	9	102	4.2%

SITE B	n =	Mean RFV	% CV
- High positive	9	1657	6.4%
- Low positive	9	975	1.7%
- Negative	9	103	6.2%

CROSS REACTIVITY

35 serum samples which were positive for one or more of the following biological markers: HSV, CMV, EBV, ANA or FR, were tested using the VIDAS Varicella-Zoster IgG. The results indicate no cross-reactivity.

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LITERATURE REFERENCES

- SCHMIDT N., Varicella-zoster Virus. In: Lennette, Balows, Hausler, and Shadomy. Manual of Clinical Microbiology. 4th Ed. Washington, D. C.: American Society for Microbiology. 1985, 720-727.
- STEELE R., Immunology of varicella-zoster virus. In: Nahmias and O'Reilly. Comprehensive Immunology, Vol. 9, Part II. New York: Plenum Publishing Corporation. 1982, 73-87.
- TAKAHASHI M., Herpesviridae: varicella-zoster virus. In: Lennette, Halonen, and Murphy. Laboratory Diagnosis of Infectious Diseases, Principles and Practices, Vol. II. New York: Springer-Verlag. 1988, 261-283.
- BRUNELL, P.A., et al. 1984. Herpesvirus Varicella. In: Ziai. Pediatrics. 3rd Ed. Boston/Toronto: Little, Brown, and Company. 522-523.
- GALLO D., SCHMIDT NJ., Comparison of anticomplement immunofluorescence and fluorescent antibody-to-membrane antigen tests for determination of immunity status to varicella-zoster virus and for serodifferentiation of varicella-zoster and herpes simplex virus infections. Journal of Clinical Microbiology, 1981, 539-543.
- WILLIAMS, V., A. GERSHON, AND P. BRUNELL. 1974. Serologic response to varicella-zoster membrane antigens measured by indirect immuno-fluorescence. J. Inf. Dis. 130:669-672.
- GERSHON A., et al.. Determination of immunity to varicella using an enzyme-linked-immunosorbent-assay. Archives of Virology. 1981, 70: 169-172.
- SHEHAB, Z., AND P. BRUNELL. 1986. Varicella-zoster virus. In: Rose, Freidman, and Fahey. Manual of Clinical Laboratory Immunology. 3rd. Ed. Washington, D.C.: American Society for Microbiology. 502-503.
- LAMPE, A.S., H.J. PEITERSE-BRUIJS, AND J.C.R. EGTERVAN Wissenderke. 1988. Wearing gloves as cause of false negative HIV tests. Lancet ii 1140-1144.
- SCHAUF V., and SALOT R. J. Varicella-zoster virus. In: Belshe. Textbook of Human Virology. Littleton, MA: PSG Publishing. 1991, 842-861.

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REVISION HISTORYChange type categories :

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
<i>Minor typographical, grammar, and formatting changes are not included in the revision history.</i>	

Release date	Part Number	Change Type	Change Summary
2015/01	07968I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS
2015/06	07968J	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE
2016/05	07968K	Technical	CONTENT OF THE KIT (60 TESTS)

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5. GALLO D., SCHMIDT NJ., Comparison of anticomplement immunofluorescence and fluorescent antibody-to-membrane antigen tests for determination of immunity status to varicella-zoster virus and for serodifferentiation of varicella-zoster and herpes simplex virus infections. *Journal of Clinical Microbiology*, 1981, 539-543.
6. WILLIAMS, V., A. GERSHON, AND P. BRUNELL. 1974. Serologic response to varicella-zoster membrane antigens measured by indirect immuno-fluorescence. *J. Inf. Dis.* 130:669-672.
7. GERSHON A., et al.. Determination of immunity to varicella using an enzyme-linked-immunosorbent-assay. *Archives of Virology*, 1981, 70 : 169-172.
8. SHEHAB, Z., AND P. BRUNELL. 1986. Varicella-zoster virus. In: Rose, Freidman, and Fahey. *Manual of Clinical Laboratory Immunology*. 3rd. Ed. Washington, D.C.: American Society for Microbiology. 502-503.
9. LAMPE, A.S., H.J. PEITERSE-BRUIJS, AND J.C.R. EGTERVAN Wissenderke. 1988. Wearing gloves as cause of false negative HIV tests. *Lancet* ii 1140-1144.
10. SCHAUF V., and SALOT R. J. Varicella-zoster virus. In : Belshe. *Textbook of Human Virology*. Littleton, MA : PSG Publishing. 1991, 842-861.

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VIDAS® EBV VCA/EA IgG (VCAG)**IVD**

VIDAS EBV VCA/EA IgG is an automated test for use on the VIDAS family instruments, for the qualitative detection of anti-VCA and anti-EA IgG in human serum using the ELFA technique (Enzyme Linked Fluorescent Assay). Detection of these specific antibodies is an aid in diagnosing infectious mononucleosis (IM).

SUMMARY AND EXPLANATION

The Epstein-Barr virus (EBV, also called Human Herpesvirus 4 (HHV4)) is ubiquitous. It was identified as the cause of infectious mononucleosis (IM). In adolescents and young adults, EBV may cause an infectious mononucleosis still known as Pfeiffer's disease (glandular fever) or the kissing disease. In young children, EBV infection is usually asymptomatic.

Transmission mainly occurs through contact with saliva. Indeed, EBV replication takes place in the oropharyngeal epithelial cells where the virions are released into the saliva by infected B-lymphocytes.

Over 95% of the adult population carries the virus. As this virus generally remains dormant in the body, reactivation may occur, particularly in immunocompromised patients. While reactivated EBV infection is asymptomatic in immunocompetent carriers, it is however associated with clinical disorders and a high morbidity and mortality rate in immunocompromised patients.

Diagnosis of IM is mainly based on clinical symptoms (sore throat, fever and swollen lymph glands). Serology is used to confirm diagnosis of IM and exclude illnesses such as lymphoma and leukemia, which may produce IM-type symptoms. Mononucleosis syndrome may also be caused by other pathogenic agents (cytomegalovirus, HHV6, adenovirus, rubella virus, mumps virus, HIV, hepatitis A virus, influenza A & B viruses and toxoplasma gondii).

Serologic diagnosis of IM includes non-specific tests such as the detection of heterophile antibodies, as well as EBV-specific tests. The latter tests are based on the detection of antibodies produced by the host in response to different antigens produced during the viral cycle. During the lytic phase, EBV early antigens (EA) are produced, then viral capsid antigens (VCA) are expressed at the same time as the viral genome. During the latent cycle, Epstein-Barr nuclear antigens (EBNA) are synthesized.

When IM occurs, heterophile antibodies appear in 60-80% of cases, anti-EA antibodies in 70-80% of cases, anti-VCA IgM antibodies in 100% of cases and anti-VCA IgG antibodies in nearly 100% of cases.

During the convalescent phase, anti-VCA IgG antibodies persist and approximately 95% of patients produce anti-EBNA IgG antibodies (1, 2, 3, 4, 5).

The VIDAS tests enable the detection of anti-EA IgG and anti-VCA IgG antibodies (VIDAS EBV VCA/EA IgG), anti-VCA IgM antibodies (VIDAS EBV VCA IgM) and anti-EBNA IgG antibodies (VIDAS EBV EBNA IgG).

Global interpretation of these assays is useful in diagnosing IM and establishing the stage of infection (see the following table):

VIDAS EBV VCA IgM	VIDAS EBV VCA/EA IgG	VIDAS EBV EBNA IgG	Interpretation of VIDAS EBV (combined results)
-	-	-	Seronegative subject (non infected)
+	-	-	Early-phase IM / primary infection (2)
+	+	-	Acute IM / primary infection
-	+	+	Seropositive subject (past infection)
-	-	+	Isolated EBNA IgG (2)
-	+	-	Isolated VCA/EA IgG (2)
+	+	+	Indeterminate profile (2)

-: Absence of antibodies.

+: Presence of antibodies.

(1) in rare cases, anti-EBNA antibodies are not detected in seropositive subjects who have had a past infection.

(2) to be checked using a specimen collected 1 to 2 weeks later.

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The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strip.

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In the first step, the sample is diluted and then cycled in and out of the SPR. The anti-VCA and/or EA antibodies in the sample will bind to the VCA P18 and EA P54 antigens coated on the interior of the SPR wall. Unbound components are eliminated during the washing steps. During the second step, mouse monoclonal anti-human IgG antibodies in Fab' form, conjugated to alkaline phosphatase, are cycled in and out of the SPR and will attach to any human IgG bound to the antigen. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibody present in the sample. At the end of the assay, the results are automatically calculated by the instrument in relation to the stored standard S1, and then printed out.

CONTENT OF THE KIT (30 TESTS):

30 EBV VCA/EA IgG Strips	STR	Ready-to-use.
30 EBV VCA/EA IgG SPRs 1 x 30	SPR	Ready-to-use. Interior of SPRs coated with VCA P18 and EA P54 antigens.
EBV VCA/EA IgG Positive Control 1 x 1,2 ml (liquid)	C1	Pool of human sera* containing anti-VCA IgG in a phosphate buffer pH 7.4 + 50 g/l BSA + preservatives. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
Negative Control 1 x 1.9 ml (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin + preservatives. MLE data indicate the index: confidence interval ("Control C2 (-) Test Value Range").
Standard EBV VCA/EA IgG 1 x 2 ml (liquid)	S1	Pool of human sera* containing anti-VCA IgG in a phosphate buffer pH 7.4 + 50 g/l BSA + preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Clip Seal		
1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib		

* This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR® is coated during production with VCA P18 and EA P54 antigens. Each SPR is identified by the code VCAG. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening using the clip seal provided.**

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Description of the VCAG strip

Wells	Reagents
1	Sample well.
2	Sample diluent: Phosphate buffered saline + Tween 20 0.05% pH 7.2 + 5 g/l BSA + preservatives (600 µl).
3 - 4 - 5 - 7 - 8	Wash buffer: Phosphate buffered saline + Tween 20 0.25% pH 7.8 + preservatives (600 µl).
6	Conjugate: alkaline phosphatase-labeled mouse monoclonal anti-human IgG antibodies in Phosphate buffered saline pH 6.1 + protein stabilizers + preservatives (400 µl).
9	Sample diluent: Phosphate buffered saline + Tween 20 0.05% + 5 g/l BSA pH 7.2 + preservatives (400 µl).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

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- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPR®s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS EBV VCA/EA IgG kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the desiccant inside using the clip seal provided, and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection

Human serum (glass or plastic plain tube with coagulation activator and tube with separation gel).

It is recommended that each laboratory validate the type of collection tube used.

Sample-related interference

None of the following factors have been found to significantly influence this assay.

- hemolysis (after spiking samples with hemoglobin: 0 to 300 µmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 30 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 610 µmol/l).

However, it is recommended not to use samples that are clearly hemolyzed, lipemic or icteric and, if possible, to collect a new sample.

Specimen preparation:

Plain tubes: wait for samples to coagulate and **centrifuge** according to the tube manufacturer's recommendations to eliminate fibrin.

Other tubes: follow the tube manufacturer's recommendations for use.

Frozen-stored samples: after thawing, these samples must be homogenized before analysis.

Sample stability

Samples can be stored in stoppered tubes for up to 5 days at 2-8°C; if longer storage is required, freeze the sera at -25 ± 6°C.

Do not exceed 3 freeze/thaw cycles.

A study performed on samples frozen for 6 months, showed that the quality of results is not affected.

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Reading VIDAS® Protocole Test Change (PTC) protocol data and MLE data

When using the assay for the first time:

With the external instrument barcode reader,

1. Scan the PTC barcode(s) at the end of the package insert. or downloadable from www.biomerieux.com/techlib. This reading allows VIDAS® PTC protocol data to be transferred to the instrument software for its update.

2. Scan the MLE data on the box label.

Note: If the MLE data have been read before the VIDAS® PTC protocol, read the MLE data again.

When opening a new lot of reagents:

Enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed before initiating the tests, the instrument will not be able to print results.

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Calibration

Calibration, using the standard provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 28 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested in **triplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "VCAG" strip and on "VCAG" SPR® for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "VCAG" code on the instrument. The standard must be identified by "S1" and tested in triplicate. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
4. Mix the standard, controls and samples using a vortex-type mixer (for serum separated from the pellet).

5. For this test, the standard, control, and sample test portion is 100 µl.

6. Insert the "VCAG" SPRs and "VCAG" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.

8. Restopper the vials and return them to 2–8°C after pipetting.
9. The results are obtained within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate biohazard receptacle.

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The RFV obtained for each sample is interpreted by the VIDAS system as follows:

Test value = patient RFV / standard RFV

This test value and the interpreted result are also included on the result sheet. The test value is interpreted as follows:

Test Value (TV)	Results
≤ 0.09	Negative
0.10 ≤ VT ≤ 0.20	Equivocal
≥ 0.21	Positive

IM diagnosis is based on the biological interpretation of test results taking into account the patient's history and the results of the other VIDAS EBV assays (VIDAS EBV EBNA IgG ref. 30 235 and VIDAS EBV VCA IgM ref. 30 237). Please refer to the VIDAS interpretation table on page 1 of this package insert.

It is advisable to check the equivocal results using a second specimen collected one or two weeks later.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS EBV VCA/EA IgG kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

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LIMITATIONS OF THE METHOD

- Biological interpretation for the diagnosis of IM must be based on at least the results of the 3 VIDAS EBV assays (VIDAS EBV EBNA IgG ref. 30 235, VIDAS EBV VCA/EA IgG ref. 30 236 and VIDAS EBV VCA IgM ref. 30 237).
- In the presence of negative results for the 3 VIDAS markers and suspicion of IM, collect a second specimen one or two weeks later and test the two samples simultaneously.
- A positive VCA/EA IgG result (non-differentiated result) does not rule out a recent EBV infection if a negative result is obtained for the VCA IgM and EBNA IgG markers.
- In some cases, at the end of the acute phase of IM (transitional phase), isolated anti-VCA/EA IgG antibodies may be encountered (drop in anti-VCA IgM antibodies and late appearance of anti-EBNA IgG antibodies).
- A negative VIDAS EBV VCA/EA IgG result does not exclude the presence of antibodies directed against VCA P23 protein.
- In rare cases, anti-EBNA antibodies are not detected in seropositive subjects who have a past infection.
- In rare cases, only anti-EBNA antibodies are detected:
- in seropositive subjects who have a past infection,
- in children under 6 months of age (passive transfer of maternal anti-EBNA IgG antibodies),
- in transfused subjects.
- The performance characteristics of this test have not been established for use within the scope of blood donation.
- The performance characteristics of this test have not been established for immunocompromised or immunosuppressed patients (transplanted patients).
- This test has not been validated for use in the diagnosis of nasopharyngeal carcinoma (NPC).
- Certain negative samples may have a test value < 0.
- Cross reactivity or interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history and the results of the two other tests: VIDAS EBV VCA IgM and VIDAS EBV EBNA IgG.

USUAL VALUES

In Western countries, only 5% of adults are never affected by an EBV infection; thus, the antibodies characteristic of a past infection (EBNA and VCA IgG) reach a prevalence rate of 95%.

PERFORMANCE

Performance studies conducted using the VIDAS EBV VCA/EA IgG assay gave the following results:

Precision

Three samples were tested for 10 days in 2 runs per day and 2 replicates per run, with 2 reagent lots on 3 instruments (N = 240).

Mean index of negative sample varies between 0.00 and 0.01.

The **intra-run precision** (repeatability) and the **total precision intra-lot** (intra-run, inter-run, inter-day, inter-calibration, inter-instrument reproducibility) were calculated using this protocol, based on the recommendations of the **NCCLS/CLSI document EP5 A2, volume 24 number 25**.

Lot 1		Repeatability		Inter-run precision		Inter-instrument precision		Total precision	
Sample	Index	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
Sample 1	0.26	0.01	5.2	0.01	3.0	0.01	5.4	0.02	8.7
Sample 2	1.21	0.04	3.7	0.05	3.9	0.03	2.8	0.08	7.0

Lot 2		Repeatability		Inter-run precision		Inter-instrument precision		Total precision	
Sample	Index	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
Sample 1	0.26	0.02	5.8	0.01	2.0	0.01	3.5	0.02	7.5
Sample 2	1.23	0.05	4.2	0.02	1.8	0.00	0.0	0.07	5.8

Sensitivity - Specificity

Specificity and sensitivity performance was established on 2 lots of VIDAS EBV EBNA IgG reagent using 621 samples characterized for their EBV status (past infections, primary infections including 53 early-phase primary infections, seronegatives) and including 93 fresh serum samples. The status of the samples was determined on the basis of the available demographic (patient age) and clinical data, and the results obtained with the routine methods used by the expert site's laboratory.

	Reference populations for sensitivity calculation	Reference population for specificity calculation
VIDAS EBV VCA/EA IgG	Acute primary infections and past infections	Seronegatives

The following consolidated results were obtained:

The 21 samples were found to be equivocal with the VIDAS EBV VCA/EA IgG assay and were not taken into account for the sensitivity and specificity calculations.

The 53 samples characterized as early-phase primary infections were not taken into account for the total sensitivity calculation.

Population	Sensitivity				Specificity
	Total sensitivity**	Early-phase primary infections**	Acute primary infections	Past infections	Seronegatives
Positive	377	25	120	257	0
Negative	14	19	9	5	165
Total*	391	44	129	262	165
%*	96.42	NA	93.02	98.09	100
95% confidence interval	94.02 – 97.88	NA	87.13 – 96.33	95.54 – 99.20	97.63 - 100

• *Equivocal results were not taken into account for the calculations.

** Early-phase primary infections are characterized by the presence of anti-VCA IgM and the absence of anti-VCA IgG with the EIA tests routinely used by the expert site.

NA: Not applicable

Concordance study

Comparison with two other EIA methods

A study was performed at an external site, using 526 samples, to establish the concordance between the VIDAS EBV VCA/EA IgG kit and two other automated EIA methods.

The discordances between the methods were analyzed on the basis of the EBV clinical status of the samples determined by the expert site, taking into consideration the patient's history and serological evolution.

		EA IgG and VCA IgG EIA			Total
		Positive	Equivocal	Negative	
VIDAS EBV VCA/EA IgG	Positive	302	1	15****	318
	Equivocal	12	0	9	21**
	Negative	6***	2	179	187
	Total*	308	3	194	
Concordance % *		95.82%			

*The equivocal results for each of the VIDAS and EIA assays were not taken into account for the calculations.

** Among the 21 samples, 1 is seronegative, 16 are primary infections and 4 are past infections.

*** The 6 samples are primary infections.

**** Among the 15 samples, 10 are primary infections and 5 are past infections.

Comparison with the clinical status of the specimen

A study was performed to establish the concordance between the EBV status defined by the expert site, and the global status obtained using the results of the 3 VIDAS EBV VCA/EA IgG, VIDAS VCA IgM and VIDAS EBNA IgG reagents. The clinical status of EBV and the global interpretations of the 3 VIDAS assays were defined by the expert site taking into consideration the patient's history and serological evolution.

		Reference serological status				
		Seronegatives	Total primary infections		Total past infections	
			Early-phase primary infections	Acute primary infections	Past infections	Isolated VCA IgG
VIDAS Serological status	Seronegatives	159	6	3	0	0
	Primary infections (early-phase and acute)	4	37	120	0	5
	Past infections	0	0	0	224	4
	Isolated VCA/EA IgG	0	5	3	7	12
	Isolated EBNA IgG	2	0	0	5	0
	Indeterminate	1	5	10	8	1
	Total*	165	48	126	236	21
Total concordance %*		92.62%				

*The indeterminate EBV statuses (4.03%) were not taken into account for the calculations.

CROSS REACTIVITY AND INTERFERENCES

Cross-reactivity

The notion of cross-reactivity is the study of samples which are negative for the test to be evaluated and positive for the potentially interfering pathology. The presence of these potentially interfering pathologies must not modify the interpretation of the evaluated test. The equivocal results were not taken into account for the analysis.

The results of the 136 samples tested are presented in the following table:

Pathology	Number of sera tested	Number of cross reactions	Number of equivocal results	Pathology	Number of sera tested	Number of cross reactions
Toxo IgG	10	0	1	HAVT	4	1
CMV IgG	12	0	0	HHV6	8	0
VZVG	30	0	0	HCV	5	0
HSV IgG	7	0	0	HIV	4	0
Rheumatoid factor	12	0	0	HBsT	5	0
Rub IgG	7	0	0	HBcT	5	0
HAMA	3	0	0	HBs Ag	5	0
Parvo B19 IgG	17	0	0	ANA	2	0

Among all the samples tested, only one sample with anti-HAV antibodies showed cross-reactivity with a weakly positive index (index = 0.21).

Interferences

The notion of interference is the study of samples which are positive for the test to be evaluated and positive for the potentially interfering pathology. The presence of these potentially interfering pathologies must not modify the interpretation of the evaluated test.

124 samples of potentially interfering pathologies (Toxo IgG, HBs Ag, Rub IgG, HBcT, HBsT, ANA+, HCV, HIV, CMV IgG, VZV IgG, HAVT, HSV IgG, FR, Parvo B19 IgG, HAMA and HHV6 IgG) with anti-VCA/EA IgG antibodies, were tested.

There was no change in the interpretation of the VIDAS EBV VCA/EA IgG test results.

No interference was observed with the VIDAS EBV VCA/EA IgG test.

WASTE DISPOSAL

Dispose of used or unused reagents, as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

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Release date	Part Number	Change Type	Change Summary
2015/01	14120E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (30 TESTS) WARNINGS AND PRECAUTIONS
2015/06	14120F	Technical	CONTENT OF THE KIT (30 TESTS) INSTRUCTIONS FOR USE
2016/05	14120G	Technical	CONTENT OF THE KIT (30 TESTS)

BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending, and/or registered trademarks belonging to bioMérieux or one of its subsidiaries or one of its companies.

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
VIDAS® EBV VCA IgM (VCAM)

VIDAS EBV VCA IgM is an automated test for use on the VIDAS family instruments, for the qualitative detection of anti-VCA IgM in human serum using the ELFA technique (Enzyme Linked Fluorescent Assay). Detection of these specific antibodies is an aid in diagnosing infectious mononucleosis (IM).

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Over 95 % of the adult population carries the virus. As this virus generally remains dormant in the body, reactivation may occur, particularly in immunocompromised patients. While reactivated EBV infection is asymptomatic in immunocompetent carriers, it is however associated with clinical disorders and a high morbidity and mortality rate in immunocompromised patients.

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The VIDAS tests enable the detection of anti-EA IgG and anti-VCA IgG antibodies (VIDAS EBV VCA/EA IgG), anti-VCA IgM antibodies (VIDAS EBV VCA IgM) and anti-EBNA IgG antibodies (VIDAS EBV EBNA IgG). Global interpretation of these assays is useful in diagnosing IM and establishing the stage of infection (see the following table):

VIDAS EBV VCA IgM	VIDAS EBV VCA/EA IgG	VIDAS EBV EBNA IgG	Interpretation of VIDAS EBV (combined results)
-	-	-	Seronegative subject (non infected)
+	-	-	Early-phase IM / primary infection (2)
+	+	-	Acute IM / primary infection
-	+	+	Seropositive subject (past infection)
-	-	+	Isolated EBNA IgG (2)
-	+	-	Isolated VCA/EA IgG (2)
+	+	+	Indeterminate profile (2)

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After a sample dilution step, the IgM are captured by the monoclonal Ab coating the interior of the SPR. The anti-EBV IgM are specifically detected by the VCA P18 antigen conjugated with alkaline phosphatase.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibody present in the sample. At the end of the assay, an index is automatically calculated by the instrument in relation to the S1 standard stored in memory, and then printed out.

CONTENT OF THE KIT (30 TESTS):

30 EBV VCA IgM Strips	STR	Ready-to-use.
30 EBV VCA IgM SPRs 1 x 30	SPR	Ready-to-use. Interior of SPRs coated with mouse monoclonal anti-human IgM antibody.
EBV VCA IgM Positive Control 1 x 0.6 ml (liquid)	C1	Ready-to-use. Pool of human plasma* containing anti-VCA IgM in a phosphate buffer Ph 7.4 + 50 g/l BSA + preservatives. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin + preservatives. MLE data indicate the index: confidence interval ("Control C2 (-) Test Value Range").
EBV VCA IgM Standard 1 x 1.6 ml (liquid)	S1	Ready-to-use. Pool of human plasma* containing anti-VCA IgM dans in a phosphate buffer pH 7.4 + 50 g/l BSA + preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR

The interior of the SPR® is coated during production with anti-human antibody (mouse). Each SPR is identified by the code VCAM. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

The Strip

The strip consists of 10 wells covered with a labeled foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the VCAM strip

1	Sample well.
2	Sample diluent: Phosphate buffered saline + Tween 20 0.05% pH 7.2 + 5 g/l BSA + preservatives (300 µl).
3 - 4 - 5 - 7 - 8	Wash buffer: TRIS buffered saline + Tween 20 0.25% pH 7.8 + preservatives (600 µl).
6	Conjugate: alkaline phosphatase-labeled VCA P18 antigen in Phosphate buffered saline pH 6.1 + Protein stabilizers + preservatives (400 µl).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

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STORAGE CONDITIONS

- Store the VIDAS EBV VCA IgM kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the desiccant inside and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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- hemolysis (after spiking samples with hemoglobin: 0 to 300 µmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 30 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 607 µmol/l).

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Sample stability

Samples can be stored in stoppered tubes for up to 2 days at 2-8°C; if longer storage is required, freeze the sera at -25 ± 6°C.

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Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "VCAM" strip and "VCAM" SPR® for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "VCAM" code on the instrument. The standard must be identified by "S1" and tested in **triplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
4. Mix the standard, controls and samples using a vortex-type mixer (for serum separated from the pellet).
5. **For this test, the standard, control, and sample test portion is 100 µl.**
5. Insert the "VCAM" SPRs and "VCAM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
6. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
7. Restopper the vials and return them to 2–8°C after pipetting.
8. The results are obtained within approximately **40 minutes**. After the assay is completed, remove the SPRs and strips from the instrument.
9. Dispose of the used SPRs and strips into an appropriate biohazard receptacle.

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IM diagnosis is based on the biological interpretation of test results taking into account the patient's history and the results of the other VIDAS EBV assays (VIDAS EBV EBNA IgG ref. 30 235 and VIDAS EBV VCA/EA IgG ref. 30 236). Please refer to the VIDAS interpretation table on page 1 of this package insert.

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QUALITY CONTROL

One positive control and one negative control are included in each VIDAS EBV VCA IgM kit.

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LIMITATIONS OF THE METHOD

- Biological interpretation for the diagnosis of IM must be based on at least the results of the 3 VIDAS EBV assays (VIDAS EBV EBNA IgG ref. 30 235, VIDAS EBV VCA/EA IgG ref. 30 236 and VIDAS EBV VCA IgM ref. 30 237).
- In the presence of negative results for the 3 VIDAS markers and suspicion of IM, collect a second specimen one or two weeks later and test the two samples simultaneously.
- A positive VCA/EA IgG result (non-differentiated result) does not rule out a recent EBV infection if a negative result is obtained for the VCA IgM and EBNA IgG markers.
- In some cases, at the end of the acute phase of IM (transitional phase), isolated anti-VCA/EA IgG antibodies may be encountered (drop in anti-VCA IgM antibodies and late appearance of anti-EBNA IgG antibodies).
- In rare cases, anti-EBNA antibodies are not detected in seropositive subjects who have a past infection.
- In rare cases, only anti-EBNA antibodies are detected:
 - in seropositive subjects who have a past infection,
 - in children under 6 months of age (passive transfer of maternal anti-EBNA IgG antibodies),
 - in transfused subjects.
- The performance characteristics of this test have not been established for use within the scope of blood donation.
- The performance characteristics of this test have not been established for immunocompromised or immunosuppressed patients (transplanted patients).
- Certain negative samples may have a test value < 0.
- Cross reactivity or interference may be encountered with certain sera containing antibodies directed against reagent components or against other potentially interfering pathologies (particularly with IgM directed against CMV and Parvovirus B19). For this reason, assay results should be interpreted taking into consideration the patient's history and the results of the two other tests: VIDAS EBV VCA/EA IgG and VIDAS EBV EBNA IgG.

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PERFORMANCE

Performance studies conducted using the VIDAS EBV VCA IgM assay gave the following results:

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Sensitivity - Specificity

Specificity and sensitivity performance was established on 2 lots of VIDAS EBV VCA IgM reagent using 621 samples characterized for their EBV status (past infections, early-phase and acute primary infections, seronegatives) and including 93 fresh serum samples. The status of the samples was determined on the basis of the available demographic (patient age) and clinical data, and the results obtained with the routine methods used by the expert site's laboratory.

	Reference population for sensitivity calculation	Reference populations for specificity calculation
VIDAS EBV VCA IgM	Early-phase and acute primary infections	Past infections and seronegatives

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The 17 samples were found to be equivocal with the VIDAS EBV VCA IgM assay and were not taken into account for the sensitivity and specificity calculations.

	Sensitivity			Specificity		
Population	Total sensitivity	Early-phase primary infections**	Acute primary infections	Total specificity	Past infections	Seronegatives
Positive	160	38	122	9	6	3
Negative	21***	13	8	414	252	162
Total*	181	51	130	423	258	165
%*	88.40	74.51	93.85	97.87	97.67	98.18
95% confidence interval	82.79 – 92.35	60.84 – 84.61	88.18 – 96.89	95.96 – 98.89	94.95 – 98.95	94.69 – 99.39

* Equivocal results were not taken into account for the calculations.

**Early-phase primary infections are characterized by the presence of anti-VCA IgM and the absence of anti-VCA IgG with the EIA tests routinely used by the expert site.

*** Among the 21 samples, additional samples were collected from 6 patients: for 4 of them, the second sample led to the detection of anti-VCA IgM antibodies. The other 2 samples, collected 2 to 4 days after the first one, remained negative.

Concordance study

Comparison with another EIA method

A study was performed at an external site, using 526 samples, to establish the concordance between the VIDAS EBV VCA IgM kit and another automated EIA method.

The discordances between the 2 methods were analyzed on the basis of the EBV clinical status of the samples defined by the expert site, taking into consideration the patient's history and serological evolution.

		VCA IgM EIA			Total
		Positive	Equivocal	Negative	
VIDAS EBV VCA IgM	Positive	153	10	2****	165
	Equivocal	3	6	5	14**
	Negative	13***	18	316	347
	Total*	166	28	318	
Concordance % *		96.90%			

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** Among the 14 samples, 1 is seronegative, 8 are primary infections and 5 are past infections.

*** Among the 13 samples, 2 are seronegatives, 9 are primary infections and 2 are past infections.

**** Among the 2 samples, 1 is seronegative and 1 is a primary infection.

Comparison with the clinical status of the sample

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		Reference serological status				
		Seronegatives	Total primary infections		Total past infections	
			Early-phase primary infections	Acute primary infections	Past infections	Isolated VCA IgG
VIDAS serological status	Seronegatives	159	6	3	0	0
	Primary infections (early-phase and acute)	4	37	120	0	5
	Past infections	0	0	0	224	4
	Isolated VCA/EA IgG	0	5	3	7	12
	Isolated EBNA IgG	2	0	0	5	0
	Indeterminate	1	5	10	8	1
	Total*	165	48	126	236	21
Total concordance %*		92.62%				

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The results of the 181 samples tested are presented in the following table:

Pathology	Number of sera tested	Number of cross reactions	Number of equivocal results
Toxo IgM	11	1	3
CMV IgM	8	2	3
VZV IgM	10	0	1
HSV IgM	18	1	0
Rheumatoid factor	12	1	0
Rub IgM	7	1	1
Parvo B19 IgM	11	0	1
Anti Pal	10	1	1

Pathology	Number of sera tested	Number of cross reactions	Number of equivocal results
HAV IgM	9	2	1
HHV6 IgM	9	0	0
HCV	12	1	0
HIV	8	0	0
HBc IgM	14	0	0
HBs Ag	23	1	0
ANA	19	0	1

Cross-reactivity may be encountered with certain sera containing antibodies directed against reagent components or against other potentially interfering pathologies. For this reason, assay results should be interpreted taking into consideration the patient's history and the results of the two other tests: VIDAS EBV VCA/EA IgG and VIDAS EBV EBNA IgG.

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115 samples of potentially interfering pathologies (Toxo IgM, HBs Ag, Rub IgM, HBc IgM, ANA+, HCV, HIV, CMV IgM, VZV IgM, HAV IgM, HSV IgM, FR, Parvo B19 IgM, anti-Pal and HHV6 IgM) with anti-VCA IgM antibodies, were tested. There was no change in the interpretation of the VIDAS EBV VCA IgM test results.

No interference was observed with the VIDAS EBV VCA IgM test.

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Release date	Part Number	Change Type	Change Summary
2015/01	14125G	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (30 TESTS) WARNINGS AND PRECAUTIONS
2015/06	14125H	Technical	CONTENT OF THE KIT (30 TESTS) WARNINGS AND PRECAUTIONS
2016/05	14125I	Technical	CONTENT OF THE KIT (30 TESTS)

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VIDAS® HBe/Anti-HBe (HBE/HBET)

IVD

VIDAS HBe/Anti-HBe is an automated qualitative test for use on the VIDAS family instruments, for the detection of the hepatitis B e-antigen (HBe Ag) or antibodies (anti-HBe) in human serum or plasma (lithium heparin, sodium citrate or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Approximately 5% of the world population is infected by the hepatitis B virus (HBV) which causes a necroinflammatory liver disease of variable duration and severity. Chronically infected patients with active liver disease carry a high risk of developing cirrhosis or hepatocellular carcinoma. The immune response to HBV-encoded antigens is responsible both for viral clearance and for disease pathogenesis during this infection. Hepatitis B can be transmitted through sexual contact, exposure to blood products, or perinatally.

Perinatal transmission can be as high as 90% in women who are chronically infected with HBV, in highly endemic areas or in regions with no systematic testing of pregnant women. The child becomes a chronic carrier of HBs antigen in 90% of cases (1).

The C gene of the HBV viral genome can express two distinct proteins: 1) core protein (HBc Ag) which forms the nucleocapsid, 2) e non-particulate protein (HBe Ag). HBe Ag can be detected in the serum of patients with hepatitis B wild-type virus during active viral replication (2). The function of this antigen in the viral replication cycle is not clearly defined. It is not indispensable to the virus, however this antigen is a major immunological target in the clearance of the virus.

Generally, HBe Ag can be detected early in an acute HBV infection. It coincides with or follows the appearance of HBs Ag. In acute cases which evolve to recovery, HBe Ag disappears after several weeks and seroconversion to anti-HBe usually follows. For chronic HBV cases, HBe Ag can persist from several months to as long as several years, indicating on-going viral replication (2). The presence of HBe Ag in serum indicates active replication of HBV and people with HBe Ag-positive results are considered highly infectious for hepatitis B (3).

HBe Ag is used to monitor chronic hepatitis and anti-viral therapy. The objective of anti-viral therapy is to prevent the progression to liver cirrhosis. HBe Ag seroconversion to anti-HBe is generally considered as an indicator of transition to a state of viral latency accompanied by normalization of aminotransferase levels. HBe Ag seroconversion to Anti-HBe reduces the risk of developing cirrhosis and decompensated liver disease (4).

A positive result for anti-HBe in patients recovering from acute hepatitis indicates normal recovery, particularly if HBsAg and HBe Ag are no longer detectable. In an HBV carrier, a positive anti-HBe result usually indicates inactivity of the virus and low infectivity of the patient (5). However a positive anti-HBe result in the presence of a positive HBV-DNA test result can indicate active viral replication and progression of liver disease in a carrier (5).

HBV mutants, unable to secrete HBe Ag, can prevail over wild-type HBVs in patients with severe acute and chronic hepatitis B and in chronic HBs Ag carriers at the time of HBe Ag/anti-HBe seroconversion (6).

The VIDAS HBe/anti-HBe assay aids in detecting the presence of HBe antigen or anti-HBe antibody which are respectively, except for infections by mutant HBe viruses, markers for a viral replication phase or an evolution towards normalization.

PRINCIPLE

HBe Ag

The principle of the assay combines an enzyme immunoassay method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®), serves as the solid phase as well as the pipetting device. The reagents for the assay, with the exception of the standard and controls, are ready-to-use and pre-dispensed in the sealed reagent strips.

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After dilution in the instrument, the sample is cycled in and out of the SPR. During this time, the HBe antigen, if present in the sample, will bind simultaneously to the specific monoclonal antibody fixed to the SPR and to another monoclonal specific antibody conjugated with biotin. Unbound components are removed by washing steps. The presence of biotin is detected by incubation with streptavidin conjugated with alkaline phosphatase. Successive washes remove unbound components. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample. At the end of the assay, results are automatically calculated by the instrument and then printed out.

anti-HBe antibodies

The test for anti-HBe uses the same protocol outlined above with the exception of a pre-assay sample incubation with added recombinant HBeAg. During this incubation, the recombinant HBe Ag is neutralized by antibodies of the sample, if present. The HBe Ag assay is then performed, to detect the residual free HBe Ag. The presence of antibodies in the sample will decrease the signal. If anti-HBe is not present in the sample, the added recombinant HBe Ag will be completely detected. The specific interpretation of the results for anti-HBe will be done automatically by the instrument, and then printed out.

CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS:

30 HBE strips	STR	Ready-to-use.
30 HBE SPRs 1 x 30	SPR®	Ready-to-use. Interior of SPRs coated with mouse monoclonal anti-HBe Ab.
HBe Ag positive control 1 x 1.5 ml (liquid)	C1	Ready-to-use. Stable protein base overloaded with recombinant HBe Ag + 0.9 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1 (+)Test Value Range").
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use negative control for both the HBe Ag and Anti-HBe assays. Phosphate buffer + protein stabilizer of animal origin + preservatives.
Anti-HBe positive control 1 x 1.5 ml (liquid)	C3	Ready-to-use. Stable protein base overloaded with mouse monoclonal anti-HBe Ab + 0.9 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C3 (+)Test Value Range").
HBe Ag standard 4 x 1 ml (lyophilized)	S1	Stable protein base overloaded with recombinant HBe Ag + protein stabilizers. Dilute the contents of a vial with 1 ml of standard diluent to reconstitute the standard. Store at 2-8°C after reconstitution for up to 6 months.
Anti-HBe standard 1 x 2 ml (liquid)	S2	Ready-to-use. Delipidated human serum* + 0.9 g/l sodium azide.
S1 Standard diluent 1 x 5 ml (liquid)	R1	Ready-to-use. Contains 0.9 g/l sodium azide
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR

The interior of the SPR is coated during production with monoclonal anti-HBe antibody. Each SPR is identified by the HBE code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

The strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The first well in the strip is for the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the HBe/Anti-HBe strip

Wells	Reagents
1	Sample well
2	Conjugate: biotin labeled monoclonal anti-HBe Ab (mouse) + 0.9 g/l sodium azide (300 µl).
3 - 4 - 6 - 7 - 8 - 9	Wash buffer: Buffered saline (0.05 mol/l) (pH 7.8) with 0.9 g/l sodium azide (600 µl).
5	Tracer: Alkaline Phosphatase-labeled streptavidin + 0.9 g/l sodium azide (400 µl).
10	Optical Cuvette with Substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine DEA* (0.62 mol/l or 6.6 %) pH 9.2 + 1 g/l sodium azide (300 µl).

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For further information, refer to the Material Safety Data Sheet

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- Pipette with disposable tip to dispense 100µl, 150µl, 1 ml.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use SPRs if the pouch is pierced.
- Do not use visibly deteriorated strips (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless gloves** as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS HBe/Anti-HBe kit at 2-8°C .

- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C, including S1 standard after reconstitution.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain the stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

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The use of heat inactivated sera has not been validated for this test. Do not heat samples.

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- hemolysis (after spiking samples with hemoglobin, from 0 to 5.6 mg/ml (monomer)),
- lipemia (after spiking samples with lipids, from 0 to 10 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin, from 0 to 490 µmol/l).

However, it is recommended not to use samples that are hemolyzed or lipemic and, if possible, to collect a new sample.

Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 1 week. If longer storage is required, freeze the sera or plasma at -25 ± 6°C.

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Calibration

Calibration, using the standards provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration data and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standards, identified by S1 (Ag HBe) and S2 (Anti-HBe) must be tested in duplicate (see User's Manual). The standard value must be within the set RFV ("Relative Fluorescence Value") and the coefficient of variation of the standard run in duplicate must be less than the normal value indicated on the MLE card. If this is not the case, recalibrate.

HBe assay procedure

1. Only remove the required reagents from the refrigerator and allow to come to room temperature for 30 minutes before use.
2. Remove one "HBE" strip and one "HBE" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBE" code on the instrument. The calibrator, which must be identified by "S1", should be tested **in duplicate** and be placed at the beginning of the run. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
4. Mix the standard, control and samples using a vortex type mixer (for serum or plasma separated from the pellet).
5. **For this test, the standard, control, and sample test portion is 150 µl.**
6. Insert the "HBE" SPRs and the "HBE" Reagent Strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Restopper the vials and return them to 2–8°C after pipetting.
9. The assay will be completed within approximately 90 minutes. After the assay is completed, remove the SPRs and the strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION OF HBe TEST

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent strip's reading cuvette for each test. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The index value is calculated by dividing the sample or control RFV by the standard RFV:

$$i = \text{index} = \text{sample RFV} / \text{S1 standard RFV}$$

This index and its interpretation appear on the result sheet:

Index	Interpretation
$i < 0.1$	Negative: absence of HBe Ag
$i \geq 0.1$	Positive: presence of HBe Ag

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As no international standard is available for the determination of HBe Ag, the VIDAS reagent is calibrated against collection sera.

Anti-HBe assay procedure

Note: The VIDAS Anti-HBe assay begins with a preliminary incubation step which can be performed in a water bath, dry incubator, or directly in the instrument.

Caution: Each time a new vial of S1 is reconstituted for the HBET test, this operation must be validated by testing the C2 and C3 controls.

Preliminary incubation

1. Mix the S1 and S2 standards, the C2 and C3 controls and the samples using a vortex type mixer.
2. For each sample standard (S2) and control to be tested: pipette 100 µl of HBe (S1) standard into a glass or plastic tube or into the sample well of the VIDAS "HBE" strip. Add 100 µl of sample, standard (S2) or control. Cover and mix the tubes using a vortex type mixer or mix by pipetting several times in the sample well.

Incubate at 37 ± 2°C for 1 hour ± 5 minutes in a water bath, dry incubator if tubes are used, or in the instrument if preliminary incubation is performed in the strips.

Testing with the instrument

3. Remove the required reagents from the refrigerator 30 minutes before the end of preliminary incubation to allow them to come to room temperature.
4. The test is identified by the "HBET" code on the instrument. The standard, must be identified by "S2", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C3". If the negative control needs to be tested, it should be identified by "C2".
5. a) If preliminary incubation is not performed in the strip, pipette 150 µl of the mixture (prepared during preliminary incubation) into the sample well of the "HBE" strip (NB: samples and controls will be tested singly and the standard in duplicate).

Insert the strips and the SPRs into the instrument.

5. b) **If preliminary incubation has been performed in the instrument, remember to insert the SPRs.**
6. Check to make sure the color labels with the assay code on the "HBE" SPRs and the "HBE" Reagent strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Restopper the vials and return them to 2–8°C after pipetting.
9. The assay will be completed within approx. 90 minutes. After the analysis is completed, remove the SPRs and the strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

Anti-HBe RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the strip's cuvette for each test. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The index value is calculated by dividing the sample or control RFV by the standard RFV:

$$i = \text{index} = \text{sample RFV} / \text{standard S2 RFV}$$

This index and its interpretation appear on the result sheet: Interpretation of results according to the index is as follows:

Index	Interpretation
$i < 0.4$	Positive: presence of anti-HBe
$0.4 \leq i < 0.5$	Equivocal
$i \geq 0.5$	Negative: absence of anti-HBe

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As no international standard is available for the determination of anti-HBe Ab, the VIDAS reagent is calibrated against collection sera.

QUALITY CONTROL

One positive Ag HBe control (C1) and one positive Anti-HBe control (C3) are included in each VIDAS HBe/Anti-HBe kit as well as one negative HBe/Anti-HBe control (C2) which can be used for both tests.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. Results cannot be validated if any of the control values deviate from the expected values.

Caution: any reconstitution of a new S1 vial used for the anti-HBe test must be validated using C2 and C3 controls.

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LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies directed against reagent

components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

HBe PERFORMANCE

Precision

Intra-assay, intra-instrument reproducibility

Five samples were tested singly in 12 runs on the same instrument over a 7-week period (calibration was performed every 14 days as described in the User's Manual).

Sample	Mean index	% CV
1	0.01	20.4
2	0.01	15
3	0.22	5.1
4	0.80	4.2
5	2.12	5.5

Inter-assay reproducibility (inter-instrument)

A panel of 3 samples was tested singly in 3 runs at 4 different sites.

Sample	Mean index	% CV
Negative	0.00	-
Weak positive	0.29	12.5
Strong positive	1.61	8.6

Analytical sensitivity

The sensitivity, determined using the PEI standard, was estimated at 0.25 PEI U/ml.

Sensitivity - Specificity

1) Random population

413 samples from blood donors were tested in parallel using VIDAS and an enzyme immunoassay technique (EIA).

		EIA	
		positive	negative
VIDAS	positive	0	0
HBe Ag	negative	3*	410

* These samples were found to be negative with a second EIA technique used to resolve discrepant results. These samples were also found to be negative for HBV-DNA.

Relative specificity after confirmation: 100%
(95% Confidence interval: 99.04% - 100%).

2) Routine diagnostic population

368 samples were tested at two different clinical sites using VIDAS HBe and an enzyme immunoassay technique. Five discrepant results (negative with VIDAS / positive with EIA) were analyzed using a second EIA technique and HBV-DNA detection. The final results are as follows:

		Final interpretation	
		positive	negative
VIDAS	positive	204	0
Ag HBe	negative	4	160

Relative sensitivity after confirmation: 98.1%

(95% Confidence interval: 95.1% - 99.3%).
Relative specificity after confirmation: 100%
(95% Confidence interval: 97.6% - 100%).

3) Study of 203 clinically negative samples

203 samples were tested using VIDAS HBe and an enzyme immunoassay technique.

Discrepant results were resolved using another EIA technique.

		Final interpretation	
		positive	negative
VIDAS	positive	1	1
Ag HBe	negative	0	201

Relative specificity after confirmation: 99.5%
(95% Confidence interval: 97.2% - 99.9%).

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69 potentially interfering samples were tested.

	VIDAS	EIA 1
	negative	negative
Antinuclear antibodies	9	9
Anti-EBV antibodies	3	3
Anti-HCV antibodies	18	18
Anti-HAV IgM	6	6
Rheumatoid factor	33	33

Anti-HBe PERFORMANCE DATA

Precision

Precision was evaluated using the positive control and the negative control tested in duplicate twice a day for eight days at three sites.

Sample	negative control C2	Positive control C3
Mean index	0.99	0.20
Intra-assay reproducibility	3.4%	6.6%
Inter-assay reproducibility	10.8%	8.2%

Analytical sensitivity

The sensitivity, determined using the PEI standard, was estimated at 0.03 PEI U/ml.

Sensitivity - Specificity

1) Routine diagnostic population

900 acute hepatitis, chronic hepatitis, and potentially interfering samples were tested at two different sites. The VIDAS results were compared to a commercially available EIA for anti-HBe. Equivocal samples were not included in the assay performance calculation. Discordant results were resolved with a second commercially available EIA for anti-HBe and taking into account the patient's history.

		Confirmed results	
		positive	negative
VIDAS	positive	292	3
Anti-HBe	negative	4	581

Relative sensitivity after confirmation 98.65%
(95% Confidence interval: 96,52% - 99.48%)
Relative specificity after confirmation: 99.49%
(95% Confidence interval: 98,47%-99.83%)

2) Study of 210 blood donor samples

210 samples were tested using VIDAS anti-HBe and an enzyme immunoassay technique.

Discrepant results were resolved using another EIA technique.

		positive	negative
VIDAS	positive	0	1
Anti- HBe	negative	0	209

Relative specificity after confirmation: 99.5%
(95% Confidence interval: 97.3% - 99.9%).

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58 potentially interfering samples were tested.

	VIDAS	EIA
	negative	negative
Anti-HBs antibodies	10	10
Antinuclear antibodies	10	10
Anti-HCV antibodies	8	8
Anti-HAV IgM	9	9
Rheumatoid factor	11	11
Anti-EBV antibodies	10	10

EXPECTED VALUES

The incidence of hepatitis B cases in Europe is approximately 20/100,000, ranging from 1/100,000 in Scandinavian countries to 60/100,000 in Central Europe. In Europe, the endemia increases North to South and West to East. In South-East Asia, in China, in Sub-Saharan Africa or in South America the prevalence can exceed 10%.

The HBe antigen only exists in HBs Ag positive subjects, its presence is an unfavorable prognosis element as it is the marker of an active viral multiplication.

The anti-HBe antibody is a favorable prognosis element, especially if it appears early on.

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INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

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Release date	Part Number	Change Type	Change Summary
2015/06	09078K	Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2016/04	09078L	Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS

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VIDAS® HAV IgM (HAVM)

VIDAS HAV IgM is an automated qualitative test for use on the VIDAS family instruments for the detection of IgM directed against the hepatitis A virus (HAV) after immunocapture, in human serum or plasma (heparin or EDTA), using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Hepatitis A is caused by the hepatitis A virus, belonging to the picornavirus family. It is transmitted by fecal-oral contact and was therefore considered to be a generally asymptomatic disease, occurring mainly in children. The fall in prevalence observed over the past few years has led to an increase in symptomatic cases in older patients (1). Fulminant hepatitis, which is often fatal, remains rare (0.1 to 0.2% of cases) (2).

Direct diagnosis is difficult since virus excretion in stools is early and generally ends before clinical signs appear (3). The viremia stage is short-lived, and the detection of viral RNA of little significance (4). Furthermore, the early and rapid development of IgG makes it difficult to detect the time of seroconversion by titrating the IgG. This explains why the only reliable method is the detection of anti-HAV IgM, especially after immunocapture. IgM can be detected as soon as the first symptoms appear, and generally persist for 2 to 4 months. In rare cases, residual IgM have been detected 6 to 12 months after the start of the infection (5).

In cases of active prevention, anti-HAV IgM can be detected in patient sera within two weeks following the first injection of the vaccine (6).

Detection of anti-HAV IgM aids in the diagnosis of acute hepatitis, but is not sufficient to validate post-vaccine immunity.

PRINCIPLE

The assay principle combines a 2-step enzyme immunoassay after immunocapture, with a final fluorescent detection (ELFA).

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After diluting the sample, the serum IgM bind with the anti-μ chain polyclonal antibodies coating the interior of the SPR. Unbound sample components are removed by washing.

The anti-HAV IgM are specifically detected by an immune complex, formed of inactivated viral antigen and an alkaline phosphatase-labeled anti-HAV monoclonal antibody. The unbound immune complex is removed by washing.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of anti-HAV IgM present in the sample. At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (30 TESTS) - RECONSTITUTION OF REAGENTS:

30 HAVM strips	STR	Ready-to-use.
30 HAVM SPRs 1 x 30	SPR	Ready-to-use. Interior of SPRs coated with anti-human μ chain polyclonal antibodies (goat).
HAVM positive control 1 x 1 ml (liquid)	C1	Ready-to-use. Human* serum with anti-HAV IgM + protein stabilizer + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1 (+)Test Value Range).
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives.
HAVM standard 1 x 1.2 ml (liquid)	S1	Ready-to-use. Human* serum with anti-HAV IgM + protein stabilizer + 1 g/l sodium azide.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

*This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR®

The SPR is coated during production with anti-human μ chain polyclonal antibody purified by affinity chromatography. Each SPR is identified by the HAVM code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the HAVM strip

Wells	Reagents
1	Sample well.
2	Sample diluent: TRIS buffer (0.05 mol/l) pH 7.4 + Tween + protein and chemical stabilizers + 1 g/l sodium azide (400 μ l).
3 - 4 - 5 - 8 - 9	Washing buffer: TRIS buffer (0.05 mol/l) pH 7.4 + protein and chemical stabilizers + 1 g/l sodium azide (600 μ l).
6	Inactivated hepatitis A antigen + protein and chemical stabilizers + 1 g/l sodium azide (300 μ l).
7	Anti-hepatitis A monoclonal antibody (mouse) conjugated with alkaline phosphatase + 1 g/l sodium azide (300 μ l).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 μ l).

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- Pipette with disposable tip to dispense 100 μ l.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS® family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.

- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS HAV IgM kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection

Serum or plasma (lithium heparinate or EDTA).
It is recommended that each laboratory checks the compatibility of collection tubes used.

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- hemolysis (after spiking samples with hemoglobin: 0 to 270 µmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 2 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples (naturally icteric plasma) with bilirubin: 0 to 500 µmol/l).

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Do not inactivate samples.

Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C.
Freeze once only.

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Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HAVM" strip and one "HAVM" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HAVM" code on the instrument. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
4. If necessary, clarify samples by centrifugation.
5. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
6. **For this test, the standard, control, and sample test portion is 100 µl.**
7. Insert the "HAVM" SPRs and "HAVM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
9. Restopper the vials and return them to 2-8°C after pipetting.
10. The assay will be completed within approximately 60 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and reagent strips into an appropriate recipient.

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$$i = \text{Test Value} = \text{Patient RFV} / \text{Standard RFV}$$

This Test Value and the result interpretation are both mentioned on the report.

The Test Value is interpreted as follows:

Test Value	Interpretation
$i < 0.4$	Negative
$i \geq .4$ and $i < 0.5$	Equivocal**
$i \geq 0.5$	Positive

** It is advisable to control equivocal results by performing a new test using a second sample.

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QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HAV IgM kit.

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PERFORMANCE

Studies performed using VIDAS HAVM gave the following results:

Precision

Within-run reproducibility:

3 samples were tested 15 times in a same run.

	Number	Mean RFV	CV (%)
Negative	15	74.5	4.05
Weak positive	15	994	1.79
Strong positive	15	2627	1.8

Between-run reproducibility:

3 samples were tested in triplicate in 3 different runs on 3 sites.

	Number	Mean Index	CV (%)
Negative	27	0.05	7.63
Weak positive	27	1.01	2.07
Strong positive	27	3.45	2.4

Specificity

Out of the 544 blood donors who were found to be negative with another EIA technique; 538 were found to be negative and 3 were found to be positive with VIDAS HAV IgM, giving a **relative specificity of 99.44%***.

*Equivocal results were not included in the specificity calculation.

Study of clinical trials

The sensitivity was tested on sera corresponding to samples from patients with clinically and biologically documented hepatitis A or from patients at risk with anti-HAV IgM, detectable using another EIA technique.

1) Acute hepatitis A or positive anti-HAV IgM antibody in patients presenting either an increase in transaminase levels or a risk factor (travel, contact with infection etc.):

Of the 205 sera tested:

- 201 were found to be positive with both VIDAS HAV IgM and the comparative EIA technique,
- 4 were found to be equivocal.

2) Monitoring of 8 clinically and biologically documented cases of acute hepatitis A:

Of the 55 sera tested:

- 23 were found to be negative with VIDAS HAV IgM. Discrepancies were found for six sera with the EIA technique which found them positive.
- 30 were found to be positive with both VIDAS HAV IgM and the comparative EIA technique,
- 2 were equivocal.

The 6 discrepant sera were from two patients who were being monitored for hepatitis A. They probably correspond to residual IgM which were later detected using the comparative EIA technique.

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* The serum that was found to be positive with VIDAS HAV IgM was confirmed to be positive with the other EIA technique, used as reference.

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Release date	Part Number	Change Type	Change Summary
2015/06	07206Q	Technical	CONTENT OF THE KIT (30 TESTS) - RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2016/05	07206R	Technical	CONTENT OF THE KIT (30 TESTS) - RECONSTITUTION OF REAGENTS

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VIDAS[®] Anti-HCV (HCV)

VIDAS Anti-HCV is an automated qualitative test for use on the instruments of the VIDAS family, for the detection of IgG antibodies to hepatitis C virus (anti-HCV) in human serum or plasma (heparin) using the ELFA technique (Enzyme Linked Fluorescent Assay). The detection of these specific antibodies, in conjunction with other clinical information, aids in the diagnosis of infection in persons with symptoms of hepatitis and in persons at risk for hepatitis C infection.

SUMMARY AND EXPLANATION

The Hepatitis C virus (HCV) discovered in 1989 using advanced molecular biology techniques, was rapidly found to account for the majority of those patients with non-A non-B hepatitis. HCV represents a major worldwide public health problem requiring global action for the diagnosis, treatment and prevention of this infection (1).

HCV is primarily parenterally transmitted through direct blood-to-blood contact between two people: use of unsterilized injection devices and transfusion of unscreened blood or blood products (2). The disease frequently progresses to chronic hepatitis C (80%), exposing patients to a greater risk of hepatic complications such as cirrhosis or hepatocellular carcinoma. (3).

The current standard of treatment for HCV is a combination of two drugs: pegylated interferon and ribavirin, but due to the high genetic variability of HCV (4), it is still only partially effective: viral eradication in less than 50% of patients infected with genotype 1 hepatitis C virus against approximately 80% of patients infected with genotype 2 or 3. New therapeutic options are under study to offer more effective and safer personalized treatments (5,6).

Diagnosis of patients infected with HCV can be performed using two categories of virological tests: indirect tests, and direct tests (7). Indirect serological tests are third-generation enzyme immunoassays that detect antibodies to HCV. The antigens used in the tests to detect antibodies are from the structural and non-structural regions of the HCV (8) (capsid, protein, cofactors, polymerase, etc.). The presence of anti-HCV antibodies indicates that an individual may have been infected with HCV in the past or may have an ongoing HCV infection. To confirm the presence of active HCV infection, a positive serological test can be completed using direct tests (e.g.: molecular assays that detect RNA genomes). The results will be used to guide patient management and determine the optimal duration of treatment.

The VIDAS Anti-HCV assay is a third-generation test using antigens corresponding to the HCV core, NS3 and NS4 proteins for the qualitative detection of anti-HCV antibodies.

PRINCIPLE

The assay principle combines a two-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR[®]) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and are pre-dispensed in the sealed reagent strips.

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During the first step, the sample is diluted, and then cycled in and out of the SPR several times. The anti-HCV antibodies present in the sample will bind to the antigens representing the HCV core, NS3 and NS4 proteins coated on the interior of the SPR. Unbound sample components are washed away.

During the second step, mouse monoclonal anti-human IgG antibodies in Fab form, conjugated to recombinant alkaline phosphatase (yeast) are cycled in and out of the SPR several times and will bind to the human Ig bound to the molecules on the solid phase. Further wash steps remove unbound components.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibody present in the sample. At the end of the assay, the results are automatically calculated by the instrument in relation to the Standard S1 stored in memory, and then printed out.

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60 Strips HCV	STR	Ready-to-use.
60 SPRs HCV (2 x 30)	SPR®	Ready-to-use. Interior of SPRs coated with antigens representing the HCV core, NS3 and NS4 proteins.
HCV Positive control 1 x 1.9 ml (liquid)	C 1	Pooled human serum or plasma* containing anti-HCV IgG in a phosphate buffer + BSA + preservatives Index: The confidence interval is indicated on the MLE card after the following mention: "Control C1 (+) Test Value Range".
HCV Negative control 1 x 1.9 ml (liquid)	C 2	Phosphate buffer + BSA + preservatives Index: The confidence interval is indicated on the MLE card after the following mention: "Control C2 (-) Test Value Range".
Standard 1 x 1.9 ml (liquid)	S 1	Pooled human serum or plasma* containing anti-HCV IgG in a phosphate buffer + BSA + preservatives The confidence interval in "Relative Fluorescence Value (RFV)" is indicated on the MLE card after the following mention: "Standard (S1) RFV Range".
1 MLE card (Master Lot Entry)		Specifications for the factory master data required to calibrate the test: to read the MLE data, please refer to the User's Manual
1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib		

* This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1 and HIV2. The product has been inactivated. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with the antigens representing the HCV core, NS3 and NS4 proteins. Each SPR is identified by the code "HCV". Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening.**

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Description of the HCV strip

Wells	Reagents
1	Sample well.
2	Sample diluent: TRIS buffered saline + Tween 20 + BSA + preservatives (600 µl)
3 – 4 – 5 – 7 – 8	Wash buffer: TRIS buffered saline + Tween 20 + preservatives (600 µl)
6	Conjugate: mouse monoclonal anti-human IgG antibodies conjugated to recombinant ALP in Phosphate buffered saline + protein stabilizer + preservatives (400 µl)
9	Reactive diluent: Phosphate buffered saline + preservative (400 µl)
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 µl).

*** IRRITANT reagent:**

- **R 36:** Irritating to eyes.
- **S 26:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

For further information, refer to the Safety Data Sheet available on request.

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- Pipette with disposable tip to dispense 100 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPR®s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the risk phrase "R" and the precautions "S" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS Anti-HCV kit at 2-8°C.
- Do not freeze reagents.
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the desiccant inside and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS**Specimen type and collection**

Human serum or plasma

Types of tubes validated:

- Plain tube,
 - Tube with lithium heparin,
 - Tube with sodium heparin,
 - Tube with lithium heparin and separation gel,
 - Plastic tube with clot activator,
 - Plastic tube with clot activator and separation gel.
- The use of heat-inactivated sera has not been validate.

Note: Depending on the manufacturer, blood collection tubes may contain materials and additives that could generate different test results.

It is the responsibility of each laboratory to validate use of these tubes in accordance with the manufacturer's recommendations for use.

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- hemolysis (after spiking samples with hemoglobin: 0 to 300 µmol/l (monomer),
- lipemia (after spiking samples with lipids: 0 to 30 mmol/l equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 220 mg/ml or 376 µmol/l).

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Serum and plasma samples separated from the clot, can be stored at 2-8 °C in stoppered tubes for 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C.

Do not exceed 3 freeze/thaw cycles.

A study performed on samples frozen for 12 months, showed that the quality of results is not affected.

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VIDAS PTC protocol data entry

When using the assay for the first time, **and before reading the MLE data**, scan the bar code(s) (at the end of the package insert) using the instrument's external bar code reader. This reading will allow VIDAS PTC protocol data to be transferred to the instrument software for its update. These data should only be read the first time the assay is used.

Master lot data entry

Note: When using the assay for the first time, enter the VIDAS PTC protocol (bar codes at the end of the package insert) before reading the MLE data. If the MLE data have been read before the VIDAS PTC protocol, read the MLE data again.

Before each new lot of reagents is used, specifications (or factory master calibration data) must be entered into the instrument using the MLE data. If this operation is not performed **before initiating the tests**, the instrument will not be able to print results. The master lot data need only be entered once for each lot.

It is possible to enter MLE data manually or automatically depending on the instrument (refer to the User's Manual).

Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 28 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by "S1", must be tested **in duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator. They can be used immediately.**
2. Use one "HCV" strip and one "HCV" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HCV" code on the instrument. The standard must be identified by "S1", and tested in duplicate. If the positive control is to be tested, it should be identified by "C1". If the negative control is to be tested, it should be identified by "C2".
4. If necessary, clarify samples by centrifugation.
5. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
6. **For this test, the standard, control, and sample test portion is 100 µl.**
7. Insert the "HCV" SPRs and the "HCV" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.

9. Restopper the vials and return them to 2-8°C after pipetting.
10. The assay will be completed within approximately 40 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

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The results are automatically calculated by the instrument.

The patient RFV is interpreted as follows:

Test value = (patient RFV / standard RFV).

This test value and the interpretation are also included on the result sheet. The interpretation depending on the test value is as follows:

Test value (TV)	Interpretation
<1.00	negative
≥1.00	positive

All positive patient results must be verified in duplicate. If at least one of the repeat values is positive, the patient result is considered as positive. In that case, additional tests should be performed (another immunoassay or HCV marker) on the same sample or on a second one.

Note: In all cases, refer to current national guidelines concerning HCV diagnosis.

Interpretation of VIDAS Anti-HCV test results should be made taking into consideration the patient history and the results of any other tests or hepatitis C markers.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS Anti-HCV kit.

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LIMITATIONS OF THE METHOD

For the diagnosis of HCV infection, the serological results should be used and interpreted taking into account the patient history, the clinical record, and further tests.

A negative test result does not exclude the possibility of exposure to HCV or infection with HCV. Anti-HCV antibodies may be undetectable in some stages of the infection (acute phase of hepatitis or presence of a serological scar) and in some clinical conditions (immunosuppression) (7,9).

Interference may be encountered with certain sera containing antibodies against reagent components.

This test has not been validated for use with any specimen matrices other than human serum or plasma.

RANGE OF EXPECTED VALUES(1)

Hepatitis C has a worldwide prevalence of 2-3% that varies according to country:

Region of the world	Anti-HCV prevalence (%)
Europe	2.3
Africa	3.2
Americas	1.5
Australia & Oceania	1.2
Asia	2.1
Middle East	4.7
Total	2.4

PERFORMANCE

The following study results demonstrate the conformity of VIDAS Anti-HCV to the Common Technical Specifications of 98/79/CE Directive:

1. Specificity for blood donor population:

5104 blood donor samples (including 2904 fresh samples with a negative status collected ≤ 24 hours previously) obtained from 2 blood transfusion centers, were tested using the VIDAS Anti-HCV assay.

VIDAS Anti-HCV	Status	
	Positive	Negative
Positive	0	20
Negative	0	5084

Diagnostic specificity of the VIDAS Anti-HCV assay on this population: 99.61%

(95% confidence interval: 99.40% - 99.76%)

2. Clinical specificity for hospitalized patients

200 samples with a negative status were tested using the VIDAS Anti-HCV assay.

Diagnostic specificity of the VIDAS Anti-HCV assay on this population: 99.50%.

(95% confidence interval: 97.25% - 99.99%)

3. Diagnostic sensitivity

439 samples with a positive status, including 102 fresh samples (collected ≤ 24 hours previously), were tested using the VIDAS Anti-HCV assay.

Genotypes 1 to 6 were tested:

Genotype	Number tested
1	21
2	21
3	23
4	22
(including non-a sub-types)	
5	6
6	2

Results on tested populations:

Population	Positive VIDAS Anti-HCV /total tested	Diagnostic sensitivity observed (95% confidence interval)
HCV/HIV-negative patient	254/254	100% [98.56% - 100%]
HCV/HIV-positive patient	60/61*	98.36% [91.20% - 99.96%]
Patient with unknown HCV/HIV status	124/124	100% [97.07% - 100%]
Total HCV population	438/439*	99.77% [98.74% - 99.99%]

* The patient who was not detected using VIDAS Anti-HCV either had a low antibody level or was not detected using equivalent methods.

4. Sensitivity for seroconversion panels

Testing of 30 seroconversion panels demonstrated the precocity of detection of the VIDAS Anti-HCV assay. The results are comparable to those obtained using the most sensitive methods.

5. Precision

The repeatability and reproducibility were determined at two sites and calculated according to the recommendations of the CLSI® documents EP5-A2 / EP12-A2.

Four human samples were tested in duplicate using two lots of reagents. Testing was performed twice a day for 10 days on three instruments at one experimental site (N=120). Each reagent lot used a single calibration curve throughout the study. Data from this study are summarized in the following table:

Sample ID / Target value	Sample 1		Sample 2		Sample 3		Sample 4	
	0.26		0.93		1.08		1.19	
	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
Repeatability	0.01	5.6	0.04	4.3	0.05	4.8	0.06	4.7
Reproducibility	0.07	26.9	0.05	5.9	0.07	6.3	0.07	5.8

6. Cross-reactivity

273 samples from patients with a physiological status that can potentially interfere with the detection of hepatitis C antibodies, were tested using VIDAS Anti-HCV. All of the samples were found to be negative with another EIA method (except one CMV IgG+ sample). No disease-related interference was observed for VIDAS Anti-HCV.

	VIDAS Anti-HCV
HSV +	0/10
VZV +	0/10
EBV +	0/10
HIV +	0/10
CMV IgG +	1/11*
LYME Ig+	0/10
HAV IgG +	0/10
HVB (HBcT +)	0/8
HVB (Ag HBs +)	0/10
Syphilis	0/10
Rubella IgG +	0/10
Toxoplasmosis IgG +	0/10
Rheumatoid factor	0/10
Anti-Nuclear Antibody	0/10
Anti-E. coli antibody	0/10
Anti-Pichia Antibody	0/10
Pregnant women**	0/114

* The reference EIA method also showed one false positive sample, but on a different sample.

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REFERENCES BIBLIOGRAPHIQUES

1. LAVANCHY D. Evolving epidemiology of hepatitis C virus - Clin Microbiol Infect 2011; 17: 107-115.
2. ALTER MJ. Epidemiology of hepatitis C virus infection, World J Gastroenterol 2007; 13(17): 2436-2441.
3. GURTSEVITCH VE. Human Oncogenic Viruses: Hepatitis B and Hepatitis C Viruses and Their Role in Hepatocarcinogenesis, Biochemistry (Moscow), 2008, 73(5): 504-513.
4. SIMMONDS P. et al. Consensus Proposals for a Unified System of Nomenclature of Hepatitis C Virus Genotypes, HEPATOLOGY 2005;42: 962-973.
5. VERMEHREN J, SARRAZIN C. New HCV therapies on the horizon, Clin Microbiol Infect 2011; 17: 122-134.
6. WEBSTER D.P, KLENERMAN P., COLLIER J., JEFFERY K. JM. Development of novel treatments for hepatitis C, Lancet Infect Dis 2009; 9: 108-17.
7. CHEVALIEZ S. Virological tools to diagnose and monitor hepatitis C virus infection, Clin Microbiol Infect 2011; 17: 116-121.
8. PENIN F, DUBUISSON, REY F-A, MORADPOUR D, ET PAWLOTSKY JM. Structural Biology of Hepatitis C Virus, HEPATOLOGY 2004; 39: 5-19.
9. THIO CL, NOLT KR, ASTEMBORSKI J, VLAHOV D, NELSON KE et THOMAS DL. Screening for Hepatitis C Virus in Human Immunodeficiency Virus-Infected Individuals, Journal of clinical microbiology, Feb. 2000, 38(2): 575-577.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests

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BIOMERIEUX, the blue logo, VIDAS and SPR are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

CLSI is a trademark belonging to Clinical and Laboratory Standards Institute Inc.

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1. LAVANCHY D. Evolving epidemiology of hepatitis C virus - Clin Microbiol Infect 2011; 17: 107-115.
2. ALTER M.J. Epidemiology of hepatitis C virus infection, World J Gastroenterol 2007; 13(17): 2436-2441.
3. GURTSEVITCH VE. Human Oncogenic Viruses: Hepatitis B and Hepatitis C Viruses and Their Role in Hepatocarcinogenesis, Biochemistry (Moscow), 2008, 73(5): 504-513.
4. SIMMONDS P. et al. Consensus Proposals for a Unified System of Nomenclature of Hepatitis C Virus Genotypes, HEPATOLOGY 2005;42: 962-973.
5. VERMEHREN J, SARRAZIN C. New HCV therapies on the horizon, Clin Microbiol Infect 2011; 17: 122-134.
6. WEBSTER D.P, KLENERMAN P., COLLIER J., JEFFERY K. JM. Development of novel treatments for hepatitis C, Lancet Infect Dis 2009; 9: 108-17.
7. CHEVALIEZ S. Virological tools to diagnose and monitor hepatitis C virus infection, Clin Microbiol Infect 2011; 17: 116-121.
8. PENIN F, DUBUISSON, REY F-A, MORADPOUR D, ET PAWLITSKY JM. Structural Biology of Hepatitis C Virus, HEPATOLOGY 2004; 39: 5-19.
9. THIO CL, NOLT KR, ASTEMBORSKI J, VLAHOV D, NELSON KE et THOMAS DL. Screening for Hepatitis C Virus in Human Immunodeficiency Virus-Infected Individuals, Journal of clinical microbiology, Feb. 2000, 38(2): 575-577.

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VIDAS® TOXO IgM (TXM)

VIDAS TOXO IgM is an automated qualitative test for use on the VIDAS family instruments, for the detection of anti-toxoplasma IgM in serum using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Toxoplasma gondii, an obligate intracellular protozoan parasite, is a significant pathogen in humans. The parasite, whose primary host is the Felidae, is scattered in nature and infects numerous other mammals (1).

Toxoplasmosis is usually benign or asymptomatic in immunocompetent subjects, but can have severe consequences if it occurs in immunodeficient subjects or fetuses. Congenital toxoplasmosis connected with maternal contamination prior to conception is exceptional and most often occurs in immunocompromised women (11). Women who are seronegative are at risk of becoming infected during gestation. The possible ensuing fetal transmission is due to the transplacental migration of toxoplasma during the acute phase of the infection. The frequency and severity of the fetal infection will depend on a number of factors including the virulence of the infecting strain, the size of the original inoculum, the immune response of the patient and when, during gestation, the mother became infected (2-7).

The diagnosis of *T. gondii* infection is most commonly made by biological examination: specific immunoglobulin detection (IgM and IgG) (1, 4, 6).

The diagnosis of an acute acquired infection during pregnancy is established by demonstration of a seroconversion, or by a significant rise in antibody titer in two sequential sera assayed concomitantly (7-10).

VIDAS TOXO IgM is intended for use as an aid in the determination of the immune status of patients.

PRINCIPLE

The assay principle combines an enzyme immunoassay method by immunocapture with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

After a sample dilution step, the IgM are captured by the polyclonal Ab coating the interior of the SPR. The anti-toxoplasma IgM are specifically detected by inactivated toxoplasma antigen (RH Sabin strain), which is itself revealed by an alkaline phosphatase-labeled murine monoclonal anti-toxoplasma antibody (anti-P30).

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antibodies present in the sample. At the end of the assay, an index is automatically calculated by the instrument in relation to the S1 standard stored in memory, and then printed out.

CONTENT OF THE KIT (60 TESTS):

60 TXM strips	STR	Ready-to-use.
60 TXM SPRs 2 x 30	SPR	Ready-to-use. Interior of SPRs coated with anti-human μ chain antibodies (goat).
Positive control TXM 1 x 2 ml (liquid)	C1	Human serum* containing anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
Negative control TXM 1 x 2 ml (liquid)	C2	Human serum* negative for anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide.
Standard TXM 1 x 1 ml (liquid)	S1	Human serum* containing anti-Toxoplasma IgM + protein stabilizer + 1 g/l sodium azide.
Specifications for the factory master data required to calibrate the test:		
- MLE data (Master Lot Entry) provided in the kit, or - MLE bar codes printed on the box label.		
1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

* This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR® is coated during production with anti-human μ chain antibodies (goat). Each SPR is identified by the TXM code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

The strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the TXM strip

1	Sample well.
2	Serum diluent: TRIS buffer (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (300 μ l).
3	Pre-wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 μ l).
4 - 5 - 7 - 8	Wash buffer: TRIS (50 mmol/l) pH 7.4 + protein and chemical stabilizers + 0.9 g/l sodium azide (600 μ l).
6	Conjugate: alkaline phosphatase-labeled immune complex (toxoplasma antigen RH Sabin strain grown in mice (12) - mouse monoclonal anti-P30 antibodies) + 0.9 g/l sodium azide + 0.02% gentamicin (400 μ l).
9	Empty well
10	Cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%, pH 9.2) + 1 g/l sodium azide (300 μ l).

* Signal Word: **DANGER**

Hazard statement

H318 : Causes serious eye damage.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information, refer to the Material Safety Data Sheet.

MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 100 μ l.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory Biosafety Manual - WHO - Geneva - latest Edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents.

It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.

- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS TXM kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR® pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMENS

Specimen type and collection:

Serum

It is recommended that each laboratory checks the compatibility of collection tubes used.

The use of cord blood or neonatal sera has not been validated. For these specimens, we recommend the use of the Toxo-ISAGA technique reference 75 361.

Since the use of hemolyzed, icteric or lipemic samples has not been validated; it is recommended to collect a new sample.

Heat inactivated sera (56°C for 30 minutes) may be tested with VIDAS TOXO IgM.

Specimen stability

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera at -25 ± 6 °C.

Avoid successive freezing and thawing.

INSTRUCTIONS FOR USE

For complete instructions, see the User's Manual.

Reading Master lot data

Before each new lot of reagents is used, enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results.

Note: the master lot data need only be entered once for each lot.

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User's Manual).

Calibration

Calibration, using the standard provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested in **duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
 2. Use one "TXM" strip and one "TXM" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
 3. The test is identified by the "TXM" code on the instrument. The standard must be identified by "S1", and tested in **duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".
 4. If necessary, clarify samples by centrifugation.
 5. Mix the standard, controls and samples using a vortex-type mixer (for serum separated from the pellet).
6. **For this test, the standard, control, and sample test portion is 100 µl.**
7. Insert the "TXM" SPRs and "TXM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
 8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
 9. Restopper the vials and return them to 2-8°C after pipetting.
 10. The assay will be completed within approximately **40 minutes**. After the assay is completed, remove the SPRs and strips from the instrument.
 11. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The instrument calculates a test value for each sample (index), which is the ratio between its RFV and that of the memorized standard.

Interpretation of test results should be made taking into consideration the patient's history, and the results of any other tests performed or other IgM assay methods.

As no international standard is available for the determination of anti-toxoplasma IgM, the VIDAS TOXO IgM reagent is calibrated against collection sera.

Thresholds and interpretation of results

Index	Interpretation
$i < 0.55$	Negative
$0.55 \leq i < 0.65$	Equivocal
$i \geq 0.65$	Positive

Equivocal samples (indices between 0.55 and 0.65) should be retested. If the interpretation remains equivocal, a new sample must be collected.

QUALITY CONTROL

A positive and a negative control are included in each VIDAS TOXO IgM kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

Results cannot be validated if the control values deviate from the expected values.

Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.

In children under 6 months old, production of IgM is sometimes impossible to detect, in which case a significant increase in specific IgG titer must be noted.

In AIDS infected patients, the presence of IgG associated with results from radiological examinations are useful in diagnosis, since the detection of an increase in IgG or the presence of IgM is unfortunately not very reliable. In nearly 100% of AIDS infected patients who develop toxoplasmosis, IgG are detected in the serum, whereas IgM are seldom detected.

Since the persistence of anti-toxoplasma IgM for several months or several years after seroconversion is a regularly described phenomenon, it is not possible to confirm recent infection, based solely on the presence of anti-toxoplasma IgM, without taking into account the serological history of the patient.

PERFORMANCE

Studies performed using VIDAS TOXO IgM gave the following results:

PrecisionWithin-run reproducibility:

3 samples were tested 20 times in 3 different runs on the same VIDAS instrument.

Sera	Mean index	Standard deviation	Intra-assay CV %
S1	0.96	0.04	3.8
C1	1.39	0.03	2.0
C2	0.09	0.01	11.5

Between-run reproducibility

3 samples were tested singly in 7 different runs on the same VIDAS instrument.

Sera	Mean index	Standard deviation	Inter-assay CV %
S4	2.09	0.10	4.8
S5	0.76	0.04	5.3
S6	0.11	0.01	8.2

Relative specificity

The relative specificity of VIDAS TOXO IgM was calculated in comparison with the Toxo-ISAGA IgM technique Ref. 75361 using unselected fresh sera samples collected from an immunized and non-immunized population, sera collected prior to contamination from pregnant women were in contact with toxoplasmosis, cord blood of neonates not contaminated with toxoplasmosis and sera from children suffering from congenital toxoplasmosis but with no anti-toxoplasma IgM.

Out of the 1,483 sera found to be negative with the Toxo-ISAGA IgM kit, 7 were equivocal and 1,465 were negative with the VIDAS TOXO IgM kit, resulting in a relative specificity of **99.25%*** [1,465/1,476] (95% confidence interval: 98.66-99.59%).

Out of the 11 samples found to be positive with the VIDAS TOXO IgM kit and negative with the Toxo-ISAGA IgM kit, 4 were positive with another anti-toxoplasma IgM detection technique and 2 corresponded to follow-up sera collected more than 8 weeks after infection (residual IgM).

*The equivocal results were not used to calculate the performance.

Relative sensitivity

The relative sensitivity of VIDAS TOXO IgM was calculated in comparison with the Toxo-ISAGA IgM technique Ref. 75361 using sera from pregnant women who were in contact with toxoplasmosis during their pregnancy and sera from neonates suffering from congenital toxoplasmosis.

A total of 77 biologically and clinically documented cases of seroconversion and 24 cases of congenital infections were studied.

Acquired toxoplasmosis:

Out of the 150 seroconversion sera found to be positive with Toxo-ISAGA IgM, 144 were positive with VIDAS TOXO IgM. No equivocal samples were found.

The relative sensitivity of sera from pregnant women who were in contact with toxoplasmosis during their pregnancy is **96.00%** (95% confidence interval: 91.43-98.18%).

Out of the 6 discrepant sera:

- Four corresponded to follow-up sera collected more than 8 weeks after the infection (residual IgM). Two of these sera were also found to be negative with two other anti-toxoplasma IgM detection techniques.

- Two corresponded to specimens collected shortly after contamination: absence or low levels of anti-toxoplasma IgG (negative or weakly positive Lyse Test) and low levels of IgM (weakly positive Toxo-ISAGA IgM result). In both these cases, seroconversion was detected with the specimens collected later.

Congenital toxoplasmosis:

Out of the 30 congenital toxoplasmosis sera found to be positive with the Toxo-ISAGA IgM technique, 27 were positive with VIDAS TOXO IgM. No equivocal samples were found.

The relative sensitivity in comparison with the Toxo-ISAGA IgM technique using sera from neonates suffering from congenital toxoplasmosis is **90.00%** (95% confidence interval: 73.47-97.89%). The positive predictive value for this population is 100% (95% confidence interval: 87.23-100%).

Out of the 3 discrepant sera, one was collected at birth from cord blood. The other two were collected at one and three months respectively. One serum showed low levels of toxoplasma IgM (weakly positive Toxo-ISAGA IgM result).

In view of the results obtained for cases of congenital toxoplasmosis, it is recommended that all sera from neonates with congenital toxoplasmosis found to be negative with VIDAS TOXO IgM, be confirmed using the Toxo-ISAGA IgM technique.

CROSS REACTIVITY AND RELEVANT INTERFERENTS

During the evaluations, 24 samples containing rheumatoid factor negative for anti-toxoplasma IgM and 10 samples containing anti-EBV IgM negative for anti-toxoplasma IgM were tested. All of these samples were found to be negative with VIDAS TOXO IgM.

PREVALENCE

T. gondii is a strict pathogen whose prevalence differs greatly from one country to another, or even from one region to another.

Contamination by *T. gondii* can vary depending on cultural customs and eating habits, resulting in a prevalence ranging from less than 10% in certain regions of Northern Europe to more than 90% in Africa.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

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INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

REVISION HISTORYChange type categories :

N/A

Correction

Technical change

Administrative

Not applicable (First publication)

Correction of documentation anomalies

Addition, revision and/or removal of information related to the product

Implementation of non-technical changes noticeable to the user

Note:*Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release date	Part Number	Change Type	Change Summary
2015/01	05933R	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE

BIOMERIEUX, the blue logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

Any other name or trademark is the property of its respective owner.

 **bioMérieux SA**
376 Chemin de l'Orme
69280 Marcy-l'Etoile - France

673 620 399 RCS LYON
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
www.biomerieux.com

VIDAS® Anti-HBc Total II (HBCT)

VIDAS Anti-HBc Total II (HBCT) is an automated qualitative test for use on the VIDAS family instruments, for the detection of total antibodies against the hepatitis B core antigen (anti-HBc) in human serum or plasma, using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

The hepatitis B virus is responsible for acute and chronic hepatitis infections. Acute hepatitis can be asymptomatic or present symptoms of varying severity which may progress to fulminant hepatitis in 0.1 to 0.5% of cases. Chronicity occurs in 5 to 10% of cases in adults, but up to 90% of cases in infants following perinatal transmission. Currently, approximately 300 million people worldwide are chronic carriers of the virus (1). Chronic hepatitis may be asymptomatic and lead to liver lesions of varying severity, possibly evolving to cirrhosis, with an evolution in 5% of cases to hepatocellular carcinoma (2). The hepatitis B virus can be transmitted by parenteral or perinatal pathways or through sexual contact. Persons most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, close friends and family of an infected subject, and newborns of an infected mother (2).

CLINICAL SIGNIFICANCE

Total antibodies against the hepatitis B core antigen (anti-HBc IgM and IgG) can be detected in the serum of patients with acute or chronic hepatitis B or in recovered patients. Anti-HBc total antibodies therefore act as an epidemiological indicator of current or former HBV infections.

In cases of acute hepatitis, anti-HBc antibodies (IgM and IgG) are generally detected 2 to 4 weeks after HBs and HBe antigens appear (3). Whereas anti-HBc IgM are transitory and progressively decrease whether infection develops towards recovery or chronicity, high anti-HBc IgG titers can be detected during the infection and after recovery (2). Anti-HBc IgG can persist for several years, and even throughout the patient's lifetime.

In cases of chronic hepatitis, only titration of anti-HBc IgM antibodies can indicate an active phase of the disease (4). Anti-HBc antibodies do not provide protection. Only the presence of anti-HBs antibodies can confirm immunity.

PRINCIPLE

The VIDAS Anti-HBc Total II assay is an enzyme-linked fluorescent immunoassay (ELFA) based on an inhibition principle.

The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed Reagent Strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

After dilution, the sample is incubated in the SPR. The anti-HBc antibodies (IgM and IgG) present in the sample bind with the recombinant HBc antigen coating the interior of the SPR. Unbound sample components are eliminated during the washing steps.

The solid phase is then incubated with the conjugate: alkaline phosphatase-labeled monoclonal anti-HBc antibody. This conjugate binds to the solid phase HBc antigen sites that have not bound to serum antibodies. Unbound conjugate is removed by washing.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm.

The intensity of the fluorescence is inversely proportional to the quantity of anti-HBc antibodies present in the sample. Results are analyzed automatically by the instrument and are expressed as an index calculated using a standard.

CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS:

60 HBCT Reagent Strips	STR	Ready-to-use.
60 HBCT SPRs 2 x 30	SPR	Ready-to-use. Interior of the SPR coated with recombinant HBc antigen.
HBCT Positive control 1 x 1.5 ml (liquid)	C1	Ready-to-use. Human* serum containing anti-HBc antibodies + 0.2 g/l gentamicin sulfate + 1 g/l sodium azide. Index strictly inferior to 1.
Negative control 1 x 1.9 ml (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives. MLE data indicate the index: confidence interval ("Control C2 (-) Test Value Range).
HBCT Standard 2 x 2 ml (lyophilized)	S1	Human* serum containing anti-HBc antibodies + 0.2 g/l gentamicin sulfate. The standard must be reconstituted with 2 ml of distilled sterile water (measured exactly). Allow to dissolve for at least 20 min and then mix using a vortex. After reconstitution, the standard can be kept for up to 1 month at 2-8°C. Aliquoted and frozen at -25 ± 6°C, it can be stored for up to 6 months. Avoid repeated freeze/thaw cycles.

Specifications for the factory master data required to calibrate the test:

- MLE data (Master Lot Entry) provided in the kit,

or

- MLE bar code printed on the box label.

1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib.

* This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1, HIV 2 and HCV. However, since no known test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with recombinant HBc antigen. Each SPR is identified by the "HBCT" code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

The strip

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the HBCT strip

Wells	Reagents
1	Sample well
2	Sample diluent: TRIS buffered saline (50 mmol/l) pH = 7.4 with protein and chemical stabilizers + 0.9 g/l sodium azide (300 µl).
3	Wash buffer: TRIS buffered saline (50 mmol/l) pH = 7.4 with protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
4 - 6 - 7 - 8 - 9	Wash buffer: TRIS buffered saline (50 mmol/l) pH = 7.4 with protein and chemical stabilizers + 1 g/l sodium azide (600 µl).
5	Conjugate: alkaline phosphatase-labeled anti-HBc monoclonal antibody (mouse) + 0.9 g/l sodium azide (400 µl).
10	Optical Cuvette with substrate: 4-Methyl-umbelliferyl phosphate + diethanolamine DEA* (0.62 mol/l or 6.6%) pH 9.2 + 1 g/l sodium azide (300 µl).

* Signal Word: **DANGER**

**Hazard statement**

H318 : Causes serious eye damage.

Precautionary statement

P280 :Wear protective gloves/protective clothing/eye protection/face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information, refer to the Material Safety Data Sheet.

MATERIALS AND DISPOSABLES REQUIRED BUT NOT PROVIDED

- Pipette with disposable tip to dispense 2 ml and 150 µl.
- Powderless, disposable latex gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument

WARNINGS AND PRECAUTIONS

- **For in vitro diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**

- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use powderless gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.

- The optical cuvette with substrate (well 10) contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the VIDAS User's Manual).

STORAGE CONDITIONS

- Store the VIDAS Anti-HBC Total II kit at 2-8°C.
- **Do not freeze SPRs and strips.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

SPECIMENS

Specimen type and collection

Use non-contaminated sera (plain tube, separator gel) or plasma (validated anticoagulants: EDTA, lithium heparin, citrate).

None of the following factors have been found to significantly influence this assay:

- hemolysis (after spiking samples with hemoglobin: 0 to 340 µmol/l (monomer)),
- lipemia (after spiking samples with lipids: 0 to 5 mg/ml equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 450 µmol/l).

Do not deplete samples.

Specimen stability

Samples containing impurities must be clarified by centrifugation prior to testing. Samples can be stored at 2-8°C in stoppered tubes for a maximum of 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C (**freeze once only**). A study performed on frozen samples over a period of 2 months, showed that the quality of results is not affected.

INSTRUCTIONS FOR USE

For complete instructions, see the User's Manual.

Reading Master lot data

Before each new lot of reagents is used, enter the specifications (or factory master data) into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results.

Note: the master lot data need only be entered once for each lot.

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User's Manual).

Calibration

Calibration, using the standard provided in the kit, must be performed upon receipt of a new lot of reagents after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested **in triplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HBCT" strip and one "HBCT" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBCT" code on the instrument. The standard must be identified by "S1", and tested **in triplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.
4. Mix the standard, control(s) and sample(s) using a Vortex-type mixer (for serum or plasma separated from the pellet).
5. **For this test, the standard, control, and sample test portion is 150 µl.**
6. Insert the "HBCT" SPRs and "HBCT" strips into the instrument. Check to make sure the color labels with the three letter assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the required temperature after pipetting.
9. The assay will be completed within approximately 90 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument and test values are expressed as an index calculated using a standard.

The patient RFV is interpreted by the VIDAS system as follows:

$$i = \text{test value} = \text{patient RFV} / \text{S1 standard RFV}$$

A report is printed which records the test value and the interpretation. The following table shows the thresholds and the interpreted results:

Index	Interpretation
$i < 1$	Presence of anti-HBc antibodies
$1 \leq i < 1.4$	Equivocal result
$i \geq 1.4$	Absence of anti-HBc antibodies

All equivocal results must be confirmed using a second sample.

Interpretation of test results should be made taking into consideration the patient history, and the results of any other tests performed.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HBCT kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2.

Results cannot be validated if the control values deviate from the expected values.

Note

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

2. Sensitivity:

a) Results obtained using the SNTS 1996 panel in 2 blood transfusion centers:

Characteristics	Panel sections	VIDAS results	Samples
Positive anti-HBc antibodies: concordant screening for the 8 EIA kits studied	7	7 pos.	14, 15, 16, 17, 18, 19 and 20
Anti-HBc antibodies: positive anti-HBc plasma dilutions, anti-HBs and/or anti-HBe positive	5	2 pos. (1 and 2), 2 equiv. (3 and 4) 1 neg. (5)	1, 2 to 5 (sample 1 - increasing dilutions)
	4	2 pos. (21 and 22), 1 equiv. (23) 1 neg. (24)	21 to 24 (same sample - increasing dilutions)
	4	2 pos. (25 and 26), 1 equiv. (27) 1 neg. (28)	25 to 28 (same sample - increasing dilutions)
	4	3 pos. (29, 30, 31), 1 neg. (32)	29 to 32 (same sample - increasing dilutions)
Anti-HBc antibodies: discrepant screening for the 8 EIA kits studied	4	1 pos. (8) 1 pos./equiv. (9) 2 neg. (10,12)	8, 9, 10, 12
Absence of anti-HBc antibodies	1	1 neg.	6

LIMITATIONS OF THE METHOD

Interference may be encountered with certain sera containing antibodies against reagent components. For this reason, assay results should be interpreted taking into consideration the patient history, and the results of any other tests performed.

PERFORMANCE

Studies performed using VIDAS Anti-HBc Total II gave the following results:

1. Specificity for blood donor population

The results obtained from 4 blood transfusion centers, where 5064 blood donor samples were tested, are as follows:

VIDAS Anti-HBc Total II	Final interpretation: other EIA techniques / other markers	
	Positive	Negative
Positive	47	3
Equivocal	0	8
Negative	0	5006

The eight samples giving an equivocal result were not included in the performance calculation.

Relative specificity with VIDAS Anti-HBc Total II for this population: 99.94%

(95% confidence interval 95%: 99.82%- 99.98%.)

b) Clinical study: 3 studies were performed:

Study 1: 228 documented positive anti-HBc samples were tested in comparison with another reference technique. Discrepant results between the 2 techniques were confirmed using a 3rd reference technique.

Clinical category	Positive VIDAS Anti-HBc Total II
Acute hepatitis B virus N = 8	8 (100%)
Chronic HBV carriers with replication marker (positive HBe Ag and/or positive viral DNA) N = 49	49 (100%)
Treated chronic HBV carriers (anti-viral treatment or anti-HBs immunoglobulins) N = 47	47 (100%)
Chronic HBV carriers without standard replication marker (negative HBe Ag, viral DNA by negative hybridization) N = 42	42 (100%)
Subjects who have acquired natural immunity against HBV N = 56	54* (100%)
"Isolated anti-HBc" (negative Ag HBs and anti-HBs) N = 26	25 (96%)

* Two uninterpretable samples after confirmation (equivocal with VIDAS Anti-HBc Total II, negative and positive with other tests) were excluded from the performance calculation.

Relative sensitivity with VIDAS Anti-HBc Total II: 99.56% (95% confidence interval: 97.46% - 99.9%).

Study 2: study performed using 130 documented positive anti-HBc samples.

Sample category	VIDAS Anti-HBc Total II		
	Positive	Equivocal	Negative
Positive	129	1*	0

* Sample excluded from sensitivity calculation.

Clinical sensitivity = **100%**

(95% confidence interval: 96.99% - 100%)

Study 3: study performed using 87 documented positive anti-HBc samples.

Sample category	VIDAS Anti-HBc Total II		
	Positive	Equivocal	Negative
Positive	85*	0	0

* Two uninterpretable samples after confirmation were excluded from the performance calculation.

Clinical sensitivity = **100%**

(95% confidence interval 95%: 95.75% - 100%)

c) Diagnostic sensitivity

33 fresh samples with a positive status (collection < 24 hours) were tested, 31 were found to be positive, 1 was found to be equivocal and 1 was found to be negative with VIDAS anti-HBc Total II. 31 fresh negative samples (collection < 24 hours) were tested and found to be negative with VIDAS anti-HBc Total II.

3. Precision

Precision was evaluated using a positive sample, an equivocal sample and a negative sample. Each sample was tested in duplicate in two different runs per day for twenty days on three sites.

Within-run reproducibility (intra-assay precision) and between-run reproducibility (total precision) were calculated according to the recommendations of the NCCLS Document EP5-T2, Volume 12-4.

The combined results for the three sites are given below:

a) Within-run reproducibility:

Sample	N	Mean index	CV (%)
Negative	240	2.569	4.36
Equivocal	240	1.180	5.65
Positive	240	0.082	7.53

b) Between-run reproducibility:

Regarded as total precision taking into account all sources of variation.

Sample	N	Mean index	CV (%)
Negative	240	2.569	7.42
Equivocal	240	1.180	8.10
Positive	240	0.082	9.25

4. Cross-reactivity

	VIDAS Anti-HBc Total II (positive)
VHC + (Anti-HBc -)	0 / 10
EBV + (Anti-HBc -)	0 / 10
HIV + (Anti-HBc -)	0 / 10
FR + (Anti-HBc -)	0 / 21
ANA + (Anti-HBc -)	0 / 20
Rubella IgG + (Anti-HBc -)	0 / 11
CMV IgG + (Anti-HBc -)	0 / 11
CMV IgM + (Anti-HBc -)	0 / 10
Pregnant women (Anti-HBc -)	0 / 24
Vaccinated: Anti-HBs + (Anti-HBc -)	1 / 20
HSV + (Anti-HBc -)	0 / 6
Patients who have received VIII (Anti-HBc -)	0 / 15

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. MILICH D.R., Immune response to the hepatitis B virus: infection, animal models, vaccination, VIRAL HEPATITIS, 1997, 3, 63 -103.
2. HOLLINGER F.B., Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. KURSTAK E., Viral Hepatitis, Current Status and Issues, Springer-Verlag Wien New York, 1993, 11, 93-104.
4. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIAONTI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, Journal of Hepatology, 1993, 19, 431-436.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

REVISION HISTORYChange type categories :

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

Note: *Minor typographical, grammar, and formatting changes are not included in the revision history.*

Release date	Part Number	Change Type	Change Summary
2015/01	09577G	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS
2015/06	09577H	Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE
2016/05	09577I	Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.
Any other name or trademark is the property of its respective owner.

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VIDAS® HBs Ag Ultra (HBS)

VIDAS HBs Ag Ultra (HBS) is an automated qualitative test for use on the VIDAS family instruments for the detection of hepatitis B surface antigen (HBs Ag) in human serum or plasma, using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

The hepatitis B virus is responsible for acute and chronic hepatitis infections. Acute hepatitis can be asymptomatic or present symptoms of varying severity which may progress to fulminant hepatitis in 0.1 to 0.5% of cases. Chronicity occurs in 5 to 10% of cases in adults, but up to 90% of cases in infants following perinatal transmission. Currently, approximately 350 million people worldwide are chronic carriers of the virus (1). Chronic hepatitis may be asymptomatic and lead to liver lesions of varying severity, possibly evolving to cirrhosis, with an evolution in 5% of cases to hepatocellular carcinoma (2). The hepatitis B virus can be transmitted by parenteral or perinatal pathways or through sexual contact. Persons most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, as well as close friends and family of an infected subject, and newborns of an infected mother (2).

The discovery of the Australia antigen in 1965 – later known as HBs antigen - combined with viral hepatitis (3, 4) was a major breakthrough in the diagnosis of hepatitis B.

HBs antigen appears several days to several weeks after contact with the virus and can persist for several months. Persistence of HBs antigen for more than 6 months serologically defines chronic HBV infection. Disappearance of the HBs antigen is normally followed by the appearance of anti-HBs antibodies, which is a sign of recovery. The anti-HBs antibody assay is performed to confirm efficacy of the HBV vaccine. Anti-HBc antibodies are normally detected at the onset of the disease (2). Anti-HBc IgM titers are elevated (> 100 UPEI/ml) during acute hepatitis, then their titer decreases or disappears. However, during chronic hepatitis, the appearance of low IgM titers, reflecting hepatic cytolysis, confirms the active phase of the disease (5). Total anti-HBc antibodies, mainly IgG, are detected during acute and chronic hepatitis and persist after recovery.

HBe antigen is a circulating protein with a sequence very similar to that of HBc antigen, but with distinct antigenicity. During acute or chronic hepatitis, the presence of HBe antigen is generally associated with intense viral replication. Its disappearance coincides with the appearance of anti-HBe antibodies. Anti-HBe seroconversion is generally an indicator of a favorable prognosis of recovery (6). Nevertheless, the possible selection of a mutant strain incapable of synthesizing Hbe antigen may also lead to anti-HBe seroconversion, thereby limiting the prognostic value of HBe and anti-HBe markers (7). In this case, only the viral DNA assay is capable of directly detecting viral replication.

PRINCIPLE

The VIDAS HBs Ag test is an enzyme-linked fluorescent immunoassay (ELFA) that is performed in the automated VIDAS system (see User's Manual). This assay can be performed according to 2 protocols: HBL long protocol (90 minutes), HBS short protocol (60 minutes).

The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready to use and pre-dispensed in the sealed reagent strips.

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After a preliminary washing step, the antigen present in the sample will bind simultaneously to the monoclonal antibody (8) coating the interior of the SPR and to the antibody conjugated with biotin. Unbound sample components are washed away. The antigen bound to the solid phase and to the biotinylated antibody is in contact with streptavidine conjugated with alkaline phosphatase, which will bind with biotin. Another wash step follows and removes unbound components.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of antigen present in the sample.

At the end of the assay, results are analyzed automatically by the instrument and are expressed as an index calculated using a standard.

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60 HBS strips	STR	Ready-to-use.
60 HBS SPRs 2 x 30	SPR [®]	Ready-to-use. Interior of SPRs coated with two monoclonal anti-HBs Ag antibodies (mouse) (5).
HBS standard 3 x 1 ml (lyophilized)	S1	Serum base* supplemented with inactivated human plasma HBs antigen + preservatives. The standard must be reconstituted with 1 ml of sterile distilled water (measured exactly). Allow to dissolve for at least 20 mins, then mix using a vortex. After reconstitution, the standard must be stored in aliquots at $-25 \pm 6^{\circ}\text{C}$ for up to 6 months, by freezing within an hour of reconstitution. Avoid successive freezing and thawing.
HBS positive control 1 x 1.5 ml (liquid)	C1	Ready-to-use. Serum base* supplemented with inactivated human plasma HBs antigen + 1 g/l sodium azide. MLE data indicate the index: confidence interval ("Control C1(+) Test Value Range)
Negative control 1 x 1.9 ml (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

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The SPR

The interior of the SPR is coated during production with monoclonal anti-HBs Ag antibody (mouse). Each SPR is identified by the HBS code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Hazard statement

H318 : Causes serious eye damage.

H373 : May cause damage to organs through prolonged or repeated exposure.

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- Pipette with disposable tip to dispense 1 ml and 150 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- **Sample to sample contamination through contact with gloves: since high concentrations of HBs Ag may be encountered, it is strongly recommended to keep tips on a stand prior to use. It is also recommended to keep one tube of sample aside for testing HBs antigen.**
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.

- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The wash buffer (wells 3, 4, 7, 8, 9) contains a harmful agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- The optical cuvette with substrate (well 10) contains an irritant agent (diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS HBs Ag Ultra kit at 2-8°C.
- **Do not freeze reagents, with the exception of S1 standard after reconstitution.**
- Store all unused reagents at 2-8°C, with the exception of S1 standard which must be kept frozen after reconstitution.
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

SPECIMENS

Specimen type and collection

Use sera (plain tube, tube with separator gel, plain tube with beads) or plasmas collected in lithium heparin.

Serum and plasma should be stored separated from the pellet. Samples containing impurities must be clarified by centrifugation.

The results obtained were not found to be influenced by: icteric samples (bilirubin concentrations up to 500 µmol/l), hemolyzed samples (hemoglobin concentrations up to 270 µmol/l of monomer) and lipemic samples (up to 30 mg/ml).

Do not inactivate specimens.

Specimen stability

Samples can be stored for 5 days in stoppered tubes at 2-8°C. If longer storage is required, freeze the sera or plasma at -25 ± 6°C. A study performed on samples frozen for 2 months, showed that the quality of results is not affected.

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The standard, identified by S1, must be tested **in duplicate** (see User's Manual). The standard value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HBS" strip and one "HBS" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBS" code on the instrument, or by the "HBL" code for the long protocol. The standard must be identified by "S1", and tested **in duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.

4. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

5. For this test, the calibrator, control, and sample test portion is 150 µl.

6. Insert the "HBS" SPRs and "HBS" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the required temperature after pipetting.
9. The assay will be completed within approximately 60 to 90 minutes, depending on the protocol selected. After the assay is completed, remove the SPRs and strips from the instrument.
10. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet. The patient RFV is interpreted by the VIDAS system as follows:

$$i = \text{test value} = \text{patient RFV} / \text{standard RFV}$$

The test value and interpretation are also indicated on the result sheet. Interpretation of results according to test values is as follows:

Test value		
Short protocol	Long protocol	Interpretation
$i < 0.13$	$i < 0.10$	Negative
$i \geq 0.13$	$i \geq 0.10$	Positive

If a positive result is obtained for a patient with no previous history, the assay must be repeated and confirmed using a neutralization test (VIDAS HBs Ag Ultra Confirmation ref. 30 317) or other tests.

A positive sample should be confirmed using the same assay protocol (short or long) as was initially used to test the sample.

Before retesting, samples should be centrifuged again so as to eliminate any interference caused by fibrin fragments or cell elements.

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VIDAS HBs Ag Ultra is calibrated against the panel of the Société Française de Transfusion Sanguine (French Blood Transfusion Society) (adw2/ayw3 mixture expressed in ng/ml).

The analytical sensitivity of VIDAS HBs Ag Ultra, which was determined using this panel, is less than 0.20 ng/ml with the short protocol (HBS), and less than 0.15 ng/ml with the long protocol (HBL).

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HBs Ag Ultra kit.

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LIMITATIONS OF THE METHOD

- As interference (i.e. anti-idiotypic antibody) may be encountered with certain sera, the test should only be declared positive after taking into account the patient's history and the results of the other hepatitis B markers.
- A negative Ag HBs result does not allow infection by the hepatitis B virus to be excluded. The Ag HBs serum concentration may be below the analytical sensitivity of the reagent. The presence of a modified HBs antigen (variant) cannot be excluded; the antigen may, in this case, have been incorrectly recognized or not recognized by the antibodies in the reagent. The results of this test must be interpreted taking into consideration the patient's history and the results of any other tests performed (confirmation by neutralization, HBV DNA etc.).
- In rare cases it is possible to detect both HBs antigen and anti-HBs antibodies.
- This assay has been validated for serum and plasma and should not be used for other biological fluids such as saliva, CSF or urine.
- This assay should not be used with specimens collected post-mortem.
- Do not use serum pools.

RANGE OF EXPECTED VALUES (9)

The world wide prevalence of hepatitis B varies depending on the countries:

Region of the world	Prevalence of HBs Ag (%)
Europe, North America, Australia	0.2 - 0.5
Eastern Europe, Mediterranean basin, Russia, southwest Asia, South America	2 - 7
Southeast Asia, tropical Africa	8 - 20

PERFORMANCE

The results of the studies, which demonstrated the conformity of VIDAS HBs Ag Ultra to the Common Technical Specifications of Directive 98/79/EC, are as follows:

1. Specificity for blood donor population

5059 blood donor samples from 2 blood transfusion centers were tested using the short protocol and the long protocol. No discrepancies were observed between the two protocols.

VIDAS HBs Ag Ultra (HBS/HBL)	Final interpretation of other EIA techniques	
	Positive	Negative
Positive	0	0
Negative	0	5059

Relative specificity with VIDAS HBs Ag Ultra for this population: 100.00%
(95% confidence interval: 99.87% - 100.00%).

2. Clinical specificity

a) for hospitalized patients:

200 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.
(95% confidence interval: 98.04% - 100.00%).

b) for outpatients in a screening center:

100 samples were tested using the short protocol, the long protocol and another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

Relative specificity of VIDAS HBs Ag Ultra, with the two protocols (short and long), for this population: 100.00%.
(95% confidence interval: 96.38% - 100.00%).

3. Analytical sensitivity

An external study performed using the panel of the French Blood Transfusion Society, showed a sensitivity of 0.12 ng/ml with the short protocol (HBS) and 0.08 ng/ml with the long protocol (HBL).

The sensitivity of VIDAS HBs Ag Ultra, determined using the international standard NIBSC 00/588, was estimated at 0.05 IU/ml with the long protocol (HBL) and 0.075 IU/ml with the short protocol (HBS).

4. Diagnostic sensitivity

30 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the short protocol (HBS). 30 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the short protocol (HBS).

31 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive using VIDAS HBs Ag Ultra with the long protocol (HBL). 31 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative using VIDAS HBs Ag Ultra with the long protocol (HBL).

5. Detectability

A study was performed using 506 samples characterized as positive, including 28 subtyped or genotyped samples and 46 samples from patients with acute hepatitis. Detectability with the long protocol (threshold of 0.08 ng/ml) and with the short protocol (threshold of 0.12 ng/ml), is 100.00 %.

(95% confidence interval: 99.22% - 100.00%).

6. Sensitivity of seroconversion panels

During a study, 32 seroconversion panels were tested using the two VIDAS HBs Ag Ultra protocols.

With the short protocol, HBs Ag was detected earlier by the VIDAS HBs Ag Ultra than with the comparison method in 13 out of 32 seroconversion panels (12 panels one specimen collection earlier, 1 panel two specimen collections earlier). HBs Ag was detected in one panel, one specimen collection later.

With the long protocol, 14 of the 32 panels tested were detected earlier than with the comparison method (9 panels one specimen collection earlier, 5 panels two specimen collections earlier).

7. Sensitivity of HBs Ag mutants

A panel of 27 recombinant proteins imitating the most significant mutations in HBs Ag amino acid sequences and 4 native samples harboring HBs Ag mutants were tested successfully (10).

8. Precision

Intra-assay and inter-assay reproducibility were determined at two sites and calculated according to the recommendations of the NCCLS EP5-T2 document, volume 12-4.

Intra-assay reproducibility

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	3.41	3.07	2.68
Weak positive	72	0.21	4.98	0.15	3.94
Negative	80	0.01	0.00 *	0.02	0.00 *

* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

Inter-assay reproducibility

The total precision takes into account all the sources of variability.

The results are expressed as an index:

Sample	N	Short protocol		Long protocol	
		Mean	CV (%)	Mean	CV (%)
Positive	72	3.42	11.58	3.07	6.17
Weak positive	72	0.21	12.41	0.15	6.75
Negative	80	0.01	0.01 *	0.02	0.01 *

* since the index value is too low to be used for the calculation of the CV, only the standard deviation is indicated.

9. Cross-reactivity

367 samples from patients whose physiological status is likely to affect the detection of HBs antigen, were tested using both the short and the long protocols. All the samples were found to be negative with another EIA technique. No discrepancies were observed between the two protocols and the technique used as reference.

	VIDAS HBs Ag Ultra positive samples
HCV +	0/20
EBV +	0/10
HIV +	0/10
CMV IgG +	0/10
HAV IgG +	0/10
HSV +	0/10
Syphilis	0/10
Rubella IgG +	0/10
Toxoplasmosis IgG +	0/10
Rheumatoid factor	0/9
Anti-nuclear antibodies	0/9
Non-viral hepatitis	0/5
Dialysis patients	0/10
Children less than 15 years old	0/10
Vaccinated: Anti-HBs +	0/10
Pregnant women*	0/214**

* including 22 multipara

** including 2 false positive results which were not repeated with the long protocol.

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LITERATURE REFERENCES

1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. Journal of Clinical Virology 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. JAMA, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, Proc Natl Acad Sci USA, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIAONTI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, Journal of Hepatology, 1993, 19, 431-436.
6. ALDERSHVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, Journal of Infectious Diseases, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in The Molecular Medicine of Viral Hepatitis, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - Journal of Medical Virology - 2001, vol. 65, p.241-249.
9. ZUCKERMAN A.J., Hepatitis Viruses. In: Baron S, eds. Medical Microbiology, 4th ed. The University of Texas M Branch at Galveston, 1996: 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, Journal of Virological Methods 135, 2006, 109-117.

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Release date	Part Number	Change Type	Change Summary
2015/01	11728I	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) - RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE
2016/05	11728J	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

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1. LAVANCHY D., Worldwide epidemiology of HBV infection, disease burden, and vaccine prevention. *Journal of Clinical Virology* 34 suppl. 1 (2005) S1-S3.
2. HOLLINGER F.B., Hepatitis B virus, in *Fields Virology*, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. BLUMBERG B.S., ALTER H.J., VISNICH S. JAMA, A "New" Antigen in Leukemia Sera, 1965, 191, 541-546.
4. PRINCE A.M., An antigen detected in blood during the incubation period of serum hepatitis, *Proc Natl Acad Sci USA*, 1968, 60, 814-821.
5. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIANTONI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, *Journal of Hepatology*, 1993, 19, 431-436.
6. ALDERSHIVILE J., FRÖSNER G.G., NEILSEN J.O. et al., Hepatitis B e antigen and antibody measured by radioimmunoassay in acute hepatitis B surface antigen-positive hepatitis, *Journal of Infectious Diseases*, 1980, 141, 293-298.
7. TONG S.T. and TREPO C., The HBe-minus mutants of hepatitis B virus in *The Molecular Medicine of Viral Hepatitis*, Harrison T.J. and Zuckerman A.J. ed., 1997, 89-104.
8. JOLIVET-REYNAUD. C., LESENECHAL. M., O' DONNELL. B. et al.- Localization of hepatitis B surface antigen epitopes present on variants and specifically recognised by anti-hepatitis B surface antigen monoclonal antibodies - *Journal of Medical Virology* - 2001, vol. 65, p.241-249.
9. ZUCKERMAN A.J., Hepatitis Viruses. In: Baron S, eds. *Medical Microbiology*, 4th ed. The University of Texas M Branch at Galveston, 1996 : 849-863.
10. WEBER B., VAN DER TAELEM-BRULE N., BERGER N. et al. Evaluation of a new automated assay for hepatitis B surface antigen (HBs Ag) detection VIDAS HBs Ag Ultra, *Journal of Virological Methods* 135, 2006, 109-117.

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VIDAS® Anti-HBs Total II (AHBS)

VIDAS Anti-HBs Total II is an automated quantitative test for use on the VIDAS family of instruments, for the immunoenzymatic detection of antibodies to hepatitis B surface antigen (Anti-HBs) in human serum or plasma using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

The hepatitis B virus (HBV) is responsible for acute and chronic hepatitis infections, possibly evolving to cirrhosis or primary liver cancer. Chronicity occurs in 5 to 10% of cases in adults, but up to 90% of cases in infants following perinatal transmission. Currently, more than 360 million people worldwide are chronic carriers of the virus (1). The hepatitis B virus can be transmitted by parenteral and perinatal pathways or through sexual contact. People most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, as well as close friends and family of an infected subject, and newborns of an infected mother (2).

HBs antigen (HBsAg) appears several days to several weeks after contact with the virus and can persist for several months: in this case, the infection is considered to be "chronic". Disappearance of the HBs antigen is normally followed by the appearance of Anti-HBs, which is a sign of recovery. In this case, the presence of Anti-HBs is associated with that of Anti-HBc.

The detection of Anti-HBs is thus performed to monitor infected patients, but also to check the efficacy of immunization against the hepatitis B virus (vaccination by HBsAg).

Strategies of vaccination against HBV depend on the epidemiological situation. It is strongly recommended for subjects at risk. In vaccinated individuals, only the Anti-HBs are positive. The persistence of Anti-HBs is related to the initial titer obtained after the complete vaccination series has been administered.

The "European Consensus Group" recommends performing Anti-HBs determination 1 to 3 months following vaccination (3). It is only possible to affirm that the immunological memory will confer long-term protection to immunocompetent adults if titers are greater than 100 mIU/mL (3, 4).

Nevertheless, numerous countries have adopted 10 mIU/mL as the lowest titer indicating protective immunity against HBV infection (4).

For the detection of Anti-HBs, there is no reference measurement method that can guarantee immunity to HBV, particularly for low titers. The nature of the Ag used in the vaccine and in Anti-HBs detection tests, may cause the results obtained to vary (5).

USE

- monitoring of patients infected with the hepatitis B virus, treated patients, and patients on the road to recovery,
- monitoring of pre- and/or post-vaccine immunity status of patients.

PRINCIPLE

The assay principle combines an enzyme immunoassay sandwich method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and are pre-dispensed in the sealed reagent strips.

The five reaction steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times. Each step is followed by a wash cycle which eliminates unbound components:

- Specific binding of antibody in the sample to the ad and ay HBs surface antigen coated on the interior of the SPR.
- Formation of an antibody/antigen complex bound to the SPR with the biotinylated antigens in the diluent (HBs sub-types ad and ay).
- Wash cycle.
- Binding of the biotin to the anti-biotin-alkaline phosphatase conjugate.
- Wash cycle.
- Detection: alkaline phosphatase catalyzes the hydrolysis of the substrate (4-Methyl-umbelliferyl phosphate) into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm.

The intensity of the fluorescence is proportional to the quantity of Anti-HBs in the sample. At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (60 TESTS):

60 AHBS Strips	STR	Ready-to-use.
60 AHBS SPRs 2 x 30	SPR	Ready-to-use. Interior of SPRs coated with recombinant hepatitis B surface antigen, subtypes ad and ay.
AHBS Positive control 1 x 2.4 mL (liquid)	C1	Ready-to-use. Human serum* containing Anti-HBs +preservatives. MLE data indicate the confidence interval in mIU/mL ("Control C1 (+) Dose Value Range").
AHBS Negative control 1 x 5.8 mL (liquid)	C2	Ready-to-use. Human serum* without Anti-HBs + preservatives. MLE data indicate the confidence interval in mIU/mL ("Control C2 (-) Dose Value Range").
Calibrator 2 x 2.4 mL (liquid)	S1	Human serum* containing Anti-HBs + preservatives.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

* This product has been tested and shown to be negative for HBs surface antigen and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR®

The interior of the SPR is coated during production with recombinant hepatitis B surface antigen subtypes ad and ay (6). Each SPR is identified by the code "AHBS". Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the AHBS strip:

Wells	Reagents.
1	Sample well.
2	Sample diluent: calf serum + Tetraborate (pH 7.8) + biotin-labeled ad and ay surface antigen + fetal calf serum + preservatives (300 µL).
3 - 4 - 6 - 7 - 8 - 9	Wash solution: buffer pH 7.8 + Tween 20 + preservatives (600 µL).
5	Detection conjugate: Tris buffer + alkaline phosphatase-labeled anti-biotin antibody + preservatives (400 µL).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine (DEA*) (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

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- Pipette with disposable tip to dispense **200 µL**.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instrument of the VIDAS family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory Biosafety Manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPR[®]s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

STORAGE CONDITIONS

- Store the VIDAS Anti-HBs Total II kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

SPECIMENS**Specimen type and collection:**

Human serum or plasma.

Frozen-stored samples: after thawing, these samples must be homogenized before analysis. Mix using a vortex-type mixer. If necessary, clarify the samples by centrifugation before analysis.

Types of tubes validated:

- Silicone coated tube
- PET tube with clot activator
- PET tube with separator gel
- Tube with lithium heparin
- Tube with sodium heparin
- Tube with lithium heparin and separator gel

Follow the tube manufacturer's recommendations for use.

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- albumin up to 60 g/L
- rheumatoid factor up to 800 IU/mL
- hemolysis, after spiking samples with purified hemoglobin up to 300 µmol/L or 500 mg/dL (monomer),
- lipemia (after spiking samples with lipids up to 30 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin up to 510 µmol/L or 30 mg/dL).

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Do not inactivate samples.**Specimen stability**

The stability of samples separated from the clot has been validated using samples in stoppered tubes:

- at 18-25°C for 8 hours
- at 2-8°C for 5 days.

If longer storage is required, freeze the sera at -25 ± 6°C.

Three freeze/thaw cycles have been validated over a period of six months.

These studies did not reveal any impact on the quality of results.

INSTRUCTIONS FOR USE

**For complete instructions, see the User Manual.
Reading VIDAS® Protocole Test Change (PTC)
protocol data and MLE data**

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Calibration

Calibration, using the calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered, and then every 28 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested in duplicate. The calibration values must be within the set RFV ("Relative Fluorescence Value"). If this is not the case, recalibrate.

Procedure

1. **Remove the required reagents from the refrigerator. They can be used immediately.**
2. Use one "AHBS" strip and one "AHBS" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "AHBS" code on the instrument. The calibrator must be identified by "S1" and tested **in duplicate**. The positive control is identified by "C1" and tested singly. The negative control is identified by "C2" and tested singly.
4. Mix the calibrator, controls, and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

5. For this test, the calibrator, control, and sample test portion is 200 µL.

6. Insert the "AHBS" SPRs and "AHBS" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
7. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
8. Reclose the vials and return them to the recommended temperature after pipetting.

9. The assay will be completed within approximately 60 minutes. After the assay is completed, remove the SPRs and strips from the instrument.

10. Dispose of the used SPRs and strips into an appropriate recipient.

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The results are automatically calculated using calibration curves which are stored by the instrument (4-parameter logistics model) and are expressed in mIU/mL. The measuring range is 3-500 mIU/mL.

VIDAS Anti-HBs Total II is calibrated against the international standard (WHO Second International Standard 07/164 for anti-hepatitis B surface antigen (Anti-HBs) immunoglobulin).

Interpretation of results is as follows:

TITER (mIU/mL)	INTERPRETATION
Titer < 8	Negative
8 ≤ Titer < 12	Equivocal
Titer ≥ 12	Positive

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Specimen dilution

Only dilute specimens which, when first tested, gave a positive result > 500 mIU/mL with VIDAS Anti-HBs Total II.

The specimens should be diluted in the C2 negative control. In most cases, a 1:10 dilution (1 volume of sample + 9 volumes of negative control) enables a result to be obtained which is within the measurement range.

If the dilution factor has not been entered when the Work List was created (see User Manual), multiply the result by the dilution factor to obtain the sample concentration.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS Anti-HBs Total II kit.

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Results cannot be validated if the control values deviate from the expected values. All the samples tested since the last valid control will have to be retested.

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LIMITATIONS OF THE METHOD

Recent studies have shown that the Anti-HBs measuring range following vaccination depends on both the vaccine and the test kit (5).

In HBV infected patients, the simultaneous presence of HBsAg and Anti-HBs can affect the detection of Anti-HBs. Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's history and the results of any other tests performed.

In the context of vaccination, samples with concentrations <10 mIU/mL are considered to be non-reactive (negative) and samples with concentrations >10 mIU/mL are considered to be reactive (positive).

Nevertheless, since the VIDAS Anti-HBs Total II assay has an equivocal zone, it is recommended to refer to the results and interpretation section of this package insert.

Given the diversity of antibodies, the results obtained can differ depending on the assay used. If assays from different manufacturers are used with the same sample, the results can vary as much as four-fold (in rare cases as much as ten-fold). If the test system is changed during vaccination follow-up, the Anti-HBs level must be determined in parallel on the two systems during the transitional period (7).

PREVALENCE

The prevalence (8) of Anti-HBs varies depending on the endemic disease and vaccination policies.

Region of the world	Anti-HBs (%)
Northern, eastern, and central Europe, North America, Australia.	4 - 6
Eastern Europe, Mediterranean basin, Russia, Southwest Asia, South America.	20 - 55
Southeast Asia, tropical Africa.	70 - 95

Vaccination produces anamnestic responses in approximately 90% of vaccinated individuals, if they are healthy adults less than 40 years of age. The factors associated with a waning antibody response are age, smoking, obesity, HIV infection or chronic disease (9).

The following table indicates the seroconversion rates which may be observed according to age (10).

Neonates	> 95%
2-19 years	~ 99%
20-29 years	~ 95%
30-39 years	~ 90%
40-49 years	~ 85%
50-59 years	~ 70%
>59 years	~ 50%
Renal failure, HIV infection, other immunosuppression	50-70%
Liver disease	60-70%

PERFORMANCE

Diagnostic Sensitivity and Specificity

The VIDAS Anti-HBs Total II assay performance was established using characterized samples. The status of samples was defined by the agreement of results obtained using three commercially available EIAs (CE marked).

The samples with indeterminate status (discordant results between the three systems) were excluded from the analysis (82 samples or 3.7% of the samples tested).

A total of 2110 samples were tested and analyzed at two laboratories: 641 blood donor samples and 1469 clinical samples. Among these samples, 1041 had a positive status (including 116 samples from naturally infected patients and 925 from patients vaccinated against hepatitis B) and 1069 had a negative status.

The following results were obtained:

VIDAS Anti-HBs Total II results	Sample Status		
	Positive	Negative	Total
Positive	1019	11	1030
Negative	10	1043	1053
Total	1029	1054	2083*

* Equivocal results obtained with VIDAS Anti-HBs Total II excluded from the analysis (34 samples or 1.6% of the samples tested).

The diagnostic specificity of VIDAS Anti-HBs Total II for this population is: **99.0%** (95% confidence interval: [98.1 – 99.5]%).

The diagnostic sensitivity of VIDAS Anti-HBs Total II for this population is: **99.0%** (95% confidence interval: [98.2 – 99.5]%).

Sensitivity for seroconversion panels

During a study, twelve commercial seroconversion panels were detected by the device and the results were comparable to those claimed by the manufacturers. They consisted of 11 follow up panels and 1 seroconversion panel.

Limits of detection and quantitation

The Limit of Detection (LoD), corresponding to the concentration of Anti-HBs in a sample likely to be distinguished from the blank sample with a 95% probability, is equal to 1.03 mIU/mL. The Limit of Quantitation (LoQ), corresponding to the lowest concentration of anti-HBs that can be quantified with an acceptable level of precision and accuracy, is equal to 2.58 mIU/mL. The study was performed as recommended in the CLSI® EP17-A document.

Linearity

The VIDAS Anti-HBs total II assay is linear over the claimed measurement range, according to the study performed following the recommendations of the CLSI® EP6-A document.

Hook effect

No hook effect was observed up to Anti-HBs concentrations of 75000 mIU/mL.

Precision

The repeatability and reproducibility of the VIDAS Anti-HBs Total II assay were evaluated according to the recommendations of the CLSI® EP5-A2 document. Four human samples were tested in duplicate using two lots of reagents (N = 240). Testing was performed twice a day for 10 days on 3 instruments. Two calibrations were used for each lot (5 test days per calibration and per lot). The combined results of this study are summarized in the following table:

Sample	Mean measured value (mIU/mL)	Repeatability		Between-lot; Between-system Reproducibility	
		Standard deviation	CV (%)	Standard deviation	CV (%)
Sample 1	5.72	0.21	3.7	0.46	8.0
Sample 2	16.82	0.55	3.3	1.37	8.1
Sample 3	241.00	10.03	4.2	16.65	6.9
Sample 4	324.80	12.88	4.0	20.95	6.5

Cross-reactivity

Cross-reactivity is based on the study of samples that are negative with the assay being evaluated and positive for the potentially interfering physiological state. The results of the potentially interfering samples tested are shown in the table below:

	Number of patients	VIDAS Anti-HBsT II positive results % cross-reactivity
Anti-nuclear antibodies	20	0 (0%)
Anti Pichia	11	0 (0%)
CMV IgG	10	0 (0%)
EBV IgG	13	0 (0%)
FR	16	0 (0%)
HAMA	10	0 (0%)
HAV T	12	0 (0%)
HBc T	16	0 (0%)
HCV	21	0 (0%)
HEV	3	0 (0%)
HIV	15	0 (0%)
HSV IgG	15	0 (0%)
Rub IgG	15	0 (0%)
Syphilis IgG	11	0 (0%)
Toxoplasma IgG	12	0 (0%)
Multiparous pregnant women	10	0 (0%)
Non-viral liver disease	3	0 (0%)

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LITERATURE REFERENCES

1. WHO : Weekly epidemiological record., 2009, October 2009. No. 40 Vol 84, p405–420. JIA-HORNG KAO., DING-SHINN CHEN. – Global control of hepatitis B virus infection – The Lancet infectious diseases July 2002 Vol 2, p.395-403.
2. HOLLINGER FB. - Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. KANE M et al. – Are booster immunisations needed for lifelong hepatitis B immunity – The Lancet 2000, vol. 355, 561-565.
4. WEST DJ., CALANDRA GB. – Vaccine induced immunologic memory for hepatitis B surface antigen : implications for policy on booster vaccination – Vaccine 1996, vol. 14, 1019-1027.
5. HEIJTINK RA., SCHNEEBERGER PM., POSTMA B et al. – Anti-HBs levels after hepatitis B immunisation depend on test reagents : routinely determined 10 and 100 IU/ml seroprotection levels unreliable -Vaccine 20, 2002, 2899-2905.
6. OTTONE S., NGUYEN X., BAZIN J, BÉRARD C, JIMENEZ S and LETOURNEUR O. Expression of hepatitis B surface antigen major subtypes in Pichia pastoris and purification for in vitro diagnosis. Protein Expr Purif. 2007 Dec;56(2):177
7. HUZLY D., SCHENK T., JILG W. and NEUMANN-HAEFELIN D. Comparison of Nine Commercially Available Assays for Quantification of Antibody Response to Hepatitis B Virus Surface Antigen. Journal of Clinical Microbiology. 2008 apr ;Vol. 46, N°4 : 1298-306
8. ZUCKERMAN AJ., Hepatitis Viruses. In: Baron S, eds. Medical Microbiology, 4th ed. The University of Texas M Branch at Galveston, 1996 : 849-863.
9. STENZEL M. - Hepatitis B : the only vaccine-preventable bloodborne pathogen – Medecine and Health, 2000, vol.83, n°7, 201-204.
10. ROBINSON WS., Hepatitis B virus and hepatitis D virus. In: Mandell GL, Bennett JE, and Dolin R. eds. Principles and Practice of Infectious Diseases, 4th ed. New York. Churchill Livingstone, 1995: 1406-1439.

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Release date	Part Number	Change Type	Change Summary
2015/01	9304018E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) WARNINGS AND PRECAUTIONS
2015/09	9304018F	Technical	CONTENT OF THE KIT (60 TESTS) INSTRUCTIONS FOR USE
2016/10	9304018G	Technical	PERFORMANCE

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1. WHO : Weekly epidemiological record., 2009, October 2009. No. 40 Vol 84, p405–420. JIA-HORNG KAO., DING-SHINN CHEN. – Global control of hepatitis B virus infection – The Lancet infectious diseases July 2002 Vol 2, p.395-403.
2. HOLLINGER FB. - Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. KANE M et al. – Are booster immunisations needed for lifelong hepatitis B immunity – The Lancet 2000, vol. 355, 561-565.
4. WEST DJ., CALANDRA GB. – Vaccine induced immunologic memory for hepatitis B surface antigen: implications for policy on booster vaccination – Vaccine 1996, vol. 14, 1019-1027.
5. HEIJTINK RA., SCHNEEBERGER PM., POSTMA B et al. – Anti-HBs levels after hepatitis B immunisation depend on test reagents : routinely determined 10 and 100 IU/ml seroprotection levels unreliable -Vaccine 20, 2002, 2899-2905.
6. OTTONE S., NGUYEN X., BAZIN J, BÉRARD C, JIMENEZ S and LETOURNEUR O. Expression of hepatitis B surface antigen major subtypes in Pichia pastoris and purification for in vitro diagnosis. Protein Expr Purif. 2007 Dec;56(2):177

7. HUZLY D., SCHENK T., JILG W. and NEUMANN-HAEFELIN D. Comparison of Nine Commercially Available Assays for Quantification of Antibody Response to Hepatitis B Virus Surface Antigen. Journal of Clinical Microbiology. 2008 apr ;Vol. 46, N°4 : 1298-306
8. ZUCKERMAN AJ., Hepatitis Viruses. In: Baron S, eds. Medical Microbiology, 4th ed. The University of Texas M Branch at Galveston, 1996 : 849-863.
9. STENZEL M. - Hepatitis B : the only vaccine-preventable bloodborne pathogen – Medecine and Health, 2000, vol.83, n°7, 201-204.
10. ROBINSON WS., Hepatitis B virus and hepatitis D virus. In: Mandell GL, Bennett JE, and Dolin R. eds. Principles and Practice of Infectious Diseases, 4th ed. New York. Churchill Livingstone, 1995: 1406-1439.

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VIDAS[®] Lyme IgM (LYM) - VIDAS[®] Lyme IgG (LYG)

IVD

The VIDAS Lyme IgM (LYM) (ref. 30 319) assay is an automated qualitative test intended for use on the VIDAS family of instruments, for the detection of IgM antibodies to *Borrelia burgdorferi sensu lato* (sl) in human **serum or plasma**, to aid in the diagnosis of Lyme disease.

The VIDAS Lyme IgG (LYG) (ref. 30 320) assay is an automated qualitative test intended for use on the VIDAS family of instruments, for the detection of IgG antibodies to *Borrelia burgdorferi sensu lato* (sl) in human **serum or plasma**, to aid in the diagnosis of Lyme disease.

The VIDAS Lyme IgG (LYGS) (ref. 30 320) assay is an automated qualitative test intended for use on the VIDAS family of instruments, for the detection of IgG antibodies to *Borrelia burgdorferi sensu lato* (sl) in **cerebrospinal fluid (CSF)**, in order to determine intrathecal antibody production, to aid in the diagnosis of neuroborreliosis.

The principle of detection of the VIDAS Lyme IgM and VIDAS Lyme IgG assays is based on the ELFA (Enzyme Linked Fluorescent Assay) technique.

SUMMARY AND EXPLANATION

Lyme disease is an infectious disease caused by different strains of *Borrelia burgdorferi* sl (*Borrelia burgdorferi sensu stricto*, *Borrelia afzelii*, *Borrelia garinii*), and it is transmitted to humans through the bite of an infected tick. The main vectors are *Ixodes ricinus* in Europe, *Ixodes pacificus* and *Ixodes scapularis* in the United States (1, 2).

The disease is seasonal (from spring to autumn). It occurs in three stages and is best described as a complex evolution from early disease to late disease (when left untreated) (2, 3):

- early localized disease (days to weeks) (stage I): essentially erythema migrans, generally accompanied with flu-like symptoms.
- early disseminated disease (weeks to months) (stage II):
 - neurological manifestations: early neuroborreliosis,
 - rheumatological manifestations: "Lyme arthritis",
 - dermatological, cardiac or ophthalmological (more rare) manifestations.
- late persistent disease (months to years) (stage III):
 - neurological manifestations: late neuroborreliosis,
 - dermatological manifestations: acrodermatitis chronica atrophicans (ACA),
 - rheumatological manifestations.

In Europe, neuroborreliosis is reported in about 15% of patients with Lyme borreliosis (2). The clinical manifestations are varied and non-specific: facial palsy, radiculoneuritis, sensory and/or motor disorders, headache, etc. The specific detection of antibodies to *Borrelia burgdorferi* in CSF is an additional argument in favor of the diagnosis of neuroborreliosis, and enables the determination of the specific intrathecal IgG antibody production index. This index is used to differentiate between passive diffusion of *Borrelia burgdorferi*-specific IgG from serum to CSF and *in situ* formation in CSF (4).

The determination of this parameter is one of the criteria retained by the Consensus Conference for the diagnosis of neuroborreliosis (2).

Diagnosis

Detection of spirochetes from biopsies (cutaneous, synovial) or biological fluids (blood, CSF or synovial fluid), either by direct staining, culture or gene amplification techniques is extremely difficult and therefore reserved for specialized laboratories.

Testing for antibodies directed against *Borrelia burgdorferi sensu lato* is the most common method for biological detection of infection.

The VIDAS Lyme IgM assay uses a Solid Phase Receptacle (SPR[®]) sensitized with a combination of recombinant chimeric proteins corresponding to DbpA and OspC protein sequences of *Borrelia burgdorferi sensu stricto* B31, *Borrelia burgdorferi afzelii* PKo and *Borrelia burgdorferi garinii* PBi and PEi and 40.

The VIDAS Lyme IgG assay uses a Solid Phase Receptacle (SPR[®]) sensitized with a combination of recombinant chimeric proteins corresponding to VlsE, DbpA and OspC protein sequences of *Borrelia burgdorferi sensu stricto* B31, *Borrelia burgdorferi afzelii* PKo and ACA-1 and *Borrelia burgdorferi garinii* PBi and Ip90.

The 16th Consensus Conference on Anti-Infective Therapy – Lyme Borreliosis: diagnosis, treatment and prevention, the European Union Concerted Action On Lyme Borreliosis (EUCALB) and the Center for Disease Control (CDC) recommend laboratories use a two-test approach for the serologic diagnosis of Lyme borreliosis (2, 4, 5):

- detection of specific antibodies should be performed first using an enzyme immunoassay screening technique (ELISA).
- positive or equivocal specimens should then be tested using an immunoblotting confirmation technique (Western blot). Sensitivity and specificity of the ELISA and Western Blot tests depend on the clinical evolution. Interpretation of results should take into account patient history and the date of infection.

Treatment

Early antibiotic treatment of Lyme disease (β-lactams, cyclines, macrolides) can resolve clinical symptoms and prevent progression of the disease to later stages (2, 3).

PRINCIPLE

The VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) assay principle combines a 2-step enzyme immunoassay sandwich method with a final fluorescent detection (ELFA) (see User Manual).

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All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR® several times.

After a preliminary wash step and a sample dilution step, the antibodies to *B. burgdorferi* present in the sample will bind to the *B. burgdorferi* *sl* specific recombinant proteins coating the interior of the SPR®.

Unbound sample components are washed away. Anti-human IgM antibodies (for the VIDAS Lyme IgM assay) or anti-human IgG antibodies (for the VIDAS Lyme IgG assay) conjugated with alkaline phosphatase, will attach to the immune complex bound to the SPR® wall.

A final wash step removes unbound conjugate. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR®. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the quantity of anti-*B. burgdorferi* IgM antibody (for the VIDAS Lyme IgM assay) or IgG antibody (for the VIDAS Lyme IgG assay) present in the sample.

At the end of the VIDAS Lyme IgM and VIDAS Lyme IgG assays, results are automatically calculated by the instrument. A test value is generated and a report is printed.

CONTENT OF THE VIDAS LYME IgM KIT (60 TESTS):

60 LYM Strips	STR	Ready-to-use.
60 LYM SPR®s 2 x 30	SPR®	Ready-to-use. SPR®s coated with a combination of recombinant chimeric proteins (<i>DbpA</i> and <i>OspC</i>)* of <i>Borrelia burgdorferi sensu lato</i> .
LYM Standard 1 x 1.4 mL (liquid)	S1	Ready-to-use. Defibrinated human plasma** containing anti- <i>Borrelia burgdorferi</i> IgM in phosphate buffer + stabilizers + bovine albumin (BSA) + MIT 1g/L pH 7.4 MLE data indicate the confidence interval in Relative Fluorescence Value (RFV) ("Standard (S1) RFV Range").
LYM Positive Control 1 x 0.7 mL (liquid)	C1	Ready-to-use. Defibrinated human plasma** containing anti- <i>Borrelia burgdorferi</i> IgM in phosphate buffer + stabilizers + BSA + MIT 1g/L pH 7.4. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
LYM Negative Control 1 x 0.7 mL (liquid)	C2	Ready-to-use. Phosphate buffer + BSA, with no detectable anti- <i>Borrelia burgdorferi</i> IgM + MIT 1g/L pH 7.4. MLE data indicate the index: confidence interval ("Control C2 (-) Test Value Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit. or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

* bioMérieux Patent Pending WO/2011/023909

**This product has been tested and shown to be negative for HBs antigen, and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

CONTENT OF THE VIDAS LYME IgG KIT (60 TESTS):

60 LYG Strips	STR	Ready-to-use.
60 LYG SPR®s 2 x 30	SPR®	Ready-to-use. SPR®s coated with a combination of recombinant chimeric proteins (<i>VisE*</i> , <i>DbpA</i> and <i>OspC**</i>) of <i>Borrelia burgdorferi sensu lato</i> .
LYG Standard 1 x 1.4 mL (liquid)	S1	Ready-to-use. Anti- <i>Borrelia burgdorferi</i> IgG in phosphate buffer + BSA + MIT 1g/L pH 7.2. MLE data indicate the confidence interval in Relative Fluorescence Value (RFV) ("Standard (S1) RFV Range").
LYG Positive Control 1 x 0.7 mL (liquid)	C1	Ready-to-use. Anti- <i>Borrelia burgdorferi</i> IgG in phosphate buffer + BSA + MIT 1g/L pH 7.2. MLE data indicate the index: confidence interval ("Control C1 (+) Test Value Range").
LYG Negative Control 1 x 0.7 mL (liquid)	C2	Ready-to-use. Phosphate buffer + BSA, with no detectable anti- <i>Borrelia burgdorferi</i> IgG + MIT 1g/L pH 7.2. MLE data indicate the index: confidence interval ("Control C2 (-) Test Value Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit. or • MLE bar code printed on the box label.		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

* bioMérieux Patent Pending WO/2011/023914

** bioMérieux Patent Pending WO/2011/023909

The SPR®

For the VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) kits, the SPR® is coated during production with a combination of 3 *B. burgdorferi* *sl* markers. Each SPR® is identified by the LYM code for the VIDAS Lyme IgM kit or the LYG code for the VIDAS Lyme IgG kit. Only remove the required number of SPR®s from the pouch and **carefully reseal the pouch after opening**.

The Reagent Strip

For the VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) kits, each reagent strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the LYM Reagent Strip:

Well	Reagents
1	Sample well.
2	Sample diluent: TRIS buffer, NaCl + neutralizer + detergent + BSA + MIT 1 g/L pH 7.4 (400 µL).
3 - 4 - 5	Wash buffer: TRIS buffer, NaCl + detergent + MIT 1g/L pH 7.8 (600 µL).
6	Conjugate: TRIS buffer, NaCl + BSA + MIT 1g/L + titrated solution of mouse monoclonal anti-human IgM antibodies conjugated to alkaline phosphatase pH 7.4 (400 µL).
7 - 8	Wash buffer: TRIS buffer, NaCl + detergent + MIT 1g/L pH 7.8 (600 µL).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/l/L) + diethanolamine DEA* (0.62 mol/l/L or 6.6%) pH 9.2 + 1 g/L sodium azide (300 µL).

* Signal Word: **DANGER****Hazard statement**

H318: Causes serious eye damage.

Precautionary statement

P280: Wear protective gloves/protective clothing/eye protection/face protection.

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For further information, please refer to the Material Safety Data Sheet.

Description of the LYG Reagent Strip:

Well	Reagents
1	Sample well.
2	Sample diluent: Phosphate buffer+ NaCl + neutralizer + detergent + BSA + MIT 1 g/L pH 7.2 (600 µL).
3 - 4 - 5	Wash buffer: TRIS buffer, NaCl + detergent + MIT 1g/L pH 7.8 (600 µL).
6	Conjugate: Phosphate buffer, NaCl + bovine proteins + titrated solution of mouse monoclonal anti-human IgG antibodies conjugated to alkaline phosphatase + MIT 1g/L pH 6.1 (400 µL).
7 - 8	Wash buffer: TRIS buffer, NaCl + detergent + MIT 1g/L pH 7.8 (600 µL).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl-phosphate (0.6 mmol/l/L) + diethanolamine DEA* (0.62 mol/l/L or 6.6%) pH 9.2 + 1 g/L sodium azide (300 µL).

* Signal Word: **DANGER****Hazard statement**

H318: Causes serious eye damage.

Precautionary statement

P280: Wear protective gloves/protective clothing/eye protection/face protection.

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For further information, please refer to the Material Safety Data Sheet.

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- Pipette with disposable tips to dispense 100 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, refer to the Instrument User Manual.
- Instrument of the VIDAS® family.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **The VIDAS Lyme IgM kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest; do not inhale).**
- The VIDAS Lyme IgM and VIDAS Lyme IgG kits contain products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest; do not inhale).

- Consider all patient specimens as potentially infectious and observe routine biosafety precautions. Dispose of all used components and other contaminated materials by acceptable procedures for potentially biohazardous human blood products (6, 7, 8).
- Although the recombinant *B. burgdorferi* proteins coating the SPR®s have been inactivated, handle the SPR®s as if they were infectious.
- Do not use the SPR®s if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the kit label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests (9).
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User Manual).

STORAGE CONDITIONS

- Store the VIDAS Lyme IgM and VIDAS Lyme IgG kits at 2–8° C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2–8°C.**
- After opening the kit, check that the SPR® pouch is correctly sealed and undamaged. If not, do not use the SPR®s.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPR®s and return the complete kit to 2–8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label.

SPECIMEN PREPARATION

Sample type and collection:

- Serum or plasma for the VIDAS Lyme IgG (LYG) and VIDAS Lyme IgM (LYM) assays.
- CSF for the VIDAS Lyme IgG (LYGS) assay.

Types of tubes tested:

- Dry tube,
- Tube with lithium heparin or sodium heparin,
- Tube with separation gel.

Do not use heat-inactivated or contaminated sera. Clarify sera or plasma containing particulate matter by centrifugation or filtration prior to use.

Pre-analytical preparation of CSF must follow the official recommendations (10).

- Collect the fluid by lumbar puncture on the same day as the blood sample.
- If necessary, clarify the samples by centrifugation prior to testing.
- Do not use grossly bloody CSF, samples containing particulate matter, or clearly contaminated samples.
- A minimum sample volume of 100 µL is required.

Sample-related interference (serum or plasma)

It is recommended not to use samples that are clearly hemolyzed, lipemic or icteric and, if possible, to collect a new sample.

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- hemolysis (after spiking samples with hemoglobin: 5 g/L (monomer)),
- lipemia (after spiking samples with lipids: 30 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0.3 g/L).

Sample stability

- Serum, plasma

If the samples cannot be tested on the day of collection, store them at 2–8°C in stoppered tubes for up to 6 days. If longer storage is required, freeze the sera at -25 ± 6°C. A study performed on samples that had been frozen for six months showed that the quality of results is not affected. Two freeze/thaw cycles have been validated.

- Cerebrospinal fluid

Prior to testing, store the samples at 2–8°C, for 48 hours maximum.

If longer storage is required, freeze the CSF samples at -25 ± 6°C or at -80°C if possible (10).

INSTRUCTIONS FOR USE

For complete instructions, see the Instrument User Manual.

VIDAS® PTC protocol data entry

When using the VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) assays for the first time and **before reading the MLE data**, scan the bar code(s) (at the end of the package insert) using the external bar code reader. This reading allows VIDAS PTC protocol data to be transferred to the instrument software for its update. These data should only be read the first time the assay is used.

Types of assay and protocols

	Lyme IgM assay	Lyme IgG assay	
Sample	Serum or plasma	Serum or plasma	CSF
Assay code	LYM	LYG	LYGS
Strips	LYM	LYG	LYG
SPR®s	LYM	LYG	LYG
Calibration	LYM	LYG	LYGS

Master lot data entry

Note: When using the VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) assays for the first time, enter the VIDAS® PTC protocol (bar codes at the end of the package insert) before reading the MLE data. If the MLE data have been read before the VIDAS PTC protocol, read the MLE data again.

Before each new lot of VIDAS Lyme IgM (ref. 30 319) and VIDAS Lyme IgG (ref. 30 320) reagents is used, specifications (or factory master data) must be entered into the instrument using the master lot entry (MLE) data.

If this operation is not performed **before initiating the tests**, the instrument will not be able to print results. The master lot data need only be entered once for each lot.

It is possible to enter MLE data manually or automatically depending on the instrument (refer to the User Manual).

Calibration of VIDAS Lyme IgM (LYM), VIDAS Lyme IgG with serum or plasma (LYG) and VIDAS Lyme IgG with CSF (LYGS)

For the VIDAS Lyme IgM or VIDAS Lyme IgG assay, calibration must be performed each time a new lot is opened using the standard provided in the kit, after the master lot data have been entered. Calibration should then be performed every 28 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The standard, identified by S1, must be tested in **duplicate** (see the User Manual). The standard value must be within the set RFV (Relative Fluorescence Value) range. If this is not the case, recalibrate.

The expected values for the standard are also stored using the MLE data. If the standard values deviate from the expected values, they will be flagged on the result sheet and the result should not be reported.

The "LYG" Standard S1 is used for the calibration of the "LYG" and "LYGS" assays.

Procedure for VIDAS Lyme IgM (LYM), VIDAS Lyme IgG with serum or plasma (LYG) and VIDAS Lyme IgG with CSF (LYGS)

1. Remove the required reagents from the refrigerator. They can be used immediately.
2. Use one "LYM" strip and one "LYM" SPR® for each serum or plasma sample, control or standard to be tested.
Use one "LYG" strip and one "LYG" SPR® for each serum or plasma sample, control or standard to be tested.
Use one "LYG" strip and one "LYG" SPR® for each CSF sample, control or standard to be tested. **In this case, the test code is LYGS.**
3. **Make sure the storage pouch has been carefully resealed after the required SPR®s have been removed.**
4. The test is identified by the "LYM", "LYG" or "LYGS" code on the instrument. If a standard is to be tested, enter "S1" for the sample ID.
5. The standard identified by "S1", must be tested in **duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by "C2".

6. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet).

7. For this test, the standard, control and sample test portion is 100 µL.

8. The **LYM and LYG assays with serum or plasma** can be performed **in the same section** in the instrument. However, the LYGS assay for the detection of IgG in CSF must be performed in a **different section** to the LYG assay for the detection of IgG in serum or plasma.
9. Insert the SPR®s and the strips into the instrument. Check to make sure the color labels with the assay code on the SPR®s and the Reagent Strips match.
10. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
11. Reclose the vials and return them to 2-8°C after pipetting.
12. The printed results are obtained within approximately **27 minutes**. After the assay is completed, remove the SPR®s and strips from the instrument.
13. Dispose of the used SPR®s and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the VIDAS Lyme IgM (ref. 30 319) or VIDAS Lyme IgG (ref. 30 320) assay is completed, the results are analyzed automatically by the instrument. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR® is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR®. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. The RFV for the assay appears on the result sheet.

The Patient RFV result is interpreted by the instrument in the following manner:

$i = \text{index} = \text{sample RFV} / \text{standard RFV}$

This index, along with the interpretation, also appears on the result sheet. Interpretation of the index is as follows:

- **Plasma, serum**

VIDAS Lyme IgM (LYM)

Index	Interpretation for LYM
$i < 0.20$	Negative
$0.20 \leq i < 0.32$	Equivocal
$i \geq 0.32$	Positive

VIDAS Lyme IgG (LYG)

Index	Interpretation for serum or plasma
$i < 0.20$	Negative
$i \geq 0.20$	Positive

The 0.20 and 0.32 cut-offs are the results obtained after the VIDAS Lyme IgG and Lyme IgM ROC curves were studied and analyzed. The clinical status of the samples used in the studies was that obtained with a commercially available immunoblotting test (Western Blot).

The EUCLAB and the CDC today recommend performing immunoblotting (Western Blotting) for all samples that give equivocal or positive results with ELISA assays. In general, these results should be interpreted taking into consideration the clinical presentation, the case history, and the results of any other biological tests performed (4, 11). Furthermore, the following table aids in the joint interpretation of VIDAS Lyme IgM and VIDAS Lyme IgG assay results.

VIDAS Lyme IgM result	VIDAS Lyme IgG result	Possible interpretation
Negative	Negative	Unlikely infection. In case of clinical suspicion (tick bite, neurological symptoms, etc.), collect another sample later.
Positive	Negative	Probable early infection.
Positive	Positive	Probable early or disseminated infection.
Negative	Positive	Probable infection at any stage (early, disseminated or persistent).
Equivocal or positive at cut-off	Positive at cut-off	Depending on context and case history, early infection or serological scar may be considered.
Equivocal or positive at cut-off	Negative	Depending on context and case history, early infection or a potentially interfering disease may be considered.

• Cerebrospinal fluid

VIDAS Lyme IgG (LYGS)

In the neurological manifestations of Lyme Borreliosis, the specific detection of antibodies to *Borrelia burgdorferi* in CSF is an additional argument in favor of the diagnosis of neuroborreliosis. It enables the determination of the specific intrathecal IgG antibody production index (IAP) through a concomitant analysis of a blood sample and a CSF sample. (2).

Indeed, the intrathecal antibody production index takes into account the detection of anti-*Borrelia burgdorferi* IgG antibodies in the patient's serum and CSF, once the protein variations between the two samples have been corrected. This index is then used to differentiate between passive diffusion of *Borrelia burgdorferi*-specific IgG from serum to CSF and *in situ* formation in CSF (4).

IAP determination is based on the VIDAS Lyme IgG indexes obtained using the patient's serum and CSF, and the albumin (Tibbling method) or total IgG (Reiber method) titers in the same two samples (12).

Tibbling method:

$$\text{IAP} = \frac{[(\text{VIDAS Lyme IgG CSF index} / \text{VIDAS Lyme IgG serum index}) / 9^*]}{[(\text{CSF albumin titer} / \text{serum albumin titer})]}$$

Reiber method:

$$\text{IAP} = \frac{[(\text{VIDAS Lyme IgG CSF index} / \text{VIDAS Lyme IgG serum index}) / 9^*]}{[(\text{CSF total IgG titer} / \text{serum total IgG titer})]}$$

* Invariable factor, specific to the test protocol.

Note: The choice of either the albumin or total IgG test method used, and the titers obtained, are the responsibility of the laboratory. The same method should be used to titer albumin or total IgG in both the blood and the CSF samples.

An index value greater than or equal to 2, points to intrathecal *Borrelia burgdorferi* IgG antibody production. An index value less than 2 can neither confirm nor exclude intrathecal antibody production.

Limitations for using the intrathecal antibody production index:

- the calculation formula is not applicable for VIDAS Lyme IgG indexes on serum ≤ 0.00 .
- when the albumin or total IgG ratio is nearing 0 (zero), the calculation formula may generate very high results.
- it is not recommended to apply the calculation formula when the VIDAS Lyme IgG index on CSF is very close to 0.00 or when the results (in index) of the patient's serum and CSF pair are very low. In fact, the calculation could lead to a result that is ≥ 2 when there are no *Borrelia burgdorferi*-specific IgGs in the test samples.
- high index values (> 6.0) with the VIDAS Lyme IgG assay on serum may be underestimated due to a level of anti-*Borrelia burgdorferi* IgG immunoglobulins that exceeds the assay's capture capability. The actual level could be greater. Therefore, it is not recommended to apply the calculation formula when indexes on serum are greater than 6.0. In this situation, if the numerator value in the formula is overestimated, the calculation could lead to a result that is ≥ 2 , which cannot confirm intrathecal antibody production.

QUALITY CONTROL

One positive control and one negative control are included in each VIDAS Lyme IgM and VIDAS Lyme IgG kit. These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2. Results cannot be validated if the control values deviate from the expected values.

Note

It is the responsibility of the user to perform Quality Control in accordance with any applicable local regulations.

LIMITATIONS OF THE METHOD

- Antibody detection methods do not provide definitive results for establishing or ruling out diagnosis of Lyme Borreliosis.
- Negative results with the VIDAS Lyme IgM and VIDAS Lyme IgG assays does not rule out the possibility of *B. burgdorferi* infection in a patient. Patients in the early stages of infection or who have undergone antibiotic therapy, may not produce measurable IgM and IgG. Patients with clinical history and/or symptoms suggestive of Lyme borreliosis, but with negative test results, should be reported as "No detectable antibodies to *B. burgdorferi*". A second specimen should be collected 4-6 weeks later.

Note: it is estimated that in 50% of subjects in the primary stage of disease, antibody levels remain below the detectable threshold (11).

- Positive results with the VIDAS Lyme IgG and IgM assays must be interpreted with caution. Cross-reactivity may be observed with certain diseases (13, 14): Refer to section "Cross-reactivity".
- Clinical symptoms, epidemiological information and other laboratory test results must all be considered when interpreting VIDAS Lyme IgM and IgG assay results.
- Interference may be encountered with certain sera containing antibodies directed against reagent components or substances that affect the reaction. For this reason, assay results should be interpreted taking into consideration the patient's history, and the results of any other tests performed.
- Testing should be done only when exposure history, epidemiology and clinical symptoms suggest Lyme disease (14).

PERFORMANCE

Studies performed using the VIDAS Lyme IgM and VIDAS Lyme IgG assays gave the following results.

1. Performance of the VIDAS Lyme IgM and VIDAS Lyme IgG assays on serum

A. Sensitivity

This retrospective study was performed by two laboratories (site 1 and site 2), using the sera of patients at different clinical stages of Lyme borreliosis: erythema migrans, neuroborreliosis, arthritis/acrodermatitis chronica atrophicans.

The definition of the clinical stage and the final diagnosis were based on a combination of several parameters, including:

- clinical symptoms,
- evaluation of the risk of tick exposure,
- IgM and IgG serological profiles, associated with Western Blot confirmation tests for any positive or equivocal results.

The serological profiles were representative of the profiles observed routinely by each of the laboratories.

For each stage of the disease, the sensitivity of the VIDAS Lyme IgM and VIDAS Lyme IgG assays, as well as the sensitivity resulting from the combination of these two tests, were expressed with a 95% confidence interval.

	Sensitivity (% [95%CI])							
	Site 1				Site 2			
	Number N	VIDAS Lyme IgM	VIDAS Lyme IgG	VIDAS Lyme IgM + IgG	Number N	VIDAS Lyme IgM	VIDAS Lyme IgG	VIDAS Lyme IgM + IgG
Erythema migrans	49	61.2 [46.2 – 74.8]	51.0 [36.3 – 65.6]	83.7 [70.3 – 92.7]	55	47.3 [33.7 – 61.2]	58.2 [44.1 – 71.4]	76.4 [63.0 – 86.8]
Neuroborreliosis	60	58.3 [44.9 – 70.9]	90.0 [79.5 – 96.2]	90.0 [79.5 – 96.2]	70	35.7 [24.6 – 48.1]	91.4 [82.3 – 96.8]	92.9 [84.1 – 97.6]
Arthritis/acrodermatitis	47*	40.0 [22.7 – 59.4]	100.0 [92.5 – 100.0]	100.0 [88.4 – 100.0]	51	37.3 [24.1 – 51.9]	96.1 [86.5 – 99.5]	96.1 [86.5 – 99.5]

Note: samples found to be equivocal with VIDAS Lyme IgM were considered to be positive.

* among the 47 samples, only 30 samples in Arthritis/acrodermatitis stage were tested using the VIDAS Lyme IgM assay.

B. Specificity

This prospective study was performed using sera from healthy volunteers from two Lyme Borreliosis endemic regions: population of an area with high prevalence (Eastern France), population of an area with low prevalence (Southern France).

The volunteers were blood donors with no known previous history of Lyme borreliosis.

The population recruited in a high prevalence area (16.1%) consisted of volunteers between 18 and 68 years of age (64% men and 36% women).

The population recruited in the low prevalence area (3.9%) consisted of volunteers between 18 and 68 years of age (52% men and 48% women).

The samples from each population were tested using commercially available IgM and IgG Enzyme Immunoassays (EIA). Samples found to be positive or equivocal were tested using an immunoblotting confirmation technique. Each positive immunoblot sample was eliminated from the specificity calculation.

The specificity of VIDAS Lyme IgM and VIDAS Lyme IgG on each population, was expressed with a 95% confidence interval.

Specificity (% [95%CI])					
Low prevalence area			High prevalence area		
Number N	VIDAS Lyme IgM	VIDAS Lyme IgG	Number N	VIDAS Lyme IgM	VIDAS Lyme IgG
273	96.3 [93.4 ; 98.2]	99.3 [97.4 ; 99.9]	302	97.7 [95.3 ; 99.1]	99.7 [98.2 ; 99.99]

Note: samples found to be equivocal with the EIA or VIDAS Lyme IgM were considered to be positive.

C. Comparison of methods

This prospective study was performed in two laboratories (site 1 and site 2), using sera from patients for whom a Lyme IgM and IgG serologic test had been prescribed.

At each laboratory, the samples were tested in parallel using commercially available IgM and IgG EIA methods (different for each site) and the VIDAS Lyme IgM and VIDAS Lyme IgG assays. The samples found to be positive with one of these methods were confirmed using an immunoblotting confirmation technique.

VIDAS Lyme IgM results:

The concordance of positive and negative results for VIDAS Lyme IgM and the EIA IgM method at each site are expressed with a 95% confidence interval.

VIDAS Lyme IgM	EIA method at Site 1			EIA method at Site 2		
	Positive	Equivocal	Negative	Positive	Equivocal	Negative
Positive	39	0	12	38	3	12
Equivocal	4	2	8	7	5	13
Negative	18	10	284	26	21	691
Number N	377			816		
Positive concordance	N = 73 61.6% [49.5% – 72.8%]			N = 100 53.0% [42.8% – 63.1%]		
Negative concordance	N = 304 93.4% [90.0% – 95.9%]			N = 716 96.5% [94.9% – 97.7%]		

Note: samples found to be equivocal with the EIA or VIDAS Lyme IgM were considered to be positive for the calculations of concordance.

An IgM status (positive or negative) was then defined for each sample in accordance with the two-test approach recommended by the 16th Consensus Conference on Anti-Infective Therapy – Lyme Borreliosis, the EUCALB and the CDC, i.e:

- a positive IgM status was reported for all samples confirmed positive by immunoblot IgM,
- a negative IgM status was reported for all samples confirmed negative by immunoblot IgM or for which all the VIDAS Lyme IgM and EIA IgM results were negative.

The following table shows the concordance of positive and negative results for VIDAS Lyme IgM and the EIA IgM at each site with the defined IgM status.

		IgM status at Site 1		IgM status at Site 2	
		Positive	Negative	Positive	Negative
VIDAS Lyme IgM	Positive	35	16	49	2
	Equivocal	4	10	16	4
	Negative	5	307	13	709
	Number N	377		793	
	Positive concordance	N = 44 88.6% [75.4% – 96.2%]		N = 78 83.3% [73.2% - 90.8%]	
	Negative concordance	N = 333 92.2% [88.8% – 94.8%]		N = 713 99.2% [98.2% - 99.7%]	
EIA methods	Positive concordance	N = 44 88.6% [75.4% - 96.2%]		N = 78 82.1% [71.7% - 89.8%]	
	Negative concordance	N = 333 89.8% [86.0% – 92.8%]		N = 713 97.2% [95.7% - 98.3%]	

Note 1: samples found to be equivocal with the EIA or VIDAS Lyme IgM were considered to be positive for the calculations of concordance.

Note 2: samples found to be equivocal with the immunoblot were not included in the study.

VIDAS Lyme IgG results:

The concordance of positive and negative results for VIDAS Lyme IgG and the EIA IgG method at each site are expressed with a 95% confidence interval.

VIDAS Lyme IgG	EIA method at Site 1			EIA method at Site 2		
	Positive	Equivocal	Negative	Positive	Equivocal	Negative
Positive	102	4	2	160	2	3
Negative	17	16	236	28	24	599
Number N	377			816		
Positive concordance	N = 139 76.3% [68.3% – 83.1%]			N = 214 75.7% [69.4% – 81.3%]		
Negative concordance	N = 238 99.2% [97.0% – 99.9%]			N = 602 99.5% [98.6% – 99.9%]		

Note: samples found to be equivocal with the EIA or VIDAS Lyme IgG were considered to be positive for the calculations of concordance.

An IgG status (positive or negative) was then defined for each sample in accordance with the two-test approach recommended by the 16th Consensus Conference on Anti-Infective Therapy – Lyme Borreliosis, the EUCALB and the CDC, i.e:

- a positive IgG status was reported for all samples confirmed positive by immunoblot IgG,
- a negative IgG status was reported for all samples confirmed negative by immunoblot IgG or for which all the VIDAS Lyme IgG and EIA IgG results were negative.

The following table shows the concordance of positive and negative results for VIDAS Lyme IgG and the EIA IgG at each site with the defined IgG status.

		IgG status at Site 1		IgG status at Site 2	
		Positive	Negative	Positive	Negative
VIDAS Lyme IgG	Positive	82	25	133	9
	Negative	6	262	21	607
	Number N	375		770	
	Positive concordance	N = 88 93.2% [85.8% - 97.5%]		N = 154 86.4% [79.9% - 91.4%]	
	Negative concordance	N = 287 91.3% [87.4% - 94.3%]		N = 616 98.5% [97.2% - 99.3%]	
EIA methods	Positive concordance	N = 88 97.7% [92.0% - 99.7%]		N = 154 99.4% [96.4% - 99.9%]	
	Negative concordance	N = 287 82.2% [77.3% - 86.5%]		N = 616 97.4% [95.8% - 98.5%]	

Note 1: samples found to be equivocal with the EIA or VIDAS Lyme IgG were considered to be positive for the calculations of concordance.

Note 2: samples found to be equivocal with the immunoblot were not included in the study.

D. Precision

For each VIDAS Lyme IgM and VIDAS Lyme IgG assay, 3 serum samples were tested in duplicate in 40 different runs (2 runs per day) with 2 reagent lots at 3 sites (n =240).

The repeatability (intra-run precision), inter-lot reproducibility and inter-site reproducibility were calculated using this protocol based on the recommendations of the CLSI® EP5-A2 document:

VIDAS Lyme IgM		Repeatability		Inter-lot reproducibility		Inter-site reproducibility	
Sample	Mean index	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
1	0.19	0.01	3.5	0.01	6.3	0.01	6.3
2	0.37	0.01	3.1	0.02	6.2	0.02	6.2
3	1.26	0.03	2.5	0.12	9.4	0.12	9.4

VIDAS Lyme IgG		Repeatability		Inter-lot reproducibility		Inter-site reproducibility	
Sample	Mean index	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
1	0.15	0.01	8.6	0.02	15.5	0.02	15.5
2	0.26	0.01	5.4	0.02	7.7	0.02	7.7
3	2.31	0.10	4.1	0.12	5.3	0.12	5.3

E. Cross-reactivity

Cross-reactivity is based on the study of samples that are negative with the test being evaluated and positive for the potentially interfering disease. The results of the samples tested according to the disease are shown in the table below:

Positive samples	Number of samples	Number of false positive VIDAS Lyme IgM results	% Cross-reactivity
Syphilis	270	9	3.33
Leptospirosis	216	9	4.16
Rickettsiosis	112	1	0.08
<i>Helicobacter pylori</i>	143	3	2.09
Epstein Barr virus	65	2	3.07
Herpes simplex virus	98	5	5.10
Hepatitis A virus	153	6	3.92
Cytomegalovirus	34	2	5.88
HIV	20	1	5.00
Rubella	19	1	5.26
Mumps	46	1	2.17
Measles	38	1	2.63
Varicella-Zoster virus	58	1	1.72
Toxoplasmosis	26	1	3.84
Rheumatoid factor	61	1	1.63
HAMA	31	1	3.22
Antinuclear antibody	60	1	1.66
C Reactive protein	61	2	3.27
Systemic lupus erythematosus	28	0	0.00

Positive samples	Number of samples	Number of false positive VIDAS Lyme IgG results	% Cross-reactivity
Syphilis	256	1	0.39
Leptospirosis	206	6	2.91
Rickettsiosis	133	3	2.3
<i>Helicobacter pylori</i>	143	2	1.4
Epstein Barr virus	34	0	0.0
Herpes simplex virus	125	1	0.8
Hepatitis A virus	150	3	2.0
Cytomegalovirus	40	0	0.0
HIV	20	1	5.0
Rubella	19	0	0.0
Mumps	46	0	0.0
Measles	38	0	0.0
Varicella-Zoster virus	58	0	0.0
Toxoplasmosis	26	1	3.85
Rheumatoid factor	28	0	0.0
HAMA	43	0	0.0
Antinuclear antibody	60	5	8.33
C Reactive protein	61	2	3.28
Systemic lupus erythematosus	28	2	7.14

2. Performance of the VIDAS Lyme IgG assay on serum/CSF pairs

This retrospective study was performed at two laboratories using serum/CSF pairs from patients with a defined intrathecal IgG antibody production index. This index was therefore calculated using the VIDAS Lyme IgG test values and the albumin (Tibbling method) or total IgG (Reiber method) titers. The calculation formulas were the following:

Tibbling Method:

$$\text{Index} = \frac{(\text{VIDAS Lyme IgG CSF TV}^*/\text{VIDAS Lyme IgG serum TV})/9^{**}}{(\text{CSF albumin/serum albumin})}$$

Reiber Method:

$$\text{Index} = \frac{(\text{VIDAS Lyme IgG CSF TV}/\text{VIDAS Lyme IgG serum TV})/9^{**}}{(\text{total CSF IgG}/\text{total serum IgG})}$$

* TV = Test Value

** "9" is a constant that is characteristic of the VIDAS Lyme IgG test protocol.

Any pairs with an index value greater than or equal to 2 were considered to be positive with the VIDAS Lyme IgG assay (presence of intrathecal IgG antibody production in the patient). Conversely, any pairs with an index value less than 2 were considered to be negative with the VIDAS Lyme IgG assay (absence of intrathecal IgG antibody production in the patient).

A. Sensitivity

This study was performed using serum/CSF pairs from patients with neuroborreliosis and whose intrathecal antibody production had previously been proven by the laboratory's techniques.

The sensitivity of the VIDAS Lyme IgG assay was expressed with a 95% confidence interval.

	Sensitivity (% [95% Confidence Interval])			
	Site 1 (Tibbling)		Site 2 (Reiber)	
	Number	Sensitivity	Number	Sensitivity
Pairs with proven intrathecal antibody production	32	90.6% [75.0% - 98.0%]	31	74.2% [55.4% - 88.1%]

B. SPECIFICITY

This study was performed using serum/CSF pairs from patients whose intrathecal antibody production had not been proven by the laboratory's techniques.

The specificity of the VIDAS Lyme IgG assay was expressed with a 95% confidence interval.

	Specificity (% [95% Confidence Interval])			
	Site 1 (Tibbling)		Site 2 (Reiber)	
	Number	Specificity	Number	Specificity
Pairs without proven intrathecal antibody production	61	83.6% [71.9% - 91.8%]	52	84.6% [71.9% - 93.1%]

C. Precision

For VIDAS Lyme IgG, ("LYGS" test), 4 CSF samples spiked with anti-*Borrelia burgdorferi* sI IgG, at different index levels, were assayed in duplicate in 40 different runs (2 runs per day) with 2 reagent lots on 3 different instruments (n=240).

The repeatability (intra-run precision), inter-instrument reproducibility and inter-lot reproducibility were calculated using this protocol based on the recommendations of the CLSI® EP5-A document:

VIDAS LYGS		Repeatability		Inter-instrument reproducibility		Inter-lot reproducibility	
Sample	Mean index	Standard deviation	CV (%)	Standard deviation	CV (%)	Standard deviation	CV (%)
1	0.11	0.005	5.0	0.009	8.7	0.009	8.7
2	0.48	0.019	4.1	0.048	10.0	0.048	10.0
3	1.72	0.069	4.0	0.161	9.3	0.161	9.3
4	5.98	0.193	3.2	0.466	7.8	0.466	7.8

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LITERATURE REFERENCES

- BURGDORFER W. Discovery of the Lyme disease spirochete and its relation to tick vectors. Yale Jnl. Biol. Med. 1984, 57 : 515-520.
- 16ème conférence de Consensus en Thérapeutique Anti-Infectieuse – Borréliose de Lyme : démarches diagnostiques, thérapeutiques et préventives – Décembre 2006.
- WORMSER GP. et al. The Clinical Assessment, Treatment, and prevention of Lyme Disease, Human Granulocytic Anaplasmosis, and Babesiosis: Clinical Practice Guidelines by the Infectious Disease Society of America. Clin Infect Dis. 2006, 43: 1089-1134.
- EUCALB. <http://meduni09.edis.at/eucalb/>: DIAGNOSIS: Serology: Diagnostic Guidelines.
- CDC. <http://www.cdc.gov>: Diagnosis: laboratory testing.
- U.S. Department of Health and Human Services. 1999. Biosafety in Microbiological and Biomedical Laboratories, 4th ed. HHS Publication (CDC) 1999. Government Printing Office, Washington D.C.
- U.S. Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Blood borne Pathogens.
- National Committee for Clinical Laboratory Standards. 1997. Approved Standard, M29-A. Protection of Laboratory Workers from Instrument Biohazard and /Infectious Disease Transmitted by Blood, Body Fluids, and Tissue. NCCLS, Wayne Pa.
- LAMPE A.S., PEITERSE BRUINS, and EGTERVAN WISSERDERKE J.C.R. Wearing gloves as cause of false negative HIV tests. Lancet ii, 1988, 1140-1144.
- TEUNISSEN C. E., PETZOLD A., BENNETT F. et al. A consensus protocol for the standardization of cerebrospinal fluid collection and biobanking, Neurology 2009; Neurology 2009; 73(22) 1914-1922.
- CDC Recommendations Standardizations and Interpretation, CDC Conference 10/28/94.
- KAISER. R. 1993. Intrathecal synthesis of specific antibodies in neuroborreliosis. Comparison of different ELISA techniques and calculation methods. Journal of the Neurological Sciences, 118, 64-72.
- MAGNARELLI L., et al. Cross reactivity in serologic tests for Lyme disease and other spirochetal infections. J. Infec. Dis. 1987, 156:183-188.
- GOLIGHTLY, M., THOMAS J., and VICIANA A. The laboratory diagnosis of Lyme borreliosis. Laboratory Medicine. 1990, 21 : 299-304.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

REVISION HISTORY

Change type categories:

N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note:	<i>Minor typographical, grammar, and formatting changes are not included in the revision history</i>

Release date	Part Number	Change Type	Change Summary
2015/01	9301296E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical change	CONTENT OF THE <u>VIDAS LYME IgM</u> KIT (60 TESTS): CONTENT OF THE <u>VIDAS LYME IgG</u> KIT (60 TESTS): WARNINGS AND PRECAUTIONS
2016/08	9301296F	Technical change	SUMMARY AND EXPLANATION RESULTS AND INTERPRETATION PERFORMANCE

BIOMERIEUX, the blue logo, VIDAS and SPR are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

The VIDAS® Lyme IgG reagent is sold under license of PCT/US97/02952 (University of Texas System) and PCT/US1998/013551 (The Administrators of the Tulane Educational Fund) patent families.

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- BURGDORFER W. Discovery of the Lyme disease spirochete and its relation to tick vectors. Yale Jnl. Biol. Med. 1984, 57 : 515-520.
- 16ème conférence de Consensus en Thérapeutique Anti-Infectieuse – Borréliose de Lyme : démarches diagnostiques, thérapeutiques et préventives – Décembre 2006.
- WORMSER GP. et al. The Clinical Assessment, Treatment, and prevention of Lyme Disease, Human Granulocytic Anaplasmosis, and Babesiosis: Clinical Practice Guidelines by the Infectious Disease Society of America. Clin Infect Dis. 2006, 43: 1089-1134.
- EUCALB. <http://meduni09.edis.at/eucalb/>: DIAGNOSIS: Serology: Diagnostic Guidelines.
- CDC. <http://www.cdc.gov>: Diagnosis: laboratory testing.
- U.S. Department of Health and Human Services. 1999. Biosafety in Microbiological and Biomedical Laboratories, 4th ed. HHS Publication (CDC) 1999. Government Printing Office, Washington D.C.
- U.S. Department of Labor, Occupational Safety and Health Administration, 29CFR Part 1910.1030, Occupational Exposure to Blood borne Pathogens.
- National Committee for Clinical Laboratory Standards. 1997. Approved Standard, M29-A. Protection of Laboratory Workers from Instrument Biohazard and /Infectious Disease Transmitted by Blood, Body Fluids, and Tissue. NCCLS, Wayne Pa.
- LAMPE A.S., PEITERSE BRUINS, and EGTERVAN WISSERDERKE J.C.R. Wearing gloves as cause of false negative HIV tests. Lancet ii, 1988, 1140-1144.
- TEUNISSEN C. E., PETZOLD A., BENNETT F. et al. A consensus protocol for the standardization of cerebrospinal fluid collection and biobanking, Neurology 2009; Neurology 2009; 73(22) 1914-1922.
- CDC Recommendations Standardizations and Interpretation, CDC Conference 10/28/94.
- KAISER. R. 1993. Intrathecal synthesis of specific antibodies in neuroborreliosis. Comparison of different ELISA techniques and calculation methods. Journal of the Neurological Sciences, 118, 64-72.
- MAGNARELLI L., et al. Cross reactivity in serologic tests for Lyme disease and other spirochetal infections. J. Infec. Dis. 1987, 156 : 183-188.
- GOLIGHTLY, M., THOMAS J., and VICIANA A. The laboratory diagnosis of Lyme borreliosis. Laboratory Medicine. 1990, 21 : 299-304.

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
VIDAS[®] HBc IgM II (HBCM)

VIDAS HBc IgM II is an automated quantitative test for use on the VIDAS family instruments for the determination of anti-hepatitis B virus core antigen IgM (anti-HBc IgM) in human serum or plasma (EDTA, citrate and lithium heparin), using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

The hepatitis B virus is responsible for acute and chronic hepatitis infections. Acute hepatitis can be asymptomatic or present symptoms of varying severity which may progress to fulminant hepatitis in 0.1 to 0.5% of cases. Chronicity occurs in 5 to 10% of cases in adults, but in up to 90% of cases in infants following perinatal transmission. Currently, approximately 300 million people worldwide are chronic carriers of the virus (1). Chronic hepatitis may be asymptomatic or lead to liver lesions of varying severity, possibly evolving to cirrhosis, with an evolution in 5% of cases to hepatocellular carcinoma (2). The hepatitis B virus can be transmitted by parenteral or perinatal pathways or through sexual contact. Persons most at risk are health workers, drug addicts, those with multiple sexual partners, multiple transfusion or hemodialysis patients, close friends and family of an infected subject, and newborns of an infected mother (2).

CLINICAL SIGNIFICANCE

Anti-HBc IgM antibodies are indirect markers of viral replication, reflecting the host's response to intra-hepatic HBc antigen (3).

Diagnosis of acute hepatitis relies on the detection of HBs antigen and anti-HBc IgM. The level of infectivity, and infection evolution are assessed using other markers, which also enable determination of patients requiring screening. In cases of acute hepatitis, anti-HBc IgM titers are generally high (> 100 PEIU/mL), then progressively decrease over a 6-8 month period whether infection develops towards recovery or chronicity. However, total anti-HBc antibodies (mainly IgG) persist after recovery. Expression of HBs antigen can be limited in cases of fulminant hepatitis. Anti-HBc IgM can then indicate active viral replication (2).

In cases of chronic hepatitis, the appearance of anti-HBc IgM reflects hepatic cytolysis, revealing an active phase of the disease (4). As with serum HBe Ag and HBV DNA, a positive result for anti-HBc IgM is of considerable assistance in the diagnosis of chronic hepatitis, but a negative result does not necessarily exclude this possibility.

The diagnosis should therefore be made using both histological and serological criteria. In patients with a significant response, titration of anti-HBc IgM enables anti-viral treatments to be monitored by a specific marker of HBV.

Quantitative and **ultra-sensitive** tests enable patients to be monitored from the time HBs Ag appear up to the appearance of anti-HBs Ab, by titering the anti-HBc IgM, and also a distinction to be made between healthy carriers and those with active chronic hepatitis. These tests are also useful for monitoring chronic hepatitis and anti-viral therapy.

PRINCIPLE

The assay principle combines a two-step enzyme immunoassay method with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR[®]) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips.

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After dilution of the sample, the serum IgM bind with the anti-μ chain antibodies coating the interior of the SPR. Unbound serum components are washed away. IgM antibodies directed against the core antigen of the hepatitis B virus (HBc) bind to an immunocomplex composed of recombinant HBc antigen and alkaline-phosphatase-labeled monoclonal anti-HBc antibody. Unbound components are eliminated by washing steps. A fluorescent substrate, 4-Methyl-umbelliferyl phosphate is introduced into the SPR and reveals the amount of bound immunocomplex. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm.

The intensity of the fluorescence is proportional to the concentration of anti-HBc IgM present in the sample.

At the end of the assay, results are automatically calculated by the instrument in relation to the calibration curve stored in memory, and then printed out.

CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS

30 HBCM strips	STR	Ready-to-use.
30 HBCM SPRs 1 x 30	SPR	Ready-to-use. Interior of the SPR coated with anti-μ antibody (goat).
HBCM positive control 1 x 2 mL (liquid)	C1	Ready-to-use. Inactivated human* serum containing anti HBc IgM + 0.2 g/L gentamicin sulfate + 1 g/L sodium azide. MLE data indicate the confidence interval in PEIU/mL (Paul Erlich Institute unit per milliliter) ("Control C1 (+) Dose Value Range").
Negative control 1 x 1.9 mL (liquid)	C2	Ready-to-use. Phosphate buffer + protein stabilizer of animal origin+ preservatives.
HBCM calibrator 2 x 1 mL (lyophilized)	S1	Inactivated human* serum containing anti HBc IgM + 0.2 g/L gentamicin sulfate + preservatives. The calibrator must be reconstituted with 1 mL of distilled sterile water (measured exactly). Allow to dissolve for at least 15 min and then mix using a vortex. After reconstitution, the calibrator can be stored at 2-8°C for up to 4 weeks. Aliquoted and frozen at -25 ± 6°C, the calibrator can be kept for up to 6 months. Avoid repeated freeze/thaw cycles.
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar codes printed on the box label.		
1 Package Insert provided in the kit or downloadable from www.biomerieux.com/techlib .		

* This product has been tested and shown to be negative for antibodies to HIV1, HIV2 and HCV. It has been heat inactivated. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious and usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with anti-human μ chain antibody. Each SPR is identified by the HBCM code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening**.

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Description of the HBCM Strip

Wells	Reagents
1	Sample well
2	Sample diluent: buffer containing Tween 20 (pH = 7.4) + protein and chemical stabilizers + 1 g/L sodium azide + products of human origin (400 μL).
3 - 4 - 5 - 8 - 9	Wash buffer: TRIS buffer (50 mmol/L) (pH = 7.4) + protein and chemical stabilizers + 0.9 g/L sodium azide (600 μL).
6	Conjugate: alkaline phosphatase-labeled anti-HBc monoclonal antibody (mouse) + 1 g/L sodium azide + inactivated negative human** serum (300 μL).
7	HBc recombinant antigen (<i>E. Coli</i>) + 1 g/L sodium azide + inactivated negative human** serum (300 μL).
10	Cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine*** (DEA) (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 μL).

** This product has been tested and shown to be negative for HBs antigen, antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

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- Pipette with disposable tip to dispense 1 mL and 100 µL.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the label.
- Do not mix reagents (or disposables) from different lots.
- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent and a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS HBc IgM II kit at 2-8°C.
- **Do not freeze reagents, with the exception of the calibrator after reconstitution.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **Carefully reseal the pouch with the desiccant inside after use to maintain stability of the SPRs and return the complete kit to 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

SPECIMENS

Specimen Type and Collection:

Use human sera (plain tube) or plasma (validated anticoagulants: EDTA, citrate, lithium heparin). Sera containing impurities should be clarified by centrifugation or filtration before analysis.

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- hemolysis (after spiking samples with hemoglobin: 0 to 320 µmol/L (monomer)),
- lipemia (after spiking samples with lipids: 0 to 5 mg/mL equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin: 0 to 560 µmol/L).

However, it is recommended not to use samples that are clearly hemolyzed or lipemic and, if possible, to collect a new sample.

Do not inactivate samples.

Specimen Stability

Samples can be stored at 2-8°C in stoppered tubes for up to 7 days; if longer storage is required, freeze the sera or plasma at -25 ± 6°C. Avoid successive freezing and thawing.

A study performed on frozen samples over a period of 2 months, showed that the quality of results is not affected.

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Calibration

Calibration, using the calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the master lot data have been entered. Calibration should then be performed every 14 days. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf-life of the kit.

The calibrator, identified by S1, must be tested **in duplicate** (see User's Manual). The calibrator value must be within the set RFV "Relative Fluorescence Value" range. If this is not the case, recalibrate.

Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HBCM" strip and one "HBCM" SPR for each sample, control or calibrator to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HBCM" code on the instrument. Indicate the number of determinations to perform. The calibrator must be identified by "S1", and tested in **duplicate**. If the positive control is to be tested, it should be identified by "C1". If the negative control needs to be tested, it should be identified by C2.
4. If necessary, clarify samples by centrifugation.
5. Mix the calibrator, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet) (5).
6. **For this test, the calibrator, control, and sample test portion is 100 µL.**
7. Insert the "HBCM" SPRs and "HBCM" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
9. Reclose the vials and return them to the required temperature after pipetting.
10. The assay will be completed within approximately 55 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each

sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after incubating the substrate with the enzyme remaining on the interior of the SPR.

The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (4-parameter logistics model). The concentrations are expressed in PEIU/mL according to the anti-HBc IgM standard of the Paul Erlich Institute (n°84).

Test Value (PEIU/mL)	Interpretation
< 5	Negative
≥ 5 and < 10	Equivocal *
≥ 10	Positive

* It is advised to follow up the evolution of anti-HBc IgM titers in this zone.

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QUALITY CONTROL

One positive control and one negative control are included in each VIDAS HBc IgM II kit.

These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1 and C2. Results cannot be validated if the control values deviate from the expected values.

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PERFORMANCE

Studies performed using VIDAS HBc IgM II gave the following results:

Measurement Range

The measurement range of the VIDAS HBc IgM II kit extends up to 200 PEIU/mL.

Sensitivity - Specificity for Clinical Population.

Sensitivity and specificity were determined by comparison with 2 other anti-HBc IgM screening tests. Final interpretation of samples was based on the "2 out of 3" rule (2 tests with the same interpretation take precedence over the other test). Sera which cannot be classed using this rule (a positive test, an equivocal test and a negative test) were not included in the performance calculation. Samples which gave an equivocal result using the VIDAS HBc IgM II kit were not included when calculating test performance.

a. Sensitivity

247 documented samples were tested:

- 97 positive samples for anti-HBc IgM, including 5 samples from patients with acute hepatitis.
- 150 sequential sera from 20 subjects with chronic hepatitis B most of whom had received anti-viral therapy.

	Sera classified according to 2/3 rule*	
	Positive	Negative
VIDAS HBc IgM II Positive	194	7**
VIDAS HBc IgM II Negative	0	15

* 24 samples gave results which could not be interpreted between the different techniques (positive with VIDAS equivocal and negative with other tests) and 7 samples gave an equivocal result with VIDAS HBc IgM II. All these results involved chronic HBV carriers with active viral replication, possibly reduced by anti-viral therapy.

** Follow-up samples for chronic HBV carriers with active viral replication, possibly reduced by anti-viral therapy. Four of these patients were positive for HBV DNA.

Relative sensitivity for VIDAS HBc IgM II: 100% (95% confidence interval: 97.98% - 100%)

Out of the 152 documented samples which were positive for HBV DNA or for HBe Ag, the positive screening percentage for anti-HBc IgM using the VIDAS test was 79%. This value was 65.8% and 50.7% using the other two kits in the comparison.

Out of the 5 acute hepatitis cases tested, 4 patients showed anti-HBc IgM titers greater than 200 PEIU/mL and one patient had a titer of 95 PEIU/mL.

In chronic hepatitis cases, the evolution of anti-HBc IgM titers is useful for patient monitoring (natural or therapeutic evolution):

Example: patient with chronic hepatitis B with HBV replication, HBe Ag negative:

Date	VIDAS HBc IgM II	Interpretation	Other viral DNA marker
03/1996	43	Positive	ADN VHB +
End of March 1996: start of anti-viral therapy			
05/1996	> 200	Positive	ADN VHB -
06/1996	115	Positive	Not tested
09/1996	13	Positive	Not tested
01/1997	6	Equivocal	ADN VHB -
07/1997	1	Negative	ADN VHB -

b. Specificity

351 serum samples sent to the laboratory for diagnosis of viral hepatitis were tested.

	Sera classified according to 2/3 rule*	
	Positive	Negative
VIDAS HBc IgM II Positive	20	14**
VIDAS HBc IgM II Negative	0	297

* 3 samples gave an uninterpretable result and 17 samples gave an equivocal result with VIDAS HBc IgM II.

** A study of other markers and case histories showed that 6 of the 14 "false-positive" samples corresponded to chronic hepatitis cases with active replication (detection of viral DNA by branched DNA or by PCR).

Relative specificity for VIDAS HBc IgM II on a clinical population: 95.50%
(95% confidence interval: 92.51% - 97.33%).

Specificity for Blood Donors:

601 unselected samples of blood donors were tested. Two samples found to be positive and one sample found to be equivocal were confirmed to be negative with the other 2 anti-HBc IgM screening tests.

The equivocal result was not included in the performance calculation.

Relative specificity of the VIDAS HBc IgM II kit for this population is 99.67% (95% confidence interval: 98.77% - 99.91%).

Sensitivity on Patients With Acute Hepatitis B:

The VIDAS HBc IgM II reagent detected anti-HBc IgM in all 40 samples tested.

Precision

Within-run reproducibility (intra-assay precision) and between-run reproducibility (total precision) were calculated according to the recommendations of the NCCLS Document EP5-T2, volume 12 number 4.

Precision was evaluated with a negative sample and two positive samples distributed on the measurement range. Each sample was tested in duplicate in two different runs per day, over a 14-day period on one site.

The combined results are given below:

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CROSS REACTIVITY AND RELEVANT INTERFERENTS**Potentially Interfering Samples**

2 equivocal samples were not taken into account in this study.

	VIDAS HBc IgM II positive
anti HIV Ab	0/19
anti EBV Ab	1*/10
anti HCV Ab	0/15
anti HAV IgM	0/15
anti CMV IgM	0/15
anti-nuclear Ab	0/15
Rheumatoid factors	1*/10
Pregnant women	0/20
anti HTLV Ab	2**/14
anti- <i>E.coli</i> Ab	0/10

*these samples were found to be positive with an HBs Ag screening test.

**these samples were found to be negative with an HBs Ag screening test.

RANGE OF EXPECTED VALUES

The reported incidence of hepatitis B cases in Europe is approximately 20/100,000, ranging from 1/100,000 in Scandinavian countries to 60/100,000 in Central Europe. In Europe, endemic cases increase North to South and West to East. In South-East Asia, China, Black Africa or South America, endemic prevalence can exceed 10%.

In cases of acute hepatitis, anti-HBc IgM titers are generally high (> 100 PEIU/mL), then progressively decrease over a 6-8 month period whether infection develops towards recovery or chronicity.

WASTE DISPOSAL

Dispose of used and unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. MILICH D.R., Immune response to the hepatitis B virus: infection, animal models, vaccination, VIRAL HEPATITIS, 1997, 3, 63-103.
2. HOLLINGER F.B., Hepatitis B virus, in Fields Virology, Third Edition, Lippincott-Raven Publishers, Philadelphia, 1996, 2739-2807.
3. SJOGREN M. A, DUSKEIKO G. M., KEW M. C. AND SONG E. - Hepatitis B virus infection and hepatocellular carcinoma: correlation between IgM antibody to hepatitis core antigen, hepatitis B e antigen, and hepatitis B DNA. Am. J. Trop. Med. Hyg., 1988, 39 (6), 582-585 (8 - 082).
4. BRUNETTO, M.R., CERENZIA M.T., OLIVERI F., PIAONTONI P., RANDONE A., CALVO P., MANZINI P., ROCCA G., GALLI C. and BONINO F., Monitoring the natural course and response to therapy of chronic hepatitis B with an automated semi-quantitative assay for IgM anti-HBc, Journal of Hepatology, 1993, 19, 431-436.
5. TREPO C., CAUSSE C. Hépatites virales: dépistage et suivi biologique - Feuilles de biologie, 1989, 30, 171: 79-89.

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Release date	Part Number	Change Type	Change Summary
2015/01	09433F	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS INSTRUCTIONS FOR USE
2016/05	09433G	Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS
2017/05	09433H	Technical	CONTENT OF THE KIT (30 TESTS) – RECONSTITUTION OF REAGENTS

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VIDAS® HIV DUO Quick (HIV6)

VIDAS HIV DUO Quick is an automated HIV infection screening test for use on the VIDAS family instruments, for the combined detection of anti-HIV-1 (groups M and O) and anti-HIV-2 total immunoglobulins and HIV-1 p24 antigen in human serum or plasma (lithium heparin or EDTA) using the ELFA technique (Enzyme Linked Fluorescent Assay).

SUMMARY AND EXPLANATION

Human immunodeficiency viruses (HIV) are RNA retroviruses transmitted via sexual contact, parenteral and perinatal pathways or the placenta. HIV-1 and HIV-2 were respectively isolated in 1983 and 1985 in patients infected by AIDS (Acquired Immunodeficiency Syndrome). Since then numerous genetic variants have been characterized. These mutations seemed to be without consequence for serological diagnosis until HIV-1 variants of group O (Outlier) were isolated, since they have only 50% homology at the *env* gene level with those of group M (Major). (1, 2, 3, 4, 5).

During 2007, 2.5 million individuals were newly infected by HIV. According to estimates, 33.2 million [30.6 – 36.1 million] people were living with HIV in 2007 (6).

Current diagnosis of HIV infection relies on the detection of anti-HIV serum antibodies using an ELISA method. However, there is a mean period of 3 weeks between contamination and the appearance of the first antibodies (7). During this period p24 antigen is present in most people infected by HIV-1 (8). Thus simultaneous detection of p24 antigenemia, anti-HIV-1 and anti-HIV-2 antibodies enables the time lapse between contamination and diagnosis of the infection to be decreased (9, 10). VIDAS HIV DUO Quick is an automated test based on the combined detection of HIV-1 p24 antigen and anti-HIV-1 and anti-HIV-2 total immunoglobulins, enabling a reduction of the seroconversion window.

Note: The use of VIDAS HIV DUO Quick does not rule out the obligation of using a second technique in countries where legislation requires the use of two different screening tests.

PRINCIPLE

The principle of the test combines 2 enzyme immunoassay reactions with a final fluorescent detection (ELFA).

The Solid Phase Receptacle (SPR®), serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and are pre-dispensed in the sealed reagent strips. All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times.

The upper part of the SPR is coated with monoclonal anti-p24 antibodies for the detection of p24 antigen. The lower part of the SPR enables the detection of anti-HIV-1 and anti-HIV-2 antibodies: it is coated with an HIV-1 gp160 protein, and HIV-1 group O and HIV-2 specific synthetic peptides.

During preliminary incubation, the sample, the biotinylated anti-p24 antibody (rabbit) in the strip and the biotinylated antigen (the same used in the solid phase) are cycled in and out of the SPR. Lysis of the virus occurs during incubation and the released p24 antigen binds to the monoclonal anti-p24 antibody coated on the SPR and is recognized by the biotinylated anti-p24 antibody. At the same time, the anti-HIV-1 and/or anti-HIV-2 antibodies bind to the gp160 and/or the peptides in the lower part of the SPR and are recognized by the biotinylated antigen. Two washing steps remove unbound components.

A second incubation with the biotinylated antigen in the strip (same used in the solid phase) is only performed in the lower part of the SPR. The biotinylated antigen binds to any anti-HIV antibody coated in the lower part of the SPR.

Any excess reagent is eliminated by washing.

A third incubation is performed with alkaline phosphatase-labeled streptavidin. During this step, the streptavidin binds to the biotinylated anti-p24 antibody, if it is present in the upper part of the SPR, and to the biotinylated antigen, if it is present in the lower part of the SPR.

Any excess reagent is eliminated by washing.

During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the presence of anti-HIV antibody and/or p24 antigen in the sample.

At the end of the assay, results are automatically calculated by the instrument in relation to the standard, and then printed out.

CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS:

60 HIV6 strips	STR	Ready-to-use.
60 HIV6 SPRs (2 x 30)	SPR	Ready-to-use. Interior of SPRs coated with gp160 (HIV-1) protein, HIV-1 group O and HIV-2 specific synthetic peptides, and monoclonal anti-p24 antibodies.
HIV6 positive antibody control (1 x 2.0 ml) (lyophilized)	C1	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C. Human serum* with anti-HIV-1 antibody + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives. MLE data indicate the confidence interval for the antibody result (Control C1 (+) Test Value Range).
HIV6 negative control (1 x 1.3 ml) (liquid)	C2	Ready-to-use. Human serum** + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives.
HIV6 positive antigen control (1 x 2.0 ml) (lyophilized)	C3	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C. Human serum** with inactivated HIV-1 viral lysate + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives. MLE data indicate the confidence interval for the antigen result (Control C3 (+) Test Value Range).-
HIV6 standard (1 x 2.0 ml) (lyophilized)	S1	Reconstitute with 2 ml of distilled water. Leave for 5 to 10 minutes and then mix. After reconstitution, stable for 2 months at 2-8°C. Human serum* with anti-HIV-1 antibody + TRIS (0.1 mol/l, pH 7.4) + protein and chemical stabilizers + preservatives. MLE data indicate the confidence interval in "Relative Fluorescence Value (RFV)" ("Standard (S1) RFV Range").
Specifications for the factory master data required to calibrate the test: • MLE data (Master Lot Entry) provided in the kit, or • MLE bar code printed on the box label.		
1 Clip seal		
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

* This product has been tested and shown to be negative for HBs surface antigen and antibodies to HCV. This product has been heat inactivated (30 minutes at 56°C). However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

** This product has been tested and shown to be negative for HBs surface antigen, and antibodies to HIV1, HIV2 and HCV. However, since no existing test method can totally guarantee their absence, this product must be treated as potentially infectious. Therefore, usual safety procedures should be observed when handling.

The SPR

The interior of the SPR is coated during production with a gp160 protein and HIV-1 group O and HIV-2 specific synthetic peptides in the lower part, and monoclonal anti-p24 antibodies in the upper part. The antigens enable the capture of anti-HIV-1 and/or anti-HIV-2 antibodies, and the monoclonal antibodies enable the capture of HIV-1 p24 antigen.

Each SPR is identified by the HIV-6 code. Only remove the required number of SPRs from the pouch and **carefully reseal the pouch after opening using the clip seal provided with the kit.**

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Description of the VIDAS HIV DUO Quick strip:

Wells	Reagents
1	Sample well.
2	HEPES buffer + biotin-labeled anti-p24 antibody (rabbit) + triton X100 + goat serum + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (300 µl).
3 - 5 - 6 - 8 - 9	Wash buffer: TRIS + protein and chemical stabilizers + 0.9 g/l sodium azide (600 µl).
4	HEPES buffer + gp160 protein and biotin-labeled synthetic peptides + skimmed milk + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (300 µl).
7	Tracer: alkaline phosphatase-labeled streptavidin + TRIS + skimmed milk + 2.5 g/l bovine albumin + triton X100 + protein and chemical stabilizers + 0.2 g/l gentamicin sulfate + 0.9 g/l sodium azide (400 µl).
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/l) + diethanolamine (DEA*) (0.62 mol/l or 6.6%) pH 9.2 + 1 g/l sodium azide (300 µl).

* Signal Word: **DANGER**Hazard statement

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- Pipette with disposable tip to dispense 2 ml and 200 µl.
- Powderless, disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User's Manual.
- VIDAS family instrument.

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- **This kit contains products of human origin. No known analysis method can totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (see Laboratory biosafety manual - WHO - Geneva - Latest edition).**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).
- Do not use the SPRs if the pouch is pierced.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.
- Do not mix reagents (or disposables) from different lots.

- Use **powderless** gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- The substrate in well 10 contains an irritant agent (6.6% diethanolamine). Refer to the hazard statements "H" and the precautionary statements "P" above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User's Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (see the User's Manual).

STORAGE CONDITIONS

- Store the VIDAS HIV DUO Quick kit at 2-8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at 2-8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPRs.
- **To maintain stability of the remaining SPRs, carefully reseal the pouch after use with the dessicant inside using the clip seal provided, and return the complete kit to 2-8°C.**

- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Refer to the kit composition table for special storage conditions.

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- hemolysis (after spiking samples with hemoglobin, up to 300 µmol/l (monomer)).
- triglyceride concentration (after spiking samples with up to 30 mg/ml of triglycerides)
- bilirubin concentration (after spiking samples with up to 200 mg/l of bilirubin).

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Specimen stability

Samples can be stored at 2-8 °C in stoppered tubes for up to 2 days. If longer storage is required, freeze the sera or plasma at -25 ± 6°C.

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A study performed on samples frozen for 2 months, showed that the quality of results is not affected.

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Procedure

1. **Only remove the required reagents from the refrigerator and allow them to come to room temperature for at least 30 minutes.**
2. Use one "HIV6" strip and one "HIV6" SPR for each sample, control or standard to be tested. **Make sure the storage pouch has been carefully resealed after the required SPRs have been removed.**
3. The test is identified by the "HIV6" code on the instrument. The standard must be identified by "S1", and tested in **duplicate**. If the positive controls are to be tested, they should be identified by "C1" and "C3". If the negative control needs to be tested, it should be identified by "C2".
4. If necessary, clarify samples by centrifugation.
5. Mix the standard, controls and samples using a vortex-type mixer (for serum or plasma separated from the pellet) in order to improve result reproducibility.
6. **For this test, the calibrator, control, and sample test portion is precisely 200 µl. The test portion is pipetted into the sample well, which is used as the reaction cuvette for this test.**
7. Insert the "HIV6" SPRs and "HIV6" strips into the instrument. Check to make sure the color labels with the assay code on the SPRs and the Reagent Strips match.
8. Initiate the assay as directed in the User's Manual. All the assay steps are performed automatically by the instrument.
9. Restopper the vials and return them to 2–8°C after pipetting.
10. The assay will be completed within approximately 80 minutes. After the assay is completed, remove the SPRs and strips from the instrument.
11. Dispose of the used SPRs and strips into an appropriate recipient.

RESULTS AND INTERPRETATION

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in each Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR is introduced into the substrate. The second reading is taken after the substrate in the SPR has been incubated. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The test value is calculated by the instrument for each sample as follows:

$$\text{Test value} = \text{patient RFV} / \text{standard RFV}$$

This test value and the interpreted result are also included on the result sheet.

The test value is interpreted as follows:

Test value	Interpreted result
< 0.25	Negative
≥ 0.25	Positive

If the message "Invalid" is generated for the sample, the result is invalid. In this case, the sample should be retested. Samples with a negative test value are considered to be negative, within the performance limitations of the reagent. In cases of suspected primary infection, values ≤ 0.25 must be interpreted with caution.

Samples with a positive test value should be retested in duplicate:

- Non repeatable positive samples (two negative reactions for three tests) are considered to be negative, within the performance limitations of the reagent.
- Repeatable positive samples (at least two positive reactions for three tests) should be confirmed using additional techniques.

The additional tests can be analyzed by Western Blotting and/or using a second screening test for detection of anti-HIV antibodies, p24 antigen and/or determination of viral load.

It is strongly recommended to test a second sample collected a few days later, particularly in the presence of clinical symptoms and/or risk factors.

Note: country-specific requirements for HIV diagnostics must be taken into account if necessary.

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QUALITY CONTROL

A positive antibody control, a positive antigen control, and a negative control are included in each VIDAS HIV DUO Quick kit. These controls must be performed immediately after opening a new kit to ensure that reagent performance has not been altered. Each calibration must also be checked using these controls. The instrument will only be able to check the control values if they are identified by C1, C2 and C3.

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LIMITATIONS OF THE METHOD

VIDAS HIV DUO Quick is an HIV infection screening test.

It should not be used as a specific test for the detection of HIV-1 p24 antigenemia. VIDAS HIV DUO Quick must not be used as a supplementary test for VIDAS HIV DUO Ultra.

This test has been validated for use with serum or plasma but not with other body fluids such as saliva, CSF or urine.

Since the use of recalcified plasma may generate false positive results, this type of sample must not be used.

Do not use a mixture of samples.

Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient history, and the results of any other tests performed.

PERFORMANCE

Studies performed using VIDAS HIV DUO Quick gave the following results:

1. Specificity on a blood donor population:

5027 blood donor samples from 3 blood banks were tested.

VIDAS HIV DUO Quick	Final interpretation	
	Positive	Negative
Positive	0	6
Negative	0	5021

Specificity of the VIDAS HIV DUO Quick reagent on this population: 99.88%

(95% confidence interval: 99.74% - 99.95%).

2. Specificity on hospitalized patients:

Two false positive results were found for the 200 samples tested.

Specificity of the VIDAS HIV DUO Quick reagent on this population: 99.00%.

(95% confidence interval: 96.35% - 99.73%)

3. Specificity on high-risk behavior patients:

No false positive results were found for the 100 samples tested.

Specificity of the VIDAS HIV DUO Quick reagent on this population: 100.00%.

(95% confidence interval: 96.38% - 100.00%)

4. Diagnostic sensitivity:

A study was performed using 440 samples presumed to be HIV-1 positive and 122 HIV-2 positive. All of the samples were found to be positive.

Sensitivity of the VIDAS HIV DUO Quick reagent on this population: 100.00% (95% confidence interval: 99.29% - 100.00%).

In order to check the sensitivity of VIDAS HIV DUO Quick on non-B subtypes, 40 group M samples (4A, 14B, 3C, 3D, 3F, 3G, 3H, 3CRF01-AE and 4 CRF02-AG) and 10 group O samples were tested. All were found to be positive.

31 fresh samples with a positive status (collection < 24 hours) were tested and found to be positive with VIDAS HIV DUO Quick. 31 fresh samples with a negative status (collection < 24 hours) were tested and found to be negative with VIDAS HIV DUO Quick.

In order to check the sensitivity of VIDAS HIV DUO Quick on p24 antigen, at least 50 culture supernatants including different HIV-1 and HIV-2 subtypes and at least 50 p24 antigen-positive samples were tested. All were found to be positive.

5. Sensitivity on seroconversion panels

During different studies, 39 seroconversion panels were tested, demonstrating the early detection capabilities of VIDAS HIV DUO Quick.

The dilutions of the 3 HIV-1 samples and the 2 HIV-2 samples of the EFS 96 panel were found to be positive.

The sensitivity of VIDAS HIV DUO Quick, determined on the Ag HIV SFTS panel (subtype B), was estimated at 16.5 pg/ml of HIV-1 antigen. The sensitivity of VIDAS HIV DUO Quick, determined using the international standard NIBSC 90/636 was estimated at 0.5 IU/ml of HIV-1 antigen.

At least 40 early HIV seroconversion samples were tested. The results are conform with the state of the art.

6. Precision

Intra-assay and inter-assay reproducibility, determined at two different sites, were calculated according to the recommendations of NCCLS document EP5-T2, volume 12-4.

Site 1

Sample	n	Mean index	Intra-assay reproducibility %	Inter-assay reproducibility %
Low positive antigen	40	2.48	2.17	3.69
Low positive antibody	40	3.96	3.53	4.14

Site 2

Sample	n	Mean index	Intra-assay reproducibility %	Inter-assay reproducibility %
Low positive antigen	40	2.23	2.07	9.88
Low positive antibody	40	3.63	3.46	7.11

Global analysis

Sample	n	Mean index	Inter-assay reproducibility %
Low positive antigen	80	2.36	10.23
Low positive antibody	80	3.79	8.20

In this analysis, the reproducibility takes into account within-run variability, and between-run and between-site dispersion.

7. Cross-reactivity

134 HIV negative samples from patients whose physiological status may interfere with the VIDAS HIV DUO Quick assay, were tested. The following results were obtained:

	Positive VIDAS HIV DUO Quick results
anti-HAV + antibody	0/10
anti-HCV + antibody	0/10
anti-HBV + antibody	0/10
anti-Gag + antibody	0/13
anti-EBV + antibody	0/10
anti-HSV + antibody	0/10
anti-CMV + antibody	0/10
anti-nuclear + antibody	0/10
anti- <i>Toxoplasma gondii</i> + antibody	0/6
anti- <i>Treponema pallidum</i> + antibody	0/5
Vaccinated against influenza	0/10
Vaccinated against hepatitis B	0/10
Rheumatoid factor +	0/10
Children under 15 years of age	0/10

Spiking of HIV positive and HIV negative samples with increasing concentrations of biotin (0 to 2 mg/l) did not affect the VIDAS HIV DUO Quick result.

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LITERATURE REFERENCES

1. F.Barré-Sinoussi, J.C. Chermann, F.Rey et al. Isolation of T-Lymphotropic Retrovirus from a Patient at Risk for acquired Immune Deficiency Syndrome (AIDS). 20 May 1983. Science.Vol.220. page 868-871.
2. F.Clavel, D.Guétard, F.Brun-Vézinet et al. Isolation of a New human Retrovirus from West African Patients with AIDS. 18 July 1986. SCIENCE Vol 233 page 343-346.
3. R.de Leys, B.Vanderborgh, M Vanden Haesevelde and al. Isolation and Partial Characterization of an Unusual Human Immunodeficiency Retrovirus from Two Persons of West-Central African Origin. Journal of Virology. Mar.1990.p1207-1216.
4. M.Guyader, M.Emerman, P.Sonigo et al. Genome organization and transactivation of the human immunodeficiency virus type 2. Nature. 16 april 1987. Vol 326 p.662-669.
5. M.Popovic, M.G.Sarngadharan, E. Read et al. Detection, Isolation, and Continuous Production of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and Pre-AIDS. Science. 4 May 1984. Vol.224. p 497-500.
6. World Health Organization, AIDS epidemic update. December 2007. p 1-60.

7. S.Kleinman, M.P.Busch, J.J.Korelitz et al. The Incidence/Window Period Model and its Use to Assess the Risk of Transfusion-Transmitted Human Immunodeficiency Virus and Hepatitis C Virus Infection. *Transfusion Medecine Reviews*. July 1997. Vol 11, N°3 p.155-172.
8. M.Von Sydow, H.Gaines, A.Sönnernborg et al. Antigen detection in primary HIV infection. *British Medical Journal* vol.296. 23 January 1988. p.238-240.
9. A.Galetto-Lacour, S.Kinloch-de Loës, b.Hirschel et al. Primo-infection VIH: un diagnostic souvent évoqué mais posé tardivement. *Schwerz.Med Wochenschr* 1995: 125 n°8. p.341-346.
10. T.Meier, E.Knoll, M.Henkes et al. Evidence for a diagnostic window in fourth generation assays for HIV. *Journal of Clinical Virology* vol.23, 2001, p.113-116.

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Note:

Minor typographical, grammar, and formatting changes are not included in the revision history.

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Release date	Part Number	Change Type	Change Summary
2015/01	12232E	Administrative	INDEX OF SYMBOLS REVISION HISTORY
		Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS WARNINGS AND PRECAUTIONS
2015/06	12232F	Technical	CONTENT OF THE KIT (60 TESTS) – RECONSTITUTION OF REAGENTS INSTRUCTIONS FOR USE

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LITERATŪRA

1. F.Barré-Sinoussi, J.C. Chermann, F.Rey et al. Isolation of T-Lymphotropic Retrovirus from a Patient at Risk for acquired Immune Deficiency Syndrome (AIDS). 20 May 1983. Science.Vol.220. page 868-871.
2. F.Clavel, D.Guébard, F.Brun-Vézinet et al. Isolation of a New human Retrovirus from West African Patients with AIDS. 18 July 1986. SCIENCE Vol 233 page 343-346.
3. R.de Leys, B.Vanderborgh, M.Vanden Haesevelde and al. Isolation and Partial Characterization of an Unusual Human Immunodeficiency Retrovirus from Two Persons of West-Central African Origin. Journal of Virology. Mar.1990.p1207-1216.
4. M.Guyader, M.Emerman, P.Sonigo et al. Genome organization and transactivation of the human immunodeficiency virus type 2. Nature. 16 april 1987. Vol 326 p.662-669.
5. M.Popovic, M.G.Sarngadharan, E. Read et al. Detection, Isolation, and Continuous Production of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and Pre-AIDS. Science. 4 May 1984. Vol.224. p 497-500.
6. World Health Organization, AIDS epidemic update. December 2007. p 1-60.
7. S.Kleinman, M.P.Busch, J.J.Korelitz et al. The Incidence/Window Period Model and its Use to Assess the Risk of Transfusion-Transmitted Human Immunodeficiency Virus and Hepatitis C Virus Infection. Transfusion Medicine Reviews. July 1997. Vol 11, N°3 p.155-172.

8. M.Von Sydow, H.Gaines, A.Sönnnerborg et al. Antigen detection in primary HIV infection. British Medical Journal vol.296. 23 January 1988. p.238-240.
9. A.Galetto-Lacour, S.Kinloch-de Loës, b.Hirschel et al. Primo-infection VIH : un diagnostic souvent évoqué mais posé tardivement. Schweiz.Med Wochenschr 1995 : 125 n°8. p.341-346.
10. T.Meier, E.Knoll, M.Henkes et al. Evidence for a diagnostic window in fourth generation assays for HIV. Journal of Clinical Virology vol.23, 2001, p.113-116.

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
Quality Control VIDAS[®] (QCV)

Automated test for use on the VIDAS[®] system to detect abnormal operation of the VIDAS and mini VIDAS instrument pipette mechanisms and optical systems.

SUMMARY AND EXPLANATION

Quality Control VIDAS (QCV) is used to detect abnormal operation of pipette mechanisms which may affect the results of biological tests. It is also intended for checking that the optical system is capable of measuring high fluorescence levels.

The QCV test must be run at least once a month in each position of the VIDAS and mini VIDAS instruments, or any time that one of the above-mentioned instrument problems is suspected.

PRINCIPLE

The test corresponds to successive aspirations/dilutions of fluorescent substrate (4-Methyl-umbelliferone) solutions which are standardized and have varying levels of concentration. The aspirations are performed at different speeds to check pump aspiration capability.

The results obtained at the end of the test correspond to:

- a fluorescence ratio "Test Value 1" (TV1 for mini VIDAS) for checking that the pipette mechanism is operating correctly,
- a fluorescence reading (Reading 3 for VIDAS and R3 for mini VIDAS) for checking that the optical system is intact.

The values obtained must be within the acceptable range defined later in this package insert.

CONTENT OF THE KIT (60 TESTS):

60 QCV Strips	STR	Ready-to-use.
60 QCV SPRs (2 x 30)	SPR	Ready-to-use.
1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib		

DESCRIPTION

The SPR[®]:

The SPR is used as a pipetting device and is identified by the code QCV. **Only remove the required number of SPRs from the pouch.**

The strip:

The strip consists of 10 wells covered with a labeled, foil seal. The label comprises a bar code which mainly indicates the assay code, kit lot number and expiration date. The last well is a cuvette in which the fluorometric reading is performed. The wells in the center of the strip contain the various reagents required for the assay.

Description of the QCV strip:

Wells	Reagents
1	Empty well: no sample is needed for this test.
2	400 µL of 4-Methyl-umbelliferone solution (100 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
3	Empty well.
4	400 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
5 - 6	600 µL of 4-Methyl-umbelliferone solution (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
7	Empty well.
8 - 9	400 µL of 4-Methyl-umbelliferone solution (26 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.
10	Optical cuvette with substrate: 300 µL of 4-Methyl-umbelliferone (0.6 µmol/L) in CHES buffer (20 mmol/L, pH 9.6) + 0.9 g/L sodium azide.

WARNINGS AND PRECAUTIONS

1. **Caution: the strips are photosensitive and must be kept in their box to be protected from light.**
2. Do not use reagents after the expiration date indicated on the label.
3. Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.

4. It is recommended to wear gloves as if handling biological tests.

STORAGE

- Store the VIDAS QCV kit at 2 - 8°C. Do not freeze reagents. **Store all unused reagents in their box at 2-8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the label. Do not use the reagents after the expiration date.

INSTRUCTIONS FOR USE

Your instrument should be checked **at least once a month** using the VIDAS QCV test, and each time abnormal operation of the pipette mechanism or optical system is suspected.

It is essential to perform the VIDAS QCV test on all the positions to ensure that the whole instrument has been checked.

It is essential to archive patient work lists and their results on a monthly basis in order to track results obtained between two QCV tests.

1. Remove necessary components from the kit and return all unused components to storage at 2-8°C.

2. Create the work list: type QCV and indicate the number of tests to be run.
3. Load the reagent strips and SPRs into the instrument (**Note:** no desiccant is needed in the SPR pouch). Be sure to close all SPR® doors and tray covers.
4. Run the test (see User's Manual). All steps will be performed automatically by the instrument. The test will be completed in approximately 20 minutes.
5. **After the test is completed, print the results and analyze them as instructed later in this package insert. Once interpretation of results has been validated, dispose of the used SPRs and strips into an appropriate recipient.**

RESULTS AND INTERPRETATION

The first reading is performed at the beginning of the assay to check that the reagent is intact:

reading 1 for VIDAS
R1 for mini VIDAS

The second reading is taken after successive aspirations which are performed at different speeds; they are used to check the ability of the pump to pipette properly:

reading 2 for VIDAS
R2 for mini VIDAS

The third reading is taken after a small quantity of highly concentrated substrate is pipetted. This operation is performed to check that the optical system is capable of measuring high fluorescence levels:

reading 3 for VIDAS
R3 for mini VIDAS


When the test is completed, the computer automatically calculates Test Value 1 and 2 (**TV1** and **TV2** for mini VIDAS) using the 3 fluorescence readings. The calculations appear on a report.

VIDAS: Test Value (TV1) = $\frac{\text{reading 2}}{\text{reading 1}}$ Test Value (TV2) = $\frac{\text{reading 3}}{\text{reading 1}}$

mini VIDAS: TV1 = $\frac{R2}{R1}$ TV2 = $\frac{R3}{R1}$

The values for reading 3 (R3 for mini VIDAS) and Test Value 1 (TV1 for mini VIDAS) must be within the acceptable ranges indicated below (see diagram):


- **Acceptable range for Test Value 1 (TV1 for mini VIDAS):** the Test Value 1 for each position must be **≥ the value indicated on the kit label**. If the result of a particular position is outside the range, two new VIDAS QCV tests must be run successively in all the positions in the section concerned (for VIDAS and mini VIDAS). If at least one other non-compliant result is produced in the same section, independent of the position, take the section offline and contact bioMérieux.
- **Acceptable range for reading 3 (R3 for mini VIDAS):** reading 3 for each position must be **≥ 4100 RFU**. If the result of a particular position is outside the range, take the instrument offline and contact bioMérieux.

LIMITATIONS OF THE METHOD

- VIDAS QCV may produce an invalid result when the reading 1 or R1 value is < 200 RFU. In this case, perform the test again using a new QCV strip from a new lot (if possible).
- VIDAS QCV should only be used to diagnose problems with the instrument pipette mechanism or optical system. It is not a calibrator.
- VIDAS QCV does not diagnose temperature problems or leakages.
- VIDAS QCV should only be used to test the instrument.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Protect from light

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IVD

Intended Use

VIDAS® PTH (1-84) (PTH) is a third generation automated quantitative assay for use on the VIDAS® family of instruments, for the quantitative measurement of the biologically active parathyroid hormone PTH (1-84) in human serum or plasma using the ELFA (Enzyme Linked Fluorescent Assay) technique.

Used in conjunction with other laboratory findings and clinical assessments, this assay is intended for use as follows:

- As an aid in the diagnosis of hyper or hypoparathyroidism.
- As an aid for the monitoring of calcium homeostasis in patients with chronic kidney disease.

Summary and Explanation

Parathyroid hormone (PTH) is the main hormone involved in calcium and phosphorus homeostasis. PTH is formed in the parathyroid glands (essentially four, located on the back of the thyroid gland) and secreted as an 84 amino-acid peptide [1-84 according to the universal numbering convention beginning with the N-terminus] when the extracellular calcium concentration decreases.

The main function of the PTH is to maintain the ionized calcium level within a narrow range by a powerful hypercalcemic effect, by stimulating calcium release from bone and distal tubular reabsorption of calcium by the kidney, and by increasing renal production of 1,25-dihydroxyvitamin D-3 (calcitriol) which stimulates intestinal calcium absorption.

PTH increases urine phosphate levels and decreases serum phosphate levels. PTH stimulates bone remodeling by mobilizing osteoblasts, thus justifying its use in osteoporosis treatment to reduce the risk of vertebral fractures.

The secretion of PTH is promoted by low calcium concentrations and high phosphate concentrations, and is inhibited by high calcium concentrations and low phosphate concentrations. It is also inhibited by calcitriol, PTH (7-84) fragment, and magnesium deficiency.

In the bloodstream, the half-life of PTH (1-84) does not exceed 4 minutes. PTH is rapidly degraded by the liver into non-(1-84) C-terminal fragments, essentially (7-84) and (53-84) with long half-lives, and which are eliminated by the kidneys. These fragments accumulate in the blood when the kidneys fail to function.

Parathyroid gland disorders lead to either hyperparathyroidism (primary, secondary or tertiary) or hypoparathyroidism. Apart from bone fragility⁽¹⁾ or nephrolithiasis, the associated symptoms may not be very specific, and asymptomatic parathyroid disorder may be unexpectedly detected during a routine blood test⁽²⁾.

This biological assay contributes in determining a patient's status regarding the function of the parathyroid glands as normal or hyper- or hypoparathyroidism, when interpreted with other biological parameters (such as calcium, phosphate, vitamin D) and clinical signs.

PTH measurement is indicated when performing a differential diagnosis for hypercalcemia or hypocalcemia. It enables the diagnosis of primary or secondary hyperparathyroidism⁽³⁾, and the diagnosis of hypoparathyroidism⁽⁴⁾ mainly after parathyroid surgery. It is also prescribed for patients with chronic kidney disease, in particular in the follow-up of dialysis patients⁽⁵⁾.

In case of kidney failure, medical monitoring is guided by the KDOQI (Kidney Disease Outcomes Quality Initiative) and KDIGO (Kidney Disease Improving Global Outcomes) recommendations.

In non-kidney disorders, PTH levels in blood are one of the results of routine biological tests investigated by a healthcare specialist (for example: Endocrinology, Rheumatology, Geriatrics, Metabolic Diseases).

Principle

The assay principle combines a 2-step sandwich enzyme immunoassay method with a final fluorescence detection (ELFA). VIDAS® PTH (1-84) is a third generation assay traced to the WHO International Standard for Parathyroid Hormone 1-84, recombinant, coded 95/646.

The single-use Solid Phase Receptacle (SPR) serves as the solid phase as well as the pipetting device. Reagents for the assay are ready-to-use and pre-dispensed in the sealed single-use reagent strips.

All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR device several times.

The sample is collected and then diluted in a suitable diluent. The mixture is cycled in and out of the SPR device several times. This operation enables PTH to bind with the capture antibodies fixed to the interior wall of the SPR device.

Unbound components are eliminated during washing steps.

The conjugate mixture is cycled in and out of the SPR device several times. This operation enables the conjugate to bind with the immune complexes that are fixed to the interior wall of the SPR device and to form a sandwich. Unbound components are eliminated during washing steps.

During the final detection step, the substrate (4-Methylumbelliferyl phosphate) is cycled in and out of the SPR device. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm. The intensity of the fluorescence is proportional to the concentration of PTH in the sample.

At the end of the assay, the results are automatically calculated and reported in pg/mL by the instrument according to the calibration curve stored in memory. The results can then be printed out.

Content of the Kit (30 tests)

30 PTH Strips ^(a)	STR	Ready-to-use.
30 PTH Solid Phase Receptacles 1 x 30	SPR	Ready-to-use. Interior of SPR device coated with mouse monoclonal anti-PTH immunoglobulins (specific to the N-terminal (1-6) fragment).
PTH Calibrator ^(b) 1 x 4.4 mL (liquid)	S1	Ready-to-use. Buffer containing recombinant PTH (1-84) + Stabilizer of animal origin + preservative. MLE (Master Lot Entry) data indicate the acceptable range in "Relative Fluorescence Value" ("Calibrator (S1) RFV Range").
PTH Control ^(b) 1 x 2.3 mL (liquid)	C1	Ready-to-use. Buffer containing recombinant PTH (1-84) + Stabilizer of animal origin + preservative. MLE data indicate the acceptable range in pg/mL ("Control C1 Dose Value Range").
Specifications for the factory master data required to calibrate the assay: MLE barcode printed on the box label.		
1 package insert provided in the kit or downloadable from www.biomerieux.com .		

^(a) DANGER
P305 + P351 + P338

WARNING

EUH208 / H317 / H318 / P261 / P280 / P302 + P352 /

^(b) WARNING

EUH208 / H317 / P261 / P280 / P302 + P352

Hazard statements

EUH208: Contains 2-methyl-2H-isothiazolin-3-one. May produce an allergic reaction.

H317: May cause an allergic skin reaction.

H318: Causes serious eye damage.

Precautionary statements

P261: Avoid breathing dust/fume/gas/mist/vapours/spray.

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P302 + P352: IF ON SKIN: Wash with plenty of soap and water.

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For further information, consult the Safety Data Sheet.

The SPR device

The interior of the SPR device is coated during production with mouse monoclonal anti-PTH N-ter (1-6) fragment immunoglobulins. Each SPR device is identified by the "PTH" code.

Only remove the required number of SPR devices from the pouch and carefully reseal the pouch after opening.

Reagent Strip

The strip consists of 10 wells covered with a labeled foil seal. The label comprises a barcode which mainly indicates the assay code, kit lot number, and expiration date. The foil of the first well is perforated to facilitate the introduction of the sample. The last well of each strip is a cuvette in which the fluorometric reading is performed. The wells in the center section of the strip contain the various reagents required for the assay.

Description of the VIDAS® PTH (1-84) (PTH) strip

Well	Reagents
1	Sample well: dispense 300 µL of calibrator, control, or sample.
2	Sample diluent: Buffer + Protein stabilizer of animal origin + preservative (300 µL).
3	Empty well.
4 - 5	Wash buffers: Buffer + preservative (600 µL).
6	Conjugate: Buffer containing anti-(7-84) PTH antibodies + Stabilizer of animal origin + preservative (400 µL).
7 - 8	Wash buffers: Buffer + preservative (600 µL).
9	Empty well.
10	Reading cuvette with substrate: 4-Methyl-umbelliferyl phosphate (0.6 mmol/L) + diethanolamine (DEA) (0.62 mol/L or 6.6%, pH 9.2) + 1 g/L sodium azide (300 µL).

Materials and Disposables Required but Not Provided

- Single-use pipette and/or micropipettes to dispense the appropriate volumes.
- Powderless disposable gloves.
- For other specific materials and disposables, please refer to the Instrument User Manual.
- Instruments of the VIDAS® family.

Warnings and Precautions

- **For *in vitro* diagnostic use only.**
- **For professional use only, by qualified laboratory personnel, in clinical laboratories.**
- This kit does not contain products of human origin.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest; do not inhale).
- Handle all patient specimens as potentially infectious material and follow the usual safety precautions.
- Although the PTH (1-84) recombinant proteins present in the S1 Calibrator and C1 Control are inactive, handle the tubes as potentially infectious products.
- Do not use the SPR devices if the pouch is pierced or if the dot sealing a SPR device has come unstuck.
- Do not use visibly deteriorated STRs (damaged foil or plastic).
- Do not use reagents after the expiration date indicated on the box label.

- Do not mix reagents (or disposables) from different lots.
- VIDAS® PTH (1-84) assay reagents are only for use with the instruments of the VIDAS® family.
- Use powderless gloves, as powder has been reported to cause false results for certain enzyme immunoassay tests.
- Kit reagents contain sodium azide which can react with lead or copper plumbing to form explosive metal azides. If any liquid containing sodium azide is disposed of in the plumbing system, drains should be flushed with water to avoid build-up.
- Refer to the hazard statements “H” and precautionary statements “P” indicated above.
- Spills should be wiped up thoroughly after treatment with liquid detergent or a solution of household bleach containing at least 0.5% sodium hypochlorite. See the User Manual for cleaning spills on or in the instrument. Do not autoclave solutions containing bleach.
- The instrument should be regularly cleaned and decontaminated (refer to the User Manual for user and preventive maintenance operations).

Storage Conditions

- Store the VIDAS® PTH (1-84) kit at +2°C/+8°C.
- **Do not freeze reagents.**
- **Store all unused reagents at +2°C/+8°C.**
- After opening the kit, check that the SPR pouch is correctly sealed and undamaged. If not, do not use the SPR devices.
- **After use, carefully reseal the pouch with the desiccant inside to maintain stability of the SPR devices, and return the complete kit to +2°C/+8°C.**
- If stored according to the recommended conditions, all components are stable until the expiration date indicated on the box label.

Samples

Sample type and collection

Human serum or plasma.

Types of tubes validated:

- Plastic tube with clot activator
- Plastic tube with clot activator and separation gel
- Plastic tube with lithium heparin
- Plastic EDTA tube

It is recommended to validate collection tubes before use as some contain substances which interfere with test results.

Note: Blood sampling tube results may vary from one manufacturer to another depending on the materials and additives used.

It is the responsibility of each laboratory to validate the type of sample tube used and to follow the manufacturer's recommendations for use.

Sample Preparation

The current WHO/DIL/LAB/99.1 document provides recommendations for sample preparation. ⁽⁶⁾

For use of sample tubes, refer to the tube manufacturer's recommendations for use.

The pre-analytical step, including the preparation of blood samples, is an essential first step when performing medical analyses. In accordance with Good Laboratory Practice, this step is performed under the responsibility of the laboratory manager.

Insufficient clot time can result in the formation of fibrin with micro-clots that are invisible to the naked eye. The presence of fibrin, red blood cells, or suspended particles can lead to erroneous results.

Samples containing suspended fibrin particles or erythrocyte stroma should be centrifuged before testing.

For serum samples, ensure that complete clot formation has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting times.

Sample stability

Samples (serum and plasma) can be stored in closed primary tubes at +18°C/+25°C for up to 8 hours and at +2°C/+8°C for up to 48 hours. If longer storage is required, freeze the serum or plasma at -19°C/-31°C.

A study performed on samples frozen for 6 months at -19°C/-31°C showed that the quality of results is not affected. Three freeze/thaw cycles were validated.

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- hemolysis (after spiking samples with hemoglobin, up to 5 g/L (monomer)),
- lipemia (after spiking samples with lipids, up to 30 g/L equivalent in triglycerides),
- bilirubinemia (after spiking samples with bilirubin, up to 0.3 g/L).

It is recommended not to use hemolyzed, lipemic, icteric samples, and, if possible, to collect a new sample.

Refer to the section **PERFORMANCE – Study of drugs and other potentially interfering substances** for the compounds tested.

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Reading VIDAS® PTC (Protocol Test Change) data and MLE data**When using the assay for the first time:**

With the external instrument barcode reader, **scan the barcodes (PTC and MLE) in the following order:**

1. According to the instrument used, scan the PTC barcode(s) downloadable from www.biomerieux.com. This reading allows VIDAS® PTC protocol data to be transferred to the instrument software for its update.
2. Scan the MLE data on the box label.

When opening a new lot of reagents:

With the external instrument barcode reader, scan the MLE data on the box label before performing the test.

Note: The master lot data need only be entered once for each lot.

It is possible to enter MLE data **manually or automatically** depending on the instrument (refer to the User Manual).

Calibration

The calibrator is traceable to the WHO International Standard 95/646 for PTH (1-84).⁽⁷⁾

Calibration, using the Calibrator provided in the kit, must be performed each time a new lot of reagents is opened, after the MLE data have been entered, and then every **84 days**. This operation provides instrument-specific calibration curves and compensates for possible minor variations in assay signal throughout the shelf life of the kit.

The Calibrator identified by S1 must be tested in **duplicate** (refer to the User Manual). The Calibrator value must be within the set RFV (Relative Fluorescence Value) range. If this is not the case, recalibrate using S1.

Caution: The assay Calibrator (S1) and Control (C1) are ready-to-use reagents made liquid by the addition of a tension-active stabilizer. When pipetting the 300 µL calibration and calibration control portion, make sure that the entire 300 µL volume is dispensed out of the pipet tip.

Kit controls

One Control is included in each VIDAS® PTH (1-84) (PTH) kit. This Control must be performed immediately after opening a new kit to ensure that reagent performance has not been altered.

Each calibration must also be checked using this Control. The instrument will only be able to check the control value if it is identified by C1.

Results cannot be validated if the control value deviates from the expected values.

Note: The aim of the VIDAS® PTH (1-84) Quality Control (C1) is to validate calibrations. The use of the VIDAS® PTH (1-84) Quality Control (C1) as an Internal Quality Control is under the customer's responsibility.

Procedure

1. **Remove the kit from storage at +2°C/+8°C and take out the required reagents. Carefully reseal the SPR pouch and return the kit to +2°C/+8°C. The reagents can be used immediately.**

2. Use one "PTH" strip and one "PTH" SPR device for each sample, control, or calibrator to be tested. Make sure the SPR pouch has been carefully resealed after the required Solid Phase Receptacles have been removed.
3. The test is identified by the "PTH" code on the instrument. The calibrator must be identified by "S1" and tested in duplicate. The control identified by "C1" must be tested singly.
4. Mix the calibrator, control, and samples using a vortex-type mixer (for serum or plasma separated from the pellet).
5. For optimal results, refer to all the paragraphs in the **SAMPLES** section.
6. Before pipetting, ensure that the samples, calibrator, control, and diluent are free of bubbles.
7. **For this test, the calibrator, control, and sample test portion is 300 µL.**
8. Insert the "PTH" Solid Phase Receptacles and "PTH" strips into the instrument. Check to make sure that the color labels with the assay code on the Solid Phase Receptacles and the Reagent Strips match.
9. Initiate the assay as directed in the User Manual. All the assay steps are performed automatically by the instrument.
10. Close the vials and return them to the required temperature after pipetting.
11. The assay will be completed within approximately **24 minutes**. After the assay is completed, remove the Solid Phase Receptacles (SPR) and strips from the instrument.
12. Dispose of the used Solid Phase Receptacles (SPR) and strips into an appropriate recipient.

Quality Control

Additional quality controls can be performed in accordance with local regulations or requirements related to accreditation, as well as requirements defined in the laboratory's quality control procedure.

Note: It is the responsibility of the user to perform Quality Control in accordance with any applicable local regulations.

Results and Interpretation

Once the assay is completed, results are analyzed automatically by the computer. Fluorescence is measured twice in the Reagent Strip's reading cuvette for each sample tested. The first reading is a background reading of the substrate cuvette before the SPR device is introduced into the substrate.

The second reading is taken after incubating the substrate with the enzyme bound to the interior of the SPR device. The RFV (Relative Fluorescence Value) is calculated by subtracting the background reading from the final result. This calculation appears on the result sheet.

The results are automatically calculated by the instrument using calibration curves which are stored by the instrument (predefined mathematical model) and are expressed as pg/mL.

Interpretation of test results should be made taking into consideration the patient's clinical history, and the results of any other tests performed.

The instrument displays the VIDAS® PTH (1-84) (PTH) assay results from 4.0 to 1500.0 pg/mL.

Sample dilution

Samples with VIDAS® PTH (1-84) concentrations greater than 1500.0 pg/mL must be retested after 1/2 dilution (dilution factor = 2) in saline solution (1 volume of sample + 1 volume of saline solution) in order to obtain a result within the measurement range.

The dilution must be performed manually on the VIDAS®, MINI VIDAS®, and VIDAS® 3 instruments.

If the dilution factor (= 2) was entered when the Work List was created, the result is calculated automatically. If the dilution factor was not entered, multiply the result by the dilution factor to obtain the sample concentration.

Limitations of the Method

- Interference may be encountered with certain sera containing antibodies directed against reagent components. For this reason, assay results should be interpreted taking into consideration the patient's clinical history and the results of any other tests performed.
- Refer to the section **PERFORMANCE – Study of drugs and other potentially interfering substances** for the compounds tested.
- Any results that do not correspond to the patient's clinical history may be due to inadequate instrument maintenance (see the Instrument User Manual).

Reference Values

These results are given as a guide. It is recommended that each laboratory establish its own reference values from a rigorously selected population.

Reference values were determined in a clinical study using 491 samples from apparently healthy individuals. The population was divided into 3 groups, depending on the vitamin D status. For each group, a 95% reference value range (or reference interval, according to CLSI EP28-A3c Guideline) was determined, as well as for the total of the apparently healthy population as described below:

	N	Median (in pg/mL)	PTH (1-84) reference interval (in pg/mL)	90% CI* of the Lower reference limit (in pg/mL)	90% CI of the Upper reference limit (in pg/mL)
Healthy population	491	20.8	9.2 - 44.6	[8.7 ; 9.7]	[42.5 ; 46.7]

Healthy population by Vitamin D level	N	Median (in pg/mL)	PTH (1-84) reference interval (in pg/mL)	90% CI of the Lower reference limit (in pg/mL)	90% CI of the Upper reference limit (in pg/mL)
Deficient ≤ 20 ng/mL	176	21.5	8.9 - 45.3	[8.0 ; 9.9]	[42.1 ; 48.8]
Insufficient > 20 and < 30 ng/mL	213	20.7	8.8 - 47.8	[7.6 ; 9.9]	[41.6 ; 57.3]
Normal ≥ 30 ng/mL	102	21.3	9.6 - 47.3	[8.7 ; 10.6]	[41.6 ; 54.1]

* CI: Confidence Interval.

According to the KDIGO guidelines, the PTH concentration of the patients treated by dialysis should be maintained within two and nine times the upper normal limit of the assay.

Performance

Studies performed using the VIDAS® PTH (1-84) (PTH) assay gave the following results.

Measurement range

The measurement range (or analytical measuring interval, according to CLSI EP36-Ed1 Guideline) is the range of values corresponding to the acceptable performance limits (precision and linearity).

The measurement range of the VIDAS® PTH (1-84) (PTH) assay is 4.0 to 1500.0 pg/mL.

Linearity

Linearity was evaluated according to CLSI EP06-A recommendations. The VIDAS® PTH (1-84) assay is linear between 4.0 and 1500.0 pg/mL.

Detection limits

Limit of Blank (LoB)	0.8 pg/mL
Limit of Detection (LoD)	2.2 pg/mL
Limit of Quantitation (LoQ)	4.0 pg/mL

LoB, LoD and LoQ values were determined according to CLSI EP17-A2 recommendations.

The Limit of Blank (LoB) is the concentration below which 95% of analyte-free samples are found.

The Limit of Detection (LoD) is the lowest concentration of analyte in a sample that can be distinguished from the analyte-free sample with a probability of 95% (observed result greater than the LoB with a 95% probability).

The Limit of Quantitation (LoQ) is the lowest concentration of analyte that can be detected and measured with an acceptable level of precision. All results under 4.0 pg/mL are printed ≤ 4.0 pg/mL. For the VIDAS® PTH (1-84) assay, the acceptable level of precision corresponds to within-lot precision fixed at 20% CV.

Hook effect

The assay uses a 2-step immunoassay sandwich method. Hook effect can therefore be excluded by design.

Precision

A precision study was performed according to CLSI EP05-A3 recommendations. A panel of human samples representing 6 concentration levels in the measurement interval was analyzed on the VIDAS® instrument to include the following main sources of variability: repeatability, run, day, calibration, and lot.

Repeatability (within-run precision), intermediate precision (total within-lot) and within-laboratory precision (between-lot within-instrument) were estimated for each sample. The values obtained are reported in the following table, as a guide.

Sample	Concentration level (pg/mL)	Repeatability		Intermediate precision		Within-laboratory	
		Standard deviation (pg/mL)	CV (%)	Standard deviation (pg/mL)	CV (%)	Standard deviation (pg/mL)	CV (%)
Sample 1	14.7	0.9	6.0	1.1	7.3	1.2	7.9
Sample 2	46.3	2.8	6.1	3.0	6.5	5.7	12.4
Sample 3	106.7	4.0	3.8	5.1	4.8	5.6	5.2
Sample 4	299.0	10.3	3.4	13.1	4.4	16.2	5.4
Sample 5	590.0	15.3	2.6	25.4	4.3	30.7	5.2
Sample 6	1130.1	42.5	3.8	51.1	4.5	62.4	5.5

These results are an indicator showing the ability to monitor chronic kidney disease patients in the routine mode.

Analytical specificity

The analytical specificity of the VIDAS® PTH (1-84) (PTH) assay was established by testing cross-reactive compounds according to CLSI document EP7-A2 recommendations. Cross-reactivity was evaluated by spiking serum samples with cross-reactive compounds (about 50 pg/mL and about 400 pg/mL of PTH).

The results of this study are reported in the following table.

Tested substance	Tested concentration (pg/mL)	Cross-reactivity %
PTH (7-84)	100 000	0.11%
PTH (1-34)	100 000	0.04%
PTH (13-34)	100 000	0.0004%
PTH (39-68)	100 000	0.0004%
PTH (44-68)	100 000	0.009%
PTH (39-84)	100 000	0.002%
PTH (53-84)	100 000	0.006%
Calcitonin	100 000	0.005%
Osteocalcin	100 000	0.014%
C-telopeptid	100 000	0.023%

Study of drugs and other potentially interfering substances

Potential interference by commonly used drugs and other substances was studied according to CLSI EP7-A2 recommendations. No significant interference was detected up to the maximum concentrations tested.

Tested drug	Maximum concentration	Tested drug	Maximum concentration
Acetaminophen	35 mg/dL	Etidronate	105 mg/dL
Acetylsalicylic acid	65 mg/dL	Ibuprofen	50 mg/dL (242.5 µmol/dL)
Alendronate	8 mg/dL	Lanthanum chloride	40 mg/dL
Alfacalcidol	2.5 µg/mL	Magnesium chloride	40 mg/dL
Aluminum sulfate	40 mg/dL	Pamidronate	18 mg/dL

Tested drug	Maximum concentration	Tested drug	Maximum concentration
Biotin	2 µg/mL	Risedronate	6 mg/dL
Calcitriol	1 ng/mL	Salicylic acid	60 mg/dL
Calcium acetate	40 mg/dL	Vitamin D2	1800 ng/mL
Calcium citrate	40 mg/dL	Vitamin D3	240 ng/mL

Tested substance	Maximum concentration
Bilirubin (unconjugated)	0.3 g/L (510 µmol/L)
Cholesterol	4.21 g/L
Hemoglobin	5 g/L (310 µmol/L)
Human Albumin	60 g/L
Human Anti Mouse Antibodies (HAMA)	2 µg/mL
Lipids (intralipids)	30 g/L
Rheumatoid factor	800 IU/mL

Method comparison

The comparison studies were performed according to CLSI EP09-A3 recommendations.

All the intended populations covering the measuring interval from 4.0 to 1500.0 pg/mL have been tested with VIDAS® PTH (1-84) (PTH) assay. Among the 215 patients enrolled, the following populations were represented:

- Adenoma
- Bone disorder
- Dialysis
- Non parathyroidian hypercalcemia
- Parathyroidectomy
- Renal failure stage 3
- Renal failure stage 4
- Renal failure stage 5
- Renal transplanted
- Apparently healthy population with $20 < \text{Vit D} < 30 \text{ ng/mL}$
- Apparently healthy population with $\text{Vit D} \geq 30 \text{ ng/mL}$

The comparison of the VIDAS® PTH (1-84) assay (Y) with another commercially available PTH (1-84) immunoassay (X) gave the following results:


Number of samples tested: 215

Samples tested between 4.1 pg/mL and 1179.4 pg/mL for the VIDAS® PTH (1-84) (PTH) assay.

Samples tested between 6.68 pg/mL and 1200.00 pg/mL for the other commercially available immunoassay.

Equation for Passing-Bablok regression: $Y = 1.0460 X - 2.2489$

Correlation Coefficient: 0.9889

Monitoring performance

A total number of 21 patients were enrolled with 3 different time points. In conclusion a total of 63 samples were collected and tested with VIDAS® PTH (1-84) (PTH) assay (Y) and another commercially available PTH (1-84) immunoassay (X).

According to the KDIGO guidelines, the PTH concentration of the patients treated by dialysis should be maintained within two and nine times the upper normal limit of the assay. For information, the highest reference value established in a clinical study was 44.6 pg/mL, meaning an interval of PTH values between 89.2 and 401.4 pg/mL. Correspondence was defined by a concordant interpretation between both methods.

The global (all subjects and time points mixed) frequency and percent of concordant interpretation between the 2 assays regarding the KDIGO criterion was calculated.

Criteria	% Concordance
Concordant interpretation between the 2 assays regarding the KDIGO criterion	88.89 %

Waste Disposal

Dispose of used or unused reagents, as well as any other contaminated disposable materials, following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced, according to their nature and degree of hazardousness, and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

Literature References

1. Camacho PM *et al.* – American Association of Clinical Endocrinologists and American College of Endocrinology Clinical Practice Guidelines for the Diagnosis and Treatment of Postmenopausal Osteoporosis - 2016. – *Endocrine Practice* – 2016 September; 22 (Suppl 4): 1-42.
2. Khan AA *et al.* – Primary hyperparathyroidism: review and recommendations on evaluation, diagnosis, and management. A Canadian and international consensus. – *Osteoporos Int.* – 2017; 28: 1-19.
3. Wilhelm SM *et al.* – The American Association of Endocrine Surgeons Guidelines for Definitive Management of Primary Hyperparathyroidism. – *JAMA Surg.* – 2016 Oct; 151 (10): 959-68.
4. Bollerslev J *et al.* – European Society of Endocrinology Clinical Guideline: Treatment of chronic hypoparathyroidism in adults. – *European Journal of Endocrinology* – 2015; 173, G1-20.

5. KDIGO 2017 Clinical Practice Guideline Update for the Diagnosis, Evaluation, Prevention, and Treatment of Chronic Kidney Disease-Mineral and Bone Disorder (CKD-MBD). – *Kidney International Supplements* – 2017 July; 7 (1): 1-59.
6. WORLD HEALTH ORGANIZATION, USE OF ANTICOAGULANTS IN DIAGNOSTIC, LABORATORY INVESTIGATIONS, 2002, WHO/DIL/LAB/99.1 Rev.2
7. WORLD HEALTH ORGANIZATION, Parathyroid Hormone 1-84, human, recombinant (1st International Standard) - WHO_BS_09.2115.

Index of Symbols

Symbol	Meaning
	Catalogue number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Date of manufacture

Limited Warranty

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the "System") other than as set forth in the IFU.

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Release Date	Part Number	Change Type	Change Summary
2018-10	050049-01	N/A	Not applicable (First publication)

Release Date	Part Number	Change Type	Change Summary
2018-12	050049-02	Correction	PERFORMANCE - Method comparison
2019-11	050049-03	Technical	Content of the Kit (30 tests) / Warnings and Precautions
2022-02	050049-04	Technical change	Materials and Disposables Required but Not Provided

BIOMERIEUX, the BIOMERIEUX logo, SPR and VIDAS are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

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CLSI yra „Clinical Laboratory and Standards Institute, Inc.“ priklausantis prekės ženklas.