Quill™ Knotless Tissue-Closure Device
Variable-Loop-Design
Comprised of Monoderm™ (PGA-PLCL)
Synthetic Absorbable
Surgical Suture Material

EN
CAUTION: Federal (USA) law restricts this device to sale and use by, on or the order of, a physician. This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device is not permitted. Reprocessing and/or reuse after expiration of the date of manufacture may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize this device.

DESCRIPTION
The Quill™ Knotless Tissue-Closure Device, Variable Loop-Design consists of barbed suture material, armed with a surgical needle on one end and a loop at the opposite end. The Quill™ Device is designed to suture with a closed loop at one end and an unidirectional barbed section on the other end. The Quill™ Device barbs are oriented in one direction to allow tissue approximation without the need for the surgical knots.
The Quill™ Device is comprised of sterile synthetic absorbable monofilament suture prepared from a copolymer of glycolide and ε-caprolactone. The suture is provided dyed violet or undyed. The pigment for the violet dyed suture is D&C Violet No. 2. This polymer has been found to be nonallergenic, nonpyrogenic, and elicits only a slight tissue reaction during absorption.

While the formation of barbs in the Quill™ Device reduces the tensile strength relative to non-barbed suture material of the same size, long term studies have shown that the active stress bearing portion of the suture, the strength of the Quill™ Device can be compared to USP knot strength of non-barbed sutures. The Quill™ Device is available in suture diameters as small as 3-0. The Quill™ Device suture material, prior to barbing, except for minor variation in suture diameter with a maximum average of 0.1 mm.

USP SUTURE SIZE	SUTURE TENSILE STRENGTH (kgf)	QUILL™ KNOTLESS TISSUE-CLOSURE DEVICE EQUIVALENT SIZE
0	3.90	1
2-0	2.68	2-0
3-0	1.77	3-0
4-0	0.95	4-0

ACTIONS
The following characteristics describe the in vivo performance of absorbable sutures: first, tensile strength retention and second, the absorption rate (loss of mass). Quill™ MONODERM™ device has been formulated to minimize the variability of these characteristics and to provide wound support through the critical wound healing period.

In vivo studies demonstrate that Quill™ MONODERM™ device retains approximately 62% of its original tensile strength at 14 days post-implantation. The Quill™ MONODERM™ device is designed for use in soft tissue approximation. Absorption of Quill™ MONODERM™ device is essentially complete between 90 and 120 days.

INDICATIONS
Quill™ MONODERM™ device is indicated for use in soft tissue approximation where the use of absorbable reported is appropriate.

CONTRAINDICATIONS
The Quill™ MONODERM™ device is not to be used where extended approximation of tissue under stress is required. The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/contraction within the loop.

WARNINGS
The Quill™ MONODERM™ device should be deployed by users familiar with surgical procedure, techniques and tissues. Physicians should consider the quantity and quality of tissue when the loop will be closed. The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/contraction within the loop.

Users should be familiar with surgical procedure and techniques involving absorbable sutures before employing Quill™ MONODERM™ device for wound closure, as risk of wound dehiscence may vary with the technique. The Quill™ MONODERM™ device should be used in accordance with the instructions for use. The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/contraction within the loop.

As with any foreign body, prolonged contact of any suture with soft tissues, such as those found in the urinary or biliary tract may result in calculus formation. As an absorbable suture, the Quill™ MONODERM™ device may not function as a foreign body. Absorbable surgical practice should be followed for the management of contracted or infected wounds.

As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the physician. The use of the device with open wound expansion, stretching or distortion or which may require additional support.

The safety and effectiveness of Quill™ MONODERM™ device have not been established for use in facial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, ligament, tendon, nerve, muscle, tissue, meningeal, skin, or for use in Small bowel obstruction (SBO) including volvulus, bowel infarction, and significant morbidity, has been reported. The use of this device in the urinary or biliary tract may result in calculus formation, as with any foreign body, prolonged contact of any suture with soft tissues, such as those found in the urinary or biliary tract may result in calculus formation.

Case should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.

PRECAUTIONS
The Quill™ Device contains a loop end and barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and the suture. The Quill™ Device should be used in accordance with the instructions for use. The use of this device may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/contraction within the loop.

When using the Quill™ Device subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.

Case should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with sharp instruments, such as forceps. Do not attempt to remove/minimize the in polymer by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound dehiscence are typical of any suture material associated with any suture and hence are also potential complications associated with Quill™ Device.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

ADVERSE REACTIONS
Adverse effects associated with the use of this device may include, wound dehiscence, failure to provide adequate wound support, infection, malnutrition or debilitated patients or in patients suffering from conditions which may delay wound healing, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calcium formation in urinary and biliary tracts when prolonged contact with self solutions occurs and transitory local infection as the wound site. Broken needles may result in extended or persistent reactions or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

APPLICATION
As required per surgical procedure. Secure nonbarbed loop portion to a robust tissue. Start by taking a bite of the designated anchoring tissue and then pass the needle through the loop and pull to tighten. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in designated sharps container.

HOW SUPPLIED
The Quill™ Device is available sterile, undyed in various lengths and USP diameter sizes 3-0 through 1 (metric sizes 2 through 4). Each Quill™ Device is single armed with various needles.

SYMBOLS USED FOR LABELING
Do Not Resterilize
Do Not Reuse
Do Not Use if package is damaged
Sterilized using ethylene oxide

CAUTION: Federal (USA) law restricts this device to sale and use by, on or the order of, a physician.

Authorized Representative in the European Community:

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