Product Specification Sheet

OptiLube Active
1160, 1161

Distributed by
optimummedical®

● General information

General description: OptiLube Active (with Lidocaine and Chlorhexidine) is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.

Intended use: OptiLube Active is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes.

OptiLube Active is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication. Additionally, OptiLube Active contains the following ancillary medicinal substances; antiseptic (Chlorhexidine) and local anaesthetic (Lidocaine) to help reduce the risk of infection and pain for the patient.

Lifetime of device: Single use

● Storage Conditions

- Store between 5-30°C (41-86°F) until expiry date.
- Keep dry and out of direct sunlight.
- Unopened packs have a 3-year shelf life.
- Do not use after expiry date.

● Packaging configuration:

1 syringe in medical grade paper-film blister pack	95.3 x 14.5 [6mL syringe] 95.7 x 18.4 [11mL syringe]
10 syringes per inner box	22.5 x 20 x 8.5 cm [1160 inner] 25 x 20 x 9.8 cm [1161 inner]
20 inner boxes per outer box	46.5 x 42.5 22.5 cm [1160 outer] 52 x 52 x 22.5 cm [1161 outer]

Product code	Format
1160	6mL syringe
1161	11mL syringe

● Physical properties:

- Visual Appearance: Colourless, transparent gel
- Density: 1.0 +/- 0.1 g/mL
- Viscosity: 950 – 3700 cP
- pH: 5.20 – 6.20
- Solubility in water: Soluble
- Sterilisation: Gamma irradiation
- Latex free

● Ingredients:

100g of OptiLube Active contains:

- 0.05g Chlorhexidine gluconate
- 2g Lidocaine hydrochloride
- Purified water, propylene glycol, hydroxyethyl cellulose, methyl hydroxybenzoate, propyl hydroxybenzoate

Syringe:

- Polyethylene/polypropylene copolymer
- Polyisoprene (piston)

Company Registration No: 6779883

optilube.TM
active

Instructions for use

PRODUCT REF/CONTENTS

1160 OptiLube Active 6ml Prefilled Syringe

1161 OptiLube Active 11ml Prefilled Syringe

This is an OptiLube Active product information leaflet for the user.

OptiLube Active has been prescribed for you by a Healthcare Professional. Do not share it with anyone else.

Read all of the product information leaflet prior to use and retain for your information. If you are unsure about how to use this product, or whether it is suitable for you, please seek medical advice.

PRODUCT DESCRIPTION

OptiLube Active (with Lidocaine and Chlorhexidine) is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.

INGREDIENTS

100g of OptiLube Active contains:

Active Ingredients

• 0.05g Chlorhexidine Gluconate; • 2g Lidocaine Hydrochloride.

Other Ingredients

• Purified Water; • Propylene Glycol; • Hydroxyethyl cellulose;
• Methyl Hydroxybenzoate; • Propyl Hydroxybenzoate.

INTENDED PURPOSE

OptiLube Active is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication. Additionally, OptiLube Active contains the following ancillary medicinal substances; antiseptic (Chlorhexidine) and local anaesthetic (Lidocaine) to help reduce the risk of infection and pain for the patient.

INDICATION FOR USE

OptiLube Active is designed to lubricate the urethra prior to insertion of a urethral catheter and other urological medical devices including cystoscopes.

CONTRAINDICATIONS

OptiLube Active must not be used in patients:

- With known allergies or hypersensitivity to any of the ingredients.
- Who have ever had a reaction to Lidocaine or Chlorhexidine.
- Who have damaged or bleeding mucous membranes because of the risk of systemic absorption of the Lidocaine and Chlorhexidine.

Do not use in children below 2 years of age.

Do not use for intravenous (IV) or intramuscular (IM) injections.

Do not use orally. If swallowed, please seek medical advice.

Do not use in eyes. If it comes into contact with the eyes, please seek medical advice.

WARNINGS AND PRECAUTIONS

OptiLube Active must only be used under the supervision of Healthcare Professionals in

accordance with local guidelines, policies and procedures.

OptiLube Active is not suitable if a patient:

- Has a damaged or bleeding urethra.

Care should be taken when using the jelly with patients who:

- Have heart problems or taking medication for treating irregular heartbeat;
- Have liver problems;
- Are epileptic.

OptiLube Active must not be used at the same time as other medical devices or medicines containing Lidocaine or other local anaesthetic.

Patients may experience a slight stinging on application of the jelly, but this usually stops once the anaesthetic starts to work. Encourage patients to report any reaction to the jelly.

Keep out of reach of children.

Ensure OptiLube Active is suitable for the intended use and compatible with other medical devices to be used in conjunction with it.

This is a single-use device. Re-use of this device may result in patient infection/ cross-contamination.

Re/sterilisation, reprocessing, cleaning and disinfection may also compromise the product characteristics, resulting in trauma or infection to the patient.

OptiLube Active does not contain natural rubber latex.

PREGNANCY AND LACTATION

Only use during pregnancy and breastfeeding under the direction of a doctor.

Always ask the patient if they are pregnant prior to using the jelly.

EFFECTS ON THE ABILITY TO DRIVE AND USE MACHINERY

The ability to drive and operate machinery may be slightly impaired after the use of lubricant jelly with Lidocaine. If affected, patients should be advised not to drive or use machinery.

INTERACTIONS OF MEDICINAL SUBSTANCES

AND SUBSTANCES THAT ARE ABSORBED

Depending on the absorption of Lidocaine, these interactions can be seen when used with the following medications:

Propranolol: Reduction in plasma clearance of Lidocaine;

Cimetidine: Reduction in plasma clearance of Lidocaine;

Antiarrhythmic products: Increase in the toxicity of Lidocaine;

Phenytoin or barbiturates: Reduction in plasma levels of Lidocaine.

Specified interactions can be seen in the long-term use and repeated high doses.

When administered as recommended, there is no clinically significant interactions reported.

UNDESIRABLE SIDE EFFECTS

Like any medication, OptiLube Active may cause side effects in some people. Side effects to this product must be documented in patient records.

In rare cases local hypersensitivity reactions may occur, such as redness, stinging, blistering or itching and/or systemic reactions to Lidocaine and/or Chlorhexidine.

There is also a risk of severe reactions including drop in blood pressure, dizziness, nausea, shortness of breath, bradycardia, convulsions and anaphylactic shock. If this occurs, stop using **OptiLube Active** and consult your Healthcare Professional.

Any serious incident or malfunction which may affect the safety of **OptiLube Active** should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

DOSAGE

Dosage recommendations:

FOR ADULTS: Max. 800 mg **Lidocaine** during a 24 hour period (x 3 syringe for 11 ml) (x 6 syringe for 6 ml),

FOR CHILDREN (between 2 and 15 years): max. 0.3ml gel/kg B.W ($\approx$ 6mg **Lidocaine**/kg) is recommended per procedure and no more than four doses should be administered within 24 hours.

In children below 2 years lubricant gel with **Lidocaine** must not be used.

OVERDOSE

Excessive dosage or short intervals between doses can result in high plasma levels and serious adverse effects.

OptiLube Active must not be used at the same time as any other medical device or medicines containing a local anaesthetic agent.

In the event of excessive absorption of **Lidocaine** into the bloodstream, symptoms may include central nervous system effects and cardiovascular reactions.

In the event of excessive use or use outside of these guidelines please seek medical advice.

HOW TO USE

The decision as to which size (6ml or 11ml) of **OptiLube Active** sterile lubricating jelly pre-filled syringe to be used is taken by the Healthcare Professional. It is recommended that this product is used as part of an aseptic technique.

- Cleanse the opening to the urethra and surrounding area prior to use;
- Evenly peel back the paper backing and remove the sterile syringe;
- Remove the cap from the end of syringe;
- Apply a drop of jelly to urethral opening to make initial insertion easier;
- Insert the nozzle into the urethral opening and press the plunger of the syringe slowly to release the necessary amount of **OptiLube Active**.

NOTE: **OptiLube Active** should be applied to the urethra, not directly onto the device;

- The lubrication characteristics of the jelly start to take effect at the time of application. The onset of anaesthetic effect takes 3-5 minutes.

Do not use if packaging is damaged or unintentionally opened prior to use.

STERILISATION

OptiLube Active is supplied sterile.

Sterilised with gamma radiation after the packaging process.

Do not resterilise.

STORAGE AND HANDLING

The scale on the syringe is for the orientation of the user.

It does not have a measurement function.

Store between 5-30°C (41-86°F) until expiry date.

Keep dry and out of direct sunlight.

Unopened packs have a 3-year shelf life. Do not use after expiry date.

DISPOSAL

OptiLube Active must be disposed of according to local policies and procedures, including any remaining fluid in the syringe.

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CH REP

Decomplic AG
Freiburgstrasse 3
3010 Bern
Switzerland

**UK Responsible Person
and Distributor:**

Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, United Kingdom
T: +44(0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

EU Importer:
Optimum Medical Solutions GmbH,
c/o Grant Thornton AG,
Senckenberganlage 19,
60325 Frankfurt am Main, Germany.
Email: enquiries@optimummedical.de

STERILE R

MD

LATEX FREE

CE 2292

1160-1161-A0140-Aug22-V2.0

Black

BACK COVER

optilube.™
active

REF 1160
REF 1161

Rev. date: 11.12.2023

CE 2292

Does not contain natural rubber latex.
No contiene látex de caucho natural.


STERILE R

MD

LATEX FREE

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug, Switzerland

Sponsor:

Optimum Medical Australasia Pty Ltd,
Level 4, 470 Church Street, North Parramatta,
New South Wales 2151, Australia

Istem Medikal Tıbbi Cihaz ve San. Tic. Ltd. Şti.,
Anadolu OSB Mah. 29. Ekim Cad.,
No: 41 Malıköy, Sincan, Ankara, Turkey
Made in Turkey

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague, The Netherlands

optimummedical®

UK Responsible Person and Distributor:
Optimum Medical Solutions Ltd
Tennant Hall, Blenheim Grove,
Leeds, LS2 9ET, UK
T: +44(0) 845 643 5479
enquiries@optimummedical.co.uk
www.optimummedical.co.uk

1160-1161-A0140-Dec23-V3.0

Black

FRONT COVER

optilube.™

active

Instructions for use

GBR INSTRUCTIONS FOR USE

PRODUCT REF/CONTENTS

1160 **OptiLube Active** 6ml Prefilled Syringe

1161 **OptiLube Active** 11ml Prefilled Syringe

This is an **OptiLube Active** product information leaflet for the user.

OptiLube Active has been prescribed for you by a Healthcare Professional. Do not share it with anyone else. Read all of the product information leaflet prior to use and retain for your information. If you are unsure about how to use this product, or whether it is suitable for you, please seek medical advice.

PRODUCT DESCRIPTION

OptiLube Active (with **Lidocaine** and **Chlorhexidine**) is a sterile, single-use, pre-filled syringe with water-soluble jelly available in 6ml and 11ml syringes.

INGREDIENTS

100g of **OptiLube Active** contains:

Active Ingredients

- 0.05g **Chlorhexidine** Gluconate; • 2g **Lidocaine** Hydrochloride.

Other Ingredients

- Purified Water; • Propylene Glycol; • Hydroxyethyl cellulose;
- Methyl Hydroxybenzoate; • Propyl Hydroxybenzoate.

INTENDED PURPOSE

OptiLube Active is a medical device that helps to prevent trauma being caused to the patient during catheterisation procedures or other urethral procedures by effective lubrication. Additionally, **OptiLube Active** contains the following ancillary medicinal substances; antiseptic (**Chlorhexidine**) and local anaesthetic (**Lidocaine**) to help reduce the risk of infection and pain for the patient.

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- With known allergies or hypersensitivity to any of the ingredients.
- Who have ever had a reaction to **Lidocaine** or **Chlorhexidine**.
- Who have damaged or bleeding mucous membranes because of the risk of systemic absorption of the **Lidocaine** and **Chlorhexidine**.

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WARNINGS AND PRECAUTIONS

OptiLube Active must only be used under the supervision of Healthcare Professionals in accordance with local guidelines, policies and procedures.

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OptiLube Active must not be used at the same time as other medical devices or medicines containing **Lidocaine** or other local anaesthetic.

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Keep out of reach of children.

Ensure **OptiLube Active** is suitable for the intended use and compatible with other medical devices to be used in conjunction with it.

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The ability to drive and operate machinery may be slightly impaired after the use of lubricant jelly with **Lidocaine**. If affected, patients should be advised not to drive or use machinery.

INTERACTIONS OF MEDICINAL SUBSTANCES AND SUBSTANCES THAT ARE ABSORBED

Depending on the absorption of **Lidocaine**, these interactions can be seen when used with the following medications:

Propranolol: Reduction in plasma clearance of **Lidocaine**;

Cimetidine: Reduction in plasma clearance of **Lidocaine**;

Antiarrhythmic products: Increase in the toxicity of **Lidocaine**;

Phenytoin or barbiturates: Reduction in plasma levels of **Lidocaine**.

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UNDESIRABLE SIDE EFFECTS

Like any medication, **OptiLube Active** may cause side effects in some people. Side effects to this product must be documented in patient records.

In rare cases local hypersensitivity reactions may occur, such as redness, stinging, blistering or itching and/or systemic reactions to **Lidocaine** and/or **Chlorhexidine**.

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Maliköy, Sincan, 41:Ekim Cad., No. 29, St., Anadolu OSB Mah. Istem Medikal Medical Devices and San Ankara, Turkey

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28 Gotthardstrasse

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NOTES

Handwritten notes area with horizontal dashed lines for writing.

[illegible]

constantly creating better

11mle
10 SYRINGES

sterile lubricating jelly with local anaesthetic (**Lidocaine**)
and antiseptic (**Chlorhexidine**) for urethral application
gel lubrifiant stérile avec anesthésiant (**Lidocaine**)
local et antiseptique (**Chlorhexidine**) pour application uréthrale
steriles gleitgel mit lokalanästhetikum (**Lidocain**) & antiseptikum
(**Chlorhexidin**) für die urethrale anwendung

11mle
10 SYRINGES
NOT TO BE USED FOR
I/M OR I/V INJECTIONS

Composition (in each 100g)	
Hydroxyethyl cellulose	0.025g
Propylene Glycol	0.060g
Purified Water	0.050g
Propyl Hydroxybenzoate	2.000g
Chlorhexidine Gluconate	
Lidocaine Hydrochloride	

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ISO 13485:2016

Scope: Design, manufacturing and sales of sterile water soluble lubricant gel, sterile lubricant gel with lidocaine, sterile lubricant gel with lidocaine (CHG Free), sterile injectable implant for vur and sui treatment, sterile 0.9% sodium chloride solution, sterile sodium hyaluronate solution, sterile sodium hyaluronate gel, sterile hydrophilic urinary catheter, sterile antibacterial hydrophilic urinary catheter, sterile nelaton catheter, sterile instillation catheter, sterile inhalation water, sterile embolization application set, sterile intravascular catheters, sterile retrievable vena cava filter, sterile amorphous hydrogel, sterile filter tur collector set, surgical mask


Walzmühler Strabe 18 19073, Dümmer, Germany

Hauptstrasse 3A 19073 Dümmer, Germany

Hereby, AKSSERT Audit and Certification Ltd. Co., certifies that the above stated company gave the appropriate management system according to the requirements of the above standard. This certificate valid for 3 years since the decision date as long as the system is effectively maintained and surveillance audits are carried out. The validity of certificate can be checked through www.akssert.com, www.jas-anz.org/register. The Certificate is property of AKSSERT Audit and Certification Ltd. Co. and shall be returned if requested.

The reference standard is ISO 13485:2016

AKSSERT Audit and
Certification Ltd. Co.

Certificate Number : 85209

Registration Date : 27.01.2016

Previous Certification Cycle Expiry Date: 26.01.2022

Recertification Audit Date: 13.01.2022

Reissue Date : 24.01.2022

Revision Date/No : -

Expiry Date : 26.01.2025

Address: Mutlukent Mahallesi 1920. Cadde No:28 Çankaya/ANKARA - TÜRKİYE

Phone: +90 312 284 99 44(pbx)

E-mail: info@akssert.com Web: www.akssert.com

EC CERTIFICATE

Full Quality Assurance System Medical Devices Directive 93/42/EEC Annex II (Excluding Section 4)

M.2021.106.14615-1 Design Examination Certificate Was Prepared for Class III Products Defined in This Certificate.

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Related Directives and Annex : 93/42/EEC Medical Devices Directive - Annex II (Excluding Section 4)

Product : Sterile Lubricant Gel with Lidocaine - Class III

GMDN : 60796

This certificate has been issued based on Ministry of Health's E-61749811-511.14-437042 numbered scientific opinion taken on 20.05.2021 according to 93/42/EEC Annex I Art .7.4

Certificate Number : M.2021.106.14615

Report Number : MD.3200.IB

Initial Assessment Date : 12.03.2021

Registration Date : 24.05.2021

Revision Date /No : -

Expiry Date : 27.05.2024

UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

UDEM hereby declares that the requirements of Annex II, excluding section 4 of the 93/42/EEC Directive have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance audits, defined by Annex II, section 5 of the forementioned directive. According to Annex II, section 4 an EC design- examination certificate is required for placing the Class III devices on the market. UDEM's responsibility for class I devices covered by the EC certificate is limited to manufacturing issues related to safeguarding and maintaining sterile conditions, if the device is sterile; and manufacturing issues related to product's conformity with metrological requirements, if it has measurement function. This certificate remains as the property of UDEM International Certification Auditing Training Centre Industry and Trade Inc. Co. to whom it must be returned upon request. The above named company and UDEM must keep a copy of this certificate for 5 years from the registration of the certificate. Usage of the CE mark is under the responsibility of the manufacturer with the completion of EC Declaration of Conformity. The above mentioned company must notify all changes related with the approved product to UDEM. If UDEM will not renew the expiry date of this certificate in question, the mentioned

Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya – Ankara – TURKEY

Phone: +90 0312 443 03 90 Fax: +90 0312 443 03 76

E-mail: info@udemltd.com.tr www.udem.com.tr

EC DESIGN EXAMINATION CERTIFICATE

93/42/EEC Directive of Medical Devices Annex II, Section 4

With the expire of the certificate M.2021.106.14615 the validity of the certificate
M.2021.106.14615-1 will also end.

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Related Directives and Annex : 93/42/EEC Medical Devices Directive - Annex II (Section 4)

Product : Sterile Lubricant Gel with Lidocaine - Class III

GMDN : 60796

This certificate has been issued based on Ministry of Health's E-61749811-511.14-437042 numbered scientific
opinion taken on 20.05.2021 according to 93/42/EEC Annex I Art .7.4

Certificate Number : M.2021.106.14615-1
Report Number : MD.3200.IB
Initial Assessment Date : 12.03.2021
Registration Date : 24.05.2021
Revision Date /No : -
Expiry Date : 27.05.2024

UDEM International Certification
Auditing Training Centre Industry
and Trade Inc. Co.

The EC desing examination certificate refers to the above mentioned product, UDEM hereby declares that the requirements of Annex II, section 4 of the 93/42/EEC Directive have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance audits, defined by Annex II, section 5 of the aforementioned directive. This certificate remains as the property of UDEM International Certification Auditing Training Centre Industry and Trade Inc. Co. to whom it must be returned upon request. The above named company and UDEM must keep a copy of this certificate for 5 years from the registration of the certificate. Usage of the CE mark is under the responsibility of the manufacturer with the completion of EC Declaration of Conformity. The above mentioned company must notify all changes related with the approved product to UDEM. If UDEM will not renew the expiry date of this certificate in question, the mentioned company should stop placing the product on the market. The validity of the certificate can be checked through www.udem.com.tr.

Address: Mutlukent Mahallesi 2073 Sokak (Eski 93 Sokak) No:10 Çankaya – Ankara – TURKEY

Phone: +90 0312 443 03 90 Fax: +90 0312 443 03 76

E-mail: info@udemltd.com.tr www.udem.com.tr

EC DECLARATION OF CONFORMITY

MANUFACTURER : İSTEM MEDİKAL TIBBİ CİHAZ ve SAN. TİC. LTD. ŞTİ.

ADDRESS : ANADOLU OSB MAH. 29 EKİM CAD. NO:41 MALİKÖY-SİNCAN/ANKARA
TURKEY

PRODUCT : Sterile Lubricant Gel with Lidocaine

MODEL : It is given in the attached list.

GMDN: 60796

CLASSIFICATION: CLASS III, ANNEX-IX, RULE 5, RULE13

CONFORMITY ASSESSMENT ROUTE: ANNEX II (including section 4)

APPLICABLE STANDARDS: EN ISO 13485:2016, EN ISO 11137-1:2015, EN ISO 11137-2:2015, EN ISO 11137-3:2017, EN 868-2:2017, EN 868-3:2017, EN ISO 11737-1:2018, EN ISO 11737-2:2020, EN ISO 14971:2019, EN ISO 10993-1:2021, EN ISO 10993-3:2015, EN ISO 10993-5:2010, EN ISO 10993-10:2014, EN ISO 15223-1:2016, EN ISO 19011:2018, EN ISO 14644-1:2016, EN ISO 14644-2:2016, EN ISO 14644-3:2019, EN ISO 14155:2020, EN ISO 11607-1:2020, EN ISO 11607-2:2020, EN ISO 7886-1:2018, EN ISO 80369-7:2017, EN 980:2011, EN 1041:2008(A1:2014), EN ISO 17050-1:2013, EN ISO 17050-2:2007, TS EN 62366-1:2015, MDD 93/42 EEC, MEDDEV 2.7.1 Rev. 4:2016, MEDDEV 2.12-1 rev8:2013, MEDDEV 2.12-2 rev2:2012, MEDDEV 2.4/1:rev9:2010, MEDDEV 2.5/5:Rev 3:1998, ASTM F1980-16

Pharmacopoeia monographs: USP 35/Page 1191, EP 07/2010:0658, EP 09/2017 2.2.1-2.2.2, EP 09/2017 2.9.17, EP 09/2017 2.2.10, EP 09/2017 2.6.1, EP 09/2017:0409, EP 09/2017:0431

WE HEREWITH DECLARE THAT THE ABOVE MENTIONED PRODUCTS MEET THE PROVISIONS OF THE COUNCIL DIRECTIVE 93/42/EEC AND 2007/47/EC FOR MEDICAL DEVICES AND THESE PRODUCTS ARE NOT USED FOR MEDICINE. ALL SUPPORTING DOCUMENTATION IS RETAINED UNDER THE PREMISES OF THE MANUFACTURER.

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Registration Date: 24.05.2021

Expiry Date: 27.05.2024

EC DESIGN EXAMINATION CERTIFICATE NUMBER: M.2021.106.14615-1

Registration Date: 24.05.2021

Expiry Date: 27.05.2024

PLACE, DATE OF ISSUE: ANKARA, 25.05.2021

SIGNATURE: LEVENT HAYYAOĞLU
GENERAL MANAGER

NAZMİYE CUMALI
QUALITY ASSURANCE MANAGER

İSTEM MEDİKAL
TIBBİ CİHAZ ve SAN. TİC. LTD. ŞTİ.
Anadolu OSB Mah. 29 Ekim Cad.
No: 41 Malıköy / Sincan / ANKARA
Tel: (0312) 344 55 62 Fax: 394 55 64
Sincan V.D. 481 008 4741

MODELS

Trade Mark	Barcode	Reference Number	Definition
OptiLube Active	5060222470604	1160	OptiLube Active Sterile Lubricant Gel with Lidocaine-6ml Prefilled Syringe (package)
	5060222470611	1160	OptiLube Active Sterile Lubricant Gel with Lidocaine-6ml Prefilled Syringe (box)
	5060222470628	1160	OptiLube Active Sterile Lubricant Gel with Lidocaine-6ml Prefilled Syringe (outer box)
	5060222470635	1161	OptiLube Active Sterile Lubricant Gel with Lidocaine-11ml Prefilled Syringe (package)
	5060222470642	1161	OptiLube Active Sterile Lubricant Gel with Lidocaine-11ml Prefilled Syringe (box)
	5060222470659	1161	OptiLube Active Sterile Lubricant Gel with Lidocaine-11ml Prefilled Syringe (outer box)