XIENCE PRO[™]S

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XIENCE PRO™S

Everolimus Eluting Coronary Stent System

INFORMATION FOR PRESCRIBERS

Table of Contents

1.0	DEVICE DESCRIPTION
	Table 1-1: Product Name and Sizes
	Table 1-2: Stent and Coating Composition of the XIENCE PRO ^S Everolimus Eluting Coronary Stents
	Table 1-3: <i>In vitro</i> Device Specifications
2.0	INDICATIONS
3.0	INTENDED USE
4.0	INTENDED USER AND USER TRAINING
5.0	INTENDED CLINICAL BENEFIT
6.0	CONTRAINDICATIONS
7.0	EXPECTED LIFETIME OF THE DEVICE
8.0	WARNINGS
9.0	LIMITATIONS
10.0	PRECAUTIONS
10.1	Stent Handling
10.2	Stent Placement
10.3	Use in Conjunction with Other Procedures
10.4	Stent / System Removal
10.5	Post-Implant
10.6	Use in Special Populations
10.6.1	Pregnancy
10.6.2	Lactation
10.6.3	Pediatric Use
10.7	Drug Interactions
10.7.1	Interactions with Drugs or Other Substances
10.7.2	Immune Suppression Potential
10.7.3	Lipid Elevation Potential
10.7.4	Pregnancy
10.7.5	Lactation
10.8	Carcinogenicity, Genotoxicity, and Reproductive Toxicity
10.9	Restricted Substances
10.10	Magnetic Resonance Imaging (MRI) Safety Information
11.0	POTENTIAL ADVERSE EVENTS
12.0	REPORTING OF ADVERSE EVENTS
13.0	SUMMARY OF SAFETY AND CLINICAL PERFORMANCE
14.0	PATIENT SELECTION AND TREATMENT
14.1	Individualization of Treatment
15.0	CLINICIAN USE INFORMATION
15.1	Inspection Prior to Use
15.2	Materials Required
15.3	Preparation
15.3.1	Packaging Removal
15.3.2	Guide Wire Lumen Flush
15.3.3	Delivery System Preparation
15.4	Delivery Procedure
15.5	Deployment Procedure
15.6	Removal Procedure
15.7	Post-Deployment Dilatation of Stent Segments
15.8	Patient Implant Card
16.0	SPIRIT AND XIENCE FAMILY OF CLINICAL TRIALS
16.1	Pre-Market Clinical Trials
	Table 16.1-1: SPIRIT Family of XIENCE V Clinical Trial Designs (Pre-Market)
	Table 16.1-2: SPIRIT Family of Clinical Trials Angiographic Results (Pre-Market)
	Table 16.1-3: SPIRIT Family of Clinical Trials Principal 1-Year Clinical Outcomes (Pre-Market)
	Table 16.1-4: SPIRIT Family of Clinical Trials Principal Clinical Outcomes from Latest Follow-up (Pre-Market)
16.2	Post-Market Clinical Trials
	Table 16.2-1: SPIRIT and XIENCE V Family Clinical Trial Designs (Post-Market)
	Table 16.2-2: SPIRIT and XIENCE V Family of Trials Principal Clinical Outcomes (Post-Market)
16.3	One-Month (as Short as 28 Days) DAPT Indication
	Table 16.3-1: Primary and Powered Secondary Endpoints Results – All Death / All MI and BARC 2–5 Bleeding from 1–6 Months in 1-Month Clear Patients
	Table 16.3-2: PS Stratified Analysis of BARC 3–5 Bleeding from 1–6 Months in 1-Month Clear Patients
	Table 16.3-3: Stent Thrombosis from 1–6 Months in 1-Month Clear Patients (Non-PS Stratified)
17.0	HOW SUPPLIED

1.0 DEVICE DESCRIPTION

The XIENCE PRO™S Everolimus Eluting Coronary Stent System (EECSS) is a medical device, which includes:

A pre-mounted L-605 cobalt chromium (CoCr) alloy XIENCE PRO^S stent which is coated with two layers, a primer layer and a drug matrix layer. The primer layer is composed of poly n-butyl methacrylate (PBMA), which is an inactive polymer to improve the adhesion between the stent and drug matrix layer. The drug matrix layer is composed of a copolymer poly vinylidene fluoride-co-hexafluoropropylene (PVDF-HFP) blended with the anti-proliferative drug everolimus.

The product family consists of:

Table 1-1: Product Name and Sizes

Product Name	Basic UDI-DI	Stent Diameter (mm)	Stent Length (mm)
XIENCE PRO ^S	8717648DES0001TC	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 4.0	8, 12, 15, 18, 23, 28, 33, 38

The L-605 cobalt chromium (CoCr) alloy stent mass, available dose of everolimus and PVDF-HFP and PBMA content on the stent vary by size as follows:

Table 1-2: Stent and Coating Composition of the XIENCE PRO^S Everolimus Eluting Coronary Stents

Stent Diameter (mm)	Stent Length (mm)	L-605 Cobalt Chromium Alloy Stent Mass (mg)	Drug Dose (µg)	PVDF-HFP (µg)	PBMA (µg)
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	8	7.8	39	189	42
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	12	11.4	58	286	64
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	15	13.8	72	350	79
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	18	16.2	85	415	93
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	23	21.0	111	543	122
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	28	24.6	131	640	144
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	33	29.4	157	769	173
2.0, 2.25, 2.5, 2.75, 3.0, 3.25	38	33.0	177	866	194
3.5, 4.0	8	11.1	53	261	59
3.5, 4.0	12	15.4	72	351	79
3.5, 4.0	15	20.6	99	485	109
3.5, 4.0	18	24.0	117	575	129
3.5, 4.0	23	29.2	145	709	159
3.5, 4.0	28	36.2	181	888	199
3.5, 4.0	33	41.3	209	1023	230
3.5, 4.0	38	46.5	236	1157	260

- Two radiopaque markers, located underneath the balloon, fluoroscopically mark the working length of the balloon and the expanded stent length.
- Two proximal delivery system shaft markers (95 cm and 105 cm proximal to the distal tip) indicate the relative position of the delivery system to the end of the brachial or femoral guiding catheter. Working catheter length is 145 cm.
- A shaft color change denotes the guide wire exit notch.

Table 1-3: *In vitro* Device Specifications

Stent Diameter (mm)	Stent Length (mm)	* Minimum Guiding Catheter Compatibility (ID)	** <i>In vitro</i> Stent Nominal Pressure		Rated Burst Pressure (RBP)		Stent Free Area (%)
			(atm)	kPa	(atm)	kPa	
2.0	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	9	912	16	1621	81
2.25	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	9	912	16	1621	83
2.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	9	912	16	1621	84
2.75	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	12	1216	16	1621	86
3.0	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	12	1216	16	1621	87
3.25	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	12	1216	16	1621	88
3.5	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	12	1216	16	1621	84
4.0	8, 12, 15, 18, 23, 28, 33, 38	5F (0.056" / 1.42 mm)	12	1216	16	1621	86

* See individual manufacturer specifications for (F) equivalent.

** Ensure full deployment of the stent (see *Section 15.5 Deployment Procedure*). Deployment pressures should be based on lesion characteristics.

2.0 INDICATIONS

The XIENCE PRO™S Everolimus Eluting Coronary Stent System is indicated for improving coronary artery luminal diameter in the following:

- Patients with symptomatic ischemic heart disease due to discrete *de novo* native coronary artery lesions.
- For restoring coronary flow in patients experiencing acute myocardial infarction who present within 12 hours of symptom onset.
- For the treatment of patients with concomitant diabetes, acute coronary syndrome, dual vessel lesions (two lesions in two different epicardial vessel(s), lesions residing within small coronary vessels; lesions where treatment results in the jailing of side branches (lesions with a side branch < 2 mm in diameter or an ostial stenosis < 50%); for the treatment of elderly patients (age ≥ 65), and for treatment of both men and women.
- For treatment of patients with high bleeding risk (HBR) under dual antiplatelet therapy (DAPT) as short as 28 days.
- For the treatment of patients presenting with in-stent restenosis in coronary artery lesions; chronic total occluded coronary artery lesions (defined as coronary artery lesions with TIMI flow 0 and lasting longer than 3 months); and coronary artery bifurcation lesions.

In all cases, the treated lesion length should be less than the nominal stent length (8 mm, 12 mm, 15 mm, 18 mm, 23 mm, 28 mm, 33 mm, or 38 mm) with a reference vessel diameter of ≥ 2.00 mm and ≤ 4.25 mm.

3.0 INTENDED USE

The XIENCE PRO™S Everolimus Eluting Coronary Stent System is intended to improve coronary artery luminal diameter.

4.0 INTENDED USER AND USER TRAINING

This device is designed to be used by an interventional cardiologist / interventional radiologist / vascular surgeon, or other physician who has been trained in cardiovascular interventional techniques, in a cardiac catheterization laboratory or a surgical suite, using imaging equipment, by means of percutaneous techniques or surgical cut-down to treat atherosclerotic narrowing in coronary arteries.

5.0 INTENDED CLINICAL BENEFIT

The following are the intended clinical benefits of the XIENCE PRO™S Everolimus Eluting Coronary Stent System in coronary indications:

- Improved outcomes with revascularization versus medical treatment
- Improved outcomes with new-generation drug-eluting stents versus other percutaneous coronary intervention (PCI) devices
- No increased ischemic risk, and decreased risk of major bleeding (Bleeding Academic Research Consortium [BARC] 3–5) with shorter DAPT for HBR population

6.0 CONTRAINDICATIONS

The XIENCE PRO™S Everolimus Eluting Coronary Stent System is contraindicated for use in:

- Patients with hypersensitivity or contraindication to everolimus or structurally related compounds, or known hypersensitivity to stent components (cobalt, chromium, nickel, tungsten, methacrylic polymer, and fluoropolymer), or with contrast hypersensitivity.
- Patients who cannot tolerate, including allergy or hypersensitivity to, procedural anticoagulation or the post-procedural antiplatelet regimen.

7.0 EXPECTED LIFETIME OF THE DEVICE

The XIENCE PRO™S Everolimus Eluting Coronary Stent System is a medical device implanted to improve the coronary artery luminal diameter in partially or totally occluded stenosed coronary vessels (*de novo* lesions or restenosed vessels).

The expected lifetime of the XIENCE PRO™S Everolimus Eluting Coronary Stent System is a minimum of 10 years.

XIENCE PRO™S Everolimus Eluting Coronary Stents are not intended to be removed after implantation.

8.0 WARNINGS

- The XIENCE PRO™S Everolimus Eluting Coronary Stent and the delivery system are for single use only. Do not reuse, reprocess, or resterilize. Note the product “Use by” date on the package. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and / or delivery system and / or lead to device failure, which may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device and / or delivery system may lead to injury, illness, or death of the patient.
 - Antiplatelet therapy should be administered post-procedure (see *Section 14.1 Individualization of Treatment*).
 - This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.
 - Judicious selection of patients is necessary, since the use of this device carries the associated risk of stent thrombosis, vascular complications, and / or bleeding events.
 - It is not recommended to treat patients having a lesion that prevents complete inflation of an angioplasty balloon.
 - The XIENCE PRO™S Everolimus Eluting Coronary Stent is coated with an everolimus and polymer coating at the full implant stent length.
 - The distal and intermediate portions of the device, the tip, and tapers of the balloon are coated with HYDROCOAT™ Hydrophilic Coating.
- Refer to *Section 15.3 Preparation* for further information on how to prepare and use this device to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

9.0 LIMITATIONS

Implantation of the stent should be performed only by physicians who have received appropriate training.

10.0 PRECAUTIONS

10.1 Stent Handling

- Stent placement should be performed at centers where emergency coronary artery bypass graft (CABG) surgery can be readily performed.
- **The foil pouch is not a sterile barrier.** The inner header bag (pouch) within the foil pouch is the sterile barrier. **Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.**
- To confirm sterility has been maintained, ensure that the inner package sterile barrier has not been opened or damaged prior to use.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and balloon withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- **Do not manipulate, touch, or handle the stent**, as this may cause coating damage, contamination, or dislodgement of the stent from the delivery balloon.
- Avoid wiping the device with dry gauze or excessive wiping of the device as this may damage the device coating.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.
- Avoid soaking the XIENCE PRO™S Everolimus Eluting Coronary Stent. See instructions in *Section 15.3.2 Guide Wire Lumen Flush*.
- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon, as this may cause uneven expansion and difficulty in deployment of the stent. If gaseous medium is used and balloon rupture occurs, there is the potential of causing air embolism and / or vessel injury.
- After use, this device, its accessories and packaging should be appropriately classified for disposal (e.g., biohazard, sharps, non-hazardous waste, etc.), and carefully disposed of in compliance with facility procedures and applicable laws and regulations.

10.2 Stent Placement

- Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system.
- Do not prepare or pre-inflate the delivery system prior to stent deployment other than as directed. Use balloon purging technique described under *Section 15.3.3 Delivery System Preparation*.
- When pre-dilatation is performed, an appropriate balloon size should be used. Failure to do so may increase the difficulty of stent placement and cause procedural complications.
- The decision to pre-dilate the lesion with an appropriately sized balloon should be based on patient and lesion characteristics. Direct stenting in less complex coronary lesions with a predicate device has been shown to be as effective and safe as stenting with pre-dilation for device lengths up to 28 mm in real-world settings. If pre-dilation is performed, limit the length of pre-dilation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the implanted stent.
- When introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.
- Do not torque the catheter more than one (1) full turn.

- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent, and may cause abrupt closure / total occlusion of the vessel, requiring additional intervention (CABG surgery, further dilatation, placement of additional stents, or other).
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the stent may be damaged or dislodged during retraction back into the guiding catheter.
- Should resistance be felt at any time during removal of the undeployed coronary stent system, please refer to the steps provided in *Section 10.4 Stent / System Removal*.
- Do not expand the stent if it is not properly positioned in the vessel. (See *Section 10.4 Stent / System Removal*).
- The inflated balloon diameter of the system used to deploy the stent should approximate the diameter of the vessel. Oversizing of the stent can result in a ruptured vessel. To ensure full expansion of the stent, the balloon should be inflated to a minimum of nominal pressure.
- Do not exceed the rated burst pressure (RBP) as indicated on the product label. Monitor balloon pressures during inflation. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and if necessary, a new dilatation catheter exchanged over the guide wire to complete the expansion of the stent.
- Do not dilate the stent beyond the limits indicated in *Section 15.5 Deployment Procedure*.
- When performed, post-dilatation should be performed at high pressure with a noncompliant balloon.
- Underexpansion of the stent may result in stent movement. Care must be taken to properly size the stent to ensure that the stent is in full contact with the arterial wall upon deflation of the balloon. All efforts should be made to ensure that the stent is not underdilated. Refer to *Section 15.0 Clinician Use Information*.
- Placement of a stent has the potential to compromise side branch patency.
- Stent retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- When treating multiple lesions within the same vessel, stent the distal lesion prior to stenting the proximal lesion. Stenting in this order obviates the need to cross the proximal stent during placement of the distal stent and reduces the chance of damaging or dislodging the proximal stent.
- Ensure the stented area covers the entire lesion / dissection site and that no gaps exist between stents.
- When multiple drug-eluting stents are required, only stent materials with similar composition (e.g., XIENCE™ Everolimus Eluting Coronary Stent-Systems with the identical cobalt-chromium stent substrate and identical drug-eluting polymer coating) should be used. Potential interaction with other drug-eluting stents or coated stents has not been evaluated and should be avoided. Placing multiple stents of different metals in contact with each other may increase the potential for corrosion *in vivo*, although *in vitro* corrosion tests using an L-605 CoCr alloy stent in combination with a 316L stainless steel alloy stent did not appear to increase corrosion.
- The extent of the patient's exposure to drug and polymer is directly related to the number of stents implanted. A patient can receive up to four XIENCE PRO™S stents or other everolimus eluting coronary stents from the XIENCE Family of Stents (i.e., XIENCE V™, XIENCE PRIME™, XIENCE Xpedition™, XIENCE Alpine™, XIENCE Sierra™ and XIENCE Skypoint™) ≤ 38 mm depending on the number of vessels treated and the lesion length. Those patients receiving bailout stenting will receive additional XIENCE Family of Stents. The use of multiple XIENCE Family of Stents will result in the patient receiving larger amounts of drug and polymer.
- The safety and effectiveness of the XIENCE PRO™S stents in patients with prior brachytherapy of the target lesion, or the use of brachytherapy for restenosis in a site treated with XIENCE PRO™S stents have not been established. Both vascular brachytherapy and the implanted XIENCE PRO™S stent alter arterial remodeling. The potential combined effect on arterial remodeling by these two treatments is not known.

10.3 Use in Conjunction with Other Procedures

- While vessel preparation in complex lesions may include the use of various mechanical atherectomy devices, the safety and effectiveness of the XIENCE PRO™S stents have not been established in clinical trials with the use of either mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters.

10.4 Stent / System Removal

Stent delivery system removal prior to stent deployment:

- If removal of a stent system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent into the guiding catheter, the stent delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.
- **Withdrawal of the stent delivery system / post-dilatation balloon from the deployed stent:**
 1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10–15 seconds longer.
 2. Position the inflation device to “negative” or “neutral” pressure.
 3. Stabilize the guiding catheter position and anchor in place. Maintain the guide wire placement across stent segment.

4. Gently remove the stent delivery system / post-dilatation balloon with slow and steady pressure.
 5. Tighten the rotating hemostatic valve.
- Notes:**
- 1) If during withdrawal of the catheter from the deployed stent, resistance is encountered, use the following steps to improve balloon rewrap:
 - o Re-inflate the balloon up to nominal pressure, deflate and change pressure to neutral.
 - o Repeat steps 1 through 5 above.
 - 2) After successful withdrawal of the balloon from the deployed stent, should any resistance be felt at any time when withdrawing the stent delivery system or post-dilatation balloon into the guiding catheter, remove the entire system as a single unit.
- Failure to follow these steps and / or applying excessive force to the delivery system can potentially result in loss or damage to the stent and / or delivery system components.
 - If it is necessary to retain guide wire position for subsequent artery / lesion access, leave the guide wire in place and remove all other system components.

10.5 Post-Implant

- If necessary to cross a newly deployed stent with a guide wire, balloon, delivery system, or imaging catheters, exercise care to avoid disrupting the stent geometry.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of stents is unknown at present.
- If the patient requires imaging, see *Section 10.10 Magnetic Resonance Imaging (MRI) Safety Information*.

10.6 Use in Special Populations

10.6.1 Pregnancy

Pregnancy Category C: See *Section 10.7.4 Pregnancy*. There are no adequate everolimus or XIENCE PRO™S stent-related studies in pregnant women. Effects on the developing foetus have not been studied¹.

This product has not been tested in pregnant women or in men intending to father children. The XIENCE PRO™S stent should be used in pregnant women only if potential benefits of the stent outweigh potential risks.

10.6.2 Lactation

See *Section 10.7.5 Lactation*. It is unknown whether everolimus is distributed in human milk. A decision should be made whether or not to discontinue nursing prior to stent implantation, considering the importance of the stent to the mother.

10.6.3 Pediatric Use

The safety and effectiveness of the XIENCE PRO™S Everolimus Eluting Coronary Stent in pediatric subjects have not been established.

10.7 Drug Interactions

10.7.1 Interactions with Drugs or Other Substances

Everolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the liver and to some extent in the intestinal wall and is a substrate for the counter transporter P-glycoprotein (P-gp). Therefore, absorption and subsequent elimination of everolimus may be influenced by drugs that also affect CYP3A4 and P-gp pathways. Several drugs are known to affect everolimus metabolism, and other drug interactions may also occur. Formal drug interaction studies have not been performed with the XIENCE PRO™S stent because of limited exposure to everolimus eluted from the stent. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall, when deciding to place the XIENCE PRO™S stent in a patient taking a drug with known interaction with everolimus, or when deciding to initiate therapy with such a drug in a patient who has recently received a XIENCE PRO™S stent.

Everolimus, when prescribed as an oral medication, may interact with the following drugs or food, including the following, but not limited to:

- CYP3A4 / P-glycoprotein isozyme inhibitors
 - o Antifungal agents (e.g., fluconazole, ketoconazole, itraconazole, posaconazole, voriconazole)
 - o Macrolide antibiotics (e.g., erythromycin, clarithromycin, telithromycin)
 - o Calcium channel blockers (e.g., verapamil, nifedipine, diltiazem)
 - o Protease inhibitors (e.g., ritonavir, atazanavir, saquinavir, darunavir, indinavir, nelfinavir, amprenavir, fosamprenavir)
 - o Other (e.g., cyclosporine, nefazodone, cisapride, metoclopramide, bromocriptine, cimetidine, danazol, sildenafil, terfenadine, astemizole, grapefruit / grapefruit juice, digoxin)
- CYP3A4 / P-glycoprotein isozyme inducers
 - o Antibiotics (e.g., rifampin, rifabutin, ciprofloxacin, ofloxacin)
 - o Anticonvulsants (e.g., carbamazepine, phenobarbital, phenytoin)
 - o Non-nucleoside reverse transcriptase inhibitors (e.g., efavirenz, nevirapine)
 - o Glucocorticoids (e.g., dexamethasone, prednisone, prednisolone)
 - o HMGCoA reductase inhibitors (e.g., simvastatin, lovastatin)
 - o Other (e.g., St. John's Wort)

For more detailed drug interaction information, please reference the most recent everolimus drug label¹.

10.7.2 Immune Suppression Potential

Everolimus, the XIENCE PRO™S stent active ingredient, is an immunosuppressive agent. Immune suppression was not observed in the SPIRIT and XIENCE family of clinical trials. However, for patients who receive several XIENCE PRO™S stents simultaneously, it may be possible for everolimus systemic concentrations to approach immunosuppressive levels temporarily, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. Therefore, consideration should be given to patients taking other immunosuppressive agents or who are at risk for immune suppression.

¹ Certican[®] UK SmPC, Afinitor[®] UK SmPC, Votubia[®] UK SmPC, Afinitor[®] US label, and Zortress[®] US label. Refer to www.MHRA.gov.uk, www.ema.europa.eu, and www.fda.gov for the most recent versions of these SmPC/labels.

10.7.3 Lipid Elevation Potential

Oral everolimus use in renal transplant and advanced renal cell carcinoma patients was associated with increased serum cholesterol and triglyceride levels, which in some cases required treatment. The effect was seen with both low and high dose prolonged oral therapy in a dose related manner. When used according to the indications for use, exposure to systemic everolimus concentrations from the XIENCE PRO[®] stent is expected to be significantly lower than concentration exposure usually obtained in transplant patients. Increased serum cholesterol and triglyceride levels were not observed in the SPIRIT and XIENCE family of clinical trials. Oral administration of everolimus in combination with cyclosporine has been associated with increased serum cholesterol and triglyceride levels.

10.7.4 Pregnancy

Pregnancy Category C: There are no adequate everolimus or XIENCE PRO[®] stent-related studies in pregnant women.

Effects on the developing foetus have not been studied¹. Effects of a similar stent, XIENCE V, on pre-natal and post-natal rat development were no different than the controls (see Section 10.8 Carcinogenicity, Genotoxicity, and Reproductive Toxicity). When administered at oral doses of 0.1 mg/kg or above to animals, everolimus has shown reproductive toxicity effects including embryotoxicity and fetotoxicity. Effective contraception is recommended to be initiated before implanting a XIENCE PRO[®] stent and continued for one year post-implantation. While there is no contraindication, the risks and reproductive effects are unknown at this time¹.

10.7.5 Lactation

It is unknown whether everolimus is distributed in human milk. Everolimus pharmacokinetic and safety profiles have not been determined in infants. Consequently, mothers should be advised of potential serious adverse reactions to everolimus in nursing infants. Prior to XIENCE PRO[®] stent implantation, decisions should be made regarding whether to discontinue nursing or conduct an alternate procedure.

10.8 Carcinogenicity, Genotoxicity, and Reproductive Toxicity

The carcinogenicity, genotoxicity, and reproductive toxicity of the XIENCE PRO[®] stent have not been evaluated; however, long-term carcinogenicity and teratology studies were performed with the similar XIENCE V stent. The test results from the XIENCE V stent, as described below, are applicable to the stent of the XIENCE PRO[®] Everolimus Eluting Coronary Stent System, due to similar stent design, delivery system materials, and identical stent coating technology.

A 26-week carcinogenicity study was conducted to evaluate the carcinogenic potential of XIENCE V stents following subcutaneous implantation in transgenic mice. During the course of the study, there were no abnormal clinical observations that suggested a carcinogenic effect of the test group (XIENCE V stent). The test group did not demonstrate an increased incidence of neoplastic lesions when compared to the negative control group. However, the positive control and the experimental positive control groups demonstrated notable increases in the incidence of neoplastic lesions compared to either the test or the negative control group. Based on the results of this study, the XIENCE V stent does not appear to be carcinogenic when implanted in transgenic mice for 26 weeks.

Genotoxicity studies were conducted on the XIENCE V stent in mammalian cells and bacteria. These studies included gene mutations in bacteria (Ames Test), gene mutations in mammalian cells (chromosomal aberration), test for clastogenicity in mammalian cells, and mammalian erythrocyte micronucleus test. Based on the results of these studies, the XIENCE V stent is not genotoxic.

In addition, a reproductive toxicity (teratology) study was conducted in female Sprague-Dawley rats. The XIENCE V stent did not affect the fertility or reproductive capability of female Sprague-Dawley rats. There was no statistical difference between the test article (XIENCE V stent) and the control system in terms of any of the evaluated parameters. The test article had no effect on litter size and caused no increase of *in utero* mortality. Additionally, the XIENCE V stent did not cause any reproductive toxicity in the offspring in this study.

10.9 Restricted Substances

One or more components of this device contains the following substance defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys and stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

10.10 Magnetic Resonance Imaging (MRI) Safety Information

Nonclinical testing has demonstrated that the XIENCE PRO[®] stent, in single and in overlapped configurations up to 71 mm in length, is MR Conditional. A patient with this device can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 or 3 Tesla
- Maximum gradient slew rate of 200 T/m/s per axis
- Maximum spatial gradient field of 3000 gauss/cm (30 T/m)
- Maximum MR System reported whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the XIENCE PRO[®] stent is expected to produce a maximum temperature rise of less than 4.5°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 6 mm from the XIENCE PRO[®] stent when imaged with a gradient echo or spin echo pulse sequence and a 3 Tesla MRI system.

The XIENCE PRO[®] stent should not migrate in this MRI environment. Nonclinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of the XIENCE PRO[®] stent.

11.0 POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with PCI treatment procedures and the use of a stent in native coronary arteries include, but are not limited to, the following:

- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device materials (cobalt, chromium, nickel, tungsten, methacrylic polymer and fluoropolymer), and drug reactions to everolimus, anticoagulation, or antiplatelet drugs
- Vascular access complications which may require transfusion or vessel repair, including:
 - o Catheter site reactions

- o Bleeding (ecchymosis, oozing, hematoma, hemorrhage, retroperitoneal hemorrhage)
- o Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
- o Embolism (air, tissue, plaque, thrombotic material or device)
- o Peripheral nerve injury
- o Peripheral ischemia
- Coronary artery complications which may require additional intervention, including:
 - o Total occlusion or abrupt closure
 - o Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - o Tissue prolapse / plaque shift
 - o Embolism (air, tissue, plaque, thrombotic material, or device)
 - o Coronary or stent thrombosis (acute, subacute, late, very late)
 - o Stenosis or restenosis
- Pericardial complications which may require additional intervention, including:
 - o Cardiac tamponade
 - o Pericardial effusion
 - o Pericarditis
- Cardiac arrhythmias (including aspecific, conduction disorders, atrial and ventricular arrhythmias)
- Cardiac ischemic conditions, including:
 - o Myocardial ischemia
 - o Myocardial infarction (including acute)
 - o Coronary artery spasm
 - o Unstable or stable angina pectoris
- Stroke / cerebrovascular accident (CVA) and transient ischemic attack (TIA)
- System organ failures:
 - o Cardio-respiratory arrest
 - o Cardiac failure
 - o Cardiopulmonary failure (including pulmonary edema)
 - o Renal insufficiency / failure
 - o Shock
- Bleeding
- Blood cell disorders (including heparin induced thrombocytopenia [HIT])
- Hypotension / hypertension
- Infection
- Nausea and vomiting
- Palpitations, dizziness, and syncope
- Chest pain
- Fever
- Pain
- Death

Adverse events associated with daily oral administration of everolimus in doses varying from 1.5 mg to 10 mg daily can be found in the Summary of Product Characteristics (SmPC) and labels for the drug¹. The risks described below include the anticipated adverse events relevant for the cardiac population referenced in the contraindications, warnings, and precaution sections of the everolimus labels / SmPCs and / or observed at incidences $\geq 10\%$ in clinical trials with oral everolimus for different indications. Please refer to the drug SmPCs and labels¹ for more detailed information and less frequent adverse events:

- Abdominal pain
- Anemia
- Angioedema (increased risk with concomitant angiotensin-converting enzyme [ACE] inhibitor use)
- Arterial thrombotic events
- Bleeding and coagulopathy (including hemolytic uremic syndrome [HUS], thrombotic thrombocytopenic purpura [TTP], and thrombotic microangiopathy – increased risk with concomitant cyclosporine use)
- Constipation
- Cough
- Diabetes mellitus
- Diarrhea
- Dyspnea
- Embryo-fetal toxicity
- Erythema
- Erythroderma
- Headache
- Hepatic artery thrombosis (HAT)
- Hepatic disorders (including hepatitis and jaundice)
- Hypersensitivity to everolimus active substance, or to other rapamycin derivatives
- Hypertension
- Infection (bacterial, fungal, viral or protozoan infections, including infections with opportunistic pathogens). Polyoma virus-associated nephropathy (PVAN) JC virus associated progressive multiple leukoencephalopathy (PML), fatal infections and sepsis have been reported in patients treated with oral everolimus
- Kidney arterial and venous thrombosis
- Laboratory test alterations (elevations of serum creatinine, proteinuria, hypokalemia, hyperkalemia; hyperglycemia, dyslipidemia including hypercholesterolemia and hypertriglyceridemia; abnormal liver function tests; decreases in hemoglobin, lymphocytes, neutrophils, and platelets)
- Lymphoma and skin cancer
- Male infertility
- Menstrual irregularities
- Nausea
- Nephrotoxicity (in combination with cyclosporine)
- Non-infectious pneumonitis (including interstitial lung disease)
- Oral ulcerations
- Pain
- Pancreatitis
- Pericardial effusion
- Peripheral edema

- Pleural effusion
- Pneumonia
- Pyrexia
- Rash
- Renal failure
- Upper respiratory tract infection
- Urinary tract infection
- Venous thromboembolism
- Vomiting
- Wound healing complications (including wound infections and lymphocele)

12.0 REPORTING OF ADVERSE EVENTS

The user and / or patient should report any serious incident that has occurred in relation to XIENCE PRO[™]S Everolimus Eluting Coronary Stent System to the manufacturer and the competent authority of the European Member State in which the user and / or patient is established.

13.0 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of the safety and clinical performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI provided in Table 1-1. This is the SSCP location after the launch of the European Database on Medical Devices (EUDAMED).

14.0 PATIENT SELECTION AND TREATMENT

14.1 Individualization of Treatment

The risks and benefits described above should be considered for each patient before using a device from the family of XIENCE[™] EECSS. Patient selection factors to be assessed should include judgment regarding risk of antiplatelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Antiplatelet drugs should be used in combination with the XIENCE EECSS, per the guidelines from the American College of Cardiology, American Heart Association, and Society for Cardiovascular Angiography and Interventions (ACC/AHA/SCAI) and European Society of Cardiology (ESC) guidelines. Physicians should use the information from the large body of clinical evidence for XIENCE stents, coupled with the current literature on drug-eluting stents, current guidelines and the specific needs of individual patients, to determine the specific antiplatelet / anticoagulation regimen to be used for their patients in general practice.

Current guidelines for dual antiplatelet therapy (DAPT) discontinuation should be followed and are recommended. The decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition. In case an unanticipated interruption or discontinuation of DAPT is required any time after one month following XIENCE coronary stent implantation, two-year data from the XIENCE coronary clinical trials show low stent thrombosis rates and no observed increased risk for stent thrombosis. Additionally, the XIENCE Family of Stents has demonstrated low stent thrombosis rates from 1–6 months in HBR patients with 1-month (as shown at 28 days) DAPT duration post-PCI.

It is very important that the patient comply with the post-procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction, or death. Prior to PCI, if a surgical or dental procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventionalist and patient should carefully consider whether a drug-eluting stent and its associated recommended antiplatelet therapy is the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature discontinuation of antiplatelet therapy.

Patients who require premature discontinuation of antiplatelet therapy secondary to significant active bleeding, should be monitored carefully for cardiac events, and once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

15.0 CLINICIAN USE INFORMATION

15.1 Inspection Prior to Use

1. Carefully inspect the sterile package before opening and check for damage to the sterile barrier. Do not use if the integrity of the sterile package has been compromised.
 2. Do not use after the "Use by" date.
 3. Tear open the foil pouch and remove the inner pouch.
- Note:** The outside of the inner pouch is NOT sterile. Open the inner pouch and pass or drop the product into the sterile field using an aseptic technique.
4. Prior to using the XIENCE PRO[™]S EECSS, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent does not extend beyond the radiopaque balloon markers. Do not use if any defects are noted. However, do not manipulate, touch, or handle the stent, which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.

Note: At any time during use of the XIENCE PRO[®] EECSS, if the stainless steel proximal shaft has been bent or kinked, do not continue to use the catheter.

15.2 Materials Required

- Appropriate arterial sheath
- Appropriate guiding catheter(s)
- 2–3 syringes (10–20 cc)
- 1,000 u / 500 cc heparinized normal saline (HepNS)
- Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
- 0.014 inch (0.36 mm) x 175 cm (minimum length) guide wire
- Torque device
- Guide wire introducer
- Contrast diluted 1:1 with heparinized normal saline
- Inflation device
- Appropriate size pre-dilatation angioplasty balloon
- Appropriate size post-dilatation noncompliant angioplasty balloon
- Three-way stopcock
- Appropriate anticoagulation and antiplatelet drugs

¹ Certican[®] UK SmPC, Afinitor[®] UK SmPC, Votubia[®] UK SmPC, Afinitor[®] US label, and Zortress[®] US label. Refer to www.MHRA.gov.uk, www.ema.europa.eu, and www.fda.gov for the most recent versions of these SmPC/labels.

15.3 Preparation

15.3.1 Packaging Removal

Note: The foil pouch is not a sterile barrier. The inner header bag (pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

1. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. When using a rapid exchange (RX) system, do not bend or kink the hypotube during removal.
2. Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product, and replace with another. Follow the product returns procedure for the unused device.

15.3.2 Guide Wire Lumen Flush

1. Flush the guide wire lumen with HepNS until fluid exits the guide wire exit notch.
- Note:** Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent on the balloon.

15.3.3 Delivery System Preparation

1. Prepare an inflation device / syringe with diluted contrast medium.
2. Attach an inflation device / syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device / syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

Note: While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

Note: If air is seen in the shaft, repeat **Section 15.3.3 Delivery System Preparation**, steps 3 through 5, to prevent uneven stent expansion.

15.4 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
2. The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. **If pre-dilatation is performed**, limit the length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the XIENCE PRO^S stent.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel.

Note: If choosing between two stent diameters for tight lesions, choose the smaller diameter stent and inflate. See the product label for compliance information.

4. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
5. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion.
6. Carefully advance the delivery system into the guiding catheter and over the guide wire to the target lesion. When using a rapid exchange (RX) system, be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent system into the coronary artery.

Note: If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and the use of excessive force may result in stent damage or dislodgement. Maintain guide wire placement across the lesion and remove the delivery system and guiding catheter as a single unit.

7. Advance the delivery system over the guide wire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm stent position. If the position of the stent is not optimal, it should be carefully repositioned or removed (see **Section 15.6 Removal Procedure**). The balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion.

Note: If removal of a stent system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guiding catheter, the stent delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

8. Tighten the rotating hemostatic valve. The stent is now ready to be deployed.

15.5 Deployment Procedure

CAUTION: Refer to Table 1-3: *In vitro* Device Specifications, for stent diameter, nominal pressure, and RBP.

1. Prior to deployment, reconfirm the correct position of the stent relative to the target lesion using the radiopaque balloon markers.
2. Deploy the stent slowly by pressurizing the delivery system in 2 atm increments, every 5 seconds, until the stent is completely expanded. Fully expand the stent by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter ratio of about 1.1 times the reference vessel diameter.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel and expand the stent to nominal pressure at minimum. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the stent to the artery wall.
4. Maintain pressure for 30 seconds for full expansion of the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter as compared to the proximal and distal native coronary artery diameters (reference vessel diameters). Optimal stent expansion and proper apposition requires that the stent be in full contact with the arterial wall.

Note: See **Section 15.6 Removal Procedure** for instruction on withdrawal of stent delivery system.

5. If necessary, the delivery system can be repressurized or further pressurized to ensure complete apposition of the stent to the artery wall.

Note: Do not exceed the labeled RBP of 16 atm (1621 kPa).

6. Fully cover the entire lesion and balloon-treated area (including dissections) with the XIENCE PRO^S stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
7. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to guiding catheter position.

Note: See **Section 15.6 Removal Procedure** for instruction on withdrawal of stent delivery system.

8. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography, intravascular ultrasound (IVUS), or optical coherence tomography (OCT).
9. If the deployed stent size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent. If the initial angiographic appearance is suboptimal, the stent may be further expanded using a low profile, high pressure, noncompliant balloon dilation catheter. If this is required, the stented segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the stent geometry. Deployed stents should not be left underdilated.

CAUTION: Do not dilate the stent beyond the following limits.

- | Nominal Stent Diameter | Dilation Limit |
|------------------------|----------------|
| 2.0 mm – 3.25 mm | 3.75 mm |
| 3.5 mm – 4.0 mm | 5.5 mm |
10. If more than one XIENCE PRO^S stent is needed to cover the lesion and balloon-treated area, it is suggested that to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE PRO^S stent should be positioned inside the deployed stent prior to expansion.
 11. Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

15.6 Removal Procedure

Withdrawal of the stent delivery catheter from the deployed stent:

1. Deflate the balloon by pulling negative pressure on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10 to 15 seconds longer.
2. Position the inflation device on "negative" or "neutral" pressure.
3. Stabilize the guiding catheter position just outside coronary ostium and anchor in place. Maintain the guide wire placement across stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Tighten the rotating hemostatic valve.

If during withdrawal of the stent delivery catheter resistance is encountered, use the following steps to improve balloon rewrap.

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

Post-stent delivery system withdrawal; stent deployment confirmation:

1. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography, intravascular ultrasound (IVUS), or optical coherence tomography (OCT).
2. If more than one XIENCE PRO^S stent is needed to cover the lesion and balloon-treated area, it is suggested that to avoid the potential for gap restenosis, the stents be adequately overlapped.
3. To ensure that there are no gaps between stents, the balloon marker bands of the second XIENCE PRO^S stent should be positioned inside the deployed stent prior to expansion.
4. Reconfirm stent position and angiographic results to assess stented area. Repeat inflations until optimal stent deployment is achieved. If post-dilation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Ensure that the stent wall is in contact with the artery wall.

15.7 Post-Deployment Dilation of Stent Segments

1. All efforts should be taken to ensure that the stent is not underdiluted.
2. If the deployed stent size is still inadequate with respect to the vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and noncompliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guide wire to avoid diluting the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the limits indicated in **Section 15.5 Deployment Procedure**.

15.8 Patient Implant Card

The physician or health care center shall complete the following steps on the patient implant card:

- Fill out patient identification (front of the card), health care center or doctor information, address and contact number, and the implant date on the implant card (back of the card).
- Apply peel-away stickers from the product label to the card.

The patient shall be provided with the **entire** patient implant card before being discharged.

16.0 SPIRIT AND XIENCE FAMILY OF CLINICAL TRIALS

The XIENCE PRO[™]S EECSS is based on the predicate devices XIENCE V[™] EECSS, XIENCE PRIME[™] EECSS, XIENCE Xpedition[™] EECSS, and XIENCE Alpine[™] EECSS.

The XIENCE PRO^S EECSS uses a similar stent platform, identical drug coating formulation, identical drug primer, similar nominal total drug content, and the identical stent contacting balloon materials as the XIENCE PRIME EECSS, XIENCE Xpedition EECSS, and XIENCE Alpine EECSS.

The XIENCE PRO^S EECSS differs from the XIENCE Alpine EECSS in the stent delivery system. The XIENCE PRO^S Stent Delivery System utilizes the same principle of operation and materials as other Abbott RX stent systems and coronary dilatation catheters.

Compared to the XIENCE V EECSS and to the XIENCE PRIME EECSS, the XIENCE PRO^S EECSS has the same stent security specification, same stent placement on the balloon between the balloon markers, a similar tip entry profile, and a similar taper length for XIENCE PRO^S stent sizes up to 28 mm in length. Based on the similar nature of the XIENCE PRO^S stent to the XIENCE PRIME and XIENCE V stents, the performance of the XIENCE PRO^S EECSS can be predicted to be similar to the performance of XIENCE V and XIENCE PRIME stents. Therefore, clinical trial data for XIENCE V and XIENCE PRIME stents are summarized in this section.

16.1 Pre-Market Clinical Trials

Principal safety and effectiveness of the XIENCE V stent have been established from a series of pre-market clinical trials. SPIRIT III RCT was the pivotal randomized clinical trial (RCT) that established the noninferiority of the XIENCE V stent to the TAXUS[®] Express[®] stent (TAXUS stent). The SPIRIT IV trial is a prospective, randomized, active-controlled, single-blinded, multicenter evaluation of the XIENCE V stent compared to the TAXUS Express stent (TAXUS stent) in the treatment of up to three *de novo* lesions ≤ 28 mm in length in native coronary arteries with RVD ≥ 2.5 mm to ≤ 4.25 mm. The SPIRIT IV trial established the superiority of the XIENCE V stent compared to the TAXUS Express stent. The SPIRIT Small Vessel (SV) Registry is a prospective, single-arm, open-label, US multicenter registry study that established the safety and effectiveness of the XIENCE V 2.25 mm diameter stent. SPIRIT PRIME is a prospective, open-label, multicenter, nonrandomized clinical trial with two study arms using the core size XIENCE PRIME and XIENCE PRIME LL stent systems that established the safety and effectiveness of the XIENCE PRIME and XIENCE PRIME LL stents. Tables 16.1-1 through 16.1-4 present the trial designs, angiographic results (for studies that required angiographic follow-up), and principal clinical outcomes at 1-year and from latest follow-up, respectively.

Table 16.1-1: SPIRIT Family of XIENCE V Clinical Trial Designs (Pre-Market)					
	SPIRIT III RCT	SPIRIT IV	SPIRIT Small Vessel Registry	SPIRIT PRIME Clinical Trial	
				Core Size Registry	Long Lesion Registry
Study Type / Design	• Multicenter • Randomized • Single-blinded • Active-control	• Multicenter • Randomized • Single-blinded • Active-control	• Multicenter • Open-label • Single-arm	• Multicenter • Open-label • Single-arm	• Multicenter • Open-label • Single-arm
Number of Subjects Enrolled	Total: 1,002 XIENCE V: 668 TAXUS Express Control: 334	Total: 3,690 XIENCE V: 2,460 TAXUS Express Control: 1,230 ¹	Total: 150 2.25 mm XIENCE V	Total: 400 XIENCE PRIME	Total: 100 XIENCE PRIME
Treatment	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to three <i>de novo</i> lesions, maximum of two lesions per epicardial vessel	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels	Up to two <i>de novo</i> lesions in different epicardial vessels
Lesion Size	RVD: ≥ 2.5 ≤ 3.75 mm Length: ≤ 28 mm	RVD: ≥ 2.5 ≤ 4.25 mm ² Length: ≤ 28 mm	RVD: ≥ 2.25 ≤ 2.50 mm Length: ≤ 28 mm	RVD: ≥ 2.25 ≤ 4.25 mm Length: ≤ 22 mm	XIENCE PRIME CS: RVD: ≥ 2.25 ≤ 4.25 mm Length: ≤ 22 mm XIENCE PRIME LL: RVD: ≥ 2.5 ≤ 4.25 mm Length: > 22 mm and ≤ 32 mm
Primary Endpoint	In-segment late loss at 240 days	Ischemia-driven target lesion failure at 1 year (composite of cardiac death, target vessel MI or ischemia-driven TLR)	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year	TLF (target lesion failure) at 1 year
Co-Primary Endpoint	TVF at 270 days	None	None	None	None
Clinical Follow-up	30, 180, 240, 270 days, 1 to 5 years	30, 180, 270 days, 1 to 3 years	30, 240 days, 1 to 3 years	30, 180 days, 1 to 3 years	30, 180 days, 1 to 3 years
Angiographic Follow-up	240 days (N = 564)	None	240 days (N = 69)	None	None

¹ In the TAXUS stent arm, there was 1 subject who received 1 TAXUS¹ Liberté² stent.
² RVD ≥ 2.5 mm up to ≤ 3.75 mm and stent sizes up to 3.5 mm until 4.0 mm TAXUS is commercially available.

Angiographic Results	SPIRIT III RCT 240 Days		SPIRIT Small Vessel 240 Days
	XIENCE V (N = 376) (M = 427)	TAXUS (N = 188) (M = 220)	2.25 mm XIENCE V (N = 69) (M = 69)
In-Stent Late Loss (mm)	0.16 ± 0.41 (342)	0.30 ± 0.53 (158)	0.20 ± 0.40 (52)
In-Segment Late Loss (mm)	0.14 ± 0.39 (343)	0.26 ± 0.46 (158)	0.16 ± 0.41 (52)
In-Stent Binary Restenosis	2.3% (8/343)	5.7% (9/158)	3.8% (2/52)
In-Segment Binary Restenosis	4.7% (16/344)	8.9% (14/158)	9.6% (5/52)

Notes:
– Data are mean (mm) ± SD or % (n/N).
– N is total number of patients. M is total number of lesions.
– SPIRIT III and SV 240-day include follow-up window (240 + 28 days).

Table 16.1-3: SPIRIT Family of Clinical Trials Principal 1-Year Clinical Outcomes (Pre-Market)							
	SPIRIT IV		SPIRIT III RCT		SPIRIT Small Vessel	SPIRIT PRIME Clinical Trial	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	4.0% (97/2416)	6.8% (81/1195)	5.3% (35/655)	9.7% (31/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
TVF	5.5% (134/2416)	7.7% (92/1195)	8.5% (56/655)	11.6% (37/319)	11.0% (15/136)	N/A	N/A
MACE	4.1% (98/2416)	6.9% (82/1195)	6.0% (39/655)	10.3% (33/319)	8.1% (11/136)	4.5% (18/399)	7.7% (8/104)
All Death	1.0% (25/2416)	1.3% (15/1195)	1.2% (8/657)	1.3% (4/320)	1.5% (2/136)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	0.4% (10/2416)	0.4% (5/1195)	0.8% (5/657)	0.9% (3/320)	1.5% (2/136)	0.3% (1/399)	0.0% (0/104)
MI	1.9% (45/2416)	3.1% (37/1195)	2.7% (18/655)	4.1% (13/319)	1.5% (2/136)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	2.2% (54/2416)	3.3% (39/1195)	3.4% (22/655)	4.7% (15/319)	2.9% (4/136)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	2.3% (56/2416)	4.6% (55/1195)	3.4% (22/655)	5.6% (18/319)	5.1% (7/136)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	2.2% (54/2416)	2.4% (29/1195)	3.2% (21/655)	4.7% (15/319)	5.9% (8/136)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.29% (7/2391)	1.10% (13/1181)	0.9% (6/650)	0.6% (2/316)	1.5% (2/136)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.3% (6/2385)	0.8% (10/1183)	0.8% (5/650)	0.3% (1/317)	0.7% (1/138)	0.5% (2/399)	0.0% (0/104)

Notes:
– All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
– 1-year includes the follow-up window (365 + 28 days) for all trials.
– TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically indicated TLR definition rather than ischemia-driven TLR.
– TVF includes cardiac death, MI, ischemia-driven TLR and TVR, nontarget lesion. SPIRIT SV and PRIME used clinically indicated TLR and TVR definition rather than ischemia-driven TLR and TVR definition, which was used for SPIRIT II, SPIRIT III, and SPIRIT IV.
– MACE includes cardiac death, MI and ischemia-driven TLR.

Table 16.1-4: SPIRIT Family of Clinical Trials
Principal Clinical Outcomes from Latest Follow-up (Pre-Market)

	SPIRIT IV 3 Years		SPIRIT III RCT 5 Years		SPIRIT Small Vessel 2 Years	SPIRIT PRIME Clinical Trial 1 Year	
	XIENCE V (N = 2458)	TAXUS (N = 1229)	XIENCE V (N = 669)	TAXUS (N = 333)	2.25 mm XIENCE V (N = 144)	Core Size Registry (N = 401)	Long Lesion Registry (N = 104)
TLF	9.5% (223/2348)	11.9% (138/1158)	13.4% (81/605)	20.6% (59/286)	8.3% (11/133)	4.5% (18/399)	7.7% (8/104)
TVF	13.3% (312/2348)	14.5% (168/1158)	20.3% (123/605)	26.6% (76/286)	12.0% (16/133)	N/A	N/A
MACE	9.8% (231/2348)	12.3% (142/1158)	14.4% (87/605)	22.0% (63/286)	8.3% (11/133)	4.5% (18/399)	7.7% (8/104)
All Death	3.4% (81/2348)	5.2% (60/1158)	6.0% (37/621)	10.3% (31/300)	1.5% (2/133)	0.8% (3/399)	1.0% (1/104)
Cardiac Death	1.4% (34/2348)	1.9% (22/1158)	2.7% (17/621)	4.3% (13/300)	1.5% (2/133)	0.3% (1/399)	0.0% (0/104)
MI	3.1% (73/2348)	4.7% (55/1158)	4.6% (28/605)	7.0% (20/286)	1.5% (2/133)	1.8% (7/399)	4.8% (5/104)
Cardiac Death or MI	4.5% (105/2348)	6.0% (70/1158)	7.1% (43/605)	11.2% (32/286)	3.0% (4/133)	2.0% (8/399)	4.8% (5/104)
Ischemia-Driven TLR	6.3% (148/2348)	7.9% (92/1158)	8.9% (54/605)	12.9% (37/286)	5.3% (7/133)	2.5% (10/399)	2.9% (3/104)
Ischemia-Driven TVR, Non TL	5.6% (132/2348)	5.4% (63/1158)	8.8% (53/605)	11.9% (34/286)	6.8% (9/133)	2.8% (11/399)	2.9% (3/104)
Stent Thrombosis							
ARC (Definite / Probable)	0.62% (14/2263)	1.73% (19/1098)	1.5% (9/582)	1.9% (5/268)	1.5% (2/132)	0.5% (2/399)	0.0% (0/104)
ARC (Definite)	0.49% (11/2263)	1.28% (14/1098)	1.2% (7/582)	0.7% (2/268)	0.8% (1/132)	0.5% (2/399)	0.0% (0/104)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- Data includes the follow-up window of + 28 days for all trials.
- TLF includes cardiac death, MI attributed to target vessel and ischemia-driven TLR. SPIRIT SV and PRIME used clinically indicated TLR definition rather than ischemia-driven TLR.
- TVF includes cardiac death, MI, ischemia-driven TLR and TVR, nontarget lesion.
- MACE includes cardiac death, MI, and ischemia-driven TLR.

16.2 Post-Market Clinical Trials

The XIENCE V USA study is a prospective, multicenter, FDA-mandated post-market study to evaluate the continued safety and effectiveness of the XIENCE V EECSS in real-world settings after it was commercialized in the U.S., and also to support the FDA DAPT initiative. The objective of the SPIRIT V Single Arm Study (SAS) was to continue the assessment of the performance of XIENCE V EECSS in the treatment of patients with *de novo* coronary artery lesions. XIENCE V India is a prospective, open-label, multicenter, observational, single-arm registry to evaluate XIENCE V EECSS's continued safety and effectiveness during commercial use in real-world settings in India. Tables 16.2-1 through 16.2-2 present clinical trial designs, principal clinical outcomes at 1 year and from latest follow-up, respectively.

The results from these post-marketing clinical trials demonstrate safety and effectiveness of the XIENCE V EECSS in real-world settings. In addition, the XIENCE V EECSS improved patient-reported outcomes (including better quality of life, reduced angina frequency, improved angina stability, and reduced physical limitation) at 6 months and the improvements were sustained for 1 year in coronary artery disease patients.

Table 16.2-1: SPIRIT and XIENCE V Family Clinical Trial Designs (Post-Market)

	XIENCE V USA Phase I Cohort	XIENCE V USA Long-Term Follow-up Cohort	XIENCE V USA AV-DAPT Cohort	SPIRIT V (SAS)	XIENCE V India
Study Type / Design	• Multicenter • Prospective • Single-arm	• Multicenter • Prospective • Single-arm	• Multicenter • Randomized • Double-blinded • Placebo control	• Multicenter • Prospective • Single-arm	• Multicenter • Prospective • Single-arm
Number of Subjects Enrolled	8040	4663	868	2663	977
Treatment	Per site standard care	Per site standard care	Patients were randomized to receive either thienopyridine or placebo treatment for additional 18 months along with aspirin	Maximum of one <i>de novo</i> , native target lesion per major epicardial vessel or side branch (no prior stent implant, no prior brachytherapy), maximum of 4 planned EES	Per site standard care
Lesion Size	No angiographic restrictions			RVD = $\geq$ 2.25 $\leq$ 4.0 mm Length $\leq$ 28 mm by visual estimation	No angiographic restrictions
Primary Endpoint	ARC definite and probable stent thrombosis up to 1 year	ARC definite and probable stent thrombosis from year 1 to 5	MACE (composite of all death, MI, and stroke) 12–33 months	Composite rate of all death, MI, TVR at 30-day	ARC ST (Definite / probable) 1 year and yearly after through 3 years
Co-Primary Endpoint	Cardiac death or any MI at 1 year	Cardiac death or any MI from year 1 to 5	ARC definite and probable ST 12–33 months	None	Cardiac death and any MI at 1 year
Clinical Follow-up	14, 30, 180 days, and 1 year	2, 3 and 4 years	15, 24, 30, and 33 months	30 days and 1 and 2 years	14, 30, 180 days and 1, 2, and 3 years
Angiographic Follow-up	None	None	None	None	None

Table 16.2-2: SPIRIT and XIENCE V Family of Trials Principal Clinical Outcomes (Post-Market)					
	XIENCE V USA Phase I 1 Year	SPIRIT V (SAS)		XIENCE V India	
		1 Year	2 Years	1 Year	2 Years
	XIENCE V (N = 8040)	XIENCE V (N = 2663)	XIENCE V (N = 2663)	XIENCE V (N = 990)	XIENCE V (N = 990)
TLF (ARC)	9.4% (707/7522)	5.25% (138/2627)	7.49% (192/2562)	2.4% (24/986)	3.4% (32/942)
TLF	6.8% (513/7505)	N/A	N/A	N/A	N/A
All Death, MI (ARC) and TVR	N/A	7.04% (185/2627)	10.34% (265/2562)	N/A	N/A
Cardiac Death or MI (ARC)	7.2% (545/7522)	4.23% (111/2627)	5.74% (147/2562)	1.9% (19/986)	2.8% (26/942)
Cardiac Death or MI	3.3% (249/7505)	N/A	N/A	N/A	N/A
TLR	4.6% (349/7522)	1.90% (50/2627)	3.04% (78/2562)	1.2% (12/986)	1.5% (14/942)
TVR, non TLR	2.3% (176/7522)	1.45% (38/2627)	2.26% (58/2562)	0.1% (1/986)	0.1% (1/942)
All Death	2.6% (194/7522)	1.71% (45/2627)	2.97% (76/2562)	0.9% (9/986)	1.7% (16/942)
Cardiac Death	1.4% (108/7522)	1.10% (29/2627)	1.87% (48/2562)	0.9% (9/986)	1.7% (16/942)
All MI (ARC)	6.3% (475/7522)	3.54% (93/2627)	4.45% (114/2562)	1.3% (13/986)	1.6% (15/942)
All MI	2.2% (162/7505)	N/A	N/A	N/A	N/A
Stent Thrombosis					
ARC (Definite / Probable)	0.81% (60/7380)	0.65% (17/2607)	0.79% (20/2523)	0.51% (5/986)	0.53% (5/939)
ARC (Definite)	0.54% (40/7380)	N/A	N/A	0.41% (4/986)	0.43% (4/939)

Notes:

- All counts presented in this table are subject counts. Subjects are counted only once for each event for each time period.
- XIENCE V USA 1-year data include 42-day window or through randomization date if occurred earlier than 407 days for the second enrollment phase. XIENCE V India 1-year data include 47-day window. SPIRIT V 1-year data include 28-day window.
- SPIRIT V 2-year data include 28-day window. XIENCE V India 2-year data include 47-day window.
- TLF (ARC) includes cardiac death, MI attributed to target vessel (per ARC definition), clinically indicated TLR. TLF includes cardiac death, MI attributed to target vessel (per protocol definition), clinically indicated TLR.

16.3 One-Month (as Short as 28 Days) DAPT Indication

The safety of the XIENCE stent in high bleeding risk (HBR) patients treated with 1-month (as short as 28 days) dual antiplatelet therapy (DAPT) post-PCI was established in the XIENCE 28 trial.

The XIENCE 28 trial is a prospective, single arm, multicenter, open label, non-randomized trial to evaluate the safety of 1-month (as short as 28 days) dual antiplatelet therapy (DAPT) in high bleeding risk (HBR) patients undergoing PCI with the XIENCE stent. Patients were eligible for inclusion in XIENCE 28 if they fulfilled clinical criteria indicative of HBR and, in the opinion of the referring physician, the risk of major bleeding with more than 1-month DAPT outweighed the antithrombotic benefit. To be considered as HBR, patients had to meet at least one of the following criteria: a) age ≥ 75 years; b) clinical indication for chronic (at least 6 months) anticoagulation therapy; c) history of major bleeding requiring medical attention within 12 months of the index procedure; d) history of stroke (either ischemic or hemorrhagic); e) renal insufficiency (creatinine ≥ 2.0 mg/dl) or failure (dialysis dependent); f) systemic conditions associated with an increased bleeding risk (e.g., hematological disorders, including a history of or current thrombocytopenia defined as a platelet count < 100,000/mm³, or any known coagulation disorder associated with increased bleeding risk); g) anemia with hemoglobin < 11 g/dl. The XIENCE stent was implanted per site standard of care. There was no angiographic restriction on lesion length, but the longest XIENCE stent allowed in the trial was 38 mm. Clinical follow-up is conducted at 1 month, 3 months, 6 months and 12 months.

The XIENCE 28 trial was comprised of 58 sites within United States and Canada (XIENCE 28 USA) and 52 sites outside of US (XIENCE 28 Global included sites in Europe and Asia). A total of 1,605 HBR patients underwent PCI with the XIENCE stent and were on DAPT (aspirin + P2Y12 receptor inhibitor) for the first month (as short as 28 days) post-index procedure. Patients were assessed as “1-month clear” and were eligible to discontinue P2Y12 at 1 month and continue on aspirin monotherapy through the end of study, if they were event-free (free from myocardial infarction, repeat coronary revascularization, stroke, or definite / probable stent thrombosis) and DAPT compliant (defined as no interruption > 7 consecutive days) during the first month post-procedure.

The objective of the trial was to show noninferiority of the primary endpoint of all death or all myocardial infarction (MI; modified ARC) from 1 to 6 months following XIENCE stent implantation in HBR patients treated with 1-month DAPT compared to a historical control, treated with up to 12-month DAPT, after propensity score (PS) stratification. The PS stratification method was used to provide a fair comparison between XIENCE 28 and the XIENCE V USA (XV USA) historical control. Patients from XIENCE 28 and XV USA were divided into five strata based on the PS calculated from their baseline characteristics. Patients with similar PS were assigned the same stratum. Event rates were then calculated for each stratum. The major secondary endpoint was to show superiority of 1-month DAPT for major bleeding (Bleeding Academic Research Consortium [BARC] type 2–5) from 1 to 6 months following XIENCE stent implantation in HBR patients, compared to the XV USA historical control after PS stratification.

Baseline Patient and Lesion Characteristics in XIENCE 28 1-Month Clear Population

Of the 1,392 1-month clear patients, 32.5% were female, and the mean age was 76 years. In terms of HBR criteria, 68.2% were 75 years or older, 47.4% had renal insufficiency, 44.4% were on oral anticoagulant, 10.4% had a history of stroke, 3.3% had a history of major bleeding, 14.5% had anemia, and 2.3% had thrombocytopenia. Risk factors having a high prevalence included hypertension (84.7%) and dyslipidemia (67.5%). In terms of cardiac status, 34.1% of patients presented with acute coronary syndrome (ACS) (non-ST segment elevation myocardial infarction [NSTEMI], or unstable angina) and 65.9% had chronic coronary syndrome. A majority of patients (80.3%) had a single lesion treated, and the mean total stent length per patient was 27.2 ± 14.4 mm.

Primary Endpoint Analysis

The primary endpoint was met demonstrating that the XIENCE stent is safe in HBR patients treated with 1-month (as short as 28 days) DAPT post-PCI. Based on the number of patients and observed rates in each stratum, the results show a noninferiority p-value of 0.0005, meeting the pre-specified significance level of 0.025. The 1–6 months PS stratified mean rate for all death / all MI in 1-month clear patients was 3.5% in XIENCE 28 vs 4.3% in XV USA historical control (Table 16.3-1). This demonstrated that the 1–6 months composite rate of all death / all MI for the 1-month DAPT regimen is noninferior to the recommended 12-month DAPT regimen in the historical control.

Major Secondary Endpoint Analysis

The PS stratification method was also used to compare the XIENCE 28 trial arm to historical control for the endpoint of BARC 2–5 bleeding. Based on the number of patients and observed rates in each stratum, the results show a superiority p-value of 0.1888, which did not meet the pre-specified significance level of 0.025. Although the significance level was not met, the 1–6 months PS stratified mean rate of BARC 2–5 bleeding, in 1-month clear patients, was numerically lower in XIENCE 28 as compared to the historical control (4.9% vs 5.9%; Table 16.3-1).

Table 16.3-1: Primary and Powered Secondary Endpoints Results – All Death / All MI and BARC 2–5 Bleeding from 1–6 Months in 1-Month Clear Patients

	XIENCE 28 1-Month Clear Patients (N = 1392)	XIENCE V USA 1-Month Clear Patients (N = 1411)	p-value
All death / All MI	3.5%	4.3%	0.0005*
BARC 2–5 bleeding	4.9%	5.9%	0.1888**

* Noninferiority p-value; noninferiority margin and one-sided significant level are 2.5% and 0.025, respectively.

** Superiority p-value; superiority significance level is 0.025.

Other Secondary Endpoints

Other secondary endpoints analyzed in XIENCE 28 included BARC 3–5 bleeding and stent thrombosis from 1–6 months in 1-month clear patients.

Although not pre-specified, the PS stratification methodology was used to provide a fair comparison of BARC 3–5 bleeding rates between XIENCE 28 and XV USA historical control. This analysis suggests better bleeding outcomes associated with 1-month DAPT, as compared to 12-months DAPT. The 1–6 months PS stratified mean rate of BARC 3–5 bleeding in 1-month clear patients was significantly lower in XIENCE 28 as compared to the historical control (2.2% vs. 4.5%; Table 16.3-2).

Table 16.3-2: PS Stratified Analysis of BARC 3–5 Bleeding from 1–6 Months in 1-Month Clear Patients

	XIENCE 28 1-Month Clear Patients (N = 1392)	XIENCE V USA 1-Month Clear Patients (N = 1411)	p-value*
BARC 3–5 bleeding (PS stratified mean rate)	2.2%	4.5%	0.0156

* Superiority p-value; superiority significance level is 0.025. Ad-hoc analysis.

The non-PS stratified definite / probable stent thrombosis rate from 1–6 months, in 1-month clear patients, was 0.3% in both XIENCE 28 and the historical control groups (Table 16.3-3).

Table 16.3-3: Stent Thrombosis from 1–6 Months in 1-Month Clear Patients (Non-PS Stratified)

	XIENCE 28 1-Month Clear Patients (N = 1392)	XIENCE V USA 1-Month Clear Patients (N = 1411)
ARC Definite / Probable stent thrombosis	0.3%	0.3%

17.0 HOW SUPPLIED

Sterile – This device is sterilized with ethylene oxide gas. Non-pyrogenic. Do not use if the package is open or damaged.

This single-use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Contents – One (1) XIENCE PRO™S Everolimus Eluting Coronary Stent System.

Storage – Store in a dry, dark, cool place. Protect from light. Do not remove from carton until ready for use. Store at 25°C (77°F); excursions permitted between 15-30°C (59-86°F).

Reference Abbott website for patent markings: www.abbott.com/patents

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

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© 2023 Abbott. All Rights Reserved.
Abbott Medical UK Limited