IMMULITE® 2000

Chemiluminescent Substrate Module

English

Intended Use: For use in the routine operation of the IMMULITE 2000 Systems, as described in the Operator's Manuals.

Catalog Number: **L2SUBM**

Summary and Explanation

Chemiluminescent Substrate is a dioxetane-based compound that is intended to be used in all IMMULITE assay reactions to generate light emission for qualitative or quantitative specific-analyte determinations.

Materials Supplied

Two bottles (L2SUBX), 205 mL each (approximately 1000 tests), of ready-to-use substrate, consisting of a phosphate ester of adamantyl dioxetane¹ in an AMP² buffer, with enhancer.³ Store at 2 – 8 °C before opening. Stable on board after opening at 15 – 28 °C for 30 days, or until expiration date.

Allow substrate to equilibrate at room temperature (about 20 minutes) prior to opening bottle and dispensing.

Before use, replace the cap on the bottle with the toggle cap (supplied).

¹ LUMIGEN® PPD: 4-methoxy-4-(3-phosphatophenyl)-spiro- (1,2-dioxetane-3,2'-adamantane). LUMIGEN® is a registered trademark of LUMIGEN, Inc., Southfield, MI. Patent numbers: 4,857,652, 4,983,779, 4,962,192, 5,386,017, and 5,707,559 (US); 254051B1 (Europe); 5-45590 (Japan).

² AMP: 2-amino-2-methyl-1-propanol.

³ Polyvinylbenzyltributylphosphonium chloride. Patent numbers: 5,393,469 and 5,439,617 (US); 2648423 (Japan).

Risk and Safety

Safety data sheets (MSDS/SDS) available on [siemens-healthineers.com](https://www.siemens-healthineers.com)

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Appearance: Precipitates or turbidity may indicate deterioration. Do *not* use if the normal appearance (clear and colorless) has been altered.

Light Sensitivity: Protect the Chemiluminescent Substrate from direct sunlight.

Avoid Contamination: Do not touch the dispensing tip on the toggle cap, introduce foreign objects into the bottle, or pool partial bottles of Substrate.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

For technical assistance, contact your National Distributor.

[siemens-healthineers.com](https://www.siemens-healthineers.com)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Shading indicates technical content that differs from the previous version.

IMMULITE® is a trademark of Siemens Healthineers.

All other trademarks and brands are the property of their respective owners.

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All rights reserved.