








## 1.5 Accessories

	REF.	DESCRIPTION	QTY
	204501	<b>Adjustable flange</b>	5
		<ul style="list-style-type: none"> <li>• To fix haemostatic catheters</li> <li>• For catheters in size Ch. 24</li> <li>• Sterile, single use</li> </ul>	
	546000	<b>Catheter stylet, made of metal</b>	1
		<ul style="list-style-type: none"> <li>• Thompson pattern (one straight and one curved end)</li> <li>• Diameter 2.5 mm</li> <li>• Approx. 55 cm long</li> <li>• <b>Non-sterile</b></li> </ul>	
	548000	<b>Catheter plug, made of blue PP</b>	100
		<ul style="list-style-type: none"> <li>• For Ch. 7–30</li> <li>• <b>Non-sterile</b></li> <li>• <b>Latex-free</b></li> </ul>	
	548500	<b>Catheter plug, made of blue PP</b>	50
		<ul style="list-style-type: none"> <li>• For Ch. 7–30</li> <li>• Individual packaging</li> <li>• Blunt-nosed</li> <li>• Sterile, single use</li> <li>• <b>Latex-free</b></li> </ul>	
	556000	<b>Silkospray</b>	10
		Silkospray prevents encrustation of rubber, latex and PVC devices and also prevents them from sticking to the mucosa. The propellant used contains no CFCs and complies with environmental regulations.	
	187501	<b>Perineal drape, made of plastic</b>	50
	187502	<ul style="list-style-type: none"> <li>• Without adhesive tape</li> <li>• With adhesive tape</li> </ul>	
	850550	<b>SilFlate, prefilled syringe<sup>1</sup></b>	10
		<ul style="list-style-type: none"> <li>• With sterile 10 % aqueous glycerine solution</li> <li>• 10 ml</li> <li>• Sterile, single use</li> </ul>	

<sup>1</sup> Attention: To inflate silicone balloons only!

When ordering, please always indicate Ref. and order size