

# Biohazard Bag

You are here: [Home](#) » [Products](#) » [General Lab Supplies](#) » [Safeguard](#) » [Biohazard Bag](#)

Share to:

Biohazard Bag 

Product Description

## Material LDPE

P/N	DIMENSION	APPROX THICKNESS	PACKING
0000-5001	14.2"x 16"	5Micrometer	100pcs/polybox , 10x100pcs/car
0000-5002	20"x 26"	5Micrometer	ton 100pcs/polybox , 10x100pcs/car ton

[Previous: Sharp Containers](#)

[Next: Biohazard Bag](#)

PRODUCT CATEGORY 

RELATED PRODUCTS

Email

Skype

poz.  
2.3



Round, PP Material



Sharp Containers



No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China



0086 513 82259348



Email

Skype

---

## Citotest Scientific Co.,Ltd

CITOTEST was established in 1996 in East China. The owner ambitiously founded CITOTEST for creating a cost-effective Manufacturer and Provider for high quality Microscope Slides and Microscope Cover Glass.

### QUICK LINK

- [Home](#)
- [Products](#)
- [About Us](#)
- [Support](#)
- [Certificates](#)
- [News](#)
- [Contact Us](#)

### PRODUCTS

- [Microscope Slides& Microscope Cover Glass](#)
- [Other Pathology Histology Products](#)
- [Microbiology/Cell Culture/Molecular Biology](#)
- [Hematology and Biochemistry](#)
- [Sampling/Transport/Storage](#)
- [Liquid Handling](#)
- [Centrifuge Tube and Bottle](#)
- [General Lab Supplies](#)
- [Mixing and Centrifugation](#)

### CONTACT US

 No. 339 Beihai West Road, Haimen, 226100 Jiangsu, P.R. China

 +86-025-86216803

 [info@citotest.com](mailto:info@citotest.com)

---

Copyright © Citotest Labware Manufacturing Co.,Ltd. All rights reserved.

---

Email

Skype

# VIRUS SPECIMEN COLLECTION AND TRANSPORT SYSTEM

## Nylon Flocked Swab

Superior specimen collection and release capability than traditional swabs, capable of more sample capacity such as Viruses, Bacteria, Chlamydia, Trichomonas vaginalis and Mycoplasma



Nasopharyngeal Swab




Oropharyngeal Swab

- Better nucleic acid adsorption capacity, improve the positive rate of virus detection effectively.



- Specimen could release to the transport medium immediately, reduce the waste of the specimen.



Flexible rods, according with ergonomic design.

Jiangsu Mole Bioscience Co.,Ltd

Zhejiang·Hangzhou  
Add: 4-5th Floor, Building No.5, Health Valley,  
No.1500 Wenyi xi Road, Yuhang District, Hangzhou,  
Zhejiang province, China  
Tel: 0086 571 87209308

Jiangsu·Taizhou  
Add: 6-7th Floor, G116 Building, No.805, Jiankang Avenue,  
Medical New&Hi-tech District, Taizhou, Jiangsu Province,  
China  
Tel: 0086 523 86086060



# VIRUS SPECIMEN COLLECTION AND TRANSPORT SYSTEM

## Sample Preservation Solution-Inactivated Medium

Guarantee the viability inactivation of specimens during transportation and testing, inhibit the activity of nuclease and maintain the stability of pathogen nucleic acid and the human genome.

## Disposable Collection Tubes of Virus Samples-Culture Medium

Hank's formula with antibiotics and phenol red, can be used for subsequent virus culture without affecting virus activity and infection ability.

## Disposable Collection Tubes of Virus Samples-PBS Medium

The main components are dibasic sodium phosphate, sodium chloride, etc. which can stabilize the virus.



### SAMPLING TUBE

- Made of medical polypropylene, resistant to low temperature (-70 °C)
- Concave screw cap and silicone O-ring ensure seal integrity of tubes, prevent the leakage of solution.

Used for transportation or transfer of pathogenic microorganism, such as virus, chlamydia mycoplasma, etc.

Three specifications of the medium meet different testing requirements: 1, 2, 3 ml.

PRODUCT DESCRIPTION	
<b>Packing size</b>	20 tests/box, 36 tests/box, 50 tests/box, 100 tests/box
<b>Tube</b>	2ml, 5ml, 10ml, 15ml
<b>Medium</b>	1ml, 2ml, 3ml, 5ml
<b>Swab</b>	Nasopharyngeal swab 8cm or 10cm breaking point; Oropharyngeal swab 3cm and/or 8cm breaking point.



# TECHNINIŲ DUOMENŲ LAPAS

<b>Produkto pavadinimas:</b>	Vienkartiniai mėgintuvėliai, mėginių, skirtų ištyrimui dėl virusų, paėmimui	<b>Kat. nr.:</b>	P045T0060100
		<b>Tūris:</b>	3 ml

## Pagrindiniai komponentai

- Natrio fosfato dvibazis dodekahidratas
- Kalio divandenilio fosfatas
- Magnio sulfato heptahidratas
- Magnio chlorido heksahidratas
- Gentamicino sulfatas
- Kalcio chloridas
- Kalio chloridas
- Natrio chloridas
- D(+)-Gliukozė
- Fenolio raudonasis

## Medžiagos

<b>Mėgintuvėlis</b>	<u>Kamštelis</u>	<u>Polipropilenas</u>	<b>Tamponėlis</b>	<u>Antgalis</u>	Suspausto nailono pluoštas
	Mėgintuvėlio korpusas	<u>Polipropilenas</u>		<u>Laikiklis</u>	Polipropilenas

## Matmenys

<u>Ilgis</u>	102.2±0.1 mm	<b>Nosiaryklės tamponas</b>	
<b>Mėgin-</b> <b>tuvėlis</b> <u>Diametras</u>	15.5±0.05 mm	<u>Ilgis</u>	149 +0/-2 mm
<u>Tūris</u>	10 mL	<b>Laikik-</b> <b>lis</b> <u>Lūžio vieta</u>	81 +0/-2 mm
<b>Kamš-</b> <b>telis</b> <u>Ilgis</u>	11.0±0.1 mm	<u>Diametras</u>	Φ1.05, Φ1.15, Φ2.6 mm
<u>Diametras</u>	20.3±0.05 mm	<b>Antga-</b> <b>lis</b> <u>Ilgis</u>	24 ±2 mm
<b>Bendras ilgis</b>	103.35±0.1 mm	<u>Diametras</u>	3 ±0.5 mm

\*\*\*Išsamūs matmenų duomenys pateikti techniniuose brėžiniuose

## Pakuotės duomenys

Rinkinio dydis: 36 vnt. rinkinyje	Dėžutės matmenys: 20*14*12 cm
Dėžės dydis: 18 rinkinių dėžėje (≈9.5 kg dėžė)	Dėžės matmenys: 42*40*45 cm



Tikslus dokumento vertimas į lietuvių kalbą

Vertėja Akvilė Gegelevičienė

Data 2021-04-23

UAB Diamedica

Gėlių g. 2, Avižieniai, Lietuva

**Jiangsu Mole BioScience Co.,Ltd.**

Add.: 6-7th Floor, G116 Building,

No.805, Jiankang Avenue,

Medical New&Hi-tech District,

Taizhou, Jiangsu Province, P.R.C

[www.molechina.com](http://www.molechina.com)

[info@molechina.com](mailto:info@molechina.com)



## 2019-nCoV Related Product List

Date: 19th April, 2021

No.	DESCRIPTION	KIT SIZE	Packaging Data	Certificate	EXP	Format/Picture
<b>1.1 Sample Collection &amp; Preservation &amp; Transportation Kit</b>						
1	<b>Sample Preservation Solution (Inactivated)</b> <i>Kit Content:</i> Inactivated solution - X Disposable Sampling Swab - X	20 sets/kit (Individual Packaging) / 36/50 sets/kit (Bulk Packaging)	18 kits/carton carton dimension: 40*44*45 cm	CE&UKCA	24 months	
2	<b>Disposable Collection Tube of Virus Samples (Culture)</b> <i>Kit Content:</i> Hank's solution - X Disposable Sampling Swab - X	300 pcs/kit (Individual & Sterile packaging)	18 kits/carton carton dimension: 40*44*45 cm	CE	36 months	
<b>CONFIGURATION OPTIONS:</b>						
<b>Disposable Sampling Swabs:</b>						
OPTION 1	1) Nylon Flocked Nasopharyngeal Swab with breakpoint at 10cm 2) Nylon Flocked Nasopharyngeal Swab with breakpoint at 8cm 3) Nylon Flocked Oropharyngeal Swab with 2 breakpoints at 3&8cm 4) Nylon Flocked Oropharyngeal Swab with breakpoint at 8cm 5) Nylon Flocked Oropharyngeal Swab with breakpoint at 3cm	300 pcs/kit (Individual & Sterile packaging)	18 kits/carton carton dimension: 40*44*45 cm	CE	36 months	
OPTION 2	<b>Solution Volume:</b> 1 mL, 2 mL, 3 mL and customization.					
OPTION 3	<b>Tube Options:</b> 2 mL, 5 mL, 10 mL.					
<b>1.2 Saliva Collection &amp; Preservation &amp; Transportation Kit</b>						
3	<b>Saliva Collection Kit</b> <i>Kit Content:</i> 1) Saliva Collector 2) 2 mL Preservation Solution in 5/7/10 mL Tube 3) Biosafety Bag 4) Label Sticker 5) IFU	1 set/kit (Individual packaging)  10 sets/Kit (Individual packed)	22 kits/carton	CE	36 months	
<b>1.3 Environmental Sample Collection &amp; Preservation &amp; Transportat</b>						
4	<b>Environmental Sample Collection Kit</b> <i>Kit Content:</i> 2 mL/5 mL Preservation Solution pre-filled in 10 ml Tube disposable Collection Swab <b>OR</b> 1 pc Disposable Rolling Swab	20 sets/kit (Individual Packaging) / 36 sets/kit (Bulk Packaging)	18 kits/carton carton dimension: 40*44*45 cm	For Laboratory Use	24 months	
<b>2. Sample Collection &amp; Preservation &amp; Transportation &amp; DNA/RNA Release Reagent</b>						
5	<b>Nucleic Acid Release Reagent</b> <i>Kit Content:</i> 10 mL Tube with 1 mL <b>Release Reagent</b> Disposable Sampling Swab	20 sets/kit (Individual Packaging) / 36/50 sets/kit (Bulk Packaging)	18 kits/carton carton dimension: 40*44*45 cm	CE	24 months	

poz.  
2.1



## 2019-nCoV Related Product List

Date: 19th April, 2021

No.	DESCRIPTION	KIT SIZE	Packaging Data	Certificate	EXP	Format/Picture
<b>3. Sample DNA/RNA Extraction Kit</b>						
6	<b>Viral DNA/RNA Nucleic Acid Extraction Reagent kit</b> (Magnetic Bead Method)_ <b>Not Pre-Packed</b> <i>Compatible with Mole 96M, Mole 32M, King Fisher 96, Allsheng Autore 32/96...</i>	64 Tests/kit	18 kits/carton carton dimension: 40*44*45 cm 13.5 kgs/carton		12 months	
7	<b>Viral DNA/RNA Nucleic Acid Extraction Reagent kit</b> (Magnetic Bead Method)_ <b>Pre-Packed</b> <i>Compatible with Mole Fast 96, Mole 96M, Mole 32M, King Fisher 96, Allsheng Autore 32/96...</i>	64/96 Tests/kit	18 kits/carton carton dimension: 40*44*45 cm 13.5 kgs/carton	CE	12 months	
8	<b>Viral DNA/RNA Nucleic Acid Extraction Reagent kit</b> (Magnetic Bead Method)_ <b>Single Test Pre-Packed</b> <i>Compatible with Mole 32M</i>	1 tests/tube, 32 tests/kit			12 months	
9	<b>Viral DNA/RNA Nucleic Acid Extraction Reagent Kit</b> (Spin column Method)	24/50/100 Tests/kit	50 tests/kit: 24 kits/carton carton dimension: 40*44*45 cm 6.6 kgs/carton		12 months	
10	<b>Automatic Nucleic Acid Extraction Instrument</b> <i>(Model: FAST 96), Extraction Time: 6 mins</i> <b>Lifetime Warranty.</b>	96 channels	1 Instrument/carton	Not available for EU market		
11	<b>Automatic Nucleic Acid Extraction Instrument</b> <i>(Model: 96 M), Extraction Time: 13 mins</i> <b>Lifetime Warranty.</b>	96 channels	1 Instrument/carton	CE, not available for EU market.		
12	<b>Automatic Nucleic Acid Extraction Instrument</b> <i>(Model: 32 M), Extraction Time: 13 mins</i> <b>Lifetime Warranty.</b>	32 channels	1 Instrument/carton	CE		

**NEW!!!**

**Note: FREE of Charge for distributors if purchasing the Extraction Kits from us.**



## 2019-nCoV Related Product List

Date: 19th April, 2021

No.	DESCRIPTION	KIT SIZE	Packaging Data	Certificate	EXP	Format/Picture
<b>4. Amplification &amp; Confirmation (PCR Test Kit)</b>						
12	<b>2019-nCoV(COVID -19) ORF1ab &amp; E&amp;N Gene Nucleic Acid RT-PCR Test Kit</b> (Fluorescent Probe-based Real Time PCR)	50 Tests/kit	2 kits/cartor	CE	6 months	
13	<b>COV+IFA+IFB Multiplex Nucleic Acid Test Kit</b> (Fluorescent Probe-based Real Time PCR)	50 Tests/kit	25 kits/cartor	CE	6 months	
14	<b>COV+FLU+RSV Multiplex Nucleic Acid Test Kit</b> (Fluorescent Probe-based Real Time PCR)	50 Tests/kit	25 kits/cartor	CE	6 months	
<b>5. Immunoassay Lateral Flow Rapid Tests</b>						
15	<b>SARS-CoV-2 IgG/IgM Rapid Test</b>	20 Tests/kit	80 kits/carton	CE	24 months	
16	<b>SARS-CoV-2 Antigen Rapid Test (Nasal Swab)</b> <i>(For Professional Use)</i>	20 Tests/kit	80 kits/carton	CE, Listed in Czech Official Exception List.	24 months	
17	<b>SARS-CoV-2 Antigen Rapid Test (Nasal Swab &amp; Saliva)</b> <i>(For Professional Use, Single Package)</i>	1 test/kit, 20 kits/box.		CE	24 months	 <i>Single Package Antigen Rapid Test</i>
18	<b>SARS-CoV-2+Influenza A+B Combo Rapid Test</b>	20 Tests/kit	80 kits/carton	CE	24 months	
19	<b>Influenza A+B Rapid Test Cassette</b>	20 Tests/kit	80 kits/carton	CE	24 months	
20	<b>Sealing Tube COVID-19 Antigen Rapid Test Strip</b>	2 tests/kit, 20 tests/kit	2 tests/kit: 690 kits (1380 tests/carton) 20 tests/kit: 108 kits (2160 tests/carton)	CE	18 months	

KEY PRODUCTS!!!

NEW!!!



## 2019-nCoV Related Product List

Date: 19th April, 2021

No.	DESCRIPTION	KIT SIZE	Packaging Data	Certificate	EXP	Format/Picture
<b>6. PCR Plastic Consumables</b>						
800108	<b>10ul filter tip, packed in rack, sterile, clear</b>	96 pcs/rack, 50 racks/carton	Carton dimension: 49*34*28 cm, 6.85 kgs/carton, 50 racks/carton	CE		
800208	<b>200ul filter tip, packed in rack, sterile, yellow</b> (Available volume: 1- 150 ul)	96 pcs/rack, 50 racks/carton	Carton dimension: 48*39*27.5 cm, 8.15 kgs/carton, 50 racks/carton	CE		
800308	<b>10ul filter tip (Extra long), packed in rack, sterile, clear</b>	96 pcs/rack, 50 racks/carton	Carton dimension: 49*34*28 cm, 6.85 kgs/carton, 50 racks/carton	CE		
800408-B	<b>1000ul filter tip, packed in rack, sterile, blue</b>	96 pcs/rack, 50 racks/carton	Carton dimension: 59*51.5*28.5 cm, 12.9 kgs/carton, 50 racks/carton	CE		
800508	<b>1250ul filter tip, packed in rack, sterile, blue</b>	96 pcs/rack, 50 racks/carton	Carton dimension: 59*48*27 cm, 14.45 kgs/carton, 50 racks/carton	CE		
800608	<b>200ul filter tip (Extra thick), packed in rack, sterile, clear</b> (Available Volume: 1-250 ul)	96 pcs/rack, 50 racks/carton	Carton dimension: 59*48*27 cm, 8.5 kgs/carton, 50 racks/carton	CE		
800708	<b>20ul filter tip, packed in rack, sterile, yellow</b>	96 pcs/rack, 50 racks/carton	Carton dimension: 48*39*27.5 cm, 8.15 kgs/carton, 50 racks/carton	CE		
	<b>Tecan 200 ul Conductive Tip (Sterile)</b> <i>Optional Specification:</i> Hanging Rack/Rack Low Retention/None With Filter/ Without Filter	96 pcs/rack, 50 racks/carton				Stock are sufficient
	<b>Tecan 1000 ul Conductive Tip (Sterile)</b> <i>Optional Specification:</i> Hanging Rack/Rack Low Retention/None With Filter/ Without Filter	96 pcs/rack, 50 racks/carton				Stock are sufficient



**NEW!!!**



\* Virus Transport Medium can be stored and transport in room temperature (2-30 C).  
\* DNA/RNA Extraction Kits can be transported in room temperature (15-30 C).  
\* PCR Test Kits need to be stored and transported under the temperature of -20 C +/- 5C.  
\* Antibody/Antigen Test: Room Temperature (2-30 C)



## 2019-nCoV Related Product List

Date: 19th April, 2021

No.	DESCRIPTION	KIT SIZE	Packaging Data	Certificate	EXP	Format/Picture
<b>7. Products in Development Pipeline</b>						
1	<b>SARS-CoV-2 IgM/IgG/IgA Combined Antibody Rapid Test Kit</b>	20 Tests/kit		CE		
2	<b>SARS-CoV-2 Total Antibody Test Kit (ELISA)</b>	96 Tests/kit		CE		<i>Samples are available for evaluation</i>
3	<b>SARS-CoV-2 Neutralization Antibody Test Kit (ELISA)</b>	96 Tests/kit		CE		<i>Samples are available for evaluation</i>
4	<b>SARS-CoV-2 Total Antibody Rapid Test Kit</b>	20 Tests/kit		CE		<b>Coming Soon...</b>
5	<b>FLU A/B+COV-2 Neutralization Antibody Rapid Test Kit</b>	20 Tests/kit		CE		

## Certification ISO 13485, CIGB Registration and CE DOC



# TECHNICAL DATA SHEET

<b>Product Name:</b>	<b>Disposable Collection Tubes of Virus Samples</b>	<b>Cat. No.:</b>	P045T0060100
		<b>Volume:</b>	3 mL

## Main Components

- Sodium phosphate dibasic dodecahydrate
- Potassium dihydrogen phosphate
- Magnesium sulfate heptahydrate
- Magnesium chloride hexahydrate
- Gentamycin sulfate
- Calcium chloride
- Potassium chloride
- Sodium chloride
- D(+)-Glucose
- Phenol Red

## Materials

<b>Tube</b>	Lid	Polypropylene	<b>Swabs</b>	Tip	Nylon Flocked
	Tube body	Polypropylene		Rod	Polypropylene

## Dimension

<b>Tube</b>	Length	102.2±0.1 mm	<b>Rod</b>	<b>Nasopharyngeal Swab</b>		
	Diameter	15.5±0.05 mm		Length	149 +0/-2 mm	
	Volume	10 mL		Breakpoint	81 +0/-2 mm	
<b>Lid</b>	Length	11.0±0.1 mm	<b>Tip</b>	Diameter	Φ1.05, Φ1.15, Φ2.6 mm	
	Diameter	20.3±0.05 mm		Length	24 ±2 mm	
<b>Total Length</b>		103.35±0.1 mm	Diameter	3 ±0.5 mm		

\*\*\*Please refer to the Technical Drawing for detailed Dimension Data.

## Packaging Data

Kit Size:	36 sets/kit	Box Dimension:	20*14*12 cm
Carton Size:	18 kits/carton (≈9.5 kgs/carton)	Carton Dimension:	42*40*45 cm



**Jiangsu Mole BioScience Co.,Ltd.**

Add.: 6-7th Floor, G116 Building, No.805, Jiankang Avenue,  
Medical New&Hi-tech District,  
Taizhou, Jiangsu Province, P.R.C  
[www.molechina.com](http://www.molechina.com)  
[info@molechina.com](mailto:info@molechina.com)