

Home / Products / Orthodontics / i-GEL Phosphoric Acid Etching Gel

## i-GEL Phosphoric Acid Etching Gel

I-DENTAL

1. P.D.  
21. P.D.



### PRODUCT DETAILS

## Overview

### i-GEL Phosphoric Acid Etching Gel

- For etching teeth and preparing surfaces
- Advantages:
  - Thixotropic
  - Optimum level of viscosity
  - 37% phosphoric acid
- Available in:
  - 4 x 4.3gm syringe + 5 tips
  - 50gm syringe + 20 tips

Pradžia » Produktai » Ėsdinimo medžiagos » i-GEL N

## i-GEL N

Fosforo rūgšties ėsdinimo gelis ir skystis, 37% fosforo rūgštis



### PASKIRTIS

- Danties paviršiaus ėsdinimui



### PRIVALUMAI

- Pasižymi optimaliu takumu
- 37% fosforo rūgštis
- Pašalina lipnų sluoksnį ir atveria dentino tubules



### PAKUOTĖ

REF	SUDĖTIS
GELIS:	
IGECP	4,3g švirkštas, 5 antgaliukai
IGEK1	4x4,3g švirkštai, 10 antgaliukų
IGET2B	12g švirkštas, 10 antgaliukų
IGES0B	42g švirkštas, 2x4,3g švirkštai, 20 antgaliukų
SKYSTIS:	
IGEGL	5ml buteliukas
IGEGC	10ml buteliukas

## ALABASTER

EN ISO 6873

Colour:	grey-white
Mixing ratio:	55 ml water / 100 g powder
Working time:	ca. 3-4 minutes
Setting time:	ca. 8 minutes
Setting expansion:	0.20 %
Compressive strength:	12 MPa after 1 h
Hardness:	12 MPa after 24 h

A general purpose dental plaster for investing prosthesis into the flask before polymerization.

REF 6624      sack      25 kg



## ARTICULATOR ARTI-LITE II

- adjustable
- inclination of Condil way 35°C
- Benet angle 0°-20°
- weight 540 g



REF 81030.2

### ADDITION

REF 81035	mounting plate
REF 81065	plastic loadstar - Calotte
REF 81070	magnet Splitcast system, á 2 pcs
REF 81062	matrix Splitcast system
REF 81060	adapter
REF 81059	socle for biting spoon

## FACEBOW

- facebow complete
- registration
- facebow bag

REF 81057



## OCCLUDATORS

4. P.D. 24. P.D.

REF 3401	occludator small
REF 3402	occludator large



[PRODUCTS](#)[BECOME A RESELLER](#)[MEDICAL PROFESSIONALS](#)[ABOUT US](#)[CONTACT](#)

5. P.D.  
25. P.D.

[< BACK TO CUPS](#)

## TOP CUPS PP | 3,000 CARTON | 200 ML



Viola

TOP (UPS PP are the cups to get, if you like a **modern looking practice**.

The cup is made of food safe **polypropylene (PP)** and is 100% recyvable.

The smooth yet very resistant material guarantees a **good grip** and **optinal stability**.

TOP (UPS come in 17 different colors and are **color-intensive** and **shiny**.

TOP (UPS fit all **standardized cup dispensers**.


Burns without residue.



### SIZE

200 ml

### COLORS



### CONSUMABLE UNITS

100 Pieces | Tube; 30 Tubes | Carton

### SALES UNITS

Carton

### CERTIFICATES



CUPDP14-AKZ-MC-L01	200	Dark Blue	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-TB-L01	200	Aqua Blue	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-MB-L01	200	Med Blue	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-RE-L01	200	Red	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-PI-L01	200	Pink	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-EG-L01	200	Mint Green	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-DG-L01	200	Dark Green	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228
CUPDP14-AKZ-X3-L01	200	Mixed	3,000	38 x 38 x 51.5	8.5	6	4	24	80 x 120 x 221	228

[ARE YOU A DENTIST?](#)

## INTERROCK NEW TYPE IV

EN ISO 6873

Colour:	ivory/gold-brown
Mixing ratio:	20 ml water / 100 g powder
Mixing time (manual):	15 seconds
Mixing time (vacuum):	60 seconds
Working time:	ca. 5 minutes
Setting time:	ca. 12 minutes
Setting expansion	0.07 % after 2 h
Compressive strength:	60 MPa after 1 h
Hardness:	250 MPa after 24 h

Synthetic extra hard stone with high resistance. It is suitable for making precise working model for partial denture and fixed prosthetics work. It is a thixotropic stone that gives very sharp, clear edges and has a very smooth surface. Colour gold-brown is suitable for CAD/CAM technique for powder-free scanning.



REF 0982	alu-sack	4 kg	ivory	REF 0981	alu-sack	4 kg	gold-brown
REF 0968	pail	5 kg	ivory	REF 0951	pail	5 kg	gold-brown
REF 0969	pail	20 kg	ivory	REF 0952	pail	20 kg	gold-brown
REF 0962	sack	25 kg	ivory	REF 0953	sack	25 kg	gold-brown

## INTERSTONE NEW TYPE IV

EN ISO 6873

Colour:	sahara
Mixing ratio:	23 ml water / 100 g powder
Mixing time (manual):	15 seconds
Mixing time (vacuum):	60 seconds
Working time:	ca. 5 minutes
Setting time:	ca. 10 minutes
Setting expansion	0.10 % after 2 h
Compressive strength:	49 MPa after 1 h
Hardness:	150 MPa after 24 h

Extra hard, thixotropic, synthetic stone with good flow properties. It is suitable for making master models, individual models in fixed prosthetics and models for partial denture.



REF 0983	alu-sack	4 kg
REF 0971	pail	5 kg
REF 0972	pail	20 kg
REF 0955	sack	25 kg

## INTERDUR TYPE III

EN ISO 6873

Colour:	blue
Mixing ratio:	26 ml water / 100 g powder
Mixing time (manual):	15 seconds
Mixing time (vacuum):	60 seconds
Working time:	ca. 4 minutes
Setting time:	ca. 9 minutes
Setting expansion	0.15 % after 2 h
Compressive strength:	40 MPa after 1 h
Hardness:	120 MPa after 24 h

High quality thixotropic hard stone for production of first class working models. Its extra surface density is ideal for use in removable prosthetics. It is also recommended for antagonistic models, when making precise prosthetic works.



REF 0956	pail	5 kg
REF 0958	pail	20 kg
REF 0957	sack	25 kg



**Polierer**

Silikon






**Polisher**

silicone

**12. P.D.**
**13. P.D.**
**32. P.D.**
**33. P.D.**


VE Stück/PU pieces

 6.000/min

Art.-Nr. montiert art.-no. mounted	Härte grade	Schaft shank Ø mm	 mm	 mm	
<b>Silikonpolierer für Prothesenkunststoffe</b> silicone polisher for prothesis plastics					
S195-1G	grob coarse	2,34	12	22	
S195-2G	grob coarse	2,34	10	20	
S195-3G	grob coarse	2,34	5,6	16	
S195-4G	grob coarse	2,34	10	24	
S195-5G	grob coarse	2,34	15	17	
S195-1M	mittel medium	2,34	12	22	
S195-2M	mittel medium	2,34	10	20	
S195-3M	mittel medium	2,34	5,6	16	
S195-4M	mittel medium	2,34	10	24	
S195-5M	mittel medium	2,34	15	17	
S195-1F	fein fine	2,34	12	22	
S195-2F	fein fine	2,34	10	20	
S195-3F	fein fine	2,34	5,6	16	
S195-4F	fein fine	2,34	10	24	
S195-5F	fein fine	2,34	15	17	



14. P.D.  
34. P.D.




**EVO ROCK**

Special plasticized dental die stone, plastic refined (formaldehyde free), conforming to ISO 6873:2013, Type 4

### Areas of application

Perfect for tooth arches, individual dies, implant, master and diagnostic models (functional diagnostics)

### Standard colours



### Advantages

- ✓ small stumps and fine cavity will be piece of cake
- ✓ very good layer capacity
- ✓ highly plasticized – very smooth surface – very high edge stability
- ✓ no splintering through the highest values

### Technical data:

Mixing ratio:	100 g : 20 ml powder / dest. water
Mixing time:	vacuum: 30 sec. manual: 60 sec.
Processing time:	5 - 6 min.
Setting time:	10 - 12 min.
Demoldability:	approx. 20 - 25 min.
Brinell hardness after 2h:	> 320 N/mm <sup>2</sup>
Compressive strength after 2h:	> 340 N/mm <sup>2</sup>
Expansion after 2h:	0,06 %
Expansion after 24h:	0,06 %





Manufacturers & Suppliers of Surgical, Dental,  
Manicure And Fishing Instruments.


**15. P.D.**  
**19. P.D.**  
**35. P.D.**  
**44. P.D.**

### **ALL THE ZONA INDUSTRIES INSTRUMENTS/PRODUCTS ARE:**

Visi Zona Industries instrumentai yra:

- **made of stainless steel** pagaminti iš nerūdijančio plieno;
- **are reusable** yra daugkartinio naudojimo;
- **are resistant to disinfection and sterilization** yra atsparūs dezinfekcijai ir sterilizacijai;
- **do not change any colors due to chemical disinfectants** nekeičia spalvų dėl cheminių dezinfekcinių priemonių
- **are sterilized at 135°C temperature** yra sterilizuojami 135 laipsnių temperatūroje

**ZONA INDUSTRIES**



EXPORT REGISTRATION NO : W - 104626  
SALES TAX NO. 09-05-9018-762-28  
N.T.N : 1320485-8

**Certified ISO 9001:2008**  
**ISO 13485:2003**  
**CE Marked.**







• Pul Aik, Aminabad Road,  
Sialkot - 51310- Pakistan.

☎ 0092 - 52 - 3522403  
☎ 0092 - 52 - 3522406

✉ zona@zona.com.pk  
✉ zona@brain.net.pk  
🌐 ww.zona.com.pk

**BANKERS: Habib Bank Ltd.**  
**Habib Metropolitan Bank Ltd.**



15. P.D.  
35. P.D.

17. P.D.  
37. P.D.

## DISPOSABLE DENTAL IMPRESSION SYRINGES

- Bend the Tip to right angle, cut the tip to adjust the flow

Ref.	Description	Pcs / Bag	Bags / Case
1A2013	white, transparent	50	10



## CURVED UTILITY SYRINGES (ECONOMICAL, MULTI-PURPOSE)

Pre-curved 80° tips !



- Vista Curved Tip Utility syringes place impression materials and liquid with ease
- Great for surgical site home irrigation
- Rigid curved tip tapers to 1/16" for tight areas
- Rubber plunger end provides an excellent seal
- Ideal for use as water syringes to eliminate any issues
- Packing: 50 pcs / bag, 24 bags / case

Ref. 1A2030

## SYRINGE CAPS

- Rigid color-coded plastic syringe caps
- Fits standard luer lock syringe design

Ref.	Description	Ref.	Description
1A2070	Colorless	1A2072	Sky blue
1A2071	Cream	1A2073	Black

## LUER-LOCK ENDO IRRIGATION SYRINGES

Clear Syringe with Scale



Ref.	Description	Pcs / Bag	Bags / Case
1A2040	Clear syringe, 1.2 CC	100	50
1A2041	Clear syringe, 3 CC	100	40
1A2042	Clear syringe, 5-6 CC	100	25
1A2043	Clear syringe, 10-12 CC	100	15

Black syringe type is available

NEW



## CUSTOM SYRINGE WITH CAP

Ref.	Description
1A2044-WVC	1ml (0.1 ml scale)



## MICRO APPLICATOR

- Bendable and disposable applicator for cavity lining, bonding, etching, fluoride varnishing. Assorted colors, designed for application of minuscule amounts



Designed bendable tip for reach inner surface of tooth wall

Ref.	Description
1A6541-S	Regular, 10.2 cm length
1A6542-S	Fine, 10.2 cm length
1A6543-S	Super fine, 8 cm length
1A6544	Extra-thin and extended, 10.2 cm length

Packing: 100 pcs/bucket, 10 buckets/box, 20 boxes/case

Ref.	Description
1A6545-R	4 in 1, Regular size, 10.2 cm length
1A6545-F	4 in 1, Fine size, 10.2 cm length
1A6545-S	4 in 1, Super Fine size, 10.2 cm length
1A6545-M	4 in 1, Mix Packing: Regular x 1, Fine x 1, Super Fine x 1, Extra thin x 1, 10.2 cm length

Packing: 100 pcs/bucket, 20 buckets/box, 10 boxes/case

NEW



## MICRO APPLICATOR DISPENSER

- New durable design, slight screw
- Faster, effective and save time
- Deliver applicator accurately, one at a time
- Helps eliminate cross-contamination

Ref.	Description	Pcs / Box
1A6546	Screw Type	1


## BRUSH APPLICATORS



Ref.	Description
1A6550	Yellow
1A6551	Green
1A6552	Blue
1A6553	Red

- Flexible and bendable to reach all areas, come with 4 different attractive colors
- Packing: 144 pcs / box, 100 boxes / case

## TWIN ENDS DISPOSABLE BRUSH




Ref. 1A6555

- Save cost and time, double-ended
- Snap-apart handle, easy for use

## COMPOSITE BRUSH AND REPLACEABLE TIP

- Made of rigid plastic to withstand repeated autoclaving
- Tip can be replaced, save the cost

Ref.	Description	Pcs / Bag	Bags / Case
1A6560	Tip	100	40
1A6561	Handle with tip	100	40



Ref. 1A6560 Ref. 1A6561