

PRODUCT LIST



ROSCO

NEO-SENSITABS™

Antifungal

REF No.	Product Name	Potency µg	A,B,C Product	QTY
81012	Amphotericin B	• 10	B	
82412	Caspofungin	• 5	B	
81112	Ciclopirox	50	C	
81212	Clotrimazole	10	B	
81312	Econazole	10	C	
82512	Fluconazole	25	B	
81512	Fluorocytosine	1	B	
81812	Itraconazole	10	B	
81912	Ketoconazole	15	B	
82012	Miconazole	10	B	
82212	Nystatin	50	B	
82612	Posaconazole	5	C	
87412	Terbinafine	30	B	
82312	Voriconazole	1	B	

INFORMATION

A Products - On stock.

B Products - Generally available on stock.

C Products - Not on stock - Special orders.

CE Marking

All products are CE marked according to the European Directive 98/79/EC for In Vitro Diagnostic Medical Devices.

Storage

Most tablets are stored at room temperature. All • marked products must be stored in a refrigerator at 2 - 8 °C.

Packaging

Neo-Sensitabs are packed in cartridges containing 1 x 50 tablets.
Sales unit size contains 5 cartridges x 50 tablets = 250 tablets.

Instructions for Use – Neo-Sensitabs™

The Instructions for use are available on our website www.rosco.dk
Instructions can be downloaded in PDF format in up to 8 different European languages.

*** Rosco reserves the right to make changes.**

poz. 95.4.

For Detection of Resistance Mechanisms

REF No.	Product Name	Detection of:	A,B,C Product	Tests	QTY
98015	KPC, MBL and OXA-48 Confirm Kit				
	Meropenem				
	Meropenem + Dipicolinic acid	• Carbapenemases	B	50	
	Meropenem + Cloxacillin	• Enterobacteriaceae			
	Meropenem + Phenylboronic acid Temocillin				poz. 95.2.
98016	Total MBL Confirm kit				
	Imipenem				
	Imipenem + EDTA	• MBL all species	B	50	
	Imipenem + Dipicolinic acid				
	Meropenem Dipicolinic acid				
98017	Detection of hGISA,GISA,VRE,Daptomycin resistance				
	Vancomycin	• Staphylococci	B	50	
	Teicoplanin	• Enterococci			
	Daptomycin				
98018	Detection of COLISTIN resistance in Enterobacteriaceae, P. aeruginosa and Acinetobacter				
	COLISTIN Prediffusion	• MDR Gram negatives	B	100	
98019	Total ESBL + AmpC Confirm Kit				
	Cefotaxime + Cloxacillin				poz 95.1
	Cefotaxime + Clavulanate				
	Ceftazidime + Cloxacillin	• ESBL + AmpC	B	50	
	Ceftazidime + Clavulanate				
	Cefotaxime + Cloxacillin + Clavulanate Ceftazidime + Cloxacillin + Clavulanate				
98022	Rapid ESBL Screen Kit				
	Cefotaxime + tazobactam + Indicator Cefotaxime + Indicator	• ESBL	B	50	
98023	OXA-48+CARB Screen kit (30-60 min)				
	Imipenem + Brthymol Blue Negative Control Blue	• Carbapenemases All species Lysis buffer Not needed	B	50	
98024	Neo-Rapid CARB Kit				
	Imipenem (x2) + Indicator (CARB) Diatabs CARB Negative Control Diatabs	• Carbapenemases Enterobacteriaceae P.aeruginosa	B	50	

For Detection of Resistance Mechanisms

REF No.	Product Name	Detection of:	A,B,C Product	Tests	QTY
	KPC/MBL in P.aeruginosa/Acinetobacter Kit Version2				
98025	Imipenem Imipenem + Phenylboronic acid Imipenem + Cloxacillin 4 mg Imipenem + Dipicolinic acid Imipenem + EDTA	Carbapenemases • P.aeruginosa/Acinetobacter	B	50	
					poz. 95.3.

INFORMATION

A Products - On stock.

B Products - Generally available on stock.

C Products - Not on stock - Special orders.

CE Marking

All products are CE marked according to the European Directive 98/79/EC for In Vitro Diagnostic Medical Devices.

Storage

Most tablets are stored at room temperature. All • marked products must be stored in a refrigerator at 2 - 8 °C.

Packaging

All Kits include 50 tests except kit no. 98012 and 98018 which include 100 tests.

Instructions for Use

The Instructions for use are available on our website www.rosco.dk

Instructions can be downloaded in PDF format in up to 8 different European languages.

*** Rosco reserves the right to make changes.**