



CHANGZHOU YUEKANG MEDICAL APPLIANCE CO., LTD

Jiaoxi Town, Wujin County 213116 Changzhou, Jiangsu China

TECHNICAL DATA SHEET (TDS)

Product	Stomach Tube					
Standard of test	Q/320412SYS012-2011 《disposable sterilized stomach tube for single use》					
Inspection Content						
Items	Technical requirement					
Appearance	The tube is transparent or translucent, uniform thickness, without sundries, internal or external surface is smooth,with funnel connector, four lateral eyes and open tip at the end					
	FR/CH	O.D	I.D	Length of tube	Length of connector	Color of connector
Size(mm)	F6	2.0±0.5	1.1±0.3	1200±20	35±0.2	light green
	F7	2.3±0.5	1.3±0.4	1200±20	35±0.2	bluish green
	F8	2.7±0.5	1.5±0.4	1200±20	35±0.2	light blue
	F10	3.3±0.5	2.2±0.4	1200±20	35±0.2	black
	F12	4.0±0.5	2.5±0.4	1200±20	35±0.2	white
	F14	4.7±0.5	3±0.4	1200±20	35±0.2	green
	F16	5.3±0.5	3.4±0.4	1200±20	35±0.2	orange
	F18	6.0±0.8	3.9±0.8	1200±20	35±0.2	red
	F20	6.7±0.8	4.3±0.8	1200±20	35±0.2	yellow
	F22	7.3±0.8	5±0.8	1200±20	35±0.2	violet
	F24	8.0±0.8	5.5±0.8	1200±20	35±0.2	dark blue
Strength	The connector should be connected with the tube firmly					
Negative pressure resistance	No depression under negative pressure 40Kpa pressure for 15 seconds					
Reliability of Structure (Minimal)	mm		Minimal Force			
	1.3~2.7		5			
	3.0~4.7		15			
	≥5		20			
Sterilizatio	Should be sterilized					
Sterilizatio n Test	It shows no bacterium growth after5-days cultivation for aerobe and anaerobe at 35°C.It shows no bacterium growth after7-days cultivation for mildew at 25°C.					
EO	Amount of residue≤10ug/g					
Photo of the Product	<div></div>					