

ROMEO[®]2

RODS

ROMEO[®]2
RODS

ONE DIAMETER 5.4 mm

TWO MATERIALS Ti & CoCr
STERILE IMPLANTS

HEXAGONAL TIPS STARTING 100 mm

With an equivalent diameter, the cobalt chromium rods are twice as stiff as titanium.



1. (1)
page number is Ti.

1.3.

TI RODS PREBENT Ø5.4

TI RODS STRAIGHT Ø5.4

steris

L30 mm	L100 mm
L35 mm	L120 mm
L40 mm	L140 mm
L45 mm	L160 mm
L50 mm	L180 mm
L55 mm	L200 mm
L60 mm	L220 mm
L70 mm	L240 mm
L80 mm	L350 mm
L90 mm	L500 mm
	L550 mm

steris

ROMEO[®]2

RODS



**RODS COBALT CHROMIUM
STRAIGHT Ø5.4**

L100 mm
L120 mm
L140 mm
L160 mm
L180 mm
L200 mm
L220 mm
L240 mm
L350 mm
L500 mm
L550 mm

1.

CONCEPT AND DESIGN

Since 2005 Spineart has been following its motto: quality, innovation, simplicity, by developing highly performant systems for the treatment of spinal pathologies.

The ROMEO^{®2} posterior fixation system incorporates smart technologies and simplicity to address spine cases.

The first system offering a complete range of spinal implants delivered sterile with an intuitive and compact instrumentation.

ROMEO^{®2} is a complete posterior fixation system with unique and efficient technologies that offers alternative solutions to the surgeons and their patients.

1(2.)

īpaļina aflikti nēpāvēlētus
stūbruro fiksācijā krūtštrūņējā
jūsāvēlētā jādabū se.

AT A GLANCE

COMPLETE TL FIXATION PLATFORM

STREAMLINED SCREW TIP

LOW PROFILE IMPLANTS

COMPACT SET

SAFETY & TRACEABILITY

INDICATIONS

ROMEO^{®2} lumbar posterior osteosynthesis system is intended to provide immobilization and stabilization of thoraco-lumbo-sacral spine in skeletally mature patients for the following indications:

DDD (degenerative disc disease), spondylolisthesis, trauma, spinal stenosis, curvatures, tumor, pseudoarthrosis and failed previous fusion.

1(4)

Naudējams, gēlānt nestabilus stūbruro sūzālojums,
pētraumies stūbruro slonksteliņā slōkecijs, āģenerā-
ciēs stūbruro ligo, stūbruro āģēfāme cijs, spīne-
liņius novīliis.



1.3)

Poly Mono Semi axiality

Streamlined tip screws



kompos aplink
ašy 50° , va-
dinasi involyvis
nuo ašies 25°

3.

TELESCOPIC *telescopic*
CRANIO CAUDAL ANGLULATION

AXIAL ROTATION



Setscrew with Torx Imprint

CROSS CONNECTOR
MULTIAXIAL

adjustable

- L33 to L36 mm
- L36 to L43 mm
- L43 to L55 mm
- L55 to L80 mm



Setscrew with Torx Imprint

CROSS CONNECTOR
MULTIAXIAL PREBENT

- L33 to L36 mm
- L36 to L43 mm
- L43 to L55 mm
- L55 to L80 mm

ROMEO®2
CONNECTORS

1(5)



ROD CONNECTOR PARALLEL

ROD CONNECTOR AXIAL

*Galvinox axis or parallelism stops
surroundings.*



ILIAC CONNECTOR

- L15 mm
- L20 mm
- L30 mm
- L40 mm
- L50 mm
- L60 mm



SETSCREW


Hexagonal Imprint Size 4

ROMEO[®]2



STREAMLINED TIP
 LOW PROFILE IMPLANTS
 STERILE IMPLANTS
 COMPLETE SCREW PLATFORM

ROMEO[®]2



1.1
 1.2

ROMEO[®]2



Diagrams

DIAMETER	SCREW SIZE			REDUCTION
	POLYAXIAL	MONOAXIAL		
Ø 4 mm	L25 to L45 mm	L25 to L45 mm	NA	
Ø 5 mm	L25 to L50 mm	L25 to L50 mm	L30 to L45 mm	
Ø 6 mm	L30 to L90 mm	L30 to L80 mm	L30 to L50 mm	
Ø 7 mm	L30 to L90 mm	L30 to L80 mm	L30 to L50 mm	
Ø 8 mm	L30 to L90 mm	L30 to L80 mm	NA	

1.5.

L	25 - 55 mm	60 - 90 mm
Inch	5	10

align

SETScrew
 Hexagonal Imprint Size 4

5

ROMEO® 2
POLYAXIAL SCREW

Ø5 L50



REF ELL-PS 05 50-N
LOT 5-0701
QTY 2



Mat: Ti6-Al 4-V / ISO 6832-3 / ASTM F136



(01)07640139345623(11)200204(10)5-0701

- VIS POLYAXIALE
- POLYAXIALSCHRAUBE
- TORNILLO POLIAXIAL
- VITE POLIASSIALE
- ΠΟΛΥΑΞΟΝΙΚΗ ΒΙΒΑ
- PARAFUSO POLIAXIAL
- ŚRUBA WIEŁOŚSIOWA
- POLIAKSIJAL VIDA

Consult instructions for use

Do not use if package is damaged

Do not reuse

MD



1.1.; 1.2.
vertebralis

Caution: Federal Law restricts this device to sale by or on the order of a physician



Spineart SA
Chemin du Pré-Flauri 3
1228 Plan-les-Cuates
Switzerland
Phone : +41 (0)22 670 12 00
Fax : +41 (0)22 684 83 08
www.spineart.com

US REP

Spineart USA INC
29332 Mill Creek Dr, Suite 100
Laguna Hills, CA 92653 USA
Phone: +1 949 313 8400
Fax: +1 949 304 1001
www.spineart.com

2020-02-04


CE 1250

6832_800000028597-1 23



PP screw range

A.b.



PP Screws *distal in*

	Ø5	Ø6	Ø7	Ø8
L35	ELL-PP 05 35-S	ELL-PP 06 35-S	ELL-PP 07 35-S	ELL-PP 08 35-S
L40	ELL-PP 05 40-S	ELL-PP 06 40-S	ELL-PP 07 40-S	ELL-PP 08 40-S
L45	ELL-PP 05 45-S	ELL-PP 06 45-S	ELL-PP 07 45-S	ELL-PP 08 45-S
L50	ELL-PP 05 50-S	ELL-PP 06 50-S	ELL-PP 07 50-S	ELL-PP 08 50-S
L55	ELL-PP 05 55-S	ELL-PP 06 55-S	ELL-PP 07 55-S	ELL-PP 08 55-S
L60	-	ELL-PP 06 60-S	ELL-PP 07 60-S	ELL-PP 08 60-S

lgw

Sterile 5 years



REDUCTION SCREWS

Diametra

LENGTH / DIAMETER	Ø5	Ø6	Ø7
L30	ELL-SS 05 30-S	ELL-SS 06 30-S	ELL-SS 07 30-S
L35	ELL-SS 05 35-S	ELL-SS 06 35-S	ELL-SS 07 35-S
L40	ELL-SS 05 40-S	ELL-SS 06 40-S	ELL-SS 07 40-S
L45	ELL-SS 05 45-S	ELL-SS 06 45-S	ELL-SS 07 45-S
L50	ELL-SS 05 50-S	ELL-SS 06 50-S	ELL-SS 07 50-S

gates



f.