

ACS-D Circular Stapler



Product Code	Color	Shaft Length	Outer Diameter (mm)	Lumen (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of Staples	Rows of Staples	QTY/BOX
ACS 21D	Gray	27 cm	21.5	12.5	3.5	1.5-1.8	16	2	3
ACS 25D	White		25.5	16	4.8	1.5-2.2	22		
ACS 28D	Blue		28.5	18.8			26		
ACS 31D	Green		31.5	22.7			30		
ACS 34D	Gold		34.5	25			32		

This product is not available for Germany

PPH 33D

Procedure Set for Hemorrhoid and Prolapse



Product Code	Cutting Area Diameter (mm)	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
PPH 33D	24.4	3.5	0.8 ~ 1.5	32	2	3

Aparatas užtaisytas viena kasete

Sterilus



RLS Linear Stapler 1 p.d. Linijinis siuvimo aparatas

Užtaisyamas 45 arba 60mm ilgio kasetėmis

Product Code	Stapler Line (mm)	Preloaded with unit	Max Firings	QTY/BOX
RLS3035S	30	RLS3035L	8 Leidžia atlikti 8 šūvius	3
RLS4535S	45	RLS4535L		
RLS6035S	60	RLS6035L		
RLS3048S	30	RLS3048L		
RLS4548S	45	RLS4548L		
RLS6048S	60	RLS6048L		

RLS Linear Stapler Cartridge

2 p.d. Kasetės siūlomam linijiniam siuvimo aparatui

Sterilus

Product Code	Stapler Line (mm)	Cartridge Color	Open Staple Height (mm)	Closed Staple Height (mm)	Number of staples	Number of staple rows	QTY/BOX
RLS3035L	30	Blue	Kabės aukštis 3,5mm 3.5	1.5	11	2	6
RLS4535L	45				15		
RLS6035L	60				21		
RLS3048L	30	Green	4.8	2.0	11		
RLS4548L	45				15		
RLS6048L	60				21		

Ilgis

Date 13th May 2022

To whom it may concern

DECLARATION

We, **Reach Surgical, Inc.** , based in **120 Xinxing Road, West Zone, TEDA, Tianjin, 300462, People's Republic of China**

Hereby declare that

Visi produktai sterilūs, sterilumas galioja 5 metus

1. **The sterility expiration and shelf time of below products are 5 years**

Endoscopic Linear Cutting Staplers; Single Use Loading Units for Endoscopic Linear Cutting Staplers

Linijinis siuvimo aparatas su kasetėmis

Linear Cutting Staplers with Single Use Loading Units; Single Use Loading Units for Linear Cutting Staplers

Linear Staplers with Single Use Loading Units; Single Use Loading Units for Linear Staplers

Circular Staplers with Staples


Procedure Sets for Prolapse and Hemorrhoids

Disposable Ultrasonic Scalpels

2. The sterility expiration and shelf time of below products are 3 years

Endoscopic Motorized Cutting Staplers

Yours faithfully,




To whom it may concern,

Jun 15, 2022

DECLARATION LETTER

We, **Reach Surgical, Inc.**, located at **120 Xinxing Road, West Zone, TEDA, Tianjin, 300462, People's Republic of China** declare as following:

The staple in the reload are made by titanium. **Kabès titaninès**



To whom it may concern,

Jun 15, 2022

DECLARATION LETTER

We, **Reach Surgical, Inc.**, located at **120 Xinxing Road, West Zone, TEDA, Tianjin, 300462, People's Republic of China** declare as following:

The REACH Staplers are forming B shape

Siuvant visos kabės įgauna taisyklingą B formą

