




Epistaxis Nasal Dressings

Osseous epistaxis nasal packs are used in the treatment of anterior and posterior nose bleeding. The soft round edges create less friction and reduce irritation during nasal passage insertion. The epistaxis packs expand to gently compress and effectively stop nose bleeding

- ✦ Epistaxis packs stop nosebleeds, safely and effectively
- ✦ Designed to treat nosebleeds
- ✦ Softness reduces irritation
- ✦ Quick and uniform expansion
- ✦ Fluid absorption causes the packing to expand producing a tamponading effect



REF Number	Product Name	Specification	Product Description	No. Of Piece/Package
E7184	Epistaxis Nasal Dressing	10.0cmX1.5cmX2.5cm	No String	10 Pcs, Gamma Sterile
E7185	Epistaxis Nasal Dressing	10.0cmX1.5cmX2.5cm	With String	10 Pcs, Gamma Sterile
E7186	Epistaxis Nasal Dressing	10.0cmX1.5cmX2.5cm	With String & Tube	10 Pcs, Gamma Sterile
E7187	Epistaxis Nasal Dressing	10.0cmX1.0cmX2.5cm	Thin	10 Pcs, Gamma Sterile
E7193	Epistaxis Nasal Dressing	5.5cmX1.0cmX2.5cm	Pediatric	10 Pcs, Gamma Sterile