

Begin forwarded message:

From: Peter <Peter@bricon.com>
Subject: AW: NEED your urgent answers
Date: 2022 June 15 15:43:08 EEST
To: "MEDICINA, Lukas" <_____@medicinalt.it>


Dear Peter,

1 pirkimo dalis
5 punktas

Thank you very much for the information!

Please could you sent us the brochure for **Trays for Flexible Reamers (7499-0-0099)**? Because we can't find it in your catalogues. We will use it as attachment to catalogue.

1 pirkimo dalis
5 punktas



I have some questions about **IM Reamers:** 1 pirkimo dalis

1. Do they are intended for drilling an intramedullary canal on a guide wire? **yes.** 1 punktas

3 pirkimo dalis

Also some questions about **Instrument Set for Universal Nail Removal System (7799-0-9999):**

1. Do **Basket with Fixation (7799-0-9020)** is perforated and with special divisions for instrument set and is suitable for the autoclave? **yes** 19 punktas

2. Do **1/1 Container - bottom non perf. (OS140.135)** is suitable for sterilization, large enough to accommodate the complete set of instruments together with the basket for the autoclave? **yes 18 punkt**
3. Do all instruments in this set (**7799-0-9999**) are made of stainless steel? **Yes, except the handles 1-19 punkt**
4. Do the handles for **Hammer (7799-0-9014)**, **T-Handle with adaptor (7799-0-9011)** and for **Handle with adaptor (7799-0-9012)** are made of non - slip composite materials? **yes 1 punkt**
2 punkt
5. Do **Diverter Insert GR1 from 6-9 mm, Diverter Insert GR2 from 9-12 mm, Diverter Insert GR3 from 12-15 mm** are suitable for intramedullary nail extractor (**7799-0-9000**)? Do the end of these diverter inserts are with conical thread? **yes 12-14 punkt**
6. Do the **Wrench (7799-0-9004)** is suitable for tightening extractor system components? **yes 15 punkt**
7. Do the **Screwdriver with guide (7799-0-9017)** is suitable for quick coupling handles? **yes 17 punkt**

4 punkt

Regards,
VAKONDA
Lukas Bubulas


lukas@vakonda.lt
www.medicinalt.lt



Support Hospitals in Ukraine



The war in Ukraine has created a grim picture around the world. We have all been horrified that civilian homes, kindergartens, schools and hospitals are being demolished. So many civilians have been killed and so many injured. In cities where hostilities have ended, civilians are trying to rebuild schools and hospitals first, because without them we cannot take a single step. The goal of our company is to support Ukrainian hospitals with surgical instruments and trauma implants. All possible funds will be used to buy surgical instruments and trauma treatment equipment for Ukrainian hospitals.



Please support. Support of any size is important. Products can also be supported. For products, please contact us by email or contact phone.

The project is organized by the Public Institution "Chirurgija" and coordinated with the Ministry of Health of Ukraine.

On 2022-06-15, at 08:53, MEDICINA, _____@medicinalt.it> wrote:

Good Day! dear Peter,

I have some questions about **IM Reamers**:

1. What is the length of **HUDSON drive connection drills**?
2. Do **HUDSON drive connection drills** are flexible and cannulated with a fixed (non-replaceable) head?
3. Do you have the original instrument sterilization container for IM reamers?
4. Do **Guide Rod with Olive Tip** are suitable for **HUDSON drive connection drills**?

Also some questions about **Instrument Set for Universal Nail Removal System (7799-0-9999)**:

1. Do **T-Handle with adaptor** and **Handle with adaptor** are quick coupling?
2. Do **Kardan Prolongation** is suitable for quick coupling handles?
3. Do **Blades SW 2,5; 3,0; 3,5; 4,0; 4,5; 5,0; 6,0** are suitable for quick coupling handles and are hexagonal shape?
4. Do **Extractor** is equipped with a sliding forging handle, clamping and tensioning system?
5. Do **Extraction Screw conical** is suitable for quick coupling handles? Do this **Extraction Screw conical** have "Left hand" (counterclockwise) thread?

Regards,
VAKONDA

+370-698-86-782
lukas@vakonda.lt
www.medicinalt.it

On 2022-06-14, at 15:38, pr@bricon.com wrote:

Dear esteemed Partner,

in few days only the EFORT congress 2022, Lisbon/Portugal is starting!

We from ORTHO SELECT / Bricon will be pleased to welcome you at our **booth # 2.83**

It will be our pleasure to show you our

- traumatology products
 - CMF implants
 - new Spine implants
 - reabsorbable implants from our sister company **INION**

We are looking forward meeting you coming week in Lisbon and wish you already today a safe journey.

Mit freundlichen Grüßen / with kind regards

Senior Manager - Director International Sales

<image004.png> <image005.jpg>

Bricon GmbH

Eisenbahnstrasse 100

78573 Wurmlingen, Germany

Phone: +49-7461-9336-663

Mobile: +49 173 340 0463

Fax: +49-7461-9336-665

E-Mail: pr@bricon.com

Homepage: <https://bricon.com>

Managing Director:
765502

• *Registered:* Stuttgart HRB

<image003.png>

MEDICINA, Dovile

From:
Sent: 15 June 2022 12:32
To:
Subject: Fwd: NEED your urgent answers
Attachments: image002.jpg; Untitled attachment 00038.htm; Quotation509356.PDF; Untitled attachment 00041.htm; gtc_10_2013.pdf; Untitled attachment 00044.htm

Begin forwarded message:

From: <[P
ler@bricon.com](mailto:Peter@bricon.com)>
Subject: AW: NEED your urgent answers
Date: 2022 June 15 11:28:40 EEST
To: "M
s" <[s@medicinalt.lt](mailto:ms@medicinalt.lt)>

Good morning as,

thanks for your Email. Please see my comments hereunder in your Email.

Best, Peter

Von: MEDICINA, L s <[s@medicinalt.lt](mailto:ms@medicinalt.lt)>
Gesendet: Mittwoch, 15. Juni 2022 07:54
An: r <[r@bricon.com](mailto:Peter@bricon.com)>
Betreff: NEED your urgent answers

Good Day! dear ,

I have some questions about **IM Reamers: 1 pirkimo dalis**

1. What is the length of **HUDSON drive connection drills? 480mm 2 punktas**
2. Do **HUDSON drive connection drills** are flexible **yes** and cannulated **yes** with a fixed (non-replaceable) head **yes? 2 punktas**
3. Do you have the original instrument sterilization container for IM reamers? **Yes, see 5 punktas attached offer**
4. Do **Guide Rod with Olive Tip** are suitable for **HUDSON drive connection drills? Use Guide 6 punktas wires 7499-0-1422 for reamers with Ø 6 – 7mm and wires 7499-1-1423 from Ø 7,5 – 18mm**

3 pirkimo dalis

Also some questions about **Instrument Set for Universal Nail Removal System (7799-0-9999):**

1. Do **T-Handle with adaptor** and **Handle with adaptor** are quick coupling? **Yes 1-2 punktai**
2. Do **Kardan Prolongation** is suitable for quick coupling handles? **Yes 3 punktas**
3. Do **Blades SW 2,5; 3,0; 3,5; 4,0; 4,5; 5,0; 6,0** are suitable for quick coupling handles and are hexagonal shape? **Yes and Hexagonal 5-10 punktai**
4. Do **Extractor** is equipped with a sliding forging handle, clamping and tensioning system? **Sliding: Yes, the wight/handle. Clamping: Inserts will be screwed in. Tensioning: Yes, the inserts /spreader do extend 11 punktas**
5. Do **Extraction Screw conical** is suitable for quick coupling handles? Do this **Extraction Screw conical** have "Left hand" (counterclockwise) thread? **Yes 16 punktas**

Regards,
VAKONDA
Lukas Bubulas
+370-698-86-782
lukas@vakonda.lt
www.medicinalt.lt