


Updated HARMONIC™ Hand Piece

Simplified setup

- Updated teardrop style connector plugs directly into the GEN11 Generator.
- HARMONIC™ Connector HGA11 will no longer be needed.



HARMONIC™ Gray and Blue Hand Pieces

Product Code	Description	Quantity per sales unit
HARHPGR	HARMONIC™ Gray Hand Piece Compatible with HAR23, HAR36, HARH23, HARH36, HARH45, HBC05, HC325, HDH05	1
HARHPBL	HARMONIC™ Blue Hand Piece Compatible with SNGCB, SNGHK, SNGHK2, HAR9F, HAR17F	1

HARHPGR and HARHPBL require GEN11 Software 2018-1 or later.



For any questions or inquiries, please contact your local Ethicon Sales Representative.

Please refer always to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

Ethicon Endo-Surgery (Europe) GmbH
Hummelsbütteler Steindamm 71
22851 Norderstedt, Germany

www.jnjmedicaldevices.com