





## 1.4 p.d. Artroskopinēs siūļu žirkļēs



2.7mm

Bestell-Nr. Ordering no.	Biegung Bending	Richtung Direction	Ø	Arbeitslänge Working length
BC-255-270-SH	45°	rechts / right	2.7mm	150mm
BC-251-270-SH	30°	rechts / right	2.7mm	150mm
BC-263-270-SH	15°	aufwärts / upwards	2.7mm	150mm
BC-259-270-SH	7°	aufwärts / upwards	2.7mm	150mm
BC-250-270-SH	0°	gerade / straight	2.7mm	150mm
BC-261-270-SH	7°	abwärts / downwards	2.7mm	150mm
BC-265-270-SH	15°	abwärts / downwards	2.7mm	150mm
BC-253-270-SH	30°	links / left	2.7mm	150mm
BC-257-270-SH	45°	links / left	2.7mm	150mm



3.5mm

Bestell-Nr. Ordering no.	Biegung Bending	Richtung Direction	Ø	Arbeitslänge Working length
BC-255-350-SH	45°	rechts / right	3.5mm	160mm
BC-251-350-SH	30°	rechts / right	3.5mm	160mm
BC-263-350-SH	15°	aufwärts / upwards	3.5mm	160mm
BC-259-350-SH	7°	aufwärts / upwards	3.5mm	160mm
BC-250-350-SH	0°	gerade / straight	3.5mm	160mm
BC-261-350-SH	7°	abwärts / downwards	3.5mm	160mm
BC-265-350-SH	15°	abwärts / downwards	3.5mm	160mm
BC-253-350-SH	30°	links / left	3.5mm	160mm
BC-257-350-SH	45°	links / left	3.5mm	160mm