

Nuo: RMED Product Support <product.support@radiometer.dk>
Išsiųsta:
Kam:
Tema:

Dear TS 13 p.

I can hereby confirm that the ABL90 FLEX internal quality control will not affect the number of activities on the Sensor Cassette.

Best regards

This Message Is From an External Sender

This message came from outside your organization. Use caution when opening.

Dear All,

Could you write an official confirmation that:

When ABL90 Flex is performing quality control (A, B or C) is not using test from Sensor Cassette, just using liquids from Solution pack?

Public tender administrator



Baltai pr. 36-11, LT-48196 Kaunas, Lithuania | Phone/fax.: 8-37 33 33 17/29
Mob. phone: +370 640 10809 | E-mail: biuras@vitrolab.lt | www.vitrolab.lt

Please be advised that this email may contain confidential information. If you are not the intended recipient, please notify us by email by replying to the sender and delete this message. The sender disclaims that the content of this email constitutes an offer to enter into, or the acceptance of, any agreement; provided that the foregoing does not invalidate the binding effect of any digital or other electronic reproduction of a manual signature that is included in any attachment.