

The facts about the ABL800 FLEX

Color touch screen

- Intuitive user interface
- Adjustable screen
- Configurable interface for custom profiles
- Parameter bar provides analyzer status at a glance
- Multimedia tutorials
- Onboard help system



Inlet

- Sample volume down to 35 μL
- High sample success rate of capillary samples
- Micromodes ensure accuracy of samples as small as 35 μL

TS 6 p.

FLEXQ

- Automatic sample mixing and measuring
- Measures up to three successive samples
- Queue function saves time and reduces risk of errors



Test menu

pH | $p\text{CO}_2$ | $p\text{O}_2$ | $s\text{O}_2$ | ctHb | FO_2Hb | FCOHb | FMetHb | FHHb
 FHbF | cK^+ | cNa^+ | cCa^{2+} | cCl^- | cGlu | cLac | cCrea | ctBil