




Displays and Control Elements



Index Numbers

- | | | | |
|----|------------------------|----|------------------------------------|
| 1 | On / Off & Standby key | 14 | Tube guides |
| 2 | Start/Stop | 15 | Stop flow mechanism (Robson clamp) |
| 3 | I (ESC) key | 16 | Air in line detector |
| 4 | II (SELECT) key | 17 | Pressure sensor |
| 5 | Rate keys | 18 | Peristaltic mechanism |
| 6 | DATA keys | 19 | Easy Clip, empty bag detector |
| 7 | Rate window | 20 | Tubing |
| 8 | DATA window | 21 | Drip chamber |
| 9 | Information window | 22 | Robson clamp |
| 10 | Alarm window | 30 | Carrying handle (see also page 59) |
| 11 | LCD | | |
| 12 | Door | | |
| 13 | Door handle | | |

CONSUMABLES
INFUSION SETS



0.1

1.1 VOLUMED INFUSION SETS PVC DEHP-FREE

Volumed DEHP-Free PVC Infusion Set for the administration of medication via gravity and infusion pump. Using high precision DEHP-free PVC tubing (TOTM) assuring high performance in all applications. All sets include vented drip chambers with 15µm fluid filter, Luer Lock connections and Robson Clamp providing additional safety to prevent free flow situations.

4.

Product Code	Product Description	Volume / Chamber	Roll Quantity	Additional Features
APPK01	Volumed Set, DEHP-Free PVC, 235cm, Robson Clamp and Male Luer Lock	16ml	100	
APPK01-270	Volumed Set, DEHP-Free PVC, 285cm, Robson Clamp and Male Luer Lock	19.5ml	100	
APDK11	Volumed Set, DEHP-Free PVC, 215cm, Robson Clamp, Finger Safe Injection Point and Male Luer Lock	14.5ml	100	
APPK11-270	Volumed Set, DEHP-Free PVC, 285cm, Robson Clamp, Finger Safe Injection Point and Male Luer Lock	19.5ml	100	
APDK15VJ	Volumed Set, DEHP-Free PVC, 215cm, Robson Clamp, Needle Free Connector and Male Luer Lock	14.5ml	100	
APC023	Volumed Set, DEHP-Free PVC, 235cm, Robson Clamp, Needle Free Connector and Rotating Male Luer Lock	15.5ml	100	
APPK01-BV	Volumed Set, DEHP-Free PVC, 235cm, Robson Clamp, Anti-Reflex Valve and Male Luer Lock	16ml	100	
APPK01F	Volumed Set, DEHP-Free PVC, 225cm, Robson Clamp, In Line 0.2µm Filter and Male Luer Lock	15ml	100	
APAUP	Volumed Set, DEHP-Free PVC, 235cm, Robson Clamp, In Line 1.2µm Filter and Male Luer Lock	16ml	100	
APPPN2	Volumed Set, DEHP-Free PVC, 285cm, Robson Clamp, In Line 1.2µm Filter and Male Luer Lock	19.5ml	100	
APPK21	Volumed Set, DEHP-Free PVC, 265cm, Robson Clamp, 2 Finger Safe Injection Points, Slide Clamp and Male Luer Lock	18ml	100	
APPK2NF	Volumed Set, DEHP-Free PVC, 265cm, Robson Clamp, 2 Needle Free Connector, Slide Clamp and Male Luer Lock	18ml	100	
AAIP	Volumed Set, DEHP-Free PVC, 285cm, Finger Safe Injection Point, In Line Filter 1.2µm and Male Luer Lock	19.5ml	100	
APPK15C	Volumed Set, DEHP-Free PVC, 250cm, Robson Clamp, Anti-Reflex Valve, 3-Way Stopcock and Male Luer Lock	16.5ml	100	
APDK15VJ-ARF	Volumed Set, DEHP-Free PVC, 220cm, Robson Clamp, Needle Free Connector, Anti-Reflex Valve and Rotating Male Luer Lock	16ml	100	
APPK21AV	Volumed Set, DEHP-Free PVC, 235cm, Robson Clamp, Anti-Reflex Valve, 2 Finger Safe Injection Points, Slide Clamp and Male Luer Lock	16.5ml	100	
APC024	Volumed Set, DEHP-Free PVC, 255cm, Robson Clamp, 2 Needle Free Connectors, Anti-Reflex Valve and Rotating Male Luer Lock	17ml	100	

PARAMETRŲ VERTIMAS IŠ ANGLŲ KALBOS

1-2 psl.

4.1. *Lašinės sistemos linijos* skirtos Volumed μ 7000 chroma tūrinei pompai.

4. (APPK(0)) *lašinės sistemos* skirtos Volumed μ 7000 chroma tūrinei pompai.

VERTIMAS TIKRAS
Direktorius
Algis Bakutis
2022-11-28

Algis Bakutis

