

158 per.

Tercross™

0.014" OTW PTA Balloon Dilatation Catheter

A workhorse PTA Balloon dedicated for Btk* vessels

Outstanding crossability

- ▶ Lowest balloon crossing profile, hydrophilic coated, provides reliable crossability
- ▶ Extreme support, OTW**, seamless, polymer shaft offers enhanced pushability, adequate kink resistance and fast deflation times

Extensive pressure working range

- ▶ The balloon can be inflated up to a Rated Burst Pressure (RBP) of 20 atm (2026 kPa) for efficient and successful recanalization

Wide size mix

- ▶ A wide choice of balloon diameters (Ø 1.25 - 4.0 mm) and lengths (↔20 - 200 mm) cover all clinical needs


PVS
Peripheral Vascular Solutions

*Btk= Below the Knee
**OTW= Over the Wire

TERUMO®

Tercross™

158 *per*



Specifications

Catheter design: Over the Wire
 Guide wire compatibility: 0.014" / 0.36 mm
 Usable shaft length: 100 cm, 148 cm
 Balloon folding: 2 folded (Ø1.25 mm, 1.5 mm)
 3 folded (Ø2.0-4.0 mm)

Coating: hydrophilic
 Shaft diameter: proximal: Ø1.25 mm, 1.5 mm: 3.2 Fr (1.07 mm)
 Ø2.0-4.0 mm: 3.6 Fr (1.21 mm)
 distal: Ø1.25 mm, 1.5 mm: 2.5 Fr (0.83 mm)
 Ø2.0-4.0 mm: 3.0 Fr (1.00 mm)

Ordering information

Shaft length: 100 cm

Balloon diameter	Balloon length / Product code					
	20 mm	40 mm	80 mm	120 mm	150 mm	200 mm
1.25	BD-T1220P4E	-	-	-	-	-
1.50	BD-T1520P4E	-	-	-	-	-
2.00	-	BD-T2040P4E	BD-T2080P4E	BD-T20120P4E	BD-T20150P4E	BD-T20200P4E
2.50	-	BD-T2540P4E	BD-T2580P4E	BD-T25120P4E	BD-T25150P4E	BD-T25200P4E
3.00	-	BD-T3040P4E	BD-T3080P4E	BD-T30120P4E	BD-T30150P4E	BD-T30200P4E
3.50	-	BD-T3540P4E	BD-T3580P4E	BD-T35120P4E	BD-T35150P4E	BD-T35200P4E
4.00	-	BD-T4040P4E	BD-T4080P4E	BD-T40120P4E	BD-T40150P4E	BD-T40200P4E

Shaft length: 148 cm

Balloon diameter	Balloon length / Product code					
	20 mm	40 mm	80 mm	120 mm	150 mm	200 mm
1.25	BD-T1220Q4E	-	-	-	-	-
1.50	BD-T1520Q4E	-	-	-	-	-
2.00	-	BD-T2040Q4E	BD-T2080Q4E	BD-T20120Q4E	BD-T20150Q4E	BD-T20200Q4E
2.50	-	BD-T2540Q4E	BD-T2580Q4E	BD-T25120Q4E	BD-T25150Q4E	BD-T25200Q4E
3.00	-	BD-T3040Q4E	BD-T3080Q4E	BD-T30120Q4E	BD-T30150Q4E	BD-T30200Q4E
3.50	-	BD-T3540Q4E	BD-T3580Q4E	BD-T35120Q4E	BD-T35150Q4E	BD-T35200Q4E
4.00	-	BD-T4040Q4E	BD-T4080Q4E	BD-T40120Q4E	BD-T40150Q4E	BD-T40200Q4E



Terumo Corporation
 44-1, 2-chome, Hatagaya
 Shibuya-ku, Tokyo 151-0072, Japan
 Tel.: +81 3 3374 8111
 Fax: +81 3 3374 8196

Terumo Europe N.V.
 Interleuvenlaan 40
 3001 Leuven, Belgium
 Tel.: +32 16 38 12 11
 Fax: +32 16 40 02 49

**Terumo Europe
 Interventional Systems
 Liaison Office**
 Immeuble Krebs
 8, rue Jean-Pierre Timbaud
 CS 10720
 78067 Saint Quentin
 en Yvelines Cedex, France
 Tel.: +33 139 30 95 75
 Fax: +33 139 30 95 70

**Terumo Europe N.V.
 Benelux Sales Division**
 Belgium:
 Tel.: +32 16 39 25 80
 The Netherlands:
 Tel.: 0800 0220396

**Terumo Europe N.V.
 Emerging Markets Division**
 Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
 Tel.: +49 6196 80 230

**Terumo Deutschland GmbH
 Zweigniederlassung Spreitenbach**
 Tel.: +41 56 419 10 10

**Terumo Deutschland GmbH
 Zweigniederlassung Österreich
 Interventional Systems**
 Tel.: +43 2236 378020

Laboratoires Terumo France S.A.
 Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
 Tel.: +39 0651 96 14 20

Terumo Europe España SL
 Tel.: +34 9021 01 298

Terumo UK Ltd
 Tel.: +44 1784 476 200

Terumo Russia LLC
 Tel.: +7 495 988 4740

Terumo Sweden AB
 Tel.: +46 3174 85 880

**Terumo Denmark
 Filial of Terumo Sweden AB**
 Tel.: +45 7020 93 80

**Terumo Corporation Dubai Branch
 (Middle East and Africa)**
 Tel.: +971 4 221 2220



* Registered Trademark
 Published by Terumo Europe N.V.