


ADP.



Terumo Interventional Oncology Division
Introduces

HydroPearl[®]

Compressible Microspheres for Embolisation

BECAUSE
EXPERIENCE
MATTERS


TERUMO
INTERVENTIONAL
SYSTEMS

HydroPearl®

Compressible Microspheres for Embolisation

170 µm

- New formulation¹
- Tighter calibration²
- Broad size range
- More resilient microspheres³



ORDERING INFORMATION

HydroPearl Sizes (µm)	Volume of Microspheres (ml)	Volume of PBS (ml)	Product Codes	Colour
75 ± 30	2	4	8HP2S75	Orange
200 ± 75	2	4	8HP2S200	Yellow
400 ± 75	2	4	8HP2S400	Blue
600 ± 75	2	4	8HP2S600	Red
800 ± 75	2	4	8HP2S800	Green
1100 ± 75	2	4	8HP2S1100	Purple

Presentation: 2 ml in Syringe | **Core Material:** Poly-ethylene Glycol (PEG) | **Storage Media:** Phosphate Buffered Saline | **Catheter Compatibility:** < 0.041" for all sizes

¹ Polyethylene Glycol (PEG) ² Data on file. Comparison with EmboSpheres Merit Medical Ref HPMCV-001 ³ Data on file. Comparison with EmboSpheres, Merit Medical Ref HPMCV-001

The HydroPearl is intended to occlude blood vessels for therapeutic and adjunctive purposes. For complete list of indications, please refer to the instructions for use provided with the product.

HydroPearl is not available for sale in all countries. Please contact your Terumo local sales representative for more information. HydroPearl has not been approved by the FDA and is therefore not available for sale in the U.S.A.

TERUMO
INTERVENTIONAL
SYSTEMS

Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8399

Terumo Europe NV
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

Terumo Europe NV
Interventional Systems
Liaison Office
Immeuble Inside
7-9 Rue Eugène et Armand Peugeot
92500 Rueil Malmaison,
France
Tel.: +33 1 47 16 09 30
Fax: +33 1 47 16 09 49

Terumo Europe NV
Benelux Sales Division
Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

Terumo Europe NV
Emerging Markets Division
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
Tel.: +49 6196 80 230

Terumo Deutschland GmbH
Zweigniederlassung
Switzerland
Tel.: +41 56 419 10 10

Terumo Deutschland GmbH
Zweigniederlassung Österreich
Tel.: +43 2236 378020

Terumo France S.A.S.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
Tel.: +39 0651 96 14 20

Terumo Europe España SL
Tel.: +34 9021 01 298

Terumo UK Ltd
Tel.: +44 1276 480440

Terumo Russia LLC
Tel.: +7 495 988 4740

Terumo Sweden AB
Tel.: +46 3174 85 880

Terumo Corporation Dubai Branch
(Middle East and Africa)
Tel.: +971 4 221 2220

Terumo BCT Tıbbi Cihazlar Dağıtım
ve Hizmetleri A.Ş.
Tel.: +90 216 590 0074



© Registered Trademark
Published by Terumo Europe NV