


## Centrifugation

---



- Two phase centrifugation
- 5 Minutes
- 24 Positions

---

## Sample Throughput

---

### •ABO Test

- Creation of Erythrocyte Suspension and Pipetting = 40 Seconds
- No Incubation
- Centrifugation = 5 Minutes
- Reading/Interpretation = 20 Seconds

### •ABS + Auto Control / Crossmatch

- Creation of Erythrocyte Suspension and Pipetting = 40 Seconds
- Incubation = 10 Minutes
- Centrifugation = 5 Minutes
- Reading/Interpretation = 20 Seconds