

BIOMONITOR III

Injectable cardiac monitor



1,2,3,4,5,6,7,8,9,10,11 techninēs charakteristikos

Ordering Information

Model	Order number
BIOMONITOR III	436066
Remote Assistant III	435292

Product Highlights

One-step injection procedure

BIOvector for reliable arrhythmia detection

High quality signal, easy rhythm classification

Excellent ECG display due to lossless signal compression

Additional clinically relevant information due to AF-ECG episodes with details regarding onset and termination

Intelligent memory management

ProgramConsult

BIOTRONIK Home Monitoring®

4 years of longevity

ProMRI

Full-body MR scan with 1.5 T & 3.0 T

BIOMONITOR III

Technical Data

Sensing parameters	
ECG quality	8 bit – 128 Hz bandwidth
R-wave sensing	Detection based on 1-vector ECG signal
Detection settings	
Atrial fibrillation - Detection	
Atrial fibrillation (AF)	OFF; ON
AF sensitivity	Low; Medium; High; Individual
Detection/termination window size	8/16; 16/24; 24/32
Detection intervals	5 ... [2] ... 23
Number of detection windows	1 ... (1) ... 4
Termination intervals	1 ... [2] ... 7
Number of termination windows	1 ... (1) ... 4
RR variability limit	6 ... [3] ... 18%
Confirmation time	1 ... (1) ... 6; 10 ... [10] ... 30 min
Bigeminy rejection	OFF; Standard; Aggressive
High ventricular rate - Detection	
High ventricular rate (HVR)	OFF; ON
HVR limit	100 ... [10] ... 200 bpm
HVR counter	8 ... [4] ... 24; 32; 48 cycles
Bradycardia - Detection	
Bradycardia	OFF; ON
Brady zone limit	30 ... (5) ... 80 bpm
Brady duration	5 ... [5] ... 30 s
Sudden rate drop - Detection	
Sudden rate drop (SRD)	OFF; ON
SRD rate decrease	20 ... [10] ... 70%
SRD sensitivity	Low; Medium; High; Individual
Baseline intervals	48; 64; 128; 256
Rate-drop intervals	8; 16; 32
Asystole detection	
Asystole	OFF; ON
Asystole duration	2 ... (1) ... 10 s
Patient trigger	OFF; ON
Indication-based programs	
ProgramConsult	Syncope; Palpitations; AF monitoring; Cryptogenic stroke
Sensing settings	
SensingConsult	Standard; Sense after large PVCs; Sense small PVCs; Sense short intervals; T-wave suppression
Sensing filter	4.5; 10; 18; 24 Hz
(SECG) Display	Normal; Inverted
(SECG) Signal filter	0.05; 0.5 Hz
Resting rate period	
Start resting period	00:00 ... [01:00] ... 23:00 hh:mm
Resting period duration	1 ... (1) ... 12 h
Storage parameters	
Total storage capacity	At least 60 min
Number and duration of episodes	60 episodes (56 automatically triggered, 4 patient-triggered)
Automatically triggered ECGs	56 episodes, at least 40 s each (at least 30 s before and 10 s after the event trigger)
Patient-triggered ECGs	4 episodes, at least 7.5 min each (7 min before and 30 s after the event trigger)
Intelligent memory management	
Patient-triggered episodes	4 overwrite-protected: 4 most recent
HVR episodes	3 overwrite-protected: oldest, most recent, longest
Bradycardia episodes	3 overwrite-protected: oldest, most recent, longest
SRD episodes	3 overwrite-protected: oldest, 2 most recent
Asystole episodes	3 overwrite-protected: oldest, 2 most recent
AF episodes	3 overwrite-protected: oldest, most recent, longest

Diagnostics

Event messages and counters	Atrial fibrillation (AF)
	High ventricular rate (HVR)
	Bradycardia
	Sudden rate drop (SRD)
	Asystole
AF details	Patient trigger
	Longest AF episode
	AF burden
	AF trend
	AF time of occurrence
Activity	AF duration
	Ventricular rate during AF
	Activity trend
	Mean and resting rate trends
	Heart rate variability
Sensing	Rate histogram
	R-wave trend
	Noise duration trend

Housing

Sensing vector length	70 mm
Dimensions (WxHxL) (overall)	77.5 mm x 8.6 mm x 4.6 mm
Dimensions (WxHxL) (housing)	47.5 mm x 8.3 mm x 4.3 mm
Weight	4 g
Volume (overall)	1.9 cm ³
Volume (housing)	1.7 cm ³
X-ray identification	



Electrode coating	Fractal iridium coating
Service time	
Average service time	4 years ¹⁾
Battery charge	Display 0 – 100%

¹⁾ 6-month storage period
Home Monitoring: ON; sensing: 60 bpm; daily device message via Home Monitoring including one automatic subcutaneous ECG per day and 2 patient-triggered subcutaneous ECGs per month

MR conditional

ProMRI	Please refer to the technical manual "ProMRI MR conditional device systems" for detailed information
Without scan exclusion zone and post-operative waiting period	1.5 T and 3.0 T

Tools

FIT OneStep injection tool	201 mm x 24 mm (Ø handle)
Incision tool	130 mm x 13 mm (Ø handle)

Remote Assistant III

Dimensions (WxHxL)	96 mm x 58 mm x 20 mm
Weight	85 g
Battery type	2 x LR1, alkaline-manganese battery; Inaccessible to the user; Not replaceable

BIOTRONIK Home Monitoring®

Programmer settings

Home Monitoring	OFF; ON
Time of transmission	STD; 00:00 ... [00:30] ... 23:30 hh:mm

Home Monitoring Service Center settings

Periodic subcutaneous ECG for	OFF; daily, 7 ... (1) ... 366 days
Home Monitoring-supported follow-up	

Transmitted data

Clinical and technical data	Battery status, programmable parameters, diagnostic data and patient data, recorded subcutaneous ECGs
ECG	Up to 6 per day