




## INSTRUMENT BASKETS

Instrumentiniai krepšeliai

# BASKET CONFIGURATOR

## ILLUSTRATIVE OVERVIEW FOR JF BASKETS

### WHICH BASKET GOES WITH WHICH CONTAINER BOTTOM?




#### JF SERIES

- Black PEEK plastic feet for stacking
- Fold-down handles on the side
- Mat surface
- 4 mm base perforation / 2 mm side perforation
- Can be combined with Aesculap InOS®
- Sloped edges

1/1 short


		Top basket								
		Basket in combination with:								
		Alone	JF211R	JF211R with lid*	JF212R	JF212R with lid*	JF213R	JF213R with lid*	JF214R	JF214R with lid*
Bottom basket	JF211R	JK440 JN440	JK441 JN441	JK441 JN441	JK441 JN442	JK442 JN442	JK442 JN444		JK444 JN444	JK444 JN444
	JF211R with lid*	JK440 JN440	JK441 JN441	JK441 JN442	JK442 JN444	JK442 JN444		JK444 JN444	JK444 JN446	JK444 JN446
	JF212R	JK440 JN441	JK441 JN442	JK442 JN442	JK442 JN444		JK444 JN444			
	JF212R with lid*	JK441 JN441	JK442 JN444	JK442 JN444		JK444 JN444		JK444 JN446		JK446 JN446
	JF213R	JK441 JN442			JK444 JN444				JK446 JN446	
	JF213R with lid*	JK442 JN444	JK444 JN444	JK444 JN444		JK444 JN446		JK446 JN446		
	JF214R	JK444 JN444	JK444 JN446		JK444 JN446	JK446 JN446	JK446 JN446			JK446
	JF214R with lid*	JK444 JN444	JK444 JN446		JK446 JN446			JK446		-

\*JF217R lid



#### 1 EXAMPLE

- Basket alone: JF212R (without lid)
- ➔ suitable container bottom: JK440 or JN441



#### 2 EXAMPLE

- Top basket: JF212R (without lid)
- Bottom basket: JF213R with lid
- ➔ suitable container bottom: JK444 or JN444

JK = Solid bottom  
(fractionated vacuum process)



JN = Perforated bottom  
(gravity process)



**NOTE:** ■ According to EN ISO 868-8, the load must not exceed a total of 10 kg in order to avoid excessive condensate formation.  
■ Caution! The inner height and therefore the appropriate basket changes depending on whether JK or JN bottoms are used.

# BASKETS PROTECT VALUABLE INSTRUMENTS



## GOOD REASONS FOR CHOOSING AESCULAP® INSTRUMENT BASKETS:

- Comprehensive product range offering solutions not only for use with the AESCULAP® sterile container system but also for use with single-use wrap
- Various handle designs for an easy aseptic presentation in the operating room
- Close-meshed side walls to prevent breakage of fine instrument tips
- Lids may be used additionally to better protect very fine instruments
- Compatible with the Aesculap InOS® System
- Rounded feed for stacking made of thermostable plastic or of stainless steel
- Baskets with large perforation (8 mm) specially suited for cleaning and disinfection of instruments with rough debris (e. g. bone chips)

## BASKET WITH FEET WITH LARGE PERFORATION











- Handles lockable in three positions
- Large base perforation: 8 mm

	Standard 1/2	Outside dim. Inside dim.	Standard 3/4	Outside dim. Inside dim.	Standard 1/1 short	Outside dim. Inside dim.	Standard 1/1 long	Outside dim. Inside dim.
56	JH112R	243 x 253 x <b>56</b> 241 x 251 x 42	JH252R	406 x 253 x <b>56</b> 404 x 251 x 42	JH212R	485 x 253 x <b>56</b> 483 x 251 x 42	JH222R	540 x 253 x <b>56</b> 538 x 251 x 42
76	JH113R	243 x 253 x <b>76</b> 241 x 251 x 62	JH253R	406 x 253 x <b>76</b> 404 x 251 x 62	JH213R	485 x 253 x <b>76</b> 483 x 251 x 62	JH223R	540 x 253 x <b>76</b> 538 x 251 x 62
106	JH114R	243 x 253 x <b>106</b> 241 x 251 x 92	JH254R	406 x 253 x <b>106</b> 404 x 251 x 92	JH214R	485 x 253 x <b>106</b> 483 x 251 x 92	JH224R	540 x 253 x <b>106</b> 538 x 251 x 92









## BASKET WITH FEET



		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 36	JF111R	243 x 253 x <b>36</b> 241 x 251 x 22	JF251R	406 x 253 x <b>36</b> 404 x 251 x 22	JF211R	485 x 253 x <b>36</b> 483 x 251 x 22	JF221R	540 x 253 x <b>36</b> 538 x 251 x 22
 56	<b>JF112R</b> 15.2	243 x 253 x <b>56</b> 241 x 251 x 42	JF252R	406 x 253 x <b>56</b> 404 x 251 x 42	<b>JF212R</b> 15.1	485 x 253 x <b>56</b> 483 x 251 x 42	JF222R	540 x 253 x <b>56</b> 538 x 251 x 42
 76	JF113R	243 x 253 x <b>76</b> 241 x 251 x 62	JF253R	406 x 253 x <b>76</b> 404 x 251 x 62	JF213R	485 x 253 x <b>76</b> 483 x 251 x 62	JF223R	540 x 253 x <b>76</b> 538 x 251 x 62
 106	JF114R	243 x 253 x <b>106</b> 241 x 251 x 92	JF254R	406 x 253 x <b>106</b> 404 x 251 x 92	JF214R	485 x 253 x <b>106</b> 483 x 251 x 92	JF224R	540 x 253 x <b>106</b> 538 x 251 x 92

## BASKET WITHOUT FEET



		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.		Outside dim. Inside dim.
 24	JG111R	243 x 253 x <b>24</b> 241 x 251 x 22	JG251R	406 x 253 x <b>24</b> 404 x 251 x 22	JG211R	485 x 253 x <b>24</b> 483 x 251 x 22	JG221R	540 x 253 x <b>24</b> 538 x 251 x 22
 44	JG112R	243 x 253 x <b>44</b> 241 x 251 x 42	JG252R	406 x 253 x <b>44</b> 404 x 251 x 42	JG212R	485 x 253 x <b>44</b> 483 x 251 x 42	JG222R	540 x 253 x <b>44</b> 538 x 251 x 42
 64	JG113R	243 x 253 x <b>64</b> 241 x 251 x 62	JG253R	406 x 253 x <b>64</b> 404 x 251 x 62	JG213R	485 x 253 x <b>64</b> 483 x 251 x 62	JG223R	540 x 253 x <b>64</b> 538 x 251 x 62
 94	JG114R	243 x 253 x <b>94</b> 241 x 251 x 92	JG254R	406 x 253 x <b>94</b> 404 x 251 x 92	JG214R	485 x 253 x <b>94</b> 483 x 251 x 92	JG224R	540 x 253 x <b>94</b> 538 x 251 x 92