



SWISS DESIGN AND QUALITY



Home

Generate a Quality Certificate

Downloads

Quality Management

Technical Information

TPP Worldwide

About TPP

Contact

Products

Tissue Culture Flasks

Tissue Culture Flasks with re-closable lid

Tissue Culture Flasks with peel-off foil

Clipmax

Tissue Culture Tubes

Cell Scraper and Cell Spatula

TubeSpin® Bioreactor

TubeSpin® Bioreactor with Septum

"easy read" measuring device

PCV Packed Cell Volume Tube

Serological Pipettes

Pipettor Turbo-Fix

Tissue Culture Test Plates

Tissue Culture Dishes

Vacuum Filtration "rapid"-Filtermax

Syringe Filters

Centrifuge Tubes

Cryo Tubes

Racks and Cryo Boxes

Triple bagged products (3-B)

Tube Racks and Cryo Boxes

Advantages of the TPP racks for centrifuge tubes 15 ml / 50 ml:

- alpha-numerical inscription to identify your samples
- a convenient form
- heat resistant autoclavable material - racks need to be disassembled when heat treated
- collapsible for storage: TPP racks with the "click" system are quickly assembled and disassembled. They show no inhibitory burrs on the insertion holes

Additional advantages for the TPP rack for tissue culture flat tubes 10:

- transparent top
- protection rim that prevents horizontally stored tissue culture flat tubes 10 from shifting

Advantages für the TPP rack for TubeSpin Bioreactor 600:

- re-usable
- stackable
- secure footprint
- fits 3 x vessels



Advantages of the TPP cryo boxes:

- intended for storage of samples in a freezer
- even after multiple usages the lids can be removed from - or replaced on - the base with ease and one-handed
- for cryo tubes larger than 4 ml boxes with raised lids are available

Advantages of the TPP cryo rack with star-shaped insertions:

- star-shaped insertion system in the grooves allows one hand operation when used with the star foot TPP cryo tubes
- excellent non-slip properties
- ready-mounted rubber feet

Good to know!



Space Saver

- The click-system allows easy assembling and disassembling of the racks. A collapsed rack saves storage space when not in use.

Sterilization

- We recommend to heat treat racks only when disassembled.






Downloads

- Instructions for Use / IFU
- Mounting/Demounting instructions
- Chemical Resistance

Type	Product No.	Capacity (pcs)	Dimensions (mm)	Material	Qty / Bag	Qty / Case	Case (kg)
	99013	3 x 600 ml (87600)	366 x 146 x 70	PP	1	6	2.4



	99016	40 x cryo tubes with star foot	109 x 209 x 24	PP	1	12	1.4
	99017	18 x 15 ml / 10 x 50 ml	209 x 98 x 62	PP	1	45	5.1
	99018	8 x 91243 / 8 x 15 ml or 8 x 5 ml	200 x 88 x 72	PP	1	28	3.0
	99019	30 x 15 ml / 20 x 50 ml	209 x 172 x 62	PP	1	30	5.5