

Invitrogen™

Barrier (Filter) Tips, 200 µL size

Catalog number: AM12655**Related applications:** [RNA Extraction](#) [PCR](#) [Real Time PCR \(qPCR\)](#)**Catalog Number****AM12655****Unit Size**

10 racks

Price (USD)[Contact Us](#) ›[Save to list](#)

Find the Essentials lab products you need

Ambion® Barrier tips prevent cross-contamination of samples by means of an aerosol filter. They are RNase- and DNase-free and are ideal for handling RNA, PCR applications, and for working with radioisotopes or toxic chemicals. Ten racks containing 8 x 12 tips each are provided. They are compatible with Eppendorf®, Finnpiquette, Oxford Benchmate®, Pipetman™, Proline®, Rainin, Socorex, and other brands of pipettors. Each lot of tips is subjected to rigorous testing for RNase contamination and is certified to be nuclease-free.

For Research Use Only. Not for use in diagnostic procedures.

Specifications

Volume (Metric)	200 µL
Filtered	Filtered
Sterility	Sterile
Autoclavable	Autoclavable
Format	Rack
Tip Style	Universal Fit
Shipping Condition	Room Temperature
For Use With (Equipment)	Universal Fit
Product Line	Invitrogen™

Contents & Storage

Store at room temperature.



[Shop now](#)






Find the Essentials lab products you need

Convenient, on-site access to the products you need. [Learn more.](#)

Documents & Downloads

Certificates

Search by lot number or partial lot number

Lot #	Certificate Type	Date	Catalog Number(s)
 00322482	Certificates of Analysis	Jun 23, 2022	AM12655
 28621482	Certificates of Analysis	Feb 23, 2022	AM12655
 16121482	Certificates of Analysis	Nov 09, 2021	AM12655
 10721482	Certificates of Analysis	Jul 21, 2021	AM12655
 20620482	Certificates of Analysis	Feb 27, 2021	AM12655

0 results displayed, search above for a specific certificate

[Request a Certificate](#)

[Shop now](#)

Find the Essentials lab products you need

Product Information



[Customer_Information_letter_Plastics_Tips_Final](#)

Manuals



[RNase-free Tips](#)