

[Home](#) ▶ [Shop All Products](#) ▶ [Microplates](#) ▶ [Nucleic Acid Reaction Assay Microplates](#) ▶

Need help choosing the right PCR/qPCR plastics? Use our on-line selection tool »



Applied Biosystems™

# MicroAmp™ 96-Well Tray/Retainer Set for Veriti™ Systems

**Catalog number:** 4381850**Related applications:** [PCR](#) [Real Time PCR \(qPCR\)](#)[Technical Support](#) | [Customer Service](#)

<b>Catalog Number</b>	<b>4381850</b>
<b>Unit Size</b>	10 trays
<b>Price (EUR)</b>	<a href="#">Contact Us</a> ›

[Save to list](#)[Product Overview](#)[Documents](#)[FAQ](#)

designed to snugly hold MicroAmp™ 8-Tube Strips (0.2 mL) and MicroAmp™ Reaction Tubes without Caps (0.2 mL) in a 96-well format for processing on the Veriti™ 96-Well Thermal Cycler.

- Engineered for efficient sample processing
- Designed specifically for MicroAmp™ 8-Tubes Strip (0.2 mL) and MicroAmp™ Reaction Tubes without Caps (0.2 mL)
- Manufactured to high specification, ensuring a quality product

**Convenience and Performance** The MicroAmp™ 96-Well Tray/Retainer Set for Veriti™ Systems allows you to organize your tubes while adding sample and mix. This handy tray reduces spillage and tube mix-ups by securing and ordering your tubes on the lab bench. When the samples are ready, simply put the tray of tubes directly on the block of the Veriti™ 96-Well Thermal Cycler. Not only does the tray eliminate the hassle of transferring the tubes, but it also provides the thermal contact necessary for thermal cycling.

**Note:** See user's manual or package insert for limited label license and trademark information.

## Specifications

<b>Certifications/Compliance</b>	Certified DNA/RNase Free
<b>Product Line</b>	MicroAmp™
<b>Color</b>	Blue
<b>Product Type</b>	Tray/Retainer Set
<b>For Use With (Equipment)</b>	SimpliAmp™ Thermal Cycler, Veriti Dx Thermal Cycler, Veriti Thermal Cycler, ProFlex™ PCR System
<b>For Use With (Application)</b>	PCR
<b>Plate Compatibility</b>	96-Well Plates
<b>Quantity</b>	10 Trays
<b>Shipping Condition</b>	Room Temperature

## Contents & Storage

Contains 10 MicroAmp™ 96-Well Tray/Retainer Sets. Store at room temperature.




## Certificates

Search by lot number or partial lot number



Search

Lot #	Certificate Type	Date	Catalog Number(s)
 <b>SG24110-5B</b>	Certificates of Analysis	Jan 29, 2019	4381850

1 result displayed, search above for a specific certificate

[Request a Certificate](#)

## Product Information

### Manuals



Engineering Diagram: MICROAMP 96 WELL TRAY/RETAINER SET

## Frequently asked questions (FAQs)

Why does the SimpliAmp Thermal Cycler starter kit contain two different types of trays?

Answer +

Why are my PCR tubes crushed/deformed after my PCR run?

Answer +

