

\ Treatment \ Post-operative Redon systems \ Low vacuum system

500ml



low vacuum system for routine use in operating theatre
 product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains
- set with 2 drains and guiding needle
- length of the drain: 800 mm
- perforation of the drain: 150 mm
- individually packed and sterilized
- special designs available on request
- UMDNS-Code: 11-307
- Klasse/Regel: IIa/7

Art.-Nr.: 21586, PRI-LOW-VAC Redon set, 500 ml with drain (500/150) Ch. 16

Packing unit (VPE): 35, Shipping unit (VSE): 35 **P.d. 39.1**



low vacuum system for routine use in operating theatre
 product features:

- starting with a vacuum of 120 mbar
- content: 500 ml
- connecting tube, length: 125 cm
- **universal stepped connector for safe connection of ch. 6 to ch. 18 redon drains**
- length of the drain: 500 mm
- perforation of the drain: 150 mm
- sets with drain
- individually packed and sterilized
- special designs available on request

UMDNS-Code: 11-307

Klasse/Regel: IIa/7

Art.-Nr.: 21587, PRI-LOW-VAC Redon set, 500 ml with drain (500/150) Ch. 18

Packing unit (VPE): 35, Shipping unit (VSE): 35