



**AIDAS KATAUSKAS**  
Key Account Manager Baltic States  
Healthcare  
M +370 620 25 085  
E [aidas.katauskas@ecolab.com](mailto:aidas.katauskas@ecolab.com)

**To whom it may concern**

**Manual reprocessing of instruments and endoscopes**

## **NOTIFICATION**

19<sup>th</sup> of October 2021

Kaunas

We do recommend using Ecolab® and Laboratoires Anios an Ecolab company products for manual reprocessing of any kind of surgical instruments and endoscopes.

For pre-cleaning stage - Sekusept™ Pure Clean (neutral detergent), Sekusept™ MultiEnzyme (enzymatic cleaner), Aniosyme Synergy 5 (multi-enzymatic cleaner), Aniosyme XL3 (tri-enzymatic cleaner/disinfectant) in combination with Sekusept™ Pulver (cleaner/disinfectant), Sekusept™ Aktiv (high level disinfectant/cleaner), Sekusept™ Easy (high level disinfectant) or Anioxyde 1000LD (high level disinfectant) for the final disinfection. Please follow recommendations from instrument manufacturers (use of wipes, brushes etc.), always read product label and SDS.

All products have great material compatibility and could be used in combinations according to hospital protocols.

Respectfully,

Aidas Katauskas  
Key Account Manager Baltic states  
Ecolab® Healthcare