



AIDAS KATAUSKAS
Key Account Manager Baltic States
Healthcare
M +370 620 25 085
E aidas.katauskas@ecolab.com

To whom it may concern

Manual reprocessing of instruments and endoscopes

NOTIFICATION

27th of February 2023

Kaunas

We do recommend using Ecolab® products for manual reprocessing of any kind of surgical instruments and endoscopes.

Sekusept™ Aktiv (high level disinfectant/cleaner), has approvals from endoscopes producers as Olympus and Fujinon.

Active substance of Sekusept™ Aktiv is peracetic acid what is formed from sodium percarbonate when granules are dissolved in water. Working solution of 2% concentration contents peracetic acid above 1 000 ppm. Efficacy of working solution can be checked with Sekusept Aktiv/Easy test strips. Following requirements for high level disinfection we do recommend using working solution fresh every day.

Please follow recommendations from manufacturers (use of wipes, brushes etc.) and always read product label and SDS.

Respectfully,

Aidas Katauskas
Key Account Manager Baltic states
Ecolab® Healthcare