

## AndraStent® Aortic

Bare Stent for Coarctation

Cobalt chromium bare metal stent  
Diameters 8 – 32 mm  
Lengths 13 – 57 mm

# AndraStent® Aortic XL




- Small introducer sheath 7 Fr – 13 Fr
- Diameters 8 – 25 mm
- Lengths 13 – 57 mm
- We recommend to use our AndraBalloon® or Andra2Balloon® as stent delivery catheter

## Order Codes AndraStent® Aortic XL

Order Code	Length in mm	Expansion Stent Diameter in mm						
Initial Values		8	12	15	18	20	22	25
ASA13XL	13	13	12	12	11	10	9	8
ASA17XL	17	17	16	16	15	14	13	11
ASA21XL	21	21	20	20	19	18	16	14
ASA26XL	26	26	24	23	22	21	19	17
ASA30XL	30	30	29	28	27	25	23	21
ASA35XL	35	35	33	31	29	27	25	24
ASA39XL	39	39	36	35	33	31	29	27
ASA43XL	43	43	40	38	36	34	32	29
ASA48XL	48	48	44	42	40	37	35	33
ASA57XL	57	56	52	50	47	44	42	40

The diameter of the introducer sheath depends on the balloon (+ 1 Fr for the stent). Packaging unit of the AndraStent® Aortic XL: 1

# AndraStent® Aortic XXL



- Small introducer sheath 10 Fr – 14 Fr
- Diameters 12 – 32 mm
- Lengths 21 – 57 mm
- Biggest stent on the market for coarctation
- We recommend to use our AndraBalloon® or Andra2Balloon® as stent delivery catheter

## Order Codes AndraStent® Aortic XXL

Order Code	Length in mm	Expansion Stent Diameter in mm							
Initial Values		12	15	18	20	22	25	28	32
ASA21XXL	21	21	20	20	20	19	18	16	13
ASA26XXL	26	26	25	24	23	22	21	19	16
ASA30XXL	30	30	29	28	27	26	24	22	18
ASA35XXL	35	35	34	33	32	31	29	26	21
ASA39XXL	39	38	37	36	35	34	33	30	25
ASA43XXL	43	41	41	40	40	39	37	34	28
ASA48XXL	48	46	45	45	44	43	40	38	32
ASA57XXL	57	55	53	52	51	49	47	43	37

The diameter of the introducer sheath depends on the balloon (+ 1 Fr for the stent). Packaging unit of the AndraStent® Aortic XXL: 1

## **Andramed GmbH**

Schiesswiesenstr. 18  
72766 Reutlingen, Germany  
Tel.: +49 (0) 7127 811 57 15  
Fax: +49 (0) 7127 811 57 11  
[info@andramed.com](mailto:info@andramed.com)  
[www.andramed.com](http://www.andramed.com)

Made in Germany  
Rev. 01/2020

