

RADIAL ARTERY INTRODUCER

Set with hemostasis valve to introduce and exchange catheters and endocavitary electrodes



■ Set components:

- sheath with a valve
- angiographic needle
- guidewire
- dilator

The set contains a guidewire of special construction. The soft distal part of the guidewire enables its introduction to the vessel, and the stiff proximal part provides an excellent scaffold for the dilator, along with preventing vessel perforation. The distal part of the sheath and dilator is finished in a way that enables an easy and atraumatic introduction of the catheter to the vessel.

Reference:	Sheath size [F]	Sheath length [mm]	Colour identification	Dilator size [F]	Angiographic needle size [G]	Guidewire size [inch]	Guidewire length [cm]	Guidewire type
INT4FKR	4	70	red	4	21	0.018	45	"J"
INT4FKRP						0.018		STRAIGHT
INT5FKR	5	110	grey	5		0.018	45	"J"
INT5FKRP						0.018		STRAIGHT
INT5F21KRP						0.021		STRAIGHT
INT5FKR23		230				60	0.018	"J"
INT5FKR23P							0.018	STRAIGHT
INT5F21KR23P							0.021	STRAIGHT
INT6FKR	6	110	green	6		0.018	45	"J"
INT6FKRP						0.018		STRAIGHT
INT6F21KRP						0.021		STRAIGHT
INT6FKR23		230				60	0.018	"J"
INT6FKR23P							0.018	STRAIGHT
INT6F21KR23P							0.021	STRAIGHT
INT7FKR	7	110	orange	7		0.018	45	"J"
INT7FKRP						0.018		STRAIGHT
INT7F21KRP						0.021		STRAIGHT
INT7FKR23		230				60	0.018	"J"
INT7FKR23P							0.018	STRAIGHT
INT7F21KR23P							0.021	STRAIGHT

When ordering use Ref. Number, e.g.: INT5FKR.

Special order can contain introducer with radiopaque marker

Model reference number: INT7FKR**M** – introducer with a marker.

Hydrophilic coating is available on request.

Model reference number: INT7FKR**H** – introducer with a hydrophilic surface protection.

