


AMPLATZER™ GUIDEWIRES



INTUITIVELY DESIGNED FOR A WIDE RANGE OF VASCULAR ANATOMY

Amplatzer Guidewires offer catheter support and tracking capabilities for a wide range of vascular anatomy. They are designed and tested to perform in conjunction with Amplatzer occluders.

Information contained herein for
DISTRIBUTION outside of the U.S. ONLY.
Always check the regulatory status of the device
in your region.



AMPLATZER GUIDEWIRES
OFFER EXCELLENT LUBRICITY
AND POSITIONING
CAPABILITIES WITH MINIMAL
INTERFERENCE
TO VASCULATURE.

INTUITIVELY DESIGNED

Stainless Steel

- Allows for enhanced positing capability*

PTFE-coating

- Provides lubricity for smooth navigation through tortuous anatomy*

J-tip Design

- Reduces the risk of internal wall damage for atraumatic positioning*

Exchange Length up to 300 cm

- Assures adequate length for exchange of a wide variety of devices over the wire*

Four Curve Properties

- Designed to reach targeted anatomy*

Two Body Types

- Super Stiff
 - Provides support to track the Amplatzer TorqVue™ Delivery System to the treatment site*
- Soft Tip, Fixed Core
 - Facilitates access to targeted defect location*

* Test(s) performed by and data on file at Abbott.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.
Always check the regulatory status of the device in your region.

ORDERING INFORMATION

AMPLATZER GUIDEWIRE SPECIFICATIONS

Model / Reorder Number	Diameter (inch)	Body	Tip Description	Usable Length (cm)
9-GW-001	0.035	Super Stiff	7.5 mm, Modified J-tip	260
9-GW-002	0.035	Super Stiff	1.5 mm, Modified J-tip	260
9-GW-003	0.035	Super Stiff	6 mm, J-tip	300
9-GW-004*	0.035	Soft Tip, Fixed Core	6 mm, J-tip	300

* The Amplatzer Noodlewire is a soft tip, flexible guidewire recommended for establishing an Arterial-Venous loop, facilitating closure of ventricular septal defects.

AMPLATZER GUIDEWIRES ARE DESIGNED TO BE USED WITH THE FOLLOWING AMPLATZER OCCLUDERS.

Model / Reorder Number	Amplatzer™ Occluder
9-GW-001	Amplatzer™ Duct Occluder
9-GW-002	Amplatzer™ Septal Occluder Amplatzer™ Cribriform Occluder
9-GW-003	Amplatzer™ Muscular VSD Occluder
9-GW-004	Amplatzer™ P.I. Muscular VSD Occluder Amplatzer™ Membranous VSD Occluder

Abbott is focused
on reducing risk by continuously
finding ways to put more control
into the hands of those who
save and enhance lives.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott

3200 Lakeside Dr, Santa Clara, CA 95054, USA

www.cardiovascular.abbott

TM Indicates a trademark of the Abbott group of companies.

© 2020 Abbott. All Rights Reserved.

MAT-2004323 v1.0 | Item approved for Global OUS use only.

