



EndoSwab technical details

Product name: EndoSwab

Usage: Innovative, single use cleaning brush for flexible endoscope channel
Test for residual protein with Pyromol-E

Swab material: Highly absorbent polyester fibres

Sizes/Versions:

32 p. d. nr.1

1.7: length 2.3m diameter 1.7mm

2.8: length 2.3m diameter 2.8mm

3.8: length 2.3m diameter 3.8mm

5.0: length 2.3m diameter 5.0mm

Endoswab	Best used for Endoscope channel
1.7 mm	1.5 – 1.7 mm
2.8 mm	2.0 – 2.5 mm
3.8 mm	2.6 – 3.0 mm
5.0 mm	3.1 – 4.5 mm

NOTE: Due to possible constriction or blockage of a channel, please always verify the right size of Endoswab with an initial test. Do not apply force when sliding the Endoswab through the channel.

Faulty material/operation cannot be ruled out 100%. After use always check if the tip of the swab is still attached and undamaged. EndoSwab is a single use product.

EndoSwabs are validated for the PEREG GmbH Pyromol-Test only, do not use in combination with unsuitable protein tests.

PEREG GmbH

PEREG GMBH
KONFIGURATION 12
7G
00
62

30.04.2018

Mar
Mar



PEREG GMBH, PORSCHESTRASSE 12, D-84478 WALDKRAIBURG

„SANOVUS“, UAB
Mr. Mindaugas Danilevicius
Daugėlišio g.32-301
LT-09300 Vilnius

PEREG GMBH
ANALYTIK. HYGIENE

PORSCHESTRASSE 12
D-84478 WALDKRAIBURG
TELEFON (0049) 84100
TELEFAX (0049) 84162
E-MAIL: info@pereg.de
<http://www.pereg.de>

32 p.d. nr. 2, 3

Waldkraiburg, 28.08.2023

Patvirtinimas dėl Pyromol-Test / Pyromol-E

Šiuo laišku patvirtinama, kad Pyromol testų rinkiniais galima aptikti baltymų likučius ant chirurginių instrumentų ir lanksčių endoskopų paviršių po rankinio arba automatinio apdorojimo procesų. Aptikimo riba - 1 µg - yra vienoda tiek šviežių, tiek denatūruotų baltymų likučiams. Rezultatas (spalvos pokytis) išryškėja per 5 minutes.

You

B H

12

URG

41 00

41 02

Mar

Mar

PEF

✓