



WE UNDERSTAND.



NEUROSURGERY

ACCESSORIES FOR MIETHKE SHUNT SYSTEMS

TITANIUM CONNECTORS











- Preferably for use with catheters of:

- internal diameter
 $d_i = 1.2 \text{ mm}$

- external diameter
 $d_o = 2.5 \text{ mm}$

- Connector diameter:
 $d_o = 1.9 \text{ mm}$


TITANIUM CONNECTORS


Straight		Art. no.	Sales unit
		FV012T	1 pc. 16 pirkimo dalis
	$\pm 1.9 \text{ mm}$	FV013T	5 pcs.
	1:1	FV014T	10 pcs.
L-shape		Art. no.	Sales unit
		FV051T	1 pc.
	$\pm 1.9 \text{ mm}$		
	1:1		
Y-shape		Art. no.	Sales unit
		FV015T	1 pc.
	$\pm 1.9 \text{ mm}$	FV016T	5 pcs.
	1:1	FV017T	10 pcs.
		18 pirkimo dalis	
T-shape		Art. no.	Sales unit
		FV018T	1 pc.
	$\pm 1.9 \text{ mm}$	FV019T	5 pcs.
	1:1	FV020T	10 pcs.
X-shape		Art. no.	Sales unit
		FV021T	1 pc.
	$\pm 1.9 \text{ mm}$		
	1:1		

TITANIUM CONNECTORS

- Step-down connectors are designed for connecting standard catheters to lumbar catheters
- Step-down connector diameters: $d_o = 1.9 \text{ mm}$ to 1.4 mm


TITANIUM CONNECTORS, STEP-DOWN

	Straight	Art. no.	Sales unit
	1.9 mm \perp 1.4 mm 1:1	FV050T 17 pirkimo dalis	1 pc.

	L-shape	Art. no.	Sales unit
	1.9 mm \perp 1.4 mm 1:1	FV052T	1 pc.

- L-P connectors are designed for connecting two lumbar catheters
- Connector diameters for L-P drainage: $d_o = 1.4 \text{ mm}$

TITANIUM CONNECTOR, L-P DRAINAGE

	L-shape	Art. no.	Sales unit
	1.4 mm \perp 1.4 mm 1:1	FV053T	1 pc.