

SOFIA® EX

Intracranial Support Catheter

8. Intermediate kateteriai

EMEA

EXcellent Trackability. EXtra Stability.

EXcellent Trackability for distal navigation and positioning.¹


- Soft distal tip designed for easy navigation to the desired target location



EXtra Stability for reliable device delivery.²

- Exhibits minimal tip movement when deploying and re-sheathing treatment devices²
- Kink resistant construction provides consistent support for device delivery in tortuous anatomy

8. Intermediate kateteriai



SOFIA® EX Intracranial Support Catheter

Catalog Number	Distal/Proximal ID		Distal OD			Proximal OD			Usable Length (cm)	Distal Soft Length (cm)	Tip Shape
	(in)	(mm)	(F)	(in)	(mm)	(F)	(in)	(mm)			
ISC5105ST	0.058	1.47	5.2	0.068	1.73	5.4	0.071	1.80	105	9	Straight
ISC5115ST									115		

Kateterio dydis 5 Fr ir 6 Fr;

Distaliai matmenys – OD ne mažiau 0,068 / 1,7; pasirinktinai 0,0815 / 2,1mm;

Proksimaliniai matmenys – OD ne mažiau 0,068 / 1,7; pasirinktinai 0,0825 / 2,1mm;

Vidinis diametras ID 0,055" ir 0,070";

Distalinio galiuko forma – tiesi;

Distalinio galiuko ilgis 17 cm ir 19 cm;

Proksimalaus galo ilgis: 96 cm, 98 cm kai kateterio darbinis ilgis 115 cm; 108 cm, 106 cm kai kateterio darbinis ilgis 125 cm; 112 cm kai kateterio darbinis ilgis 131 cm.

Darbinis kateterio ilgis: 115 cm; 125 cm; 131cm;

1. TR19-145A; n=3
2. TR19-073A; n=3

INDICATIONS FOR USE

The SOFIA® EX Catheter is indicated for general intravascular use, including the neuro and peripheral vasculature. The SOFIA® EX Catheter can be used to facilitate introduction of diagnostic agents or therapeutic devices. The SOFIA® EX Catheter is not intended for use in coronary arteries.

For Professional Use Only

MICROVENTION and SOFIA® are registered trademarks of MicroVention, Inc.
© 2019 MicroVention, Inc. MM872(i) EMEA 06/19

MicroVention Worldwide Innovation Center

PH +1.714.247.8000

35 Enterprise

Aliso Viejo, CA 92656, USA

MicroVention UK Limited

PH +44 (0) 191 258 6777

MicroVention Europe, S.A.R.L.

PH +33 (1) 39 21 77 46

MicroVention Deutschland GmbH

PH +49 211 210 798-0

Web

microvention.com