

▶ **Wide Range**

- For single use or reusable
- With manometer or pressure indicator
- Available in following sizes:  
500 ml, 1000 ml, 3000 ml

▶ **Patient Safety**


- Reliable pressure control
- Latex free

▶ **Hygiene**

- Easy cleaning and disinfection
- Infection control

▶ **Design**

- Fits most fluid bag configurations
- Durable construction



# Pressure Infusion Cuffs

Manual systems for rapid pressure infusion and irrigation purposes

# Pressure Infusion Cuffs – Du

Thanks to a large choice of options, the VBM pressure infusion cuffs cover all requirements for a safe management of fluid infusions or irrigations.

Available as reusable and single-use versions, all VBM options incorporate a high quality design of cuffs with durable materials and reliable pressure control mechanisms. The different configurations will fulfil the expectation of use in the OR, ICU, EMS and outdoor emergency.

## Pressure Infusion Cuff

- ▶ Reusable
- ▶ Easy cleaning & disinfection thanks to water-repellent coating
- ▶ Cost effective



**Hand Inflator**  
- Large volume allows quick inflation  
- Ergonomic & latex free

**Three-Way Stop Cock**  
- For pressure control  
- Quick deflation

### Versions / Features



**Precision Manometer**  
Robust & durable



**Hand inflator with colour coded pressure indicator**

**Double Hook**  
Easy fixation of most fluid bag configurations

**Translucent Netting**  
Clear view on the fluid bag and the fluid level

## Order information

**Pressure Infusion Cuff**  
Pocket design with double hook, reusable

Version	500 ml	1000 ml	Box
<b>A</b> with hand inflator, stop cock and precision manometer	58-08-050	58-08-100	1
<b>B</b> with hand inflator, stop cock and colour coded pressure indicator	58-02-050	58-02-100	1

# able – Hygienic – Reliable – Ease of use

The special water-repellent material of the reusable systems makes it easy to clean or disinfect and contributes to reduce the risk of cross contamination.

A large choice of pressure control versions gives the user the possibility to set up their priority for precision of measurement and cost effectiveness.

## Infusor

- ▶ Reusable
- ▶ Easy cleaning & disinfection thanks to smooth surface
- ▶ 360° see-through design



**Hook**  
Easy fixation of most fluid bag configurations

**Transparent**  
Clear 360° view on the fluid bag and the fluid level


**Hand Inflator**

- Large volume allows quick inflation
- Ergonomic & latex free

**Three-Way Stop Cock**

- For pressure control
- Quick Deflation

### Versions / Features



### Wrap-Around Design

Accommodates most fluid bag configurations thanks to wrap-around design and provides an excellent pressure distribution



## Order information

### Infusor

Transparent wrap-around design with hook, reusable

Version	500 ml	1000 ml	3000 ml	Box
<b>A</b> with hand inflator, stop cock and precision manometer	57-08-050	57-08-100	57-08-300	1
<b>B</b> with hand inflator, stop cock and colour coded pressure indicator	57-02-050	57-02-100	57-02-300	1
<b>C</b> with Luer Lock connection only	-	-	57-01-300	1

# Dispo Infusor

- ▶ For single use
- ▶ Infection control



## Hand Inflator

- Large volume allows quick inflation
- Ergonomic & latex free

## Three-Way Stop Cock

- Pressure control
- Quick deflation

## Versions / Features



Hand inflator with colour coded pressure indicator



Hand Inflator and stop cock only



Luer Lock connection



## Hook Fastener

Easy fixation of most fluid bag configurations

## Translucent Netting

Clear view on the fluid bag and the fluid level

## Order information

### Dispo Infusor

Pocket design with hook fastener (500 ml & 1000 ml) and hook (3000 ml), for single use

Version	500 ml	1000 ml	Box	3000 ml	Box
<b>A</b> with hand inflator, stop cock and colour coded pressure indicator	56-02-050	56-02-100	5	56-02-300	1
<b>B</b> with hand inflator and stop cock	56-03-050	56-03-100	5	-	-
<b>C</b> with Luer Lock connection only	56-01-050	56-01-100	5	56-01-300	1

**VBM** Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany

Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

CE  
0123

PRO006\_1.0\_EN\_02.19