


ESVIsensor ECG electrode

ES GS40 Electrode

Application: monitoring



- 1- One side siliconized PET release liner.
- 2- Ring of gel cured by UV lamp.
- 3- PE foam 1 mm thick
- 4- Round $\varnothing 24$ mm PE label.
- 5- Sensor Ag/AgCl
- 6- Stainless sensor stud
- 7- Paper finger tip

Packing:

- pouch 50 pcs.
- box 2 000 pcs.

ESVICO Sp.z o.o.

84-214 Bożepole Wielkie, ul. Chmieleniecka 15, Poland

NIP: 588-237-42-30, REGON: 221198703, KRS: 0000420978

ALIOR BANK S.A. nr konta: 96 2490 0005 0000 4500 3009 0184