


66, 67, 68, 69 pishimo delys

Misago®

RX Self-expanding Peripheral Stent

For a lifetime of flexibility




- + Extended durability
- + Precise deployment
- + Easy to use RX design

 **TERUMO**

New Generation RX Self-expanding stent

- + The special zig-zag cell design provides **optimal flexibility and durability** for the best long term patency for the patient
- + The unique delivery system allows **precise and easy stent deployment**
- + The 0.035"/0,89 mm guide wire compatible Rapid Exchange technology **allows for single operator handling**



✓ 3 radiopaque gold markers
at both ends

✓ 8 cells for 6-8 mm
9 cells for 9-10 mm

✓ Nitinol fabrication




Misago

Application:
Occlusive/stenotic
lesions in
iliac/femoro-popliteal
arteries





Ergonomic handle

Thumb wheel



Durability compared with compression




In-vitro test at Terumo Corporation, Japan

Ordering information

+
new length
120/150mm
+

Stent Diameter	Stent length	40 mm	60 mm	80 mm	100 mm	120 mm	150 mm
6 mm		SX-VXA0640LN	SX-VXA0660LN	SX-VXA0680LN	SX-VXA06X0LN	SX-VXA06X2LN	SX-VXA06X5LN
7 mm		SX-VXA0740LN	SX-VXA0760LN	SX-VXA0780LN	SX-VXA07X0LN	SX-VXA07X2LN	SX-VXA07X5LN
8 mm		SX-VXA0840LN	SX-VXA0860LN	SX-VXA0880LN	SX-VXA08X0LN	-	-
9 mm		SX-VXA0940LN	SX-VXA0960LN	-	-	-	-
10 mm		SX-VXA1040LN	SX-VXA1060LN	-	-	-	-
Coaxial part length		46 cm	48 cm	50 cm	52 cm	64 cm	67 cm

Effective catheter length: 135 cm for all sizes



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-8572
Japan
Tel: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel: +32 16 38 12 11
Fax: +32 16 40 02 49

**Terumo Europe
Interventional Systems
Liaison Office**
Immeuble Krebs
8, rue Jean-Pierre Timbaud
CS 10720
78067 Saint Quentin
en Yvelines Cedex,
France
Tel: +33 1 39 30 95 75
Fax: +33 1 39 30 95 70

European Offices

**Terumo Europe N.V.
Benelux Sales Division**
Interleuvenlaan 40
3001 Leuven,
Belgium
Tel: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands
Tel: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.
Parc Ariane - Bât. Uranus
1, rue Helene Boucher
78284 Guyancourt Cedex
France
Tel: +33 1 30 96 13 00
Fax: +33 1 30 43 60 85

Terumo Italia Srl
Via Simone Martini, 143/145
00142 Roma RM,
Italy
Tel: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo UK Ltd
Tamesis, The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel: +44 1784 476 200
Fax: +44 1784 476 201

Terumo Deutschland GmbH
Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Europe España SL
Avda. Juan Carlos I, N° 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcala de Henares (Madrid),
Spain
Tel: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Sweden AB
Nya Varvet, Byggnad 90
Sven Kattefells Gata 18
426 71 Västra Frölunda, Sweden
Tel: +46 3174 85 880
Fax: +46 3174 85 890

**Terumo Denmark
Filial of Terumo Sweden AB**
Tel: +45 7020 93 80
Fax: +45 7020 94 80

**Terumo Europe N.V.
Moscow Representative Office**
JSC Northern Tower
13F, 10 Testovskaya Street
Moscow 123317, Russia
Tel: +7 495 988 4740
Fax: +7 495 988 4739

**Addresses subject to change,
please consult our Website.**

© Registered Trademark
Published by Terumo Europe N.V.

CE 0197